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R2 18.8 mL

M 16.0 mL

COBAS, COBAS E and
ELECSYS are trademarks
of Roche.

Cal1 0.67 mL

Cal2 0.67 mL

LOT 50226801

 2021-08-31

 Roche Diagnostics GmbH
Sandhofer Str. 116
D-68305 Mannheim

IVD  0123  2-8 °C

Distribution in USA by:
Roche Diagnostics,
Indianapolis, IN
Made in Germany

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
REF 07028024190

For USA: **CONTENT**

- R1 Toxoplasma-antigen labeled with ruthenium complex > 1 mg/L
- R2 Biotinylated monoclonal anti-h-IgM antibody (mouse) > 500 µg/L
- M Streptavidin-coated microparticles 0.72 mg/mL
- Cal1 Human serum, negative for anti-Toxo IgM
- Cal2 Human serum, reactive for anti-Toxo IgM

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
Precautions and warnings

For in vitro diagnostic use.

Exercise the normal precautions required for handling all laboratory reagents. Disposal of all waste material should be in accordance with local guidelines. Safety data sheet available for professional user on request.

For USA: Caution: Federal law restricts this device to sale by or on the order of a physician.

This kit contains components classified as follows in accordance with the Regulation (EC) No. 1272/2008:

2-methyl-2H-isothiazol-3-one hydrochloride

EUH 208 May produce an allergic reaction

Product safety labeling follows EU GHS guidance

All human material should be considered potentially infectious.

All products derived from human blood (TOXOIGM Calt, TOXOIGM Calt) are prepared exclusively from the blood of donors tested individually and should be free from HBsAg and antibodies to HCV and HIV.

The serum containing anti-Toxo IgM (TOXOIGM Calt) is not sterile.

The testing methods used assays approved by the FDA or cleared in compliance with the European Directive 98/79/EC, Annex II, List A.

However, as no testing method can rule out the potential risk of infection with absolute certainty, the material should be handled with the same level of care as a patient specimen. In the event of exposure, the directives of the responsible health authorities should be followed.

Avoid foam formation in all reagents and sample types (specimens, calibrators and controls).

Method Sheet download

The product code together with the respective document version (if applicable) is stated in the download information boxes listed above. Please go to the website listed in the relevant box to download this document. The website displays either the version listed here or, if available, any later version that is also valid for this product code. If you have no internet access, please contact your local Roche affiliate to obtain the document free of charge.

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2-Methyl-2H-isothiazol-3-on-Hydrochlorid

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