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APOSTILLE

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1. COUNTRY: UNITED STATES OF AMERICA
2. THIS PUBLIC DOCUMENT HAS BEEN SIGNED BY:
JANICE M LOPEZ
3. ACTING IN THE CAPACITY OF:
NOTARY PUBLIC OF NEW JERSEY
4. BEARS THE SEAL/STAMP OF:
JANICE M LOPEZ, NOTARY

CERTIFIED

5. AT TRENTON, NEW JERSEY
6. THE 15TH DAY OF NOVEMBER 2022
7. BY: Elizabeth Maher Muoio
State Treasurer
8. NO: A824841
9. SEAL/STAMP:
10. SIGNATURE

A handwritten signature in black ink, appearing to read "Elizabeth Maher Muoio".


Elizabeth Maher Muoio
State Treasurer

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
Ailin Ortiz – Coordinator International Regulatory Affairs

Integra LifeSciences Corporation

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
JANICE M. LOPEZ
NOTARY PUBLIC OF NEW JERSEY
Commission # 50048048
My Commission Expires 10/19/2026

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
monopolar and bipolar cables, electrodes, and instruments, which are cleared and available separately.

The Generator is part of an electrosurgical system to be used with electrosurgical instruments such as a pair of bipolar forceps and can be used standalone or with the Codman Electrosurgical Irrigator when irrigation is needed. International power cords will be supplied with the Generator also, depending on which country the Generator is going to. Power cord is supplied with the Irrigator also.


Codman Electrosurgical Generator Front View:

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
Codman Electrosurgical Generator Back View:

Codman Electrosurgical Generator Components Descriptions

#	Component	Description
External Features		
1	Power Button	Controls the primary AC power to the Generator. When on, the display

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(2.2) RF generator

The RF generator is placed on the left front side at the main board. It has been configured as a push pull amplifier and works at a fixed frequency of 4 MHz. The microcontrollers generate a pulse control signals for the RF generator. The RF generator produces different signals for the chosen kind of current and dependent on the current mode.

There is a matching network between the RF generator and the output transformer. RF power is uncoupled through separate output transformers for the monopolar and bipolar RF modes and guide to the respective connector socket via the RF output board.

The output transformers are mounted at the front of the main board and have plugs at the winding ends coming out of the transformer. The plugs are connected with the RF output board.

(2.3) Power supply The DC power supply provides the RF generator with a DC voltage from 0V to 60V. The DC voltage value is a setting value to regulate the RF output power.

The voltage is generated by a primary synchronized switching supply unit in voltage fed half bridge circuits located on the right side at the rear of the main board. The switch mode power supply is addressed by a switching power supply controller on the secondary side. The switch transistors on the primary side are addressed by pulse transformers with a working frequency of about 100 kHz. The power supply is protected by a current limiting circuit in a case of failure, such as by a short circuit in the RF generator. The main fuses also serve as a primary side protection for the switching power supply. The fuse is a T3,15A/ 250V protecting the entire system for voltages from 100V to 240V.

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The RF current can be activated either by the foot pedal or a finger switch. The activation signals CUT and COAG from the finger switch of the handpiece are transmitted to the controller system via two optical fibers. The auxiliary energy is generated by galvanic isolated power supply.

(4) Balun (monopolar and bipolar)

Means for suppressing high frequency leakage currents.


(5) Front board

The front board includes all keys for selecting the current waveform, special function and program mode. The activation LEDs, NE LED and the seven-segment displays are also on the front board. When the unit is switched on, the keys are illuminated and are visible even in a dark environment.

(6) Rotary encoder

The set points can be adjusted via the rotary encoder.

Power cord supply together with Electrosurgical Generator. The imagine is provide below:

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




ISOCOOL bipolar forceps tip (REF 8135025S; 8135050S; 8135100S; 8135200S; 8145025S; 8145050S; 8145100S; 8145200S; 8155025S; 8155050S; 8155100S; 8155200S; 823510GL).

ISOCOOL tip should be used together with ISOCOOL Handle as bipolar forceps for Codman Electrosurgical Generator for coagulation and manipulating of tissue. ISOCOOL tips are available in a variety of lengths and diameters, as provided below:

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
In addition, the tip design helps ensure properly- aligned forceps each time. These features of the ISOCOOL Bipolar Forceps can help save procedure time by minimizing or eliminating the need to clean or replace the forceps.


The unique AHT (Active Heat Transfer) Technology employed in the forceps tips is designed to eliminate sticking of tissue. In addition, the tip design helps ensure properly-aligned forceps each time. These features of the ISOCOOL bipolar forceps can help save procedure time by minimizing or eliminating the need to clean or replace the forceps. By drawing heat away from the tips, the electrode maintains a low enough tip temperature to avoid sticking but high enough to coagulate effectively.

Schematic Diagram of AHT Technology in the ISOCOOL Tips:


Enlarged wicking structure

Controlling heat at the tip by continuously evaporating and condensing the internal working fluid

Isocool Bipolar Forceps with AHT Technology


Competitive forceps without AHT Technology

The ISOCOOL Tip contain a heat pipe, which consists of a sealed, vacuum tight copper pipe, containing a wick structure and a working fluid (water) enclosed

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

Isocool Bipolar Forceps with AHT Technology

Competitive forceps without AHT Technology

within the heat pipe. The heat pipe is evacuated and then back filled with a small quantity of water, in order to saturate the wick before being sealed. The atmosphere inside the heat pipe is set by equilibrium of liquid and vapor. As heat enters the pipe at the evaporator end (in contact with the tissue), this equilibrium is upset, generating vapor at a slightly higher pressure. The higher-pressure vapor travels to the condenser end where the slightly lower temperature causes the vapor to condense, giving up its latent heat of vaporization. The condensed liquid is then transported back to the evaporator by capillary forces present in the wick structure. This continuous cycle can transfer a large quantity of heat with very low thermal gradients. The heat pipe operates passively, driven only by the heat that is transferred.

Electrical energy is delivered from the electrosurgical generator through the electrode and ultimately to the tissue, causing the tissue to heat. Heat conducts back to the tip of the electrode and heats the tip. The heat pipe transfers the heat from the tip back towards the handle of the device, where it is released into surrounding air by convection. By drawing enough heat away from the electrode tip, tissue sticking is eliminated. The effectiveness of the device's ability to coagulate is not altered by the heat pipe technology because coagulation is accomplished through heat generation in the tissue rather than heat generation at the tip. The heat pipe extracts a sufficient amount of heat from the electrode tip to keep it cool, but not enough to significantly cool the tissue.

ISOCOOL Tip is provided sterile, 1 pcs per package.

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ISOCOOL bipolar forceps handle (REF 8010010H):


ISOCOOL bipolar forceps handle design allows the disposable ISOCOOL tips to be attached and detached from the reusable handles. ISOCOOL handle should be used together with ISOCOOL tip as bipolar forceps for Codman Electrosurgical Generator.

ISOCOOL bipolar forceps handle is designed for use with bipolar electrosurgical generator for coagulating and manipulating tissue. The handle accommodates an assortment of tips, which are available in a variety of lengths and diameters. Use only ISOCOOL tips with the ISOCOOL handle.

Electrical Insulation:

The ISOCOOL handle receives two coatings for electrical insulation: a PARYLENE (clear) coating applied by a urethane (blue) coating. Please be aware that the hybrid coating might begin to peel over time through normal usage, cleaning, and sterilization cycles. The peeling or flaking appearance of the coating does not constitute product failure. This potential condition is cosmetic in nature and does not compromise product performance or the safety of the user or the patient. As required by the standards of the International Electrotechnical Commission (IEC), the handle has been validated through simulated use testing and was demonstrated to be safe through 20 cleaning and sterilization cycles.

Pictures of the ISOCOOL bipolar forceps handle is provided below:

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Bipolar forceps (REF 80-2900; 80-2901; 80-2913; 80-2917; 80-2918; 80-2919):

The irrigating insulated bipolar forceps are designed with a non-stick mirror finish tip surface, to provide the surgeon with forceps that will incur less coagulum build up when used under normal conditions. The product line includes a broad selection of forceps size and types. The variety of patterns and sizes provides bipolar cutting and coagulating capabilities for a wide range of neurosurgical procedures. The bayonet style forceps enhance visibility by keeping the surgeon's hand farther from the line of sight. Insulated forceps shafts protect the surrounding tissue when performing deep procedures. The unique feature of the irrigating forceps is that the irrigation pipe is completely recessed into the right blade of the forceps. This offers the surgeon increased visibility at the forceps tips. Irrigating type forceps are provided with a stylet to help keep the irrigation tube free of dust and debris during storage.

The bipolar forceps are provided non-sterile and must be cleaned and sterilized before each use, including initial use. For multiple use.

Picture of bipolar forceps is provided below:

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Disposable Bipolar Cord, Unitized Plug (REF 80-1188):

For use with Codman non-irrigating disposable and reusable bipolar forceps. This cord is designed for use with bipolar electrosurgical units only.

Supply sterile, for single use

Precautions:

Inspect the sterile package carefully. Do not use if:

- the package or seal appears damaged,
- contents appear damaged, or
- the expiry date has passed.

Do not reuse. This product is for single use only. Use aseptic technique in all phases of handling this product. Leave cord loose. Do not secure cord.

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Codman Electrosurgical Generator is part of an electrosurgical system, it could use with electrosurgical irrigator, foot pedal. Other accessories, available separately, include a cord/tubing set, interconnecting cable, foot pedal, and floor stand.

1. Codman Electrosurgical Irrigator (REF 901001IRRO):

Irrigation module for delivering irrigation fluid to bipolar forceps during an electrosurgical procedure. Codman Electrosurgical Irrigator consists of a power supply board and a control board which includes a microcontroller and motor controller. The motor controller operates the Irrigator by use of a peristaltic pump controlled by a stepper motor. The Irrigator functions when connected to an electrosurgical generator, driven by a foot pedal connected to the generator. The tubing portion of a cord and tubing set is spiked into an IV bag at one end and inserted into the irrigation plug of a pair of bipolar irrigating forceps. The appropriate silicone portion of the tubing is placed onto the pump of the Irrigator and delivers flow from the IV bag to the tips of the bipolar irrigating forceps when the Irrigator is activated. The Irrigator has 21 settings from 0 – 20. These settings correspond to a flow profile of a slow drip to constant stream. The Irrigator is operational in either cut or coagulation mode of the electrosurgical generator and the Irrigator's flow profile is not impacted by either mode.

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



Rear and Bottom Irrigator View:

Irrigator Components Descriptions		
#	Component	Description
External Features		

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
Codman Electrosurgical System Interconnecting Cable was designed through reverse engineering to replicate the form, fit, and function of Codman devices. It function as an interconnecting cable between an irrigator and a generator system. Table below shows the features that support design equivalence.

REF	Overall Length	Connector Diameter	Connector Type	Pin to Pin continuity	Label Placement ($\pm 3\%$)
901001CBL1	1574 mm $\pm 3\%$	19 mm	One 4 pin and one 3 pin	< 1 Oh m	3.81 mm from connector

Schematic picture is provide below:

Shell Contact Schematic for Codman Electrosurgical System Interconnecting Cable:

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
Note: numbering at contact socket and connector pin areas optional

3. Codman Electrosurgical System Foot pedal (REF 901001FPED):

Codman Electrosurgical Foot Pedal (REF 901001FPED) used to drive the functionality of the Generator, and Irrigator when the Irrigator is connected and powered on. Foot Pedal with two step controls (one for activating COAG and one for CUT) and cable to connect to RF generator.

Schematic picture is provide below:

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

20. Adjustment Bar When raised, the adjustment bar allows the user to reposition the Foot Pedal with his/her foot mid-procedure.
21. Cord management feature Allows the user to wrap the Foot Pedal cord for storage.
22. CUT pedal (yellow) Depressing this pedal activates cutting power to the instrument.
23. COAG pedal (blue) Depressing this pedal activates coagulation power and irrigation to the instrument.

4. CODMAN Integrated Bipolar Cord and Tubing Set, Fixed Leads, Gravity Feed (REF 9190001GF)
and
5. CODMAN Integrated Bipolar Cord and Tubing Set, Fixed Leads, Rotary Pump (REF 9190001RP)

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
Both are for use with irrigating bipolar forceps. The CODMAN Integrated Bipolar Cord and Tubing Sets are intended to provide irrigation and energy simultaneously to bipolar forceps specifically designed for irrigation.

CODMAN Integrated Bipolar Cord and Tubing Set, Fixed Leads, Rotary Pump (REF 9190001RP) and CODMAN Integrated Bipolar Cord and Tubing Set, Fixed Leads, Gravity Feed (REF 9190001GF) will use with Irrigator.

CODMAN Integrated Bipolar Cord and Tubing Set, Fixed Leads, Gravity Feed (REF 9190001GF) use if surgeon want to use a gravity feed system instead of irrigator. 9190001GF will be connected to an pole IV. Gravity will push the solution into the forceps. The IV pole must therefore be suspended higher than the forceps. There is a part of the tubing which is special to add a clamp system to prevent the liquid from flowing unnecessarily.

Bipolar Cord and Tubing Sets consists of an inlet spike with protective cap, an on/off clamp, silicone tubing and fittings, a male luer connector, two plugs, and a female socket connector. The set is packaged sterile for single use. The components of devices are provided below:

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
Both supply sterile, for single use, packed 10 pcs/box.

The pictures of device and its component is provided below:

Medical devices to be used in combination with medical device submitted for registration:

Bipolar forceps named "CODMAN VersaTru" (manufacturer by Integra LifeSciences Switzerland Sarl) could be used together with Codman Electrosurgical Generator also.

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
Codman Electrosurgical Generator (main unit):

Mains Supply	100-240 VAC ; 50/60 Hz
Power Consumption	Without HF output: approximately 50 VA At maximum output power: approximately 500 VA
Protection Class	I
IP protection	IP2X
LF and RF Leakage Currents	According to IEC 60601-2-2
Application Type	CF; defibrillator safe (protection of the product against defibrillator discharge depends on the cables used)
Operating Mode ($\pm 5\%$)	Intermittent 10 seconds ON / 30 seconds OFF
Startup time	Max. 10 sec
Mains Fuse	2 x T 3.15 AH 250 V 5x20 mm
Signal Level	RF indicator: 40-65 dB(A) adjustable. Alarm: 65 dB(A)
Weight ($\pm 5\%$)	5.0 kg
Dimensions ($\pm 5\%$)	Width 403 mm Height 170 mm Depth 313 mm
Rated Duty Cycle	< 30 sec on / < 30 sec off
Temperature sensor	Yes
Contact Quality Monitor (CQM)	Yes
Power Activation Control	Finger Switch and/or Foot Pedal
AC Powered	yes
Battery Operated	No
Light Source	No
Sterile	no

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Power Cord:

Rating	10A 250V
Insulation Resistance	Not less 100M Ω
CONDUCTOR COLORS	BLUE / BROWN / GREEN-YELLOW

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CABLE JACKET COLOR	black
OVERALL LENGTH, in mm (±6%)	1800

Codman Electrosurgical Irrigator:

Supply Voltage	60 VA at 100–240 VAC
Supply Frequency	50/60 Hz
IPX PROTECTION	IPX1
Fuses	5 x 20 mm, F 2H, 250 V (fast blow)
Weight (±5%)	Approximately 1.78 kg
Dimensions (±5%)	173 x 149,86 x 247,14 mm
Power cord 2	Length is 2000 mm Single medical grade power cord with the following detachable connectors: Australian UK European Schuko

Codman Electrosurgical System Foot pedal:

Supply	25V AC/60V DC
IPX Protection Class	IPX8
Weight (±6%)	2,2 kg
Dimensions (±6%)	Width 403 x height 170 x depth 313 mm

Codman Electrosurgical System Interconnecting Cable:

Length, mm (±3%)	1574 mm
------------------	---------

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Reusable Bipolar Cord, Unitized Plug:

Length, mm ($\pm 2\%$)	3660 mm
Weight, g ($\pm 6\%$)	60 g

Disposable Bipolar Cord, Unitized Plug:

Length, mm ($\pm 2\%$)	3660 mm
Weight, g ($\pm 6\%$)	60 g

CODMAN Integrated Bipolar Cord and Tubing Set, Fixed Leads, Gravity Feed:

Length, mm ($\pm 6\%$)	4792 mm
Weight, g ($\pm 6\%$)	36 g
Wall thickness, mm ($\pm 6\%$)	0,7239 mm
Overall internal surface area, cm ² ($\pm 6\%$)	297,9 cm ²

CODMAN Integrated Bipolar Cord and Tubing Set, Fixed Leads, Rotary Pump:

Length, mm ($\pm 6\%$)	4699 mm
Weight, g ($\pm 6\%$)	36 g
Wall thickness, mm ($\pm 6\%$)	1,905 mm
Overall internal surface area, cm ² ($\pm 6\%$)	283,6 cm ²

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The maximum T of tips during the procedure of coagulation is 80 °C.
The technical specification of all REF of ISOCOOL tips are provided below.
Tolerance is $\pm 1.3\%$

Name	REF	Overall Length of Forceps, then in complete with handle	Tip Length (working length)	Tip Size	Tip Properties
Low profile					
ISOCOOL bipolar forceps tip, working length 89 mm, tip 0,25 mm	8135025 S	8 ½ (216 mm)	3 ½ (88,9 mm)	0,25 mm	Non-Stick, Non-irrigating, Gold Tips
ISOCOOL bipolar forceps tip, working length 89 mm, tip 0,50 mm	8135050 S	8 ½ (216 mm)	3 ½ (88,9 mm)	0,50 mm	Non-Stick, Non-irrigating, Gold Tips
ISOCOOL bipolar forceps tip, working length 89 mm, tip 1,00 mm	8135100 S	8 ½ (216 mm)	3 ½ (88,9 mm)	1,00 mm	Non-Stick, Non-irrigating, Gold Tips

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Bipolar forceps.

The maximum T of tips during the procedure of coagulation is 80 °C.
The technical specification of all REF of bipolar forceps are provided below.
General tolerance is $\pm 2\%$, if not specified for a specific amount.

Name	REF	Type	Overall length without pin	Working length	Tip diameter
Bipolar forceps, irrigating, bayonet, overall length 216 mm, tip 0,75 mm	80-2900	Insulated irrigating, bayonet, non-stick, mirror finish	8 ½ inches 213,5 $\pm$ 3 mm	4 ½ inches 114 mm	0,75 mm
Bipolar forceps, irrigating, bayonet, overall length 197 mm, tip 1,5 mm	80-2901	Insulated irrigating, bayonet, non-stick, mirror finish	7 ¾ inches 195.07 $\pm$ 3,05 mm	3 inches 76 mm	1,5 mm
Bipolar forceps, irrigating, straight, overall length 130 mm, tip 0,75 mm	80-2913	Insulated irrigating, straight, non-stick, mirror finish	5 inches 127,76 $\pm$ 3,0 mm	1 ¾ inches 41 mm	0,75 mm

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Generator; Irrigator; Connecting Cable; Codman Electrosurgical Foot Pedal	MD in contact with intact skin and hands protected by sterile gloves	category A - MD of short-term contact - the total duration of a single, repeated or continuous exposure of which is not more than 24 hours
Reusable Bipolar Cord, Unitized Plug; Disposable Bipolar Cord, Unitized Plug	MD contacting the skin and mucous membranes	category A - MD of short-term contact - the total duration of a single, repeated or continuous exposure of which is not more than 24 hours
CODMAN Integrated Bipolar Cord and Tubing Set, Fixed Leads, Gravity Feed; CODMAN Integrated Bipolar Cord and Tubing Set, Fixed Leads, Rotary Pump	MD, attached from the outside, contact with soft tissues, bone tissues, serve as conductors for the delivery of fluids	category A - MD of short-term contact - the total duration of a single, repeated or continuous exposure of which is not more than 24 hours
ISOCOOL tips; ISOCOOL handle; Bipolar Forceps	MD attached from the outside, contact with soft tissues, bone tissues.	category A - MD of short-term contact - the total duration of a single, repeated or continuous exposure of which is not more than 24 hours

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For Codman Electrosurgical Generator, foot pedal, ISOCOOL tips and ISOCOOL handle and Bipolar Forceps:

Operating Temperature Range	+10° C to +40° C
Operating Humidity Range	30% to 75% relative humidity, non-condensing
Operating Pressure Range	700 hPa to 1050 hPa

Precautions:

- The Generator must be separated from other electrosurgical devices and their cords. Close proximity to the Generator and its cords to electrosurgical units that produce excessive RF current radiation could cause the Generator to change output power. Surgical electrode cords must be positioned to avoid contact with the patient or other leads. Temporarily unused active electrodes should be stored in a location that is isolated from the patient.
- Do not operate the Generator at temperatures below 50°F (10°C). Allow the Generator to warm up to at least 50°F (10°C) before use.
- Visually inspect the power cord, bipolar cord and Foot Pedal for damage prior to each use. Connect Foot Pedal, adjust Generator to desired power setting, and activate Foot Pedal to verify proper function. Replace accessories if necessary.
- Non-insulated bipolar forceps are not recommended for use with higher power settings. When higher power settings are employed, insulated bipolar forceps are recommended.
- Keep bipolar instrument tips clean. Build-up of eschar may reduce the instrument's effectiveness. Do not activate the instrument while cleaning, as injury to operating personnel may result.
- The Generator is designed for use without a neutral electrode while operated in bipolar modes.
- The Generator does not generate waste. Disposable products used with the Generator must be disposed of in accordance with hospital procedures and regulations.

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- Always position Generator so that the power cord can be easily unplugged from the rear panel.
- The output power should be set as low as is necessary to achieve the desired effect.
- Anyone who connects additional units to the input or output section thereby configures the system and is thus responsible for compliance with section IEC 60601-1:16.
- In order to be able to completely separate the unit from the grid in the event of danger, either the power socket of the unit or the outlet into which the mains cable is plugged should remain accessible.
- Do not immerse the Foot Pedal in any liquids.
- Do not use ISOCOOL Bipolar Forceps at generator power settings that exceed peak- to-peak RF voltage of 500 Volts.
- For generators, refer to the operator's manual to determine the maximum allowable power settings.
- Activating the electrosurgical unit simultaneously with the aspiration/irrigation might alter the path of electrical energy away from the tissue.
- Use only with bipolar cords available from Integra.

For Irrigator, Interconnecting Cable and Bipolar Cords:

Operating Temperature Range	+10° C to +40° C
-----------------------------	------------------

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Operating Humidity Range	20% to 80% relative humidity, non-condensing
Operating Pressure Range	700 hPa to 1050 hPa

Precautions:

- Use only Codman Integrated Bipolar Cord and Tubing Set, Rotary Pump, either Fixed Leads or Flying Leads with the Irrigator (Ref 9190001RP and 9190002RP).
- The Irrigator does not generate waste. Disposable products used with the Irrigator must be disposed of in accordance with hospital procedures and regulations.
- Portable and mobile RF communications equipment can affect electrical medical equipment.
- Refer to individual Instructions for Use for additional information regarding accessories or use with other devices.

For CODMAN Integrated Bipolar Cord and Tubing Set, Fixed Leads, Rotary Pump and CODMAN Integrated Bipolar Cord and Tubing Set, Fixed Leads, Gravity Feed:

Operating Temperature Range	+10° C to +40° C
Operating Humidity Range	20% to 80% relative humidity, non-condensing
Operating Pressure Range	700 hPa to 1050 hPa

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
- Silicone has a low cut and tear resistance. Do not allow the silicone tubing portion of the set, between the plastic fittings, to come in contact with sharp objects.
- Do not peel apart the tubing and electrical cord as this can cause the tubing portion to leak.
- This product is designed for use only with collapsible solution containers.
- This product is designed for use with bipolar electrosurgical units only.
- Do not use this bipolar cord set at generator power settings that exceed peak-to-peak RF voltage of 500 volts. Refer to the operator's manual supplied with the generator to determine the maximum allowable power settings.

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Generator Setup.


- 1) Place the Generator on a stand, shelf, or other flat surface capable of securely supporting the device.

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

- a. Plug the power cord into the rear panel of the Generator.
- b. Plug the Generator into a grounded electrical outlet of the appropriate voltage rating

Foot Pedal & Electrosurgical Irrigator Setup.

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

- a. Insert the plug of the Foot Pedal into the Foot Pedal Receptacle on the rear of the Generator. Place the Foot Pedal in an area accessible to the user.
 - b. Twist the cable collar 1/4 turn to lock the Foot Pedal cable in place.
- Note: Verify activation of the Foot Pedal before each use.
- c. To disconnect the Foot Pedal, twist the cable collar counterclockwise 1/4 turn to unlock the cable; grasp the plastic housing of the Foot Pedal cable and pull. Always grasp and pull the plastic housing; do not pull the cord.

2) Codman Electrosurgical Irrigator Setup (Optional):


If you are using the Electrosurgical Irrigator with the Generator:

- a. Follow the Instructions provided in the next point.
- b. Using the Electrosurgical Interconnecting Cable, connect one end of the cable to the Generator receptacle on the rear of the Irrigator. This cable allows you to control irrigation while coagulating when using the Foot Pedal to activate the system.
- c. Connect other end of the cable to the receptacle on the rear of the generator.

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Power on the Generator by pushing the power button. The power button will illuminate blue. An error in the start-up indicates a device defect. Refer to the Error Diagnosis/Troubleshooting section.

Self-Test:

The Generator performs a self-test after it is powered on. All display elements on the front panel are turned on during this process. The following light up briefly in the order they are listed: The MONOCUT and BICUT mode buttons, then MONOCOAG and BICOAG mode buttons, followed by the red indicator light above the neutral electrode receptacle. The two Output Power displays light up with "88" and the white CUT and COAG Indicator lights illuminate sequentially; the respective activation tone must sound simultaneously. The default start-up of the device is BICUT and BICOAG with both displays lit up with "00".

An error in the self-test sequence indicates a device defect. Refer to the Error Diagnosis/Troubleshooting section.

Codman Electrosurgical Irrigator Setup.

NOTE: To avoid unintentional irrigation, ensure that the Irrigator is switched off during set-up.

1. Place the Irrigator on an appropriate surface.
2. Insert the appropriate corresponding Interconnecting Cable into the receptacle on the rear panel of the Irrigator:
Codman Electrosurgical Interconnecting Cable (REF 901001CBL1)
Codman Electrosurgical Interconnecting Cable for CMC-V (REF 901001CBL2)
3. Insert the other end of the Interconnecting Cable into the receptacle on the rear panel of the Generator.

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4. Insert the Power Cord into the power receptacle on the rear panel of the Irrigator and into an appropriately grounded 100-240 VAC 50/60 Hz power supply.

WARNING: Do not attempt to bypass the grounding prong on the Power Cord by using an adapter. The Irrigator must be grounded properly to ensure operator and patient safety. Grounding reliability can be achieved only when connected to a power supply marked "Hospital Only" or "Hospital Grade."

NOTE: Always position the Irrigator so that the Power Cord can be easily unplugged from the rear panel.

5. Depress the power button located on the front panel to power on the Irrigator. The power button will illuminate. The screen above the flow rate knob will also illuminate and display the Codman logo and software version. It will then display the initial flow rate setting of "0."

NOTE: If the service icon is displayed, or the service indicator light is activated, contact your Integra sales representative.

NOTE: Ensure that the pump hood is closed prior to operation.

NOTE: To verify the above-mentioned operation of the Irrigator before loading the tubing set or connecting the forceps cable, you can do the following:

- Turn the flow rate knob clockwise until the display screen indicates a flow rate setting of "10."
- Depress the foot pedal connected to the Generator and verify that the pump head is operational by hearing an audible whirring sound. The pump head rollers should rotate when the foot pedal is depressed and stop rotating when the foot pedal is released.

If the Irrigator does not operate as noted above, check all connections (including Power Cord) and confirm that the unit is powered on and not displaying an indicator to service the unit.

Irrigation Operation:

1. Rotate the flow rate knob to select the desired irrigation flow rate setting. The Irrigator display will indicate the selected flow rate setting.

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NOTE: The irrigation flow rate setting can be increased by rotating the flow rate knob clockwise. Conversely, the irrigation flow rate setting can be decreased by rotating the knob counterclockwise.

NOTE: Always adjust the flow rate setting with the bipolar forceps attached to ensure that sufficient irrigation fluid is delivered into the operative field to meet clinical needs.

2. Depress and hold the foot pedal to deliver the desired irrigation fluid to the surgical field. Release the foot pedal to stop flow.

Additional Operating Features of Irrigator:

The irrigation flow rate can be adjusted while the foot pedal is depressed and irrigation fluid is delivered. Irrigation fluid delivery can be suspended by pressing the Pause button on the front of the Irrigator. The Pause button will illuminate and the display will indicate that flow has been suspended. Pressing the Pause button again will permit irrigation fluid delivery while the foot pedal is depressed. NOTE: The Irrigate Button on the Electrosurgical Generator will not control the suspension of irrigation fluid delivery. Flow suspension can only be controlled by the Pause button on the Irrigator.

Loading the CODMAN Integrated Bipolar Cord and Tubing Set, Fixed Leads, Rotary Pump:

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
NOTE: To ensure proper operation, use only the Codman Integrated Bipolar Cord and Tubing Set, Rotary Pump, either Flying Leads or Fixed Leads.

1. Using sterile technique, remove the Cord and Tubing Set from the peel pouch and place in the sterile field. Secure the male luer and the socket end in the sterile field.

2. Transfer the Generator plugs and irrigation inlet spike from the sterile field to permit subsequent connection to the Generator and Irrigator. Remove the protective cap from the irrigation tubing inlet spike. Prepare the irrigation solution bag to receive the spike. While holding the bag vertically, push the spike into the center of the port until secure. DO NOT TWIST.

3. Engage the irrigation tubing clamp.


4. Suspend the bag, port end down, on an IV pole.

Note: To ensure sufficient flow, position the irrigation container so the set's inlet spike is above the level of the forceps. Increasing the height of the irrigation solution container increases the pressure of the flow.

5. Identify the silicone tubing located between the two plastic fittings on the tubing set.

6. Open the irrigator module door and lift the pump hood.

7. Starting from the top, thread the tubing through the slot, around the rollers, and through the exit opening at the bottom of the pump head, as shown below:

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

The plastic fitting nearest the irrigation solution bag must be at the top of the pump head and the fitting nearest the forceps must be at the base.

8. To prime the system, release the irrigation tubing clamp and allow fluid to flow freely through the tubing. Remove as many air bubbles as possible and ensure there is fluid throughout the irrigation tubing. Close the pump hood to stop the fluid flow.

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

WARNING: Ensure that the silicone tubing is not pinched when closing the pump hood.

9. Insert the Cord and Tubing Set plug(s) into the Generator. The plug(s) must fit securely in the Generator's isolated bipolar output jack(s).

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10. Connect the forceps to the Cord and Tubing Set socket.
WARNING: Ensure the connection is completely seated. Failure to observe warning may result in electrical shock or injury to the patient or operating room personnel.

11. Connect the Cord and Tubing Set and forceps designed to provide irrigation by inserting the male luer connector into the female forceps connector.

12. Priming can be achieved by performing one of the two following techniques:

a. Turn the flow rate knob clockwise until the display screen indicates a flow rate setting of "20." Depress and hold the foot pedal until irrigation fluid is flowing freely through the forceps tips and there are no air bubbles visible.

b. Manually prime the system by releasing the irrigation tubing clamp and allowing irrigation fluid to flow through the tubing. Remove as many air bubbles as possible from the tubing before use. Verify flow in the system by observing irrigation through the forceps tips. When the tubing set is fully primed, close the pump hood to stop the fluid flow. Avoid pinching the silicone tubing when closing the pump hood.

Loading the CODMAN Integrated Bipolar Cord and Tubing Set, Fixed Leads, Gravity Feed:

1. Using aseptic technique, remove the CODMAN Integrated Bipolar Cord and Tubing Set from the peel pouch and place in the sterile field. Secure the male luer and the socket end in the sterile field.

2. Transfer the generator plugs and irrigation inlet spike from the sterile field to permit subsequent connection to the generator and irrigation module. Remove the protective cap from the irrigation tubing inlet spike. Prepare the irrigation solution bag to receive the spike. While holding the bag vertically push the spike into the center of the port until secure. **DO NOT TWIST.**

3. Engage the irrigation tubing clamp.

4. Suspend the bag, port end down, on an IV pole.

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
5. Insert the silicone tubing, located between the two plastic fittings on the tubing set, into the keyhole slots at the top and base of the irrigation module.

The plastic fitting nearest the irrigation solution container must be at the top of the irrigation module and the fitting nearest the forceps must be at the base, as shown below:

6. Press the load switch on the irrigation module to retract the occluding arm. With the other hand, position the silicone tubing so that it is occluded by the arm when the load button is released.

7. To prime the system, release the irrigation tubing clamp and press the load switch on the irrigation module. Remove as many air bubbles as possible and ensure there is fluid throughout the irrigation tubing before releasing the load switch.

8. Insert the two plugs into the generator. The plugs must fit securely in the generator's isolated bipolar output jack.

9. Connect the forceps to the cord set socket.

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
WARNING: Ensure the connection is completely seated. Failure to observe warning may result in electrical shock or injury to the patient or operating room personnel.

10. Turn on power to the irrigation module and the bipolar generator. Set the irrigation module and generator to allow for irrigation to be delivered through the tubing set.

11. Verify flow through the irrigation tubing by observing irrigation to the luer connector. Connect to the forceps' irrigator receptacle.

Disconnect the power cord:

To remove power cord from the device, unplug the power cord from the power receptacle or wall socket. Do not use extension cords with three-prong to two-prong power plug adaptors, or extra-length Power Cords with the Irrigator. Before each use, visually inspect the Power Cord and plug for frayed or broken insulation. If necessary, replace the Power Cord with the same type, length, gauge, and insulation

ISOCOOL tip and handle assembly:

CAUTION:

- Make sure the generator is turned off (or is in standby mode) before connecting the forceps handle or changing the tips. User injury could result if power is applied to the forceps during assembly.
- Handle the working end of the tips with care; they can be sharp and can cause injury.

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
1. Remove the tips from the sterile package. Leave the protective tubes in place on the tips to protect hands and gloves.
2. Hold one tip so the flat surface faces the midline of the handle, as shown below:

- A. Tip
- B. Handle
- C. Protective tube
- D. Flat surface


This will align the bend in the tip to the bend in the handle.

3. While depressing the tip release button on the handle blade, insert the shaft of the tip into the receptacle of the handle, as shown below:

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

A. Tip release button

4. Push the tip completely into the handle, as shown below:

A. Silicone seal


Release the tip release button.

5. Check the security of the tip by attempting to pull it out of the handle without depressing the tip release button.

6. Check to ensure that the distal seal is covering the junction of the tip and the handle. If it is folded back on itself, depress the tip release button, pull the tip out and reassemble.


7. Follow Steps 2 through 6 with the other tip.

8. Remove the protective tubes from the tips.

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9. Close the forceps and check the tip alignment. To adjust the tip alignment, grasp the tip and swivel it as needed. Do not depress the tip release button while performing this step. Use a piece of gauze to protect your gloves from the sharp edges of the tips.

10. Connect the handle to the bipolar cord (Disposable Bipolar Cord REF 80-1188) and the cord to the bipolar generator in the usual manner.

Disassemble the tips from the handle by depressing the tip release button on the handle while pulling the tips out. Use forceps, if desired, to assist with disassembly.

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Connecting the Neutral Electrode:


Connect the Neutral Electrode to the receptacle on the front of the device.

NOTE:

- If the neutral electrode is NOT connected during monopolar mode, the "Neutral electrode not connected" indicator above the receptacle will light up red and flash.
- If the finger switch or Foot Pedal is activated while in monopolar mode without the neutral electrode connected, an audible warning sounds. RF current cannot be activated. This has no influence on the bipolar current.
- When using a multi-surface electrode, the "Neutral electrode not connected" indicator will turn off after the safe application state is reached.
- The entire area of the neutral electrode should be reliably attached to a suitably prepared and appropriate area of the patient's body as defined by the manufacturer.

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
Connecting Monopolar Accessories:


1. Connect a 3-pronged connection monopolar cord or bovie connection monopolar cord to the appropriate location in the monopolar receptacle. Ensure the 3-pronged connection is in the correct orientation (see illustration below).
2. Connect the desired Monopolar Instruments to this cord.

NOTE:

Connect a Monopolar device with a finger switch for CUT and COAG activation **OR** Connect a Monopolar device without a finger switch, in combination with the Foot Pedal for CUT and COAG activation.

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
Connecting Bipolar Accessories:


1. Connect a Codman Bipolar cord to the receptacle.
2. Connect the desired bipolar instrument to the cord. The Foot Pedal is used for CUT and COAG activation.

NOTE:

- The neutral electrode is NOT required for bipolar applications.

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For Codman Electrosurgical Generator:

- Read entire instruction manual before attempting to operate the Generator.
- Use of monopolar operating modes should be restricted to extracranial cutting/coagulation during neurosurgical procedures (i.e., the opening and closing of neurosurgical cases) and to tissues contained in the extradural compartment. The monopolar operating modes should NOT be used in the cutting/coagulation of tissues in the subdural compartment or the neural parenchyma.
- Do not remove the Generator cover. There are no user serviceable components in this system.
- Replacement fuses must meet the specifications outlined in this manual. See Technical Specifications.

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- This Generator must not be modified in any way. Unauthorized modifications could cause the Generator to malfunction or fail and will void the warranty.
- The operator of the Generator must be a healthcare professional. The operator must be experienced in the use of electrosurgical devices and generators in surgical procedures.
- Failure of the Generator could result in unintentional increase of output power.
- The Generator performs a self-test upon initial start-up. An indicator error in any of the self-test sequence indicates a device defect. See the Error Diagnosis/Troubleshooting section to determine if the error can be resolved. If not, discontinue using the device and contact your Integra sales representative.

For Electrosurgical Irrigator:

- Use only Codman Integrated Bipolar Cord and Tubing Set manufacturer by Integra
- The Irrigator does not generate waste. Disposable products used with the Irrigator must be disposed of in accordance with hospital procedures and regulations.
- Portable and mobile RF communications equipment can affect electrical medical equipment.
- Refer to individual Instructions for Use for additional information regarding accessories or use with other devices.

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For CODMAN Integrated Bipolar Cord and Tubing Set, Fixed Leads, Rotary Pump and

For CODMAN Integrated Bipolar Cord and Tubing Set, Fixed Leads, Gravity Feed:

Failure to follow these warnings could result in fire or patient injury.

- This product contains a chemical (DEHP) known to the State of California to cause cancer, or birth defects, or other reproductive harm.
- Do not use in the presence of flammable anesthetics or oxidizing gases or in close proximity to volatile solvents, as explosion may occur.
- Do not use in patients who have cardiac pacemakers or other active implants without first consulting a qualified professional (e.g., cardiologist). A possible hazard exists because interference with the action of the electronic implant may occur, or the implant may be damaged.
- Do not place instruments directly on or near the patient when not in use. Place the active forceps and cord and tubing set on a visible, non-conductive area away from the patient or in a non-conductive instrument holster when not in use. Inadvertent contact with the patient may cause burns.
- Do not place active forceps and cord and tubing set near or in contact with flammable materials such as gauze or surgical drapes. Instruments that are activated or hot from use may cause a fire.
- Do not activate the forceps connected to the cord and tubing set when not in contact with target tissue, as this may cause injuries due to capacitive coupling with other surgical equipment.
- Do not exceed 500 Volts peak to peak in COAG mode. Failure to observe mode and setting warnings may reduce efficacy of non-stick performance, result in unintended burns, and/or cause failure of device insulations.
- Inspect instruments for damage prior to use. This may be done visually under magnification or with a high voltage insulation testing device. Insulation failures may result in burns or other injuries to the patient or operator.
- Ensure the bipolar electrosurgical generator is in the "off" position when connecting the cord and tubing set. Failure to observe warning may result in electrical shock or injury to the patient or operating room personnel.
- The surface of the active electrode may remain hot enough to cause burns after the RF current is deactivated.
- Select the minimum COAG output power possible to achieve coagulation.

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- Position electrosurgical cords to avoid contact with patient or other leads.
- Observe all additional warnings provided with the bipolar electrosurgical generator.
- Note Instructions for Use below for additional warnings.

For ISOCOOL bipolar forceps handle and tips:

- Electrocoagulation is a user-controlled surgical task that is based on clinical judgment and surgical context.
- When using instruments in electrosurgery, keep the power setting as low as possible to achieve the desired clinical effect.
- Activate the device only for as long as necessary to achieve the desired clinical effect.
- Thoroughly read the instructions for use that is packaged with the bipolar electrosurgical generator before connecting and using the ISOCOOL handle and tips.
- Do not apply bipolar electrosurgical current directly to staples or clips.
- Coagulation of tissue upon application of current that raises tissue temperature above 55-60° Celsius.

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- Activating the electrosurgical unit simultaneously with the aspiration/irrigation might alter the path of electrical energy away from the tissue.

For bipolar forceps:

- Thoroughly read the IFU that are packaged with the bipolar generator before attaching and using the forceps.
- It is recommended that the output power be set as low as possible for the intended purpose.
- Do not use Insulated Bipolar Forceps at generator power settings that exceed peak-to-peak RF voltage of 500 Volts.
- Non-insulated forceps do not provide protection against electrical current. They are intended only for use in low power electrosurgical procedures.
- Activating the electrosurgical unit simultaneously with aspiration/irrigation may alter the path of electrical energy away from the tissue.
- The use of non-insulated forceps in confined spaces (e.g., tonsillectomy) may result in unintended contact with and possible injury to non-target tissue during the application of power to the forceps. The use of insulated forceps is recommended for such procedures. The use of non-insulated forceps may also result in possible injury to the user during application of power to the forceps.
- To minimize coagulum deposits, use the lowest possible power setting for the intended use.
- Remove coagulum deposits as often as necessary to keep working surfaces clean.

For safe and convenient in-use removal of coagulated blood and debris from the tips of the bipolar forceps, Integra recommends the use of any nonabrasive cotton type material.

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Any serious incident that has occurred in relation to the device for the user and/or the patient should be reported to the manufacturer and the competent authority of the member state in which the user and/or patient is established.

The Generator conducts extensive functional tests during the power-on self-test. All safety-related functions are also checked continuously during operation. In case of errors, the RF energy output is stopped and the error is indicated by an error message.

Table below describes an example of an error diagnosis, possible cause, and solutions.

The messages where resolution by the user is possible are listed in the Possible Solutions column.

Note: After an error message is displayed, the device can only be reset by turning it off and back on again.

Note: Contact your manufacturer's authorized representative if a displayed message is not described below.

Error Messages	Possible Causes	Solutions
Apparent low output or failure of the HF surgical equipment to function correctly at the normal operating settings	Faulty application of the neutral electrode or poor contact in its connections.	Application of the neutral electrode and its connections should be checked before selecting a higher output power.
Visual elements on the front are dark; device does not function.	Power cord not plugged into the outlet	Plug power cord into the outlet
	Power cord is not fully connected to Generator	Assure that the power cord is completely plugged into both the device and the outlet

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


- Foot pedal: Connectors are undamaged; connector pins are not bent.
- Connect the Foot Pedal to the Generator and activate Foot Pedal to verify functionality.
- If the Foot Pedal does not comply with these conditions or has any other problems that might affect function, do not use the Foot Pedal.
 - Handle of Forceps the working end of the tips with care; they can be sharp and can cause injury.
- Disassemble the tips from the handle by depressing the tip release button on the handle while pulling the tips out. Use forceps, if desired, to assist with disassembly

Service and Repair:

Annual safety and functional checks are recommended. After this check (or after the product is out of order), if any damage / unfitness for use is found, the product should be handed over to an authorized representative of the manufacturer to assess the possibility of repair. Disposable items cannot be repaired.

The manufacturer will provide schematics, parts list, and other technical information to authorized representative. For service, repair or replacements, contact your local Integra representative. Equipment under warranty should be returned to Integra LifeSciences Production Corporation. It will be repaired or replaced without charge to the purchaser

Irrigation Tube of Bipolar Forceps Care:

1. Remove stylet from irrigation tube prior to inserting the irrigation tubing.
2. Upon completion of the procedure, flush the tube with nonsaline, sterile water. Use compressed air (using a dry syringe) to remove the excess moisture from the inside of the tube. The stylet should be inserted into place and remain there during handling and nonsterile storage.

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Medical electrical equipment needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided in the points below.

NOTE: Portable and mobile RF communications equipment can affect medical electrical equipment.

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For Generator:

- To avoid risk of electric shock, the Generator must only be connected to supply mains with protective ground.
- Disconnect the power cord to isolate the Generator from the supply mains. Serious shock hazards may be present even with the power switch off.
- Do not bypass the grounding prong on the Generator by using an adapter. The Generator must be properly grounded to ensure user and patient safety. Grounding reliability can only be achieved when connected to a receptacle marked "Hospital Only" or "Hospital Grade."
- The output from the Generator must not be applied directly to the heart.
- Do not use in patients who have electronic implants such as cardiac pacemakers without first consulting a qualified professional (e.g., cardiologist). A possible hazard exists because interference with the action of the electronic implant may occur, or the implant may be damaged.
- Do not connect multiple portable socket outlets (MPSO) or extension cords to the Generator.
- Connect adapters and accessories to the Generator only when the energy is off. Failure to do so may result in an injury or electrical shock to the patient or operating room personnel.
- Never simultaneously touch the patient and the external fuse.
- Inspect all cords, monopolar and bipolar accessory devices, and Foot Pedal for proper function and intact insulation prior to each use.

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- Use in typical commercial or hospital environment, except in RF (radiofrequency) shielded room for magnetic resonance imaging, where the intensity of EM (electromagnetic) disturbances is high.
- Avoid HF (high-frequency) output settings where maximum output voltage may exceed rated accessory voltage.
- Activating the RF current is not permitted while inserting the electrode and during electrode replacement due to risk of burns.
- Do not use in the presence of flammable anesthetics or oxidizing gases (such as nitrous oxide (N₂O) and oxygen) or in close proximity to volatile solvents (such as ether or alcohol), as explosion may occur.
- Do not immerse the Generator in any type of liquid.
- Non-flammable agents should be used for cleaning and disinfection whenever possible.
- Flammable agents, including disinfectants, cleaning fluids, and solvents, must fully evaporate prior to HF surgical procedures. There is a risk of pooling of flammable solutions under the patient or in body depressions such as the umbilicus, and in body cavities such as the vagina. Any fluid pooled in these areas should be mopped up before HF surgical equipment is used. Attention should be called to the danger of ignition of endogenous gases. Some materials, for example cotton and gauze, when saturated with oxygen may be ignited by sparks produced in normal use of the HF surgical equipment.
- Unless a compatible monitoring neutral electrode is used with a contact quality monitor, loss of safe contact between the neutral electrode and the patient will not result in an auditory alarm.
- Other equipment may be adversely influenced by the interference produced by the operation of HF surgical equipment. Portable and mobile RF communications equipment can affect medical electrical equipment.
- Prevent patient contact with metal parts that are grounded or that have an appreciable capacitance to ground; e.g., operating table supports. Use antistatic sheeting.
- Prevent patient skin-to-skin contact (e.g., between the arm and the body) when the Generator is in use. Use dry gauze between body parts that might touch during surgery.
- When HF surgical equipment and physiological monitoring equipment are used simultaneously on a patient, any monitoring electrodes must be placed as far as

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possible from the surgical electrodes. Needle monitoring electrodes are not recommended. If needle electrodes are used, do not place in direct contact with forceps tips. Monitoring systems incorporating high frequency current-limiting devices are recommended.

- During surgical procedures where the HF current could flow through parts of the body having a relatively small area, bipolar techniques can prevent unwanted tissue damage.

- To prevent injury to the patient and/or user, do not connect the Generator to any item not specified as part of the system or identified as a recommended accessory. Use of unspecified items with the system could result in excessive leakage current. Refer to individual Instructions for Use for additional information regarding accessories or use with other devices. The system should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the system should be observed to verify normal operation in the configuration in which it is used.

- Refer to the Guidance and Manufacturer's Declaration section for EMC (Electromagnetic Compatibility) information.

- In order for the Generator to be able to fulfill its intended clinical function effectively, a maximum cord and accessory length of 4.5 meters between the Generator and accessories is defined for Bipolar Accessories, 3.2 meters is defined for Monopolar Accessories and 3.0 meters for Neutral Electrodes.

- The use of accessories, transducers and/or cords other than those specified, except for those sold by the manufacturer as replacement parts for internal components, could result in increased emissions or decreased immunity of the equipment or system.

- The application of a current with high voltage, and in particular of a high voltage monopolar coagulation current, may cause neuromuscular stimulation in the patient.

- When vaporizing tissue, do not inhale the combustion fumes that form by design over a longer period in concentrated form. During regular use of the Generator, no pollutants other than these combustion products will be formed. For removal of such combustion products, a fume suction system can be used.

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- To prevent injury to the patient and/or user, do not connect the Irrigator to any item not specified as a recommended accessory.
- Use of unspecified items with the Irrigator could result in excessive leakage current.
- Do not immerse the Irrigator in any type of liquid.
- Check the Power Cord and Interconnecting Cable for proper function prior to use. Replace if necessary.
- The operator of this equipment must be a healthcare professional. The operator must have experience in the use of surgical devices and procedures.
- Use in typical commercial or hospital environment, except in RF shielded room for magnetic resonance imaging, where the intensity of EM DISTURBANCES is high.
- Refer to the Guidance and Manufacturer's Declaration section for EMC information.

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Emissions Test	IEC 60601-1-2, 3rd edition Compliance Level	IEC 60601-1-2, 4th edition Compliance Level	Electromagnetic Environment-Guidance
Radiated RF Emissions CISPR 11	Group 2	Group 2	The Generator must emit electromagnetic energy in order

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			to perform its intended function. Nearby electronic equipment may be affected.
Conducted RF Emissions CISPR 11	Class A	Class A	The Generator is suitable for use in all establishments, other than domestic and those directly connected to the public low voltage power supply network that supplies buildings used for domestic purposes.
Harmonics IEC 61000-3-2	Class A	Class A	
Emissions of voltage fluctuations/flicker IEC 61000-3-3	Complies	Complies	

NOTE: The Generator must be installed and operated according to the EMC regulations laid out in this document.

NOTE: The EMISSIONS characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 Class A). If it is used in a residential environment (for which CISPR 11 class B is normally required), this equipment might not offer adequate protection to radiofrequency communication services. The user might need to take mitigation measures, such as relocating or re-orienting the equipment.

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The Codman Irrigator is intended for use in the electromagnetic environment specified below. The customer or user of the Irrigator should ensure that it is used in such an environment.

Emissions Test	IEC 60601-1-2, 3rd edition Compliance Level	IEC 60601-1-2, 4th edition Compliance Level	Electromagnetic Environment-Guidance
RF Emissions CISPR 11	Group 1		The Codman Irrigator creates/uses electromagnetic energy for its internal functioning only; therefore, its electromagnetic interferences and disturbances are very low and are not likely to cause any interference in nearby electronic equipment.
RF Emissions CISPR 11	Class A		
Harmonics IEC 61000-3-2	Class A		The Codman Irrigator is suitable for use in all establishments, other than domestic and those directly connected to the public low voltage power supply network that supplies buildings used for domestic purposes.
Flicker IEC 61000-3-3	Complies		

WARNING: The Irrigator should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the System should be observed to verify normal operation in the configuration in which it is used.

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NOTE: The EMISSIONS characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 Class A). If it is used in a residential environment (for which CISPR 11 class B is normally required) this equipment might not offer adequate protection to radio-frequency communication services. The user might need to take mitigation measures, such as relocating or re-orienting the equipment.”

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Guidance and Manufacturer's Declaration—Immunity All ME Equipment and ME Systems:

The Generator is intended for use in the electromagnetic environment specified below. The customer or user of the Generator should ensure that it is used in such an environment.

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The **Codman Irrigator** is intended for use in the electromagnetic environment specified below. The customer or user of the Irrigator should ensure that it is used in such an environment.

Immunity Test	IEC 60601-1-2, 3rd edition Compliance Level	IEC 60601-1-2, 4th edition Compliance Level	Electromagnetic Environment-Guidance
ESD IEC 61000-4-2	±6 kV Contact ±8 kV Air	±8 kV Contact ±15 kV Air	Floors should be wood, concrete or ceramic tile. If floors are synthetic, the relative humidity should be at least 30%.
EFT IEC	±2 kV Mains ±1 kV I/Os	±2 kV Mains ±1 kV I/Os	Mains power quality should be that of a typical commercial or

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61000-4-4			hospital environment.
Surge IEC 61000-4-5	± 1 kV Differential ± 2 kV Common	± 1 kV Differential ± 2 kV Common	Mains power quality should be that of a typical commercial or hospital environment.
Voltage Dips/ Dropout IEC 61000-4-11	<5% UT for 0.5 Cycle 40% UT for 5 Cycles 70% UT for 25 Cycles <5% UT for 5 Seconds	0% UT for 0.5 Cycle 0% UT for 1 Cycle 70% UT for 25/30 Cycles 0% UT for 250/300 Cycles	Mains power quality should be that of a typical commercial or hospital environment. If the user of the Irrigator requires continued operation during power mains interruptions, it is recommended that the Irrigator be powered from an uninterruptible power supply or battery.
Power Frequency Magnetic Field IEC 61000-4-8	3 A/m, 50/60 Hz	30 A/m, 60 Hz	Power frequency magnetic fields should be that of a typical commercial or hospital environment.

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Recommended Separation Distances for the Generator, per IEC 60601-1-2, 3rd edition:

The Generator is intended for use in the electromagnetic environment specified below. The customer or user of the Generator should ensure that it is used in such an environment.

Maximum Output Power (Watts)	Separation (m) 150 kHz to 80 MHz $d = 1.2\sqrt{P}$	Separation (m) 80 to 800 MHz $d = 1.2\sqrt{P}$	Separation (m) 800 MHz to 2.5 GHz $d = 2.3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

Recommended separation distances between portable and mobile RF telecommunications equipment and the Generator.

Requirements for high-frequency wireless communication equipment

Frequency band (MHz)	Test frequency (MHz)	Modulation	Compliance level (V/m)
380-390	385	Pulse ^a -18 Hz	27
430-470	450	FM ± 5 kHz Hub or Pulse ^a -18 Hz	28

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704-787	710, 745, 780	Pulse ^a –217 Hz	9
800-960	810, 870, 930	Pulse ^a –18 Hz	28
1,700-1,990	1,720, 1,845, 1,970	Pulse ^a –217 Hz	28
2,400-2,570	2,450	Pulse ^a –217 Hz	28
5,100-5,800	5,240, 5,500, 5,785	Pulse ^a –217 Hz	9
Note	A minimum safety distance of 30 cm (12") between portable HF communication equipment emitting in the respective frequency band and the Generator should be observed. This includes e.g. cell phones and WIFI, RFID and Bluetooth devices. Non-observance may result in reduction of the performance features of the device.		
a) Pulse modulation is defined as a square wave signal with 50% duty cycle.			

WARNING: Per IEC 60601-1-2, 4th edition, portable RF communications should be used no closer than 30 cm (12 inches) to any part of the Generator, including cords specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

Recommended Separation Distances between portable and mobile RF communications equipment and the Codman Irrigation System, ME Equipment and ME Systems that are not life supporting:

Recommended Separation Distances for the Codman Irrigator, per IEC 60601-1-2, 3rd edition:

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The Codman Irrigator is intended for use in the electromagnetic environment in which radiated disturbances are controlled. The customer or user of the Irrigator can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment and the Irrigator as recommended below, according to the maximum output power of the communications equipment.

Maximum Output Power (Watts)	Separation (m) 150 kHz to 80 MHz $d = 1,2\sqrt{P}$	Separation (m) 80 to 800 MHz $d = 1,2\sqrt{P}$	Separation (m) 800 MHz to 2.5 GHz $d = 2,3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

WARNING: Per IEC 60601-1-2, 4th edition, portable RF communications should be used no closer than 30 cm (12 inches) to any part of the Irrigator, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

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RF Output:

The below power charts show typical output characteristics of the Generator for all available modes. For each mode, the power output at the maximum (1) and half (2) power level as a function of the resistive load is depicted.

Actual measurements of the output power may vary depending on the chosen test accessory and the measurement setup. According to IEC 60601-2-2:2017

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
paragraph 201.12.1.101, the power output shall not deviate by more than ± 20 %.

Monopolar:

Curve 1: Output setting "maximum output"

Curve 2: Output setting "half output"

Curve 1: Output setting "maximum output"

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




Curve 1: Output setting "maximum output"
 Curve 2: Output setting "half output"

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RF Voltage:



The below circuit voltage charts show typical output characteristics of the Generator for all available modes. For each mode, the peak voltage (Vp) output as a function of the selected power level is depicted.

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





Control Characteristic Curve:

The control characteristics show that the screen settings are equivalent to the actual power output in all modes (no gear ratio).

The vertical lines, labeled with their respective mode, indicate the maximum power output of each mode.

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







Malis Units to Watts Conversion:

Values below are approximate and the actual conversion rate may vary.

BIPOLAR COAG Tested at a 50 Ω Load		BIPOLAR CUT Tested at a 400 Ω Load	
Output Power Setting (Watts)	Malis Units	Output Power Setting (Watts)	Malis Units
0	0	0	0
1	10	1	10
3	15	2	15
5	20	4	20

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Do not allow the devices to remain in environmental extremes. After exposure to environmental extremes, such as those found in transport or storage, allow the generator and irrigator to come to room temperature before operating.
Store sterile packaged instruments in a manner that provides protection from dust, moisture, insects, vermin, and extremes of temperature and humidity.

For Generator, Bipolar Cords, Foot pedal, Bipolar Forceps, ISOCOOL tips and handle:

Transport & Storage Temperature Range	- 40° C to + 60° C
Transport & Storage Humidity Range	15% to 85% relative humidity, non-condensing
Transport & Storage Pressure Range	500 hPa to 1050 hPa

For Irrigator and cables:

Transport & Storage Temperature Range	- 40° C to + 50° C
Transport & Storage Humidity Range	20% to 85% relative humidity, non-condensing
Transport & Storage Pressure Range	500 hPa to 1050 hPa

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The Generator, Foot Pedal, Irrigator, Interconnecting Cable and Power Cord are not distributed as sterile devices and are not intended to be sterilized at any time.

CODMAN Integrated Bipolar Cord and Tubing Set Flying Leads, Gravity Feed and CODMAN Integrated Bipolar Cord and Tubing Set Flying Leads, Rotary Pump are supply sterile. This product is intended for SINGLE USE ONLY; DO NOT RESTERILIZE. Use aseptic technique in all phases of handling. Integra Single Use devices have not been designed to undergo or withstand any form of alteration, such as disassembly, cleaning or resterilization, after a single patient use. Reuse can potentially compromise device performance and any usage beyond the design intent of this single use device may result in unpredictable loss of functionality. Integra will not be responsible for any product that is resterilized, nor will we accept for credit or exchange any product that has been opened but not used. As long as the individual package of the product is not opened or damaged, the product is sterile with a nonpyrogenic fluid path.

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The Codman Integrated Bipolar Cord and Tubing Sets are sterilized via ethylene oxide. The Cord and Tubing Sets have been validated to ensure a sterility assurance level (SAL) of 10^{-6} in accordance with ISO 11135:2014 Sterilization of health care products – Ethylene Oxide - Requirements for the Development, Validation and Routine Control of a Sterilization Process for Medical Devices and ISO 11737-1:2006, “Sterilization of medical devices - Microbiological methods - Part 1: Determination of a population of microorganisms on products.” Sterilization testing also demonstrated that ethylene oxide residuals can be reduced to an acceptable level in accordance with ISO 10993-7:2008 Biological Evaluation of Medical Devices – Ethylene Oxide Sterilization Residuals and the proposed device can be successfully adopted into Codman’s existing sterilization cycle.

Disposable Bipolar Cord is supply sterile. This product is for SINGLE USE ONLY; DO NOT RESTERILIZE. Integra LifeSciences will not be responsible for product that is resterilized, nor accept for credit or exchange product that has been opened but not used. Integra Single Use devices have not been designed to undergo or withstand any form of alteration, such as disassembly, cleaning or resterilization, after a single patient use. Reuse can potentially compromise device performance and any usage beyond the design intent of this single use device may result in unpredictable loss of functionality. As long as the inner unit is not opened or damaged, the product is sterile.

Reusable Bipolar Cord, Unitized Plug is supply non-sterile and must be cleaned and sterilized before use, including initial use.

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WARNING: Under certain classifications of risk, the World Health Organization (WHO), or local regulatory authorities recommend special CJD

inactivation processing procedures. Consult site WHO and local regulations for further information.

Usual cleaning and sterilization could not be enough.

ISOCOOL Handle, Bipolar Forceps and Reusable bipolar cord, Unitized Plug supply non-sterile and must be cleaned and sterilized before use (including the first use).

Cleaning for Reusable bipolar cord:

After each use, wipe the bipolar cord with a damp, soft cloth

Cleaning for Bipolar Forceps:

1. Clean separately from other instruments.
2. Prepare an enzymatic cleaning solution in accordance to the manufacturer's instructions.
3. Soak soiled instruments for 5 minutes in the protein solubilizing detergent.
4. Brush with a soft non-metallic bristle brush to remove all traces of blood and debris.
5. Rinse the instrument thoroughly with warm tap water.
6. Ultrasonically clean the instrument for 10 minutes in neutral pH detergent, prepared in accordance with the manufacturer's instructions.
7. Rinse the instrument thoroughly with warm tap water.
8. Dry the instrument immediately after final rinse.

Warning: For effective sterilization of irrigating type forceps, remove the stylet from the microtube.

Caution: Repeated exposure to hydrogen peroxide sterilization will cause the instrument to deteriorate.

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Cleaning for ISOCOOL Handle:

1. Remove the disposable tips (if necessary) and inspect the handle for damage.
2. Cracked or bent instruments should be set aside and sent for service or repair.
3. Submerge handle in an enzymatic detergent safe for use with metal instruments. Prepare the detergent solution per the manufacturer's instructions.
4. Soak the soiled handle for five (5) minutes in the protein solubilizing detergent.
5. Scrub the submerged handle with a soft bristle brush.
6. Agitate the handle in the solution while scrubbing. Actuate the device.
7. Rinse the handle with warm (38 – 49 ° C) tap water. Place the handle into a bath containing warm (38 – 49 ° C) tap water. Agitate the handle by hand for at least one (1) minute.
8. Repeat this rinse process two (2) additional times.
9. Ultrasonically clean the handle for ten (10) minutes in a neutral pH detergent. Prepare the detergent solution per the manufacturer's instructions.
10. Rinse the handle with warm (38 – 49 ° C) tap water for at least one (1) minute.
11. Dry the exterior of the handle with clean, lint free cloth.

Cleaning Inspection:

- Inspect all instruments before sterilization or storage to ensure the complete removal of soil from surfaces, tubes and holes, moveable parts.
- If areas are difficult to inspect visually, check for blood by immersing or flushing the instrument in a 3% hydrogen peroxide solution. If bubbling is observed, blood is present. Rinse instruments thoroughly after using hydrogen peroxide solution.
- If soil is still present, reclean the instrument.

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Sterilization:

These devices are provided nonsterile and must be cleaned and sterilized before use, including initial use.

In a properly functioning, calibrated autoclave steam sterilizer, effective sterilization can be achieved with the following parameters:

Method	Cycle type	Temperature	Exposure time
Steam sterilization	Prevacuum	132–134°C	4 minutes

CAUTION: Do not sterilize the handle by hydrogen peroxide. The handle design is not compatible with this method of sterilization. Repeated use of hydrogen peroxide sterilization will cause the handle to deteriorate and eventually to fail.

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Medical device and accessories are manufactured in accordance with the Therapeutic Goods Act 1989 and the Therapeutics Goods (Medical Devices) Regulations 2002, Essential Requirements of the Medical Device Directive 2007/47/EC effective March 21, 2010 amending the directive 93/42/EEC, Medical Devices Regulations SOR/98-282, United States Code of Federal Regulations 21 CFR Part 820 Quality System Regulation for medical devices and to the following International Standards:

Standard Number	Year	Title
EN ISO 13485	2016	Medical devices - Quality management systems - Requirements for regulatory purposes

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



Symbol	Title
	Alternating Current
	Authorized Representative in the European Community
	Catalog Number
	Serial Number
	Caution
	Consult Instructions for Use/Operating Instructions

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





	
NONSTERILE	Nonsterile
STERILE EO	Sterilized by Ethylene Oxide
STERILE R	Sterilized by gamma radiation
	Do not reuse
	Do not re-sterilize
	Contains Phthalate: bis (2-ethylhexil) phthalate (DEHP)
MADE IN	Made in
QTY	Quantity
IPX1	Protection against vertically falling water drops
	Monopolar Connection
	Bipolar Connection
	Neutral Electrode Indicator Light

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	Bipolar Tip
	Low profile Tip
	STERILE only if package is undamaged
	Dangerous Voltage
	Non-ionizing Electromagnetic Radiation
	Irrigator Connection
	Electrosurgical Generator
	Service Indicator
	Do not reach in while pump is in operating mode
	To identify the control or the indicator which stops operation intermittently and keeps the equipment in operating mode
IP2X	Protection against solid foreign objects of 12.5 mm diameter and greater

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Medical device consists of electrical components. Consider recycling at the end of service life. The product can be returned to manufacturer for recycling. Dispose of the product in accordance with local ordinances.

Class of waste:

Directive 2012/19/EU – Waste Electrical and Electronic Equipment (WEEE directive)

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The Codman Electrosurgical Generator with accessories - is considered an EEE and the requirements outlined in this directive apply. Integra LifeSciences Production Corporation has suitable processes and means in place to ensure that medical device is disposed of and treated appropriately and customer is informed accordingly. To show that this device is compliant with the directive, the symbol as outlined in Annex IX of the directive is displayed in its labeling as well as the indication that the system must be disposed of in accordance with regional hospital protocols and in accordance with WEEE regulations.

The ISOCOOL tips on its own is not considered an EEE as defined by the directive, hence the device under consideration does not have to comply with the applicable requirements. Dispose of the tips for forceps in accordance with the facility's guidelines for potentially biohazardous sharps materials.

Disposable materials and materials with an intrinsic biological risk must be disposed of in accordance with current local regulations governing hospital waste

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The warranty period of storage for Generator and Irrigator (in accordance with GOST R 50444-2020 at least 3 months) is 24 months and exploitation is 24 months

The warranty period of storage for following devices (in accordance with GOST R 50444-2020 at least 3 months) is 12 months and exploitation is 12 months: Codman Electrosurgical System Interconnecting Cable, Codman Electrosurgical System Foot pedal, Reusable Bipolar Cord, Unitized Plug, Disposable Bipolar Cord, Unitized Plug, Codman Integrated Bipolar Cord and Tubing Set, Fixed Leads, Gravity Feed, Codman Integrated Bipolar Cord and Tubing Set, Fixed Leads, Rotary Pump, ISOCOOL bipolar forceps tips, ISOCOOL bipolar forceps handle, Bipolar forceps.

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PRODUCT INFORMATION DISCLOSURE:

INTEGRA LIFESCIENCES CORPORATION ("INTEGRA") HAS EXERCISED REASONABLE CARE IN THE SELECTION OF MATERIALS AND THE MANUFACTURE OF THESE PRODUCTS. INTEGRA WARRANTS THAT THESE PRODUCTS SHALL CONFORM TO THE PRODUCT LIMITED WARRANTY AS PROVIDED IN THE PRODUCT LABELING OR APPLICABLE PRODUCT CATALOG. THIS WARRANTY IS EXCLUSIVE, AND INTEGRA DISCLAIMS ALL OTHER WARRANTIES, WHETHER EXPRESSED OR IMPLIED, INCLUDING BUT NOT LIMITED TO, ANY IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. INTEGRA SHALL NOT BE LIABLE FOR ANY INCIDENTAL OR CONSEQUENTIAL LOSS, DAMAGE, OR EXPENSE, DIRECTLY OR INDIRECTLY ARISING FROM USE OF THESE PRODUCTS. INTEGRA NEITHER ASSUMES NOR AUTHORIZES ANY PERSON TO ASSUME ANY OTHER OR ADDITIONAL LIABILITY OR RESPONSIBILITY IN CONNECTION WITH THESE PRODUCTS.

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Shelf life is not applicable for the **Codman Electrosurgical Generator, Irrigator Codman Electrosurgical Foot Pedal, Reusable Bipolar Cord, Unitized Plug and Interconnecting Cable** as they are non-sterile, reusable devices. Annual safety and functional checks are recommended.

The product lifetime of the **Codman Electrosurgical Generator** and **Codman Electrosurgical Foot Pedal** are eight (8) years. The product lifetime of the **Irrigator** are seven (7) years.

Disposable Bipolar Cord, Unitized Plug – shelf life is 5 years from the date of manufacturing.

CODMAN Integrated Bipolar Cord and Tubing Set, Fixed Leads, Rotary Pump - shelf life is 3 years from the date of manufacturing.

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CODMAN Integrated Bipolar Cord and Tubing Set, Fixed Leads, Gravity Feed - shelf life is 3 years from the date of manufacturing.

ISOCOOL bipolar forceps tips - shelf life is 5 years from the date of manufacturing. The lifetime for the ISOCOOL tips starts when the sterile barrier packaging has been opened and ends with the proper disposal of the tips after the procedure.

ISOCOOL bipolar forceps handle - shelf life is not applicable. The lifetime of the handles have been validated to withstand 20 reprocessing cycles without a loss of performance.

Bipolar Forceps are intended for short term use during surgical procedures. Procedural times may vary for numerous reasons such as the medical condition to be treated as well as surgical approach and techniques. The lifetime of the forceps has not been defined at this time. The time between reprocessing cycles may vary based on the frequency of use often determined by case load of the facility where the devices are in use, as well as surgeon preferences. Cleaning, sterilization and handling may impact the useful life of the devices if the provided manufacturer's instructions for use are not followed.

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For additional information related to the use of the medical product, please contact an authorized representative of the manufacturer:
 "Neuroproject LLC",
 127486, Moscow, Korovinskoe shosse, 10, bldg. 2, floor 1, of.101
 tel. +7(495) 660 81 73
 www. neurop.ru

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