Abbott

Abbott Ireland Diagnostics
Division,
Lisnamuck,
Longford
Co. Longford
Ireland

Tel. +353 43 3331000
Fax +353 43 3331001

I, Nora Fitzgerald, Specialist Regulatory Affairs, Abbott Ireland Diagnostics Division, Longford, Ireland, hereby certify that the attached documents that keep the number of pages indicated below are the Insert for Use. These documents are in a full accordance with actual manuals of manufacturer.

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Signed: _____

Nora Fitzgerald
Specialist Regulatory Affairs
Abbott Ireland Diagnostics Division
Lisnamuck, Longford
Co Longford
Ireland

Date: _____

18 - Dec - 19

ABBOTT DIAGNOSTICS DIVISION LONGFORD

SLIGO CHAMBER OF COMMERCE (INC)
CERTIFIED CORRECT IN ALL GOOD FAITH

Mary Harty 19/12/2019 **Authorised Official**

ARCHITECT

STAT High Sensitive Troponin-I Controls

REF 3P25-11

en

STAT High Sensitive Troponin-I
3P25
G56529R02
C3P2U0

Read Highlighted Changes: Revised March 2017.

INTENDED USE

The ARCHITECT STAT High Sensitive Troponin-I Controls are for the estimation of test precision and the detection of systematic analytical deviations of the ARCHITECT I System with STAT protocol capability when used for the quantitative determination of cardiac troponin I (cTnI) in human serum and plasma with the ARCHITECT STAT High Sensitive Troponin-I reagent kit.

Refer to the ARCHITECT STAT High Sensitive Troponin-I reagent package insert and the ARCHITECT System Operations Manual for additional information.

CONTENTS

3 Bottles (8.0 mL each) of ARCHITECT STAT High Sensitive Troponin-I Controls. Low Control, Medium Control, and High Control contain recombinant human cardiac troponin IC complex in phosphate buffer with protein (bovine) stabilizer.

Preservative: ProClin 300.

The controls are at the following concentrations:

Controls	Concentration (pg/mL)	Range (pg/mL)
CONTROL L	20	12.0 - 28.0
CONTROL M	200	120.0 - 280.0
CONTROL H	15,000	9,000.0 - 21,000.0

Each laboratory should establish its own concentration ranges for new control lots at each control level.

Note: The insert ranges for the controls are not lot specific and represent the total range of values which may be generated throughout the life of the product. It is recommended that each laboratory establish its own means and acceptable ranges which should fall within the package insert ranges. Sources of variation that can be expected include:

- Calibration
- Calibrator lot
- Control lot
- Reagent lot
- Instrument

PRECAUTIONS

- **IVD**
- For In Vitro Diagnostic Use

The following warnings and precautions apply to: CONTROL L / CONTROL M / CONTROL H

WARNING	Contains methylisothiazolones, H317
	May cause an allergic skin reaction.
Prevention	
P261	Avoid breathing mist / vapors / spray.
P272	Contaminated work clothing should not be allowed out of the workplace.
P280	Wear protective gloves / protective clothing / eye protection.

Response	
P302+P352	IF ON SKIN: Wash with plenty of water.
P333+P313	If skin irritation or rash occurs: Get medical advice / attention.
P362+P364	Take off contaminated clothing and wash it before reuse.
Disposal	
P501	Dispose of contents / container in accordance with local regulations.

Safety Data Sheets are available at www.abbottdiagnostics.com or contact your local representative.


For a detailed discussion of safety precautions during system operation, refer to the ARCHITECT System Operations Manual, Section 8.

Handling Precautions

- **Caution:** Due to the sensitivity of the ARCHITECT STAT High Sensitive Troponin-I assay, use caution when handling control bottles and caps to prevent cross-contamination.

STORAGE

- Controls are shipped frozen and must be stored at -10°C or colder upon receipt until first usage.
- Store thawed controls at 2-8°C. Do not refreeze controls after thaw.
- Controls are stable until the expiration date when stored and handled as directed.
- Do not use past expiration date.

PREPARATION FOR USE

- Remove from carton and allow controls to stand at room temperature (15-30°C) until completely thawed (approximately 120 to 150 minutes).
- Thawed controls should be thoroughly mixed by gentle inversion before use.
- After each use, immediately close the caps tightly and return the thawed controls to 2-8°C.
- Thawed controls may be stored refrigerated at 2-8°C for up to 3 months, not to exceed the expiration date printed on the bottle.
- To obtain the recommended volume requirements for the controls, hold the bottles vertically and dispense 10 drops of each control into each respective sample cup.
- For information on ordering controls, refer to the ARCHITECT System Operations Manual, Section 5.

QUALITY CONTROL PROCEDURES

The recommended control requirement for the assay is that a single sample of each control level be tested once every 24 hours each day of use. If your laboratory quality control procedures require more frequent use of controls to verify test results, follow those procedures.

Additional controls may be tested in conformance with local, state, and/or country regulations or accreditation requirements and your laboratory's quality control policy.

The control values must be within the acceptable ranges specified in the control package insert. If a control is out of its specified range, the associated test results are invalid and must be retested. Recalibration may be indicated.

Refer to the ARCHITECT STAT High Sensitive Troponin-I assay reagent package insert and ARCHITECT System Operations Manual for additional information.

Key to Symbols

	Consult instructions for use
	Manufacturer
	Temperature limitation
	Upper limit of temperature
	Use by/Expiration date
AFTER THAW	After thaw store at
CONC	Concentration
CONTROL L	Control Low, Medium, High (L,M,H)
IVD	In Vitro Diagnostic Medical Device
LOT	Lot Number
PRODUCT OF IRELAND	Product of Ireland
RANGE	Range
REF	List Number
UNTIL FIRST USE	Until first use store at
WARNING: SENSITIZER	Warning: May cause an allergic reaction.

ARCHITECT is a trademark of Abbott Laboratories in various jurisdictions. All other trademarks are property of their respective owners.

Abbott Ireland
Diagnostics Division
Lisnamuck, Longford
Co. Longford
Ireland
+353-43-3331000

Customer Service: Contact your local representative
or find country-specific contact information on
www.abbottdiagnostics.com

Revised March 2017.
©2015, 2017 Abbott Laboratories

ADDENDUM

to the Instruction for Use ARCHITECT STAT High Sensitive Troponin I Controls (for Russia)

SHIPPING AND STORAGE CONDITIONS

Recommended shipping conditions for device is - 10°C or colder in single-use thermo containers containing coolant elements, or in reusable thermo containers with automatically maintained temperature or in autorefrigerators using thermal indicators. Device is stable until the expiration date indicated on the package. Device is no longer suitable for use if shipping temperature limits have been violated. Ship in an enclosed carrier vehicle of any kind in accordance with shipping rules, applicable to the respective carrier.

Medical device must be stored at a temperature - 10°C or colder in refrigerators/refrigerator chambers to ensure regulated temperature conditions with a daily temperature record.

WARRANTY

Manufacturer guarantees that device remains stable until the expiration date when stored and handled as directed. The maximum shelf life of the medical device is 6 months from date of manufacturing. Device must be discarded if stored with deviations in manufacturer's requirements.

MAINTENANCE

Not required.

REPAIRS

Not required.

DESTRUCTION AND DISPOSAL OF A MEDICAL DEVICE

Destruct and dispose of a medical device in accordance with the requirements of SanPin 2.1.7.2790-10 dated December 12, 2010 "Sanitary-epidemiological Requirements for the Management of Medical Waste".

Transmit and dispose of a product and its components and packaging in accordance with local regulations.

LIST OF APPLICABLE STANDARDS

ISO 9001, ISO 13485, ISO 13485 (including CMDCAS)

POTENTIAL USER

ARCHITECT STAT High Sensitive Troponin I Controls are intended to be used in a variety of clinical and laboratory conditions, including, without limitation, hospitals, diagnostic centers, reference laboratories and blood transfusion facilities around the world, whose personnel have been trained in the use of the ARCHITECT i analyzer.

INFORMATION ABOUT PRESENCE OF DRUGS IN COMPOSITION OF MEDICAL DEVICE

Does not contain any drugs

STERILIZATION

Medical device is not required to be sterile.

INFORMATION ABOUT POSSIBILITY OF MULTIPLE APPLICATION

Medical device cannot be re-applied.

MATERIALS OF ANIMAL ORIGIN

Materials of animal origin are tested and risks of using MD are indicated in SDS. Safety Data Sheets are available at www.abbottdiagnostics.com or contact your local representative.

INDICATIONS AND CONTRA-INDICATIONS FOR USE

Not applicable.

**ARCHITECT****STAT High Sensitive Troponin-I Controls**

REF 3P25-11

STAT High Sensitive Troponin-I
3P25
G54998R02
C3P2UR

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Mary Harte 19/12/2019 Authorised Official

Read Highlighted Changes: Revised June 2017.

INTENDED USE

The ARCHITECT STAT High Sensitive Troponin-I Calibrators are for the calibration of the ARCHITECT i System with STAT protocol capability when used for the quantitative determination of cardiac troponin I (cTnI) in human serum and plasma with the ARCHITECT STAT High Sensitive Troponin-I reagent kit.

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CONTENTS

6 Bottles (4.0 mL each) of ARCHITECT STAT High Sensitive Troponin-I Calibrators. Calibrator A contains phosphate buffer with protein (bovine) stabilizer. Calibrators B - F contain different concentrations of recombinant human cardiac troponin IC complex in phosphate buffer with protein (bovine) stabilizer.

Preservative: ProClin 300.

The calibrators are at the following concentrations:

Calibrators	Concentration (pg/mL)
CAL A	0
CAL B	20
CAL C	100
CAL D	500
CAL E	10,000
CAL F	50,000

STANDARDIZATION

The calibrators are manufactured gravimetrically using recombinant human cardiac troponin IC complex. The calibrators are matched to internal reference standards prepared using recombinant human cardiac troponin IC complex. Calibrator B is traceable to NIST SRM 292 at 20 pg/mL. Calibrators may not be commutable with other assays.

PRECAUTIONS

- For In Vitro Diagnostic Use

The following warnings and precautions apply to: CAL A - CAL F	
WARNING	Contains methylisothiazolones.
NOT	May cause an allergic skin reaction.
Prevention	
200	Avoid breathing mist / vapors / spray.
202	Contaminated work clothing should not be allowed out of the workplace.
200	Wear protective gloves / protective clothing / eye protection.

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Safety Data Sheets are available at www.abbottiagnostics.com or contact your local representative.

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Handling Precautions

- Caution: Due to the sensitivity of the ARCHITECT STAT High Sensitive Troponin-I assay, use caution when handling calibrator bottles and caps to prevent cross-contamination.

STORAGE

- Calibrators are shipped frozen and must be stored at -10°C or colder upon receipt.
- Calibrators are stable until the expiration date when stored and handled as directed.
- Do not use past expiration date.
- Calibrators must be stored at -10°C or colder.

PREPARATION FOR USE

- Remove from carton and allow calibrators to stand at room temperature (15-30°C) until completely thawed (approximately 90 to 120 minutes).
- Thawed calibrators should be thoroughly mixed by gentle inversion before use and used immediately.
- After each use, immediately close the caps tightly and return the calibrators to -10°C or colder. Incorrect handling of the calibrators could result in inaccurate calibration of the instrument.
- Discard calibrators after 3 freeze-thaw cycles, regardless of remaining expiration.
- To obtain the recommended volume requirements for the calibrators, hold the bottles vertically and dispense 15 drops of each calibrator into each respective sample cup.
- To perform a calibration, test the calibrators in duplicate. The calibrators should be priority loaded.
- For information on ordering calibrations, refer to the ARCHITECT System Operations Manual, Section 6.

QUALITY CONTROL PROCEDURES


A single sample of each control level must be tested to evaluate the assay calibration. For information on ordering controls, refer to the ARCHITECT System Operations Manual, Section 5.

- Ensure that assay control values are within the ranges specified in the control package insert.

Once an ARCHITECT STAT High Sensitive Troponin-I calibration is accepted and stored, all subsequent samples may be tested without further calibration unless:

- A reagent kit with a new lot number is used.
- Controls are out of range.

Key to Symbols

	Consult instructions for use
	Manufacturer
	Upper limit of temperature
	Use by/Expiration date
	Calibrator (A,B,C,D,E or F)
	In Vitro Diagnostic Medical Device
	Lot Number
	Product of Ireland
	List Number
	Warning: May cause an allergic reaction.

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Revised June 2017.
©2015, 2017 Abbott Laboratories

ADDENDUM

to the Instruction for Use ARCHITECT STAT High Sensitive Troponin I Calibrators (for Russia)

SHIPPING AND STORAGE CONDITIONS

Recommended shipping conditions for device is - 10°C or colder in single-use thermo containers containing coolant elements, or in reusable thermo containers with automatically maintained temperature or in auto refrigerators using thermal indicators. Device is stable until the expiration date indicated on the package. Device is no longer suitable for use if shipping temperature limits have been violated. Ship in an enclosed carrier vehicle of any kind in accordance with shipping rules, applicable to the respective carrier.

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WARRANTY

Manufacturer guarantees that device remains stable until the expiration date when stored and handled as directed. The maximum shelf life of the medical device is 9 months from date of manufacturing. Device must be discarded if stored with deviations in manufacturer's requirements.

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ISO 9001, ISO 13485, ISO 13485 (including CMDCAS)

POTENTIAL USER

ARCHITECT STAT High Sensitive Troponin I Calibrators are intended to be used in a variety of clinical and laboratory conditions, including, without limitation, hospitals, diagnostic centers, reference laboratories and blood transfusion facilities around the world, whose personnel have been trained in the use of the ARCHITECT i analyzer.

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