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Director, Quality Make, Mitek, Spine & Trauma platform Lead

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
Senior Director, DePuy Synthes SC Collaboration Management

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
Cerenovus Enterprise 2 intracranial vessel stenting device on the delivery system in variants:

1. Cerenovus Enterprise 2 intracranial vessel stenting device on the delivery system with radiopaque distal guidewire tip 12 mm, working diameter of stent 4 mm, minimal working (effective) length of stent 16 mm in composition:
 - 1.1. Self-expanding stent – 1 pc.
 - 1.2. Delivery system in composition: guidewire with radiopaque distal tip – 1 pc., introducer 1 pc.
 - 1.3. Instructions for use – 1 pc.
 - 1.4. Implant card – 1 pc.
2. Cerenovus Enterprise 2 intracranial vessel stenting device on the delivery system with radiopaque distal guidewire tip 12 mm, working diameter of stent 4 mm, minimal working (effective) length of stent 23 mm in composition:
 - 2.1. Self-expanding stent – 1 pc.
 - 2.2. Delivery system in composition: guidewire with radiopaque distal tip – 1 pc., introducer 1 pc.
 - 2.3. Instructions for use – 1 pc.
 - 2.4. Implant card – 1 pc.
3. Cerenovus Enterprise 2 intracranial vessel stenting device on the delivery system with radiopaque distal guidewire tip 12 mm, working diameter of stent 4 mm, minimal working (effective) length of stent 30 mm in composition:
 - 3.1. Self-expanding stent – 1 pc.
 - 3.2. Delivery system in composition: guidewire with radiopaque distal tip – 1 pc., introducer 1 pc.
 - 3.3. Instructions for use – 1 pc.

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3.4. Implant card – 1 pc.

4. Cerenovus Enterprise 2 intracranial vessel stenting device on the delivery system with radiopaque distal guidewire tip 12 mm, working diameter of stent 4 mm, minimal working (effective) length of stent 39 mm in composition:

- 4.1. Self-expanding stent – 1 pc.
- 4.2. Delivery system in composition: guidewire with radiopaque distal tip – 1 pc., introducer 1 pc.
- 4.3. Instructions for use – 1 pc.
- 4.4. Implant card – 1 pc.

5. Cerenovus Enterprise 2 intracranial vessel stenting device on the delivery system without radiopaque distal guidewire tip, working diameter of stent 4 mm, minimal working (effective) length of stent 16 mm in composition:

- 5.1. Self-expanding stent – 1 pc.
- 5.2. Delivery system in composition: guidewire with radiopaque distal tip – 1 pc., introducer 1 pc.
- 5.3. Instructions for use – 1 pc.
- 5.4. Implant card – 1 pc.

6. Cerenovus Enterprise 2 intracranial vessel stenting device on the delivery system without radiopaque distal guidewire tip, working diameter of stent 4 mm, minimal working (effective) length of stent 23 mm in composition:

- 6.1. Self-expanding stent – 1 pc.
- 6.2. Delivery system in composition: guidewire with radiopaque distal tip – 1 pc., introducer 1 pc.
- 6.3. Instructions for use – 1 pc.
- 6.4. Implant card – 1 pc.

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7. Cerenovus Enterprise 2 intracranial vessel stenting device on the delivery system without radiopaque distal guidewire tip, working diameter of stent 4 mm, minimal working (effective) length of stent 30 mm in composition:

7.1. Self-expanding stent – 1 pc.

7.2. Delivery system in composition: guidewire with radiopaque distal tip – 1 pc., introducer 1 pc.

7.3. Instructions for use – 1 pc.

7.4. Implant card – 1 pc.

8. Cerenovus Enterprise 2 intracranial vessel stenting device on the delivery system without radiopaque distal guidewire tip, working diameter of stent 4 mm, minimal working (effective) length of stent 39 mm in composition:

8.1. Self-expanding stent – 1 pc.

8.2. Delivery system in composition: guidewire with radiopaque distal tip – 1 pc., introducer 1 pc.

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Figure 1— Cerenovus Enterprise 2 Vascular Reconstruction Device

- A. Reference marker (150 cm from delivery wire tip)
- B. Proximal delivery wire
- C. Introducer
- D. Proximal marker of delivery wire
- E. Proximal stent markers
- F. Stent positioning marker of delivery wire
- G. Stent body
- H. Distal stent markers
- I. Distal marker of delivery wire (if applicable)
- J. Introducer Flush Holes

The CERENOVUS ENTERPRISE 2 Vascular Reconstruction Device is designed for use under fluoroscopy with the PROWLER SELECT Plus Microcatheter (Infusion Catheter), (a 0.021" inner diameter, 5 cm distal length microcatheter (infusion catheter) manufactured by Cerenovus).

CAUTION: Compatibility with other microcatheters (infusion catheters) has not been established.

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Specifications are represented in tables 2-6.

Table 2 Nominal characteristic of stents

Recommended Parent Vessel Diameter (mm)	Stent Length @ 4.0 mm Vessel Diameter (mm)	Length Foreshortening @ 4.0 mm Vessel Diameter		% Open Area		Unconstrained Stent Diameter (mm)
		%	mm	4.0 mm	2.0 mm	
2.0-4.0	16	7.3	1.2	96.3	92.9	5.0
2.0-4.0	23	7.7	2.0	96.1	92.7	5.0
2.0-4.0	30	8.5	2.8	96.1	92.6	5.0
2.0-4.0	39	9.9	4.3	96.0	92.4	5.0

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Table 3 – Stents Specifications (for all codes)

Parameter	Value
Unconstrained Stent OD	5.0 mm (+ 0.5 / - 0.3) mm
Unconstrained Stent Flare OD	≥ 6.8 mm
Markerband	4 proximal markerbands and 4 distal ¹ markerbands
Markerband OD	0.2 +/- 0.038 mm
Markerband Length	1.00 +/- 0.05 mm
Coating thickness	2.5 ± 0.5 μm

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Indications for use

The CERENOVUS ENTERPRISE 2 Vascular Reconstruction Device is intended for use with occlusive devices in the treatment of intracranial aneurysm.

Note – The final decision on the application of Cerenovus Enterprise 2 must be accepted by user on the basis of existing knowledge and risks associated with a particular case.

Contraindications

Intracranial artery stenting is generally contraindicated in the following patient types:

- Patients in whom antiplatelet and/or anticoagulation therapy is contraindicated.
- Patients in whom the angiography demonstrates the anatomy is not appropriate for endovascular treatment, due to conditions such as:
 - Severe intracranial vessel tortuosity or stenosis
 - Intracranial vasospasm not responsive to medical therapy.

Note – The final decision on the application of Cerenovus Enterprise 2 must be accepted by user on the basis of existing knowledge and risks associated with a particular case.

Warnings

- The stenting procedure should be carried out under the direction of personnel with the requisite interventional training, especially intracranial stent procedures. Appropriate facilities should be available for managing the potential complications of the procedure.
- The device is designed to be manipulated while under high-quality fluoroscopic observation. If resistance is met during manipulation, determine the cause of resistance before proceeding.

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- This product is equipped with a temperature indicator label, found on the inner pouch. The acceptance criterion for the product is delineated by the graphic below the indicator label with the green check mark. The reject criterion for the label is delineated by the graphics marked with a red "X". Do not use if the temperature exposure indicator mirrors the reject criterion graphics, as the constrained stent diameter may have been compromised by exposure to high temperature.

- Do not use if the inner package is opened or damaged.

- Persons with allergic reactions to nickel titanium (Nitinol) may suffer an allergic response to this implant.

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- Select a stent length that is at least 10 mm longer than the aneurysm neck to maintain a minimum of 5 mm on either side of the aneurysm neck. A large differential in diameter between the proximal and distal segments of the parent vessel may increase the risk of stent migration.

- Adverse events may occur without warning. At all times a fully equipped emergency cart and resuscitation equipment should be readily available, and personnel competent in recognizing and treating adverse events of any severity should be on hand.

- In parent vessels with a diameter of less than 2.5 mm, only the Microcatheter Jailing Technique should be used. At diameters of less than 2.5 mm, the stent cell opening will not be able to accommodate a pass-through access of a Cerenovus microcatheter.

- During the Microcatheter Jailing Technique, loss of microcatheter position of the jailed microcatheter may lead to inability to re-access the aneurysm or coil protrusion out of aneurysm.

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- Selection of improper delivery system configuration may cause stent instability, inaccurate placement during deployment, or vessel damage.

Precautions

- Experience with stent implants indicates that there is a risk of stenosis. Stenosis may require dilatation of the vessel segment containing the stent. The risks and long-term outcome following dilatation of endothelialized stents is unknown at present.

- The CERENOVUS ENTERPRISE 2 Vascular Reconstruction Device is not intended for use as a stand-alone device, i.e. without subsequent coil embolization of the aneurysm.

- Do not use the CERENOVUS ENTERPRISE 2 Vascular Reconstruction Device if any component appears damaged or missing.

- Confirm that the device labeling clearly indicates the size of the stent to be used.

- Do not expose the system to organic solvents (e.g., alcohol).

- For single use only. Do not resterilize or reuse.

- Use product prior to the "Use By" date.

- Store system in a cool, dark, dry place.

- Dispose of all used devices in accordance with hospital policy for biohazardous materials.

- Coil protrusion during embolization may not be visualized fluoroscopically because of the superimposition of the stent and coil mass. Intermittent angiograms in multiple views may be

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necessary to ensure there are no coil loops protruding into the parent artery.

- Do not recapture the stent more than once.
- During deployment, the vascular reconstruction device may foreshorten. Refer to the “Device Description” section to examine foreshortening values for each of the vascular reconstruction device sizes.
- The performance and safety of two or more overlapped stents has not been established. The ability of the stent to withstand post balloon dilatation has not been established.
- Use caution when crossing the deployed stent with guidewires or accessory devices.
- The use of the CERENOVUS ENTERPRISE 2 Vascular Reconstruction Device in patients in whom antiplatelet and/or anticoagulant therapy is contraindicated could result in a higher risk of thrombosis.

Potential adverse events

Adverse events that may be associated with the use of Cerenovus Enterprise 2 in the intracranial arteries include:

- Aneurysm recanalization
- Myocardial infarction
- Allergic reaction including, contrast, Nitinol metal, and medications
- Neurologic deficit
- Pain and/or infection (at insertion site)
- Arrhythmia
- Perforation
- Arteriovenous fistula

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- Clopidogrel: A loading dose of 300 mg Clopidogrel given at time of catheterization prior to the procedure.

Elective intra-procedure medications

It was recommended that a weight-adjusted bolus of intravenous heparin be given to maintain the activated clotting time near 250–300 seconds throughout the procedure. If a GP IIb-IIIa Inhibitor was given, the study protocol recommended that the ACT be maintained at less than 220 seconds to minimize the risk of bleeding.

Post-procedure medications

- ASA (75–325 mg/day) should be continued for at least one year.
- Clopidogrel (75 mg/day) for at least 30 days.

Table 7 – Adjudicated Serious Adverse Events Related to Procedure or Device

Device or Procedure Related ¹	All Subjects (N = 31)	95% Confidence Interval
Any Procedure or Device Related Event	16.1% (5/31)	[5.5%, 33.7%]
Stroke/CVA	6.5% (2/31)	[0.8%, 21.4%]
Stent thrombosis	3.2% (1/31)	[0.1%, 16.7%]
Coil prolapse or migration	3.2% (1/31)	[0.1%, 16.7%]
Neurological deterioration	3.2% (1/31)	[0.1%, 16.7%]
Extracranial/intracranial bypass	3.2% (1/31)	[0.1%, 16.7%]
Incomplete Embolization of aneurysm	3.2% (1/31)	[0.1%, 16.7%]
Infection of femoral vein	3.2% (1/31)	[0.1%, 16.7%]
Leg Hematoma	3.2% (1/31)	[0.1%, 16.7%]
Retroperitoneal hematoma	3.2% (1/31)	[0.1%, 16.7%]

¹ One subject experienced five Serious Adverse Events

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Clinical Study

The Cerenovus Enterprise 2 Vascular Reconstruction Device and Delivery System Study was conducted at three investigative centers in Europe. Subjects were eligible if they presented with a ruptured (Hunt and Hess Grade I – III) or non-ruptured intracranial wide neck saccular aneurysm that was deemed by the attending interventional neuroradiologist to be an acceptable candidate for endovascular coil embolization. Wide neck was defined as having a neck width >4 mm or a dome-to-neck ratio <2.

Thirty-one subjects were enrolled into the study. The mean age of the 31 subjects enrolled in this study was 54.5 years and a majority of the subjects were female (87.1%). A medical history of hypertension was reported in 38.7% of the subjects enrolled and 25.8% had an endocrine/metabolic disease. The most commonly reported neurological event was subarachnoid hemorrhage (38.7%) with 16.1% reporting a history of cranial nerve palsy and 16.1% reporting of pounding/pulsatile headaches. At baseline, 45.2% of the subjects had undergone an intracranial embolization procedure.

Of the 31 subjects who were entered into the clinical study and treated with a Cerenovus Enterprise 2 Vascular Reconstruction Device, the aneurysms were most commonly located in the internal

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carotid/ophthalmic location (32.3%). Table 8 summarizes the locations of the treated aneurysms.

Table 8 - Aneurysm Location

Aneurysm	Total (All Subjects N = 31)	
	N	%
Carotid-Ophthalmic	10/31	32.3
Posterior Communicating Artery	6/31	19.4
Basilar Artery	5/31	16.1
Middle Cerebral Artery	4/31	12.9
Anterior Communicating Artery	3/31	9.7
Carotid, Other	2/31	6.5
Vertebral	1/31	3.2

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Aneurysm and Parent Artery Dimensions

An independent angiographic core laboratory was used to assess the aneurysm and parent artery dimensions. Table 9 contains the pre-procedure, post-procedure and six months aneurysm dimensions and Table 10 contains the parent artery dimensions. At pre-procedure, the mean aneurysm dome height and width were 6.3 mm and 7.0 mm, respectively. The mean neckwidth was 4.8mm, and mean dome width-to-neck ratio was 1.47. The mean proximal/distal diameter of the parent vessel pre-procedure was 3.1/2.8, and 2.8/2.5 at six months.

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Numbers are % (numerator/denominator) or mean $\pm$ 1 SD

1 Sample sizes vary due to observations classified as N/D = Not Determinable.

2 One subject was lost to follow-up and did not complete the six-month assessment.

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Procedure Success Measures

Of the 31 subjects entered into the study who received the VRD, 30 were followed for six months (one subject was lost to follow-up and did not complete the six-month assessment). Technical success was judged based on angiographic core lab assessment. Successful VRD placement with satisfactory coil mass position immediately post-procedure, and without occurrence of any procedural serious adverse events reported by the site was 90.0%. Successful VRD placement was defined as stable VRD placement with complete coverage across the aneurysm neck and parent artery patency, while satisfactory coil mass position was defined as VRD maintains coil position within the sac with parent artery patency. Maintenance of coil mass position was 96.7% at six months. The mean post-procedure and six month percent aneurysm occlusion were 93.0% and 93.9% at six months. The protocol-defined procedural success measures are provided in Table 11. Two subjects treated by the same investigator are known to have had the VRD placed, and then returned at 30 days for coiling of the aneurysm when the VRD was fixed due to endothelialization. This method of treatment did not affect the study's procedural success measures.

Table 11 - Procedure Success Measures

Parameters Measured	Post-Procedure (N = 31) ¹	95% Confidence Interval	Six months (N = 31) ^{1, 2}	95% Confidence Interval
Procedure Success Measures				

Proximal Parent Artery Coverage (mm)				
Mean $\pm$ SD (n)	9.7 $\pm$ 4.4 (31)	[8.1, 11.3]	- $\pm$ - (0)	[-, -]
Range (min., med., max.)	(3.5, 9.4, 18.6)		(-, -, -)	
Distal Parent Artery Coverage (mm)				
Mean $\pm$ SD (n)	6.6 $\pm$ 3.0 (31)	[5.5, 7.7]	- $\pm$ - (0)	[-, -]
Range (min., med., max.)	(0.0, 6.7, 12.6)		(-, -, -)	
Stable VRD Placement with Complete Neck Coverage	96.8% (30/31)	[83.3%, 99.9%]	- (- / -)	[-, -]
Maintenance of Coil Mass Position	96.7% (29/30)	[82.8%, 99.9%]	96.7% (29/30)	[82.8%, 99.9%]
Percent Aneurysm Occlusion				
Mean $\pm$ SD (n)	93.0% $\pm$ 5.6 (31)	[90.9%, 95.0%]	93.9% $\pm$ 10.8 (30)	[89.9%, 98.0%]
Range (min., med., max.)	(77.3%, 95.0%, 100.0%)		(50.0%, 99.0%, 100.0%)	
Procedure Success 3	90.0% (27/30)	[73.5%, 97.9%]	- (- / -)	[-, -]

Numbers are % (numerator/denominator) or mean $\pm$ 1 SD.

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2One subject was lost to follow-up and did not complete the six-month assessment.

3 Successful VRD placement with satisfactory coil mass position angiographically assessed immediately post procedure, and without the occurrence of procedural serious adverse events. The duration of the procedure was defined as the time the stent delivery system was introduced into the infusion catheter until the time the guiding catheter was removed from the subject

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MRI Safety Information

MR Conditional

The device is conditionally MR compatible. Non-clinical testing demonstrated that the Enterprise 2 Stent is MR Conditional. A patient with this device can be scanned safely, immediately after placement under the following conditions:

- Static magnetic field of 1.5-Tesla and 3-Tesla, only
- Maximum spatial gradient magnetic field of 4,000 Gauss/cm (40T/m)
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 4-W/kg with the MR system operating in the First Level Controlled Operating Mode

MRI-Related Heating

Under the scan conditions defined above, the Enterprise 2 Stent is expected to produce a maximum temperature rise of 2.5°C after 15 minutes of continuous scanning (per pulse sequence). The effect of MRI related heating for overlapping stents or stents with fractured struts is unknown.

Artifact Information

The maximum artifact size extends approximately 5-mm from the Enterprise 2 stent when imaged using a gradient echo pulse sequence on a 3.0 T MRI System. The lumen of this stent cannot be visualized.

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Preparations for Use

In addition to the CERENOVUS ENTERPRISE 2 Vascular Reconstruction Device, the following items are recommended:

- Appropriately sized Sheath Introducer of any Manufacturer.
- A **PROWLER SELECT Plus** Microcatheter (Infusion Catheter) (0.021" inner diameter, 5 cm distal length microcatheter (infusion catheter) manufactured under Cerenovus brand) (RC FSZ 2010/07085 Devices for Angiography and embolization).
- **TRUFILL** family of Detachable Coils (RC FSZ 2010/07084 Angioplasty and stenting devices with accessories) and **TRUFILL** DCS Syringe (RC FSZ 2010/07085 Devices for Angiography and embolization).
- Appropriately sized exchange length guidewire of any Manufacturer for the selected infusion catheter.
- **ENVOY** Guiding Catheter (RC FSZ 2010/07085 Devices for Angiography and embolization and RC RZN 2013/676 Intravascular catheters, variations: PROWLER, ENVOY, REVIVE) for the selected microcatheter (infusion catheter).
- Two or more Y connector/Rotating Hemostasis Valves (RHV) with >0.076" lumen diameter of any Manufacturer.
- Sterile heparinized saline solution of any Manufacturer.

Vascular Reconstruction Device Selection

Appropriate selection of the **CERENOVUS ENTERPRISE 2 Vascular** Reconstruction Device is important for patient safety. In order to choose the optimal **CERENOVUS ENTERPRISE 2 Vascular** Reconstruction Device for any given lesion, examine pre-stent procedure angiograms.

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Directions for Use

1. Gain vascular access according to standard practice. **Note:** If using the Microcatheter Jailing Technique, select a Guiding Catheter of sufficient lumen (ID) to accept the Microcatheter (Infusion Catheter) required for delivery of the **CERENOVUS ENTERPRISE 2 Vascular Reconstruction Device (0.021" PROWLER SELECT Plus** Microcatheter (Infusion Catheter), 5 cm distal length) as well as the Microcatheter (Infusion Catheter) required to deliver embolic coils, while allowing adequate contrast infusion around these two Microcatheters (Infusion Catheters) for fluoroscopy mapping.

2. Navigate the microcatheter (infusion catheter) for stent deployment (0.021" **PROWLER SELECT Plus** Microcatheter (Infusion Catheter), 5 cm distal length) over a guidewire at least 1.2 cm distal to the aneurysm neck. **Note:** If using the Microcatheter Jailing Technique, navigate the microcatheter (infusion catheter) delivering the embolic coils, using standard endovascular practice, until the aneurysm site has been accessed.

3. Remove the guidewire from the microcatheter (infusion catheter) being used for stent deployment.

4. Maintain flush through all microcatheters (infusion catheters) per standard endovascular practice.

5. Select the appropriate **CERENOVUS ENTERPRISE 2 Vascular Reconstruction Device**.

6. Carefully inspect the Vascular Reconstruction Device package for damage to the sterile barrier.

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7. Peel open the pouch using aseptic technique.

8. Carefully place the dispenser hoop into the sterile field.

9. Remove the delivery wire from the clip on the dispenser hoop. Grasp the proximal end of the introducer and the delivery wire at the point where it exits the introducer. Hold the delivery wire and introducer together to prevent stent movement. Remove the **CERENOVUS ENTERPRISE 2 Vascular** Reconstruction Device from the dispenser hoop.

- Do not partially deploy the stent from the introducer.

- Confirm that the delivery wire does not move relative to the introducer during removal of the **CERENOVUS ENTERPRISE 2 Vascular** Reconstruction Device from the dispenser hoop.

- Confirm that the tip of the delivery wire is entirely within the introducer.

- Confirm that the delivery wire is not kinked and that the introducer tip is not damaged. **DO NOT CONTINUE** if either defect is observed; return the device to Manufacturer.

- **WARNING:** Do not shape the tip of the delivery wire if a 12 mm tip configuration is selected. This could lead to damage or separation of delivery wire.

10. Advance the distal end of the introducer into the RHV connected to the microcatheter (infusion catheter) until it is fully engaged with the microcatheter (infusion catheter) hub (0.021" **PROWLER SELECT Plus** Microcatheter (Infusion Catheter)), then tighten the RHV locking ring. Flush the y-connector of the RHV with sterile saline and verify that fluid exits the proximal end of the introducer.

- **WARNING:** Purge the device carefully to avoid the accidental introduction of air into the system.

- **WARNING:** Confirm that there are no air bubbles trapped anywhere in the system.

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CAUTION: The introducer must be properly engaged with the microcatheter (infusion catheter) hub to enable stent introduction into the microcatheter (infusion catheter).

11. Advance the delivery wire to transfer the stent from the introducer into the microcatheter (infusion catheter).

WARNING: Delivery wire should not be torqued to gain access into the aneurysm.

12. Continue advancing the delivery wire into the microcatheter (infusion catheter) until the distal edge of the delivery wire reference marker (150 cm from the delivery wire distal tip) enters the introducer. Loosen the RHV locking ring, remove the introducer, and set it aside. **Note:** Fluoroscopy may be used up to this point at the physician's discretion.

WARNING: Do not apply undue force if resistance is encountered at any point during stent manipulation. Withdraw the unit and advance a new one.

13. Track the stent through the microcatheter (infusion catheter) to the tip.

14. Position the stent for deployment by aligning the stent positioning marker of the delivery wire with the target site (Figure 2).

15. If stent positioning is satisfactory, carefully retract the microcatheter (infusion catheter), while maintaining the position of the delivery wire, to allow the stent to deploy across the neck of the aneurysm. The stent will expand as it exits the microcatheter (infusion catheter).

WARNING: Do not deploy the stent if it is not properly positioned in the vessel. **Note:** If using the Microcatheter Jailing Technique, verify that the position of the microcatheter (infusion catheter) delivering the embolic coils has been maintained. If necessary, reposition the microcatheter (infusion catheter) using standard endovascular practice.

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WARNING: Do not deploy the stent if the position of the microcatheter (infusion catheter) delivering the embolic coils has been lost.

16. If stent positioning is not satisfactory, the stent may be recaptured and repositioned. The stent may be recaptured until the point where the proximal end of the stent positioning marker is aligned with the microcatheter (infusion catheter) distal markerband (recapturability limit) (Figure 2). If stent repositioning is required, gently advance the microcatheter (infusion catheter) over the deployed stent (do not pull the stent back into the microcatheter (infusion catheter)), reposition the system, and re-deploy the stent in the new location. **Note:** When advancing the microcatheter (infusion catheter) over the stent during recapture, it may be necessary to keep the stent stable with tension on the delivery wire.

CAUTION: If resistance is felt while recapturing the stent, do not continue to recapture the device. Withdraw the microcatheter (infusion catheter) slightly to unsheath the stent (without exceeding the recapture limit), and then attempt to recapture the stent again.

CAUTION: The stent may be fully recaptured once. **Note:** Be careful to maintain delivery wire access through the detached stent to facilitate access distal to the deployed stent, if necessary. **Note:** Maintain adequate stent length (approximately 5 mm) on each side of the aneurysm neck to ensure appropriate neck coverage.

17. Prior to removing the delivery wire, position the microcatheter distal to the stent to maintain access through the stent. Slowly remove and discard the **CERENOVUS ENTERPRISE** delivery wire.

18. Advance an exchange-length guidewire through the 0.021" **PROWLER SELECT Plus** Microcatheter (Infusion Catheter).

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19. When using the Microcatheter Jailing Technique remove the 0.021" **PROWLER SELECT Plus** Microcatheter (Infusion Catheter) and skip to step 25.

20. Remove the 0.021" **PROWLER SELECT Plus** Microcatheter (Infusion Catheter) and advance a 0.014" (according to the **TRUFILL** family of Detachable Coils Instructions for Use) inner diameter microcatheter (infusion catheter).

21. Use the guidewire and microcatheter (infusion catheter) to access the aneurysm through the stent cells (Figure 3). **Note:** Access to the aneurysm may be facilitated by the use of an microcatheter (infusion catheter) that has been shaped.

22. After the microcatheter (infusion catheter) is positioned within the aneurysm, detachable coils may be delivered into the aneurysm according to conventional methods.

CAUTION: Compatibility of the **CERENOVUS ENTERPRISE 2 Vascular** Reconstruction Device has only been established with the **TRUFILL** family of Detachable Coils.

WARNING: Observe stent marker position during the coiling procedure to ensure that the stent does not migrate from its deployed position.

23. After placement of the last coil, verify that the stent has remained patent and properly positioned. Carefully remove the microcatheter through the stent cells.

24. After completing the procedure, withdraw and discard all applicable accessory devices.

The following steps apply to the Microcatheter Jailing Technique only:

25. Deliver detachable coils using the microcatheter (infusion catheter) previously placed within the aneurysm according to conventional methods.

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WARNING: Observe markers on the microcatheter (infusion catheter) used for coil delivery to ensure catheter position has been maintained.

26. After placing the last coil, verify that the stent has remained patent and properly positioned. Carefully remove the catheter used for coil delivery.

WARNING: Remove microcatheter (infusion catheter) used for coil delivery carefully to prevent movement of stent.

27. Complete procedure by withdrawing and discarding all applicable accessory devices.

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
Figure 2 – Vascular Reconstruction Device deployed to recapturability limit (Proximal end of the stent positioning marker is even with the distal microcatheter (infusion catheter) marker band).

LEFT: vessel and microcatheter (infusion catheter) shown as transparent to indicate Vascular Reconstruction Device and Delivery Wire landmarks. MIDDLE: schematic of fluoroscopic view with 12 mm Tip Delivery Wire. RIGHT: schematic of fluoroscopic view with No Distal Tip Delivery Wire.

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




Warranty



Manufacturer warrants that this medical device is free from defects in both materials and workmanship. Any other express or implied warranties, including warranties of merchantability or fitness, are hereby disclaimed. Suitability for use of this medical device for any particular surgical procedure should be determined by the user in conformance with the manufacturer's instructions for use. There are no warranties that extend beyond the description on the face hereof.

Warranty shelf life is 3 years from the manufacturing date.

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Delivery Wire	Corewire	Nitinol
	Distal Marker (coil)	Platinum/Tungsten
	Stent Positioning Marker (coil)	Platinum/Tungsten
	Proximal Marker (coil)	Platinum/Tungsten
	Solder/Flux (for Delivery Wire Marker attachment)	Solder: EutecRod 157 Flux: Indalloy flux #2

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Package contents

Cerenovus Enterprise 2 Vascular Reconstruction Device in variants:

1. Cerenovus Enterprise 2 Vascular Reconstruction Device consisting of:

1) Cerenovus Enterprise 2 Vascular Reconstruction Device – 1 pc:

- Cerenovus Enterprise 2 Vascular Reconstruction Device 4.0mmx16mm.
- Cerenovus Enterprise 2 Vascular Reconstruction Device 4.0mmx23mm.
- Cerenovus Enterprise 2 Vascular Reconstruction Device 4.0mmx30mm.

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- Cerenovus Enterprise 2 Vascular Reconstruction Device 4.0mmx39mm.

2) Instructions for use – 1 pc.

3) Stent implant card – 1 pc.

2. Cerenovus Enterprise 2 Vascular Reconstruction Device NO TIP
consisting of:

1) Cerenovus Enterprise 2 Vascular Reconstruction Device NO TIP – 1 pcs:

- Cerenovus Enterprise 2 Vascular Reconstruction Device 4.0mmx16mm NO TIP.
- Cerenovus Enterprise 2 Vascular Reconstruction Device 4.0 mm x 23 mm NO TIP.
- Cerenovus Enterprise 2 Vascular Reconstruction Device 4.0 mm x 30 mm NO TIP.
- Cerenovus Enterprise 2 Vascular Reconstruction Device 4.0 mm x 39 mm NO TIP.

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Medos International S.p.A. - Chemnitz - Blanc 38, case postale, CH-2400 Le Locle - T : +41 32 934 8000