State of Florida

Department of State

APOSTILLE

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1. Country: **United States of America**

This public document

2. has been signed by

Katherine Perez

3. acting in the capacity of

Notary Public of Florida

4. bears the seal/stamp of

Notary Public, State of Florida

Certified

5. at **Tallahassee, Florida**

6. the **Fourth day of September, A.D., 2013**

7. by **Secretary of State, State of Florida**

8. No. **2013-112161**

9. Seal/Stamp:

10. Signature:

Ken Retzner

Secretary of State

DSDE 99 (2/12)

This document contains a true watermark. Hold up to light to see "SAFE" and "VERIFY FIRST."


"State of Florida" appears in small letters across the face of this 8 1/2 x 11" document

ATTESTATION

STATE OF FLORIDA

COUNTY OF DUVAL

On this 27th day of August 2013, I attest that the preceding or attached documents are true, exact, complete, and unaltered photocopies made by me of the following Instruction for Use No. M726750B765 A for ALAR™ Nasal Valve Stent by Medtronic Xomed, Inc., presented to me by the document's custodian, Julia Ondruska, and, to the best of my knowledge, that the photocopied document is neither a public record nor a publicly recordable document, certified copies of which are available from an official source other than a notary public.

Official Notary Signature


United States of America
State of Florida
County of Duval

Subscribed and sworn to before me
this 27th day of August, 2013

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English

Indications for Use

The ALAR Nasal Valve Stent supports and immobilizes intranasal tissues and cartilage during the post-operative healing period following intranasal and nasal valve surgery.

Device Description

The elliptically shaped fluoroplastic stent is approximately 20 mm long and 16 mm wide. Slits symmetrically spaced along the stent's lateral edges allow for left and right nostril usage. The suture and needle attached to the stent provide stent anchoring and hold cartilage plates stable during post-operative healing.

Contraindications

Contraindications for the ALAR Nasal Valve Stent include use in patients with insufficient or over-resected alar cartilages from prior surgery or non-operated, weak or malpositioned lateral crura.

Warnings and Precautions

- Do not use cauteries and other surgical equipment capable of generating high temperatures close to the stents. Contact with such devices may cause damage to the device and result in burns to the patient and/or medical personnel.
- Do not modify the devices or malfunctions may occur and the device can become displaced and inhaled or ingested or efficacy can be compromised.
- Use caution when opening the pouch which is supplied with an anchoring suture and suture needle attached, so as not to puncture gloved hands and compromise the sterile field.
- Use caution when securing with the suture to avoid pressure necrosis of tissue.

Possible Adverse Effects

- Possible adverse effects associated with the use of the ALAR Nasal Valve Stent include possible infection, recurrent or persistent bleeding, dislodgement of the packing, trauma to intranasal tissues, and/or tissue necrosis.
- Toxic shock syndrome is possible with certain fluid absorbing nasal packing materials. Monitor patients closely for indications of this condition.
- Altered respiratory function and reduced arterial oxygen (PaO₂) levels may result when using intranasal packing. Monitor patients closely for signs of respiratory distress and decreased oxygenation, particularly those at risk for impaired respiratory status.

Instructions for Use

It is the responsibility of the physician to decide which techniques and procedures are appropriate for the application of the ALAR Nasal Valve Stent based on medical training and the specific needs of the patient.

1. Following dissection and/or preparation of tissues as required per patient

English

- needs, secure stent to tissues with supplied needle/suture.
2. Using needle driver, pass needle through the tissues then
3. Pass needle through the center hole in the stent and clench suture in slit at desired tension to achieve desired tissue result.
4. Complete additional anchoring by passing suture through lateral cleats.
5. Tie suture to free end of suture to secure anchoring.
6. Cut suture tails and discard needle appropriately.

Limited Warranty

A. This LIMITED WARRANTY provides assurance for the customer who purchases a Medtronic Xomed Product (hereinafter the "Product") that should the Product fail to function to Medtronic Xomed's published specifications during the term of this **LIMITED WARRANTY** (one year from the date of shipment for new Product, 90 days from date of shipment for refurbished or used Product), Medtronic Xomed will either replace, repair, or issue a credit (adjusted to reflect the age of the Product) for the Product or any portion thereof. This **LIMITED WARRANTY** is extended only to the buyer purchasing the Product directly from Medtronic Xomed or from its affiliate or its authorized distributor or representative.

B. To qualify for this LIMITED WARRANTY, the following conditions must be met:

- (1) The Product must be used on or before its "Use By" or "Use Before" date, if applicable.
- (2) The Product must be used in accordance with its labeling and may not be altered or subjected to misuse, abuse, accident or improper handling.
- (3) Medtronic Xomed must be notified in writing within thirty (30) days following discovery of a defect.
- (4) The Product must be returned to Medtronic Xomed within thirty (30) days of Medtronic Xomed receiving notice as provided for in (3) above.
- (5) Upon examination of the Product by Medtronic Xomed, Medtronic Xomed shall have determined that: (i) the Product was not repaired or altered by anyone other than Medtronic Xomed or its authorized representative, (ii) the Product was not operated under conditions other than normal use, and (iii) the prescribed periodic maintenance and services have been performed on the Product.

C. This LIMITED WARRANTY is limited to its express terms. **THIS LIMITED WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED WHETHER STATUTORY OR OTHERWISE, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.** In no event shall Medtronic Xomed be liable for any consequential, incidental, prospective or other similar damage resulting from a defect, failure, or malfunction of the Product, whether a claim for such damage is based upon the warranty, contract, negligence or otherwise.

D. The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. Users may benefit from statutory warranty rights under legislation governing the sale of consumer goods. If any part or term of this **LIMITED WARRANTY** is held by any court of competent jurisdiction to be illegal, unenforceable, or in conflict with applicable law, the validity of the remaining portion of the **LIMITED WARRANTY** shall not be affected, and all rights and obligations shall be construed and enforced as if this **LIMITED WARRANTY** did not contain the particular part or term held to be invalid.

International Customer Service

International customers should contact their local Medtronic Xomed office.

Note: When contacting Customer Service and Technical Support, please have the appropriate product number, date of purchase, and nature of inquiry available.

The following is a trademark or registered trademark of Medtronic, Inc. in the United States and other countries: ALAR™. All other trademarks, service marks, registered trademarks or registered service marks are the property of their respective owners in the United States and other countries.

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For Russia

Appendix to Instruction for Use ALAR Nasal Valve Stent

1. Packaging	The elliptically shaped fluoroplastic stent is approximately 20 mm long and 16 mm wide. The stent packaged one stent per pouch with two pouches per box and with the addition of the Tyvek envelope.
2. Shelf life	4 years (1460 days)
3. Utilization	Utilize in accordance with the requirements of the national legislation on environment and health.
4. Methods and conditions of disinfection and sterilization	The devices are single use devices, sterilization method is 100% EtO. They should not be resterilized or disinfected after use.
5. Operating conditions	The product can be used under the standard conditions. Operation room temperature range from 15 °C to 30 °C. During operation at nominal temperature values from 32°C to 42°C, humidity range from 80% to 100% (+2%). Other requirements set by customer request.
6. Transportation conditions	Packed by the manufacturer sheaths are transported by all models of transport in closed vehicles in accordance with the rules of carriage of goods, operating in each mode of transport.
7. Storage conditions	Storage in its original packaging including shelf carton in a clean, cool and dry area to protect the product and minimize the potential for contamination. Temperature: no more 27°C (80° F) Rotate inventory so that devices are used prior to the expiration date on package label.
8. Warranty	Limited Warranty A. This LIMITED WARRANTY provides assurance for the customer who purchases a Medtronic Xomed Product (here in after the “Product”) that should the Product fail to function to Medtronic Xomed’s published specifications during the term of this LIMITED WARRANTY (one year from the date of shipment for new Product, 90 days from date of shipment for refurbished or used Product), Medtronic Xomed will either replace, repair, or issue a credit (adjusted to reflect the age of the Product) for the Product or any portion thereof. This LIMITED WARRANTY is extended only to the buyer purchasing the Product directly from Medtronic Xomed or from its affiliate or its authorized distributor or representative. B. To qualify for this LIMITED WARRANTY, the following conditions must be met: (1) The Product must be used on or before its “Use By” or “Use Before” date, if applicable.

	(2) The Product must be used in accordance with its labeling and may not be altered or subjected to misuse, abuse, accident or improper handling. (3) Medtronic Xomed must be notified in writing within thirty (30) days following discovery of a defect. (4) The Product must be returned to Medtronic Xomed within thirty (30) days of Medtronic Xomed receiving notice as provided for in (3) above. (5) Upon examination of the Product by Medtronic Xomed, Medtronic Xomed shall have determined that: (i) the Product was not repaired or altered by anyone other than Medtronic Xomed or its authorized representative, (ii) the Product was not operated under conditions other than normal use, and (iii) the prescribed periodic maintenance and services have been performed on the Product. C. This LIMITED WARRANTY is limited to its express terms. THIS LIMITED WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED WHETHER STATUTORY OR OTHERWISE, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. In no event shall Medtronic Xomed be liable for any consequential, incidental, prospective or other similar damage resulting from a defect, failure, or malfunction of the Product, whether a claim for such damage is based upon the warranty, contract, negligence or otherwise. Français ALAR Nasal Valve Stent 7 D. The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. Users may benefit from statutory warranty rights under legislation governing the sale of consumer goods. If any part or term of this LIMITED WARRANTY is held by any court of competent jurisdiction to be illegal, unenforceable, or in conflict with applicable law, the validity of the remaining portion of the LIMITED WARRANTY shall not be affected, and all rights and obligations shall be construed and enforced as if this LIMITED WARRANTY did not contain the particular part or term held to be invalid.
9. Customer Service	Customer Service: If you have any device returns or questions about Medtronic Xomed, Inc., USA device please call US Medtronic Customer Service at (+1-800) 468-9710 or Medtronic Russia tel.: (+7495) 580-7377.
	International customers should contact their local Medtronic Xomed office.

For Russia

	Note: When contacting Customer Service and Technical Support, please have the appropriate product number, date of purchase, and nature of inquiry available. The following is a trademark or registered trademark of Medtronic, Inc. in the United States and other countries: ALAR. All other trademarks, service marks, registered trademarks or registered service marks are the property of their respective owners in the United States and other countries.
--	---

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©Medtronic, Inc. 2013