ANITA Dr. Helbig GmbH
Grafenstrasse 23
83098, Brannenburg

Johannes Weber-Unger
Member of the executive board

19.12.2022

User manual

Version 1.1. Change date 19.12.2022

Exoprosthesis of the mammary gland,
models 1019X, 1050X, 1051X,
1054X, 1066X, 1080L ,
1080 R , 1081 L ,
1081R , 1086X

Manufacturer: Anita Dr. _ Helbig GmbH
Grafenstrasse 23, 83098, Brannenburg, Germany

User manual
Version 1.1. Change date 19.12.2022

CURRENT VERSION.

**Exoprosthesis of the mammary gland,
models 1019X, 1050X, 1051X,
1054X, 1066X, 1080L ,
1080R, 1081L,
1081R, 1086 X**
(hereinafter referred to as the exoprosthesis(s), product(s))

**Please read this instruction manual carefully before using this product and
keep it for future reference!**

SAFETY INSTRUCTIONS

READ THE OPERATING MANUAL BEFORE USING THE PRODUCT

WARNINGS AND PRECAUTIONS FOR USE:

- Do not use the product if it is damaged.
- For the correct selection of the product, you need to contact a specialist.

CONTENT

Introduction	5
1 Name, purpose and contraindications	6
2 General description of the exoprosthesis	7
3 Design features and appearance of the exoprosthesis	8
4 Specifications	14
5 Materials used in manufacturing	18
6 Use of an exoprosthesis	19
7 Rules for care and cleaning	20
8 Marking	21
9 Packaging	22
10 Safety requirements and environmental protection	23
1 1 Transport and storage	24
1 2 Maintenance and repair	25
1 3 Disposal information	26
1 4 Manufacturer's warranties	27
1 5 Claim information	29
List of regulatory documents used by the manufacturer to confirm compliance	30

INTRODUCTION

The main purpose of the instructions for use is to provide the user with instructions for using the exoprosthesis, as well as information regarding its cleaning. The manual provides exoprosthesis performance, cleaning, and safety information. Before use, the user must carefully read the instructions for use in order to use the exoprosthesis correctly and without error, while achieving the specified safety and performance standards.

1 PURPOSE AND CONTRAINDICATIONS

Breast exoprosthesis models 1019X, 1050X, 1051X, 1054X, 1066X, 1080 L , 1080 R , 1081 L , 1081 R , 1086X are designed to simulate the breast and maintain body symmetry in everyday life using a bra .

Scope : in hospitals and clinics, as well as at home (domestic) conditions .

Intended User: Patients after mastectomy, removal of breast implants, if necessary, breast augmentation.

Indications for use:

- after mastectomy;
- after removal of breast implants;
- if necessary, to replenish the volume of the breast;
- wearing until the healing of scars after surgery, but no more than 6 months. from the beginning of application (model 1019X).

Contraindications for use:

- individual intolerance to materials used in the manufacture of exoprotheses.

Possible complications:

- itching and irritation at the site of an allergic reaction to the materials used in the manufacture of exoprotheses.

Possible side effects:

- allergic reaction to materials used in the manufacture of exoprotheses.

Risk assessment: The product is designed and manufactured in such a way that, when used under conditions and for purposes consistent with the intended use, determined by the user as the intended user, in strict accordance with the indicated indications, contraindications, side effects, acts as intended and is safe for the user and third parties. persons (Risk management file F MDR - A - II -5 dated 25.05.2021).

Any risk associated with their use is acceptable when weighed against the benefit to the user (Risk Management File F MDR - A - II -5 from 25.05.2021)

2 GENERAL DESCRIPTION OF THE EXOPROSTHESIS

The breast exoprosthesis is an external prosthesis for use with a bra.

Depending on the model, the design of the exoprosthesis may include volume control, air circulation, and a fold on the back surface to provide similar behavior to a natural breast during use.

Depending on the design features, operational and other characteristics, exoprostheses are produced in the following models:

- 1019X ;
- 1050X ;
- 1051X ;
- 1054X ;
- 1066X;
- 1080 L ;
- 1080R ;
- 1081L ;
- 1081R ;
- 1086 H.

3 DESIGN FEATURES AND APPEARANCE OF THE EXOPROSTHESIS

DESIGN FEATURES AND APPEARANCE OF THE EXOPROSTHESIS

3.1 Breast exoprosthesis model 1019X


Figure 1. Appearance of the model 1019X breast exoprosthesis

Used on both sides.

Microfiber back.

The pre-formed soft shell ensures the optimal shape of the breast.

3.2 Breast exoprosthesis model 1050X


Figure 2. Appearance of the model 1050X breast exoprosthesis

Used on both sides.

The shape of the exoprosthesis imitates the shape of the mammary gland.

The rear surface of the exoprosthesis is covered with a layer of soft silicone film, which allows the exoprosthesis to fit perfectly.

A layer of soft silicone film on the front surface of the product allows you to make the product up to 35% lighter than without it.

The shape of the thinned edge of the exoprosthesis allows you to hide the transitions from the product to the body.

The walking behavior of the exoprosthesis is similar to the natural mammary gland, due to the presence of a fold on the back surface.

3.3 Breast Exoprosthesis Model 1051X


Figure 3. Appearance of the model 1051X breast exoprosthesis

Used on both sides.

The shape of the exoprosthesis imitates the shape of the mammary gland.

The walking behavior of the exoprosthesis is similar to the natural mammary gland, due to the presence of a fold on the back surface.

3.4 Breast Exoprosthesis Model 1054X


Figure 4. Appearance of the model 1054X breast exoprosthesis

Used on both sides.

The shape of the exoprosthesis imitates the shape of the mammary gland.

The design of the rear surface of the exoprosthesis is ribbed to ensure air circulation.

When swimming, water drains faster than other prostheses.

Exoprosthesis Model 1066X


Figure 5. Appearance of the model 1066X breast exoprosthesis

Used on both sides.


The shape of the exoprosthesis imitates the shape of the mammary gland.

The rear surface of the exoprosthesis is covered with a layer of soft silicone film, which allows the exoprosthesis to fit perfectly.

A layer of soft silicone film on the front surface of the product allows you to make the product up to 35% lighter than without it.

The design of the rear surface of the exoprosthesis is made of testurized, which allows for optimal air circulation and maintaining the temperature of the product, in connection with which the exoprosthesis is as comfortable as possible to use and creates a massage effect.

3.6 Breast exoprosthesis model 1080 L

L breast exoprosthesis

Used on the left side.

The shape of the exoprosthesis imitates the shape of the mammary gland on the left.


The back surface of the exoprosthesis is covered with a layer of soft silicone film, which allows the exoprosthesis to fit perfectly.

A layer of soft silicone film on the front surface of the product allows you to make the product up to 35% lighter than without it.

The shape of the thinned edge of the exoprosthesis allows you to hide the transitions from the product to the body.

The walking behavior of the exoprosthesis is similar to the natural mammary gland, due to the presence of a fold on the back surface.

3.7 Breast Exoprosthesis Model 1080 R

R breast exoprosthesis

Used on the right side.

The shape of the exoprosthesis imitates the shape of the breast on the right. The back surface of the exoprosthesis is covered with a layer of soft silicone film, which allows the exoprosthesis to fit perfectly.

A layer of soft silicone film on the front surface of the product allows you to make the product up to 35% lighter than without it.

The shape of the thinned edge of the exoprosthesis allows you to hide the transitions from the product to the body.

The walking behavior of the exoprosthesis is similar to the natural mammary gland, due to the presence of a fold on the back surface.

3.8 Breast Exoprosthesis Model 1081 L


Figure 8. Appearance of the exoprosthesis of the mammary gland of the model 1081 L
Used on the left side.

The shape of the exoprosthesis imitates the shape of the mammary gland on the left.

The walking behavior of the exoprosthesis is similar to the natural mammary gland, due to the presence of a fold on the back surface.

3.9 Breast Exoprosthesis Model 1081 R


Figure 9. Appearance of the breast exoprosthesis model e 1081 R

Used on the right side.

The shape of the exoprosthesis imitates the shape of the breast on the right.

The walking behavior of the exoprosthesis is similar to the natural mammary gland, due to the presence of a fold on the back surface.

3.10 Breast exoprosthesis model 1086X


Figure 10. Appearance of the model 1086X breast exoprosthesis

Used on both sides.

The shape of the exoprosthesis imitates the shape of the mammary gland.

The design of the rear surface of the exoprosthesis is ribbed to ensure air circulation.

When swimming, water drains faster than other prostheses. Resistant to chlorine and salt water. Can be used while taking a thermal bath.

EQUIPMENT:

Breast exoprosthesis, model 1019X, consisting of:

- 1) Breast exoprosthesis, model 1019X, size: 60/65/70/75/80/85/90/95/100/105/110/115 - 1 pc;
- 2) Instruction manual - 1 pc.

Breast exoprosthesis, model 1050X, consisting of:

- 1) Breast exoprosthesis, model 1050X, size: 65/70/75/80/85/90/95/100/105/110/115 - 1 pc;
- 2) Instruction manual - 1 pc.

Breast exoprosthesis, model 1051X, consisting of:

- 1) Breast exoprosthesis, model 1051X, size: 60/65/70/75/80/85/90/95/100/105/110/115/120/125 - 1 pc;
- 2) Instruction manual - 1 pc.

Breast exoprosthesis, model Active 1054X, consisting of:

- 1) Breast exoprosthesis, model 1054X, size: 55/60/65/70/75/80/85/90/95/100/105/110/115 - 1 pc;
- 2) Instruction manual - 1 pc.

Breast exoprosthesis, model 1066X, consisting of:

- 1) Breast exoprosthesis, model 1066X, size: 65/70/75/80/85/90/95/100/105/110/115 - 1 pc;
- 2) Instruction manual - 1 pc.

Breast exoprosthesis, model 1080 L , consisting of:

- 1) Breast exoprosthesis, model 1080 L , size: 65/70/75/80/85/90/95/100/105/110/115 - 1 pc;
- 2) Instruction manual - 1 pc.

Breast exoprosthesis, model 1080 R , consisting of:

- 1) Breast exoprosthesis, model 1080 R , size: 65/70/75/80/85/90/95/100/105/110/115 - 1 pc;
- 2) Instruction manual - 1 pc.

Breast exoprosthesis, model 1081 L , consisting of:

- 1) Breast exoprosthesis, model 1081 L , size: 55/60/65/70/75/80/85/90/95/100/105/110/115/120 - 1 pc;
- 2) Instruction manual - 1 pc.

Breast exoprosthesis, model 1081 R , consisting of:

- 1) Breast exoprosthesis, model 1081 R , size: 55/60/65/70/75/80/85/90/95/100/105/110/115/120 - 1 pc;
- 2) Instruction manual - 1 pc.

Breast exoprosthesis model 1086 X consisting of:

- 1) Breast exoprosthesis, model 1086 X , size: 55/60/65/70/75/80/85/90/95/100/105/110/115 - 1 pc;
- 2) Instruction manual - 1 pc.

4 TECHNICAL DATA




Operational and technical characteristics of exoprostheses correspond to table 1.

The exoprostheses are suitable for use at temperatures from -30°C to +50°C, relative air humidity up to 85% and atmospheric pressure from 70 to 106 kPa.

Exoroses are resistant to repeated manual cleaning with soap and water.

Service life - 10 years.

Table 1

Name of the exoprosthesis model	For bras sizes according to bust, cm	The size is similar to the size of the bra cup	Volume control yes/no	Air circulation yes/no	Color options (code)	Weight, g	Structure, single layer / double layer	The form
1	2		3	4	5	6	7	8
1019X	60/65-110/115	AA (0), A (1), B (2), C (3), D (4), E (5), F (6), G (7)	Yes	Not	007	20	single layer	
1050X	65-115	AA (0), A (1), B (2), C (3), D (4), E (5), F (6), G (7)	Not	Not	007	227	double layer	
1051X	60-125	AA (0), A (1), B (2), C (3), D (4), E (5), F (6), G (7)	Not	Not	007	227	single layer	

1054X	55-115	AA (0), A (1) , B (2), C (3), D (4), E (5), F (6), G (7)	Not	Yes	007 736	227	single layer	
1066X	65-115	AA (0), A (1) , B (2), C (3), D (4), E (5), F (6), G (7)	Not	Yes	007	227	double layer	
1080 L	65-115	AA (0), A (1) , B (2), C (3), D (4), E (5), F (6), G (7)	Not	Not	007	227	double layer	
1080 R	65-115	AA (0), A (1) , B (2), C (3), D (4), E (5), F (6), G (7)	Not	Not	007	227	double layer	
1081L	55-120	AA (0), A (1) , B (2), C (3), D (4), E (5), F (6), G (7)	Not	Not	007	180	single layer	

1081R	55-120	AA (0), A (1) , B (2), C (3) , D (4), E (5), F (6), G (7)	Not	Not	007	180	single layer	 
1086 X	55-115	AA (0), A (1) , B (2), C (3) , D (4), E (5), F (6), G (7)	Not	Yes	300	227	single layer	 

Note - Tolerances for all characteristics are $\pm 10\%$, unless otherwise indicated.

1. For models 1050X, 1051X, 1044X, 1066X, 1080L, 1080R, 1081L, 1081R, 1086X: The tensile strength of the shell at break must be at least 1 MPa. The relative elongation at rupture of the shell must be at least 150%, the residual elongation must not exceed 40%.
2. For model 1019X: The breaking load of the seam should be: along 140 ± 10 N / across 90 ± 10 N.
3. The exoprosthesis should not collapse under compression under the influence of a load of 5 kg for 10 minutes: there should not be a rupture of the shell, the formation of cracks, delamination of the connecting elements of the shell and the sealing lining.

5 MATERIALS USED IN MANUFACTURING

Exoprostheses are made from the materials presented in Table 2.

Table 2

No.	Name of the exoprosthesis model	Material
1	1019X	OEKO -TEX® polyester , 35% OEKO -TEX® polyurethane , 5% OEKO -TEX® elastane
2	1050X	100% OEKO brand silicone - TEX ®
3	1051X	100% OEKO brand silicone - TEX ®
4	1054X	100% OEKO brand silicone - TEX ®
5	1066X	100% OEKO brand silicone - TEX ®
6	1080 L	100% OEKO brand silicone - TEX ®
7	1080 R	100% OEKO brand silicone - TEX ®
8	1081L	100% OEKO brand silicone - TEX ®
9	1081R	100% OEKO brand silicone - TEX ®
10	1086 X	100% OEKO brand silicone - TEX ®

To fill the exoprosthesis model 1019X, OEKO -TEX ® cotton is used

In the manufacture of exoprostheses, flesh dyes, brown dyes, all dyes of the OEKO -TEKS ® trademark are used.

All color options are provided by mixing basic dyes.

6 USING THE EXOPROSTHESIS

6.1 Open the package and check the exoprosthesis for any mechanical damage or contamination.

If any problems are found, contact the manufacturer's authorized representative immediately.

Warning

• **Disposal of packaging materials must be carried out in accordance with your local regulations.**

6.2 The product is intended for use by women who have undergone breast surgery.

6.3 In addition to imitation and optical balance, the exoprosthesis maintains body symmetry (prevention of shoulder displacement) due to its own weight.

6.4 How to use:

- Before use, check the product for dirt and mechanical damage.
- Fix the exoprosthesis in a special bra.
- Make sure the exoprosthesis is fixed.
- Put on a bra with the exoprosthesis(s).

7 CARE AND CLEANING

If the product becomes dirty, it must be cleaned manually using water and soap. It is forbidden to use chemicals when cleaning.

After cleaning, dry the product on a flat surface.

For model 1019 X, remove the filler before cleaning, and after cleaning, blot the wet surface with a paper or ordinary towel, and then dry on a flat surface .

8 MARKING


8.1 Marking applied to the exoprosthesis, except for models 1019X and 1086 X (marking is not applied to exoprostheses models 1019X and 1086 X):

- manufacturer's logo;
- exoprosthesis model;
- the size;
- name and address of the manufacturer;
- CE marking.

8.2 Marking applied to the packaging of the exoprosthesis:

- product name;
- exoprosthesis model;
- manufacturer's logo;
- serial number;
- date of manufacture;
- composition of the material;
- the size;
- name and address of the manufacturer;
- CE marking;
- Disposal instructions;
- Barcode;
- QR code;
- Information about the authorized representative;
- number and date of the registration certificate.

8.3 Decoding of all marking symbols:

Symbol	Meaning	Description
	Serial number	Product serial number
	Manufacturer	Name of the manufacturer of the product
	date of manufacture	Product production date
	C E marking	C E marking

9 PACKAGING

The exoprosthesis is packed in a consumer box. The consumer box has a plastic liner with a recess for the exoprosthesis. The exoprosthesis is placed in a plastic insert.

The consumer box is sealed with a polyethylene sticker with an adhesive layer so that the package cannot be opened without breaking the integrity.

10 SAFETY REQUIREMENTS AND ENVIRONMENTAL PROTECTION

10.1 Safety requirements

Serviceability in accordance with IEC 62366.

The exoprosthesis materials must be biologically safe and comply with the requirements of the ISO 10993 group of standards.

10.2 Environmental protection

The materials used in the manufacture of the exoprosthesis should not pose a danger to life, human health and the environment both during operation and after the end of its service life.

The main type of possible hazardous impact on the environment is pollution of the atmospheric air of populated areas, soils and waters as a result of:

- emergency leaks (placers) of production materials;
- unorganized waste disposal on the territory of the manufacturer or outside it;
- random dumping of waste in places not intended for this purpose.

11 TRANSPORT AND STORAGE

11.1 Transport conditions

Transportation of the exoprosthesis can be carried out by all types of covered vehicles and the rules for the carriage of goods in force for each type of transport.

Transportation of the exoprosthesis by sea should be carried out in accordance with the "Safety Rules for the Sea Transportation of General Cargoes". Type of shipment - containers and small shipment.

The exoprosthesis, during transportation, must be resistant to climatic influences, at temperatures from -40°C to $+50^{\circ}\text{C}$, relative air humidity up to 90% and atmospheric pressure from 70 to 106 kPa.

The exoprosthesis, during transportation, must be resistant to mechanical stress, with the following parameters:

- vibration loads - frequency range 10-55 Hz with a displacement amplitude of 0.15 mm;
- shock loads - peak shock acceleration 100 ms⁻² with a duration of 16 ms.

11.2 Storage conditions

The exoprosthesis, during storage, must be resistant to climatic influences, at temperatures from -40°C to $+50^{\circ}\text{C}$, relative air humidity up to 90% and atmospheric pressure from 70 to 106 kPa.

12 MAINTENANCE AND REPAIR

The exoprosthesis is not subject to maintenance, except for cleaning, and repair.

13 DISPOSAL INFORMATION

Disposal of the exoprosthesis must be carried out in accordance with the disposal rules in force on the territory of the Russian Federation, according to waste class A, belongs to class A - epidemiological safe waste.

The exoprosthesis must be disposed of if:

- end of service life;
- confirmation of facts and circumstances that pose a threat to the life and health of the user and indicate failure to fulfill the intended purpose.

Recycling is also subject to all packaging, including transport. Paper, wood, polyethylene and plastic must be disposed of separately in accordance with the rules for the collection, accounting and disposal of packaging.

14 MANUFACTURER WARRANTY

The exoprosthesis is guaranteed for a period of 12 months from the date of sale.

It is guaranteed that the exoprosthesis is free from defects associated with the production and quality of production materials.

The manufacturer or an authorized representative of the manufacturer in the territory of the Russian Federation shall not be liable for any personal / material damage caused by erroneous use, operation for purposes other than the intended purpose, and negligent cleaning.

Authorized representative of the manufacturer in the territory of the Russian Federation: Limited Liability Company "MEDPARTNER" (LLC "MEDPARTNER") at the address 194100, St. Petersburg, Bolshoy Sampsonievsky Prospekt, 68, letter N, office 268; phone +7 812 200 94 15; email hi @ medpartner . pro .

WARRANTY CARD

Product category _____
 Manufacturer _____
 Version _____ Stamp _____
 Serial number/art. _____ about sale
 Date of sale _____
 Organization of the seller _____
 Seller's signature _____
 The product has been tested and has no damage.
 I have read and agree to the terms of use.
 Buyer's signature _____
 The warranty period is 12 months from the date of sale.

For warranty issues, please contact: Limited Liability Company "MEDPARTNER" (LLC "MEDPARTNER") at 194100, St. Petersburg, Bolshoi Sampsonievsky Prospekt, 68, letter H, office 268; phone +7 812 200 94 15; email hi @ medpartner . pro .

ATTENTION! Make sure that the warranty card is filled out correctly. The serial number of the product / article must exactly match the number in the coupon. All columns of the coupon (with the exception of the columns for the guarantee center) must be filled in.

In case of inaccurate or incomplete filling of the coupon, the warranty for the product is considered lost.

The manufacturer's warranty is carried out in accordance with the Civil Code of the Russian Federation on a general basis.

The product must not bear traces of mechanical or thermal effects, as well as traces of exposure to various chemical agents. The product may only be operated in accordance with the instruction manual for the product.

WARRANTY STATEMENT

THE DATE	DESCRIPTION	WARRANTY CENTER

15 COMPLAINT DETAILS

In case of complaints about the product, contact: Limited Liability Company "MEDPARTNER" (LLC "MEDPARTNER") at 194100, St. Petersburg, Bolshoy Sampsonievsky Prospekt, 68, letter N, office 268; phone +7 812 200 94 15; email hi @ medpartner . pro .

The list of regulatory documentation used by the manufacturer to confirm compliance:

EN 1041 , Information provided by the manufacturer of medical devices
ISO 13485 , Medical devices - Quality management systems - Requirements for regulatory purposes
ISO 14971 , Medical devices - Application of risk management to medical devices
ISO 15223-1 , Medical devices - Symbols to be used on medical device labels, markings and information to be provided - Part 1: General requirements
IEC 62366 , Medical devices - Serviceability applications
ISO 10993-1 , Biological evaluation of medical devices - Part 1: Evaluation and testing as part of the risk management process
ISO 10993-5 , Biological evaluation of medical devices - Part 5: Tests for cytotoxicity by the in method vitro
ISO 10993-10 , Biological evaluation of medical devices - Part 10: Skin irritation and sensitization tests
EN 980 - Symbols used on the labeling of medical devices

File.Nr. 0136 / 2023

City: Rosenheim
Country: Germany/Bavaria

I hereby certify that the document affixed and sealed hereto was signed in my presence by

Mr. Johannes Weber-Unger,
date of birth: 16. August 1982,
place of residence: Pienzenauerstr. 15, 6330 Kufstein, Austria,
personally known to me, Notary out of service,


acting on behalf of

ANITA Dr. Helbig GmbH
with registered office at Grafenstr. 23, 83098 Brannenburg, Germany.

as this company's managing director (Geschäftsführer).

Upon my inspection of the Commercial Register on today, I further certify that ANITA Dr. Helbig GmbH is registered with the Municipal Court of Traunstein - Commercial Register - under No. HRB 746 and that Mr. Johannes Weber-Unger is entitled to act individually as this company's legal representative.

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