80

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CERTIFICATE

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China Council for the Promotion of International Trade is China Chamber of International Commerce

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China Council for the Promotion of International Trade
China Chamber of International Commerce

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THIS IS TO CERTIFY THAT: the seal of CHIFENG TIANBO MEDICAL DEVICE CO., LTD on the annexed DOCUMENT is genuine.

China Council for the Promotion
of International Trade

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
Authorized
Signature: Xiao Guang

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(Date: Dec. 13, 2022)

证明书验证网址 Website for verifying the certificate: <http://www.rzccpit.com/validate.html>

赤峰天波医疗器材有限责任公司
CHIFENG TIANBO MEDICAL DEVICE CO., LTD.
ADD: Zhongjing Industrial park, Economic Development Area,
Ningcheng District, 024200 Chifeng, Inner Mongolia, China
Tel/ Fax: 010-476 4243071; E-mail: 1215460574@qq.com

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I affirm
Name: Wang Yu Fen
Director
CHIFENG TIANBO MEDICAL DEVICE CO.,
LTD.
Date: November 7, 2022

INSTRUCTIONS FOR USE OF MEDICAL DEVICE

«Sterile syringe for single use with needle SALUVIT®»

Document number: IFU-02. Ver. No.: 2

«CHIFENG TIANBO MEDICAL DEVICE CO., LTD.»

NOTIFICATION: The information contained here is CONFIDENTIAL and PROTECTED BY PROPRIETARY RIGHTS. The contents of this document may not be reproduced or made available to third parties without the prior consent of CHIFENG TIANBO MEDICAL DEVICE CO., LTD., and should not be used in any unauthorized way.

This technical document has been prepared to provide information for the purpose of registering a medical device on the territory of the Russian Federation. The information is provided to the extent required in accordance with the Rules for Registration of Medical devices (Decree of the Government of the Russian Federation No. 1416 of 27.12.2012), Order of the Ministry of Health of the Russian Federation No. 11n of 19.01.2017 "On approval of requirements for the content of technical and operational documentation of the manufacturer (manufacturer) of a medical device".

**INSTRUCTIONS FOR USE
OF MEDICAL DEVICE**
Sterile syringe for single use with needle SALUVIT®

PRODUCT NAME

Sterile syringe for single use with needle SALUVIT®.
(hereinafter referred to as "syringe" or "product").

VARIANTS OF THE MEDICAL DEVICE AND DELIVERY SET:

Sterile syringe for single use with needle SALUVIT®. Variants:

1. Sterile syringe for single use with needle SALUVIT® tuberculin volume 1 ml: Type 1, needle size: 0,40 x 12 mm (27G x 1/2") – 1 pcs., brief instructions for use - 1 pc. or 180 pcs. / pack., brief instructions for use - 1 pc. (applied to the group packaging).
2. Sterile syringe for single use with needle SALUVIT® tuberculin volume 1 ml: Type 1, needle size: 0,45 x 12 mm (26G x 1/2") – 1 pcs., brief instructions for use - 1 pc. or 180 pcs. / pack., brief instructions for use - 1 pc. (applied to the group packaging).
3. Sterile syringe for single use with needle SALUVIT® tuberculin volume 1 ml: Type 2, needle size: 0,40 x 12 mm (27G x 1/2") – 1 pcs., brief instructions for use - 1 pc. or 180 pcs. / pack., brief instructions for use - 1 pc. (applied to the group packaging).
4. Sterile syringe for single use with needle SALUVIT® tuberculin volume 1 ml: Type 2, needle size: 0,45 x 12 mm (26G x 1/2") – 1 pcs., brief instructions for use - 1 pc. or 180 pcs. / pack., brief instructions for use - 1 pc. (applied to the group packaging).
5. Sterile syringe for single use with needle SALUVIT® volume 2 ml: Type 1, needle size: 0,6 x 25 mm (23G x 1") – 1 pcs., brief instructions for use - 1 pc. or 120 pcs. / pack., brief instructions for use - 1 pc. (applied to the group packaging).
6. Sterile syringe for single use with needle SALUVIT® volume 2 ml: Type 1, needle size: 0,6 x 30 mm (23G x 1 1/4") – 1 pcs., brief instructions for use - 1 pc. or 120 pcs. / pack., brief instructions for use - 1 pc. (applied to the group packaging).
7. Sterile syringe for single use with needle SALUVIT® volume 2 ml: Type 2, needle size: 0,6 x 25 mm (23G x 1") – 1 pcs., brief instructions for use - 1 pc. or 120 pcs. / pack., brief instructions for use - 1 pc. (applied to the group packaging).
8. Sterile syringe for single use with needle SALUVIT® volume 2 ml: Type 2, needle size: 0,6 x 30 mm (23G x 1 1/4") – 1 pcs., brief instructions for use - 1 pc. or 120 pcs. / pack., brief instructions for use - 1 pc. (applied to the group packaging).
9. Sterile syringe for single use with needle SALUVIT® volume 3 ml: Type 1, needle size: 0,6 x 25 mm (23G x 1") – 1 pcs., brief instructions for use - 1 pc. or 120 pcs. / pack., brief instructions for use - 1 pc. (applied to the group packaging).
10. Sterile syringe for single use with needle SALUVIT® volume 3 ml: Type 1, needle size: 0,6 x 30 mm (23G x 1 1/4") – 1 pcs., brief instructions for use - 1 pc. or 120 pcs. / pack., brief instructions for use - 1 pc. (applied to the group packaging).
11. Sterile syringe for single use with needle SALUVIT® volume 3 ml: Type 2, needle size: 0,6 x 25 mm (23G x 1") – 1 pcs., brief instructions for use - 1 pc. or 120 pcs. / pack., brief instructions for use - 1 pc. (applied to the group packaging).
12. Sterile syringe for single use with needle SALUVIT® volume 3 ml: Type 2, needle size: 0,6 x 30 mm (23G x 1 1/4") – 1 pcs., brief instructions for use - 1 pc. or 120 pcs. / pack., brief instructions for use - 1 pc. (applied to the group packaging).
13. Sterile syringe for single use with needle SALUVIT® volume 5 ml: Type 1, needle size: 0,7 x 38 mm (22G x 1 1/2") – 1 pcs., brief instructions for use - 1 pc. or 90 pcs. / pack., brief instructions for use - 1 pc. (applied to the group packaging).
14. Sterile syringe for single use with needle SALUVIT® volume 5 ml: Type 1, needle size: 0,8 x 38 mm (21G x 1 1/2") – 1 pcs., brief instructions for use - 1 pc. or 90 pcs. / pack., brief instructions for use - 1 pc. (applied to the group packaging).
15. Sterile syringe for single use with needle SALUVIT® volume 5 ml: Type 2, needle size: 0,7 x 38 mm (22G x 1 1/2") – 1 pcs., brief instructions for use - 1 pc. or 90 pcs. / pack., brief instructions for use - 1 pc. (applied to the group packaging).
16. Sterile syringe for single use with needle SALUVIT® volume 5 ml: Type 2, needle size: 0,8 x 38 mm (21G x 1 1/2") – 1 pcs., brief instructions for use - 1 pc. or 90 pcs. / pack., brief instructions for use - 1 pc. (applied to the group packaging).
17. Sterile syringe for single use with needle SALUVIT® volume 10 ml: Type 1, needle size: 0,8 x 38 mm (21G x 1 1/2") – 1 pcs., brief instructions for use - 1 pc. or 60 pcs. / pack., brief instructions for use - 1

- pc. (applied to the group packaging).
18. Sterile syringe for single use with needle SALUVIT® volume 10 ml: Type 2, needle size: 0,8 x 38 mm (21G x 1 1/2") – 1 pcs., brief instructions for use - 1 pc. or 60 pcs. / pack., brief instructions for use - 1 pc. (applied to the group packaging).
 19. Sterile syringe for single use with needle SALUVIT® volume 20 ml: Type 1, needle size: 0,8 x 38 mm (21G x 1 1/2") – 1 pcs., brief instructions for use - 1 pc. or 36 pcs. / pack., brief instructions for use - 1 pc. (applied to the group packaging).
 20. Sterile syringe for single use with needle SALUVIT® volume 20 ml: Type 2, needle size: 0,8 x 38 mm (21G x 1 1/2") – 1 pcs., brief instructions for use - 1 pc. or 36 pcs. / pack., brief instructions for use - 1 pc. (applied to the group packaging).
 21. Sterile syringe for single use with needle SALUVIT® volume 30 ml: Type 1, needle size: 0,8 x 38 mm (21G x 1 1/2") – 1 pcs., brief instructions for use - 1 pc. or 30 pcs. / pack., brief instructions for use - 1 pc. (applied to the group packaging).
 22. Sterile syringe for single use with needle SALUVIT® volume 30 ml: Type 2, needle size: 0,8 x 38 mm (21G x 1 1/2") – 1 pcs., brief instructions for use - 1 pc. or 30 pcs. / pack., brief instructions for use - 1 pc. (applied to the group packaging).
 23. Sterile syringe for single use with needle SALUVIT® volume 50 ml: Type 1, needle size: 1,2 x 38 mm (18G x 1 1/2") – 1 pcs., brief instructions for use - 1 pc. or 15 pcs. / pack., brief instructions for use - 1 pc. (applied to the group packaging).
 24. Sterile syringe for single use with needle SALUVIT® volume 50 ml: Type 2, needle size: 1,2 x 38 mm (18G x 1 1/2") – 1 pcs., brief instructions for use - 1 pc. or 15 pcs. / pack., brief instructions for use - 1 pc. (applied to the group packaging).
 25. Sterile syringe for single use with needle SALUVIT® volume 60 ml: Type 1, needle size: 1,6 x 38 mm (16G x 1 1/2") – 1 pcs., brief instructions for use - 1 pc. or 15 pcs. / pack., brief instructions for use - 1 pc. (applied to the group packaging).
 26. Sterile syringe for single use with needle SALUVIT® volume 60 ml: Type 2, needle size: 1,6 x 38 mm (16G x 1 1/2") – 1 pcs., brief instructions for use - 1 pc. or 15 pcs. / pack., brief instructions for use - 1 pc. (applied to the group packaging).

3. MANUFACTURER

Company name: CHIFENG TIANBO MEDICAL DEVICE CO., LTD.

Location address: Zhongjing Industrial Park, Economic Development Area, Ningcheng District 024200 Chifeng, Inner Mongolia, China

4. INTENDED USE

Intended use: Syringe is intended for the introduction of liquids/gases (for example, medicines) into a medical device or body or the extraction of liquids /gases from a medical devices /body (i.e., used for both administration and suction) and/or conducting a sensitivity test to allergen.

5. INDICATIONS FOR USE

- Injection of drugs to the patient
- Carrying out medical and diagnostic procedures

6. PRECAUTIONS WHEN USING

- 1) The product usage must comply with the applicable requirements of the medical sector and relevant regulations, limited to use only by qualified doctors or medical personnel, or independently, if necessary, by trained patients;
- 2) The products are delivered in a sterile form. Repeated sterilization of products is prohibited.
- 3) Do not use if the package is damaged.
- 4) This product is "disposable", subject to destruction after use in accordance with the rules for the disposal of medical waste to be disposed of, reuse or other use is prohibited.

7. CONTRADICTIONS FOR THE PRODUCT USAGE AND POSSIBLE SIDE EFFECTS:

No

Warnings!

A single-use product, reuse is prohibited!

Do not use if the tightness of the individual packaging is broken!

8. QUALIFICATION OF PERSONNEL

The product is applicable in nursing, general medical practice of medical organizations in the provision of inpatient, outpatient, emergency, emergency medical care by medical organizations and at home during injections by medical personnel or independently, if necessary, by trained patients.

10. TIME OF CONTACT WITH THE BODY

The product is intended for injecting drugs into the human body, as well as aspiration of liquid contents and flushing of the patient's body cavities.

11. CATEGORIES OF PATIENTS

- Age: children and adults.
- Body weight: not applicable.
- Health status: not applicable.
- Nationality: different.
- Patient's condition: The patient is not a user of the product, if necessary could be used independently by trained patients.

11. BASIC CONDITIONS OF USE


It is used for subcutaneous, intradermal, intramuscular and intravenous injections with precise dosing of injectable agents, as well as for aspiration of liquid contents and flushing of body cavities.

The recommended application temperature of the product is from $+10^{\circ}\text{C}$ to $+40^{\circ}\text{C}$, and the relative humidity of the air is not more than 80%.

12. PRODUCT DESCRIPTION

12.1 Product design

12.1.1. The syringe consists of the following parts: a transparent cylinder, a plunger with a sealing ring (cuff) and a needle with a protective cap. The designation of the components of the syringes are shown in Figure 1.

1 - protective cap, 2- needle tube, 3 - needle head, 4 - Syringe tip, 5- Cuff (sealing ring), 6 - Scale, 7 - cylinder, 8 - Fingers stoppers, 9 - plunger, 10 plunger stopper

Figure 1. Syringe scheme

12.1.2. The syringe graduation is contrasting, clear, indelible, applied in a contrasting color. Provides the most accurate dosage.

12.1.3. Syringes are divided into 2 types: Type 1 - syringes with normal graduation, Type 2 - syringes with additional graduation and inscription "DISCARD AFTER USE".


Figure 2. Type 1 - Scale without additional graduation


Figure 3. Type 2 - Scale with additional graduation

12.1.4. The plunger on the end part has a rubber sealing ring (cuff), which provides a well-visible and clear edge that acts as a reference line at the end of the piston due to close contact with the inner surface of the cylinder, ensures its sliding without jerks and provides smooth injection of the injection agent. The design of the plunger allows to activate it without difficulty and, while holding the syringe with one hand, down it with the thumb of the same hand. The thrust of the plunger prevents the fingers from slipping during the operation of the syringe.

12.1.5. Atraumatic injection needle is strong, straight, with variable cross-section and wall thickness, with a long cut and triangular laser sharpening of the tip at an angle. The needle head is identified by a color code according to the international standard ISO 6009.

12.1.6 Actual length of needle tube

The actual length of the needle tube should be equal to the nominal length within the tolerances specified in Table 1.

Table 1 (to the figure 4)

Nominal length of needle tube, mm (13)	Tolerance
12, 25	+1 -2
30, 38	+1,5 -2,5


Figure 4. Schematic representation of the needle

12.1.7 The nominal outer diameter of the needles are indicated in Table 2

Table 2

Indicated metric size, mm	Caliber	ODmin, mm	ODmax, mm	Wall	IDmin	IDmax	Color coding
0,40	27	0.400	0.420	RW	0.184	0.240	Medium grey
				TW	0.241	0.350	
0,45	26	0.440	0.470	RW	0.232	0.291	Brown
				TW	0.292	0.390	
0,60	23	0.600	0.673	TW	0.370	0.459	Deep blue
0,70	22	0.698	0.730	TW	0.440	0.521	Black
0,80	21	0.800	0.830	TW	0.547	0.609	Deep green
1,20	18	1.200	1.300	TW	0.910	1.040	Pink
1,60	16	1.600	1.690	TW	1.283	1.389	White

The main technical parameters of the product

Table 3

No.	Quality indicator name	Standard
STERILE SYRINGE FOR SINGLE USE WITH NEEDLE		
1	Purity	The surface of the syringe in contact with the injection fluid is clean during normal use and has no foreign particles during visual inspection with normal or adjusted to normal vision, without the use of magnification and with illumination from 300 lux to 700 lux.
2	Limits of acidity or alkalinity	The pH level of the extract should correspond to the pH of the control liquid within one unit when measured by a laboratory pH meter with a standard electrode.
3	Limit content of extracted metals	The syringe extract should contain no more than 5 mg/l of the total amount of lead, tin, zinc and iron when tested by a generally accepted microanalytical method, for example, by atomic absorption. The cadmium content in the extract should not exceed 0.1 mg/l.
4	Lubrication	Lubrication of the inner surface of the syringe cylinder, including the plunger, should not be noticeable in the form of drops when observed with normal or adjusted to normal vision.
5	Finger stoppers	The open end of the cylinder must have finger stoppers, Finger stoppers must have the appropriate size, shape and strength that meet their purpose and ensure reliable retention of the syringe when it is used. Finger stoppers should not have burrs and sharp edges
6	Tip	The tip of the syringe should have a 6% conical connection (Luer)
7	Tip hole	The diameter of the central hole of the tip should be at least 1.2 mm.
8	Lubrication of the needle tube	The amount of lubricant does not form liquid droplets on the outer and inner surfaces of the needle when viewed with normal, or adjusted to normal, vision.
9	Needle point	The needle tip should be sharp, without burrs when viewed with a 2.5-fold magnification. Needle sharpening angle - $11^{\circ} \pm 2^{\circ}$.

12.3. The nominal capacity of the syringe, the value of the "dead" space, the scale dimensions are shown in Table 4.

Table 4

Nominal capacity of syringe V, ml	Maximum "dead" space, ml	The total length of the scale up to the mark of the nominal capacity, not less than, mm	Scale interval, ml	Capacity between scale lines to be numbered, ml
$V < 2$	0,07	57	0,05	0.1
$2 \leq V < 5$	0,07	27	0,2	0,5 or 1
$5 \leq V < 10$	0,075	36	0,5	1
$10 \leq V < 20$	0,1	44	1,0	5
$20 \leq V < 30$	0,15	52	2,0	10
$30 \leq V < 50$	0,17	67	2,0	10

Nominal capacity of syringe V, ml	Maximum "dead" space, ml	The total length of the scale up to the mark of the nominal capacity, not less than, mm	Scale interval, ml	Capacity between scale lines to be numbered, ml
50 ≤ V	0,2	75	5,0	10

12.4. Table 5 contains information about the weight of the product

Table 5

Nominal capacity of the syringe, ml	The weight of the syringe without a needle, not more than, g	Weight of a syringe with a needle, not more than, g
1	2.5	3.2
2	3.3	4
3	3.4	4.1
5	4.8	5.6
10	7.6	8.4
20	12	13
30	19	20
50	35	36
60	36	37

13. OPERATION PROCEDURES

1. Prepare the injection site by treating it with a cloth with a disinfectant.
2. Open the package and take out the product.
3. Inject into the prepared place, observing the injection technique and in accordance with the injection standard, depending on the type of injection (subcutaneous, intradermal, intramuscular, intravenous).
4. Remove the syringe with the needle. Lightly press the injection site with a cloth with a disinfectant solution. It is permissible for a drop of blood to appear at the injection site of the needle.
5. If it is necessary to aspirate the liquid contents from the body cavities, proceed to suction after removing the syringe from the package.
6. If it is necessary to flush the body cavities, pre-fill the syringe with a solution.

14. PACKAGING

Individual packaging

Each syringe is hermetically packed in an individual polyethylene package. The product is packed in an individual sealed package with heat-sealable seams. The materials used for packaging do not have harmful effects on the contents.

The material and design of the packaging of the product ensure the sterility of the contents during transportation, storage and minimizes the risk of contamination of the contents during the opening process.

Middle packaging

Syringes in individual packaging can be packed in a middle package that protects the contents during transportation, storage and handling. Each group package has a brief instruction for use.

Transport packaging

Syringes in individual packaging or in a middle package must be packed in a transport package for storage and transportation. The transport package is a box made of corrugated cardboard.

15. STORAGE AND TRANSPORTATION CONDITIONS

Storage conditions: The product must be stored closed at a temperature from 0 ° C to + 40 ° C, with a relative humidity of no more than 80%, excluding the ingress of corrosive gas and direct sunlight.

Transportation conditions: The product must be transported at temperatures from -25 ° C to +40 ° C in covered vehicles of all types.

Fragile, handle carefully. Avoid exposure to sunlight. Keep away from moisture.

16. SHELF LIFE

The shelf life of the product is 5 years from the date of manufacture, subject to storage and transportation conditions. It is forbidden to use the product after the expiration date.

17. WARRANTY OBLIGATIONS

The manufacturer guarantees that the products comply with technical requirements during the entire shelf life (5 years), subject to the consumer's storage and transportation conditions.

CHIFENG TIANBO MEDICAL DEVICE CO., LTD. guarantees the high quality of the product and the materials used for its manufacture. CHIFENG TIANBO MEDICAL DEVICE CO., LTD. is not responsible for any incidental, special or indirect costs, damage or additional costs associated, directly or indirectly, with the use of this product.

CHIFENG TIANBO MEDICAL DEVICE CO., LTD. expressly declares the inadmissibility of reuse of this product, in which case the warranty and other implied or direct obligations (including the warranty of merchantability and the possibility of use for special purposes) are considered invalid.

CHIFENG TIANBO MEDICAL DEVICE CO., LTD. also, it is not responsible for any accidental or indirect damage that may occur as a result of repeated use of the product. For any questions regarding the described products, please contact directly CHIFENG TIANBO MEDICAL DEVICE CO., LTD., or an official representative of the manufacturer.

18. AUTHORIZED REPRESENTATIVE OF THE MANUFACTURER IN THE RUSSIAN FEDERATION

Limited liability company «RAFEL» («RAFEL» LLC)

Location address: Russia, 420108, Republic of Tatarstan, Kazan, Mazita Gafuri, 71, manufacturing bld. no. 1

Telephone: 8(843)211-12-21

E-mail: info@rafel.ru

19. COMPLIANCE STANDARDS IN THE RUSSIAN FEDERATION

GOST R ISO 15223-1, GOST ISO 7886-1, GOST ISO 11607-1, GOST ISO 7864, GOST R ISO 6009, GOST EN 556-1, standards series GOST ISO 10993.

20. STERILITY

Disposable sterile syringes with a needle sterilized with ethylene oxide (sterility level (SAL) < 10⁻⁶).

The product has certain storage requirements, if the storage conditions are met, the product remains sterile for five years.

The finished product is sterile, non-pyrogenic, non-toxic. Repeated sterilization of the product is prohibited. Repeated use is prohibited.

21. CLEANING AND DISINFECTION

This product does not need cleaning and disinfection.

22. MAINTENANCE AND CALIBRATION

The product is used by trained doctors or independently by trained patients and does not require maintenance and calibration.

23. DISPOSAL

This product is "disposable". After use, the product is subject to destruction as epidemiologically hazardous waste, in accordance with the national rules for the disposal of medical waste (SanPiN 2.1.3684: Class B). Reuse or other use of the product is prohibited.

Products that have become unusable or expired are subject to destruction in accordance with SanPiN 2.1.3684 as waste belonging to class "A" (epidemiologically safe waste, close to MSW), in a way that prevents their reuse.

24. SYMBOLS ON THE MARKING OF A MEDICAL DEVICE

Prohibition of reuse

Fragile, handle carefully

Date of manufacture

Avoid exposure to sunlight

Use before

Protect from moisture

LOT

Batch code

Temperature range

Do not use if the
package is damaged

Sterilization with
ethylene oxide

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