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Instruction for use
of medical device

**blueflow venous self-expanding stent with the
delivery system**

«APPROVED»

Contract Medical International GmbH

Renate Westermann
Senior Regulatory Affairs Manager

Ver. 02, 2023

Medical device name

blueflow venous self-expanding stent with the delivery system:

1. blueflow venous self-expanding stent with the delivery system, standard sizes: 12 mm x 60 mm, 12 mm x 100 mm, 14 mm x 60 mm, 14 mm x 100 mm, 14 mm x 150 mm, 16 mm x 60 mm, 16 mm x 100 mm, 16 mm x 150 mm, 18 mm x 60 mm, 18 mm x 100 mm – 1 pc.;
2. Implant card – 1 pc.;
3. Instructions for use – 1 pc.

Developer/Manufacturer/Manufacturing facility

Contract Medical International GmbH

Lauensteiner Str. 37, 01277 Dresden, Germany

+49 351 213-8888

Email: hmdd.info@heraeus.com

Authorized representative on the territory of the Russian Federation

"SP-CARDIO" LLC

107076, Russia, Moscow,

Matrosskaya Tishina str., h. 23., bld. 1

Tel.: +7 (903) 616-15-71

e-mail: sp_cardio@rambler.ru

Intended Use

The blueflow Venous Stent is intended to treat symptomatic, dilatable obstructions, such as those caused by May-Thurner syndrome (MTS) and other clinical conditions, deep venous thrombosis (DVT), prevent post thrombotic syndrome (PTS) and other forms of iliac vein compression (tumor, radiation) in the common iliac, external iliac or femoral veins.

Contact with human body

The blueflow Venous Stent is an implantable device that has permanent (more than 30 days) contact with blood.

The delivery system is an attaching from the outside, has a short-term (less than 24 hours) contact with circulating blood.

Terms of use, potential consumer

This medical device can only be used by qualified doctors who have the necessary knowledge for the implantation of venous stents in medical and preventive institutions in an operating unit equipped with X-ray equipment and the necessary tools for interventional manipulations.

Indications for use

- symptomatic, dilatable obstructions, such as those caused by May-Thurner syndrome (MTS) and other clinical conditions,
- deep venous thrombosis (DVT),
- prevent post thrombotic syndrome (PTS),
- and other forms of iliac vein compression (tumor, radiation) in the common iliac, external iliac or femoral veins.

Contraindications

The blueflow Venous Stent is contraindicated for use:

- in patients with complete venous occlusion without sufficient dilation to allow insertion of the delivery system.
- in the arterial vasculature.

- for patients with uncorrected bleeding disorders, including patients with conditions that preclude adequate hemostasis at the access site
 - for patients with contraindications to thrombolytic medications and anticoagulants
- Note - Common contraindications concerning the PTA of vascular lesions also apply to vascular stent implantations, especially in case of blood coagulation disorders.*

Warnings

Read the instructions prior to use.

- Flushing needed prior to procedure.
- Nitinol may cause allergic reactions as Nitinol is a nickel titanium alloy. Therefore, a physician should be consulted if the patient has a nickel allergy.
- The blueflow Venous Stent is made of Nitinol. Its compatibility with other stent materials for overlapping/extension has not been tested. Therefore, contact with other stent materials except for Nitinol should be avoided.
- The IP medical Stent is only designed for use in the common and external iliac veins and femoral veins. Use in arterial vasculature is not recommended.
- Do not use the device if the package or contents appear damaged in any way.
- Do not use disinfectant to clean product prior to usage!
- This device is intended for single use only. Do not reuse, re-sterilize or autoclave. Reuse of single-use devices creates a potential risk of patient or user infections. Contamination of the device may lead to injury, illness or death of the patient. Cleaning, disinfection and sterilization may compromise essential material and design characteristics leading to device failure. Contract Medical International GmbH (CMI) will not be responsible for any direct, incidental or consequential damages resulting from re-sterilization or reuse.
- The system is intended to be used by physicians trained in placing venous stents only in a hospital and surgical environment equipped with fluoroscopy and support for interventional procedures.
- Do not attempt to advance the device if resistance is felt. Always use fluoroscopy to determine the cause. If the cause cannot be determined and corrected, discontinue the procedure and withdraw the device. Continued advancement or retraction against resistance may result in serious injury, and/or breakage of the interventional medical device.
- Do not alter this device. Alterations may impair function. Be aware that attempting to alter the shape of this device by twisting or similar methods may compromise the integrity and cause damage to the stent resulting in improper placement of the stent in the vasculature, which may require additional intervention.
- Verify blueflow Venous Stent device and accessory components compatibility prior to use.
- Do not attempt to use a guide wire with a maximum diameter greater than specified on package label. Device damage or breakage may occur.
- This device should only be used by or under the direction of physicians thoroughly trained in the technique of vein stent delivery systems.
- Observe sterile technique at all times when handling and inserting the blueflow Venous Stent.
- Adequate heparinization and coagulation must be maintained during the procedure. The risk of systemic anticoagulation must be weighed against the benefits of the stent placement procedure. The safety and effectiveness of the blueflow Venous Stent has not been established for use in patients for whom antiplatelet, anticoagulation therapy, or thrombolytic drugs are contraindicated or who exhibit coagulopathy.
- The medical techniques and procedures described in these instructions are presented as an example only, and do not represent ALL medically acceptable protocols. They are

not intended as a substitute for the physician's experience and judgment in treating any specific patient.

- Only the stent of the IP Medical system is MR Conditional. Do not attempt to use the stent delivery system under MRI conditions.

Precautions

Depending on the morphological structure of the vein lumen, an angioplasty may be performed prior to stent implantation.

Do not attempt to remove stent from the delivery system before use. Do not use organic solvents (e.g. alcohol) with the delivery system.

Do not advance the delivery system without the guidewire extending distally from the tip.

Potential complications and adverse events

Only physicians familiar with possible complications may implant this stent. Complications may occur both during and following the procedure. Possible complications include:

- Failure of stent deployment
- Failure of stent application system
- Stent mal-implantation – immediate (inadequate anchorage, dislocation due to insufficient stent diameter)
- Stent migration or dislodgement - later
- Stent thrombosis, stent occlusion
- Stent fracture
- Hemorrhage/bleeding
- Allergic reactions or hypersensitive reactions
- Hypotension, shock
- Hematoma (stent site, access site, retroperitoneal)
- Renal failure
- Pain
- Vascular damage (including dissection, rupture, hematoma, pseudo aneurysm) at stent site and vascular access site
- Formation of arterio-venous fistula
- Tissue necrosis
- Organ failure
- Infection
- Fever
- Sepsis
- Limb ischemia
- Hypertension
- Pulmonary embolism, air embolism
- Cerebrovascular dysfunction (stroke, TIA)
- Cardiac and cardiovascular dysfunction (hypertension, myocardial infarction, cardiac arrest)
- Death

DEVICE DESCRIPTION

The blueflow Venous Stent consists of two main components: a blueflow braided stent and a delivery system.

Vascular access is gained via standard percutaneous technique (Seldinger technique) with a 10 French vascular sheath and a 0.035 inch guide wire (maximum 0.95 mm). The minimum

or diameter of the vascular sheath that may be used is 3.38 mm. Fluoroscopy is needed for anatomically correct stent placement.

The blueflow Venous Stent is a self-expanding stent loaded in a stent application system for insertion and deployment in various sizes (Table 1).

Table 1: Size specifications of blueflow Venous Stent

Length (mm)	Diameter (mm)			
	12 mm	14 mm	16 mm	18 mm
60 mm	X	X	X	X
100 mm	X	X	X	X
150 mm		X	X	

The self-expanding stent is made of biocompatible Nitinol, a 55 weight-% Nickel and balanced Titanium alloy with a long safety history in medical device manufacturing and stent design.

Delivery system


The application system (Figure1) is a 10 Fr caliber catheter with a rod length of 100 cm, into which a venous stent is preloaded.

Figure 1: Delivery system

Delivery system of the blueflow Venous Stent medical device consists of following parts (figure 6).

Figure 2: Stent application system with individual components labelled

- (10) White radiopaque catheter tip made of Pebax/Bariumsulfate
- (20) Guide wire lumen made of Polyimide and Pebax
- (30) Blue, stainless steel wire enforced, braided outer shaft made of Pebax
- (40) Stent retriever made of Pebax
- (50) Stent pusher, overmolded, made of Pebax with embedded Pt/Ir marker band
- (60) Natural colored plastic pusher tube made of Pebax
- (70) Molded distal handle made of Pebax
- (80) Stent size, printed with white ink
- (90) Clear Y-Tuohy borst valve with Luer side arm, made of Polycarbonate, silicone and PTFE
- (100) Marker for stent retriever stop, laser-marked
- (110, 120) Stainless steel, inner and outer metal pusher tubes
- (130) Molded proximal handle with Luer, made of Pebax

Magnetic Resonance Imaging (MRI) Safety Information:

Magnetic Resonance Conditional

Non-clinical testing has demonstrated that the blueflow Venous Stents are MR Conditional. Only the stent of the blueflow Venous Stent System is MR Conditional. Do not attempt to use the stent delivery system under MRI conditions.

A patient with this stent implanted can be safely scanned in an MR system meeting the following conditions:

- Static magnetic field of 1.5-Tesla (1.5 T) or 3-Tesla (3 T)
- Maximum spatial field gradient of:
 - o 5,180 G/cm (51.80 T/m) for 1.5 T systems
 - o 2,590 G/cm (25.90 T/m) for 3 T systems
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of:
 - o 1,0 W/kg at 1.5 T and 3 T

RF Heating

Under the scan conditions defined above, blueflow Venous Stent's common iliac, external, iliac and femoral Vein Stents are expected to produce a maximum temperature rise of less than 7.0°C after 15 minutes of continuous scanning.

Caution: The RF heating behaviour does not scale with static field strength. Devices that do not exhibit detectable heating at one field strength may exhibit high values of localized heating at another field strength.

MR Artefact

In non-clinical testing, the image artefact caused by the device extends approximately 0.7 cm from blueflow Venous Stent's common iliac, external iliac and femoral vein stents when imaged with a gradient-echo pulse sequence in a 3 T MRI system.

TECHNICAL SPECIFICATION

blueflow Venous Stent

Overall dimensions of the blueflow Venous Stent

The table (Table 2) shows the overall dimensions of the stent for all versions (types and sizes).

Table 2. Overall dimensions of the blueflow Venous Stent

Stent size	Outer Ø of expanded stent, mm			Length of expanded stent, mm			The ratio of the actual length and diameter of the opened stent
	Nominal	Min	Max	Nominal	Min	Max	
12x60	12,5	12,0	13,0	60,0	55,0	65,0	4,80 (+0,62/-0,57)
12x100	12,5	12,0	13,0	100,0	95,0	105,0	8,00 (+0,75/-0,69)
14x60	14,5	14,0	15,0	60,0	55,0	65,0	4,14 (+0,50/-0,47)
14x100	14,5	14,0	15,0	100,0	95,0	105,0	6,90 (+0,60/-0,56)
14x150	14,5	14,0	15,0	150,0	145,0	155,0	10,34 (+0,73/-0,68)
16x60	16,5	16,0	17,0	60,0	55,0	65,0	3,64 (+0,43/-0,40)
16x100	16,5	16,0	17,0	100,0	95,0	105,0	6,06 (+0,50/-0,47)
16x150	16,5	16,0	17,0	150,0	145,0	155,0	9,09 (+0,60/-0,56)
18x60	18,5	18,0	19,0	60,0	55,0	65,0	3,24 (+0,37/-0,35)
18x100	18,5	18,0	19,0	100,0	95,0	105,0	5,41 (+0,43/-0,41)

Functional characteristics of the blueflow Venous Stent.

Table 3. Functional characteristics of the blueflow Venous Stent

Characteristic	Requirements
Minimum bending radius of the stent, mm	10


Fracture resistance under a plane-parallel load	The load required to cause compression of the stent along the entire length by 50% of the initial distance between the walls should be at least 6 N.
Local compression	The stent must withstand constant deformation under local (point) load. After local compression of the stent applied perpendicular to the longitudinal axis of the device, up to 50% of the diameter, a complete restoration of the geometry of the shape of the stent should be observed.
Radial force	The radial force (expansion and compression) caused by the self-expanding stent is 0,004 – 0,046 N/mm.

Delivery system

Basic dimensions

Depending on the size of the loaded venous stent, delivery systems differ in total length.

Figure 3. Delivery system with size symbols

Letter	Size
A	Total length
B	Effective length (shaft length)
C	Shaft outer diameter
D	Shaft inner diameter

Table 4. Overall and basic dimensions of the delivery system

Stent size, mm	A (mm)	B (mm)	C (mm)	D (mm)
12×60	1252	1000,0	3,3	1,0
12×100	1310			
14×60	1261			
14×100	1321			
14×150	1398			
16×60	1270			
16×100	1333			
16×150	1413			
18×60	1278			
18×100	1349			

Functional characteristics


Table 5. Functional characteristics of the delivery system

Characteristic	Requirements
Joint strength	All connections between individual parts of the delivery system must remain intact after applying a longitudinal force of 35 Nm.
Tensile strength of the tubes	The tubular parts of the delivery system must withstand a force of 35 N.
Resistance to bending/kinking	The delivery system must provide the ability to bend/kink, which may be required when performing access and delivery. With a minimum bending radius of 10 mm, there should be no deterioration in the appearance and operational properties of the delivery system.
Expanding force	No more than 30 N

Luer Connector

There are two Luer connectors in the delivery system: for flushing the lumen of the guidewire and for flushing the lumen of the stent (Y-shaped Tuohy connector).

Figure 4. Luer connectors of the delivery system

Luer connectors comply with the requirements of the ISO 80369-7 standard.

X - ray transparency

Table 6. X-ray transparency of the components of the device

Component	Description
Stent	Visible with conventional radiography
Delivery system	Marker strip is visible with conventional radiography
Tip	Visible with conventional radiography

MATERIALS OF ANIMAL AND (OR) HUMAN ORIGIN

The blueflow Venous Stent and the Delivery System do not contain materials of animal and/or human origin.

DATA FROM THE STATE REGISTER OF MEDICINES

The blueflow Venous Stent and the Delivery System do not contain drugs and pharmaceutical substances.

INSTRUCTIONS FOR USE

Recommended Materials (not included in the set):

- 10Fr introducer sheath
- 0.035" (0,89mm) guidewire
- Sterile syringes with Luer-Lock for flushing
- Sterile, heparinized physiological saline solution

Recommended Materials (not included in the set):

- 10Fr introducer sheath
- 0.035" (0,89mm) guidewire
- Sterile syringes with Luer-Lock for flushing
- Sterile, heparinized physiological saline solution

Ensure compatibility with other medical products used during procedure or already implanted.

Note - Indications, contraindications and application of the above product must always comply with current medical standards. The guidelines and recommendations of the relevant specialist medical associations must be observed.

Stent Selection:

Proper stent sizing is critical for correct positioning and to avoid stent migration, and thrombosis. When choosing the appropriate stent size, take into account the foreshortening of the stent based on the vessel diameter as referenced in the table below.

- Prior to stent implantation, the diameter of the vessel and the lesion length has to be measured.
- The diameter of the stent should be 10-20% larger than the vein diameter.
- The stent should be at least 1cm longer than the obstructive venous lesion (0.5 cm distally; 0.5 cm proximally).
- Overlapping of stents has not been tested. We recommend to use a single stent with the appropriate size.

Table 7. Appropriate stent size based on the vessel diameter

Catheter Size	Dimension of Fully Expanded Stent			Nominal Vessel Diameter and Approximate Implanted Stent Length			
	OD Stent	Nominal length	Actual length	Vessel Diameter	Stent Length	Vessel Diameter	Stent Length
Fr (mm)	mm	mm	mm	mm	mm	mm	mm
10 (3,30)	12	60	61	11,5	64	10	70
		100	104		109		119
		60	64		69		76
	14	100	105	13	113	11,5	124
		150	156		169		185
		60	65		72		81
	16	100	104	15	115	13	129
		150	156		169		190
		60	64		74		81
	18	100	106	16,5	121	15	133

Preparation:

The procedure requires fluoroscopy at all times. The blueflow Venous Stent is aseptically

removed from its packaging immediately prior to the procedure.
Check for damage of package or product: A damaged product must not be used!
In case of vein stenosis or scarring (with decrease in vein diameter) a balloon angioplasty ought to be performed up to normal vessel diameter prior to releasing the self-expanding stent.

Vascular access and introduction of the delivery system:

Vascular access is gained via standard percutaneous technique (Seldinger technique) with a 10 French (minimum 3,5 mm) vascular sheath and a 0.035 inch guidewire (maximum 0.95 mm). The minimum inner diameter of the vascular sheath that may be used is 3.38 mm. During manipulation of the application system the guidewire must always remain in the guidewire lumen of the application system.

Implantation of blueflow Venous Stent:

Before stent release, the deployment lock has to be inactivated by opening the Tuohy-Borst Valve (Figure 1, No. 90). The lock prevents unintended or premature stent deployment during device handling. Do not unlock the valve before reaching the implantation site!
To release the stent, a characteristic two-handed-maneuver is performed by the operator, described here for a right-handed person:

1. The correct stent position and stent size – extending beyond the distal and proximal end of the occlusion – is fluoroscopically confirmed. The radiopaque catheter tip and the radiopaque marker band visualize the distal and proximal end of the loaded stent.
2. Deployment is a two-handed maneuver. The left hand on the distal grip pulls toward a fixed right hand, which holds the proximal grip still. It is important that the right hand stays completely fixed; otherwise it would change the position of the stent within the occlusion. This movement releases the loaded stent out of the blue outer shaft of the application system.
3. Before the full release of the stent, the physician checks the correct placement of the stent. If needed, the stent can still be pulled back into the system, as long as the marker on the metal pusher tube (Figure 1, No. 100) is still visible.
4. To repositioning or remove the stent during the procedure, gently advance the distal handle while holding the proximal handle in place. The application system will retrieve the stent back into the stent pocket.
Do not retrieve the stent more than three times!
5. Release the stent fully by pulling the distal grip fully back, holding the right hand completely fixed.

Removing delivery system:

After the stent is fully released and expanded in the venous lumen, the application device is closed and retracted from the vessel over the guidewire, which is removed last. Following stent deployment, dilatation of the stent via balloon angioplasty might occasionally be required for full stent expansion. The balloon chosen for this additional expansion should have a nominal diameter not larger than the diameter of the stent.

Stent removal or replacement:

The blueflow Venous Stent System is designed for implantation and cannot be removed or replaced after implantation procedure is finished.

Registration:

A patient card is inside the package. Stent details, indication and implantation areas are to be entered. The removable labels on the package can also be used for this purpose.

INSTALLATION, ASSEMBLY AND MAINTENANCE

Installation and assembly is not required for the use of this medical device.
The product does not require maintenance and should not be repaired.

VISUALIZATION AND DEIMPLANTATION

The installation of the stent should be carried out under X-ray control.

Previously implanted stents (presence and position) are determined by standard procedures using X-ray examination.

The blueflow Venous Stent System is designed for implantation and cannot be removed or replaced after implantation procedure is finished.

STERILITY

The blueflow stent on the delivery system is delivered sterile in a sealed package, sterilization is carried out using ethylene oxide gas (EO) used as a sterilizing agent.

Do not use if the tightness of the package is broken.

For one-time use only. Repeated sterilization and reuse is prohibited.

SHELF LIFE

The shelf life (guaranteed period of preservation of sterility) is 13 months.

The exact expiration date is indicated on the package.

This medical device cannot be used after the expiration date.

CLEANING AND DISINFECTION

The product is not subject to cleaning and/or disinfection.

This is a single-use product; it cannot be reused, because after the first application it loses the necessary functional properties. Changes in mechanical, physical and (or) chemical properties due to repeated use, cleaning, disinfection and (or) repeated sterilization may lead to a violation of the integrity of the product and (or) destruction of the materials from which it is made, which, in turn, contributes to contamination (infection) of the product.

CONDITIONS OF TRANSPORTATION, STORAGE and OPERATION

The products can be transported by all modes of transport in covered vehicles in accordance with the rules of cargo transportation applicable to each type of transport.

A medical device should be handled with caution. Store in the original packaging in a dry place, protect from direct sunlight.

During transportation, storage and operation, it is necessary to adhere to the following modes:

Table 8. Transportation, storage and operation conditions

Condition	During transportation	During storage	During operation
Temperature	from -20 to +55°C	from +20 to +25°C	from +32 to +42°C
Humidity, special conditions	from 30 to 75%	from 30 to 65%	the product is resistant to biological fluids (blood)
Atmospheric pressure	101 ±10 kPa	101 ±10 kPa	101 ±10 kPa

SCOPE OF DELIVERY

The product delivery package includes (for all versions, the delivery package is the same):

- 1) blueflow self-expanding venous stent on the delivery system - 1 pc.;
- 2) Patient card - 1 pc.;
- 3) Instruction for use - 1 pc.

PACKAGING

One delivery system, with a loaded self-expanding blueflow stent, is fixed on a substrate, which is packed in a Tyvek package (individual packaging).




The Tyvek package is placed in a cardboard box (consumer packaging).

Each package and packaging are labeled. Each package is provided with instructions for use and patient cards.

5 (five) consumer packages are placed in a corrugated cardboard transport box for transportation and labeled accordingly.

LABELING

Table 9. Symbols on the label

Symbol	Description
	Catalog number
	Batch code
	Used-by date
	Sterilized Using Ethylene Oxide
	Keep away from direct sunlight
	Keep dry
	CE Mark
	Caution: Please read instructions for use
	Do not use if package is damaged
	Do not re-sterilize
	Do not re-use
	Manufacturer

	Magnetic Resonance Conditional (implant only)
	Non-pyrogenic

WARRANTY OBLIGATIONS

The manufacturer guarantees the compliance of this medical device with the declared characteristics within the limits of preserving its sterility period, provided that the user complies with the requirements specified in the manufacturer's operational documentation.

DISPOSAL OR DESTRUCTION

Products that have come into contact with the blood and (or) other biological fluids of patients should be disposed of in accordance with the requirements for the treatment of medical waste established by SanPiN 2.1.3684-21 for Class B waste (epidemiologically hazardous waste), and other regulations of the Russian Federation.

Unused products, including expired products that have not come into contact with the biological fluids of patients, should be disposed of in accordance with the requirements for the treatment of medical waste established by SanPiN 2.1.3684-21 for Class A waste (epidemiologically safe waste, composition close to MSW).

ENVIRONMENTAL PROTECTION REQUIREMENTS

Taking into account the fulfillment of the requirements of the operational documentation, the product does not have a negative impact on humans and the environment during use, transportation and storage.

LIST OF INTERNATIONAL NORMATIVE DOCUMENTS/STANDARDS THAT MEDICAL DEVICE COMPLIES WITH

Standard	Name
Council Directive 93/42/EEC of June 14, 1993 concerning medical devices	Council Directive 93/42/EEC of June 14, 1993 concerning medical devices
Directive 2007/47/EC of 5 September 2007, amending Council Directive 93/42/EEC concerning medical devices	Directive 2007/47/EC of 5 September 2007, amending Council Directive 93/42/EEC concerning medical devices
EN 1707:1996	Conical fittings with 6% (Luer) taper for syringes, needles and certain other medical equipment
EN ISO 10555-1	Intravascular catheters – Sterile and single-use catheters – Part 1: General requirements
ISO 14630	Inactive surgical implants - General requirements catheters
EN ISO 10993-1	Biological evaluation of medical devices -- Part 1: Evaluation and testing within a risk management process
EN ISO 10993-3	Tests for Genotoxicity, Carcinogenicity, and Reproductive Toxicity
EN ISO 10993-4	Biological evaluation of medical devices – Part 4:

EN ISO 10993-5	Selection of tests for interactions with blood
ISO 10993-6	Biological evaluation of medical devices – Part 5: Tests for invitro cytotoxicity
EN ISO 10993-7	Biological evaluation of medical devices – Part 6: Tests for local effects after implantation
ISO 10993-10	Biological evaluation of medical devices – Ethylene oxide sterilization residuals
ISO 10993-11	Biological evaluation of medical devices – Part 10: tests for irritation and skin sensitization
EN ISO 10993-18	Biological evaluation of medical devices – Part 11: Tests for system toxicity
ISO 11135-1	Biological evaluation of medical devices – Part 18: Chemical characterization of materials
ISO 11138-2:2017	Sterilization of health care products – Ethylene oxide - Part 1: Requirements for development, validation and routine control of a sterilization process for medical devices
EN ISO 11607-1	Sterilization of health care products – Biological indicators – Part 2: Biological indicators for ethylene oxide sterilization processes
EN ISO 11607-2	Packaging for terminally sterilized medical devices – Part 1: Requirements for materials, sterile barrier systems and packaging systems
EN ISO 11737-1	Packaging for terminally sterilized medical devices – Part 2: Validation requirements for forming, sealing and assembly processes
ISO 11737-2	Sterilization of medical devices – Microbiological methods – Part 1: Determination of a population of microorganisms on products
EN ISO 13485	Sterilization of medical devices – Microbiological methods – Part 2: Tests of sterility performed in the definition, validation and maintenance of a sterilization process
EN ISO 14155	Medical devices – Quality Management Systems – Requirements for regulatory purposes
EN ISO 14644-1	Clinical investigations of medical devices for human subjects – Good clinical practice
EN ISO 14971	Cleanrooms and associated controlled environments – Part 1: Classification of air cleanliness by particle concentration
EN ISO 15223-1	Medical devices – Application of risk management to medical devices
DIN EN ISO 25539-2	Symbols for use in the labelling of medical devices
	Cardiovascular implants – Endovascular devices – Part 2: Vascular stents

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