APPROVAL

Ethicon, LLC
Organization

Highway 183 Km 8.3, San Lorenzo, 00754, Puerto-Rico, USA
Address

Associate Director – Regulatory Affairs

Occupation

Debbie Fazen

Surname, Name

08.02.2023


Date DD.MM.YYYY

Debbie Fazen

Signature

Subscribed and sworn to me before me
this February 8th, 2023

Keith Petersen
Notary Public

Instruction for use for the medical device

Surgicel Powder, sterile resorbable hemostatic based on oxidized regenerated cellulose, in the applicator

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tissue reaction. Absorption depends upon several factors, including the amount used, degree of saturation with blood, and the tissue bed. In addition to its local haemostatic properties, SURGICEL Powder is bactericidal *in vitro* against a wide range of gram-positive and gram-negative organisms including aerobes and anaerobes. SURGICEL Powder is bactericidal *in vitro* against strains of species including those of:

methicillin-resistant <i>Staphylococcus aureus</i> (MRSA)	<i>Shigella sonnei</i>
penicillin-resistant <i>Streptococcus pneumoniae</i> (PRSP)	<i>Serratia marcescens</i>
vancomycin-resistant <i>Enterococcus</i> (VRE)	<i>Bacillus subtilis</i>
methicillin-resistant <i>Staphylococcus epidermidis</i> (MRSE)	<i>Proteus vulgaris</i>
<i>Staphylococcus aureus</i>	<i>Corynebacterium xerosis</i>
<i>Staphylococcus epidermidis</i>	<i>Mycobacterium phlei</i>
<i>Micrococcus luteus</i>	<i>Clostridium tetani</i>
<i>Streptococcus Group A</i>	<i>Clostridium perfringens</i>
<i>Streptococcus Group B</i>	<i>Bacteroides fragilis</i>
<i>Streptococcus salivarius</i>	<i>Enterococcus sp.</i>
<i>Branhamella catarrhalis</i>	<i>Enterobacter cloacae</i>
<i>Escherichia coli</i>	<i>Pseudomonas aeruginosa</i>
<i>Klebsiella pneumoniae</i>	<i>Pseudomonas stutzeri</i>
<i>Lactobacillus sp.</i>	<i>Proteus mirabilis</i>
<i>Salmonella enteritidis</i>	

SURGICEL Powder is not, however, an alternative to systemically applied therapeutic or prophylactic antimicrobial agents.

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INDICATIONS AND CONTRAINDICATIONS FOR MD USE AND POSSIBLE ADVERSE EFFECTS

Indications

SURGICEL Powder is used adjunctively in surgical procedures to assist in the control of capillary, venous, and small arterial haemorrhage when ligation or other conventional methods of control are impractical or ineffective.

SURGICEL Powder can also be applied in laparoscopic or other endoscopic procedures when used with the SURGICEL Endoscopic Applicator.

Contraindications

- Do not inject or place SURGICEL Powder into an open blood vessel. Do not use to treat bleeding from large defects in arteries or veins.
- SURGICEL Powder should not be used to control haemorrhage from large arteries or veins.
- The SURGICEL Powder device was not designed for intraluminal procedures.
- When SURGICEL Powder is used to help achieve haemostasis in, around, or in proximity to foramina in bone, areas of bony confine, the spinal cord, or the optic nerve and chiasm, it must always be removed after haemostasis is achieved since it will swell and could exert unwanted pressure.
- Unlike other SURGICEL products, SURGICEL Powder particles are integrated with blood clots and complete removal of the device (hemostatic agent) may disrupt the clot and increase the risk of re-bleeding. In cases where complete removal is required, SURGICEL fabrics should be applied.
- SURGICEL Powder should not be used for implantation in bone defects, such as fractures, since there is a possibility of interference with callus formation and a theoretical chance of cyst formation.
- SURGICEL Powder should not be used on non-haemorrhagic serous oozing surfaces, since body fluids other than whole blood, such as serum, do not react with SURGICEL Powder to produce satisfactory haemostatic effect.

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• SURGICEL Powder is an absorbable haemostat, and should not be used as an adhesion prevention product.

Warnings

• Do not pump device intravascularly. Life-threatening embolic events may occur if the product is applied intravascularly.

• SURGICEL Powder is not intended for use on dry (non-bleeding) surfaces or for prevention of bleeding.

• SURGICEL Powder is supplied sterile and as the material is not compatible with autoclaving or ethylene oxide sterilization, SURGICEL Powder should not be resterilized.

• SURGICEL Powder is not intended as a substitute for careful surgery and the proper use of sutures and ligatures.

• Closing with SURGICEL Powder in a contaminated wound without drainage may lead to complications and should be avoided.

• The haemostatic effect of SURGICEL Powder is greater when it is applied dry; therefore, it should not be moistened with water or saline prior to application.

• SURGICEL Powder should not be impregnated with anti-infective agents or with other materials such as buffering or haemostatic substances. Its haemostatic effect is not enhanced by the addition of thrombin, the activity of which is destroyed by the low pH of the product.

• Although SURGICEL Powder may be left *in situ* when necessary, it is advisable to remove excess powder with irrigation and aspiration once haemostasis is achieved, without disturbing the clot.

• SURGICEL Powder is dry and there may be difficulties in precise delivery under certain circumstances. Unintentional device placement may result in powder scattering and device migration that may increase the risk of adhesion formation. In preclinical *in vivo* animal studies it was demonstrated that SURGICEL Powder does not increase the incidence of remote adhesions in laparoscopic procedures.

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- Dislodgement of SURGICEL Powder could possibly occur by intraoperative manipulation, lavage, exaggerated respiration, etc. With other SURGICEL products there have been reports that in procedures such as lobectomy, laminectomy, and repair of a frontal skull fracture and lacerated lobe, when the product was left in the patient after closure it migrated from the site of application into foramina in bone around the spinal cord, resulting in paralysis and, in one case, the product migrated into the left orbit of the eye, causing blindness. While these reports cannot be confirmed to be related to SURGICEL products, special care must be taken by physicians, regardless of the type of surgical procedure. Consider removing SURGICEL Powder in these applications (procedures) after haemostasis is achieved.

- Although SURGICEL Powder is bactericidal against a wide range of pathogenic microorganisms, it is not intended as a substitute for systemically administered therapeutic or prophylactic antimicrobial agents to control or to prevent postoperative infections.

- Do not resterilize/reuse. Reuse of this device (or portions of this device) may create a risk of product degradation, which may result in device failure and/or cross-contamination, which may lead to infection or transmission of blood-borne pathogens to patients and users.

- Do not attempt to trim the applicator tip.


Precautions

- SURGICEL Powder should not be used in conjunction with autologous blood salvage circuits, because its fragments may pass through the transfusion filters of blood-scavenging systems.

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
- Use only as much SURGICEL Powder (oxidised regenerated cellulose) as is necessary and apply only where needed for haemostasis. Remove any excess before surgical closure in order to facilitate absorption and to minimize the possibility of foreign body reaction.
- Use minimal amount of SURGICEL Powder required to achieve haemostasis, and remove excess powder in the area of drains to prevent clogging. In urological procedures, minimal amounts of SURGICEL Powder should be used and care must be exercised to prevent plugging of the urethra, ureter, or a catheter by dislodged portions of the product.
- Since absorption of SURGICEL Powder could be prevented in chemically cauterized areas, its use should not be preceded by application of silver nitrate or any other escharotic chemicals.
- If SURGICEL Powder is used temporarily to line the cavity of open wounds, it should be removed by irrigation with sterile water or saline solution after bleeding has stopped.
- Precautions should be taken in otorhinolaryngologic surgery to ensure that none of the material is aspirated by the patient (e.g., controlling haemorrhage after tonsillectomy and controlling epistaxis).
- The applicator tip provided on the SURGICEL Powder device is not intended for laparoscopic or other endoscopic use. If laparoscopic or other endoscopic use is desired, remove the existing applicator tip from the SURGICEL Powder device, and replace with the SURGICEL Endoscopic Applicator tip (supplied separately). In laparoscopic or other endoscopic procedures, SURGICEL Powder should only be applied using the SURGICEL™ Endoscopic Applicator. Consult the SURGICEL Endoscopic Applicator Instructions for Use (IFU) for proper assembly and directions for use with the SURGICEL Powder device.

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




2. By holding the body of the applicator, remove SURGICEL Powder delivery device from inner card (Fig. 2).

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
4. SURGICEL Powder delivery device is now ready for use. Compress bellows using light pressure to apply powder to treatment site (Fig. 5).

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
5. As an approximate guide, one SURGICEL Powder delivery device can provide coverage of 213 cm² (14.6 cm x 14.6 cm area) when applied in one layer at a distance of approximately 5 cm from the bleeding surface. Depending on bleeding intensity and anatomical location, more than one layer may be needed to achieve complete haemostasis at the targeted bleeding surface.

6. If necessary, powder may be held firmly against the tissues until haemostasis is obtained.

7. To prevent clogging, do not touch the tip to wet surface. Do not disassemble device bellows. Do not trim the applicator tip (see WARNINGS).

8. In the event of clogging, the tip can be wiped off using dry sterile surgical gauze to remove clog. If clog cannot be removed by dry gauze, dispose of the clogged device and use a new SURGICEL Powder delivery device.

Sterile technique should be observed in removing the SURGICEL Endoscopic Applicator from its sterile packaging. Apply adequate amount of SURGICEL Powder to cover bleeding area. Use of a non-adhering substrate to apply pressure may prevent adhesion of the formed clot to the surgical glove or other instrumentation. Unused SURGICEL Powder should be discarded.

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The maximum single dose is 2 packs of powder, each pack contains 3 g.

Sterile technique should be observed in removing the applicator for SURGICEL™ Powder (oxidised regenerated cellulose) from its sterile packaging. Apply adequate amount of SURGICEL Powder to cover bleeding area. Use of a non-adhering substrate to apply pressure may prevent adhesion of the formed clot to the surgical glove or other instrumentation. Unused SURGICEL Powder should be discarded. The maximum single dose is 2 packs of powder, each pack contains 3 g.

INFORMATION ON POTENTIAL USERS OF MD AND DIRECTIONS FOR USE

Directions for use

Devices are intended for use in healthcare facilities.

Devices are recommended for use under operation conditions indicated below.

Use/operation conditions of MD

The device is used in the internal environment of the body.

Operation conditions

Humidity	100%
Temperature	от +32°C до +42°C
Pressure	No special conditions

Potential users of MD


SURGICEL Powder is intended for use during surgical procedures, by health care professionals.

DESCRIPTION OF MAIN FUNCTIONAL ELEMENTS OF MD AND PRINCIPLE OF OPERATION

General description of functional elements of MD

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The SURGICEL Powder form of the product is made of the same Oxidized Regenerated Cellulose (ORC) as the rest of the SURGICEL Family. General view of the medical device is shown in Figure 6.

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
The powder is prefilled in an applicator (1) to dispense on a target bleeding site. The SURGICEL Powder applicator is provided with 3 grams of powder, which is equivalent in weight to approximately two 4 x 8 inch pieces of SURGICEL Original.

The applicator enables surgeons to topically apply powder to a desired anatomical site by manually compressing the bellows (3). The tip is flexible such that the surgeon can angle the tip in the direction of bleeding. The applicator is provided in the “closed” or “off” position so powder cannot be expressed (see Figure 6).

To use, the Trap is twisted to the “open” or “on” position by rotating 90 degrees (see Figure 7). This rotation exposes an output orifice to the powder.

Trap rotated to the on or open position

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
The applicator is pumped by depressing the bellows, generating air flow through the powder and expressing the powder through the now exposed output orifice. The bellows are repeatedly pumped to express the desired amount of powder to the target anatomy.

The applicator is designed so that the tip can be removed in order to attach an alternative tip (not included with SURGICEL Powder device) if required. See Figure 8.

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The possibility and ways of integrating with other medical devices

SURGICEL Powder can be applied through a trocar for endoscopic or laparoscopic use when the SURGICEL Powder Device is used with the SURGICEL Endoscopic Applicator.

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*for the users convenience, the labeling and other materials may indicate the rounded value of the mass of powder in the applicator - 3 g.

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Table 2 – MD materials

MD part	Material of MD
1. Applicator (Reservoir)	Polypropylene
2. Open Flex Tube	Polyurethane
3. Plunger	Polypropylene
4. Filter	Polyolefin Fibers
5. Hub	Polycarbonate, colorant
6. Trap	Polypropylene, colorant
7. Bellows	Ethylene Vinyl Acetate, Colorant, Slip Agent
8. Powder	Oxidized Regenerated Cellulose

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DL-JNJRUMedicalQA@ITS.JNJ.com

Highway 183 Km 8.3, San Lorenzo, 00754, Puerto-Rico, USA

**AUTHORIZED REPRESENTATIVE OF THE MANUFACTURER
IN RUSSIA:**

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