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CERTIFICATE

中国国际贸易促进委员会暨中国国际商会
China Council for the Promotion of International Trade is China Chamber of International Commerce

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Translation from English into Russian

From ZHEJIANG BIOREAR BIOTECHNOLOGY CO., LTD. ("ZHEJIANG BIOREAR BIOTECHNOLOGY CO., LTD."), I certify that the document provided is accurate, correct and valid.

Approved: October 28, 2022

Signature/Name:  (Bian Li)

Position: Vice President

INSTRUCTIONS FOR USE OF THE MEDICAL DEVICE

Tip comb for IDW magnets, clear, polypropylene, non-sterile, compatible
with automated nucleic acid extraction systems, in the assortment

Version 28102022 dated October 28, 2022

1. Name and composition of the medical device:

Tip comb for DW magnets, clear, polypropylene, non-sterile, compatible with automated nucleic acid extraction systems, including

- comb of tips - 5 packs of 2 pieces per pack, 10 packs per box;
- instructions for use - 1 pc.

2. Indications and contraindications for the use of a medical device.

Indications for use

"The tip comb for DW magnets, colorless, polypropylene, non-sterile, compatible with automated nucleic acid extraction systems, in the assortment" is to be used once when performing automatic extraction of nucleic acids from any type of human biological material, regardless of the analyte being determined.

Contraindications for use, possible side effects

There are no contraindications for use, possible side effects within the established purpose.

3. Information about potential consumers of a medical product.

Potential consumers are the personnel of clinical diagnostic laboratories (specialists in the field of clinical laboratory diagnostics). For professional use only.

4. Specifications. Description of the main functional elements of the comb of tips.

Production material	polypropylene
Color	colorless
Number of tips, pcs.	96
non-sterile	

The medical product is biologically inert. Free from DNA, RNases, DNases, PCR inhibitors.

The dimensions of the comb of the tips are shown in Figure 1.


Figure 1 - Dimensions of the comb of tips, mm

5. Purpose of the medical device, principles of operation and procedure for the use of medical devices

Intended use

"Comb of tips for DW magnets, colorless, polypropylene, non-sterile, compatible with automated systems for the extraction of nucleic acids, as a part" is intended for placement on the magnetic head of compatible automated systems for purification of nucleic acids and retention of magnetic particles in the process of simultaneous multichannel processing of biomaterial (magnetic separation) when using deep well plates. The product is intended for single use by a professional user when performing diagnostic procedures in vitro.

Functional use

Tip comb for DW magnets, colorless, polypropylene, non-sterile, compatible with automated nucleic acid extraction systems, as a part, used in clinical diagnostic laboratories as a consumable (auxiliary) for diagnostics in vitro.

Operating principle

The success of the subsequent analysis depends on the qualitative and reproducible purification of nucleic acids, proteins and cells. Nucleic acid extraction systems provide quality results with minimal manual preparation time, helping to automate much of your workflow.

Tip combs made of polypropylene are resistant to biomolecules and are ideal for working with magnetic particles. When using tip combs specially designed for nucleic acid extraction systems, the quality of protein and nucleic acid extraction and the yield of the target product are significantly improved (see Figure 2).

In nucleic acid extraction systems, magnetic rods with disposable tips collect, transfer and mix magnetic particles:

- 1) When the magnetic rod located inside the tip is lowered into the solution, the magnetic balls are collected at the end of the tip.
- 2) When the tip is placed in the wells of the plate, particles are released as the magnetic bars move out of the tip.
- 3) The tip makes it easier to mix reagents with particles as the magnetic head moves up and down.



Figure 2 - Principle of operation

The order of application of tip comb.

Refer to the instructions for use of the product for compatible nucleic acid extraction systems.

Before installation in the system holders, inspect the product for any damage.

6. Operational / functional characteristics and methods of control.

6.1 Compatibility of tip combs with corresponding nucleic acid extraction system holders :

The comb of tips should be symmetrically captured from both sides by the holder of the nucleic acid extraction system and fixed in it during the entire cycle of nucleic acid extraction carried out in accordance with the specified extraction protocol.

Compliance is assessed visually in 8-10 repetitions when running on a compatible nucleic acid isolation system.

6.2 Mutual geometric compatibility of tip combs and deep well plates:

The comb tips should be free to dip into the appropriate wells of the plate. The comb of tips placed in the tablet should be located in it without a visible slope and displacement. The comb tips should be free to remove from the wells of the plate.

Compliance is assessed visually by manually aligning plates and combs of tips by testing 8-10 pieces of each product.

The tip combs are geometrically compatible with deep well plates having the dimensions shown in Figure 3:

Figure 3 - Dimensional characteristics of compatible deep well plates

6.3 Geometric compatibility of tip combs with magnetic heads of the nucleic acid extraction system :

The rods of the magnetic head of the system for extracting nucleic acids must freely, to the depth specified in the operational documentation of the system, be immersed in the combs of the tips.

Compliance is assessed visually in 8-10 repetitions when running on a compatible nucleic acid isolation system.

6.4 Retention of magnetic particles during operation in a nucleic acid extraction system :

The underside of all comb tips, when removed from the wells of the plate, which contained the solution with magnetic particles, should be visible to the naked eye evenly distributed coating, having the color of the magnetic particles included in the compatible nucleic acid extraction system used .

In this case, the liquid that remains in the wells of the plate should have the color and transparency of the original (before the introduction of magnetic particles) solution. A slight (hardly visible to the eye) sediment is allowed in single (1-10) wells of the tablet *.

Compliance is assessed visually by testing 10-20 samples in 3 repetitions. The presence of sediment in 1-2 wells is allowed.

Note:

** If this indicator does not meet the specified requirements, it is necessary to exclude the malfunction of the functional modules of the magnetic head of the nucleic acid extraction system.*

6.5 Preservation of the integrity of the walls of the combs of the tips during the extraction of nucleic acids:

At the end of the extraction cycle, the tips should be free of visible traces of liquid.

Compliance is assessed visually in 3 repetitions.

7. Marking.

On the labeling of packages (primary and transport), symbols are used, the decoding of which is presented in table 1.

Table 1

Symbol name	Symbol
Batch code	
Use up to	
Manufacturer (name and address of manufacturer)	
Date of manufacture	
Medical device for diagnostics in vitro	
Reuse prohibition	

Symbol name	Symbol
Do not use if packaging is broken	
Protect from moisture / Keep dry	
Avoid exposure to sunlight / Keep away from heat sources	
Temperature limit	

8. Expiry date.

The manufacturer guarantees the conformity of products to the requirements of technical documentation. Warranty period - 36 months from the date of manufacture.

9. Transportation and storage.

9.1 The product is stored in a transport package in closed, clean and dry warehouses protected from moisture and direct sunlight, at a distance of at least 1 m from heating devices on racks at least 5 cm from the floor in conditions that exclude the impact of aggressive environments (acidic, alkaline, etc.), as well as flammable and combustible liquids at a temperature of 5 to 30 ° C in a dry place.

9.2 The product is transported by all means of transport in covered vehicles in accordance with the rules for the carriage of goods in force for this type of transport.

9.3 During storage and transportation, it is not recommended to stack transport packages on top of each other.

9.4 It is not allowed to transport the product without transport packaging.

9.5 Environmental conditions during transportation:

- air temperature from minus 40°C to 60°C,
- relative air humidity ≤ 95 %;
- atmospheric pressure 84.0-106.7 kPa.

10. Disposal.

Used products are classified as epidemiologically hazardous waste. Unused products are epidemiologically safe waste. Disposal must be carried out in

accordance with the rules in force in the country of the consumer, observing all established restrictions at the time of disposal.

In the Russian Federation, the product is disposed of in accordance with the requirements of SanPiN 2.1.3684-21 "Sanitary and epidemiological requirements for the maintenance of territories of urban and rural settlements, for water bodies, drinking water and drinking water supply, atmospheric air, soils, residential premises, operation of industrial, public premises, organization and implementation of sanitary and anti-epidemic (preventive) measures. Waste should be disposed of through organizations licensed for this type of activity.

11. Instructions for use.

The use of the product must be carried out in accordance with the instructions for use of the product itself and compatible systems for the extraction of nucleic acids, for example, the KingFisher Magnetic Particle Processor for Purification of Nucleic Acids, Cells and Proteins (FSZ 2009/05562).

Permissible operating conditions:

Air temperature from 5 to 30°C.

The order of transportation, storage and use, set out in the instructions for use, must be observed.

12. Requirements for the maintenance and repair of a medical device.

Medical polymer consumable products compatible with automated nucleic acid extraction systems are single-use products that are not subject to maintenance.

13. Applicable standards.

Standard	Name of the standard
ISO 13485	Medical products. Quality management systems. Requirements for Regulatory Purposes
ISO 14971	Medical products. Application of risk management to medical devices
ISO 15223-1	Medical products. Symbols used with information provided by the manufacturer. Part 1. General requirements
EN ISO 23640	Medical products for in vitro diagnostics. Evaluation of the stability of reagents for in vitro diagnostics
EC Directive 98/79/EC	Directive 98/79/EC of the European Parliament and of the Council of 27 October 1998 on medical devices for in -diagnosis vitro

14. Medical device manufacturer:

ZHEJIANG BIOREAR BIOTECHNOLOGY CO., LTD.
("ZHEJIANG BIOREAR BIOTECHNOLOGY CO., LTD.")

Address of the manufacturer of the medical device:

1st Floor , Building 2, No.3535 , Linggongtang Road , Nanhu district , Jiaxing ,
Zhejiang , China / 1 Floor, Building 2, No .3535, Linggongtan Road, Nanhu
District, Jiaxing, Zhejiang, China

Phone: + 86 057382618891

Email: sales@biorear.com

15. Authorized representative of the manufacturer:

Cosmofarm Limited Liability Company (Cosmofarm LLC)

Address of the authorized representative of the manufacturer of the medical device:

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