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
Ashley Luloff

Ashley Luloff

Associate Regulatory Affairs Specialist

Medtronic Inc

«08» July 2014 г.

Spinal needle for intrathecal catheter insertion Ascenda

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Spinal needle for intrathecal catheter insertion Ascenda (8786, 8787)

Component	Material	Material contacts Human tissue	Material contacts drug
Introducer needle	Stainless steel 304	Yes	No
Needle hub	Plastic Acrylic / polycarbonate	Yes	No
Stylet	Stainless steel 304	Yes	No

Instructions for use

Implanting physicians should be experienced in intrathecal-access procedures and should be thoroughly familiar with all product labeling.

Cautions:

To prevent damage to surrounding tissue and injury to the patient, exercise care when using electrocautery near the introducer needle.

Positioning the patient

1. Place the patient in a lateral recumbent position with the lumbar region slightly flexed. Drape to allow fluoroscopic visualization of the

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spine in regions where the catheter will be inserted and where the catheter tip will be placed. Position the patient's hips, legs, and arms to avoid pressure points and interference with fluoroscopic imaging of the spine.

2. Administer the appropriate anesthesia.

Inserting the introducer needle

Select the introducer needle of the desired length for the patient's needs.

Note: A selection of needles is available from Medtronic needle kits.

2. Make certain that the key (Figure 2) of the needle stylet cap is properly seated in the slot of the needle hub.

3. Orient the bevel of the 16 T-gauge introducer needle cephalad or parallel

to the dural fibers and insert while using fluoroscopy. Use a paramedian oblique approach to insert the needle into the intrathecal space at a shallow angle (approximately equal to 30°) (Figure 3). The entry point of

the needle into the skin (or fascia, if needle insertion is performed through

an open incision) should be parallel to the vertebral pedicle, approximately

1 to 1.5 vertebral levels below the interlaminar space selected for dural puncture and 1 to 2 cm lateral to the midline on the side of the intended pump pocket (Figures 2 and 3).

Note: In retrospective analysis, low catheter complication rates were associated with a shallow paramedian oblique technique.

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
Figure 2. Shallow paramedian oblique insertion top view.

Note: The orientation indicator is the dot, arrow, color, or other indicator on the needle hub that depicts which direction the bevel of the needle is positioned.


Figure 3. Shallow paramedian oblique insertion side view.

4. Advance the needle until the dura is penetrated.
5. Remove the needle stylet and confirm needle location by observing cerebrospinal fluid (CSF) backflow.
6. Replace the needle stylet to stop CSF flow. Make certain the key of

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the needle stylet cap is properly seated in the slot in the needle hub (Figure 2).

Placing the catheter

1. Verify that the bevel of the needle has a cephalad orientation by making sure that the key on the stylet cap and the orientation indicator (dot, arrow, color, or other indicator) on the needle hub face up (Figure 2).
2. Remove the stylet.
3. Confirm that the catheter is seated within the channel of the guide wire handle (Figure 4).


Figure 4. Catheter seated within the guide wire handle.

4. Thread the spinal segment catheter tip (the end with fluid-outlet holes)

through the needle and advance the catheter to the desired location. **Note:** When the catheter tip reaches the curved point of the needle, a slight

increase in advancement pressure will be noted.

#Caution: If the catheter must be retracted during positioning, do not withdraw the catheter through the introducer needle. The needle tip can damage the catheter, requiring additional surgery to repair or replace the catheter. Withdraw the introducer needle and catheter together. Then, carefully retract the catheter from the introducer needle. Begin the procedure again.

5. Verify the catheter position using fluoroscopy.

6. Prepare the anchor site by making an incision at the needle site to expose

an area of the fascia that is large enough to place an anchor. The incision

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should be long enough to avoid introducing kinks or sharp bends when anchoring the catheter.

Note: To protect the catheter, keep the introducer needle in place while the incision is made. Alternatively, this incision can be made prior to catheter placement.

7. Dissect the subcutaneous tissue to develop a smooth fascial plane upon which to suture the anchor beneath the edges of the incision.

8. Using heavy nonabsorbable suture thread, place two tie-down sutures adjacent to, and on opposite sides of, the needle entry point. Leave 10-12 cm of suture thread exposed.

9. Place a purse-string suture on the fascia around the catheter but do not close the fascia at this time.

Caution: To prevent catheter damage or dislodgement during introducer needle and guide wire removal:

- hold the catheter straight.
- hold the catheter securely at the exit site (do not squeeze the catheter because excessive compression will increase the difficulty of guide wire removal).
- slowly remove the needle and guide wire.

10. Carefully remove the needle from the fascia and grasp the catheter near the exit site to prevent catheter dislodgement.

11. Being careful not to damage the catheter with the needle tip, slide the needle to the guide wire handle until the needle hub touches the guide wire handle.

12. Tilt the guide wire handle to bring the end of the catheter fully away from the channel of the guide wire handle.

13. Hold the catheter securely at the fascia exit site and remove the needle and guide wire from the catheter simultaneously (Figure 5).

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
Figure 5. Removing the needle and guide wire.

Products are supplied in a package labeled with the following information:

- name and number of products in the catalogue;
- lot number;
- trade mark and address of the manufacturer;
- single use;
- ethylene oxide sterilization, sterile;
- the expiry date of the product;
- the certificate CE number;
- date of production;

Refer to the appropriate product for symbols that apply.

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Caution, consult accompanying documents

Use by

Manufacturer

Date of manufacture

Temperature limitation

Lot number

Authorized representative in the European community

Conformité Européenne (European Conformity).
This symbol means that the device fully complies with
European Directive AIMD 90/385/EEC.

General Information:

Sterile, EO

Single-Use, Do Not Reuse

Caution: Read Instructions For Use (IFU) before using this device.

Conditions of transportation and storage :

Storage temperature -30°F to 135°F -34°C to 57°C

Temperature transport -30°F to 135°F -34°C to 57°C

The device is transported by all modes of transport in covered vehicles in conformity with the rules of transportation of goods, operating in each mode.

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Disposal:

Dispose of used devices as a biohazard.

Shelf life:

The shelf life of the product is 2 years.

Warranty Limited Warranty-Areas outside the United States should contact their local Medtronic representative for exact terms of the Limited Warranty.

Customer Service

If you have any device returns or questions about Medtronic's device please call US Medtronic Customer Service at 1-800-328-0810 or Medtronic Russia Tel.: (+7495) 580-7377.

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