20 1

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CERTIFICATE

中国国际贸易促进委员会暨中国国际商会
China Council for the Promotion of International Trade is China Chamber of International Commerce

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China Council for the Promotion of International Trade
China Chamber of International Commerce

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THIS IS TO CERTIFY THAT: the seal of YAFHO BIO-TECHNOLOGY CO., LTD. on the annexed DOCUMENT is genuine.

China Council for the Promotion
of International Trade

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Authorized

Signature:

Lu Xiaohui

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(Date: Mar. 24, 2022)

网站 Website for verifying the certificate: <http://www.rzccpit.com/validate.html>

Yafho Bio-Technology Co., Ltd

Second Floor Room 202 and Third Floor No. 81, Junfeng Road, Huangpu
District Guangzhou 510760 Guangdong China
Tel: 020-85295649 Fax: 020-85295797 Website: www.yafho.com

APPROVE

GENERAL MANAGER
Yafho Bio-Technology Co.,Ltd
HANMING-MA

February 01, 2021

M.P.

INSTRUCTIONS FOR USE FOR MEDICAL DEVICE

Polymer-based dressing

OMEGASOFT PU (OMEGASOFT PU)

2021

TABLE OF CONTENTS

1	Description of the medical device	2
2	Technical description of the medical device	4
3	Completeness	6
4	Medical device packaging	7
5	Information about the sterilization of the medical device	8
6	Instructions for use	8
7	Storage and transport conditions	10
8	Warranty	10
9	List of foreign standards used by the manufacturer (manufacturer) of medical devices	10
APPENDIX A		12
	Explanation of symbols used in marking	12

1 Description of the medical device

Table 1



1.1	Name of the medical device	OMEGASOFT PU polymer-based dressing (OMEGASOFT PU), in versions (hereinafter referred to as the dressing)
1.2	Versions	1. OMEGASOFT PU polymer-based dressing (OMEGASOFT PU), size, cm: 5x7.2, consisting of: - OMEGASOFT PU polymer-based dressing (OMEGASOFT PU), size, cm: 5x7.2 -1 pc. - Individual packaging with instructions printed on top - 1 pc. 2. OMEGASOFT PU polymer-based dressing (OMEGASOFT PU), size, cm: 5x8, consisting of: - OMEGASOFT PU polymer-based dressing (OMEGASOFT PU), size, cm: 5x8 - 1 pc. - Individual packaging with instructions printed on top - 1 pc. 3. OMEGASOFT PU polymer-based dressing (OMEGASOFT PU), size, cm: 6x7, consisting of: - OMEGASOFT PU polymer-based dressing (OMEGASOFT PU), size, cm: 6x7 - 1 pc. - Individual packaging with instructions printed on top - 1 pc. 4. "OMEGASOFT PU polymer-based dressing (OMEGASOFT PU)" size, cm: 6x10, consisting of: - OMEGASOFT PU polymer-based dressing (OMEGASOFT PU), size, cm: 6x10 -1 pc. - Individual packaging with instructions printed on top - 1 pc. 5. OMEGASOFT PU polymer-based dressing (OMEGASOFT PU), size, cm: 8x10, consisting of: - OMEGASOFT PU polymer-based dressing (OMEGASOFT PU), size, cm: 8x10 -1 pc. - Individual packaging with instructions printed on top - 1 pc. 6. OMEGASOFT PU polymer-based dressing (OMEGASOFT PU), size, cm: 9x10, consisting of: - OMEGASOFT PU polymer-based dressing (OMEGASOFT PU), size, cm: 9x10 - 1 pc. - Individual packaging with instructions printed on top - 1 pc. 7. OMEGASOFT PU polymer-based dressing (OMEGASOFT PU), size, cm: 9x15, consisting of: - OMEGASOFT PU polymer-based dressing (OMEGASOFT PU), size, cm: 9x15 -1 pc. - Individual packaging with instructions printed on top - 1 pc. 8. OMEGASOFT PU polymer-based dressing (OMEGASOFT PU), size, cm: 10x20, consisting of: - OMEGASOFT PU polymer-based dressing (OMEGASOFT PU), size, cm: 10x20 -1 pc. - Individual packaging with instructions printed on top - 1 pc. 9. OMEGASOFT PU polymer-based dressing (OMEGASOFT PU), size, cm: 10x25, consisting of: - OMEGASOFT PU polymer-based dressing (OMEGASOFT PU), size, cm: 10x25 -1 pc. - Individual packaging with instructions printed on top - 1 pc. 10. OMEGASOFT PU polymer-based dressing

		(OMEGASOFT PU), size, cm: 10x30, consisting of: - OMEGASOFT PU polymer-based dressing (OMEGASOFT PU), size, cm: 10x30 - 1 pc. - Individual packaging with instructions printed on top - 1 pc. 11. OMEGASOFT PU polymer-based dressing (OMEGASOFT PU), size, cm: 10x35, consisting of: - OMEGASOFT PU polymer-based dressing (OMEGASOFT PU), size, cm: 10x35 - 1 pc. - Individual packaging with instructions printed on top - 1 pc. 12. OMEGASOFT PU polymer-based dressing (OMEGASOFT PU), size, cm: 6x8, consisting of: - OMEGASOFT PU polymer-based dressing (OMEGASOFT PU), size, cm: 6x8 - 1 pc. - Individual packaging with instructions printed on top - 1 pc.
1.3	Purpose of the medical device	The dressing is intended for closing superficial wounds, for postoperative treatment of wounds, protection against secondary infection, for fixing catheters on the skin surface.
1.4	Scope of the medical device	The dressing is applied by medical personnel in clinical conditions in medical institutions, and is also applied independently at home.
1.5	Information about potential consumers of a medical device	- Patients after surgery - patients with cuts, abrasions, acute wounds, lacerations, - Patients who have undergone venipuncture
1.6	Indications for use	The use of a dressing is indicated for: - the need to apply an absorbent dressing over the suture on the postoperative wound; - the need to fix medical devices and catheters; - for cuts, burns, abrasions, acute wounds, lacerations.
1.7	Contraindications	Individual intolerance to materials
1.8	Possible side effects	If the precautionary measures, instructions for use, transportation and storage are violated, the following side effects are possible: - Cross-contamination when reusing the dressing after another person; - Allergic reactions in case of individual intolerance to the materials used in the manufacture of the dressing, or in case of violation of the conditions of transportation / storage and subsequent changes in the physicochemical properties of the dressing; - Contamination of the wound in case of repeated use or in case of violation of the conditions of transportation / storage, precautions and subsequent violation of the sterility of the dressing.
1.9	Information about restrictions on the sharing of medical devices	There are no restrictions on the sharing of medical devices


1.10	Potential risk class for the use of a medical device	2a
1.11	Medical device developer	Yafho Bio - Technology Co. , Ltd
1.12	Manufacturer and place of production of the medical device	Yafho Bio - Technology Co. , Ltd The address: Second Floor Room 202 and Third floor No.81 Junfeng Rd, Huangpu District Guangzhou 510760 Guangdong, China Contact details: Allen Lu Email the address : allen @ yafho . com Tel. 86-13416187329
1.13	Authorized representative of the manufacturer in the Russian Federation	Limited Liability Company " Omegamed " Abbreviated name of the enterprise: LLC " Omegamed " Address: 620014, Sverdlovsk region, Yekaterinburg, Lenina avenue, building 5b, office 6 Contact details: +7 (343) 372-09-72 www . omegasoft.ru

2 Technical description of the medical device

Table 2

2.1	The composition of the medical device	Version 1: Plaster - polymer-based dressing OMEGASOFT PU, dimensions, cm: 5x7.2; 5x8; 6x7; 6x10; 8x10; 9x10; 9x15; 10x20; 10x25; 10x30; 10x35  1. Polymer (polyurethane) film 2. Absorbent pad (hereinafter referred to as the pad) 3. Protective paper 4. External protective polymer (polyurethane) layer A - dressing width B - dressing length a - pad width b - pad length Version 2: Plaster - a dressing on a polymer basis OMEGASOFT PU, 6 cm x 8 cm  1. Polymer (polyurethane) film 2. Absorbent pad (hereinafter referred to as the pad) 3. Protective paper 4. External protective polymer (polyurethane) layer A - dressing width B - dressing length a - pad width b - pad length
-----	---------------------------------------	--

2.2	Overall dimensions of the medical device	Version 1: Plaster - polymer-based dressing OMEGASOFT PU, dimensions, cm: 5x7.2; 5x8; 6x7; 6x10; 8x10; 9x10; 9x15; 10x20; 10x25; 10x30; 10x35 A - dressing width B - dressing length a - pad width b - pad length Dimensions - Dressing size A=5 cm: B= 7.2 cm, pad size a = 3 cm, b = 4 cm - Dressing size A=5 cm: B= 8 cm, pad size a = 2.2 cm, b = 4.5 cm - Dressing size A=6 cm: B= 7 cm, pad size a = 3 cm, b = 4 cm, - Dressing size A= 6 cm, B= 10 cm, pad size a = 3 cm, b = 7 cm, - Dressing size A= 8 cm, B= 10 cm pad size a = 4 cm, b = 6.5 cm, - Dressing size A= 9 cm, B= 10 cm pad size a = 5 cm, b = 6 cm, - Dressing size A= 9 cm, B= 15 cm, pad size a = 5 cm, b = 10 cm - Dressing size A= 10 cm, B= 20 cm, pad size a = 5 cm, b = 15 cm - Dressing size A= 10 cm, B= 25 cm, pad size a = 5 cm, b = 20 cm, - Dressing size A= 10 cm, B= 30 cm, pad size a = 5 cm, b = 25 cm - Dressing size A= 10 cm, B= 35 cm, pad size a = 5 cm, b = 30 cm Permissible deviation from the specified dimensions $\pm 10\%$. Version 2: Plaster - a dressing on a polymer basis OMEGASOFT PU, 6 cm x 8 cm
-----	--	---

		 Dimensions Dressing size A= 6 cm, B= 8 cm, pad size a = 2 cm, b = 2 cm , Permissible deviation from the specified dimensions is $\pm 10\%$.
2.3	General description of the main functional elements	Sterile small narrow flexible strip on a polymer basis, coated on one side with a self-adhesive layer and having an absorbent pad. The rounded edges of the dressing prevent it from peeling off, prolonging the stay of the dressing on the wound bed. Absorbent pad - rectangular shape, covered with atraumatic mesh.
2.4	General description of the medical device	The dressing is hypoallergenic, does not cause skin irritation. The dressing does not contain latex, materials of human or animal origin, drugs, antimicrobial substances and does not have antimicrobial properties. Released without a doctor's prescription.
2.5	Information about sterility	Sterile; ethylene oxide sterilization
2.6	Information about single use	The dressing is a single-use medical product
2.7	Materials from which the medical device is made	Basis: Polymeric (polyurethane) film. Absorbent pad: 60% viscose, 40% terylene , covered with an atraumatic polyethylene mesh. Protective paper: Glassine. External protective polymeric (polyurethane) layer: Polymeric (polyurethane) film. Glue: Acrylic glue.
2.8	Classification according to ГОСТ-Р 53498-2019	Functional properties: fixing and protective Purpose - fixing, therapeutic and prophylactic Form - application in plates, polymer-based, with an adhesive layer, with a functional pad, with a protective coating

3 Completeness

The OMEGASOFT PU polymer-based dressing is supplied in the following composition:

Dressing	1pc;
Individual packaging with instructions printed on top	1 pc.

4 Medical device packaging

4.1 Individual packing

Each dressing is packed in an individual sterile package.

Individual packaging material: Dialysis paper; Polyethylene composite film; Dye.

The product is sold in individual packaging. In the quantity confirmed by the customer, the product is completed in group and transport packaging. The maximum number of products in group and transport packaging is given in tables 4.5.

Material:

dialysis paper; the brand is not applied; density 70 mg/m²;

Polyethylene composite film; the brand is not applied; density 100 g/m²;

Permissible deviation from the specified parameters is $\pm 10\%$.

Dye: the brand is not applied; manufacturer Ningbo Ltd. jianxiang Packing Material Technology".

Manufacturer: Yafho LLC Bio Technology Co., Ltd., China.

Overall dimensions and weight of individual packaging are given in table 3.

Table 3

Dressing size, cm	Overall dimensions of individual packaging, cm, $\pm 10\%$	Weight, g, $\pm 10\%$
Plaster - polymer-based dressing OMEGASOFT PU, dimensions, cm: 5x7.2; 5x8; 6x7; 6x10; 8x10; 9x10; 9x15; 10x20; 10x25; 10x30; 10x35		
5x7.2	12x9	2.2
5x8	13x9	2.0
6x7	12x9	2.1
6x10	15x9	3.0
8x10	15x12	4.0
9x10	15x12	4.5
9x15	20x12	7.1
10x20	25x12.8	10.0
10x25	30x12.8	12.5
10x30	35x12.8	15.0
10x35	40x12.8	17.5
Plaster - polymer-based dressing OMEGASOFT PU, 6 cm x 8 cm		
6 x 8	13x9	2.4

Group packaging

Dressings in individual packaging are placed in group packaging in the amount in accordance with table 4.

Table 4

Dressing size, cm	Number of products in a group package, pcs.	Maximum Number of products in transport packaging, pcs.
Plaster - polymer-based dressing OMEGASOFT PU, dimensions, cm: 5x7.2; 5x8; 6x7; 6x10; 8x10; 9x10; 9x15; 10x20; 10x25; 10x30; 10x35		
5x7.2	one hundred	2000
5x8	one hundred	2000
6x7	one hundred	2000
6x10	50	1000
8x10	thirty	720
9x10	thirty	720
9x15	thirty	720
10x20	25	450
10x25	25	450
10x30	25	450
10x35	25	450
Plaster - polymer-based dressing OMEGASOFT PU, 6 cm x 8 cm		
6 x 8	one hundred	2400

Group packaging material: White cardboard; Dye.

4.2 Transport packaging

Dressings in a group package are placed in a transport package in the amount in accordance with Table 3.

Transport packaging material: reinforced multilayer corrugated cardboard; Dye.

5 Information about the sterilization of a medical device

Sterilization is carried out with ethylene oxide in accordance with the requirements of GB 18279.1. Sterilization parameters are shown in Table 5.




Table 5

Sterilization temperature	Humidity	Prevacuum	Sterilization duration	Temperature maintenance time	Pressure maintenance time	Ethylene oxide concentration
+54±4°C	30%-85%	-70kPa	480 min	60 min	5-10 min	0.6 kg/m ³

6 Operating instructions

Table 6

1	Precautions and Warnings	1. Single use only. Reuse is prohibited. 2. Do not use if any signs of packaging damage are visible. Dressing in intact packaging is sterile for 3 years.
---	--------------------------	--

		3. If irritation occurs, discontinue use of the dressing and seek medical advice. 4. After use, dispose of the dressing in accordance with the disposal instructions.
2	Operating conditions	Dressing operating conditions: temperature from +10°C to +35°C; relative humidity is not more than 90%.
3	User manual	Step 1 : Before applying the dressing, make sure the skin around the wound is clean and dry after cleaning with an antiseptic liquid. Open the package and take the dressing out of the package.  Step 2: Remove the protective paper backing from the dressing, exposing the sticky surface.  Step 3: Apply a dressing to the wound, pressing the absorbent pad from the center to the edges. Press down on the rest of the dressing to ensure full and comfortable contact with the skin without wrinkling. 
4	Maximum duration of skin contact	Prolonged contact (up to 72 hours) with intact skin and wound surface

14

5	Maintenance Information	The dressing is a single-use medical device and is not subject to maintenance.
6	Disposal instructions	Recycling and disposal of unused products, after the expiration date, is carried out in accordance with the requirements for class A waste (СанПин 2.1.3684-21). Used products are disposed of as class B waste (СанПин 2.1.3684-21).

7 Storage and transportation conditions

Table 7

Storage conditions	Store in transport packaging at temperatures from +10°C to +35°C; relative humidity not more than 90% in a dry place, protected from direct sunlight, away from heat sources.
Transportation conditions	Transportation is carried out in transport packaging by all means of transport in covered vehicles in accordance with the rules for the carriage of goods in force for this type of transport. Transportation conditions: Temperature from -10 °C to +49°C; relative humidity is not more than 90%.

8 Warranty

The manufacturer guarantees that the dressing complies with the manufacturer's documentation, subject to the precautions and the conditions of operation, storage and transportation.

The dressing has a warranty period of 3 years from the date of manufacture. The production date is the same as the sterilization date.

Questions about the quality of dressings are accepted by an authorized representative of the manufacturer:

Omegamed Limited Liability Company (Omegamed LLC)

Abbreviated name of the enterprise: Omegamed LLC

Address: 620014, Sverdlovsk region, Yekaterinburg, Lenina avenue, building 5b, office 6

Contact details:

Phone: +7 (343) 372-09-72

Email: www.omegasoft.ru

ATTENTION! If an undesirable event occurs during the use of dressings, which has signs of an adverse event (incident), it is necessary to send a message about this event to an authorized representative of the manufacturer.

9 List of foreign standards used by the manufacturer (manufacturer) of medical devices

GB 18279.1-2015 "Sterilization of healthcare products - Ethylene oxide - Part 1: Requirements for the development, validation and routine control of sterilization for medical devices";

ISO 11135:2014 "Sterilization of healthcare products - Ethylene oxide - Part 1: Requirements for the development, validation and routine control of a sterilization process for medical devices";

ISO 15223-1:2016 "Medical devices - Symbols used for marking on medical devices, on labels and in accompanying documentation - Part 1: General requirements";

EN 556-1:2001 Sterilization of medical devices. Requirements for medical products of the "sterile" category. Part 1. Requirements for the sterilization of packaged medical devices";

EN 556-2:2015 Sterilization of medical devices. Requirements for medical products of the "sterile" category. Part 2. Requirements for medical devices manufactured under aseptic production conditions";

ISO 11607-1:2019 "Packaging of medical devices subject to terminal sterilization - Part 1: Requirements for materials, sterile barrier systems and packaging systems";

ISO 11607-2:2019 "Packaging of terminally sterilized medical devices - Part 2: Requirements for the validation of forming, packaging and assembly processes";

ISO 14971:2019 "Medical devices - Application of risk management to medical devices";

ISO 13485:2016 "Medical devices - Quality management systems - Regulatory requirements".



ГОСТ R 53498-2019 "Patch type medical products. General technical requirements. Test Methods»

ГОСТ EN ISO 11737-1-2006 Sterilization of medical products. microbiological methods. Part 1. Determination of the population of microorganisms on products.

APPENDIX A

Deciphering the symbols used in marking

Table A.1

Symbol	Description
	Avoid exposure to sunlight
	Keep away from moisture
	Reuse prohibited
	Do not use if packaging is damaged
	Latex free
	Ethylene oxide sterilization
	EU Compliant
	Batch code
	Use until ... (yy.mm.dd)
	Date of manufacture

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