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CERTIFICATE

01325800

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INSTRUCTIONS FOR USE
OF THE MEDICAL DEVICE FOR IN VITRO DIAGNOSTICS

**Test container for the detection of psychoactive substances and their
metabolites in urine by immunochromatographic method for in vitro
diagnostics**

("NARCOLAB test system"), in versions

Operation manual VER 07.06.2022

2022

Table of contents

1. Description of the medical device	4
1.1. Name of medical device	4
1.2. Purpose of the medical device	5
1.3. Scope and potential users	10
1.4. Indications and contraindications for use of the medical device.....	10
1.5. Diagnostic significance.....	11
1.6. Description of target analytes	11
1.6.1. Amphetamine (AMP)	11
1.6.2. Barbiturates (BAR).....	12
1.6.3. Benzodiazepines (BZO)	12
1.6.4. Cocaine (COC)	13
1.6.5. Oxycodone (OXY).....	13
1.6.6. Methamphetamine (MET)	13
1.6.7. Methadone (MTD).....	14
1.6.8. Opiates (OPI).....	14
1.6.9. Synthetic cannabinoids (K3).....	14
1.6.10. Marijuana (THC)	15
1.6.11. Synthetic cathinones (Methylenedioxypropylone – MDPV)	15
1.6.12. Phencyclidine (PCP).....	15
1.6.13. Ethylglucuronide (ETG).....	16
1.6.14. Methylenedioxymethamphetamine (MDMA).....	16
2. Method principle.....	17
3. Characteristics of a medical device	18
3.1. Specifications of the medical device	18
3.2. Diagnostic sensitivity and specificity	20
3.3. Analytical sensitivity	30
3.4. Analytical specificity	36
3.4.1. Interference	36
3.4.2. Cross-reactivity	36
3.5. Repeatability	49
4. Precautions and biosafety data	56

4.1. Precautions.....	56
4.2 Biosafety data for users, patients and others	56
5. Preparing the device for use	57
6. Collection and storage of urine samples.....	57
6.1. Methods of collecting biological material and precautions.....	58
7. Equipment and materials not included in the test set, but required for the analysis.....	58
8. Methodology for analysis	58
9. Interpretation of results.....	60
9.1. Results of tests to detect psychoactive substances	60
9.2. Results of a test to establish adulteration	61
10. Quality Control	62
11. Limitations of the test system	63
12. Supply package.....	63
13. Conditions of use, storage and transport	65
13.1. Operating conditions.....	65
13.2. Storage conditions	65
13.3. Transport conditions	65
14. Shelf life and stability.....	65
14.1. Real-time stability.....	66
14.2. Accelerated stability	66
14.3. Stability after opening the packaging	67
14.4. Stability during transport.....	68
15. Disposal and destruction of the test system.....	69
16. Information on the developer and manufacturer of the medical device.....	70
17. Manufacturer's authorised representative in the Russian Federation	70
18. Warranty obligations	70
19. List of normative documents	71
20. Literature references	71
Appendix A List of control materials used for external quality control of the test system procedure	73
Appendix B: Symbols used in labelling	74

1. Description of the medical device

1.1. Name of medical device

Test container for the detection of psychoactive substances and their metabolites in urine by immunochromatographic method for in vitro diagnostics ("NARCOLAB test system"), in versions:

Version 1.

Test container for the determination of 10 types of psychoactive substances and their metabolites with adulterated urine test strip for creatinine, pH, specific gravity and a thermal strip for sample temperature control.

The composition includes:

- Disposable container with key and test strip for the determination of 10 types of psychoactive substances and their metabolites: amphetamine, barbiturates, benzodiazepines, cocaine, methamphetamine, methadone, opiates, phencyclidine, synthetic cathinones, marijuana with adulteration control strip for creatinine, pH, specific gravity and thermal strip for sample temperature control - 25 pieces
- Tamper-evident tape - 25 pieces.
- Colour chart for determining adulteration - 1 pc.
- Instructions for use - 1pc.

Version 2.

Test container for the determination of 10 types of psychoactive substances and their metabolites with adulterated urine test strip for creatinine, pH, oxidants and a thermal strip for sample temperature control

The composition includes:

- Disposable container with key and test strip for the determination of 10 types of psychoactive substances and their metabolites: amphetamine, barbiturates, benzodiazepines, cocaine, methamphetamine, methadone, opiates, oxycodone, synthetic cannabinoids (K3), marijuana with adulteration control strip for

creatinine, pH, oxidants and thermal strip for sample temperature control - 25 pieces

- Tamper-evident tape - 25 pieces.
- Colour chart for determining adulteration - 1 pc.
- Instructions for use - 1pc.

Version 3.

Test container for the determination of 12 types of psychoactive substances and their metabolites with adulterated urine test strip for creatinine, pH, specific gravity and a thermal strip for sample temperature control.

The composition includes:

- Disposable container with key and test strip for the determination of 12 types of psychoactive substances and their metabolites: amphetamine, barbiturates, benzodiazepines, cocaine, methamphetamine, methadone, opiates, phencyclidine, synthetic cathinones, marijuana, ethylglucuronide, methylenedioxymethamphetamine with adulteration control strip for creatinine, pH, specific gravity and thermal strip for sample temperature control - 25 pieces
- Tamper-evident tape - 25 pieces.
- Colour chart for determining adulteration - 1 pc.
- Instructions for use - 1pc.

1.2. Purpose of the medical device

"The test container for detection of psychoactive substances and their metabolites in urine by immunochromatographic method for in vitro diagnosis ("NARCOLAB test system"), in versions" (hereinafter referred to as the test system) is intended for collection of biomaterial and simultaneous qualitative determination of psychoactive substances and their metabolites in human urine as well as for control of urine adulteration.

Functional purpose is to screen for psychoactive substances. Screening analysis of various psychoactive substances in urine ranges from simple immunochromatographic tests to complex analytical procedures. The speed and sensitivity of immunochromatographic tests ensure their widespread use in the analysis of urine for various psychoactive substances. The test system is a device for immunochromatographic analysis for the qualitative detection of substances, without the need for additional devices. The threshold concentrations and analytes for the psychoactive substances to be determined are given in Tables 1, 2 and 3.

Table 1. Threshold concentrations and calibrator for version 1

Parameter	Parameter name	Calibrator	Threshold value (ng/ml)	Detection limit* (ng/ml)
AMP	Amphetamine	D-amphetamine	300	375
BAR	Barbiturates	Secobarbital	300	375
BZO	Benzodiazines	Oxazepam	100	125
COC	Cocaine	Benzoyllecgonine	100	125
MET	Methamphetamine	Methamphetamine	300	375
MTD	Methadone	Methadone	300	375
OPI	Opiates	Morphine	1000	1250
THC	Marijuana	11 nor- Δ^9 -THC-9-COOH	25	31
PCP	Phencyclidine	Phencyclidine	25	31
MDPV	Synthetic cathinones	Methylenedioxypropylone	500	625

Table 2. Threshold concentrations and calibrator for version 2

Parameter	Parameter name	Calibrator	Threshold value (ng/ml)	Detection limit* (ng/ml)
AMP	Amphetamine	D-amphetamine	500	625
BAR	Barbiturates	Secobarbital	300	375
BZO	Benzodiazines	Oxazepam	300	375
COC	Cocaine	Benzoyllecgonine	150	188
OXY	Oxycodone	Oxycodone	100	125
MET	Methamphetamine	Methamphetamine	500	625
MTD	Methadone	Methadone	300	375
OPI	Opiates	Morphine	300	375
K3	Synthetic cannabinoids	AB- PINACA	20	25
THC	Marijuana	11 nor- Δ^9 -THC-9-COOH	20	25

Table 3. Threshold concentrations and calibrator for version 3

Parameter	Parameter name	Calibrator	Threshold value (ng/ml)	Detection limit* (ng/ml)
AMP	Amphetamine	D-amphetamine	300	375
BAR	Barbiturates	Secobarbital	300	375
BZO	Benzodiazines	Oxazepam	100	125
COC	Cocaine	Benzoyllecgonine	100	125
MET	Methamphetamine	Methamphetamine	300	375
MTD	Methadone	Methadone	300	375
OPI	Opiates	Morphine	1000	1250

THC	Marijuana	11 nor- Δ^9 -THC-9-COOH	25	31
PCP	Phencyclidine	Phencyclidine	25	31
MDPV	Synthetic cathinones	Methylenedioxypyrovalerone	500	625
ETG	Ethylglucuronide	Ethylglucuronide	500	625
MDMA	Methylenedioxymethamphetamine	3,4-Methylenedioxymethamphetamine	500	625

***The threshold value (ng/ml)** is the concentration of an analyte at which the test system can detect both positive and negative results.

The detection limit (ng/ml) is the concentration of an analyte at which the test system detects a result as positive with a probability of 95% or higher.

The test system is an immunochromatographic rapid assay for the qualitative preliminary detection of psychoactive substances, their metabolites and combinations specified in Tables 1, 2 and 3 in a human urine sample. Each test uses monoclonal antibodies to selectively detect elevated preparation levels in the urine. At a concentration equal to or higher than the detection limit (threshold value +25%), the test system detects a positive result in more than 95% of cases. This analysis is recommended by the United States Substance Abuse and Mental Health Services Administration (SAMHSA).

Specific pathology, condition or risk factor:

The test system is used for the medical examination of persons for suspected substance use, pre-trip and post-trip examinations, admission to hospital and other high-risk cases.

Type of sample to be analysed:

- freshly collected urine stored at 18 to 25°C for up to 4 hours;
- urine stored in the refrigerator at 2 to 8°C for up to 48 hours;
- frozen urine stored at minus 20 to minus 30°C for up to 30 days.

The test results are not influenced by population and demographic aspects of the test system.

The test container is equipped with a test strip to detect sample adulteration and a thermal strip to control the temperature of the sample.

Adulteration is the tampering with a urine sample in order to alter the test results. The use of impurities can lead to false negatives in tests for psychoactive substances, or interfere with the screening test, and/or destroy the psychoactive substances present in the urine. Dilution can also be used to obtain false negative results for substance tests.

One of the best ways to check for adulteration or dilution is to determine certain urine characteristics such as creatinine, pH, oxidants, specific gravity and sample temperature.

Creatinine (CRE): test for sample dilution. Creatinine is the metabolic product of creatine as an amino acid. It is found in muscle tissue and is also present in the urine. In an attempt to obtain false test results for the banned preparations, the test subject may drink excessive amounts of water or use diuretics (diuretics), such as herbal tea, to flush the preparation out of the body.

Creatinine determinations and specific gravity values are methods of testing for dilution and attempted drug excretion are two of the most common methods used to give false positive test results for the banned preparations. Low creatinine levels and low specific gravity may indicate diluted urine. Lack of creatinine (<5 mg/dL) indicates that the sample is not human urine, or is highly diluted.

Specific gravity (SG): test for sample dilution. The normal value should be between 1.003 and 1.030. Values outside this range may be the result of sample dilution or adulteration.

pH: test for the presence of acidic or alkaline substances added to urine for adulteration purposes. The normal pH level should be between 4.0 and 9.0. Values outside this range may indicate an adulterated sample.

Oxidants (OXI): test for the presence of oxidising reagents such as bleach, hydrogen peroxide, pyridinium chlorochromate. Oxidants are widely used

substances for adulteration of urine samples. Oxidants should not be found in normal human urine. Determination with a test strip of oxidant concentrations:

- Sodium hypochlorite ≥ 100 ppm;
- pyridinium chlorochromate ≥ 10 mg/dl;
- hydrogen peroxide ≥ 20 ppm;

is not the norm and is seen as a fact of adulteration.

Sample temperature

The container is fitted with a thermal strip to determine the temperature of the sample. The appearance of a green mark under a certain number (32 to 38) indicates the temperature of the urine sample. The proper temperature range for the sample is 32°C to 38°C.

1.3. Scope and potential users

Scope: Clinical laboratory diagnostics. The test container is intended to be used in a hospital laboratory for the rapid detection of psychoactive substances and their metabolites in human urine.

Potential user: a laboratory technician with a medical degree as well as a laboratory technician with a specialized medical degree.

For in vitro diagnostics only.

The test container is intended for single use.

1.4. Indications and contraindications for use of the medical device

Indications: Detection of psychoactive substances and their metabolites in human urine.

Version 1 is designed to detect 10 types of psychoactive substances and their metabolites: amphetamine, barbiturates, benzodiazepines, cocaine, methamphetamine, methadone, opiates, phencyclidine, synthetic cathinones, marijuana as well as to control urine adulteration by basic parameters: sample temperature, pH, creatinine, specific gravity.

Version 2 is designed to detect 10 types of psychoactive substances and their metabolites: amphetamine, barbiturates, benzodiazepines, cocaine, methamphetamine, methadone, opiates, oxycodone, synthetic cannabinoids (K3),

marijuana as well as to detect urine adulteration by basic parameters: sample temperature, pH, creatinine oxidants.

Version 3 is designed to detect 12 types of psychoactive substances and their metabolites: amphetamine, barbiturates, benzodiazepines, cocaine, methamphetamine, methadone, opiates, marijuana, phencyclidine, synthetic cathinones, ethyl glucuronide, methylenedioxymethamphetamine as well as to control urine adulteration by basic parameters: sample temperature, pH, creatinine, specific gravity.

Contraindications: No contraindications within the scope of the manufacturer's prescription.

1.5. Diagnostic significance

The test container is used to detect psychoactive substances and their metabolites in urine and can be used in investigations wherever a result of this test is required.

This analysis provides only a preliminary analytical result. A more precise alternative chemical method should be used to confirm the analytical result. Gas chromatography with mass spectrometer (GC/MS) or liquid chromatography with mass spectrometer (HPLC) are the preferred confirmatory methods.

The urine adulterant test strip is a semi-quantitative test based on colour comparison to determine creatinine, oxidants, pH and specific gravity of human urine. This test provides only preliminary data. A more precise alternative chemical method must be used to obtain a validated analytical result. The results should be sent to the laboratory for confirmation.

1.6. Description of target analytes

1.6.1. Amphetamine (AMP)

Amphetamines belong to the category of potent adrenomimetic drugs used for therapeutic purposes. Amphetamines are chemically related to human natural catecholamines: adrenaline and norepinephrine. High doses of amphetamine lead to increased central nervous system stimulation and induce a state of euphoria and heightened alertness, reduce appetite and cause a feeling of increased energy and

strength. Cardiovascular reactions to amphetamines include increased blood pressure and cardiac arrhythmia. More acute reactions consist of anxiety disorders, paranoia, hallucinations and psychotic behaviour. The effects of amphetamines last for 2-4 hours after consumption, with a half-life of 4-24 hours. About 30% of amphetamines are excreted unchanged in the urine with residues of hydroxylation products and desamidised derivatives.

1.6.2. Barbiturates (BAR)

Barbiturates weaken the central nervous system. They are used therapeutically as sedatives, hypnotics and anticonvulsants. Barbiturates are usually taken orally as capsules or tablets. The effect is similar to that of alcohol intoxication. Long-term use of barbiturates leads to immunological tolerance and physical dependence. Rapid-acting barbiturates taken at 400 mg/day for 2-3 months form a clinically significant degree of physical dependence. Withdrawal symptoms can be severe enough to be fatal. Very small amounts (less than 5%) of barbiturates are excreted unchanged in the urine. The detection period for barbiturates in urine is 4-7 days.

1.6.3. Benzodiazepines (BZO)

Benzodiazepines are medications often prescribed for symptomatic treatment of anxiety and sleep disorders. They act through specific receptors using gamma-aminobutyric acid (GABA). Being safer and more effective, benzodiazepines have replaced barbiturates in the treatment of anxiety and sleep disorders. Benzodiazepines are also used as sedatives before surgical and medical procedures and to treat epilepsy and alcohol withdrawal.

Taking benzodiazepines regularly (e.g. daily) for several months increases the risk of physical dependence, especially if large doses are taken. Abrupt discontinuation can lead to symptoms such as sleep disorders, gastrointestinal distress, impaired well-being, loss of appetite, sweating, shivering, weakness, anxiety disorders and altered perception. Only tiny amounts of most benzodiazepines (less than 1%) are excreted unchanged in the urine; concentrations in the urine are formed by paired preparations. The detection period

for benzodiazepines in urine is 3-7 days.

1.6.4. Cocaine (COC)

Cocaine is a potent central nervous system stimulant and local anaesthetic. It initially causes excessive energy and agitation, gradually leading to tremor, hypersensitivity and spasms. In large quantities, cocaine causes feverishness, insensibility, difficulty breathing and loss of consciousness.

Cocaine is often taken on its own by inhaling through the nose, by intravenous injections and by smoking pure cocaine. Rapidly excreted in the urine as benzoylecgonine. Benzoylecgonine, the main metabolite of cocaine, has a longer biological half-life (5-8 hours) than cocaine (0.5-1.5 hours) and can be detected within 24-48 hours of cocaine exposure.

1.6.5. Oxycodone (OXY)

Oxycodone is a semi-synthetic opioid derived from thebaine. It is a powerful pain reliever. The preparation has a destructive effect on the body. Dependence on oxycodone develops rapidly. Very often heroin addicts use it under the illusion that by using it they can get rid of their addiction, but instead of relief they get a new addiction of the same strength and severity as the old one.

Once in the body, oxycodone targets opioid receptors that are located in both the brain and spinal cord. This is followed by reactions that can not only quench the pain but also produce a feeling of euphoria. Oxycodone can overdose if used uncontrollably, it has a depressant effect on the respiratory centre and can be fatal if not treated in time.

1.6.6. Methamphetamine (MET)

Methamphetamine is an addictive, stimulant drug that intensely excites certain brain systems. Methamphetamine is chemically very similar to amphetamine, but its effects on the central nervous system are much stronger. Methamphetamine produced in illegal laboratories is highly addictive. The preparation is taken orally, by injection or inhalation. High doses cause increased central nervous system stimulation and induce a state of euphoria and heightened alertness, reduce appetite and induce a feeling of increased energy and strength.

Cardiovascular reactions to methamphetamine include increased blood pressure and cardiac arrhythmia. More acute reactions consist of anxiety disorders, paranoia, hallucinations, psychotic behaviour, depression and exhaustion.

The effects of methamphetamine last for 2-4 hours after consumption, with a half-life of 9-24 hours. Methamphetamine is excreted in the urine as amphetamine, oxidised and deamidised derivatives. However, 10-20% of methamphetamine is excreted unchanged in the urine. Thus, the presence of the parent compound in urine indicates the use of methamphetamine. Methamphetamine can usually be detected in the urine within 3-5 days, depending on the pH of the urine.

1.6.7. Methadone (MTD)

Methadone is a long-acting painkiller with an effect lasting from twelve to forty-eight hours. Methadone relieves consumers from the strain of seeking illegal heroin, the dangers of injecting and the mood swings caused by most opiates. Prolonged use of high doses of methadone can lead to very long withdrawal. Compared to other opiates, which require a withdrawal period of a week to ten days, methadone users can go through a withdrawal period of five to six weeks.

1.6.8. Opiates (OPI)

Opiates include any drug derived from the opium poppy, including natural substances, morphine and codeine and semi-synthetic preparations such as heroin. Opioid is a more general name that refers to any preparation that acts on opioid receptors.

Opioid analgesics make up a large group of substances that control pain by inhibiting the central nervous system. Large doses of morphine can cause high levels of tolerance and physical dependence and lead to addiction and abuse. Morphine is excreted unmetabolised and is the main metabolic product of codeine and heroin. Morphine can be detected in urine for several days after an opioid dose.

1.6.9. Synthetic cannabinoids (K3)

Synthetic cannabinoids are smoking mixtures sold as weed with the chemical AB- PINACA on them. They have a psychoactive effect somewhat

similar to that of marijuana, but with more severe psychopathological symptoms. Like marijuana, K3 affects cannabinoid receptors. In acute intoxication with cannabinoids, spatial perception and movement disorders may occur. The world is perceived as unfamiliar, alienated, psychomotor agitation occurs, accompanied by feelings of anxiety and fear. Along with JWH-073 and JWH-018, AB- PINACA is one of the most common synthetic cannabinoids among new psychoactive substances.

1.6.10. Marijuana (THC)

THC (Δ^9 - tetrahydrocannabinol) is the main active component of cannabinoids (marijuana).

When smoked or taken orally it induces a state of euphoria. Consumers have impaired memory and delayed comprehension. They may also experience short-term mental confusion and fear. Prolonged use of large doses can lead to behavioural disorders. The peak effect of marijuana smoking occurs after 20-30 minutes and lasts 90-120 minutes after smoking one cigarette. Elevated metabolic levels in urine are observed over time of exposure and can be detected within 3-10 days of smoking. The main metabolite excreted in the urine is 11-nor- Δ^9 - tetrahydrocannabinol - 9-carboxylic acid (Δ^9 - THC-COOH).

1.6.11. Synthetic cathinones (Methylenedioxypropylone – MDPV)

Methylenedioxypropylone is a cathinone-class psychostimulant that acts as an inhibitor of noradrenaline dopamine reuptake. MDPV is not used for medicinal purposes, but has been sold as a legal psychostimulant since 2007. It has similar features to amphetamine and cocaine. Causes effects such as: palpitations, increased blood pressure, vasoconstriction, sweating, restlessness, seemingly decreased need for food and sleep. Tachycardia and hypertension last 3 to 4 hours, mild stimulation lasts 6 to 8 hours. High doses can cause prolonged panic attacks, psychosis, lack of sleep and hallucinations in people who cannot tolerate stimulants. The detection period for MDPV in urine is 2-3 days.

1.6.12. Phencyclidine (PCP)

A preparation for general anaesthesia. It began to be used for pain relief in

the middle of the last century, and at that time it was used exclusively for medical purposes. Phencyclidine is currently listed as a psychoactive substance and belongs to the group of dissociative anaesthetic hallucinogens. The effects of the drug can be felt after half an hour or an hour if smoked, and after a couple of minutes if injected. The main effect of the drug on the addict's psyche is to develop a state similar to that of heavy alcoholic intoxication combined with euphoria. But this effect is only seen in small doses that are individual for each body. When the dose is increased, the person feels negative emotions, detachment from the world, a feeling of brokenness, sleepiness and apathy. On average, phencyclidine lasts up to 6-10 hours, but higher doses can lead to health problems for several weeks.

1.6.13. Ethylglucuronide (ETG)

Ethylglucuronide (ETG) is a minor non-oxidised metabolite of ethyl alcohol formed in vivo by the conjugation of ethanol to glucuronic acid by UDP-glucuronosyltransferase. ETG, which is a metabolic product of rapidly metabolised alcohol (ethanol) in the body, is excreted in the blood, hair and urine. An ETG rapid test panel can detect ethylglucuronide in urine, which will confirm alcohol consumption. The ETG metabolite stays in the body longer and therefore has a longer detection time (8 to 80 hours). The ETG test can be used in alcohol withdrawal or rehabilitation programmes.

1.6.14. Methylenedioxymethamphetamine (MDMA)

Methylenedioxymethamphetamine (Ecstasy) is a designer drug first synthesised in 1914 by the German Pharmaceutical Company to treat obesity. Those who took the drug often reported side effects such as increased muscle tension and sweating. Although methylenedioxymethamphetamine is clearly not a stimulant, like amphetamine it can increase blood pressure and heart rate. Methylenedioxymethamphetamine causes changes in perception, with increased sensitivity to light, difficulty in focusing and blurred vision in some people who use the substance. The mechanism of action of methylenedioxymethamphetamine is considered to be through the release of Serotonin, a neurotransmitter. Methylenedioxymethamphetamine can also release Dopamine, but the general

consensus is that this is a minor effect of the drug (Nichols and Oberlender, 1990). The most common effect of methylenedioxymethamphetamine, occurring in almost all people who took a moderate dose of the drug, was jaw clenching.

2. Method principle

The test container for the detection of drug substances and their metabolites in urine by immunochromatographic method works according to the principle of one-step immunochemical analysis based on the principle of competitive binding. Drug substances that can be detected in the urine sample compete with labelled drug substances (drug-protein conjugate) pre-applied to the membranes of the test strips in the area of each test strip for binding zones on specific antibodies.

During the test, a sample of urine is lifted upwards by capillary force. If the concentration of the psychoactive substance in the urine sample is below the threshold value, it does not saturate the binding sites of the specific antibody. The antibody then reacts with the drug-protein conjugate and a coloured line appears in the test area (T) of the strip for that psychoactive substance. The presence of a psychoactive substance at a concentration equal to or higher than the threshold value of determination leads to saturation of all antibody binding sites, and no coloured line is formed in the test area.

A urine sample positive for the psychoactive substance does not form a coloured line in the test area of the strip due to competition of the psychoactive substance for binding sites, whereas a urine sample negative for the preparation forms a line in the test area due to lack of competition. The coloured line is always displayed in the control area, acting as an internal (procedural) control and indicating that sufficient sample volume has been added and that the urine is flowing evenly across the membrane.

The test container is also equipped with strips to detect sample adulteration. The sample adulteration test is designed to detect tampering or any other manipulation of a urine sample in order to alter the test results. The use of extraneous substances for adulteration can lead to false negative test results for

psychoactive substances and their metabolites in urine by influencing the testing process and/or by diluting the psychoactive substances present in the urine.

3. Characteristics of a medical device

3.1. Specifications of the medical device


The test system is a container which has:

- Test material container with a nominal volume of 85 ± 8 ml. Container dimensions: height (84.2 ± 1) mm, length (65 ± 0.5) mm, width (62.6 ± 0.5) mm. The label on the container wall has 40ml and 60ml divisions, a minimum fill level mark (dosing accuracy of $\pm 10\%$ of the indicated volume) and a scale with a temperature indication from 32 to 38°C;
- a cylinder with a nominal volume of 2 ± 0.2 ml;
- the panel built into the container, in version 1, contains 6 test strips: 5 for the detection of psychoactive substances and their metabolites and 1 for the detection of adulterated urine. In version 2 and 3 the panel contains 7 test strips: 6 for the detection of psychoactive substances and their metabolites and 1 for the detection of adulterated urine. Panel dimensions for test strips: height (54 ± 0.5) mm, width (48 ± 0.5) mm;
- the swivel key to ensure that the test sample flows from the storage tank into the tank. Key dimensions: height (32 ± 0.5) mm, width (16 ± 0.5) mm, length (27 ± 0.5) mm.

Types of test strips for the detection of psychoactive substances in version 1:

- Test strip for the detection of BAR/MDPV/MTD;
- Test strip for the detection of THC/COC/PCP;
- Test strip for the detection of AMP/BZO;
- Test strip for the detection of MET;
- Test strip for the detection of OPI.

Types of test strips for the detection of psychoactive substances in version 2:

- Test strip for the detection of BZO/AMP/MOR(OPI);
- Test strip for the detection of BAR/MTD;
- Test strip for the detection of THC/COC;
- Test strip for the detection of K3;
- Test strip for the detection of OXY;
- Test strip for the detection of MET.

Types of test strips for the detection of psychoactive substances in E on 3:

- Test strip for the detection of BAR/MDMA/MTD;
- Test strip for the detection of BZO/AMP/OPI;
- Test strip for the detection of COC/MET/THC;
- Test strip for the detection of PCP;
- Test strip for the detection of MDPV;
- Test strip for the detection of ETG.

The structure of the test strip for the detection of psychoactive substances and their metabolites is shown in Figure 1:

- A - Sample area substrate;
- B - Conjugate application area substrate;
- C - Nitrocellulose membrane;
- D - Absorbent plate;
- E - Lower substrate.

Figure 1 – Structure of the test strip for the detection of psychoactive substances and their metabolites

The total mass of the product in versions 1 and 2 without packaging is 66 ± 0.5 g, the mass of the test container in version 3 is 67 ± 0.5 g.

3.2. Diagnostic sensitivity and specificity

Sample matching was performed on urine samples. All samples were confirmed by gas chromatography-mass spectrometry (GC/MS) or high performance liquid chromatography (HPLC). Samples were randomly selected and analysed using a test system. The results were classified as positive or negative after 5 minutes.

The test procedure was executed exactly according to the instructions provided by Assure Tech (the manufacturer of this product). The collected clinical samples were tested by rapid tests and then confirmed by GC/MS (Agilent 7890-5975) or HPLC-MS (ABI / API 3200 Q TRAP LC / MS / MS). This type of test is convenient and accessible.

The results are shown in Table 4.

Table 4. Diagnostic sensitivity and specificity

Version 1:

AMP300:

GC/MS

NARCOLAB test system		+	-
	+	125	0
	-	5	130

Diagnostic sensitivity: $125/(125+5) = 96.1\%$

Diagnostic specificity: $130/(130+0) = 100\%$

Diagnostic significance of the method: $(125+130)/(125+5+0+130) = 98.1\%$

BAR300:

GC/MS

NARCOLAB test		+	-
	+	132	3

system	-	3	162
--------	---	---	-----

Diagnostic sensitivity: $132/(132+3) = 97.8\%$

Diagnostic specificity: $162/(162+3) = 98.1\%$

Diagnostic significance of the method: $(132+162)/(132+3+3+162) = 98\%$

BZO100:

GC/MS

NARCOLAB test system		+	-
	+	47	1
	-	2	50

Diagnostic sensitivity: $47/(47+2) = 95.9\%$

Diagnostic specificity: $50/(50+1) = 98\%$

Diagnostic significance of the method: $(47+50)/(47+2+1+50) = 97\%$

COC100:

GC/MS

NARCOLAB test system		+	-
	+	165	5
	-	3	260

Diagnostic sensitivity: $165/(165+3) = 98.2\%$

Diagnostic specificity: $260/(260+5) = 98.1\%$

Diagnostic significance of the method: $(165+260)/(165+3+5+260) = 98.2\%$

MDPV 500:

GC/MS

NARCOLAB test system		+	-
	+	85	0
	-	0	151

Diagnostic sensitivity: $85/(85+0) = 100\%$

Diagnostic specificity: $151/(151+0) = 100\%$

Diagnostic significance of the method: $(85+151)/(85+0+0+151) = 100\%$

MET300:

GC/MS

NARCOLAB test system		+	-
	+	61	0
	-	2	64

Diagnostic sensitivity: $61/(61+2) = 96.8 \%$

Diagnostic specificity: $64/(64+0) = 100 \%$

Diagnostic significance of the method: $(61+64)/(61+2+0+64) = 98.4 \%$

MTD300:

GC/MS

NARCOLAB test system		+	-
	+	99	0
	-	4	110

Diagnostic sensitivity: $99/(99+4) = 96.1 \%$

Diagnostic specificity: $110/(110+0) = 100 \%$

Diagnostic significance of the method: $(99+110)/(99+4+0+110) = 98.1 \%$

OPI 1000:

GC/MS

NARCOLAB test system		+	-
	+	55	2
	-	2	48

Diagnostic sensitivity: $55/(55+2) = 96.5 \%$

Diagnostic specificity: $48/(48+2) = 96 \%$

Diagnostic significance of the method: $(55+48)/(55+2+2+48) = 96.3 \%$

PCP 25:

HPLC

NARCOLAB test system		+	-
	+	55	0
	-	0	168

Diagnostic sensitivity: $55/(55+0) = 100\%$

Diagnostic specificity: $168/(168+0) = 100 \%$

Diagnostic significance of the method: $(55+168)/(55+0+0+168) = 100 \%$

THC25:

GC/MS

NARCOLAB test system		+	-
	+	60	1
	-	2	59

Diagnostic sensitivity: $60/(60+2) = 96.8 \%$

Diagnostic specificity: $59/(59+1) = 98.3 \%$

Diagnostic significance of the method: $(60+59)/(60+2+1+59) = 97.5 \%$

Version 2:

AMP500:

GC/MS

NARCOLAB test system		+	-
	+	116	0
	-	5	150

Diagnostic sensitivity: $116/(116+5) = 95.9\%$

Diagnostic specificity: $150/(150+0) = 100 \%$

Diagnostic significance of the method: $(116+150)/(116+5+0+150) = 98.1 \%$

BAR300:

GC/MS

NARCOLAB test system		+	-
	+	132	3
	-	3	162

Diagnostic sensitivity: $132/(132+3) = 97.8\%$

Diagnostic specificity: $162/(162+3) = 98.1 \%$

Diagnostic significance of the method: $(132+162)/(132+3+3+162) = 98 \%$

BZO300:

GC/MS

NARCOLAB test system		+	-
	+	41	4
	-	2	52

Diagnostic sensitivity: $41/(41+2) = 95.3\%$

Diagnostic specificity: $52/(52+4) = 92.9\%$

Diagnostic significance of the method: $(41+52)/(41+2+4+52) = 93.9\%$

COC150:

GC/MS

NARCOLAB test system		+	-
	+	53	2
	-	2	57

Diagnostic sensitivity: $53/(53+2) = 96.4\%$

Diagnostic specificity: $57/(57+2) = 96.6\%$

Diagnostic significance of the method: $(53+57)/(53+2+2+57) = 96.5\%$

OXY100:

GC/MS

NARCOLAB test system		+	-
	+	49	0
	-	2	54

Diagnostic sensitivity: $49/(49+2) = 96.1\%$

Diagnostic specificity: $54/(54+0) = 100\%$

Diagnostic significance of the method: $(49+54)/(49+2+0+54) = 98.1\%$

MET500:

GC/MS

NARCOLAB test system		+	-
	+	63	0
	-	2	54

Diagnostic sensitivity: $63/(63+2) = 96.8\%$

Diagnostic specificity: $54/(54+0) = 100\%$

Diagnostic significance of the method: $(63+54)/(63+2+0+54) = 98.3\%$

MTD300:

GC/MS

NARCOLAB test system		+	-
	+	99	0
	-	4	110

Diagnostic sensitivity: $99/(99+4) = 96.1 \%$

Diagnostic specificity: $110/(110+0) = 100 \%$

Diagnostic significance of the method: $(99+110)/(99+4+0+110) = 98.1 \%$

OPI 300:

GC/MS

NARCOLAB test system		+	-
	+	61	1
	-	2	47

Diagnostic sensitivity: $61/(61+2) = 96.8 \%$

Diagnostic specificity: $47/(47+1) = 97.9 \%$

Diagnostic significance of the method: $(61+47)/(61+2+1+47) = 97.3 \%$

K3 20:

GC/MS

NARCOLAB test system		+	-
	+	46	1
	-	1	59

Diagnostic sensitivity: $46/(46+1) = 97.9 \%$

Diagnostic specificity: $59/(59+1) = 98.3 \%$

Diagnostic significance of the method: $(46+59)/(46+1+1+59) = 98.1 \%$

THC20:

GC/MS

NARCOLAB test system		+	-
	+	60	1
	-	2	59

Diagnostic sensitivity: $60/(60+2) = 96.8 \%$

Diagnostic specificity: $59/(59+1) = 98.3 \%$

Diagnostic significance of the method: $(60+59)/(60+2+1+59) = 97.5 \%$

Version 3:

AMP300:

GC/MS

NARCOLAB test system		+	-
	+	125	0
	-	5	130

Diagnostic sensitivity: $125/(125+5) = 96.1\%$

Diagnostic specificity: $130/(130+0) = 100 \%$

Diagnostic significance of the method: $(125+130)/(125+5+0+130) = 98.1\%$

BAR300:

GC/MS

NARCOLAB test system		+	-
	+	132	3
	-	3	162

Diagnostic sensitivity: $132/(132+3) = 97.8\%$

Diagnostic specificity: $162/(162+3) = 98.1 \%$

Diagnostic significance of the method: $(132+162)/(132+3+3+162) = 98 \%$

BZO100:

GC/MS

NARCOLAB test system		+	-
	+	47	1
	-	2	50

Diagnostic sensitivity: $47/(47+2) = 95.9\%$

Diagnostic specificity: $50/(50+1) = 98 \%$

Diagnostic significance of the method: $(47+50)/(47+2+1+50) = 97 \%$

COC100:

GC/MS

NARCOLAB test system		+	-
	+	165	5
	-	3	260

Diagnostic sensitivity: $165/(165+3) = 98.2 \%$

Diagnostic specificity: $260/(260+5) = 98.1 \%$

Diagnostic significance of the method: $(165+260)/(165+3+5+260) = 98.2 \%$

MDPV 500:

GC/MS

NARCOLAB test system		+	-
	+	85	0
	-	0	151

Diagnostic sensitivity: $85/(85+0) = 100 \%$

Diagnostic specificity: $151/(151+0) = 100 \%$

Diagnostic significance of the method: $(85+151)/(85+0+0+151) = 100 \%$

MET300:

GC/MS

NARCOLAB test system		+	-
	+	61	0
	-	2	64

Diagnostic sensitivity: $61/(61+2) = 96.8 \%$

Diagnostic specificity: $64/(64+0) = 100 \%$

Diagnostic significance of the method: $(61+64)/(61+2+0+64) = 98.4 \%$

MTD300:

GC/MS

NARCOLAB test system		+	-
	+	99	0
	-	4	110

Diagnostic sensitivity: $99/(99+4) = 96.1 \%$

Diagnostic specificity: $110/(110+0) = 100 \%$

Diagnostic significance of the method: $(99+110)/(99+4+0+110) = 98.1 \%$

OPI 1000:

GC/MS

NARCOLAB test system		+	-
	+	55	2
	-	2	48

Diagnostic sensitivity: $55/(55+2) = 96.5 \%$

Diagnostic specificity: $48/(48+2) = 96 \%$

Diagnostic significance of the method: $(55+48)/(55+2+2+48) = 96.3 \%$

PCP 25:

HPLC

NARCOLAB test system		+	-
	+	55	0
	-	0	168

Diagnostic sensitivity: $55/(55+0) = 100\%$

Diagnostic specificity: $168/(168+0) = 100 \%$

Diagnostic significance of the method: $(55+168)/(55+0+0+168) = 100 \%$

THC25:

GC/MS

NARCOLAB test system		+	-
	+	60	1
	-	2	59

Diagnostic sensitivity: $60/(60+2) = 96.8 \%$

Diagnostic specificity: $59/(59+1) = 98.3 \%$

Diagnostic significance of the method: $(60+59)/(60+2+1+59) = 97.5 \%$

ETG500:

GC/MS

NARCOLAB test		+	-
	+	50	10

system	-	12	55
--------	---	----	----

Diagnostic sensitivity: $50/(50+12) = 80.6\%$

Diagnostic specificity: $55/(55+10) = 84.6\%$

Diagnostic significance of the method: $(55+50)/(55+12+10+50) = 82.6\%$

MDMA500:

GC/MS

NARCOLAB test system		+	-
	+	85	0
	-	0	151

Diagnostic sensitivity: $85/(85+0) = 100\%$

Diagnostic specificity: $151/(151+0) = 100\%$

Diagnostic significance of the method: $(85+151)/(85+0+0+151) = 100\%$

Conclusion: Based on the results of the study, the diagnostic performance of the test system has been determined. The diagnostic sensitivity and diagnostic specificity of the tests is more than 90% for all of the stated analytes except ethylglucuronide. The diagnostic performance in determining ethyl glucuronide is over 80%.

3.3. Analytical sensitivity

The analytical sensitivity was determined by adding a control sample of the analyte to a negative urine sample at $\pm 50\%$ threshold $\pm 25\%$ threshold and threshold. Samples of the preparation were randomly selected and coded. The results were confirmed by GC/MS. Each preparation sample was analysed 50 times in accordance with the instructions for use of the drug. The analysis results shown in Table 5 were considered after 5 minutes of adding the sample.

Table 5. Results of the analytical sensitivity study

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Conclusion: Samples containing an analyte at a concentration below the threshold concentration show a negative result, samples containing an analyte at a concentration equal to the threshold concentration (cut-off) are defined as positive or negative. All samples containing the preparation at a concentration at or above the limit of detection (LoD) show a positive result. All samples not containing an

analyte are defined as negative.

3.4. Analytical specificity

3.4.1. Interference

Interference of the Ketorolac preparation at a concentration of 5×10^6 ng/ml was detected in the determination of ethyl glucuronide (ETG) and benzodiazepines (BZO). Paracetamol preparation at a concentration of 3.3×10^6 ng/ml causes interference in the determination of benzodiazepines (BZO) and ethylglucuronide (ETG). Levofloxacin preparation at a concentration of 8.3×10^5 ng/ml cause interference in the determination of benzodiazepines (BZO) and ethylglucuronide (ETG).

3.4.2. Cross-reactivity

The aim of the study was to determine the cross-reactivity of the test with substances in non-analyte urine. The following substances showed no cross-reactivity when tested with the test system at a concentration of 100 µg/ml. Table 6 lists compounds that do not exhibit cross-reactivity.

Table 6. Compounds that do not show cross-reactivity

No. pp.	Compound name	No. pp.	Compound name	No. pp.	Compound name
1	(-)- Ephedrine (except MET)	17	Chlorpheniramine	33	Pheniramine
2	(+)-Naproxen	18	Creatine	34	Phenothiazine
3	(+/-)- Ephedrine (except MET)	19	Dextromethorphan	35	Procaine
4	4-Dimethylaminoantipyrine	20	Dextrorfan tartrate	36	Protonix
5	Acetaminophen	21	Dopamine	37	Pseudoephedrine
6	Acetone	22	Erythromycin	38	Quinidine
7	Albumin	23	Ethanol	39	Ranitidine
8	Amitriptyline	24	Furosemide	40	Sertralin
9	Ampicillin	25	Glucose	41	Tyramine

10	Aspartame	26	Guaiacol glyceryl ether	42	Vitamin C (ascorbic acid)
11	Aspirin	27	Haemoglobin	43	Trimepravne
12	Benzocaine	28	Ibuprofen	44	Venlafaxine
13	Bilirubin	29	Imipramine	45	Ibuprofen
14	b-Phenylethylamine	30	Isoproterenol		
15	Coffein	31	Oxalic acid		
16	Chloroquine (excluding MET)	32	Penicillin - G		

The compounds listed in Table 7 relating to each analyte listed, respectively, exhibit cross-reactivity towards it at the indicated concentration.

Table 7. Compounds that show cross-reactivity

Version 1:	
<u>Components related to Amphetamine 20</u>	<u>ng/ml</u>
D-Amphetamine	300
L-Amphetamine	50.000
Hemisulphate salt of mephentermine	>100.000
3,4-methylenedioxyamphetamine (MDA)	625
Phentermine	625
Para-methoxyamphetamine (PMA)	625
Para-methoxymethamphetamine (PMMA)	>100.000
Tyramine	>100.000
<u>Components related to Barbiturates 20</u>	
Secobarbital	300
Allobarbital	1.250
Alphenal	625
Amobarbital	625
Aprobarbital	188
Butabarbital	94
Butalbital	2.500
Butetal	200
Cyclopentobarbital	400
Pentobarbital	1.000
Phenobarbital	300
<u>Components related to Benzodiazepines 100</u>	
Oxazepam	100
Alprazolam	42
Bromazepam	208

Chlordiazepoxide	833
Clobazam	21
Clonazepam	833
Clorazepate	1.110
Desalcflurazepam	83
Diazepam	83
Estazolam	1.667
Phentanyl	>100.000
Flunitrazepam	125
Flurazepam	>100.000
Lorazepam	417
Lormetazepam	417
Medazepam	>100.000
Midazolam	>100.000
Nitrazepam	8.333
Norchloridiazepoxide	83
Nordiazepam	167
Prazepam	>100.000
Temazepam	21
Triazolam	1.667
<u>Components related to Cocaine 100</u>	
Benzoylecgonine	100
<u>Components related to Synthetic cathinones 500</u>	
MDPV	500
<u>Components related to Methamphetamine 300</u>	
D- Methamphetamine	300

Chloroquin	7.500
Fenfluramine	12.500
l-Methamphetamine	10.000
Hemisulphate salt of mephentermine	31.250
Methyldiethylamine (MDEA)	50.000
Methylenedioxymethamphetamine (MDMA)	313
Polymethylmethacrylate (PMMA)	625
(-)-Ephedrine	2.000
<u>Components related to Methadone 300</u>	
Methadone	300
(-)-alpha-metadol	2.000
<u>Components related to Opiates 1000</u>	
Morphine	1.000
Acetylcodeine	1.000
Buprenorphine	>10.000
Codeine	1.000
Diacetylmorphine (heroin)	3.000
Dihydrocodeine	1.000
Ethylmorphine	200
Hydromorphone	25.000
Hydrocodone	50.000
Merperidine	>100.000
6-monoacetylmorphine (6-MAM)	3.000
Morphine-3-β-D-glucoronide	10.000
Nalorphine hydrochloride	>100.000
Oxycodone	>100.000
Oxymorphone	>100.000

Rifampicin	>100.000
Tebain	50.000
<u>Components related to Phencyclidine 25</u>	
Phencyclidine	25
Hydrocodone	>100.000
Hydromorphone	>100.000
4-hydroxyphencyclidine	75
<u>Components related to Marijuana 25</u>	
11 nor- Δ^9 -THC-9-COOH	25
11 nor- Δ^8 -THC-9-COOH	15
Δ^8 – Tetrahydrocannabinol	7.500
Δ^9 – Tetrahydrocannabinol	7.500
Cannabinol	10.000
Version 2:	
<u>Components related to Amphetamine 500</u>	ng/ml
D-Amphetamine	500
L-Amphetamine	50.000
3,4-methylenedioxyamphetamine (MDA)	625
Phentermine	1.250
Para-methoxyamphetamine (PMA)	625
Tyramine	100.000
<u>Components related to Barbiturates 20</u>	
Secobarbital	300
Allobarbital	1.250
Alphenal	625

Amobarbital	625
Aprobarbital	188
Butabarbital	94
Butalbital	2.500
Butetal	200
Cyclopentobarbital	400
Pentobarbital	1.000
Phenobarbital	300
<u>Components related to Benzodiazepines 300</u>	
Oxazepam	300
Alprazolam	125
Bromazepam	625
Chlordiazepoxide	2.500
Clobazam	63
Clonazepam	2.500
Clorazepate	3.330
Desalcflurazepam	250
Diazepam	250
Estazolam	5.000
Phentanyl	>100.000
Flunitrazepam	375
Flurazepam	>100.000
Lorazepam	1.250
Lormetazepam	1.250
Medazepam	>100.000
Midazolam	>100.000
Nitrazepam	25.000
Norchloridiazepoxide	250

Nordiazepam	500
Prazepam	>100.000
Temazepam	63
Triazolam	5000
<u>Components related to Cocaine 150</u>	
Benzoylecgonine	150
Cocaine	125
Ecgonine	10.000
Ecgonine methyl ester	>10.000
<u>Components related to Oxycodone 100</u>	
Oxycodone	100
Hydrocodone	6.250
Hydromorphone	50.000
Naloxone	50.000
Oxymorphone	250
<u>Components related to Methamphetamine 500</u>	
D- Methamphetamine	500
Chloroquin	12.500
Fenfluramine	12.500
l-Methamphetamine	3.125
Hemisulphate salt of mephentermine	25.000
Methyldiethylamine (MDEA)	12.500
Methylenedioxyamphetamine (MDMA)	1.875
Polymethylmethacrylate (PMMA)	625
(-)-Ephedrine	2.000

<u>Components related to Methadone 300</u>	
Methadone	300
(-)-alpha-metadol	2.000
<u>Components related to Opiates 300</u>	
Morphine	300
Acetylcodeine	150
Buprenorphine	>10.000
Codeine	250
Diacetylmorphine (heroin)	250
Dihydrocodeine	586
Ethylmorphine	200
Hydromorphone	12.500
Hydrocodone	12.500
6-monoacetylmorphine (6-MAM)	250
Morphine-3- β -D-glucoronide	2.500
Nalorphine hydrochloride	25.000
Tebain	25.000
<u>Components related to Synthetic cannabinoids 20</u>	
AB-PINACA	20
AB-PINACA-5-Pentanoic acid	20
AB-PINACA-5-hydroxypentyl	20
AB-FUBINACA	30
AB-PINACA-4-hydroxypentyl	>10.000
UR-144-5-Pentanoic acid	5.000
UR-144	>10.000
UR-144-5-hydroxypentyl	>10.000
UR-144-4-hydroxypentyl	>10.000

APINACA	>10.000
APINACA-5-hydroxypentyl	>10.000
ADB-PINACA N-(5-hydroxypentyl)	40
ADB-PINACA-Pentanoic acid	20
5-fluoro- AB-PINACA N-(4-hydroxypentyl)	40
<u>Components related to Marijuana 20</u>	
11 nor- Δ 9-THC-9-COOH	20
11 nor- Δ 8-THC-9-COOH	15
Δ 8 – Tetrahydrocannabinol	7.500
Δ 9 – Tetrahydrocannabinol	7.500
Cannabinol	10.000
Version 3:	
<u>Components related to Amphetamine 20</u>	<u>ng/ml</u>
D-Amphetamine	300
L-Amphetamine	50.000
Hemisulphate salt of mephentermine	>100.000
3,4-methylenedioxyamphetamine (MDA)	625
Phentermine	625
Para-methoxyamphetamine (PMA)	625
Para-methoxymethamphetamine (PMMA)	>100.000
Tyramine	>100.000
<u>Components related to Barbiturates 20</u>	
Secobarbital	300
Allobarbital	1.250
Alphenal	625
Amobarbital	625

Aprobarbital	188
Butabarbital	94
Butalbital	2.500
Butetal	200
Cyclopentobarbital	400
Pentobarbital	1.000
Phenobarbital	300
<u>Components related to Benzodiazepines 100</u>	
Oxazepam	100
Alprazolam	42
Bromazepam	208
Chlordiazepoxide	833
Clobazam	21
Clonazepam	833
Clorazepate	1.110
Desalcflurazepam	83
Diazepam	83
Estazolam	1.667
Phentanyl	>100.000
Flunitrazepam	125
Flurazepam	>100.000
Lorazepam	417
Lormetazepam	417
Medazepam	>100.000
Midazolam	>100.000
Nitrazepam	8.333
Norchloridiazepoxide	83
Nordiazepam	167

Prazepam	>100.000
Temazepam	21
Triazolam	1.667
<u>Components related to Cocaine 100</u>	
Benzoylecgonine	100
<u>Components related to Synthetic cathinones 500</u>	
MDPV	500
<u>Components related to Methamphetamine 300</u>	
D- Methamphetamine	300
Chloroquin	7.500
Fenfluramine	12.500
l-Methamphetamine	10.000
Hemisulphate salt of mephentermine	31.250
Methyldiethylamine (MDEA)	50.000
Methylenedioxymethamphetamine (MDMA)	313
Polymethylmethacrylate (PMMA)	625
(-)-Ephedrine	2.000
<u>Components related to Methadone 300</u>	
Methadone	300
(-)-alpha-metadol	2.000
<u>Components related to Opiates 1000</u>	
Morphine	1.000
Acetylcodeine	1.000
Buprenorphine	>10.000

Codeine	1.000
Diacetylmorphine (heroin)	3.000
Dihydrocodeine	1.000
Ethylmorphine	200
Hydromorphone	25.000
Hydrocodone	50.000
Merperidine	>100.000
6-monoacetylmorphine (6-MAM)	3.000
Morphine-3- β -D-glucoronide	10.000
Nalorphine hydrochloride	>100.000
Oxycodone	>100.000
Oxymorphone	>100.000
Rifampicin	>100.000
Tebain	50.000
<u>Components related to Phencyclidine 25</u>	
Phencyclidine	25
Hydrocodone	>100.000
Hydromorphone	>100.000
4-hydroxyphencyclidine	75
<u>Components related to Marijuana 25</u>	
11 nor- Δ^9 -THC-9-COOH	25
11 nor- Δ^8 -THC-9-COOH	15
Δ^8 – Tetrahydrocannabinol	7.500
Δ^9 – Tetrahydrocannabinol	7.500
Cannabinol	10.000

<u>Components related to Methylenedioxymethamphetamine 500</u>	
3,4-Methylenedioxymethamphetamine	500
D-amphetamine	>100.000
l-Amphetamine	>100.000
d-Methamphetamine	>100.000
l-Methamphetamine	>100.000
3,4- Methylenedioxyamphetamine	2.500
3,4- Methylenedioxyethylamphetamine	156
Paramethoxyamphetamine	50.000
Paramethoxymethamphetamine	>100.000
<u>Components related to Ethylglucuronide 500</u>	
Ethylglucuronide	500
Ethanol	>100.000
D-Glucuronic acid	>100.000
Morphine-3- β -D-glucuronide	>100.000

Conclusion: The compounds listed in Table 6 showed no cross-reactivity when tested with the test system at a concentration of 100 μ g/ml. A number of compounds listed in Table 7 exhibit cross-reactivity towards the test analytes at the indicated concentration, which is a desirable property for psychoactive substance testing techniques.

3.5. Repeatability

This study was conducted in 2 series per day and lasted 25 days for each group of three batches. It was conducted by three specialists unaware of the sample identification system. Each of the three operators analysed 2 aliquots of urine samples with different concentrations of the psychoactive substance (2 series per day). A total of 50 analyses were performed by each specialist for each concentration. The results are shown in Table 8.

Version 1:

Drug	n	AMP 300	BAR 300	BZO 100	COC 100
------	---	---------	---------	---------	---------

conc.									
Result		-	+	-	+	-	+	-	+
Negative sample	50	50	0	50	0	50	0	50	0
50% threshold value	50	50	0	50	0	50	0	50	0
75% threshold value	50	50	0	50	0	50	0	50	0
Threshold value	50	20	30	11	39	12	38	23	27
125% threshold value	50	0	50	0	50	0	50	0	50
150% threshold value	50	0	50	0	50	0	50	0	50
3X threshold value	50	0	50	0	50	0	50	0	50

Drug conc.	n	MDPV 500		MET 300		MTD 300		OPI 1000	
Result		-	+	-	+	-	+	-	+
Negative sample	50	50	0	50	0	50	0	50	0
50% threshold value	50	50	0	50	0	50	0	50	0
75% threshold value	50	50	0	50	0	50	0	50	0
Threshold value	50	17	33	15	35	6	44	13	37
125% threshold value	50	0	50	0	50	0	50	0	50

150% threshold value	50	0	50	0	50	0	50	0	50
3x threshold value	50	0	50	0	50	0	50	0	50

Drug conc.	n	PCP 25		THC 25	
Result		-	+	-	+
Negative sample	50	50	0	50	0
50% threshold value	50	50	0	50	0
75% threshold value	50	50	0	50	0
Threshold value	50	9	41	11	39
125% threshold value	50	0	50	0	50
150% threshold value	50	0	50	0	50
3x threshold value	50	0	50	0	50

Version 2:

Drug conc.	n	AMP 500		BAR 300		BZO 300		COC 150	
Result		-	+	-	+	-	+	-	+
Negative sample	50	50	0	50	0	50	0	50	0
50% threshold value	50	50	0	50	0	50	0	50	0
75% threshold value	50	50	0	50	0	50	0	50	0
Threshold value	50	14	36	11	39	17	33	24	26
125% threshold value	50	0	50	0	50	0	50	0	50
150% threshold value	50	0	50	0	50	0	50	0	50
3x threshold value	50	0	50	0	50	0	50	0	50

Drug conc.	n	OXY 100		MET 500		MTD 300		OPI 300	
Result		-	+	-	+	-	+	-	+
Negative sample	50	50	0	50	0	50	0	50	0
50% threshold value	50	50	0	50	0	50	0	50	0
75% threshold value	50	50	0	50	0	50	0	50	0
Threshold value	50	19	31	10	40	6	44	24	26

125% threshold value	50	0	50	0	50	0	50	0	50
150% threshold value	50	0	50	0	50	0	50	0	50
3x threshold value	50	0	50	0	50	0	50	0	50

Drug conc.	n	K3 20		THC 20	
Result		-	+	-	+
Negative sample	50	50	0	50	0
50% threshold value	50	50	0	50	0
75% threshold value	50	50	0	50	0
Threshold value	50	16	34	11	39
125% threshold value	50	0	50	0	50
150% threshold value	50	0	50	0	50
3x threshold value	50	0	50	0	50

Drug conc.	n	AMP 300		BAR 300		BZO 100		COC 100	
Result		-	+	-	+	-	+	-	+
Negative sample	50	50	0	50	0	50	0	50	0
50% threshold value	50	50	0	50	0	50	0	50	0
75% threshold value	50	50	0	50	0	50	0	50	0
Threshold value	50	20	30	11	39	12	38	23	27
125% threshold value	50	0	50	0	50	0	50	0	50
150% threshold value	50	0	50	0	50	0	50	0	50
3X threshold value	50	0	50	0	50	0	50	0	50

Drug conc.	n	MDPV 500		MET 300		MTD 300		OPI 1000	
Result		-	+	-	+	-	+	-	+
Negative sample	50	50	0	50	0	50	0	50	0
50% threshold value	50	50	0	50	0	50	0	50	0
75% threshold value	50	50	0	50	0	50	0	50	0
Threshold value	50	17	33	15	35	6	44	13	37

125% threshold value	50	0	50	0	50	0	50	0	50
150% threshold value	50	0	50	0	50	0	50	0	50
3x threshold value	50	0	50	0	50	0	50	0	50

Drug conc.	n	PCP 25		THC 25		MDMA 500		ETG 500	
Result		-	+	-	+	-	+	-	+
Negative sample	50	50	0	50	0	50	0	50	0
50% threshold value	50	50	0	50	0	50	0	50	0
75% threshold value	50	50	0	50	0	50	0	50	0
Threshold value	50	9	41	11	39	13	37	19	31
125% threshold value	50	0	50	0	50	0	50	0	50
150% threshold value	50	0	50	0	50	0	50	0	50
3x threshold value	50	0	50	0	50	0	50	0	50

Conclusion: The test container complies with the declared characteristics of Repeatability.

4. Precautions and biosafety data

4.1. Precautions

- For *in vitro* diagnosis only.
- Do not use test containers in case of:
 - after the expiry date printed on the packaging;
 - damage to the packaging;
- Do not re-use the test container.
- Avoid cross-contamination of samples by using a new test container for each new sample.
- Read the instructions for use carefully before carrying out the test.
- The test container contains products of animal origin. A certificate of origin and/or animal health status cannot fully guarantee the absence of infectious pathogenic agents. It is therefore advisable to take the usual precautions (e.g. do not swallow or inhale).
- Do not eat, drink or smoke in the area where samples and test containers are stored. Established precautions should be observed throughout the procedure and standard procedures should be followed for the correct disposal of the samples. Wear protective clothing such as a lab coat, disposable gloves and eye protection when analysing samples.
- Humidity and temperature can have an adverse effect on the results.
- Used sets must be disposed of in accordance with local regulations.

4.2 Biosafety data for users, patients and others

The test container contains materials of animal origin. The test containers are made from animal materials, namely bovine serum albumin (BSA) fractions, goat antibodies and mouse monoclonal antibodies. The materials are non-toxic, non-poisonous and do not come into direct contact with the human body. The materials are not a pharmacological substance or a medicinal preparation. The procedures for processing, preservation, testing and handling of tissues, cells and substances of animal origin are performed in such a way as to ensure the highest level of biological safety for users, patients and others. Tissues, cells and biological

substances of animal origin used in the manufacture of the test system are obtained from animals that have been veterinarily inspected and supervised, taking into account the intended use of these tissues.

Veterinary certificates are available for all materials of animal origin. All antibodies and antigens in the test system have certificates of analysis (see risk management file). The waste produced by the use of the test is classified as class B according to Sanitary Standard 2.1.3684-21.

5. Preparing the device for use

The container is ready to use and requires no special preparation steps.

After opening the package, the test container must be used within 1 hour at a room temperature of 18°C to 25°C and a relative humidity of no more than 80%.

6. Collection and storage of urine samples

Testing can be performed immediately after the sample has been collected. Urine samples can be stored at 2-8°C for no more than 48 hours. For long-term storage, urine samples should be frozen and stored at minus 20 to minus 30°C for up to 30 days (as per Ministry of Health Order No. 40 of 27/01/2006). Frozen samples must be brought to room temperature and mixed well before testing. Do not leave urine samples prepared for analysis for more than 2 hours at room temperature and more than 4 hours in the refrigerator.

Avoid freezing and thawing the samples repeatedly.

If samples need to be sent to a laboratory for further analysis, they should be packed in accordance with all applicable regulations for the transport of biological samples.

Samples for transport can be collected in a test container or additional clean, dry biomaterial collection containers. Urine samples collected at any time of day can be used. Centrifuge turbid samples for 10 min at 1500 g, filter through an ash-free paper filter or allow to stand until clear supernatant is obtained. The resulting clear sample is pipetted into the test container in the required quantity. Only clear samples should be used for testing.

6.1. Methods of collecting biological material and precautions

The urine is collected under conditions that prevent the possibility of substitution or adulteration of the biological object.

Urine is collected by the subject into a container not exceeding the 60 ml mark, but above the filling level marked as "minimum level" on the test container label. The subject hands the container to the medical staff, then the container is closed with the lid.

The test system contains no parts that come into contact with unprotected body surfaces. Biomaterial should be collected without direct contact with unprotected areas of the body.

The biomaterial is taken with gloves on.

7. Equipment and materials not included in the test set, but required for the analysis

- A container for the collection, storage and transport of urine, duly registered;
- a stopwatch or timer;
- laboratory centrifuge (centrifugal force is 1500g);
- an ash-free paper filter;
- Medical non-sterile gloves;
- Aseptic wipes (if necessary);
- a refrigerator with a temperature between 2 and 8 °C;
- a freezer compartment with a temperature between minus 20 °C and minus 30 °C;
- pipette for transferring liquids.

8. Methodology for analysis

In order to obtain reliable results, the Instructions for Use of the test system must be strictly adhered to.

Before use, check the integrity of the test system package by visual inspection.

Do not use the product if the packaging is damaged or poorly sealed.

The test container, if stored between 2 and 8°C, must be kept at room temperature between 18 and 25°C for at least 10 minutes before analysis. The test container, if stored between 8°C and 18°C, must be kept at room temperature between 18°C and 25°C for at least 5 minutes before analysis. If the test container is stored at a temperature between 25 and 30 °C, the test can be started immediately.

The test container can be used to collect biomaterial. A duly registered urine collection and transport container other than the one supplied can also be used (see point 7 of these instructions for use).

Open the individual packaging of the test system shortly before the test.

Remove the test device from its sealed packaging and use it as soon as possible.

After opening the individual sealed package the test system has to be analysed within 1 hour at the room temperature of 18°C to 25°C and humidity of no more than 80%.

After the subject has provided the urine sample, check the lid for leaks and seal it with test tape. The volume of the urine sample provided must reach the minimum fill level mark.


Check the temperature indicator strip. The appearance of a green colour indicates the temperature of the urine sample. The temperature of the sample is assessed no later than 2 minutes after the start of analysis. The proper temperature range for a unadulterated sample is 32-38°C. Temperature monitoring only applies to freshly collected urine samples and takes place immediately after the patient has submitted the biomaterial.

Take the key out of the cover, insert it into the hole, turn it 180° in the direction of the arrow and time the start of analysis (figure 2).

Remove the label from the test result panel.

For tests for adulteration, visually compare the colour of the indicator strip with the colour chart for interpreting the results of the adulteration test (figure 6). Do not interpret the results later than 5 minutes. If the test shows an adulterated

urine sample, it is recommended not to assess the content of psychoactive

substances and their metabolites, but to repeat the test with a different test container or to repeat the urine collection.

The test results for psychoactive substances and their metabolites are determined by the presence or absence of a coloured line/lines in the result area. The result must be determined within 5 minutes of the start of the analysis.

Positive test results must be confirmed by another test method. This involves sending a test container with a urine sample to a toxicology laboratory to confirm the results.

Figure 2– Conducting an analysis

9. Interpretation of results

9.1. Results of tests to detect psychoactive substances

A positive result (figure 3): only one coloured strip appears in the control area (C). No visible colour strip in the test area (T).

Negative result (figure 4): one or more coloured stripes appear on the membrane. One strip appears in the control area (C) and another and/or a third appears in the test area (T).

Invalid result (figure 5): no control strip is displayed. The results of any test which does not yield a reference strip within the specified reading time must be rejected. Please review the procedure and repeat the test with a new test system. If the problem persists, stop using the test system immediately and contact the manufacturer's authorised representative in the RF.

Figure 3 - Positive result Figure 4 - Negative result Figure 5 - Invalid result

NOTE:

The colour intensity in the test area (T) may vary depending on the concentration of the analytes present in the sample. Any shade of colour in the test area should therefore be considered negative.

Insufficient sample volume, non-compliance with methodology or expired test system are the most likely causes of a missing test strip.

9.2. Results of a test to establish adulteration

To interpret the results, use the colour chart (Figure 6) included in the test system.

NOTE: The adulterant detection test strips (urine) are designed to identify samples with abnormalities in basic parameters: specific gravity, creatinine, pH, oxidants.

Temperature: The temperature of the sample is determined by the thermal strip on the container. The appearance of a green marker below a specific temperature number (between 32 and 38) indicates this temperature value. The temperature of the sample is assessed no later than 2 minutes after the start of analysis.

Specific gravity: Elevated levels of protein in urine can lead to abnormally high specific gravity values.

Creatinine: The normal level of creatinine is between 20 and 350 mg/dl. In rare cases, in certain kidney diseases, its amount may be reduced.

pH: Normal pH levels should be between 4.0 and 9.0.

Oxidants: normal human urine should not contain oxidants such as bleach, hydrogen peroxide or pyridinium chlorochromate. The presence of high levels of antioxidants, such as ascorbic acid, in a sample can lead to false-negative results.

Indicators	Normal Range				Deviation from the Rated Value			
Oxidants								
Creatinine	 100 mg/dl		 20 mg/dl		 0 mg/dl		 10 mg/dl	
pH	 4,0	 7,0	 9,0	 2,0	 3,0	 10,0	 11,0	 12,0
Specific Density	 1,003	 1,005	 1,015	 1,025	 1,000		 1,035	

Figure 6– Colour table for interpreting adulteration test results (oxidants, creatinine, pH, specific gravity)

10. Quality Control

The test system provides internal control of the assay procedure with a different antigen/antibody response in the control area (C) on each strip. This control line should always appear regardless of the presence or absence of the psychoactive substance or its metabolite in the sample. If the control line does not appear, the test system should be discarded. The presence of this control strip in the control area serves to confirm sufficient sample volume, the correct technique for the analysis procedures and to verify that optimum chromatographic conditions have been achieved.

It is possible to use external control materials to ensure proper operation of the test system. Samples of test materials are commercially available. When performing quality control by means of external control material, the analysis procedure of point 8 of these instructions for use shall be followed. If the reference

values do not fall within the set limits, then the analysis results are invalid. A list of checklists can be found in Appendix A.

11. Limitations of the test system

The test system is designed for in vitro laboratory diagnosis and should only be used for the qualitative determination of psychoactive substances.

The test system is designed to perform the test using human urine only.

The test system is intended for single use only.

A positive result on either test only indicates the presence of a psychoactive substance/metabolite and does not establish the degree of intoxication.

There is a possibility that a technical or procedural error, associated substances, or factors not listed may interfere with the test and lead to false results.

If a psychoactive substance/metabolite is detected in a urine sample, the test will not show the frequency of use of these substances, the effect of food and medication intake on the test result.

12. Supply package

Test container for the detection of psychoactive substances and their metabolites in urine by immunochromatographic method for in vitro diagnosis ("NARCOLAB test system"), in versions:

Version 1.

Test container for the determination of 10 types of psychoactive substances and their metabolites with adulterated urine test strip for creatinine, pH, specific gravity and a thermal strip for sample temperature control.

The composition includes:

- Disposable container with key and test strip for the determination of 10 types of psychoactive substances and their metabolites: amphetamine, barbiturates, benzodiazepines, cocaine, methamphetamine, methadone, opiates, phencyclidine, synthetic cathinones, marijuana with adulteration control strip for creatinine, pH, specific gravity and thermal strip for sample temperature control - 25 pieces
- Tamper-evident tape - 25 pieces.

- Colour chart for determining adulteration - 1 pc.
- Instructions for use - 1pc.

Version 2.

Test container for the determination of 10 types of psychoactive substances and their metabolites with adulterated urine test strip for creatinine, pH, oxidants and a thermal strip for sample temperature control

The composition includes:

- Disposable container with key and test strip for the determination of 10 types of psychoactive substances and their metabolites: amphetamine, barbiturates, benzodiazepines, cocaine, methamphetamine, methadone, opiates, oxycodone, synthetic cannabinoids (K3), marijuana with adulteration control strip for creatinine, pH, oxidants and thermal strip for sample temperature control - 25 pieces
- Tamper-evident tape - 25 pieces.
- Colour chart for determining adulteration - 1 pc.
- Instructions for use - 1pc.

Version 3.

Test container for the determination of 12 types of psychoactive substances and their metabolites with adulterated urine test strip for creatinine, pH, specific gravity and a thermal strip for sample temperature control.

The composition includes:

- Disposable container with key and test strip for the determination of 12 types of psychoactive substances and their metabolites: amphetamine, barbiturates, benzodiazepines, cocaine, methamphetamine, methadone, opiates, phencyclidine, synthetic cathinones, marijuana, ethylglucuronide, methylenedioxymethamphetamine with adulteration control strip for creatinine, pH, specific gravity and thermal strip for sample temperature control - 25 pieces

- Tamper-evident tape - 25 pieces.
- Colour chart for determining adulteration - 1 pc.
- Instructions for use - 1 pc.

13. Conditions of use, storage and transport

13.1. Operating conditions

After opening the package, the test system should be used within 1 hour at a room temperature of 18°C to 25°C and a relative humidity of no more than 80%.

The test system, if stored between 2 and 8°C, must be kept at room temperature between 18 and 25°C for at least 10 minutes before analysis.

In order to obtain reliable results, the Instructions for Use must be strictly adhered to.

13.2. Storage conditions

The test system should be stored in the manufacturer's packaging at 2 to 30°C and humidity not exceeding 60% until the expiry date stated on the packaging.

Do not freeze. The product is moisture-sensitive and must be used immediately after opening. Any test in an incorrectly sealed package should be discarded.

13.3. Transport conditions

Transport the test system by all means of covered transport in accordance with the transport regulations applicable to the particular means of transport. Transport at a temperature between 2 °C and 30 °C.

Handle with care and avoid excessive pressure on the packaging during transport and storage

Items that have not been transported in the correct temperature range must not be used.

Symbols used in marking are shown in Appendix B.

14. Shelf life and stability

The system has a shelf life of 24 months from the manufacture date.

Stability tests have shown that the test system is stable for 24 months if properly stored in sealed aluminium foil packaging with moisture absorbers.

Based on the stability of the MD after opening, a test is recommended within 1 hour of opening the package.

14.1. Real-time stability

A real-time stability analysis of the test containers was conducted using samples from three different batches (3 lots) for each version. These samples were placed in an incubator at a temperature set at 2-8°C and 30±3°C and a relative humidity of 60%. Stability test series were performed after 0, 3, 6, 9, 12, 15, 18, 21, 24 and 27 months. Test containers were analysed using urine samples containing no detectable analyte as well as those containing psychoactive substances at a concentration of ±50% of the threshold value. Ten repeated experiments with the sample were held on day 0 and three repeated experiments with the sample at other times. The results were reviewed after 5 minutes of analysis. The tests were conducted in accordance with the instructions for use for the medical device.

Stability analysis of the test systems for real-time urine adulteration parameters was performed using samples from three different batches (3 lots) for each version. These samples were placed in an incubator at a temperature set at 2-8°C and 30±3°C and a relative humidity of 60%. A series of stability tests were performed at 0, 3, 6, 9, 12, 15, 18, 21, 24 and 25 months. Each 25 samples were enclosed in one container with 7g of desiccant. The sampling frequency follows the ASSURE stability protocol. Five repeat sample experiments were held on day 0 and one repeat sample experiment at other times. The result was assessed according to a colour table (Fig. 6)

no later than 5 minutes from the start of the analysis.

Conclusion: The test system is stable at 2-30°C for 24 months.

14.2. Accelerated stability

An accelerated stability study of the NARCOLAB test system was conducted using samples from three different batches (three lots) for each version. The test containers were placed in an incubator at a temperature set at 45°C with a relative humidity of 60%. Stability test series were held after 0, 7, 14, 21, 28, 35, 42, 56, 70, and 77 days. The test systems were analysed using urine samples

containing no analyte of interest as well as those containing an analyte at a concentration of $\pm 50\%$ of the threshold value. Three repeated tests were carried out for each sample. The tests were performed in accordance with the instructions for use.

An accelerated stability study of the urine adulteration test systems was conducted using container samples from three different batches (3 lots) for each version. These samples were placed in a thermostat with a temperature of 45°C and a relative humidity of 60%. After 0, 7, 14, 21, 28, 35, 42, 56, 70 and 77 days a series of stability measurements were carried out on the samples. Measurements were made using standard solutions of creatinine, oxidants (sodium hypochlorite, hydrogen peroxide and pyridinium chlorochromate), solutions with known pH and specific gravity values.

Five strips of each sample were taken for the initial examination (day 0) and one strip was taken for the following days.

The data obtained are presented in detail in the manufacturer's documentation "Accelerated Stability Study Report".

Conclusion: The test system is stable at 45°C for 77 days. The data from the study were plotted on an Arrhenius graph. The product has been determined to have a minimum shelf life of 24 months from the manufacture date.

14.3. Stability after opening the packaging

The study was conducted on three batches of a medical device (3 lots) for each version. Sealed aluminium foil bags were opened and the test sets were stored at $10 \pm 2^{\circ}\text{C}$ and $30 \pm 2^{\circ}\text{C}$.

The study was performed at each temperature regime, in pairs, combining it with relative humidity of 15% and 80%. Evaluation was done after 0, 1, 2, 3 and 4 hours using urine containing no detectable analytes as well as those containing analytes at concentrations of $\pm 50\%$ of the threshold. At each point in time, testing was performed in 3 repetitions for each sample. The tests were carried out according to the instructions for use. The results were examined after 5 minutes.

A post-opening stability study for adulteration parameters (pH, creatinine,

specific gravity, oxidants) was carried out on samples of three batches (3 lots) for each version. The study was performed at each temperature regime, in pairs, combining it with relative humidity of 15% and 80%. Evaluation was done after 0, 1, 2, 3 and 4 hours using standard solutions.

At each point in time, testing was performed in 3 repetitions for each sample of each lot. The tests were carried out according to the instructions for use. The results were reviewed no later than 5 minutes after the start of the analysis.

The data obtained are detailed in the manufacturer's documentation "Post-opening stability report".

Conclusion: Acceptable test results were demonstrated at all time points for all test containers stored at 10 °C and 30 °C at 15 % and 80 % humidity for four hours after opening the package. It is recommended that the test be carried out within one hour of opening the package.

14.4. Stability during transport

A study of medical device stability during transport was conducted to assess the effects of extreme temperature and humidity on the real-time performance of the test containers using simulated transport conditions. The problem of transport preceded the real-time determination of shelf life. This study found that the transport conditions do not reduce the shelf life of the set.

Stability was assessed with a human biomaterial sample containing no detectable analyte and a sample containing a detectable analyte at a concentration of $\pm 50\%$ of the threshold value. A sufficient sample of each biomaterial sample variant (urine containing no analyte and urine containing a detectable analyte at a concentration of $\pm 50\%$ of the threshold) was prepared for the transport stability tests. Three batches of sets (3 lots) of each design were tested in the study. After simulating the transport conditions, the sets were divided into two groups: one group was stored at $2-8\text{ }^{\circ}\text{C}$, the other at $30\pm 2\text{ }^{\circ}\text{C}$. Four groups of stress tests under extreme temperatures, humidity or reduced atmospheric pressure were carried out to simulate harsh delivery conditions:

3 freeze-thaw cycles: 3 freeze-thaw cycles were performed and at the last

thaw the samples were placed at normal storage temperature until the end of the study.

2 days at 55°C: the test samples were placed in storage conditions at 55°C for 2 days and then transferred to a room at normal storage temperature until the end of the test.

Low/High Humidity: the test samples were placed at low (20%) and high (80%) humidity respectively for 2 days and then moved to normal storage temperature until the end of the test.

Reduced atmospheric pressure: The test samples were placed in a vacuum to 500 mTorr for 24 hours and then transferred to normal storage conditions until the end of the study.


The results were recorded at 0 and after 3, 6, 9, 12, 15, 18, 21, 24 and 27 months. At each time of recording the results, the test system was brought to room temperature and the stability was assessed. Stability assessment for each sample was carried out in three repetitions. The test procedure and interpretation of the results were conducted in accordance with the instructions for use. Under simulated transport conditions, three batches of samples of each design showed the expected results given in the manufacturer's documentation "Report on stability study during transport".

Conclusion: after exposure to stressful environmental conditions, stability data are available for real-time evaluation for 27 months. Acceptable medical product performance was demonstrated at all time points for sets stored at 2-8 °C or 30 °C. These results show that the product in all versions is stable for at least 24 months at 2-30 °C after transport.

15. Disposal and destruction of the test system

The test system, unused after opening or due to tampering with the individual packaging, or after use in an environment outside the clinical laboratory, the sealed container and packaging must be disposed of in a waste bin.

The test system, unused after opening or tampering with the individual

packaging, must be disposed of as Class A medical waste when professionally used in clinical laboratories.

For professional use in clinical laboratories, all test system components must be disposed of after analysis in accordance with the Sanitary Regulations 2.1.3684-21. The used test container is classified as class B - epidemiologically hazardous waste - and requires appropriate precautions for its disposal.

Packaging is classified as epidemiologically safe waste, approximating in composition to municipal solid waste.

16. Information on the developer and manufacturer of the medical device

Assure Tech. (Hangzhou) Co., Ltd.

Address: Building 4, No.1418-50, Moganshan Road, Gongshu District, Hangzhou, 310011 Zhejiang, P.R. China

phone: +86-571-81022698, fax: +86571-8886920, email: contact@diareagent.com

17. Manufacturer's authorised representative in the Russian Federation

Limited Liability Company "NARCOLAB", OOO "NARCOLAB", legal address: 129515, Moscow, Akademika Koroleva Street, 13, building 1, office 831, website: narkolab.ru, email info@narkolab.ru; phone: +7(499) 758-44-44, INN 9717028246.

18. Warranty obligations

The manufacturer guarantees the compliance of the test system with the specified requirements of the standards and manufacturer's documentation. The test system is stable throughout its shelf life if the conditions for transport, storage and use are observed.

The warranty period for the test system is 24 months from the manufacture date.

Complaints, feedback on the quality of the test system should be addressed to Limited Liability Company "NARCOLAB", OOO "NARCOLAB", legal address: 129515, Moscow, Akademika Koroleva Street, 13, building 1, office 831, website: narkolab.ru, email info@narkolab.ru; phone: +7(499) 758-44-44, INN

9717028246.

19. List of normative documents

- GOST R ISO 18113-1-2015 - Medical devices for in vitro diagnosis . Information provided by the manufacturer (labelling). Part 1. Terms, definitions and general requirements;
- GOST R ISO 18113-2015 - Medical devices for in vitro diagnosis. Information provided by the manufacturer (labelling). Part 2. In vitro diagnosis reagents for professional use;
- GOST R ISO 23640-2015 - Medical devices for in vitro diagnosis. Assessment of the stability of reagents for in vitro diagnosis;
- GOST R ISO 51352-2013 - Medical devices for in vitro diagnosis. Test methodology;
- GOST ISO 15223-1-2020 -- Medical devices. Symbols to be used with medical device labels, labelling, and information to be supplied. Part 1. General requirements;
- GOST R EN 14254-2010 (sections 4-8) - Disposable containers for collection of human samples (except blood).

20. Literature references



1. Baselt RC. Распределение токсических веществ в организме человека. 2-е издание Davis: Biomedical Publications; 1982.
2. Hawks RL, Chiang CN, eds. Urine tests for psychoactive substances. Rockville: Department of Health and Human Services, National Institute on Drug Abuse; 1986.
3. Substance Abuse and Mental Health Services Administration. Mandatory Guidelines for Federal Workplace Drug Testing Programs. 53 Federal Register; 1988.
4. McBay AJ. Drug-analysis technology--pitfalls and problems of drug testing. Clin Chem. 1987 Oct; 33 (11 Suppl): 33B-40B.

5. Gilman AG, Goodman LS, Gilman A, eds. Goodman and Gilman's The Pharmacological Basis of Therapeutics. 6th ed. New York: Macmillan; 1980.

Appendix A List of control materials used for external quality control of the test system procedure

- Control material "Amphetamine sulphate - analytical sample" 1 ampoule (1 ml), batch 011120, State Institution of Health "MNPC Narcology of the DZM";
- Control material "Ethyl- β -D-glucuronide" - 1 ampoule (1 ml), REF FN12091902, State Institution of Health "MNPC Narcology of the DZM";
- Control material "DOA MC I Control Set" - 4 vials (5 ml each), REF 61732FD01, State Institution of Health "MNPC Narcology of the DZM";
- Control material "DOA MC I Cal 1" - 1 vial (10 ml), REF 29216FD01, State Institution of Health "MNPC Narcology of the DZM";
- Control material "DOA MC I Cal 2" - 1 vial (10 ml), REF 29217FD01, State Institution of Health "MNPC Narcology of the DZM";
- Control material "DOA MC I Cal 3" - 1 vial (10 ml), REF 29227FD01, State Institution of Health "MNPC Narcology of the DZM";
- Control material "DOA MC I Cal 4" - 1 vial (10 ml), REF 29232FD01, State Institution of Health "MNPC Narcology of the DZM";
- Control material Cannabinoids 200 Calibrator - 1 vial (5 ml), REF 58774FD01, State Institution of Health "MNPC Narcology of the DZM";
- Control material "Oxycodone hydrochloride"
- Control material "K3 AB-PINACA"

Appendix B: Symbols used in labelling

	Manufacturer		Refer to the instructions for use
	Use until		Careful!
	Batch code		Medical device for in vitro diagnosis
	Non-sterile		Enough content for n tests
	Do not use if the package is damaged		Compliance with Directive 98/79/EC
	Temperature limit		Keep out of moisture
	No re-use		Keep out of sunlight
	Up		Fragile, handle with care.
	Humidity range is within 60%		Manufacture date

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