Tokuyama Dental

38-9, Taitou 1-chome, Taitou-ku, Tokyo 110-0016 JAPAN

Phone: +81 3-5817-8665 Fax: +81 3-3835-2265

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Date: 20th of October 2021

«Approved by»

Hideki Kazama

President

Tokuyama Dental Corporation

Instruction for use for the medical device

Dental adhesive material Tokuyama Universal Bond II in kit

Version 1

1. Medical device name

Dental adhesive material Tokuyama Universal Bond II KIT:

1. BOND A 5 ml (bottle) - 1 pc.
2. BOND B 4 ml (bottle) - 1 pc.
3. MIXING WELL - 1 pc.
4. 1 –Dimple Disposable Mixing Well - 5 pcs.
5. DISPOSABLE APPLICATOR (FINE) - 25 pcs.
6. Instructions for Use - 1 pc.
7. Clinical guide - 1 pc.

Hereinafter referred to as material, device, Bond A, Bond B.

2. Information about the manufacturer, designer and authorized representative

Designer

Tokuyama Dental Corporation
38-9, Taitou 1-chome, Taitou-ku, Tokyo 110-0016, Japan
Tel.: +81-3-3835-7201
URL: <http://www.tokuyama-dental.com>

Manufacturer

Tokuyama Dental Corporation
38-9, Taitou 1-chome, Taitou-ku, Tokyo 110-0016, Japan
Tel.: +81-3-3835-7201
URL: <http://www.tokuyama-dental.com>

Manufacturing site

Tokuyama Dental Corporation Kashima Factory
26, Sunayama, Kamisu-shi, Ibaraki 314-0255, Japan

Authorized representative of the manufacturer in the Russian Federation

PROTECO JSC

Legal address: 192102, St. Petersburg, 57 Salova str., bldg. 1, lit. I.

Tel.: +7 (812) 635-88-90;

Email: info@protecodent.ru

3. Intended use

The product is intended for direct and indirect restorations, as well as for fixing composite materials such as polymer composite and adhesive polymer cement to the dental cavity formed by caries, to increase the adhesion strength of the polymerizable composite material with all materials used for indirect restoration.

4. Potential consumer, scope of application

Potential consumers of a medical device: dentists.

Scope: dentistry.

Instruction for use for the medical device:
Dental adhesive material Tokuyama Universal Bond II in kit

5. Indications, contraindications and possible side effects associated with the use of the medical device as directed

Indications

- Direct anterior and posterior restorations with light-curing, dual-curing, and self-curing composite materials;
- Intraoral repair of composite restorations, porcelain fused to metal, metal, and all-ceramic restorations without an additional primer;
- Cementation of indirect restorations and veneers when combined with light-cure, dual-cure, and self-curing resin cements;
- Bonding of core build-ups made of core build-up materials;
- Bonding of denture resin to metal base, clasp or attachment;
- Repair of denture with metal base, clasp or attachment;
- Bonding of opaque resin to a metal base in the fabrication of resin-faced crowns.

Contraindications:

This device contains methacrylic monomers, peroxide, organic solvents and acids. DO NOT use this device for patients allergic to or hypersensitive to methacrylic monomers, related monomers, peroxide, organic solvents and acids.

Possible side effects: allergic reactions to methacrylic and similar monomers, peroxide, organic solvents and acids or any other ingredients.

6. Device description, principle of operation of the medical device

The kit consists of two-component adhesive (primer) Tokuyama Universal Bond II (BOND A and BOND B), mixing well, 1-dimple disposable mixing wells and disposable applicators (fine).

Tokuyama Universal Bond II (primer) is a two-component (Bond A and Bond B) self-cured dental universal adhesive system for both direct and indirect restorations that can be used with self-etch, selective-enamel-etch and total-etch techniques. As a universal adhesive, Tokuyama Universal Bond II has been designed to be fully compatible with light-cured, self-cured and dual-cured composite materials. Tokuyama Universal Bond II improves the bond strength of the polymerizable resin material (adhesive resin cement, acrylic resin and composite resin) to indirect restorative materials such as glass ceramics (porcelain), oxide-ceramics (zirconia and alumina), metals (precious and non-precious) and resin materials including inorganic filler. The used mixture of Bond A and Bond B of Tokuyama Universal Bond II polymerizes, forming a "bonding layer" on the adhesion surface after air drying without light irradiation (by self-curing), after which this "bonding layer" is bonded to the polymerizable composite material. Tokuyama Universal Bond II bonds the entire adhesive surface, including the tooth structure (enamel and dentin), as well as materials used in indirect restorations such as glass ceramics, oxide ceramics, metals and composites without additional primer.

Other components included in the kit: disposable applicators (Fine), mixing well, 1-Dimple disposable mixing wells for dispensing, mixing and applying materials.

Mixing well is a plastic platform with contoured cells and rubber cup for material mixing, in which bottles with Tokuyama Universal Bond II, 1-Dimple disposable mixing well (5 pcs) are placed.

7. Precautions

- DO NOT use Tokuyama Universal Bond II for any purpose other than those listed in these instructions. Use Tokuyama Universal Bond II only as directed herein.
- Tokuyama Universal Bond II is designed for sale and use by licensed dental care professionals only. It is not designed for sale nor is it suitable for use by non-dental care professionals.
- DO NOT USE Tokuyama Universal Bond II if the safety seals are broken or appear to have been tampered with.
- If Tokuyama Universal Bond II causes an allergic reaction or oversensitivity, discontinue its use immediately.
- Use examination gloves (plastic, vinyl or latex) at all times when handling Tokuyama Universal Bond II to avoid the possibility of allergic reactions from methacrylic monomers. Note: Certain substances/materials may penetrate through examination gloves. If Tokuyama Universal Bond II comes in contact with the examination gloves, remove and dispose of the gloves, and wash hands thoroughly with water as soon as possible.
- Avoid contact of Tokuyama Universal Bond II with eyes, mucosal membrane, skin and clothing.
 - If Tokuyama Universal Bond II comes in contact with eyes, thoroughly flush eyes with water and immediately contact an ophthalmologist.
 - If Tokuyama Universal Bond II comes in contact with mucosal membrane, wipe the affected area immediately, and thoroughly flush with water after completing the restoration. Affected areas may whiten from protein coagulation; however, whitening should disappear within 24 hours. If whitening does not disappear within 24 hours, immediately contact a physician, and the patient should be so advised.
 - If Tokuyama Universal Bond II comes into contact with the skin or clothing, immediately saturate the area with an alcohol-soaked (not less than 70%) cotton swab or gauze and flush with water.
 - Instruct the patient to rinse his mouth immediately after treatment.
- Tokuyama Universal Bond II should not be ingested or aspirated. Ingestion or aspiration may cause serious injury.
- To avoid the unintentional ingestion of Tokuyama Universal Bond II, do not leave it unsupervised within the reach of patients and children.
- DO NOT expose Tokuyama Universal Bond II or its vapor to open flame because Tokuyama Universal Bond II is flammable.
- To avoid cross infection or lowering bond strength, DO NOT reuse the disposable applicator and the disposable mixing well. Clean the mixing well thoroughly with alcohol (not less than 70%) after each use.
- DO NOT use Tokuyama Universal Bond II directly to the cavity being in close proximity to the pulp. Pulp protection with glass ionomer lining or calcium hydroxide is recommended. DO NOT use EUGENOL BASED MATERIALS to protect the pulp as these materials will inhibit curing of Tokuyama Universal Bond II.
- Be aware that when a restoration is chipped from misaligned occlusion or bruxism (clenching, grinding or tapping), the repaired restoration may chip again.
- DO NOT use mixing wells supplied by other manufacturers.
- Use the mixing well or disposable mixing well provided when dispensing and mixing of Tokuyama Universal Bond II. If Tokuyama Universal Bond II is not dispensed into the rubber it may be difficult to remove the cured Tokuyama Universal Bond II from mixing well. The cured Tokuyama Universal Bond II can be removed by inverting the rubber.
- DO NOT MIX Tokuyama Universal Bond II with other brands of primers or adhesives.

Instruction for use for the medical device:
Dental adhesive material Tokuyama Universal Bond II in kit

- To avoid splashing the liquid, the bottle of Tokuyama Universal Bond II should not be opened immediately after taking out of the refrigerator. Tokuyama Universal Bond II should be brought to room temperature before using. After removing from the refrigerator, allow the product to sit for 20 minutes or until it reaches room temperature.
- Do not use Tokuyama Universal Bond II with resin materials for which special adhesives or primers are specified.
- Keep away from heat, direct sunlight, sparks and open flames.
- DO NOT USE after the indicated expiration date on the bottle/package.

PRECAUTIONS FOR MEDICATIONS AND MATERIALS

- Some materials and medications (hemostatic material) inhibit setting/adhesion of Tokuyama Universal Bond II.

DO NOT USE products that contain:

- eugenol,
 - aluminum chloride,
 - ferric sulfate,
 - calcium sulfate,
 - aluminum sulfate,
 - diamine Silver fluoride [molecular formula: $\text{Ag}(\text{NH}_3)_2\text{F}$].
- Tokuyama Universal Bond II does not adhere to tooth structure immediately after applications of materials that contain:
 - hydrogen peroxide (oxidol),
 - sodium hypochlorite.

We advise avoiding the use of these materials on the day of cementation, however, these materials can be used for root canal treatment as the inhibition caused by these materials usually evaporates within five days.

8. Application procedure

Clinical application for direct restorations and intraoral repair of restorations with composite resin

1. Cleaning

Thoroughly clean the tooth surface with a rubber cup and a fluoride-free paste then rinse with water.

2. Isolation

The rubber dam is the preferred method of isolation.

3. Cavity preparation

Prepare the cavity and rinse with water. Add bevels to the enamel margins of anterior preparations (class III, IV, V), as well as chamfers to the margins of posterior preparations (class I, II) because bevels and chamfers assist in erasing demarcations between the cavity margins and the restoration, thereby enhancing both esthetics and retention.

- If the adherent surface involves uncut enamel, apply etching agent to the uncut enamel. See section "Selective-enamel-etch technique".

4. Drying

Dry the cavity by using a blotting technique or an air syringe technique.

- DO NOT desiccate the vital tooth. Desiccation can lead to post-operative sensitivity.
- The substances listed below, which inhibit curing of Tokuyama Universal Bond II, should be removed from the tooth surface by thoroughly cleaning the tooth surface with alcohol (not less than 70%), citric acid, or the application of phosphoric acid for 2 to 3 seconds before application:
 - 1) Oil mist from handpiece,
 - 2) Saliva, blood and exudates.

Instruction for use for the medical device:

Dental adhesive material Tokuyama Universal Bond II in kit

5. Etching of the tooth structure (optional)

(a.) Selective-enamel-etch technique

Apply phosphoric acid gel only to the enamel surrounding the margin of prepared cavity and leave phosphoric acid gel in place for 10-15 seconds. Rinse the etched surface thoroughly (at least 15 seconds) with water, and then dry with mild air.

(b.) Total etch technique

Apply phosphoric acid gel onto the prepared enamel first, and then the dentin. The etchant should be left to react on the enamel and dentin for 10-15 seconds. Rinse thoroughly (at least 15 seconds) with water and dry with mild air or with cotton pellets; do not overdry.

6. Preparing the adhesive surface during repair

6.1. Ceramics other than porcelain, and composite materials

Bevel the margin if necessary. Roughen the area using air abrasion or grinding with a diamond bur using a slow speed to prepare the surface for adhesion. Thoroughly rinse and dry the surface.

6.2. Metal restorations

Roughen the area using air abrasion or grinding with a diamond bur using a slow speed to prepare the surface for adhesion. Thoroughly rinse and dry the surface.

6.3. Porcelain

Bevel the margin if necessary. DO NOT roughen the surface; Polish the area to prepare a clean and fresh surface for adhesion. Thoroughly rinse and dry the surface.

- When using air abrasion, use at least 30 to 50 micron grid of alumina powder.
- The substances listed below inhibit adhesion and must be removed from the adhesive surface by mechanical cleaning or alcohol soaking (not less than 70%).
 - 1) Oil mist from handpiece,
 - 2) Plaque, tartar, saliva, blood and exudate fluids.
- If necessary, protect the surface not being bonded with metal separator, paraffin film, Vaseline, etc.

7. Dispensing of Tokuyama Universal Bond II

Allow Tokuyama Universal Bond II to reach room temperature prior to use. Open the bottle by removing the cap and dispense one drop each of Tokuyama Universal Bond II A and B into the mixing well or disposable mixing well. Mix thoroughly with disposable applicator until the mixed bonding agent turns green.

- Dispense Tokuyama Universal Bond II B IMMEDIATELY after dispensing A. Tokuyama Universal Bond II A may turn pink due to solvent volatilization if allowed to sit too long. If Bond A turns pink, dispose of it and re-dispense fresh Tokuyama Universal Bond II A and B.
- Complete the application within 1 minute of dispensing when using the mixing well since Tokuyama Universal Bond II contains volatile solvents. After mixing, the color of Tokuyama Universal Bond II changes gradually based on the reaction of the catalyst of Tokuyama Universal Bond II, this is normal and does not indicate a problem with Tokuyama Universal Bond II.
- Complete the application within 3 minutes of dispensing when using disposable mixing well.
- Be sure to hold the bottles of Tokuyama Universal Bond II A and B vertically when dispensing. Do not hold the bottles diagonally or horizontally while dispensing; this could cause the Tokuyama Universal Bond II flow back into the nozzle and become contaminated.
- DO NOT touch the nozzle to the mixing well while dispensing. DO NOT interchange the bottle caps of Tokuyama Universal Bond II A and B. If Tokuyama Universal Bond II A and B come into contact while dispensing, Tokuyama Universal Bond II may cure in the nozzle.

- Wipe off excess Tokuyama Universal Bond II from the tip of the nozzle and the inside of the cap prior to closing.
- Immediately after dispensing, tightly replace the bottle caps of Tokuyama Universal Bond II A and B.
- DO NOT over tighten the bottle caps.

8. Application of Tokuyama Universal Bond II

Using the disposable applicator, apply the mixed Tokuyama Universal Bond II to the entire adherent surface(s). Take care not to omit any surfaces where the Tokuyama Universal Bond II should be applied.

- In the case of intraoral repair of restorations, if the adherent surface includes tooth structure and restorations (ceramics, porcelain, metals or composite materials), apply the mixed Tokuyama Universal Bond II to the entire adherent surface at once. No extra primers are needed.
- If the applied Tokuyama Universal Bond II pools on the cavity floor or a cavosurface angle and is too thick to air thin, remove the excess with a paper point or a new disposable applicator.
- Be careful not to allow excess Tokuyama Universal Bond II to flow into the subgingival. If Tokuyama Universal Bond II flows into the subgingival, thoroughly rinse and dry the area.
- If saliva, blood, or other fluids contaminate the applied Tokuyama Universal Bond II, thoroughly rinse the cavity with water, dry and re-apply fresh Tokuyama Universal Bond II.
- DO NOT RINSE the applied Tokuyama Universal Bond II with water except in unintentional contamination.
- DO NOT wait too long after the application of Tokuyama Universal Bond II. Air dry within 30 seconds after application, to ensure the appropriate film thickness.

9. Air Dry

Using an oil-free air/water syringe, apply weak air continuously to the Tokuyama Universal Bond II surface until the runny Tokuyama Universal Bond II stays in the same position without any movement. Finish by using mild air to the surface. Use a vacuum aspirator to prevent spatter of the Tokuyama Universal Bond II.

- If saliva, blood, or other fluids contaminate the applied Tokuyama Universal Bond II, thoroughly rinse with water, dry and re-apply fresh Tokuyama Universal Bond II.
- DO NOT RINSE the applied Tokuyama Universal Bond II with water except in unintentional contamination.
- If accidental spattering occurs, it may cause the tissue to whiten or possible allergic reaction.

10. Composite resin

Restore with a light-cured, self-cured or dual-cured composite resin according to its manufacturer's instructions. The overfilled composite resin should be thoroughly finished and polished.

- Tokuyama Universal Bond II layer may appear fluorescent under the dental curing light irradiation.

Clinical application for cementation of indirect restorations

1. Preparing the adhesive surface of the restoration

1.1. Ceramics other than Porcelain, and Composite materials:

Roughen the interior of the restoration by sandblasting (0.1 to 0.2 MPa), air abrasion or grinding with a diamond bur using a slow speed to prepare the surface for adhesion. Thoroughly rinse and dry the surface.

1.2. Metal restorations

Indtruction for use for the medical device:

Dental adhesive material Tokuyama Universal Bond II in kit

Roughen the area by sandblasting (0.3 to 0.5 MPa), air abrasion or grinding with a diamond bur using a slow speed to prepare the surface for adhesion. Thoroughly rinse and dry the surface.

1.3. Porcelain

DO NOT roughen the surface. Polish the area to prepare a clean and fresh surface for adhesion. Thoroughly rinse and dry the surface.

- When using the sandblaster or air abrasion, use at least 30 to 50 micron grid of alumina powder.
- If unable to roughen or polish the adhesive surface, etch the area for adhesion with phosphoric acid then thoroughly rinse and dry.
- The substances listed below inhibit adhesion and must be removed from the adhesive surface by mechanical cleaning or alcohol (not less than 70%) soaking. Clean the restoration after adjustment is complete:
 - Gypsum and investment
 - Silicone oil from fit check material
 - Oil mist from headpiece
 - Plaque, tartar, saliva, blood and exudate fluids
 - Temporary restoration and temporary cement

If necessary, protect the surface not being bonded with metal separator, paraffin film, Vaseline, etc.

2. Dispensing of Tokuyama Universal Bond II

Allow Tokuyama Universal Bond II to reach room temperature prior to use. Open the bottle by removing the cap and dispense one drop each of Tokuyama Universal Bond II A and B into a mixing well or disposable mixing well. Mix thoroughly with a disposable applicator until the mixed bonding agent turns green.

- Dispense Tokuyama Universal Bond II B IMMEDIATELY after dispensing A. Tokuyama Universal Bond II A may turn pink due to solvent volatilization if allowed to sit too long. If Bond A turns pink, dispose of it and re-dispense fresh Tokuyama Universal Bond II A and B.
- Complete the application within 1 minute of dispensing when using the mixing well since Tokuyama Universal Bond II contains volatile solvents. After mixing, the color of Tokuyama Universal Bond II changes gradually based on the reaction of the catalyst of Tokuyama Universal Bond II, this is normal and does not indicate a problem with Tokuyama Universal Bond II.
- Complete the application within 3 minutes of dispensing when using the disposable mixing well.
- Be sure to hold the bottles of Tokuyama Universal Bond II A and B vertically when dispensing. Do not hold the bottles diagonally or horizontally while dispensing; this could cause the Tokuyama Universal Bond II flow back into the nozzle and become contaminated.
- DO NOT touch the nozzle to the mixing well while dispensing. DO NOT interchange the bottle caps of Tokuyama Universal Bond II A and B. If Tokuyama Universal Bond II A and B come into contact while dispensing, Tokuyama Universal Bond II may cure in the nozzle.
- Wipe off excess Tokuyama Universal Bond II from the tip of the nozzle and the inside of the cap prior to closing.
- Immediately after dispensing, tightly replace the bottle caps of Tokuyama Universal Bond II A and B.
- Do not over tighten the bottle caps.

3. Application of Tokuyama Universal Bond II

Instruction for use for the medical device:

Dental adhesive material Tokuyama Universal Bond II in kit

Using the disposable applicator, apply the mixed Tokuyama Universal Bond II on the the surface to be bonded.

- Ensure Tokuyama Universal Bond-II covers all surfaces where the cement to be applied.
- DO NOT over apply. One application of Tokuyama Universal Bond II is sufficient. Too thick of a layer may compromise the fit of the restoration.
- DO NOT allow saliva, blood, oil, water, etc. t come in contact with the surface to be bonded. If the trated surface is contaminated, thoroughly rinse with water, dry and re-apply fresh Tokuyama Universal Bond II.
- DO NOT wait too long after the application of Tokuyama Universal Bond II. Air dry within 30 seconds after application, to ensure the appropriate film thickness.

4. Air dry

Using an oil-free air/water syringe, apply mild air to the surface. Use a vacuum aspirator to prevent spatter of the Tokuyama Universal Bond II.

5. Cementing

Cement with a resin cement according to its manufacturer's instructions.

- DO NOT apply resin cement directly to the prepared tooth or the cavity treated with Tokuyama Universal Bond II to avoid improper seating of the restoration.
- Curing may be accelerated depending upon the resin cement used. To keep the appropriate film thickness, seat the restoration as soon as possible after the application of the resin cement (within 1 minute after dispensing), if not there is a possibility that the restoration will not seat properly with cementation.
- The Tokuyama Universal Bond II layer may appear fluorescent under the dental curing light irradiation.

Clinical application for bonding of acrylic resin to metal base, clasp or attachment

1. Preparing the adhesive surface

Roughen the area by sandblasting (0.1 to 0.2 MPa), or by grinding with a diamond bur using a slow speed to prepare the surface for adhesion. Thoroughly rinse and dry the surface.

2. Dispensing, Application and Air dry

Use the same procedures as described in Clinical Application for Cementation of Indirect Restorations.

3. Application, curing, trimming and finishing of acrylic resin

Proceed with the procedure according to the manufacturer's instructions for the use of polymerizable acrylic resin.

IMPORTANT NOTE: The manufacturer is not responsible for damage or injury caused by improper use of this product. It is the personal responsibility of the user to ensure that this product is suitable for an appropriate application before use.

Specifications are subject to change without notice. When the product specifications change, the instructions and precautions may also change.

9. Functional characteristics

Table 1

Test item	Specification	
Appearance	BOND A	Translucent liquid, yellowish - orangish
	BOND B	Transparent - translucent liquid, bluish – blue greenish
Working time	Minimum 90 seconds	
Bond strength	To the tooth:	
	To enamel: minimum 12 MPa	
	To dentine: minimum 10 MPa	
	To the prosthesis:	
	To metal: minimum 12 MPa	

Indtruction for use for the medical device:
Dental adhesive material Tokuyama Universal Bond II in kit

	To ceramics (porcelain, ceramics and composite resin): minimum 10 MPa
Coating property	Form uniform thickness film

10. Devices incorporating medicinal products or substances absorbed by or dispersed in human body

Not applicable. This device is not a device incorporating a substance considered to be medicinal product, and is not a device that is composed of substances that are absorbed by or locally dispersed in the human body.

11. Devices incorporating materials of biological origin

Not applicable. This device is not a device incorporating a material of biological origin.

12. List of standards that the medical device comply with

EN ISO 10993-1, EN ISO 10993-18, ISO 7405, EN ISO 13485, EN ISO 14971, EN 1041, EN 1641, EN ISO 15223-1, ISO/TS 11405

The latest versions of reference standards should be used: the organization should maintain and update the list of reference standards as a record of its own management system in compliance with ISO 13485.

13. Interpretation of symbols used in the labeling of consumer and transport packaging

Table 2. Symbols and Description.

Symbol	Description	Symbol	Description
	Manufacturer		Lot number
	Catalog number		Expiration date
	Possibility for the occurrence of skin allergic reaction		"Caution"
	Consult instructions for use		Highly flammable liquid and vapour

14. Requirements for the maintenance and repair of a medical device

Devices do not need maintenance and repair.

15. Sterilization, cleaning and disinfection

The medical device is supplied non-sterile and is not intended to be sterilized before use.

Measures during use (Miscellaneous goods)

Instruction for use for the medical device:
Dental adhesive material Tokuyama Universal Bond II in kit

Instruments that have high probability for causing infection are mixing well, a disposable mixing well, and a disposable applicator. Disposable mixing well and disposable applicator are disposable goods and are specified in the instructions for not being reused.

Mixing well is required to be cleaned after use with an antimicrobial liquid such as an alcohol (not less than 70%). Alcohol is effective for disinfection of bacteria and fungi.

16. Conditions of transportation, storage and operation

Transportation of devices packed in shipping containers should be carried out by all types of closed vehicles in compliance with the precautions indicated on the shipping container and in accordance with the rules for the carriage of goods in force for this type of transport.

Kit shipping conditions:

- Temperature - from 0 to +25 °C;
- Relative humidity: 10 – 100 %;
- Atmospheric pressure: 800 ~ 1,060 hPa.

Store in the refrigerator. Keep away from heat, direct sunlight, sparks and open flames.

Kit storage conditions:

- Temperature - from 0 to +25 °C;
- Relative humidity: 10 – 70%;
- Atmospheric pressure: 900 ~ 1,060 hPa.

Indoor conditions when using a medical device:

- Temperature from +18 °C to +30 °C.

17. Shelf life

The shelf life of Tokuyama Universal Bond II is 2 years at the temperature 0 °C to 25 °C.

18. Procedure for the disposal and destruction of a medical device

The medical device **Dental adhesive material Tokuyama Universal Bond II in kit** depending on the degree of their epidemiological, toxicological and radiation hazard, is classified as epidemiologically safe waste, similar in composition to solid domestic waste (SDW). However, it is recommended to dispose of them separately as industrial waste through specialized organizations licensed for this type of activity.

Unused Tokuyama Universal Bond II should be absorbed into a neutral material (gauze, cotton wool) for safe disposal of device residues.

The devices should be collected and handed over to a specially licensed medical waste disposal organization. Waste transportation from institutions to the place of disposal is carried out by specialized companies.

Collection of non-hazardous waste is carried out in reusable containers or disposable bags (with the exception of bags in the color accepted for hazardous waste). Disposable bags are placed inside reusable containers. Containers for collecting waste should be marked and installed on a special site (indoor). Reusable containers after emptying should be washed and disinfected. The procedure for washing and disinfecting reusable containers is determined in accordance with the waste management scheme in each specific organization.

Medical devices should be disposed of in accordance with local medical procedures and local, state, and federal regulations. The choice of methods for the safe disposal and destruction of products depends on the capacity and profile of the medical organization, the availability of facilities for the disinfection/disposal of waste, the method of neutralization/destruction of waste adopted in the administrative territory (incineration, removal to landfills, disposal).


Registration No. 561 of 2021

Certificate of Acknowledgment of Notary

On this 21st day of December 2021 , Tomoko MIFUNE personally appeared before me, KATO Tomohiro , a notary in and for Tokyo Legal Affairs Bureau , with satisfactory evidence of her identification and of her being the agent of Hideki KAZAMA, President of Tokuyama Dental Corporation, in accordance with the Japanese Notary Act, stated that the principal , said Hideki KAZAMA , had acknowledged to have affixed his signature and seal to the attached document.

In witness whereof , I set my hand and seal .

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