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
INSTRUCTION FOR USE

For the medical device

Prosthetic Components for Neobiotech Implants

1. Medical device name

Prosthetic components for Neobiotech implants,

1. IS System Cemented Abutment: ISAH3537P, ISAH3547P, ISAH3557P, ISAH3577P.
2. IS System Cemented Abutment with screw ISCS16S: ISAH3537, ISAH3547, ISAH3557, ISAH3577.
3. IS System Cemented Abutment non-hex type: ISAN3537P, ISAN3547P, ISAN3557P, ISAN3577P.
4. IS System Cemented Abutment non-hex type with screw ISCS16S: ISAN3537, ISAN3547, ISAN3557, ISAN3577.
5. IS System A Type Angled abutment: ISAHA1739P.
6. IS System A Type Angled abutment with screw ISCS16S: ISAHA1739.
7. IS System B Type Angled abutment: ISAHB1739P.
8. IS System B Type Angled abutment with screw ISCS16S: ISAHB1739.
9. IS System temporary abutment hex type: ISAHT300P.
10. IS System temporary abutment hex type with screw ISCS16S: ISAHT300.
11. IS System temporary abutment non-hex type: ISANT300P.
12. IS System temporary abutment non-hex type with screw ISCS16S: ISANT300.
13. IS System Abutment Screw: ISCS16S.

hereinafter referred to as the device, abutment.

2. Information on the developer, manufacturer and authorized representative

Developer

Neobiotech Co., Ltd.

#104,24 / #102, 103, 104-1, 104-2, 105, 106, 205, 212, 306, 312, 509, 510, 511, 608, 704, 708, 710 & 10F, 36, Digital-ro 27gil, Guro-gu, Seoul, 08381, Korea

Tel: +82 2 582 2885

E-mail: info@neobiotech.com

Manufacturer

Neobiotech Co., Ltd.

#104,24 / #102, 103, 104-1, 104-2, 105, 106, 205, 212, 306, 312, 509, 510, 511, 608, 704, 708, 710 & 10F, 36, Digital-ro 27gil, Guro-gu, Seoul, 08381, Korea

Tel: +82 2 582 2885

E-mail: info@neobiotech.com

Manufacturing sites

Neobiotech Co., Ltd.

1. #104,24 / #102, 103, 104-1, 104-2, 105, 106, 205, 212, 306, 312, 509, 510, 511, 608, 704, 708, 710 & 10F, 36, Digital-ro 27gil, Guro-gu, Seoul, 08381, Korea
2. 1F, 60, Gasan digital 2-ro, Geumcheon-gu, Seoul, 08589, Korea

Authorized representative of the Manufacturer in the Russian Federation

GLAV DENTAL Limited Liability Company

Legal address: 49 Vyatskaya Street, Bldg. 2, Floor 2, Facility I, Room 60, Moscow, 127015, Russia

INSTRUCTION FOR USE

For the medical device

Prosthetic Components for Neobiotech Implants

3. Intended purpose of the device with indication of the potential consumer, scope of application

3.1. Intended purpose of the device

Products are designed for connection (retention) of implants with prosthetic restorations (single crown, bridgework or overdenture).

3.2. Potential consumers and scope of application

Dental surgeons, maxillofacial surgeons and dental implant surgeons.

Scope of application: dental implantology, maxillofacial surgery.

4. Indications, contraindications and possible side effects associated with the use of the medical device as intended

Indications:

It is used in dentistry for restorations with the loss of one tooth, with partial or complete loss of anterior teeth or already existing edentulism in the upper and lower jaws by means of making a relevant bone bed for the jawbone.

In bones of any density and in the following clinical situations:

- The absence of one or more teeth as well as complete adentia;
- Post-traumatic and postoperative defects of jaw bones;
- With appropriate bone characteristics (density, width and height of the bone):
- Single prosthetic dentures at the intervals from incisors to premolars.
- Bridgework on the completely and partially edentulous jaw.
- Rigid attachment of prostheses by means of dental bars and systems of ball-type fixation.

Contraindications:

- Undergoing radiation therapy to treat a malignant tumor;
- Patients with problems with jaw alignment and occlusion;
- Patients with oral infection or inflammation (inappropriate oral hygiene, bruxism, xerostomia);
- Pregnant, smoking patients, patients with mental disorders, for example, alcohol abusers, drug abusers, patients with diseases of the central and peripheral nervous system;
- Diseases of the blood and blood-forming organs or decompensated diseases of the cardiovascular system;
- Allergy to device components;
- General periodontal and surgical grounds for refusing any intervention.

Age restrictions exist only for young people under 18 years of age. This is due to the fact that they have not yet completed the growth of bone tissue and the formation of a bite.

Possible side effects:

- Complications such as malocclusion, paresthesia caused by neurological damage, infection, injury or damage to other teeth or blood vessels, edema, hypodermal bleeding, discomfort and pain, soft tissue ulcer can arise.
- Loosening of the implant screw structure, occurrence of a fracture in the implant or prosthetic.
- Bone defect around the implant and mucositis.

5. Precautions

- Consider patient condition with bone disease (osteoporosis, osteomalacia) or metabolic bone disorder in advance.
- Do not reuse devices.
- Dispose of opened devices, even if they have not been used.
- Dispose of the packaging of the used product.
- The implantation procedure can only be performed by a dentist who has undergone special training and practice. The decision on the possibility or impossibility of performing the dental implantation operation is made only by the doctor.
- The doctor must fully master the technique to be applied before the implantation. In addition, he must choose the most suitable type of implants. The wrong approach to the patient and the wrong application of methods can be the reasons for failed dental implantation or loss of bone strength around the implants.
- Patients with excessive tooth bite force need to have an implant of sufficient thickness and length to prevent fracture of the implant.
- MR compatibility statement: This device has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration, or image artifact in the MR environment. The safety of the device in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.
- This product is a non-sterile medical device. It must be sterilized before use.
- The devices are recommended to use in association with Neobiotech implants: IS-III active Fixture S-Narrow and Instruments Neobiotech for dental implantology, manufactured by "Neobiotech Co., Ltd.", Korea, registered in the appropriate manner in the Russian Federation.

6. Application procedure

1. Check the state of osseointegration between bone and dental implant.
2. Select the appropriate abutment size and angle to take into account the degree of retention and anatomical conditions.
3. Select a suitable screwdriver according to the shape of the screw.
4. Before use, sterilize medical devices with high-pressure steam at 132 °C for at least 4 minutes, drying time: 30 minutes.
5. The sterile device should be placed in the surgical tray to avoid any contact with foreign matter.
6. Use a screwdriver of an appropriate size and check the connection between abutment and screwdriver.
7. Install the device on the implant, observing the recommended tightening torque of 30Ncm. A low speed fixation is recommended for the handpiece.
8. Mark the gingival and vertical area to be removed.

9. Shape the gingiva and vertical space with a drill with sufficient irrigation.
10. Check for faults with another screwdriver.
11. Further orthopedic work.

7. Description of the device, principle of operation and characteristics of the medical device

7.1. General description

The devices are the upper parts of dental implant systems, made from various raw materials. Abutment, screws, etc. are used as a base or a support for restoration after the implant substructure has been placed. Cementing or screw retention is used as a technique for attaching the abutment and prosthesis retainers. The design of the abutment interface elements corresponds to the implants for the exact fit and the easy "implant-abutment" connection.

7.2. Connection.

IS Connection

Hex Type

Non-hex type

- ① The prosthesis is cemented to a cross-sectional abutment to prevent rotation of the prosthesis. Has a convex and concave shape to expand the adhesive surface.
- ② Hex type: This type has a hexagonal design that can be connected to the implant.
Non-hex type: This type has a non-hexagonal design that can be connected to the implant.

8. General dimensions*

- 8.1. IS System Cemented Abutment: ISAH3537P, ISAH3547P, ISAH3557P, ISAH3577P;
IS System Cemented Abutment with screw ISCS16S: ISAH3537, ISAH3547, ISAH3557, ISAH3577.

Table 1. Cemented abutment IS System: dimensions

No.	Item No.	Dimensions, mm		
		Ø max	Length	Hex
1.	ISAH3537P	3.5	12.1	2.08
2.	ISAH3547P	3.5	13.1	2.08
3.	ISAH3557P	3.5	14.1	2.08
4.	ISAH3577P	3.5	16.1	2.08
5.	ISAH3537	3.5	12.1	2.08
6.	ISAH3547	3.5	13.1	2.08
7.	ISAH3557	3.5	14.1	2.08
8.	ISAH3577	3.5	16.1	2.08

* - Tolerance max. ± 1 %


8.2.IS System Cemented Abutment non-hex type: ISAN3537P, ISAN3547P, ISAN3557P, ISAN3577P.

IS System Cemented Abutment non-hex type with screw ISCS1 6S: ISAN3537, ISAN3547, ISAN3557, ISAN3577.

Table 2. IS System Cemented Abutment non-hex type: dimensions

No.	Item No.	Dimensions, mm	
		Ø max	Length
1.	ISAN3537P	3.5	11.5
2.	ISAN3547P	3.5	12.5
3.	ISAN3557P	3.5	13.5
4.	ISAN3577P	3.5	15.5
5.	ISAN3537	3.5	11.5
6.	ISAN3547	3.5	12.5
7.	ISAN3557	3.5	13.5
8.	ISAN3577	3.5	15.5

* - Tolerance max. ± 1 %

8.3.IS System A Type Angled abutment: ISAHA1739P.

IS System A Type Angled abutment with screw ISCS16S: ISAHA1739.

Table 3. IS System A Type Angled abutment: dimensions

No.	Item No.	Dimensions, mm				
		Ø max	Length	Length of a crown part	Hex	Angle (°)
1	ISAHA1739P	3.2	14.1	9.0	2.08	17
2	ISAHA1739	3.2	14.1	9.0	2.08	17

* - Tolerance max. ± 1 %

8.4.IS System B Type Angled abutment: ISAHB1739P.
IS System B Type Angled abutment with screw ISCS16S: ISAHB1739.

Table 4. IS System B Type Angled abutment: dimensions

No.	Item No.	Dimensions, mm				
		Ø max	Length	Length of a crown part	Hex	Angle (°)
1	ISAHB1739P	3.2	14.1	9.0	2.08	17
2	ISAHB1739	3.2	14.1	9.0	2.08	17

* - Tolerance max. $\pm 1\%$

8.5.IS System temporary abutment hex type: ISAHT300P.
IS System temporary abutment hex type, with screw ISCS16S: ISAHT300.

Table 5. IS System Temporary abutment: dimensions

No.	Item No.	Dimensions, mm		
		Ø max	Length	Hex
1	ISAHT300P	3.0	14.1	2.08
2	ISAHT300	3.0	14.1	2.08

* - Tolerance max. $\pm 1\%$

8.6. IS System temporary abutment non-hex type: ISANT300P.
IS System temporary abutment non-hex type, with screw ISCS16S: ISANT300.

Table 6. IS System Temporary abutment: dimensions

No.	Item No.	Dimensions, mm	
		Ø max	Length
1	ISANT300P	3.0	13.5
2	ISANT300	3.0	13.5

* - Tolerance max. $\pm 1\%$

8.7.IS System Abutment Screw: ISCS16S.

Table 7. IS System Abutment Screw: dimensions

No.	Item No.	Dimensions, mm		
		Ø max	Length	Carving Ø
1	ISCS16S	2.0	10.2	1.55

* - Tolerance max. $\pm 1\%$

Table 8. Physical characteristics.

No.	Measure	Regulatory requirements
1	Compressive strength in shear	≥ 400 N
2	Fatigue strength	≥ 110 N


9. List of standards that the medical device complies with.

Directive 93/42/EEC as amended by 2007/47/EC; ISO 13485; EN ISO 14971; EN 1041; ISO 15223-1; ISO/IEC 17050-1; EN 1642 ISO 5832-3; ASTM F 136-13; ASTM F67-13; ASTM F1537; ISO14801; ISO 10993-1; ISO 10993-5; ISO 10993-10; MEDDEV 2.7.1 Rev.4; KS I ISO 14644-1; KS Q ISO2859-1.

Reference standards are used in the latest edition: the organization maintains and updates the list of reference standards as a record of its own management system in accordance with ISO 13485.

10. Interpretation of symbols used in the marking of consumer and transport package.

Table 9. Symbols and Interpretation.

Symbol	Interpretation	Symbol	Interpretation
	Manufacturer		Lot number
	Catalogue number (item number)		Best before
	Date of manufacture		Caution!
	Do not reuse		Do not use if the packaging is damaged
	Temperature limits		"Rx only" text – prescription only (in English)
	Non-sterile		Refer to the instructions for use
	Authorized representative in EU		Keep away from sunlight

11. Requirements for maintenance and repair of the medical device

The devices do not require maintenance and repair.

12. Sterilization, cleaning and disinfection

This product is a non-sterile medical device. It must be sterilized at 134 °C, sterilization time: min. 4 minutes, drying time: 30 minutes.

13. Transportation, storage and operation conditions

Devices packed into shipping containers should be transported by all types of covered vehicles considering the precautions indicated on the shipping container and in accordance with the guidelines for goods transportation by this mode of transportation.

Transportation conditions:

- Temperature: +1 ~ +30 °C;
- Relative humidity: no special requirements;
- Atmospheric pressure: no special requirements.

Storage conditions: in a dry place, away from direct sunlight:

- Temperature: +1 ~ +30 °C;
- Relative humidity: no special requirements;
- Atmospheric pressure: no special requirements.

Operation conditions:

- Operation temperature: +32 ~ +42 °C.

Indoor operation conditions:

- Temperature: +18 ~ +30 °C.

14. Shelf life

The shelf life of prosthetic components for Neobiotech implants is not limited.

15. Procedure for the disposal and destruction of a medical device

Medical devices "Prosthetic Components for Neobiotech Implants" after their intended use are classified as epidemiologically hazardous waste (products contaminated with blood and/or other bodily fluids of patients, as well as products contacting with patients having infectious diseases).

In this case, a medical organization performs only collection and temporary storage of the above waste. Waste is collected to single-use bags or containers with color marking adopted for hazardous medical waste in a room or in places for temporary storage of medical waste before further transportation to a destruction or disposal site.

Transportation of medical waste (used components) from medical facilities and their disposal is provided by specialized organizations licensed to treat epidemiologically hazardous medical waste.

The consumer packaging, rejected products and unused products are classified as epidemiologically safe waste (non-hazardous waste not contacting bodily fluids of patients, patients having infectious diseases), which is close by composition to municipal solid waste. However, it is recommended to have them disposed of separately as industrial waste by specialized

organizations licensed to carry out this kind of activity.

Non-hazardous waste is collected in reusable containers or single-use bags (with the exception of bags of a color used for hazardous waste). Single-use bags are placed inside reusable containers. Waste receptacles should be labeled and placed on the special site (indoors). After emptying, reusable containers should be cleaned and disinfected. The procedure for cleaning and disinfecting reusable containers is determined in accordance with the waste management plan in each specific organization.

Transportation and disposal are carried out with regard to requirements adopted at the given territory, as per the sanitary legislation of a particular country applied to the maintenance of populated areas and treatment of production and consumption waste.

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appeared before me and admitted
said principal`s subscription to
the attached INSTRUCTION FOR USE
For the medical device. _____

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This is hereby attested on this
17th day of Nov. 2022 at this office.

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Belong to Seoul Southern
District Prosecutor's Office

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Address of the office
E & C Venture Dream Tower 3rd Rm#509, 38-21,
Digital-ro 31-gil, Guro-gu, Seoul, Korea

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Jeong Kyeong Tae
Signature of the Notary Public
Jeong Kyeong Tae

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This office has been authorized by the
Minister of Justice, the Republic of
Korea, to act as Notary Public Since
7, Feb. 2020 Under Law No.15150.

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