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"APPROVED"

Zambelli Flavio Augusto

CEO STERYLAB S.r.l., Italy

"23" March 2023

stamp

STERYLAB S.r.l.

Via Magenta, 77/6
20017 RHO (MI) Italy
Part. IVA 10124260153
Tel. 02 93508427 r.a.
Fax 02 93508420

INSTRUCTIONS FOR USE

for a medical device

Huber needle PORT-JECT disposable, sterile

production: STERYLAB S.r.l., Italy

MDD 93/42/EEC Medical Devices Directive (MD Directive, MDD)

ISO 14971:2016 Medical devices. Application of risk management to medical devices

NS-EN ISO 13485:2012/AC:2012 Medical devices. Quality management systems. Requirements for regulatory purposes.

ISO 11135-1:2014 Sterilization of medical devices. Ethylene oxide. Part 1. Requirements for the development, validation and routine control of the sterilization process of medical devices

EN 556-1:2002 Sterilization of medical devices. Requirements for a medical device to have the designation "STERILE". Part 1: Requirements for terminally sterilized medical devices

EN 980-2008 Symbols for use in the labeling of medical devices

EN ISO 1041 + A1:2013 Information provided by the medical device manufacturer

ISO 15223-1:2016 Medical devices. Symbols to be used with medical device labels, labeling and information to be provided. Part 1. General requirements

ISO 11607-1:2019 Packaging for terminally sterilized medical devices. Part 1: Requirements for materials, sterile barrier systems and packaging systems

ISO 11607-2:2019 Packaging for terminally sterilized medical devices. Part 2: Validation requirements for forming, sealing and assembly processes

ISO 14644-1:2015 Clean rooms and associated controlled environments. Part 1. Classification of air purity by particle concentration

ISO 14644-2:2015 Clean rooms and associated controlled environments. Part 2: Monitoring to confirm the effectiveness of indoor cleaning related to clean air by particle concentration

ISO 11737-1:2018 Sterilization of medical devices. microbiological methods. Part 1. Evaluation of the population of microorganisms on products

ISO 11737-2:2009 Sterilization of medical devices. microbiological methods. Part 2: Tests of sterility to be performed in the definition, validation and maintenance of the sterilization process

ISO 10993-1:2018 Biological evaluation of medical devices. Part 1. Evaluation and testing in the risk management process

ISO 10993-5:2009 Biological evaluation of medical devices. Part 5: In vitro cytotoxicity tests

ISO 10993-7:2016 Biological evaluation of medical devices. Part 7: Ethylene Oxide Sterilization Residues

ISO 10993-10:2010 Biological evaluation of medical devices. Part 10: Skin irritation and sensitization tests

ASTM F1980-2016 Standard Guide for Accelerated Aging of Sterile Medical Device Packaging

MEDDEV 2.7/1 revision 4 (2016) Clinical evaluation: guidance for manufacturers and notified bodies in accordance with directives 93/42/EEC and 90/385/EEC

EN 62366-1:2015 Medical devices. Part 1: Designing medical devices for serviceability.

INTRODUCTION

Huber Needle PORT-JECT, versions: PORT-JECT A, PORT-JECT B, PORT-JECT C, PORT-JECT D, PORT-JECT E, PORT-JECT F has been designed for safe and easy insertion into vascular access devices and can be used for continuous and/or bolus infusions of pain medications, chemotherapy, antibiotic therapy, or parenteral nutrition.

All needles according to their device can be divided into two groups. The first includes straight and curved needles, which consist of a needle and a colorless/colored cap-connector. They can be connected directly to any standard infusion set as long as they are equipped with a Luer connector or extension lines.


The second group includes products with a more complex design. This is a curved Huber needle with wings ("butterfly") with an extension, which includes the needle itself, the wings that serve to fix it on the skin, in some variations there is a Y-connector for additional administration of drugs, and connecting tubes. Suitable for all implantable port systems.

The first group of needles is used for short and medium-term infusion, and the second - for long-term (up to 7 days). Butterfly wings are elastic and low-profile, which, on the one hand, facilitates the puncture, and on the other hand, does not interfere with the application of a fixing bandage over the needle (this is necessary for long-term infusions). The patient's skin is not injured, and the removal of the needle is quick and painless.

The straight Huber needle is suitable for bolus infusions, i.e. drugs weighing up to 0.5 g. and for flushing the port; curved Huber needle, suitable for infusions lasting up to 30-60 minutes in stationary conditions; the "butterfly" needle, thanks to a special clip and wings, is securely fixed, which allows you to administer drugs for up to 7 days.

NEEDLE POINT DESCRIPTION

Huber's needle is a hollow medical needle with a special sharpening of the cut and the shape of the tip. The tip of the needle is somewhat bent relative to its shaft so that the plane of the cut and the lumen of the needle are practically parallel to each other.

Huber needles are available in different lengths and different sizes - their caliber is measured on the Gage (G) scale. This scale is used to determine the outer diameter of tubular needles, which include Huber's needles.

The length and size of the needle must match the parameters of the implantable port system. If you choose a too short needle, then the introduction of the drug will be difficult due to the incorrect location of the needle tip in the silicone membrane, and in the case of blood sampling, it will be difficult to produce it with proper quality. If you take an excessively long needle, then it will not fit snugly against the skin, which will lead to discomfort during the infusion.

The unique non-coring design of the Huber needle tip promotes ease of penetration and maximum life of the Peart septum.

Available in commonly used needle gauges and lengths, standard kits include an extension line with or without a Y connector, a single acting clamp and a Luer connector.

1 PURPOSE AND SCOPE OF APPLICATION

Huber needle PORT-JECT disposable, sterile (hereinafter referred to as the needle, needles) is intended for chemotherapy, antibiotic therapy, long-term anesthesia, parenteral nutrition.

SCOPE OF APPLICATION – Medical institutions, medical offices, rehabilitation centers, clinical laboratories, medical vehicles.

2 LIST OF VERSIONS

Huber needle PORT -JECT disposable, sterile:

- ❖ 1. PORT-JECT A (Huber needle straight, Luer connector), size:
19G (1,10) x 15mm; 19G (1,10) x 20mm; 19G (1,10) x 25mm; 19G (1,10) x 32mm;
20G (0,90) x 15mm; 20G (0,90) x 20mm; 20G (0,90) x 25mm; 20G (0,90) x 32mm;
22G (0,70) x 15mm; 22G (0,70) x 20mm; 22G (0,70) x 25mm; 22G (0,70) x 32mm;
- 2. PORT-JECT B (Huber needle curved at right angle, Luer connector), size:
19G (1,10) x 15mm; 19G (1,10) x 20mm; 19G (1,10) x 25mm; 19G (1,10) x 32mm;
20G (0,90) x 15mm; 20G (0,90) x 20mm; 20G (0,90) x 25mm; 20G (0,90) x 32mm;
22G (0,70) x 15mm; 22G (0,70) x 20mm; 22G (0,70) x 25mm; 22G (0,70) x 32mm;
- 3. PORT-JECT C (Huber needle curved at a right angle with wings («butterfly»), assembled with an extension line, Luer connector, single-action clamp), size:
19G (1,10) x 15mm; 19G (1,10) x 20mm; 19G (1,10) x 25mm; 19G (1,10) x 32mm;
20G (0,90) x 15mm; 20G (0,90) x 20mm; 20G (0,90) x 25mm; 20G (0,90) x 32mm;
22G (0,70) x 15mm; 22G (0,70) x 20mm; 22G (0,70) x 25mm; 22G (0,70) x 32mm;
- 4. PORT-JECT D (Huber needle curved at a right angle with wings («butterfly»), assembled with an extension line, Luer connector, Y- connector for the 2nd infusion line, single-action clamp), size:
19G (1,10) x 15mm; 19G (1,10) x 20mm; 19G (1,10) x 25mm; 19G (1,10) x 32mm;
20G (0,90) x 15mm; 20G (0,90) x 20mm; 20G (0,90) x 25mm; 20G (0,90) x 32mm;
22G (0,70) x 15mm; 22G (0,70) x 20mm; 22G (0,70) x 25mm; 22G (0,70) x 32mm;
- 5. PORT-JECT E (Huber needle curved at a right angle with a gripper system (acupuncture), assembled with an extension line, Luer connector, single-action clamp), size:
19G (1,10) x 15mm; 19G (1,10) x 20mm; 19G (1,10) x 25mm; 19G (1,10) x 32mm;
20G (0,90) x 15mm; 20G (0,90) x 20mm; 20G (0,90) x 25mm; 20G (0,90) x 32mm;
22G (0,70) x 15mm; 22G (0,70) x 20mm; 22G (0,70) x 25mm; 22G (0,70) x 32mm;
- 6. PORT-JECT F (Huber needle curved at right angles with a gripper system (acupuncture), assembled with an extension line, Luer connector, Y-connector for the 2nd infusion line, single-action clamp), size:
19G (1,10) x 15mm; 19G (1,10) x 20mm; 19G (1,10) x 25mm; 19G (1,10) x 32mm;
20G (0,90) x 15mm; 20G (0,90) x 20mm; 20G (0,90) x 25mm; 20G (0,90) x 32mm;
22G (0,70) x 15mm; 22G (0,70) x 20mm; 22G (0,70) x 25mm; 22G (0,70) x 32mm;
- ❖ Sterile, single use: reduced risk of infection.

INDICATIONS FOR USE

The PORT-JECT Huber needle has been designed for safe and easy insertion into vascular access devices and can be used for continuous and/or bolus infusions of pain medication, chemotherapy, antibiotic therapy, or parenteral nutrition.

CONTRAINDICATIONS - none found

SIDE EFFECTS - Allergic reactions to the composition of materials are possible.

TYPE OF CONTACT – prolonged (from 24 hours to 30 days) contact with the bloodstream, indirectly.

3 SUPPLY

Huber needle PORT -JECT disposable, sterile:

1. PORT-JECT A (straight Huber needle, Luer connector), (50 pcs/pack), size:
19G (1,10) x 15mm; 19G (1,10) x 20mm; 19G (1,10) x 25mm; 19G (1,10) x 32mm;
20G (0,90) x 15mm; 20G (0,90) x 20mm; 20G (0,90) x 25mm; 20G (0,90) x 32mm;
22G (0,70) x 15mm; 22G (0,70) x 20mm; 22G (0,70) x 25mm; 22G (0,70) x 32mm;
- Instructions for use.
2. PORT-JECT B (right-angled Huber needle, Luer connector), (50 pcs/pack), size:
19G (1,10) x 15mm; 19G (1,10) x 20mm; 19G (1,10) x 25mm; 19G (1,10) x 32mm;
20G (0,90) x 15mm; 20G (0,90) x 20mm; 20G (0,90) x 25mm; 20G (0,90) x 32mm;
22G (0,70) x 15mm; 22G (0,70) x 20mm; 22G (0,70) x 25mm; 22G (0,70) x 32mm;
- Instructions for use.
3. PORT-JECT C (Huber needle bent at right angle, complete with wings ("butterfly"), extension line, Luer connector, single action clamp), (20pcs/pack), size:
19G (1,10) x 15mm; 19G (1,10) x 20mm; 19G (1,10) x 25mm; 19G (1,10) x 32mm;
20G (0,90) x 15mm; 20G (0,90) x 20mm; 20G (0,90) x 25mm; 20G (0,90) x 32mm;
22G (0,70) x 15mm; 22G (0,70) x 20mm; 22G (0,70) x 25mm; 22G (0,70) x 32mm;
- Instructions for use.
4. PORT-JECT D (right angled Huber needle, complete with wings ("butterfly"), extension line, Luer connector, Y-connector for the 2nd infusion line, single action clamp), (20pcs/pack), size:
19G (1,10) x 15mm; 19G (1,10) x 20mm; 19G (1,10) x 25mm; 19G (1,10) x 32mm;
20G (0,90) x 15mm; 20G (0,90) x 20mm; 20G (0,90) x 25mm; 20G (0,90) x 32mm;
22G (0,70) x 15mm; 22G (0,70) x 20mm; 22G (0,70) x 25mm; 22G (0,70) x 32mm;
- Instructions for use.
5. PORT-JECT E (right-angled Huber needle, complete with gripping system (needle clip), extension line, Luer connector, single action clamp), (20pcs/pack), size:
19G (1,10) x 15mm; 19G (1,10) x 20mm; 19G (1,10) x 25mm; 19G (1,10) x 32mm;
20G (0,90) x 15mm; 20G (0,90) x 20mm; 20G (0,90) x 25mm; 20G (0,90) x 32mm;
22G (0,70) x 15mm; 22G (0,70) x 20mm; 22G (0,70) x 25mm; 22G (0,70) x 32mm;
- Instructions for use.
6. PORT-JEC F (right angled Huber needle, complete with gripping system (needle clamp), extension line, Luer connector, Y-connector for the 2nd infusion line, single action clamp), (20pcs/pack), size:
19G (1,10) x 15mm; 19G (1,10) x 20mm; 19G (1,10) x 25mm; 19G (1,10) x 32mm;
20G (0,90) x 15mm; 20G (0,90) x 20mm; 20G (0,90) x 25mm; 20G (0,90) x 32mm;
22G (0,70) x 15mm; 22G (0,70) x 20mm; 22G (0,70) x 25mm; 22G (0,70) x 32mm;
- Instructions for use.

4 OPERATION INFORMATION

4.1 Needles must be used for the purposes specified in this technical document, in strict accordance with the instructions for use.

4.1.1 Huber's needles are straight and curved.

- ✓ Straight Huber needle is designed for the introduction of chemotherapy drugs in a short period of time and flushing of implanted port systems.
- ✓ The curved Huber needle is used for long-term administration of chemotherapeutic drugs, antibacterial agents, parenteral nutrition, bolus administration of drugs and for transfusion of components.
- ✓ Huber's needle, with a gripping system (needle clamp) or with wings ("butterfly"). The most convenient to use, it is perfect for treatments that last more than 24 hours, because, securely fixed in one position and reduces the risk of accidental change of its position, and therefore increase the safety of the procedure.

Curved are divided into needles for short-term infusions (30-60 min, not for ambulatory patients), which can be used both for washing port systems and for long-term infusions (up to 7 days) with wings ("butterfly"), with a gripping system (needle clip) and an extension cord.

Needles with an extension line come with a Y-connector, which allows you to additionally inject some drug in a jet without disconnecting the needle from the port, and without it. Needles are also divided by diameter and length.

Which needle to use is decided by the doctor, based on the preparation used and the complexion patient.

Needles with wings ("butterfly") are used up to 7 days.

4.1.2 Huber needles are for single use only.

4.1.3 Needles are not subject to reuse.

4.2 Instructions before use

- ✓ Put on sterile medical gloves;
- ✓ Wash hands with aseptic solution;
- ✓ Disinfect the skin area where the port is installed and wait until the disinfection site dries;
- ✓ Prepare a syringe with saline solution - a volume of at least 10 ml, since excessive pressure may be created in syringes of small volume;
- ✓ Use two fingers to fix the port system and carefully, so as not to bend the tip of the needle, insert it at a right angle to the titanium chamber located at the bottom of the device;
- ✓ After that, check the correct placement of the needle and the patency of the catheter by performing a reverse flow (reflux) of blood.
- ✓ Before the introduction of medicines, it is necessary, without fail, to aspirate (2 ml of blood into a syringe), and also inject Sodium Chloride 0.9% (5 ml) into the system;
- ✓ To avoid malfunction of the port system, the recommended pressure $2,549 \text{ kg/cm}^2 / 0,25 \text{ MPa}$ (36,26 psi/2,5 Bar) must not be exceeded.
- ✓ In accordance with the manufacturer's recommendations, contrast agents must be heated to 37C (98.6F) before use. Failure to follow these recommendations may result in a malfunction of the injection system or a decrease in consumption by 50%.
- ✓ For pediatric and shoulder port systems, it is recommended to use 20 or 22 G needles. The use of 19 G needles may lead to leakage of contrast medium.
- ✓ Before starting a high pressure injection, make sure that the needle is correctly inserted into the port and securely fixed.
- ✓ The use of catheters longer than 20 cm will result in reduced flow rates.
- ✓ For additional stabilization of the "butterfly" needle during infusion, a sterile swab can be placed under its wings;
- ✓ To fix the needle, it is recommended to use a patch that is glued to the wings of the "butterfly" and a breathable bandage over the needle and port.

4.2.1 Huber Port-Ject needle insertion procedure:

- 1) Wash hands with aseptic solution.
- 2) Locate the implantable port and disinfect the skin at that location.
- 3) For Port-Ject A and Port-Ject B: Attach an extension line with clamp (not included) to the Port-Ject Luer connector.
- 4) Close the clip.
- 5) Place the port in the triangle formed by the thumb and forefinger, and determine the center of this triangle, which will thus correspond to the center of the port itself.
- 6) Insert the Port-Ject needle perpendicularly until the needle touches the solid base of the port.
- 7) Open the clamp.
- 8) Check the correct placement of the needle by aspirating blood.
- 9) Introduce saline solution to flush the system.
- 10) In case of difficult aspiration or injection of solution, remove the needle and start from step 2).
- 11) Perform the activities for which the needle was inserted (infusion, withdrawal, heparinization)

4.2.2 Patency test after needle installation:

- 1) connect a 20 ml syringe with sterile saline solution to the free end of the needle. Then pull the plunger towards you until 1-2 ml of venous blood appears in the lumen of the syringe.

Important! If there is no reverse blood flow, then the needle is installed incorrectly and it should be replaced with a new one.

- 2) After that, slowly press the piston to flush the port. The introduction should take place unhindered and not cause any unpleasant sensations.
- 3) The port is ready for use. Connect the infusion system to the extension line from the needle and begin the procedure of administering medications or other medications.

4.2.3 Procedure for removing the Port-Ject Hub needle:

- 1) Wash your hands with aseptic solution.
- 2) Close the clamp and remove the infusion or sampling kit.
- 3) Attach a syringe with saline solution.
- 4) Open the clamp and rinse the port with saline solution.
- 5) Close the clamp and attach the syringe with the heparinized solution.
- 6) Inject the heparinized solution by closing the clamp or removing the Port-Ject needle, exerting constant pressure.
- 7) Remove the Port-Ject needle and disinfect the skin.

When working with the port system, only Goober needles are used, because ordinary needles pierce the membrane, while the Goober needle only pushes its surface without damage.

4.3 Features

The Huber needle is straight and curved, with wings ("butterfly")

→ It is not recommended to use on outpatient patients!

→ Directly connected to the infusion line or 3-way tap.


→ High-pressure contrast agent injections can be performed for CECT (Computed tomography with enhanced contrast) up to 8 ml/s (19G) at 2,549 kg/cm² (36,26 psi). inch.

The Huber needle is curved, with a gripping system (needle clamp), with an extension line and a Y-connector

- Equipped with a soft coating of polyethylene foam, which provides a high level of comfort for the patient's skin, and also provides good fixation.
- The needles have a unique device- "clutches", which provides ease of installation, facilitates port handling, and prevents the possibility of infection.

5 INSTRUCTIONS FOR USE

1) INSTRUCTIONS FOR USE HUBER NEEDLE PORT-JECT A (Huber needle straight, Luer connector)

WARNINGS:

- ☐ store in a dry, clean place, away from direct heat sources, at a temperature of +5 to +35°C
- ☐ before use, check the integrity of the sterile packaging: in case of damage or if the packaging is open, discard the product.
- ☐ check the expiration date and conformity of the product.
- ☐ open the package by pulling the special opening tab.
- ☐ be careful not to violate the sterility of the product before use.
- ☐ side effects: hematomas may occur after a biopsy.
- ☐ in case of defects or improper functioning, fill out the claim form (yellow form) placed in a cardboard box and save the defective product.
- ☐ after several punctures, it may be necessary to replace the port partition: check the integrity of the port after each puncture.

Note:

The product is intended for use by qualified medical personnel with knowledge of its application.

The product is sterile (sterilized with ethylene oxide) and it is intended for one-time use. After use, it must be disposed of. (SanPiN 2.1.3684-21 waste class B. Expired and unused products are disposed of as class A household waste).

Sterylab rejects any kind of liability arising from improper or repeated use of the product or its repeated sterilization.

In case of repeated use of the product, the following risks may arise: cross-infection, lack of sterility of the product, damage to the sharpening of the mandrel or cannula, damage to the connection between the plastic and metal parts of the product.


PROCEDURE FOR INSERTION OF HUBER NEEDLE PORT-JECT A:

- 1) Wash your hands with aseptic solution.
- 2) Determine the position of the implantable port and disinfect the skin in this place.
- 3) Connect the extension line with a clip (not included) to the Port-Ject Luer connector.
- 4) Close the clamp.
- 5) Place the port in the triangle formed by the thumb and index fingers, and determine the center of this triangle, which will correspond in this way to the center of the port itself.
- 6) Insert the Port-Ject needle perpendicular until the needle touches the solid base of the port.
- 7) Open the clip.
- 8) Check the correct location of the needle with blood aspiration.
- 9) Enter a saline solution to flush the system.
- 10) In case of difficult aspiration or injection of the solution, remove the needle and start from point 2).
- 11) Perform the actions for which the needle was inserted (infusion, sampling, heparinization)

PROCEDURE FOR REMOVING THE HUBER NEEDLE PORT-JECT A:

- 1) Wash your hands with aseptic solution.
- 2) Close the clamp and remove the infusion or intake kit.
- 3) Attach a syringe with saline solution.
- 4) Open the clamp and rinse the port with saline solution.
- 5) Close the clamp and attach the syringe with the heparinized solution.
- 6) Inject the heparinized solution by closing the clamp or removing the Port-Ject needle, exerting constant pressure.
- 7) Remove the Port-Ject needle and disinfect the skin.

2) INSTRUCTIONS FOR USE HUBER NEEDLE PORT-JECT B (Huber needle curved at right angles, Luer connector)

WARNINGS:

- ☐ store in a dry, clean place, away from direct heat sources, at a temperature of +5 to +35°C
- ☐ before use, check the integrity of the sterile packaging: in case of damage or if the packaging is open, discard the product.
- ☐ check the expiration date and conformity of the product.
- ☐ open the package by pulling the special opening tab.
- ☐ be careful not to violate the sterility of the product before use.
- ☐ side effects: hematomas may occur after a biopsy.
- ☐ in case of defects or improper functioning, fill out the claim form (yellow form) placed in a cardboard box and save the defective product.
- ☐ after several punctures, it may be necessary to replace the port partition: check the integrity of the port after each puncture.

Note:

The product is intended for use by qualified medical personnel with knowledge of its application.

The product is sterile (sterilized with ethylene oxide) and it is intended for one-time use. After use, it must be disposed of. (SanPiN 2.1.3684-21 waste class B. Expired and unused products are disposed of as class A household waste).

Sterylab rejects any kind of liability arising from improper or repeated use of the product or its repeated sterilization.

In case of repeated use of the product, the following risks may arise: cross-infection, lack of sterility of the product, damage to the sharpening of the mandrel or cannula, damage to the connection between the plastic and metal parts of the product.

PROCEDURE FOR INSERTING HUBER NEEDLE PORT-JECT B:

Wash your hands with aseptic solution.

Determine the position of the implantable port and disinfect the skin in this place.

Connect the extension line with a clip (not included) to the Port-Ject Luer connector.

Close the clamp.

Place the port in the triangle formed by the thumb and index fingers, and determine the center of this triangle, which will correspond in this way to the center of the port itself.

Insert the Port-Ject needle perpendicular until the needle touches the solid base of the port.

Open the clip.

Check the correct location of the needle with blood aspiration.

Enter a saline solution to flush the system.


In case of difficult aspiration or injection of the solution, remove the needle and start from point 2).

Perform the actions for which the needle was inserted (infusion, sampling, heparinization)

PROCEDURE FOR REMOVING HUBER NEEDLE PORT-JECT B:

- ✓ Wash your hands with aseptic solution.
- ✓ Close the clamp and remove the infusion or intake kit.
- ✓ Attach a syringe with saline solution.
- ✓ Open the clamp and rinse the port with saline solution.
- ✓ Close the clamp and attach the syringe with the heparinized solution.
- ✓ Inject the heparinized solution by closing the clamp or removing the Port-Ject needle, exerting constant pressure.
- ✓ Remove the Port-Ject needle and disinfect the skin.

3) INSTRUCTIONS FOR USE HUBER NEEDLE PORT-JECT C (Huber needle curved at right angles, assembled with wings ("butterfly"), extension line, Luer connector, single-action clamp)

WARNINGS:

- store in a dry, clean place, away from direct heat sources at a temperature of +5 to +35°C
- before use, check the integrity of the sterile packaging: in case of damage or if the packaging is open, discard the product.
- check the expiration date and conformity of the product.
- open the package by pulling the special opening tab.
- be careful not to violate the sterility of the product before use.
- side effects: hematomas may occur after a biopsy.
- in case of defects or improper functioning, fill out the claim form (yellow form) placed in a cardboard box and save the defective product.
- after several punctures, it may be necessary to replace the port partition: check the integrity of the port after each puncture.

Note:

The product is intended for use by qualified medical personnel with knowledge of its application.

The product is sterile (sterilized with ethylene oxide) and it is intended for one-time use. After use, it must be disposed of. (SanPiN 2.1.3684-21 waste class B. Expired and unused products are disposed of as class A household waste).

Sterylab rejects any kind of liability arising from improper or repeated use of the product or its repeated sterilization.

In case of repeated use of the product, the following risks may arise: cross-infection, lack of sterility of the product, damage to the sharpening of the mandrel or cannula, damage to the connection between the plastic and metal parts of the product.


PROCEDURE FOR INSERTING HUBER NEEDLE PORT-JECT C:

- 1) Wash your hands with aseptic solution.
- 2) Determine the position of the implantable port and disinfect the skin in this place.
- 3) Close the clamp.
- 4) Place the port in the triangle formed by the thumb and index fingers, and determine the center of this triangle, which will correspond in this way to the center of the port itself.
- 5) Insert the Port-Ject needle perpendicular until the needle touches the solid base of the port.
- 6) Open the clip.
- 7) Check the correct location of the needle with blood aspiration.
- 8) Enter a saline solution to flush the system.
- 9) In case of difficult aspiration or injection of the solution, remove the needle and start from point 2).
- 10) Perform the actions for which the needle was inserted (infusion, sampling, heparinization)

PROCEDURE FOR REMOVING HUBER NEEDLE PORT-JECT C:

- 1) Wash your hands with aseptic solution.
- 2) Close the clamp and remove the infusion or intake kit.
- 3) Attach a syringe with saline solution.
- 4) Open the clamp and rinse the port with saline solution.
- 5) Close the clamp and attach the syringe with the heparinized solution.
- 6) Inject the heparinized solution by closing the clamp or removing the Port-Ject needle, exerting constant pressure.
- 7) Remove the Port-Ject needle and disinfect the skin

4) INSTRUCTIONS FOR USE HUBER NEEDLE PORT-JECT D (Huber needle curved at right angles, assembled with wings ("butterfly"), extension line, Luer connector, Y-connector for the 2nd infusion line, single-action clamp)

WARNINGS:

- ☐ store in a dry, clean place, away from direct heat sources, at a temperature of +5 to +35°C.
- ☐ before use, check the integrity of the sterile packaging: in case of damage or if the packaging is open, discard the product.
- ☐ check the expiration date and conformity of the product.
- ☐ open the package by pulling the special opening tab.
- ☐ be careful not to violate the sterility of the product before use.
- ☐ side effects: hematomas may occur after a biopsy.
- ☐ in case of defects or improper functioning, fill out the claim form (yellow form) placed in a cardboard box and save the defective product.
- ☐ after several punctures, it may be necessary to replace the port partition: check the integrity of the port after each puncture.

Note:

The product is intended for use by qualified medical personnel with knowledge of its application.

The product is sterile (sterilized with ethylene oxide) and it is intended for one-time use. After use, it must be disposed of. (SanPiN 2.1.3684-21 waste class B. Expired and unused products are disposed of as class A household waste)

Sterylab rejects any kind of liability arising from improper or repeated use of the product or its repeated sterilization.

In case of repeated use of the product, the following risks may arise: cross-infection, lack of sterility of the product, damage to the sharpening of the mandrel or cannula, damage to the connection between the plastic and metal parts of the product.


PROCEDURE FOR INSERTION OF HUBER NEEDLE PORT-JECT D:

- 1) Wash your hands with aseptic solution.
- 2) Determine the position of the implantable port and disinfect the skin in this place.
- 3) Close the clamp.
- 4) Place the port in the triangle formed by the thumb and index fingers, and determine the center of this triangle, which will correspond in this way to the center of the port itself.
- 5) Insert the Port-Ject needle perpendicular until the needle touches the solid base of the port.
- 6) Open the clip.
- 7) Check the correct location of the needle with blood aspiration.
- 8) Enter a saline solution to flush the system.
- 9) In case of difficult aspiration or injection of the solution, remove the needle and start from point 2).
- 10) Perform the actions for which the needle was inserted (infusion, sampling, heparinization)

PROCEDURE FOR REMOVING THE HUBER NEEDLE PORT-JECT D:

- 1) Wash your hands with aseptic solution.
- 2) Close the clamp and remove the infusion or intake kit.
- 3) Attach a syringe with saline solution.
- 4) Open the clamp and rinse the port with saline solution.
- 5) Close the clamp and attach the syringe with the heparinized solution.
- 6) Inject the heparinized solution by closing the clamp or removing the Port-Ject needle, exerting constant pressure.
- 7) Remove the Port-Ject needle and disinfect the skin.

- 5) INSTRUCTIONS FOR USE HUBER NEEDLE PROJECT E (Huber needle curved at right angles, assembled with a gripping system (needle clamp), extension line, Luer connector, single-action clamp)

WARNINGS:

- ☐ store in a dry, clean place, away from direct heat sources at a temperature of +5 to +35°C.
- ☐ before use, check the integrity of the sterile packaging: in case of damage or if the packaging is open, discard the product.
- ☐ check the expiration date and conformity of the product.
- ☐ open the package by pulling the special opening tab.
- ☐ be careful not to violate the sterility of the product before use.
- ☐ side effects: hematomas may occur after a biopsy.
- ☐ in case of defects or improper functioning, fill out the claim form (yellow form) placed in a cardboard box and save the defective product.
- ☐ after several punctures, it may be necessary to replace the port partition: check the integrity of the port after each puncture.

Note:

The product is intended for use by qualified medical personnel with knowledge of its application.

The product is sterile (sterilized with ethylene oxide) and it is intended for one-time use. After use, it must be disposed of. (SanPiN 2.1.3684-21 waste class B. Expired and unused products are disposed of as class A household waste).

Sterylab rejects any kind of liability arising from improper or repeated use of the product or its repeated sterilization.

In case of repeated use of the product, the following risks may arise: cross-infection, lack of sterility of the product, damage to the sharpening of the mandrel or cannula, damage to the connection between the plastic and metal parts of the product.


PROCEDURE FOR INSERTION OF HUBER NEEDLE PORT-JECT E:

- 1) Wash your hands with aseptic solution.
- 2) Determine the position of the implantable port and disinfect the skin in this place.
- 3) Close the clamp.
- 4) Place the port in the triangle formed by the thumb and index fingers, and determine the center of this triangle, which will correspond in this way to the center of the port itself.
- 5) Insert the Port-Ject needle perpendicular until the needle touches the solid base of the port.
- 6) Open the clip.
- 7) Check the correct location of the needle with blood aspiration.
- 8) Enter a saline solution to flush the system.
- 9) In case of difficult aspiration or injection of the solution, remove the needle and start from point 2).
- 10) Perform the actions for which the needle was inserted (infusion, sampling, heparinization)

PROCEDURE FOR REMOVING HUBER NEEDLE PORT-JECT E:

- 1) Wash your hands with aseptic solution.
- 2) Close the clamp and remove the infusion or intake kit.
- 3) Attach a syringe with saline solution.
- 4) Open the clamp and rinse the port with saline solution.
- 5) Close the clamp and attach the syringe with the heparinized solution.
- 6) Inject the heparinized solution by closing the clamp or removing the Port-Ject needle, exerting constant pressure.
- 7) Remove the Port-Ject needle and disinfect the skin.

- 6) INSTRUCTIONS FOR USE HUBER NEEDLE PORT-JECT F (Huber needle curved at right angles, assembled with a gripping system (needle clamp), extension line, Luer connector, Y-connector for the 2nd infusion line, single-action clamp)

WARNINGS:

- ☐ store in a dry, clean place, away from direct heat sources, at a temperature of +5 to +35°C.
- ☐ before use, check the integrity of the sterile packaging: in case of damage or if the packaging is open, discard the product.
- ☐ check the expiration date and conformity of the product.
- ☐ open the package by pulling the special opening tab.
- ☐ be careful not to violate the sterility of the product before use.
- ☐ side effects: hematomas may occur after a biopsy.
- ☐ in case of defects or improper functioning, fill out the claim form (yellow form) placed in a cardboard box and save the defective product.
- ☐ after several punctures, it may be necessary to replace the port partition: check the integrity of the port after each puncture.

Note:

The product is intended for use by qualified medical personnel with knowledge of its application.

The product is sterile (sterilized with ethylene oxide) and it is intended for one-time use. After use, it must be disposed of. (SanPiN 2.1.3684-21 waste class B. Expired and unused products are disposed of as class A household waste).

Sterylab rejects any kind of liability arising from improper or repeated use of the product or its repeated sterilization.

In case of repeated use of the product, the following risks may arise: cross-infection, lack of sterility of the product, damage to the sharpening of the mandrel or cannula, damage to the connection between the plastic and metal parts of the product.

PROCEDURE FOR INSERTION OF HUBER NEEDLE PORT-JECT F:

- 1) Wash your hands with aseptic solution.
- 2) Determine the position of the implantable port and disinfect the skin in this place.
- 3) Close the clamp.
- 4) Place the port in the triangle formed by the thumb and index fingers, and determine the center of this triangle, which will correspond in this way to the center of the port itself.
- 5) Insert the Port-Ject needle perpendicular until the needle touches the solid base of the port.
- 6) Open the clip.
- 7) Check the correct location of the needle with blood aspiration.
- 8) Enter a saline solution to flush the system.
- 9) In case of difficult aspiration or injection of the solution, remove the needle and start from point 2).
- 10) Perform the actions for which the needle was inserted (infusion, sampling, heparinization)

PROCEDURE FOR REMOVING THE HUBER NEEDLE PORT-JECT F:

- 1) Wash your hands with aseptic solution.
- 2) Close the clamp and remove the infusion or intake kit.
- 3) Attach a syringe with saline solution.
- 4) Open the clamp and rinse the port with saline solution.
- 5) Close the clamp and attach the syringe with the heparinized solution.
- 6) Inject the heparinized solution by closing the clamp or removing the Port-Ject needle, exerting constant pressure.
- 7) Remove the Port-Ject needle and disinfect the skin.

6 SAFETY AND ENVIRONMENTAL REQUIREMENTS

- ☐ The products supplied to the consumer are non-toxic, do not contain substances that pose a danger to the environment and human health.
- ☐ The products do not have a sensitizing or irritating effect.
- ☐ The products do not contain components and substances that require particularly strict compliance with the requirements of industrial and fire safety during production.
- ☐ Fire safety requirements must be observed in the manufacture of products.
- ☐ Production facilities must be equipped with all necessary fire extinguishing equipment.

Environmental protection requirements

The main types of possible hazardous environmental impact are pollution of the atmospheric air of populated areas, soil and water as a result of unorganized incineration and disposal of waste materials on the territory of the manufacturer or outside it, as well as their arbitrary dumping in places not intended for these purposes.


- ☐ Production waste must be disposed of.
- ☐ It is allowed to dispose of waste materials and chemicals in the production process on a contractual basis with an organization licensed for waste disposal.

7 CONSTRUCTION

7.1 Huber needle PORT -JECT is disposable, sterile:

1. PRO-JECT A (Straight Guber needle, Luer connector), size:

19G (1,10) x 15mm; 19G (1,10) x 20mm; 19G (1,10) x 25mm; 19G (1,10) x 32mm;
 20G (0,90) x 15mm; 20G (0,90) x 20mm; 20G (0,90) x 25mm; 20G (0,90) x 32mm;
 22G (0,70) x 15mm; 22G (0,70) x 20mm; 22G (0,70) x 25mm; 22G (0,70) x 32mm;

L = 15 - 20 - 25 - 32 mm


D = 0.70 - 0.90 - 1.10 mm

PORT-JECT A:

CODE	PART	MATERIAL
1	NEEDLE	STAINLESS STEEL AISI 304
2	LUER CONNECTOR WITH SCREW CONNECTION	POLYCARBONATE (PC), MAKROLON 2458 Covestro AG brand, Germany
	PRIMARY PACKAGING	Packaging paper for medical products +plastic film for gas sterilization ETO (ETHYLENE OXIDE), Vereinigte Papierwaren Fabriken GmbH, Germany

2. PORT-JECT B (right-angled Huber needle, Luer connector), size:

19G (1,10) x 15mm; 19G (1,10) x 20mm; 19G (1,10) x 25mm; 19G (1,10) x 32mm;
 20G (0,90) x 15mm; 20G (0,90) x 20mm; 20G (0,90) x 25mm; 20G (0,90) x 32mm;
 22G (0,70) x 15mm; 22G (0,70) x 20mm; 22G (0,70) x 25mm; 22G (0,70) x 32mm;

L = 15 - 20 - 25 - 32 mm

D = 0,70 - 0,90 - 1,10 mm

PORT-JECT B:

CODE	PART	MATERIAL
1	NEEDLE 90°	STAINLESS STEEL AISI 304
2	LUER CONNECTOR WITH SCREW CONNECTION	POLYCARBONATE (PC), MAKROLON 2458 Covestro AG brand, Germany
	PRIMARY PACKAGING	Packaging paper for medical products +plastic film for gas sterilization ETO (ETHYLENE OXIDE), Vereinigte Papierwaren Fabriken GmbH, Germany


16

3. PORT-JECT C (Huber needle bent at right angle, complete with wings ("butterfly"), extension line, Luer connector, single action clamp), size:

19G (1,10) x 15mm; 19G (1,10) x 20mm; 19G (1,10) x 25mm; 19G (1,10) x 32mm;

20G (0,90) x 15mm; 20G (0,90) x 20mm; 20G (0,90) x 25mm; 20G (0,90) x 32mm;

22G (0,70) x 15mm; 22G (0,70) x 20mm; 22G (0,70) x 25mm; 22G (0,70) x 32mm;

PORT-JECT C:


CODE	PART	MATERIAL
1	LUER CONNECTOR	Polyvinyl chloride (PVC), brand rigid PVC (Nakan), Carro luigi, Italy
2	CLAMP	Polyvinyl chloride (PVC), brand rigid PVC (Nakan), Carro luigi, Italy
3	NEEDLE 90° with wings	STAINLESS STEEL AISI 304 / wings Polyvinyl Chloride (PVC (Nakan), brand FM 801 R, Carra luigi, Italy)
4	Wings	Polyvinyl chloride (PVC (Nakan)), brand FM 801 R, Carra luigi, Italy
5	EXTENSION LINE	Polyvinyl chloride (PVC (Nakan), brand FM 801 R, Carra luigi, Italy)
	PRIMARY PACKAGING	Packaging paper for medical devices and PVC blister (Polyvinyl chloride) for sterilization with ETHYLENE oxide gas, (Vereinigte Papierwaren Fabriken GmbH, Germany)

4. PORT-JECT D (right angled Huber needle, complete with wings ("butterfly"), extension line, Luer connector, Y-connector for the 2nd infusion line, single action clamp), size:

19G (1,10) x 15mm; 19G (1,10) x 20mm; 19G (1,10) x 25mm; 19G (1,10) x 32mm;

20G (0,90) x 15mm; 20G (0,90) x 20mm; 20G (0,90) x 25mm; 20G (0,90) x 32mm;

22G (0,70) x 15mm; 22G (0,70) x 20mm; 22G (0,70) x 25mm; 22G (0,70) x 32mm;

PORT-JECT D:


CODE	PART	MATERIAL
1	LUER CONNECTOR	Polyvinyl chloride (PVC), brand rigid PVC (Nakan), Carro luigi, Italy
2	CLAMP	Polyvinyl chloride (PVC), brand rigid PVC (Nakan), Carro luigi, Italy
3	Y-CONNECTOR	ABS plastic, Terluxe brand, BASF, Germany
4	NEEDLE 90° with wings	STAINLESS STEEL AISI 304 / wings Polyvinyl Chloride (PVC (Nakan), brand FM 801 R, Carra luigi, Italy)
5	Wings	Polyvinyl chloride (PVC (Nakan)), brand FM 801 R, Carra luigi, Italy
6	EXTENSION LINE	Polyvinyl chloride (PVC (Nakan), brand FM 801 R, Carra luigi, Italy)
	PRIMARY PACKAGING	Packaging paper for medical products+plastic film for gas sterilization ETO (ETHYLENE OXIDE), Vereinigte Papierwaren Fabriken GmbH, Germany

5. PORT-JECT E (right-angled Huber needle, complete with gripping system (needle clamp), extension line, Luer connector, single action clamp), size:

19G (1,10) x 15mm; 19G (1,10) x 20mm; 19G (1,10) x 25mm; 19G (1,10) x 32mm;

20G (0,90) x 15mm; 20G (0,90) x 20mm; 20G (0,90) x 25mm; 20G (0,90) x 32mm;

22G (0,70) x 15mm; 22G (0,70) x 20mm; 22G (0,70) x 25mm; 22G (0,70) x 32mm;

L = 15 - 20 - 25 - 32 mm

D = 0.70 - 0.90 - 1.10 mm

PORT-JECT E:


CODE	PART	MATERIAL
1	LUER CONNECTOR	Polyvinyl Chloride (PVC), rigid PVC brand (Nakan), Carro luigi, Italy
2	CLAMP	Polyvinyl Chloride (PVC), rigid PVC brand (Nakan), Carro luigi, Italy
3	GRIPPING SYSTEM (ACUPUNCTURE)	Polyvinyl Chloride (PVC), rigid PVC brand (Nakan), Carro luigi, Italy
4	NEEDLE 90°	STAINLESS STEEL AISI 304
5	EXTENSION LINE	Polyvinyl chloride (PVC (Nakan), brand FM 801 R, Carra luigi, Italy
	PRIMARY PACKAGING	Packaging paper for medical products and PVC blister (polyvinyl chloride) for gas sterilization ETO (ETHYLENE OXIDE), Vereinigte Papierwaren Fabriken GmbH, Germany

6. PORT-JECT F (right-angled Huber needle, complete with gripping system (needle clamp), extension line, Luer connector, Y-connector for the 2nd infusion line, single action clamp), size:

19G (1,10) x 15mm; 19G (1,10) x 20mm; 19G (1,10) x 25mm; 19G (1,10) x 32mm;

20G (0,90) x 15mm; 20G (0,90) x 20mm; 20G (0,90) x 25mm; 20G (0,90) x 32mm;

22G (0,70) x 15mm; 22G (0,70) x 20mm; 22G (0,70) x 25mm; 22G (0,70) x 32mm;

L = 15 – 20 – 25 – 32 mm



D = 0.70 – 0.90 – 1.10 mm

PORT-JECT F:



CODE	PART	MATERIAL
1	LUER CONNECTOR	Polyvinyl Chloride (PVC), rigid PVC brand (Nakan), Carro luigi, Italy
2	CLAMP	Polyvinyl Chloride (PVC), rigid PVC brand (Nakan), Carro luigi, Italy
3	Y-CONNECTOR	ABS plastic, Terluxe brand, BASF, Germany
4	GRIPPING SYSTEM (ACUPUNCTURE)	Polyvinyl Chloride (PVC), rigid PVC brand (Nakan), Carro luigi, Italy
5	NEEDLE 90°	STAINLESS STEEL AISI 304
6	EXTENSION LINE	Polyvinyl Chloride (PVC (Nakan), brand FM 801 R, Carro luigi, Italy
	PRIMARY PACKAGING	Packaging paper for medical products+plastic film for gas sterilization ETO (ETHYLENE OXIDE), Vereinigte Papierwaren Fabriken GmbH, Germany

7.2 Needle tip design (for all types of needles)

The tip of the needle is slightly curved relative to its rod so that the plane of the cut and the lumen of the needle are almost parallel to each other.

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
α - primary cut angle; β - needle point angle.

Table

Size G	Primary cut angle, (α), °	Needle point angle (β), °
19 (1,10)	10	20
20 (0,90)	10	20
22 (0,70)	10	20

8 TECHNICAL REQUIREMENTS

8.1 Dimensions and weights should be in accordance with Table 1.

Table 1.

VIEW	Size (G), mm	External diameter, mm	Needle length, mm	Weight, g	Extension length lines, mm ± 10 mm
PORT-JECT A	19	1.10	15/20/25/32	0,53/0,57/0,61/0,67	-
	20	0,90	15/20/25/32	0,54/0,59/0,64/0,71	-
	22	0,70	15/20/25/32	0,56/0,63/0,70/0,80	-
PORT-JECT B	19	1.10	15/20/25/32	0,48/0,52/0,56/0,62	-
	20	0,90	15/20/25/32	0,49/0,54/0,59/0,66	-
	22	0,70	15/20/25/32	0,51/0,58/0,65/0,74	-
PORT-JECT C	19	1.105	15/20/25/32	3,62/3,66/3,70/3,78	200mm
	20	0,90	15/20/25/32	3,63/3,68/3,73/3,80	200mm
	22	0,70	15/20/25/32	3,65/3,72/3,79/3,89	200mm
PORT-JECT D	19	1.105	15/20/25/32	6,21/6,25/6,29/6,35	220mm
	20	0,90	15/20/25/32	6,22/6,27/6,32/6,39	220mm
	22	0,70	15/20/25/32	6,24/6,31/6,38/6,48	220mm
PORT-JECT E	19	1.10	15/20/25/32	4,90/4,94/4,98/5,04	200mm
	20	0,90	15/20/25/32	4,91/4,96/5,01/5,08	200mm
	22	0,70	15/20/25/32	4,93/5,00/5,07/5,17	200mm
PORT-JECT F	19	1.10	15/20/25/32	7,08/7,12/7,16/7,22	200mm
	20	0,90	15/20/25/32	7,09/7,14/7,19/7,26	200mm
	22	0,70	15/20/25/32	7,11/7,18/7,25/7,35	200mm

! The length and size of the needle must match the parameters of the implantable port system. If you take a too short needle, then the introduction of the drug will be difficult due to the incorrect location of the needle tip in the silicone membrane, and in the case of blood sampling, it will be difficult to produce it with proper quality. If you choose an excessively long needle, then it will not fit snugly against the skin, which will lead to discomfort during the infusion.

8.2 Dimensional tolerances shall be $\pm 2\%$.

Needle tube length tolerances:

Length: up to 25mm: $+1/-2$ mm; 25 to 39mm: $+1.5/-2.5$ mm.

8.3 The cannula (connector) of the needle must be made of a colored or colorless material. If the material of the needle cannula, wings ("butterfly"), gripping system is colored, then its color must comply Table 2. Color coding for quick product identification (color of plastic cannula, wings ("butterfly"), gripping system depends on the diameter of the needle).

The color code is used for types A, B, C, D. It is not used for types E, F.

Table 2

Size (G), mm	Color of the cannula needle
22 (0,70)	Black
20 (0,90)	Yellow
19 (1,10)	White

8.4 Needle cannula (PORT-JECT C/D/E/F), extension line and protective cap: Polyvinyl chloride (PVC).

The needle cap must be made of a colorless material.

Attaching the cannula to the sleeve, connecting the parts of the product: epoxy resin.

8.5 Sharpening the needle. The tip of the needle has a trihedral atraumatic sharpening of the point, while each facet is ultrasonic polished for less painful insertion.

8.6 Lubricant – silicone, brand LIVEO MDX4-4159, DUPONT, USA.

Active silicone is an amino-functional dimethylsiloxane copolymer. If the needle tube is lubricated, then the amount of lubrication should be minimal and not form droplets of liquid on the outer and inner surface of the needle when viewed with normal or corrected to normal vision.

8.7 On the surface of the needle there should be no cracks, shells, nicks, scratches, chipped places, burrs, delaminations, burns and other contaminants (dross, grinding materials, polishing and traces of lubrication).

Wall of the Needle. The needle has a normal wall thickness (RW).

8.8 The blunting radii of the working parts of the needles should be no more than 0.03 mm.

8.9 The needle must be corrosion-resistant under the conditions of operation, transportation and storage.

8.10 The surface roughness Ra of the needles should not exceed 0.80 μm .

8.11 The length of the extension line is given in Table 1. The tolerance for the length of the extension line is ± 10 mm.

8.12 Needles must be sterile. Sterilization of needles is carried out by a gas method (ethylene oxide).

8.13 Needles must be suitable for operation at ambient temperature from plus 10 °C to plus 35 °C, relative humidity not higher than 80% at a temperature of 25 °C and atmospheric pressure from 86 to 106.7 kPa.

8.14 Needles must be resistant to vibration loads with the following parameters: frequency range from 10 to 55 Hz, displacement amplitude up to 0.15 mm.

8.15 For the manufacture of needles, the following materials should be used:

Component name	Material	Brand/manufacturer
the Huber needle is straight, curved	stainless steel	Brand AISI 304 и 302, Act-Med AG, Freiensteinau, Germany
luer connector (PORT-JECT A/B)	Polycarbonate (PC)	brand MAKROLON 2458 Covestro AG, Germany
single action clamp	Polyvinyl chloride (PVC (nakan))	brand rigid PVC (Nakan), Carro luigi, Italy
gripping system (needle forceps)		
luer connector, (PORT-JECT C/D/E/F)		
extension line		brand FTM 801 R, Carro luigi, Italy
protective cap for the needle wings ("butterfly")		
Y-connector	ABS plastic	brand Terluxe, BASF, Germany
lubricant	Silicone	brand LIVEO MDX4-4159, DUPONT, USA

Needle cannula (connector), Luer connector, single action clamp, butterfly wings: Dyes 100470 (White), 200840 (Black), 30021 (Yellow) manufactured by Elian S.A., Italy. All color options are provided by mixing basic dyes.

8.16 Individual packages must be made from:

for PORT-JECT A/B/D/F: DuPont™ Tyvek® LDPE (Medical Packaging Paper + Polyethylene Film), Vereinigte Papierwaren Fabriken GmbH, Germany.

for PORT-JECT C/E: DuPont™ Tyvek® brand packaging material (SEH brand paper and PVC blister (Polyvinyl chloride)), Vereinigte Papierwaren Fabriken GmbH, Germany.

8.17 Schemes and dimensions (linear, angular) of the needle point (see section 7).

8.18 Extension tube diameter:

for PORT-JECT C/E/F:

Outer diameter: 2.40 ± 0.10 mm

Inner diameter (ID): 1.50 ± 0.10 mm

for PORT-JECT D:

Outer Diameter: $2.5\text{mm} \pm 0.1\text{mm}$;

Inner diameter (ID): $1.3 \text{ mm} \pm 0.1 \text{ mm}$.

8.19 Connecting dimensions of Luer connector see section 7.

8.20 The shape and overall dimensions of the wings, see section 3. Tolerance ± 0.2 mm.

8.21 Maximum fluid flow $8 \text{ ml/s} \pm 0.5 \text{ ml/s}$ at a pressure of 2.549 kg/cm^2 .

8.22 The strength of the connection of parts of the needles (tensile strength) must be at least: 15N.

The connection of the tube with the needle head must be strong and tight.

8.23 The bending resistance of the needles should be:

Table 3

Size (G)	Bending force, N, ± 0.1	Maximum deviation, not more than, mm
19 (1,10)	9,7	0,71
20 (0,90)	9,0	0,48
22 (0,70)	6,7	0,42

8.24 Needles shall not break (resistance to fracture of the needle tube) when bent through an angle of 25° under a load applied at the distance between the rigid support and the point of application of the bending force given in Table 4:

Table 4

Size G	The distance between the rigid support and the point of application of the bending force, mm ± 0.1
19 (1,10)	27,50
20 (0,90)	25,0
22 (0,70)	17,5

9 MARKING

9.1 On the individual packaging of the needle must be applied:

- ✓ name and (or) trademark of the manufacturer;
- ✓ address of the manufacturer;
- ✓ authorized representative of the manufacturer (name, address, telephone number);
- ✓ description of the contents of the package;
- ✓ dimensions of the needle, extension line (if available);
- ✓ symbols corresponding to international standards "Sterile", "sterilized with ethylene oxide", "single use", "read instructions for use before use";
- ✓ batch number (series);
- ✓ date of sterilization (production);
- ✓ expiry date;
- ✓ do not reuse;
- ✓ sterile, apyrogenic, non-toxic;
- ✓ do not use if the integrity of the package is broken;
- ✓ dispose of after use;
- ✓ symbol CE;
- ✓ registration certificate number.

Symbols used for marking:

Symbol "Manufacturer" (name and address of the manufacturer)

YYYY-MM

Symbol "Use until"

Symbol "Batch number"

Symbol "Catalog number"

YYYY-MM-DD

Symbol "Date of manufacture"

Symbol "Ethylene oxide sterilization"

PYROGEN FREE NON TOXIC

Symbol "Caution! Please refer to the instructions for use"

Symbol "Do Not Reuse", "Single Use", "Use Only Once", or "Single Use Only"

XXXX (0426)

Symbol CE marking

Special disposal symbol (dispose of in the safest way)

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9.2 Each packaging box (multipack) must be labeled with a printed label:

- ✓ name and (or) trademark of the manufacturer;
- ✓ address of the manufacturer;
- ✓ authorized representative of the manufacturer (name, address, telephone number);
- ✓ description of the contents of the package;
- ✓ dimensions of the needle, extension line (if available);
- ✓ Catalog number (REF);
- ✓ symbols corresponding to international standards “Sterile”, “sterilized with ethylene oxide”, “single use”, “read instructions for use before use”;
- ✓ batch number (series);
- ✓ shelf life;
- ✓ do not reuse;

- ✓ sterile, apyrogenic, non-toxic;
- ✓ do not use if the integrity of the package is broken;
- ✓ dispose of after use;
- ✓ CE marking conforming to international standards;
- ✓ registration certificate number;
- ✓ quantity in a packing box;
- ✓ barcode.

Symbols used for marking:

Symbol "Manufacturer" (name and address of manufacturer)

YYYY-MM-

Symbol "Date of Manufacture"

DD

YYYY-MM


Symbol "Use until"

Symbol "Batch Number"

Symbol "Catalog number"

Symbol "Ethylene Oxide Sterilization"

Symbol "Carefully! Please refer to the instructions for use"

Symbol "Do Not Reuse", "Single Use", "Use Only Once", or "Single Use Only"

XXXX
(0426)

Symbol "CE marking"

Symbol "Sign of processing of cellulose products": 20 - cardboard
(located on the packing box)

Contains or contains phthalates (DEHP) (for PORT-JECT C, PORT-JECT D, PORT-JECT E, PORT-JEC F)

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9.3 Transport marking - depending on the types of vehicles used. Marking is applied directly to containers or labels: plywood, metal, paper or fiberboard. At the same time, the handling signs applied on the container must correspond to the values: "Fragile. Caution, Keep away from moisture."

Shipping container labeling must contain:

- ✓ name and address of the manufacturer and trademark (if any);
- ✓ authorized representative of the manufacturer (name, address, telephone number);
- ✓ name of the medical product;
- ✓ expiration date;
- ✓ date of manufacture;
- ✓ storage conditions;
- ✓ transportation conditions;
- ✓ batch number;
- ✓ number of units of consumer packaging of products in group packaging or shipping containers;
- ✓ net weight and gross weight on the group packaging of products or shipping container;
- ✓ CE mark;
- ✓ registration certificate number.

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10 PACKAGING, STERILIZATION and PRODUCTION PROCESS

10.1 PACKAGING

Individual packing

10.1.1 All products must be packed individually.

10.1.2 The design of the individual packaging must ensure the tightness of the product and prevent the penetration of microorganisms. The packaging of the product must not have any damage, holes/tears.

Individual packaging material:

for PORT-JECT A/B/D/F: DuPont™ Tyvek® LDPE (Medical Packaging Paper + Polyethylene Film), Vereinigte Papierwaren Fabriken GmbH, Germany.

for PORT-JECT C/E: DuPont™ Tyvek® brand packaging material (SEH brand paper and PVC blister (Polyvinyl chloride)), Vereinigte Papierwaren Fabriken GmbH, Germany.

10.1.3 The dimensions of the individual packaging differ depending on the version of the product.

10.1.3.1 The tensile strength in the longitudinal direction shall be at least 4.4 kN/m, in the transverse direction - at least 2.2 kN/m.

10.1.3.2 The thickness shall be 52 µm with a tolerance of ± 2 µm;

10.1.3.3 The tensile strength of the seams shall be > 1.2 N/15mm.

10.1.3.4 The width of the joints shall be: 6 ± 1 mm;

10.1.3.5 Opening force of the package: 5 ± 1 N.

10.1.3.6 Air permeability must be at least 0.2 µm/Pa s;

Group packaging

10.1.4 Products in individual packages, together with instructions for use, must be placed in a group package, which should provide protection from external mechanical damage.

Number of products in a group (consumer) package:

PORT-JECT A/B: 50 pieces/pack

PORT-JECT C/D/E/F: 20 pieces/pack

Secondary packaging material - two-layer cardboard, trademark GR Cartotecnica, Cartotecnica Sangiorgio, Italy.

The size of the group (consumer) package (LxHxW), cm:

PORT-JECT A/B: 36.5x10x10 ± 0.5 cm.

PORT-JECT C/D/E/F: 42.5x10x14 0.5 cm.

Shipping container

10.1.5 Each group package with needles, together with instructions for use, must be placed in a shipping container (corrugated cardboard box).

Corrugated cardboard box, trademark GR Cartotecnica, Cartotecnica Sangiorgio, Italy.

The box must be pasted over with paper tape, paper-based adhesive tape or polyethylene tape with an adhesive layer so that the package cannot be opened without breaking the integrity.

Quantity in a transport package - from 1 pack. up to 1000 pack.

10.2 STERILIZATION

INFORMATION ABOUT THE STERILIZATION METHOD, DESCRIPTION OF THE STERILIZATION AND PACKAGING PROCESS. – Sterilization method – gas, sterilizing agent: ethylene oxide. (See the document "Validation of the sterilization process").

10.3 PRODUCTION PROCESS

DESCRIPTION OF THE PRODUCTION PROCESS. – see the Document "Description of the Production Process".

11 REPAIR DETAILS

The product cannot be repaired and must be disposed of after the first use.

12 STORAGE AND TRANSPORTATION IN THE TERRITORY OF THE RUSSIAN FEDERATION

Transportation

12.1 Transportation of products can be carried out by all types of covered vehicles in accordance with the rules for the carriage of goods in force for each type of transport.

Transportation of cargo by sea should be carried out in accordance with the "Safety Rules for the Sea Transportation of General Cargoes". Type of shipment - containers and small shipment.

12.2 Products must be suitable for transportation at ambient temperature from minus 50 °C to plus 50 °C, relative humidity not higher than 80% at a temperature of 25 °C and atmospheric pressure from 86 to 106.7 kPa.

12.3 Products in shipping containers must be resistant to vibration loads with the following parameters: frequency range from 10 to 55 Hz, displacement amplitude up to 0.35 mm and mechanical shocks with parameters: peak shock acceleration up to 100 ms⁻², duration of impact acceleration 16 ms, beat frequency up to 120 min⁻¹.

Storage

12.4 Storage of products in transport packaging in the warehouses of the manufacturer (consumer), excluding the possibility of exposure to moisture and aggressive environments, at a distance of at least 1 m from heating devices, should be carried out at an ambient temperature of +5 °C to +35 °C and relative humidity air up to 80%.

12.5 When loading, unloading, storing and transporting, measures must be taken to protect the needle from damage and contamination.

12.6 Guaranteed shelf life of mirrors is 5 (five) years from the date of sterilization.

13 INFORMATION ON DISPOSAL IN THE TERRITORY OF THE RUSSIAN FEDERATION

13.1 Products must be disposed of in accordance with local regulations.

The product after use is subject to destruction in accordance with the norms and laws in force in the territory of the Russian Federation that regulate the disposal of potentially infectious medical waste (waste class B (SanPiN 2.1.3684-21)).

Products with expired shelf life and not used are disposed of as household waste of Class A.

13.2 The needle cannot be reused.

13.3 Recycling is also subject to individual and transport packaging. Paper, wood, polyethylene and plastic must be disposed of separately in accordance with the rules for the collection, accounting and disposal of packaging.

13.4 To work with used needles, personnel must be provided with personal protective equipment (dressing gown, scarf, gloves, mask).

13.5 Disposal is carried out in accordance with the current legislation of the Russian Federation.

14 MANUFACTURER WARRANTY

14.1 The manufacturer guarantees that the product complies with the requirements of this technical file subject to the rules of operation, transportation and storage.

14.2 Warranty period of storage (term of preservation of sterility) - 5 years from the date of sterilization.

14.3 A product with an expired sterilization period cannot be reprocessed.

14.4 In the event of defects or malfunctions, complete the claims form enclosed in the carton and keep the defective product.

Authorized representative of the manufacturer in the Russian Federation: MEDPARTNER LLC, legal address: 194100, St. Petersburg, Bolshoi Sampsonievsky Prospect, 68, lit. N, office 268.

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