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CERTIFICATE

中国国际贸易促进委员会暨中国国际商会
China Council for the Promotion of International Trade is China Chamber of International Commerce

01556917

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China Council for the Promotion of International Trade
China Chamber of International Commerce

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THIS IS TO CERTIFY THAT: the seal of DONGGUAN OUWEN MEDICAL DEVICE TECHNOLOGY CO., LTD. on the annexed DOCUMENT is genuine.

China Council for the Promotion
of International Trade

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查询网址 Website for verifying the certificate: <http://www.rzccpit.com/validate.html>

Dongguan Ouwen Medical Device Technology Building 3, No. 216, Huanchang North Changping Town, Dongguan City, Guang- dong Province, China	INSTRUCTION FOR USE	Version	01
		Date	03/2023

APPROVED BY

Enterprise Director of Dongguan Ouwen
Medical Device Technology Co., Ltd.

Hou Jianwei
(real signature)

«31» March 2023

INSTRUCTION FOR USE OF THE MEDICAL DEVICE: Microapplicator DENTAL FORMULA®

2023

Table of contents

Information on the medical device	3
Name of the medical device.....	3
Manufacturer of the medical device	3
1.3 Designation of the medical device.....	3
1.4 Scope of application of the medical device	3
1.5 Expected users of the medical device	3
1.6 Indications for use of the medical device	3
1.7 Contraindications for the use of the medical device	4
1.8 Adverse effects	4
1.9 Principle of operation of the medical device	4
1.10 Instruction for use (application) of the medical device.....	6
1.11 Conditions of use (application) of the medical device.....	6
1.12 Warnings and precautions.....	6
1.13 Shelf life of the medical device.....	6
2 Conditions for transportation, storage and disposal of the medical device	6
2.1 Conditions for transportation of the medical device.....	6
2.2 Conditions for storage of the medical device	6
2.3 Conditions for disposal of the medical device.....	7
3 Methods of cleaning and disinfecting the medical device.....	7
4 Technical description of the medical device	7
5 Levels of control and acceptable quality level	8
6 Information about packaging and scope of delivery of the medical device	9
6.1 Packaging of the Medical Device	9
6.2 Scope of Delivery of the Medical Device.....	10
7 List of National Standards Applied by the Manufacturer of the Medical Device	11
8 Circumstances under which a consumer should consult a healthcare professional	11
9 Warranties of the manufacturer	11
10 Information about the authorized representative of the manufacturer	11

1. Information on the Medical Device

1.1 Name of the Medical Device

Microapplicator DENTAL FORMULA®

Designs:

- 1) Microapplicator DENTAL FORMULA® Ultrafine, color: blue, orange, green, yellow, purple, mauve, white, black, in tubes of 10/50/100 pcs., 4/8/10/20 tubes per pack.
- 2) Microapplicator DENTAL FORMULA® Fine, color: blue, orange, green, yellow, purple, mauve, white, black, in tubes of 10/50/100 pcs., 4/8/10/20 tubes per pack.
- 3) Microapplicator DENTAL FORMULA® Regular, color: blue, orange, green, yellow, purple, mauve, white, black, in tubes of 10/50/100 pcs., 4/8/10/20 tubes per pack.
- 4) Microapplicator DENTAL FORMULA® Cylinder, color: blue, orange, green, yellow, purple, mauve, white, black, in tubes of 10/50/100 pcs., 4/8/10/20 tubes per pack.

The following terms and designations "medical device", "device" and "microapplicator" may be used to refer to the *DENTAL FORMULA® Microapplicator* medical device.

1.2 Manufacturer of the Medical Device

Manufacturer: Dongguan Ouwen Medical Device Technology Co., Ltd, China.

Manufacturing site: Room 403, Building 3, No. 216, Huanchang North Road, Changping Town, Dongguan City, Guangdong Province, China.

1.3 Designation of the Medical Device

The microapplicator is designed for cost-effective and convenient application of different medical gels, liquid or paste-like materials and other dental materials to the tooth/teeth of the patient.

1.4 Scope of Application of the Medical Device

The microapplicator is used in dental practice.

1.5 Expected Users of the Medical Device

Medical personnel involved in dental procedures.

1.6 Indications for Use of the Medical Device

Indications for the use of the medical device is the need to apply (insert) different medical gels, liquid or paste-like and other dental materials on the tooth/teeth of the patient.

1.7 Contraindications for the Use of the Medical Device

1.8 Adverse Effects

Not detected

1.9 Principle of Operation of the Medical Device

The microapplicator is designed to apply different compounds and materials to the surfaces to be treated. The device is comprised of a handle with a tip brush that holds the material to be applied on its surface/bristles. Non-slip surface and optimal rigidity of plastic rod allows easy fixation of the required tip angle and application of gels, varnishes, pharmaceuticals and bonding systems in any spatial position and in hard-to-reach places. If required, you can change the length of the handle of the microapplicator by easily breaking off the lower part (see fig.No.1, No.2).

In addition, you can bend the microapplicator to the required angle (see fig.No.3, No.4) (it is not recommended to bend the microapplicator many times).

The medical device does not adsorb liquids, solutions, varnishes, completely giving the material when applied to the working surface.

Different color marking of microapplicators minimizes the risk of erroneous application of medication, allowing to exclude cross-infection.

The plastic packaging of microapplicators in the form of a tube acts as a dispenser (the bottom cover of the tube has an outlet hole and a flap, thereby simplifying the work of staff with the device).

Figure No.1.


Figure No. 1

Figure No. 2

Figure No. 3

1.10 Instruction for Use (Application) of the Medical Device

Choose a suitable size of the microapplicator.

2. Carefully turn the bottom cap of the tube and take the microapplicator out.
3. Take out the required amount of dental material using the microapplicator.
4. Apply the material to the desired area of the oral cavity strictly according to the instructions for use and safety procedures for the compounds used.
5. Dispose of the microapplicator according to the rules of medical waste disposal.

1.11 Conditions of Use (Application) of the Medical Device

The device is used in dental or preventive medical institutions in an enclosed space in a dry place, away from direct sunlight, at temperature from $+15^{\circ}\text{C}$ to $+35^{\circ}\text{C}$ and relative humidity from 45% to 85%, in accordance with instruction for use of the medical device.

1.12 Warnings and Precautions

- The device is intended for single use;
- For individual use only;
- The medical device is non-sterile and should not be sterilized;
- Do not wet and keep away from moisture;
- Avoid exposure to direct sunlight.

1.13 Shelf Life of the Medical Device

The shelf life of microapplicators is 5 years.

2 Conditions for Transportation, Storage and Disposal of the Medical Device

2.1 Conditions for Transportation of the Medical Device

The device can be transported by all modes of transport (sea, road, rail or air transport) in accordance with the transportation rules in force for this type of transport. It is not allowed to violate the storage conditions during transportation.

2.2 Conditions for Storage of the Medical Device

Store device in a dry place protected against direct sunlight, heat sources, ignition, ultraviolet radiation, as well as ozone, at temperatures from 0°C to $+35^{\circ}\text{C}$ and relative humidity from 0% to 8

2.3 Conditions for Disposal of the Medical Device

The medical device shall be disposed of in the Russian Federation by an organization carrying out medical activities as class "B" waste in accordance with SanPiN 2.1.3684-21. Any unused microapplicators that have come into disrepair, or microapplicators with expired shelf life, shall be disposed of as "A" class waste in accordance with SanPiN 2.1.3684-21.

3 Methods of Cleaning and Disinfecting the Medical Device

The microapplicator is a non-sterile single-use medical device and shall not be cleaned. It is prohibited to reuse the medical device.

4 Technical Description of the Medical Device


Description and dimensions of the medical device. The microapplicator is a portable hand-held dental instrument in the form of a rod (handle) with a tip/brush that holds the material to be applied on its surface brush bristles (see Figure No.5).

Dimensions of the medical device, with regard to its design option, are shown in Table No.1.

Table No.1 – Dimensions of the medical device:

Design option	Rod length	Rod diameter	Tip/brush length	Tip/brush diameter
Ultrafine	10.2 ± 0.2 cm	0.30 ± 0.05 cm	0.20 ± 0.05 cm	0.15 ± 0.05 cm
Fine			0.30 ± 0.05 cm	0.20 ± 0.05 cm
Regular			0.35 ± 0.10 cm	0.27 ± 0.05 cm
Cylinder			0.40 ± 0.05 cm	0.20 ± 0.05 cm

Figure No.5


gn and physical form of the medical device tip/brush depending on the design option.

Figure No.6 –
Ultrafine Small
tip brush.

Figure No.7 - Fine
Medium tip/standard
shaped brush

Figure No.8 - Regular
Round tip/brush.

Figure No.9 -
Cylinder.
Elongated tip brush
shape.

Technical characteristics of the structure. The angle of flexion of the microapplicator and the force to break off the lower part of the medical device are presented in Table 2.

Table No.2 – Technical characteristics of the structure:

Characteristic	Indicator
Angle of flexion	50° - 75°
Force to break off the lower part	35 N - 37 N

5 Levels of Control and Acceptable Quality Levels

For quality control purposes, the medical device is sampled according to ISO 2859-1. The control levels and acceptable quality level (AQL) of the microapplicators are shown in Table No. 3.

Levels of control and acceptable quality levels (AQL):

Characteristic		Control level	Acceptable quality level (AQL)
Dimensions (length, diameter)		S-2	4.0
Visual defects	Large visual defects	C-1	2.5
	Small visual defects	C-1	4.0

6 Information on the Packaging and Scope of Delivery of the Medical Device

6.1 Packaging of the Medical Device

Tube. Devices of one design and color are packaged in plastic tubes in the form of a rounded cylinder, which provide protection of microapplicators against external mechanical damage and climatic factors during transportation and storage. It is unacceptable to pack microapplicators of different colors and or designs in a tube.

The packaging in the form of a tube acts as a dispenser (the bottom cap of the tube has an outlet hole and a flap), due to which it can be easily opened and closed by simply turning the cap (see Figure No.10, No.11).

The tube and the bottom cap of the tube are made of polypropylene (CAS No.9003-07-0), manufactured by Dongguan Ouwen Medical Equipment Technology Co.

The volume of the tube is up to 0.5 dm³. The tube withstands one fall to the bottom without damage and leakage, as well as the compression force in the axial direction: with a capacity of up to 0.5 dm³ incl - no less than 49 N.

Figure No 10



Figure No 11


Consumer packaging. The tubes with microapplicators are packed in boxes of lined cardboard of especially durable material, manufactured by Lee & Meng Paper Mfg. Ltd., Hong Kong. The consumer packaging provides additional protection of the microapplicators against external mechanical damage and climatic factors during transportation and storage.

Depending on the scope of delivery, the consumer packaging contains a certain number of tubes with microapplicators of one or more colors within one design option. It is not allowed to pack in the consumer package microapplicators of different designs.

Shipping box. The microapplicators in consumer packaging are packaged in cardboard shipping boxes, manufactured by Lee & Meng Paper Mfg. Ltd., Hong Kong. The cardboard boxes are sturdy enough to keep devices safe during transportation and storage.

6.2 Scope of Delivery of the Medical Device

Microapplicators are supplied in the following configurations: in tubes of 10/50/100 pcs., 4/8/10/20 tubes per pack..

List of National Standards Applied by the Manufacturer of the Medical Device

No 4 - List of standards used by the manufacturer

Designation	Name
ISO 13485:2016	Medical devices. Quality management systems. Requirements for regulatory purposes.
ISO 14971:1	Medical devices. Application of risk management to medical devices.
ISO 2859-1	Statistical methods. Procedures for sampling control by Alternative attribute. Part 1. Sampling control plans. Consecutive batches based on an acceptable level of quality.
ISO 15223-1	Medical device. Symbols, used in marking of the medical devices, labels and accompanying documents. Part 1. Basic requirements.

8 Circumstances under which a consumer should consult a healthcare professional

Not detected.

9 Warranties of the Manufacturer

The manufacturer guarantees the preservation of all parameters and characteristics of the microapplicators for 5 years from the date of manufacture, provided that the conditions for transportation, storage and use are in accordance with the instructions.

10 Information about the Authorized Representative of the Manufacturer

In case of any questions related to the circulation of the medical device in the Russian Federation, please, contact the Authorized Representative of the manufacturer:

"BE Plus" Company Limited (Be Plus Co. Ltd.), 6th Verkhny per. 12, liter A, part of the office, 5/3 1, Sergievskoe municipality, Saint Petersburg, Russia, 194292. +7 (812) 309-36-62, benovydenal.ru, beplus@beplus.spb.ru

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