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
## Intended Use

The **PerforMax EE Total Face Mask** is intended to provide an interface for application of CPAP or bi-level therapy to patients. This mask is for single patient use in the home or multi-patient use in the hospital/institutional environment. The small size mask is to be used on patients 7 years or older (>20 kg) for whom CPAP or bi-level therapy has been prescribed. The large and extra large size masks are to be used on patients (>30 kg) for whom CPAP or bi-level therapy has been prescribed.

 **Important Note:** An exhalation port is built into the PerforMax EE Total Face Mask so a separate exhalation device is not required.

The **PerforMax SE Total Face Mask** is intended to provide a patient interface for application of noninvasive ventilation. The mask is to be used as an accessory to ventilators which have adequate alarms and safety systems for ventilation failure, and which are intended to administer CPAP or positive pressure ventilation for treatment of respiratory failure, respiratory insufficiency or obstructive sleep apnea. This mask is for multi-patient use in the hospital/institutional environment only. The small size mask is to be used on patients 7 years or older (>20 kg) who are appropriate candidates for noninvasive ventilation. The large and extra large size masks are to be used on patients (>30 kg) who are appropriate candidates for noninvasive ventilation.

 **Important Note:** An exhalation port is not built into the PerforMax SE Total Face Mask. A separate exhalation device must be used with this mask.

**Caution:** U.S. federal law restricts this device to sale by or on the order of a physician.

 **Note:** This mask does not contain natural rubber latex or DEHP.

## ⚠ Warnings:

- This mask is not suitable for providing life support ventilation.
- Hand wash prior to first use. Inspect the mask for damage or wear (cracking, crazing, tears, etc). Discard and replace any components as necessary.
- The **PerforMax SE Total Face Mask** requires a separate exhalation device.
- The **PerforMax SE Total Face Mask** does not include an anti-asphyxia valve that would allow a patient to breathe if the ventilator failed.

- This mask is designed for use with CPAP or bi-level systems recommended by your health care professional or respiratory therapist. Do not wear this mask unless the CPAP or bi-level system is turned on and operating properly. **Do not block or try to seal the exhalation port.** **Explanation of the Warning:** CPAP systems are intended to be used with special masks with connectors which have vent holes to allow continuous flow of air out of the mask. When the CPAP machine is turned on and functioning properly, new air from the CPAP machine flushes the exhaled air out through the attached mask exhalation port. However, when the CPAP machine is not operating, enough fresh air will not be provided through the mask, and exhaled air may be rebreathed. This warning applies to most models of CPAP systems. Rebreathing of exhaled air for longer than several minutes can, in some circumstances, lead to suffocation.
- A minimum of 4 cm H<sub>2</sub>O (hPa) must be maintained when using the Performax EE Total Face mask.
- Do not block or try to seal the entrainment valve on the Performax EE Total Face Mask.
- If the Performax EE Total Face Mask is dropped, verify that the exhalation vents have not been obstructed prior to use.
- If oxygen is used with the device, the oxygen flow must be turned off when the device is not operating. **Explanation of the Warning:** When the device is not in operation, and the oxygen flow is left on, oxygen delivered into the ventilator tubing may accumulate within the device enclosure. Oxygen accumulated in the device enclosure will create a risk of fire.
- Oxygen supports combustion. Oxygen should not be used while smoking or in the presence of an open flame.
- At a fixed flow rate of supplemental oxygen flow, the inhaled oxygen concentration will vary, depending on the pressure settings, patient breathing pattern, mask selection, and the leak rate. This warning applies to most types of CPAP and bi-level machines.
- System pressure and oxygen therapy may be different for a Respirationics Performax mask than for a nasal or nasal/oral interface. When changing mask types, determine the appropriate pressure level and oxygen therapy. When using the elbow with built-in exhalation and air entrainment valve, continuous intentional leak exists at the exhalation ports, you must ensure that the pressure and oxygen therapy are maintained against this higher leak level.
- At low CPAP or EPAP pressures, the flow through the exhalation port may be inadequate to clear all exhaled gas from the tubing. Some rebreathing may occur.
- If an additional exhalation device is added to the patient circuit, you may need to adjust the pressure level to compensate for the additional leak of the exhalation device.
- Some users may experience skin redness, irritation, or discomfort. If this happens, discontinue use and contact your healthcare professional.
- The patient's physician should be contacted if the patient experiences the following symptoms while using the masks or after removing it: Unusual chest discomfort, shortness of breath, stomach distension, belching, or severe headache; drying of the eyes, eye pain, or eye infections; blurred vision. (Consult an ophthalmologist if symptoms persist.)
- Use of a nasal or full face mask may cause tooth, gum, or jaw soreness or aggravate an existing dental condition. Consult your physician or dentist if symptoms occur.
- This mask should not be used for patients who are uncooperative, obtunded, unresponsive, or unable to remove the mask.
- This mask is not recommended for patients who are taking a prescription drug that may cause vomiting.
- Appropriate patient monitoring should be used as medically necessary.
- The mask contains small parts which could result in a choking hazard.

#### Contraindications

This mask may not be suitable for use on patients with the following conditions: glaucoma, recent eye surgery or dry eyes, hiatal hernia, impaired cardiac sphincter function, excessive reflux, impaired cough reflex; or on patients unable to remove the mask by themselves.

#### Symbols

Warning or Caution

Note

Tip

Consult instructions for use

Leak Symbol and Value

Does Not Contain Natural Rubber Latex

## Before Use

- Read and understand the instructions completely.
- Hand wash the mask.
- Clean the patient's face.
- Verify that the mask and headgear are the correct size.
- If a nasogastric tube or similar device is in place, use the optional nasogastric tube sealing pad. Position the pad so that its flat surface is against the patient's face and the C-shaped opening surrounds the tube.
- Verify that the therapy device, i.e., ventilator, including the alarms and safety systems, has been validated prior to use.
- For the EE, verify that the air entrainment valve functions correctly.
- For the SE, check for the blue elbow on the mask, and note that the mask does not include an anti-asphyxia valve that would allow a patient to breathe if the ventilator failed.
- Inspect the mask and replace it if the cushion has hardened or is torn, or if any parts are broken.
- Verify therapy device pressure(s).

## Cleaning Instructions

1. Hand wash mask in warm water with a mild liquid dishwashing detergent.

**⚠ Caution:** Do not use bleach, alcohol, cleaning solutions containing bleach or alcohol, or cleaning solutions containing conditioners or moisturizers.

2. Rinse thoroughly. Air dry completely before use.

**⚠ Caution:** Do not attempt to clean the headstrap. Replace if the headstrap becomes soiled.

**Institutional Disinfection:** For multi-patient use in the hospital/institutional environment, use the Disinfection Guide to reprocess the mask between patients. These instructions can be obtained by visiting us online at [www.philips.com/performax](http://www.philips.com/performax) or by contacting Philips Respironics Customer Service at 1-800-345-6443 (USA or Canada) or 1-724-387-4000 or Respironics Deutschland ☎ +49 8152 93060.

## Sizing the Mask

To select the appropriate size mask, the mask cushion should surround the patient's face without obstructing the patient's sight.

**Figure 1-A:** The top of the cushion should rest comfortably against the middle of the forehead above the eyebrows.

**Figure 1-B:** The bottom of the cushion should rest comfortably under the mouth and above the chin.

## Connecting the Headgear and Mask

**Disconnecting:** Undo the headgear tabs and pull the straps through the slots on the top of the faceplate. Disengage the bottom clips from the mask with a gentle tug.

**Connecting:** Thread the headgear tabs through the slots on the top of the faceplate and fold them back to catch onto the headgear straps. Engage the bottom clips to the mask with gentle push.

## Mask Elbows

**Exhalation Elbow (EE):** The amber EE has built-in exhalation and an air entrainment valve.

**Standard Elbow (SE):** The blue SE does not have built-in exhalation or an air entrainment valve.


## Leak Symbol for the EE Elbow

Some ventilators may incorporate the use of a leak symbol and value in the mask selection setup procedures. The leak characteristics of this mask with the EE elbow is leak symbol (1 2). The leak symbol and value represents the intentional leak characteristics of the interface. On ventilators equipped with a Mask Selection control, enter the leak symbol value (1 2) that corresponds with the leak symbol value on the mask.

## Achieving the Right Fit


**⚠ Warning:** Verify the correct mask elbow is in place.

1. Install the elbow that is appropriate for the non-invasive ventilator with which the PerformMax will be used.

**Figure 1**

2. If using the EE elbow with built-in exhalation/air entrainment valve, verify that the air entrainment valve is functioning correctly. Locate the air entrainment valve flapper inside the mask elbow and verify it functions as follows:
  - With the airflow turned off, verify that the entrainment valve flapper is lying flat so that room air can flow through the large opening in the valve. Next, with the airflow on, the flapper should now cover the large opening, and air from the CPAP or bi-level device should flow into the mask. If the flapper does not close or does not function properly, replace.
  - **⚠ Warning:** Do not block the fresh air inlet on the entrainment valve. Ensure that the valve is not blocked with secretions and that the flapper is dry. Do not block or seal the exhalation vents.
3. Open the headgear to a large size setting and disconnect one (or both) of the bottom clips (Figure 2-D). Slide the headgear and mask on over the patient's head. Hold the mask over the patient's face by the faceplate protrusion. Do not hold the mask by the elbow.
4. Slide the headgear into position. The top strap will lay at the crown of the patient's head. The cross-strap will sit low on the back of the patient's head. Connect the bottom clips and tighten the straps using the headgear tabs. Do not overtighten. Final adjustments will be made when the patient is lying down.
5. The top of the cushion will contact the face just above the eyebrows (Figure 2-A). Adjust the top strap, but do not overtighten. The bottom of the cushion will contact the face just above the chin. The entire cushion should fit snugly, but comfortably against the patient's face. Re-evaluate the mask size if you cannot achieve facial contact points described in the Sizing the Mask section.
6. The patient should be lying down or the bed should be adjusted to a reclining position.
7. If using the EE elbow with the built-in Exhalation/Air Entrainment valve, connect the flexible tubing from the therapy device to the female 22 mm port. If using the SE elbow, connect the patient circuit to the SE elbow.

**Figure 2**  
 A = Mask Cushion  
 B = Headgear Tabs  
 C = Mask Elbow (EE shown)  
 D = Bottom Clip

8. Turn on the airflow. Verify there is only minor lifting of the mask from the patient's face (6.35 mm to 12.7 mm) during a full respiratory cycle. Have the patient assume various sleeping positions (supine, lateral, prone) with the mask in place.
9. Make any final adjustments using the headgear tabs (Figure 2-B).

### Removing the Headgear and Mask

If the patient is lying flat or against the bed, first disconnect one (or both) bottom clips from the mask and slide the mask upward. Gently pull the headgear free of the patient. If the patient can sit up, disconnect one (or both) of the bottom clips and slide the entire mask and headgear over the patient's head.


### Replacing the Mask Elbow

**To remove:** Hold the faceplate hub and gently pull and twist the elbow from the faceplate hub.

**To replace or change:** Press the mask elbow into the faceplate hub and slightly twist while holding the faceplate hub. Do not use excessive force or any tools.

**Figure 3**


**Figure 4**  
 A = Slot for Top Headgear Strap  
 B = Mask Cushion  
 C = Post for Bottom Headgear Clip  
 D = Faceplate Hub  
 E = Exhalation Port  
 F = Elbow with built-in Exhalation Ports and Air Entrainment Valve (EE)  
 G = Flapper (inside)  
 H = Pressure Pick-off Ports  
 I = Blue Standard Elbow (SE)  
 J = Accessory Swivel (connects elbows to tubing)

## Specifications

**Warning:** The technical specifications of the **PerforMax EE Total Face Mask** are provided for your healthcare professional to determine if it is compatible with your CPAP or bi-level therapy device. If used outside these specifications, or if used with incompatible devices, the mask may be uncomfortable, the seal of the mask may not be effective, optimum therapy may not be achieved, and leak, or variation in the rate of leak, may affect device function.

### Intentional Leak for EE Elbow

### pressure Drop cm H<sub>2</sub>O (hPa) for EE Elbow

|                         | 50 SLPM | 100 SLPM |
|-------------------------|---------|----------|
| Small/Large/Extra Large | 0.52    | 1.07     |

### Deadspace

Small = 375 mL

Large = 550 mL

Extra Large = 717 mL

### Storage Conditions

Temperature: -20° to +60°C

Relative Humidity: 15% to 95% non-condensing

### Disposal

Dispose of in accordance with local regulations.

### Reorder Numbers

Visit [www.philips.com/performax](http://www.philips.com/performax) for a complete list of products and reorder information, or contact Philips Resironics Customer Service at 1-800-345-6443 (USA or Canada) or 1-724-387-4000 or Resironics Deutschland at +49 8152 93060 for additional information.

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©2014 Koninklijke Philips N.V.