Instructions for Use
Transseptal Sheath
Dated June 12, 2023
Consisting of 18 pages

ACKNOWLEDGEMENT


State of Utah)


§

County of Salt Lake)

On this 12th day of June, 2023, before me, Jennifer Fordham, a notary public, personally appeared Tom Haueter, Director, Regulatory Affairs, proved on the basis of satisfactory evidence to be the person whose name is subscribed to the attached instrument, and acknowledge he executed the same.

Witness my hand an official seal.

Notary Public

My Commission Expires: June 27, 2025

APPROVAL

MERIT MEDICAL SYSTEMS, INC.

Organization

1600 West Merit Parkway, South Jordan, UT 84095, U.S.A.

Address

Director, Regulatory Affairs

Occupation

Haueter, Tom

Surname, Name

12.06.2023

Date DD.MM.YYYY

Signature

INSTRUCTIONS FOR USE

Transseptal Sheath for access to the heart CardiaGuide

In the event of any discrepancy between the English language text set forth in this document and the Russian language text, the English language text shall be controlling in all respects.

Carefully read all instructions prior to use. Observe all contraindications, warnings and precautions noted in these directions. Failure to do so may result in patient complications. Merit Medical Systems, Inc. relies on the physician to determine, assess, and communicate to each patient all foreseeable risks of the procedure.

MEDICAL DEVICE (MD) NAME AND CLASSIFICATION

MD name, variations of MD

Transseptal Sheath for access to the heart CardiaGuide

In variants:

1. CardiaGuide fixed curve transseptal sheath in composition:

1.1. CardiaGuide fixed curve transseptal sheath, effective length 63 cm, inner diameter 8,5F (2.7 mm), curve shape 50° – 1 pc.;

1.2. Dilator, effective length 67.4 cm – 1 pc.;

1.3. Guidewire – 1 pc.;

1.4. Instructions for use – 1 pc.;

2. CardiaGuide fixed curve transseptal sheath in composition:

2.1. CardiaGuide fixed curve transseptal sheath, effective length 63 cm, inner diameter 8,5F (2.7 mm), curve shape ML1– 1 pc.;

2.2. Dilator, effective length 67.4 cm – 1 pc.;

2.3. Guidewire – 1 pc.;

2.4. Instructions for use – 1 pc.;

3. CardiaGuide fixed curve transseptal sheath in composition:

3.1. CardiaGuide fixed curve transseptal sheath, effective length 81 cm, inner diameter 8,5F (2.7 mm), curve shape 50° – 1 pc.;

3.2. Dilator, effective length 85.4 cm – 1 pc.;

3.3. Guidewire – 1 pc.;

3.4. Instructions for use – 1 pc.;

4. CardiaGuide fixed curve transseptal sheath in composition:

4.1. CardiaGuide fixed curve transseptal sheath, effective length 81 cm, inner diameter 8,5F (2.7 mm), curve shape ML1– 1 pc.;

4.2. Dilator, effective length 85.4 cm – 1 pc.;

4.3. Guidewire – 1 pc.;

4.4. Instructions for use – 1 pc.;

Further in the text the following names of the medical device can be used: the Braided Transseptal Sheath, Transseptal Sheath, CardiaGuide, Sheath, CardiaGuide Transseptal Sheath and other shortened names that facilitate the perception of information.

MD – medical device.

HOW SUPPLIED:

Sterile: Sterilized with ethylene oxide gas. Non-pyrogenic.

Contents: One (1) Radiopaque Sheath, one (1) Radiopaque Dilator, one (1) 0.032" (0.81 mm) x 180 cm Guidewire

NOTE: Sheath Length, Diameter, and Curve configurations are indicated on the product label.

DESCRIPTION:

The Braided Transseptal Sheath is designed to provide a conduit to deliver diagnostic and therapeutic catheters to specific heart chambers and locations. It provides support for positioning and maintaining the position of catheters at specific locations in the heart. The sheath may be used for percutaneous entry.

The kit consists of three components: a sheath, a dilator, and a J-tipped guidewire.

The sheath has a radiopaque marker band which aids in defining the location of the tip, an atraumatic soft tip, and a lubricious coating on the inner and outer surface.

The dilator is designed to conform to the inner diameter and curve of the sheath, and has a tapered tip.

To facilitate access to a variety of cardiac structures and sites, the sheaths are available in various sizes, lengths and tip curve configurations.

WARNINGS:

1. Contents supplied STERILE using an ethylene oxide (EO) process. Do not use if sterile barrier is damaged.
2. For single patient use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.
3. The device(s) should be used by physicians engaged in the practice of specialized invasive cardiology techniques. Use of the device should be restricted to those physicians specifically trained in the approach to be used.
4. When the sheath is left in the vessel, a continuous heparinized infusion under pressure is strongly recommended through the sheath sideport.
5. Infusion through the sideport should only be done after all air is removed from the unit.
6. Dilators and catheters should be removed slowly from the sheath. Rapid removal may damage the valve components resulting in blood flow through the valve, as well as cause a vacuum which may allow air to enter the sheath.
7. Aspiration of the side port is recommended when withdrawing the catheter, probe, or dilator to remove any fibrin deposition which may have accumulated in or on the tip of the sheath.
8. Careful sheath manipulation must be performed in the presence of an implantable cardiac device of any kind to minimize the potential to displace or dislodge lead placement.
9. Direct percutaneous insertion of the sheath requires the use of the dilator to minimize the potential risk of vessel injury due to a flared tip.
10. Fluoroscopic monitoring of the location of the distal tip of the sheath using the radiopaque marker, especially when used in a transseptal approach, is recommended.

PRECAUTIONS:

1. Aspiration and flushing of the sheath, dilator, and catheter should be performed frequently to help minimize the potential for air embolism.
2. Indwelling sheaths should be internally supported by a catheter, electrode, or dilator.
3. Never advance, torque, or withdraw guidewire or sheath when resistance is met. Determine cause by fluoroscopy and take remedial action.
4. Use the sideport for injection or aspiration of sheath and sideport assembly. Assure that stopcock is in the closed position after flushing, to prevent back-bleeding.
5. The following conditions require that special care be taken when using this product involving the transseptal approach.
 - enlarged aortic root
 - marked right atrial enlargement
 - small left atrium
 - marked distortion of the thoracic configuration (e.g. kyphosis or scoliosis)

6. Care should be taken to avoid excessive bending of the sheath and/ or dilator before and during use.
7. Fluoroscopic procedures involve exposure to ionizing radiation by the patient and staff. Precautions to minimize exposure should be taken and protective equipment should be used.
8. Fluoroscopic guidance should be used when advancing the Braided Transseptal Sheath and/or dilator. When advancing the sheath and/ or dilator across a valve, a guidewire or pigtail should be used.
9. The sheath, dilator, and guidewire are designed for single use only. Reuse may expose the patient to communicable disease and/or injury.
10. Arrhythmias may occur during the use of any intracardiac device. Careful monitoring and availability of emergency equipment are mandatory.
11. When using the Braided Transseptal Sheath in the presence of radio frequency ablation, care must be taken to assure all ablating elements are outside the sheath.

ADVERSE REACTIONS:

Adverse reactions to cannulation of the peripheral vasculature and intracardiac placement of the sheath and dilator may include:

- infection
- local nerve damage
- perforation
- dissection
- AV fistula formation
- pseudoaneurysm formation
- arrhythmias
- hematoma
- hemorrhage
- thromboembolic events
- catheter entrapment
- valve damage
- pacemaker/defibrillator lead displacement
- air embolus
- vasovagal reaction
- vessel trauma
- vessel spasm
- atrial septal defect
- aortic puncture
- perforation and/or tamponade
- coronary artery spasm and/or damage
- stroke
- myocardial infarction
- pericardial/pleural effusion
- pulmonary edema

REQUIRED EQUIPMENT (See details in section "The possibility and ways of integrating with other medical devices" below):

Carefully read the instructions for each accessory before use.

- Needle of any manufacturer
- Heparinized Normal Saline of any manufacturer

INSPECTION PRIOR TO USE:

Carefully inspect the package for any breach of the sterile barrier or damage to the contents prior to use.

USE STERILE TECHNIQUE:

Suggested Procedure

1. Open package and place contents on sterile field.
2. Prep skin and drape in area of anticipated vein-puncture as desired.
3. Perform skin wheal using 25 gauge needle (not supplied).
4. Locate vessel using small gauge needle and syringe.
5. Insert Introducer Needle (not supplied) into vessel. Look for blood return to confirm proper needle position.
6. Insert soft tip of guidewire through Introducer Needle into vessel. Advance guidewire to required depth. At no time should the guidewire be advanced or withdrawn when resistance is met. Determine the cause of resistance before proceeding.
7. Hold guidewire in place and remove Introducer Needle. Do not withdraw the guidewire back into the cannula as this may result in separation of the guidewire. The cannula should be removed first.
8. Enlarge cutaneous puncture site with scalpel.
9. Assemble dilator and sheath together until the dilator hub locks into the sheath hub. Thread the dilator/sheath assembly over the guidewire, using a slight twisting motion. The Braided Transseptal Sheath may now be positioned to deliver catheters to desired locations.
10. Aspirate all air from the sheath valve assembly by using a syringe connected to the sideport. Flush the sheath through the sideport. If transseptal puncture is required, continue with steps below.

Suggested Procedure

1. Advance the Braided Transseptal Sheath and Dilator assembly into the superior vena cava (SVC) just above the right atrium.
2. Separate the dilator and sheath hub approximately 1cm while slowly advancing the sheath over the dilator. This will aid the introduction of the curved transseptal needle (not supplied).
3. While the sheath and dilator hubs are separated, slowly remove the guidewire from the dilator. Remove all air from the dilator by slowly aspirating blood. After making certain no air is in the dilator, flush the dilator.
4. Completely flush the transseptal needle.
5. Introduce the needle into the dilator hub. Carefully advance the curved section of the needle into the dilator, being certain not to restrict the movement of the needle.
6. Withdraw the sheath approximately one centimeter while maintaining the dilator position. Re-attach the dilator and sheath hubs.
7. While maintaining the sheath's position, slowly advance the curved section of the needle until it is about to protrude from the dilator tip.
8. Monitor the right atrial pressure by connecting the needle hub to pressure monitoring equipment. Good right atrial pressure should be observed before proceeding.
9. Position the needle and sheath set in the right atrium. Verify position using fluoroscopy.
10. Position the unit (dilator and needle point) against the atrial septum in the region of the fossa ovale by gradually rotating the needle posteriorly and toward the left scapula during withdraw. Use continual pressure monitoring and repeated anterior-posterior and lateral visualization of the tip under fluoroscopy during all positioning procedures.

11. After confirming the position of the dilator tip needle point against the atrial septum, advance needle and complete the transseptal puncture. Successful needle entry into the left atrium is confirmed by pressure monitoring and a sudden reduction in resistance. It is critical that acceptable left atrial pressure monitoring immediately after the needle penetration is felt through the interatrial septum. Do not advance the dilator if acceptable pressure is not seen. Disconnect the pressure monitoring line from the needle. This will show the location of the needle.

Reconnect the pressure monitoring line to the needle.

12. Advance the dilator with the needle in place through the septum. Acceptable left atrial pressure should be continually observed. A sequential increase in resistance to movement followed by a sharp decrease in resistance will indicate location of the dilator in the left atrium.

13. Withdraw the needle point even with the dilator tip. The dilator with the needle point within should be freely located in the left atrium. Verify with fluoroscopy.

14. Advance the sheath slowly over the dilator-needle combination until it is in the left atrium. A slow rotating motion of the sheath as firm pressure is applied will aid in this procedure. The sheath will be in position when a sharp reduction of resistance is felt.

15. Advance the sheath approximately 2 cm into the left atrium while maintaining the dilator needle position.

16. Disconnect the pressure monitoring line from the needle.

17. Slowly remove the needle from the dilator.

18. Slowly remove the dilator from the sheath.

19. Attach the sheath sideport to the monitoring line. Gently aspirate blood through the side arm for sampling and to be sure the sheath is clear of air.

Caution: Remove the dilator slowly to reduce the possibility of creating a vacuum in the sheath. Blood should aspirate freely through the sideport. If not, withdraw the sheath 0.5-1.0 cm (sheath tip may be resting against the wall of the atrium or a pulmonary vein). **Note:** Do not apply strong vacuum.

20. In order to maintain the location of the sheath in the left atrium, monitor the location of the radiopaque tip marker frequently under the fluoroscopy.

21. Introduce the properly prepared catheter through the hemostatic valve into the left atrium. Improved catheter manipulation may be obtained by withdrawing the sheath into the right atrium. The

sheath should be returned to the left atrium over the catheter prior to removing the catheter. Location of sheath can be confirmed by comparison to position established in step 20.

22. After removal of the sheath, use standard technique to achieve hemostasis.

INTERPRETATION OF SYMBOLS USED ON THE PACKAGE

	Caution!
	Inner diameter
	Outer diameter
	Non-pyrogenic
	Sterilized using Ethylene Oxide
	Do not resterilize

	Do not use if package is damage..
	Do Not Re-Use
	Keep away from sunlight
	Keep dry
	Contents: 1
	Distributor
REF	Catalog Number
	Manufacturer
LOT	Batch Code
	Date of Manufacture
	Use-By Date
MD	Medical Device
Sterile Package	Sterile Package
EC REP	Authorized Representative in the European Community

Additional information for Russian Federation

Further in the text the following names of the medical device can be used: the Braided Transseptal Sheath, Transseptal Sheath, CardiaGuide, Sheath, CardiaGuide Transseptal Sheath and other shortened names that facilitate the perception of information.

MD – medical device.

MD name, variations of MD	REF
Transseptal Sheath for access to the heart CardiaGuide	
Transseptal Sheath for access to the heart CardiaGuide: variant: CardiaGuide fixed curve transseptal sheath in composition: CardiaGuide fixed curve transseptal sheath, effective length 63 cm, inner diameter 8,5F (2.7 mm), curve shape 50° – 1 pc.; Dilator, effective length 67.4 cm – 1 pc.; Guidewire – 1 pc.; Instructions for use – 1 pc.	FCL20100

Transseptal Sheath for access to the heart CardiaGuide: variant: CardiaGuide fixed curve transseptal sheath in composition: CardiaGuide fixed curve transseptal sheath, effective length 63 cm, inner diameter 8,5F (2.7 mm), curve shape ML1 – 1 pc.; Dilator, effective length 67.4 cm – 1 pc.; Guidewire – 1 pc.; Instructions for use – 1 pc.	FCL20101								
Transseptal Sheath for access to the heart CardiaGuide: variant: CardiaGuide fixed curve transseptal sheath in composition: CardiaGuide fixed curve transseptal sheath, effective length 81 cm, inner diameter 8,5F (2.7 mm), curve shape 50° – 1 pc.; Dilator, effective length 85.4 cm – 1 pc.; Guidewire – 1 pc.; Instructions for use – 1 pc.	FCL20102								
Transseptal Sheath for access to the heart CardiaGuide: variant: CardiaGuide fixed curve transseptal sheath in composition: CardiaGuide fixed curve transseptal sheath, effective length 81 cm, inner diameter 8,5F (2.7 mm), curve shape ML1 – 1 pc.; Dilator, effective length 85.4 cm – 1 pc.; Guidewire – 1 pc.; Instructions for use – 1 pc.	FCL20103								
Intended use Transseptal Sheath for access to the heart CardiaGuide is intended for the percutaneous introduction of various types of cardiovascular catheters to all heart chambers, including the left atrium via transseptal puncture. This is a single-use device.									
Indications for use The CardiaGuide Fixed Curve Braided Transseptal Sheath is indicated for use in patients who require transseptal puncture for the treatment of cardiac conditions including ablation of cardiac arrhythmias.									
Contraindication There are no known contraindications for the CardiaGuide Fixed Curve Braided Transseptal Sheath. Contraindications are determined by the contraindications of the transseptal needle and cardiovascular catheters (not included) used in conjunction with the device in the treatment of cardiac diseases, including ablation of cardiac arrhythmias. Read the instructions for each compatible product carefully before use.									
Potential users of MD The device should be used in medical facilities by physicians engaged in the practice of specialized invasive cardiology techniques. This device should be used only by physicians thoroughly trained in percutaneous procedures.									
Use conditions of MD This medical device should only be used by specially trained personnel in hospitals and healthcare facilities.									
Operation conditions of MD Devices are recommended for use under operation conditions indicated in Table 1. Table 1 – Operation conditions									
<table border="1"> <tr> <th colspan="2">Operation conditions</th></tr> <tr> <td>Humidity</td><td>80-100%</td></tr> <tr> <td>Temperature inside human body</td><td>From +32°C to +42°C</td></tr> <tr> <td>Atmospheric pressure</td><td>From 700 hPa to 1060 hPa</td></tr> </table>		Operation conditions		Humidity	80-100%	Temperature inside human body	From +32°C to +42°C	Atmospheric pressure	From 700 hPa to 1060 hPa
Operation conditions									
Humidity	80-100%								
Temperature inside human body	From +32°C to +42°C								
Atmospheric pressure	From 700 hPa to 1060 hPa								
Type and duration of contact of a medical device with the human body - external communicating – circulating blood, limited contact (< 24 h)									

DESCRIPTION OF MAIN FUNCTIONAL ELEMENTS OF MD

General description of functional elements of MD

The CardiaGuide Fixed Curve Braided Transseptal Sheath Introducer is comprised of:

- A radiopaque sheath, which includes a curved, hubbed braided tube. The sheath has a radiopaque marker band which aids in defining the location of the tip, and an atraumatic soft tip. On the proximal end is a valve for hemostasis and a sidearm with a three-way stopcock.
- A radiopaque dilator, which consists of a curved plastic tube with a proximal plastic hub. The dilator is designed to conform to the inner diameter and curve of the sheath and has a tapered tip.
- A guide wire assembly, which consists of a soft tip guidewire, a guidewire dispenser and a J-tipped straightener.

The CardiaGuide Fixed Curve Braided Transseptal Sheath Introducer is available in two lengths: 63 and 81 cm, as well as two ML1 and 50° tip curves.

There are no differences in indications for use: the type of introducer is determined by the doctor depending on the anatomical features of the patient.

A pictorial representation of the device is provided in Figure 1


Figure 1. CardiaGuide Curve Braided Transseptal Sheath Introducer


The possibility and ways of integrating with other medical devices


The device can be used only with the compatible dilator and guidewire included in its composition, and a Transseptal Needle 18 or 19 Gauge and Tip Gauge 21 or 22 of any manufacturer, cardiovascular catheters of compatible outside diameter maximum 8.5F (2.83 mm) of any manufacturer, Heparinized Normal Saline of any manufacturer, 25 gauge needle of any manufacturer.


Technical characteristics of MD


All the values in this section have $\pm 5\%$ Tolerance if not otherwise indicated

Sheath (Drawings in Order FCL20100, FCL20101, FCL20102, FCL20103):


Dimensions in Order FCL20100, FCL20101, FCL20102, FCL20103

Specifications	Value
Usable Length, for FCL20100	63 cm
Usable Length, for FCL20101	63 cm
Usable Length, for FCL20102	81 cm
Usable Length, for FCL20103	81 cm
Sheath tube Inner Diameter	0.119 in/ 3,02 mm
Sheath tube Outer Diameter	0.144 in/ 3,66 mm
Overall Length, for FCL20100	84.39 cm
Overall Length, for FCL20101	84.39 cm
Overall Length, for FCL20102	102.39 cm
Overall Length, for FCL20103	102.39 cm
Sideport tube, Tube and Stopcock Joint, and Tip Pull Force	3.4 lbs minimum (2.34 N/cm ²)

Side Hole location	degrees from Sideport
	0.4" +/- 0.025" / 10,16 ± 0,64 mm from the tip edge
Side Hole Diameter	0.031 – 0.034" / 0,79 – 0,86 mm
Marker band location	6.4 mm maximum from the tip
Minimum inner diameter	.107" / 2,72 mm
Mass, g:	
For FCL20100, FCL20101	7.8
For FCL20102, FCL20103	8.7

Dilator

Specifications	Value
Outer diameter	0.112 in/ 2.84 mm
Maximum outer diameter	.1176"/ 2.9 mm
Minimum inner diameter	.033"/ 0.84 mm
Effective length (A)	67.4 cm and 85.4 cm
Mass, g:	
For FCL20100, FCL20101	5.9
For FCL20102, FCL20103	7.1
Guidewire:	
Maximum outer diameter	0.032 in/ 0.81 mm
Total Length	180 cm

Mass, g	5
Curve Type	J-curve

LIST AND DESCRIPTION OF MD MATERIALS
The brands/grades of medical device materials are listed in the extract from the manufacturer's technical documentation and can be provided to the user upon request

Device's components	Material name
Sheath	
Valve Housing	Polyurethane with colorant
Tube Assembly	Polyether Block Amide with BaSO ₄
Soft Tip	Polyether Block Amide with BaSO ₄
Braid wire	Stainless Steel
Marker Band	Platinum 90% / Iridium 10%
Sidearm Tubing	Pellethane
3-Way Stopcock	Polycarbonate / HDPE / Silicone
Valve Cap	Polyurethane with colorant
Hemostasis Valve	Silicone
Valve lubrication	Silicone
Sheath lubrication	Silicone 5% / Hexane 95%
Solvent Cement	50%/50% Cyclohexanone/Tetrahydrofuron
Adhesive	Cyanoacrylate
Dilator	
Tubing	Thermoplastic Polyester Elastomer with BaSO ₄
Hub	Thermoplastic Polyester Elastomer with colorant
Lubricant	Silicone 5% / 95% Hexane
Tipping lubricant - processing aid	Silicone 5% / 95% HEXANE
Guidewire	
Guidewire	Stainless Steel

Shelf life, use by date
Shelf life of the product is 3 years.
After expiration date this medical device is to be disposed.

LIST OF MATERIALS OF ANIMAL AND (OR) HUMAN ORIGIN
Not applicable. This medical device does not contain materials of animal and / or human origin.

LIST OF APPLICABLE STANDARDS FOR MD

The medical device has passed the verification and validation tests and meet the requirements of the following standards: ISO 13485: 2016, ISO 14971: 2012, EN 1041:2008, ISO 15223-1: 2016, ISO 11070: 2014, ISO 14644-1: 2015, ISO 14644-2: 2015, ISO 11135: 2014, AAMI TIR 28: 2016, ASTM D4169-09: 2009, ISO 10993-1: 2009, ISO 2233: 2000, ISO 10993-7: 2008, ISO 11607-1:2017, ISO 11737-1: 2006, ISO 8601: 2004, ISO 594-1: 1986, ISO 594-2: 1998, GOST 20790-93, GOST R ISO 11070-2010, GOST EN 556-1-2011, GOST ISO 14971-2021, GOST R ISO 15223-1-2020, GOST ISO 11607-1-2018, GOST ISO 11607-2-2018, GOST R 52770-2016, GOST 31209-2003, GOST 4011-72, GOST 31870-2012, MG 4.1.3166-14, GOST R 55227-2012, I-880-71, MVL.MN 1924-2003, MR 1436-76, MR 2413-81, MG 1704-77, GOST ISO 10993-7-2016, GOST 31214-2016, GOST ISO 10993-1-2021, GOST R ISO 10993-2-2009, GOST ISO 10993-4-2020, GOST ISO 10993-5-2011, GOST ISO 10993-10-2011, GOST ISO 10993-11-2021, GOST ISO 10993-12-2015, GPM.1.2.4.0005.15, GPM. 1.2.4.0003.15

INFORMATION ABOUT MEDICINES AND (OR) PHARMACEUTICAL SUBSTANCES INCLUDED IN MD

This section is not applicable to the product as the devices do not contain medicinal substances and are not used in combination with medicinal substances.

DESCRIPTION OF STERILIZATION METHOD

Merit Medical Inc.'s CardiaGuide Fixed Curve Braided Transseptal Sheath Introducer is sterilized with a 100% Ethylene Oxide (EO) gas cycle and was validated to Sterility Assurance Level (SAL) of 10^{-6} .

REQUIREMENTS FOR TECHNICAL SUPPORT AND REPAIR OF MD

The product is single-use, does not require maintenance and is not subject to repair. In case of damage, the product must be disposed of in accordance with the requirements of the legislation in the country of application.

Description of transport packaging

All the values in this section have $\pm 5\%$ Tolerance if not otherwise indicated.

Table 9 - Weight and dimensions of transport packaging

Variant	Packaging overall dimensions, (Length x Width x Height)	Weight
1	Length x Width x Height 78,74 x 15,88 x 19,37 cm	635 g
2	Length x Width x Height 96,52 x 15,88 x 19,37 cm	753 g
3	Length x Width x Height 105 x 17 x 18,5 cm	592 g

CONDITIONS OF TRANSPORTATION AND STORAGE

Storage conditions

Temperature	+5 - +25 °C
-------------	-------------

Humidity	20 – 60%
Atmospheric pressure	1010 +/- 10 hPa

Transportation conditions	
Temperature	0 - +40 °C
Humidity	20 – 60%
Atmospheric pressure	1010 +/- 10 hPa

MD DISPOSAL RULES AND CONDITIONS

The product has no harmful effect on the environment to the conditions of operation, transportation, storage and disposal.

Used products, as well as packaging, must be disposed of in accordance with the actual standards - SanPiN 2.1.3684-21 provided for the disposal of potentially polluting medical waste.

The hazard class of medical waste for used device is based on the characteristics of the morphological composition: class B (epidemiologically hazardous waste). Medical waste class of unused product – A.

In the case of expiration date or damage to the packaging, the products are also subject to disposal following the applicable legislation, as class A waste.

WARRANTY

Manufacturer warrants that this medical device is free from defects in both materials and workmanship. Any other express or implied warranties, including warranties of merchantability or fitness, are hereby disclaimed. Suitability for use of this medical device for any particular surgical procedure should be determined by the user in conformance with the manufacturer's instructions for use. There are no warranties that extend beyond the description on the face hereof.

Guaranteed shelf life of the product is 3 years.

NAME AND LEGAL ADDRESS OF MANUFACTURER AND AUTHORIZED REPRESENTATIVE IN RUSSIAN FEDERATION

DEVELOPER: MERIT MEDICAL SYSTEMS, INC.

1600 West Merit Parkway, South Jordan, UT 84095, U.S.A.

MANUFACTURER: MERIT MEDICAL SYSTEMS, INC.

1600 West Merit Parkway, South Jordan, UT 84095, U.S.A.

AUTHORIZED REPRESENTATIVE OF THE MANUFACTURER IN RUSSIA: Johnson & Johnson LLC, Krylatskaya street, 17 building 2, Moscow 121614, Russia

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