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CERTIFICATE

中国国际贸易促进委员会暨中国国际商会
China Council for the Promotion of International Trade is China Chamber of International Commerce

01502693

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China Council for the Promotion of International Trade
China Chamber of International Commerce

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THIS IS TO CERTIFY THAT: the seal of TAIAN DALU MEDICAL INSTRUMENT CO., LTD. on the annexed OPERATING INSTRUCTIONS is genuine.

China Council for the Promotion
of International Trade

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Authorized

Signature:

Yu Yaning

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(Date: Jan. 30, 2023)


证明书查询网址 Website for verifying the certificate: <http://www.rzccpit.com/validate.html>

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APPROVED BY
Taian Dalu Medical Instrument Co., Ltd
Director
Sherry Li
(signature)
28th January, 2023

OPERATING INSTRUCTIONS

Single-dose inhaler for delivery of powdered medicines in variations

Produced by
Taian Dalu Medical Instrument Co., Ltd
West Part of Yitianmen Street, Hi-tech Zone, Taian, 271000, Shandong, China

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CONTENT

1.	DEVICE INFORMATION	3
2.	INTENDED USE	4
3.	TECHNICAL SPECIFICATION	7
4.	SCOPE OF SUPPLY	9
5.	PACKAGING	9
6.	STERILITY	10
7.	MAINTENANCE	10
8.	LABELING	10
9.	COMPLIANCE WITH STANDARDS	13
10.	OPERATION CONDITIONS	13
11.	TRANSPORTATION AND STORAGE	13
12.	INFORMATION ON DISPOSAL	13
13.	SERVICE LIFE	14
14.	WARRANTY LIABILITIES	14
15.	DIRECTIONS FOR USE	14

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1. DEVICE INFORMATION

Name:

Single-dose inhaler for delivery of powdered medicines in variations (hereinafter referred to as the inhaler, device, product).

Single-dose inhaler for delivery of powdered medicines in variations:

1. Single-dose inhaler for the delivery of powder medicines DL-D01 comprised of:

- protective cap - 1 pc.
- mouthpiece - 1 pc.
- capsule camera - 1 pc.
- puncture button - 2 pcs.
- needle - 2 pcs.
- spring - 2 pcs.
- direction for use - 1 pc.

2. Single-dose inhaler for the delivery of powder medicines DL-D01 A, comprised of:

- protective cap - 1 pc.
- mouthpiece - 1 pc.
- capsule camera - 1 pc.
- puncture button - 2 pcs.
- needle - 4 pcs.
- spring - 2 pcs.
- direction for use - 1 pc.

Manufacturer:

Taian Dalu Medical Instrument Co., Ltd.

West Part of Yitianmen Street, Hi-tech Zone, Taian, 271000, Shandong, China

Phone: +86 538 8931968;

E-mail: sales@charamedical.com

Website: <http://www.charamedical.com>

Manufacturing site:

West Part of Yitianmen Street, Hi-tech Zone, Taian, 271000, Shandong, China

Authorized representative of the manufacturer in the Russian Federation:

PSK Pharma Limited Liability Company (PSK PHARMA Ltd.),

5 Programmistov Str., bld. 1, 141983, Dubna, Dubna Urban District, Moscow Region

Phone: 8 (499) 400 16 99, 8 (800) 234 16 99

E-mail: r4@rusbiopharm.ru

Class depending on the potential risk of use of medical application in accordance with the Nomenclature Classifier of Medical Devices approved by Order No. 4н of the Ministry of Healthcare of the Russian Federation dated June 6, 2012 — 2a.

Type of contact with the human body:

See Table 2 — Materials, type of contact with the human body.

Conditions of use: The device is intended for use at such medical institutions as hospitals, polyclinics, and doctors' consulting rooms, in conventional accommodation facilities, and outdoors.

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2. INTENDED USE

Single-dose inhaler for delivery of powdered medicines in variations is a portable device intended for the administration of powder medicines into the bronchi through the respiratory tract, for multiple use by one patient.

2.1. Information on potential consumers of the medical device

Patients with a wide range of respiratory diseases (adults, children, infants), using an inhaler as a means to deliver powder medicines to the bronchi through the respiratory tract, or healthcare professionals performing inhalation procedures in patients.

2.2 Scope

Anaesthetic and respiratory medical devices.

2.3 Precautions

The device should only be used dry.

Do not insert capsules into the mouthpiece hole.

Do not use excessive force when opening and closing the device to avoid its breakdown.

Do not use the device if damaged.

Do not push the buttons and do not open the air intake holes when inhaling.

Do not exhale via the mouthpiece.

Do not push the buttons more than once per capsule.

The device is intended for single-patient use. Do not allow other people to use the device.

2.4 Side effects

There is a possibility of occurrence of adverse events described in the instructions for medical use of medicines to be used in this device.

To reduce the risk of adverse events during the device use, it is required to adhere strictly to:

- instruction for use of medicines prescribed or recommended by physicians and permitted for use with this inhaler.

- instruction for use of this device.

2.5 Contraindications

Do not use the inhaler with medicines for the indications and in the dosages that differ from the indications and dosages described in the instruction for medical use of such medicines or from the indications and dosages prescribed by the attending physician.

The inhaler should not be used by patients with severe asthma and breathing difficulties.

2.6 Specific features of the device

Single-dose inhaler for delivery of powdered medicines in variations is a portable device that the patient can carry in their pocket or bag for convenience. This device is intended for multiple use by one patient.

The inhaler is used with a wide range of various medicines. The spectrum is based on either fragmented or spray-dried powders filled into solid gelatin capsules or hydroxypropyl methicellulose (HPMC) capsules.

Capsule-based delivery system is a well-established and convenient method for the delivery of medicines into the lungs. Patient's compliance with the dosage regimen is satisfactory since the dose is pre-filled into capsules that are stored in an air-tight package. This medicine administration form is an alternative to aerosol inhalation.

Capsule inhaler use is a simple procedure during which the patient's respiratory effort provides energy for capsule emptying. It eliminates the need to coordinate inhalation with the

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dosage delivered and also helps the patient adhere to the dosage regimen.

If necessary, the device can be disassembled for disposal with no medicine residues subject to proper cleaning. Therefore, it can be considered an environmentally-friendly delivery system.

The device contains no radiation sources. The device is not a source of electromagnetic interference.

Principle of operation

The principle of operation of the device consists in the use of the air flow exhaled by the patient as a source of energy for the dosage system emptying and delivery of minor medicine particles into the respiratory tract.

The general view of the inhaler is shown in Figure 1.


Figure 1 — General view of the inhaler

Description of operation principles

The functioning of the medical device is intuitive for the consumer since its operation principles are similar to the functions of existing devices: remove the protective cap, open the mouthpiece to access the capsule compartment, insert the capsule into the capsule compartment, close the mouthpiece, press the push buttons to puncture the capsule, and inhale. The device does not need calibration.

The user's inhale raises the capsule from the capsule compartment and makes it rotate at a high speed around its main axis, thus promoting the powder disaggregation.

The powder comes out of the capsule via two holes made by the puncture system and is inhaled by the patient via the mouthpiece.

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The lower part of the mouthpiece includes a mesh that prevents the capsule from entering the mouthpiece channel.

The product has very good basic performance in terms of total delivered dose and fine particle fraction.

The dose constancy that is common to most single-dose devices is further enhanced by the improved puncture system inherent to this inhaler.


Figure 2 — Description of the principle of inhaler operation

There are no limitations to the combined use of the medical devices.

Troubleshooting

Possible signs of failure, their causes and ways of elimination are given in Table 1.


Table 1

Failure	Signs	Elimination measures
Needle surface damage/contamination	Insufficient capsule puncture.	Clean the needle surface using a dry cloth.
Clogging of the mesh in the lower part of the mouthpiece with fine particles	Insufficient flow of fine particles up the mouthpiece	Rinse the mouthpiece with warm water (it must be completely dry before using the inhaler again). The capsule compartment can be cleaned with a small brush.

2.7 Device structure description of model designations

Variations DL-D01 and DL-D01 A have the same functional characteristics and vary only in technical characteristics, namely dimensions and shape.

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The inhaler consists of a protective cap, mouthpiece, puncture buttons, and capsule compartment, as shown below:


Figure 3 — Inhaler components

The service life of the inhaler is 36 months (3 years) from the manufacturing date specified on the device package.

CAUTION!

The medical device does not include any medicines. The inhaler is intended for the administration of medicines registered in the Russian Federation in accordance with the established procedure in the prescribed dosage.

2.8 Information on the circumstances requiring consultation of the consumer with a health care professional


If you experience pain, discomfort, or other unpleasant sensations while using the device, stop using it and contact your consulting physician.

3. TECHNICAL SPECIFICATION

3.1 Materials used in the manufacture of the device




Materials used in the inhaler components and type of their contact with the body are specified in Table 2.

Table 2

Image of an option component		Quantity of parts per assembly	Component name	Material of construction, grade	Manufacturer of the material of construction	Type of contact with the human body
DL-D01	DL-D01 A					
		1 pc.	Protective cap	Acrylonitrile butadiene styrene (ABS), PA-757 (+) (granules)	ZHENJIANG CHIMEI CHEMICAL CO., LTD	Surface contact with the body
		1 pc.	Mouthpiece	Acrylonitrile butadiene styrene (ABS), PA-757 (+) (granules)	ZHENJIANG CHIMEI CHEMICAL CO., LTD	Contact with the mucous membrane

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		1 pc.	Bottom cover of the device	Acrylonitrile butadiene styrene (ABS), PA-757 (+) (granules)	ZHENJIANG CHIMEI CHEMICAL CO., LTD	Surface contact with the body
		1 pc.	Capsule compartment cavity	Acrylonitrile butadiene styrene (ABS), PA-757 (+) (granules)	ZHENJIANG CHIMEI CHEMICAL CO., LTD	Surface contact with the body
		2 pcs.	Puncture button	Acrylonitrile butadiene styrene (ABS), PA-757 (+) (granules)	ZHENJIANG CHIMEI CHEMICAL CO., LTD	Surface contact with the body
			Bright blue dye	Panton SZ94502, food grade dye under 2011/65/EU	Jingliang Technology (Suzhou) Co., Ltd	
2 pcs.	4 pcs.	-	Needle	Medical grade steel, SUS 304	Yingtian Xinfu Hardware Technology Co., Ltd	-
		2 pcs.	Spring	Medical grade steel, SUS 304	Yingtian Xinfu Hardware Technology Co., Ltd	-
			Polyethylene pouch (for single package of the inhaler)	Non-oriented polypropylene film and polymer film GB/T 16886.5	Xiongxiang Jinglian plastic packaging Co., Ltd	-
			Shipper box	Corrugated cardboard	Taian Yiye industry and Trade Co., Ltd	-

3.2 Basic parameters

The dimensions of Single-dose inhaler for delivery of powdered medicines DL-D01 are specified in Table 3.

Table 3

Component	Length (mm)	Diameter (mm)	Width (mm)	Loops (qty)	Weight (g)
Protective cap	43.4 ± 0.2	-	30.68 ± 0.2	-	14.5 ± 1
Inhaler housing	Upper part height 32.9 ± 0.2	-	-	-	-
	Lower part height 20.4 ± 0.2		30.7 ± 0.2		
Needle	10.2 ± 0.2	1.2	-	-	-
Spring	11.6 ± 0.1	5.8	-	9	

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The dimensions of Single-dose inhaler for delivery of powdered medicines DL-D01 A are specified in Table 4.

Table 4

Component	Length (mm)	Diameter (mm)	Width (mm)	Loops (qty)	Weight (g)
Protective cap	55.5 ± 0.2	-	32.9 ± 0.2	-	18.7
Inhaler housing	Upper part height 30.5 ± 0.2		17.69 ± 0.2	-	
	Lower part height 17.5 ± 0.2		30.43 ± 0.2	-	
Needle	10.0 ± 0.2	1.2	-	-	
Spring	11.6 ± 0.1	5.8	-	7	

The package dimensions are specified in Table 5.

Table 5

Packaging	Length, cm	Height, cm	Width, cm
Polyethylene pouch	12.5	-	5.7
Shipper box	40.5	32	34

4. SCOPE OF SUPPLY

Supplied components include:

Single-dose inhaler for delivery of powdered medicines in variations:

1. Single-dose inhaler for delivery of powdered medicines DL-D01 comprised of:

- protective cap — 1 pc.
- mouthpiece — 1 pc.
- capsule compartment — 1 pc.
- puncture button — 2 pcs.
- needle — 2 pcs.
- spring — 2 pcs.
- direction for use — 1 pc.

2. Single-dose inhaler for delivery of powdered medicines DL-D01 A comprised of:

- protective cap — 1 pc.
- mouthpiece — 1 pc.
- capsule compartment — 1 pc.
- puncture button — 2 pcs.
- needle — 4 pcs.
- spring — 2 pcs.
- direction for use — 1 pc.

5. PACKAGING

The inhaler is packed in an individual plastic pouch for medical purposes and then in a multiple package in the quantity of 500 pcs.

CAUTION!

The product is supplied only packed in bulk, in a multiple package, by 500 pcs. (cardboard box).

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6. STERILITY

The device is not sterile. Not subject to sterilization prior to use.

7. MAINTENANCE

The inhaler should be cleaned every week.

Cleaning methods: Rinse the mouthpiece with warm water. Before using the inhaler again, it must be completely dry. The capsule compartment can be cleaned with a small brush. Press the buttons several times so that the perforation needles could also be cleaned. Do not use water to clean the inhaler base and/or perforation needles. Allow the inhaler to air dry.

The device cannot be repaired.

8. LABELING

The label for the multiple (shipping) package (cardboard box) contains the following information:

- name of the medical device;
- name, address, phone, e-mail of the manufacturer's authorized representative;
- manufacturer's name, logo, address;
- registration certificate number;
- part number;
- information on the medical device color;
- batch number;
- date of medical device manufacturing;
- use by date;
- quantity of devices per package;
- net and gross weight;
- Symbols: "CE", "Manufacturer", "Date of manufacture", "Use by..."

West Part of Yitannan Street,
Hi-tech Zone, Taian, Shandong,
China 271000

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_____ **CE** _____

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Variation labeling

Single-dose inhaler for delivery of powdered medicines DL-D01

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West Part of Yitianshen Street,
 Hi tech Zone, Taian, Shandong,
 China 271000

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Variation labeling

Single-dose inhaler for delivery - of powdered medicines DL-D01A

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Single-dose inhaler for delivery of powdered medicines DL-D01A

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

8.1 Interpretation of labeling symbols

Labeling symbols comply with GOST R ISO 15223-1. The interpretation of the symbols is given in Table 6.

Table 6

Symbol	Description
	Manufacturer's logo
	Manufacturer
	CE mark
	Date of manufacture
	Use by date

Dalu Medical

Taian Dalu Medical Instrument Co., Ltd.**泰安大陆医疗器械有限公司****9. COMPLIANCE WITH STANDARDS**

No.	Document No.	Document Title
1	2017/745/EU	Medical device regulation
2	MEDDEV 2.12/1, rev. 8	Supervisory control report form for the field safety corrective action report. Manufacturer's form. Field safety corrective action report
3	MEDDEV. 2.7.1 rev. 4	Clinical evaluation: a guide for manufacturers and notified bodies
4	EN ISO 14971:2019	Medical devices. Application of risk management to medical devices
5	EN ISO 15223-1:2021	Medical devices. Symbols to be used with medical device labels, labelling, and information to be supplied. Part 1. General requirements
6	EN ISO 20417:2021	Medical devices. Information provided by medical device manufacturers
7	EN ISO 10993-1:2020	Biological evaluation of medical devices — Part 1: Evaluation and testing within a risk management process
8	EN ISO 10993-5:2009	Biological evaluation of medical devices — Part 5: Tests for <i>in vitro</i> cytotoxicity
9	EN ISO 10993-10:2021	Biological evaluation of medical devices. Part 10. Tests for irritation and skin sensitization
10	EN ISO 10993-12:2021	Biological evaluation of medical devices — Part 12: Sample preparation and reference materials
11	ISO/IEC 17025:2017	General requirements for the competence of testing and calibration laboratories
12	EN 62366-1:2015/AMD 1:2020	Medical devices — Part 1: Application of usability engineering to medical devices / Amendment 1

10. OPERATION CONDITIONS

The device shall be operated at temperatures from +5°C to +35°C and relative humidity up to 70 %. It can be used in proximity with other medical devices.

11. TRANSPORTATION AND STORAGE

11.1 Inhalers in the manufacturer's multiple package are transported by all types of covered vehicles in accordance with the shipping rules applicable to each type of transport, at temperatures from +5°C to +35°C, relative humidity of up to 70 %.

11.2 The device shall be stored at temperatures from +5°C to +35°C and relative humidity of up to 70 %. Excessive heat must be avoided, and the device should be protected from direct sunlight and moisture.

The claimed shelf life is 3 years from the date of manufacture.

12. INFORMATION ON DISPOSAL

According to SanPiN 2.1.3684-21, Sanitary and epidemiological requirements for the maintenance of the territories of urban and rural settlements, for water bodies, drinking water and drinking water supply to the population, atmospheric air, soils, residential premises, operation of industrial, public premises, organization and implementation of sanitary and anti-epidemic (preventive) measures, in terms of waste hazard class depending on the degree of its epidemiological, toxicological and radiation hazard, as well as the negative impact on the environment, the medical device Single-dose inhaler for delivery of powdered medicines in variations belongs to Class A (municipal solid waste).

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13. SERVICE LIFE

36 months from the date of manufacture.

14. WARRANTY LIABILITIES

Manufacturer:

Taian Dalu Medical Instrument Co., Ltd.

West Part of Yitianmen Street, Hi-tech Zone, Taian, 271000, Shandong, China

Phone: +86 538 8931968;

E-mail: sales@charamedical.com

Website: <http://www.charamedical.com>

The manufacturer guarantees the compliance of this medical device with all the parameters and properties specified in this document, provided that the consumer observes the conditions of operation, transportation and storage throughout the entire guaranteed period.


The guaranteed service life is 3 years from the date of manufacture subject to compliance with the storage conditions.

15. DIRECTIONS FOR USE

1. Remove the protective cap from the inhaler.

2. Grab the middle part of the mouthpiece and the inhaler base, turn the mouthpiece in the direction indicated by the arrow on the housing, and open the capsule slot in the inhaler base.

3. With dry hands, take one capsule from the vial or blister and place the capsule in the capsule slot.

Caution! Do not insert capsules straight into the mouthpiece hole.

Do not swallow!

Do not use the device with capsules that are not intended for use with this device.

Do not store any capsules in the device.

4. Close the mouthpiece by turning it in the direction opposite from that indicated by the arrow, and hold the inhaler upright.

5. While holding the inhaler upright, push the buttons at the inhaler base and puncture the capsule.

Release the buttons. Exhale fully.

Caution! Due to the capsule destruction, its particles can enter the mouth or throat as a result of inhalation. In order to minimize such risk, do not pierce the capsule more than once.

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6. Place the inhaler into your mouth, clasp the mouthpiece tightly with your lips. Take a deep breath.

Caution! Do not push the buttons and do not open the air intake holes.

The capsule can rotate freely inside the compartment, creating extraneous noise.

7. Remove the mouthpiece from your mouth and hold your breath for as long as possible. Then exhale.

Caution! Do not exhale via the mouthpiece!

8. Open the capsule slot: grab the middle part of the mouthpiece and the inhaler base, turn the mouthpiece in the direction indicated by the arrow on the housing.

9. Remove the capsule shell and check it for punctures and powder.

If the capsule is not punctured and is still full, repeat the procedure from step 3.

If the capsule is punctured but not empty, reinsert the capsule into the capsule slot and repeat the procedure from step 6. Do not push the bottom more than once!

If the capsule is empty, discard it.

10. After use, if necessary, wipe the mouthpiece with a dry cloth, close the mouthpiece and put the protective cap back.

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