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



Lifotronic  **Wheisman**

Production name: Auto 5-part Hematology Analyzer

Model: AC 610

Power: a.c. 200V—240V, 50/60Hz, 300VA

Wheisman Medical Technology Co., Ltd.

District 1101A, 16th Bldg., Tongfu Industrial Park, Houdikeng, Shijia Community,
Matian Street, Guangming District, Shenzhen City, Guangdong Province, China

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

AC 610 Hematology Analyzer Operator's Manual

For this Operator's Manual, the issuing time is September, 2018 (Version: 1.5)
Whelshman Medical Technology Co., Ltd
11F, 16th Bldg, 8Jiang Technology Park, Zhoujia Blvd., Matian, Guangming District,
Shenzhen 518106, P. R. China
Email: info@whelshman.com
Tel/Fax: +86 755 86542379

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Concise Operation Procedure and Maintenance Guidance of The Auto 5-Part Hematology Analyzer

1 Installation

1. Please check whether the packaging of analyzer has any physical damage before unpacking. If there is any damage, please contact the service department or the regional agent timely.
2. After unpacking, check whether the items on the packing list are complete.
3. Open the right-side door of analyzer, remove the strapped blocks (two in total) in the horizontal-vertical direction of sampling probe.
4. Connect the analyzer with power supply and reagents in accordance with the installation requirements in user manual.

2 Preparation before startup

Check whether the reagent tubes are connected correctly and the remaining reagent is enough or not.

(I) Category and specification of reagent.

SN	Reagent	Specification	Corresponding joint on analyzer
1	DILUENT	20 L	DILUENT
2	LILYSE	250/500ml	LILYSE
3	LILYSE	250/500ml	LILYSE
4	PROBI. CLEANSER	50 ml	

(II) Waste Liquid Container: check whether the waste liquid container is full or not. If it is, it shall be replaced with empty waste liquid container.

3 Daily operations

3.1 Startup:

Power on analyzer, run the software, input the correct username and password (Initial user name is admin and password is also admin) and click the "Login" button, and then the analyzer will clean the tubes automatically.

3.2 Daily quality control

Enter L-J quality control mode to conduct daily quality control analysis of analyzer. If the measured result of each parameter is in the specified reference range of the quality control, sample analysis then can be conducted.

3.3 Sample analysis

1. Sampling and preparation

Whole blood mode: use vacuum tube with anticoagulant EDTAK₂ or EDTAK₁ to take blood samples, the blood volume must not be less than 1ml, the blood stored in cold environment must be placed at room temperature for at least 30 minutes and then can be tested.

Pre-diluted mode: add 20μL blood into the centrifuge tube with 400μL diluent and mix properly (diluent can be assigned by adding diluents function of analyzer or pipettor), then perform sample analysis.

Capillary whole blood: use 0.5ml centrifuge tube (with clot activator) to collect at least 50ul capillary blood and mix properly, then perform sample analysis.

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


Installation Instruction

The Hematology Analyzer Installation Instruction

V 2.0 (2018.01)

Preparation

1. Installation Requirements

a. Space

To ensure the needed space for maintenance, heat dissipation and avoid the liquid tubing be squeezed, please keep the analyzer:

- With ≥ 100 cm distance from wall in both Left and Right side
- With ≥ 50 cm distance from wall to back side

And make sure there is enough space under the workbench to place the Diluent and Waste barrel.

b. Power

- The Analyzer must be well grounded
- UPS is highly suggested using for the Hematology Analyzer
- Make sure the input power is as below required

Power Requirement	A.C. 100V ~ 240V	300 VA	50/60 Hz
-------------------	------------------	--------	----------


- Please use the power cord provided together with the analyzer, other power cord might be influence to the analyzer's performance and test results

c. Environment

- Temperature: $15^{\circ}\text{C} \sim 30^{\circ}\text{C}$
- Humidity: $30\% \sim 85\%$
- Air pressure: $70\text{ kPa} \sim 106\text{ kPa}$
- Avoid from direct sun light, direct wind from air conditioner, dust, vibration and electromagnetic interference etc.
- The analyzer must be put on a steady workbench

d. Tools

These tools are needed for the installation:

Pliers (Open the Wooden case and Cut the Strap inside the Analyzer)

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