TO THE
FEDERAL SERVICE ON
SURVEILLANCE IN HEALTHCARE
(ROSZDRAVNADZOR)

Contact: Dr. Daniela Schul
Fon: +4915114384269
Fax: -
Email: daniela.schul@bbraun.com
Internet: <http://www.bbraun.de>
Date: March 29, 2023

COVERING LETTER

With this letter we

B. Braun Melsungen AG
Carl-Braun-Str. 1
34212 Melsungen
Germany

as responsible/legal manufacturer according to the European Council Directive 93/42/EEC of medical devices of the product group

Infusomat® Plus Line and Infusomat® plus Line SafeSet for infusion by Infusomat® compact plus pump

herewith declare that we have taken into account the notes in accordance with letter no. 10-5475/23 issued on February, 02 2023 by Roszdravnadzor. Given that we are accepting the remarks raised by Roszdravnadzor and we submit for your attention and consideration the following updated document:

Instructions for Use (specific version of instructions for use (IFU) in English and Russian languages)
Infusomat® Plus Lines (Group: Infusomat® Plus Line SafeSet, Infusomat® Plus Line, Type IV-Standard)
Date of issue: 03.2023

for the medical devices registration and submission of the above document to the Federal Service for Surveillance in Healthcare (Roszdravnadzor).

Chairwoman of the Supervisory Board:
Dr. Annette Beller

Executive Board:
Markus Strotmann
(Chairman)
Priv.-Doz. Dr. Stefan Ruppert
Jurgen Stuhl

Corporate Office: Melsungen
Register Court
Local Court Fritzlar
HRB 11 000
WEEE-Reg.-No. DE 42690900


Address:
B. Braun Melsungen AG
Carl-Braun-Straße 1
34212 Melsungen
Germany

In case of any discrepancy between the English and Russian translations of this IFU, the English version shall prevail.


For and on behalf of

B. Braun Melsungen AG

i. A.

Aida Sadykova
Manager Regulatory Affairs
CoE IV-Systems

i. A.

Olga Siegmund
Administrator Regulatory Affairs
CoE IV-Systems

Instruction for use

Name of product

Infusomat® Plus Line and Infusomat® plus Line SafeSet for infusion by Infusomat® compact plus pump, variants:





1. Infusomat® Plus Line for infusion by Infusomat® compact plus pump, infusion (art.№: 8700310, 8700310-20) 100 pc./box, with instruction for use (1 pc.)
2. Infusomat® Plus Line for infusion by Infusomat® compact plus pump, infusion (art.№: 8700320) 100 pc./box, with instruction for use (1 pc.)
3. Infusomat® Plus Line for infusion by Infusomat® compact plus pump, infusion (art.№: 8700330) 100 pc./box, with instruction for use (1 pc.)
4. Infusomat® Plus Line for infusion by Infusomat® compact plus pump, infusion (art.№: 8700340) 100 pc./box, with instruction for use (1 pc.)
5. Infusomat® Plus Line SafeSet for infusion by Infusomat® compact plus pump, infusion (art.№: 8700200, 8700200-20) 100 pc./box, with instruction for use (1 pc.)
6. Infusomat® Plus Line SafeSet for infusion by Infusomat® compact plus pump, infusion (art.№: 8700210) 100 pc./box, with instruction for use (1 pc.)
7. Infusomat® Plus Line SafeSet for infusion by Infusomat® compact plus pump, infusion (art.№: 8700220) 100 pc./box, with instruction for use (1 pc.)
8. Infusomat® Plus Line SafeSet for infusion by Infusomat® compact plus pump, infusion (art.№: 8700240, 8700240-20) 100 pc./box, with instruction for use (1 pc.)
9. Infusomat® Plus Line SafeSet for infusion by Infusomat® compact plus pump, infusion (art.№: 8700250) 100 pc./box, with instruction for use (1 pc.)
10. Infusomat® Plus Line SafeSet for infusion by Infusomat® compact plus pump, infusion (art.№: 8700300) 20 pc./box, with instruction for use (1 pc.)
11. Infusomat® Plus Line SafeSet for infusion by Infusomat® compact plus pump, infusion, UV protect (art.№: 8700260, 8700260-20) 100 pc./box, with instruction for use (1 pc.)
12. Infusomat® Plus Line SafeSet for infusion by Infusomat® compact plus pump, infusion, UV protect (art.№: 8700270) 100 pc./box, with instruction for use (1 pc.)

(further in the text «Infusomat® Plus Line and Infusomat® Plus Line SafeSet, infusion»)

Intended Use

Intended for continuous or intermittent IV administration by a compatible pump (primary method) and by gravity (secondary method). Intended for IV administration of infusion solutions, liquid drugs and fat emulsions, as well as low-soluble drugs and cytotoxic drugs.

Device description

			
8700310, 8700310-20, 8700320, 8700200, 8700200-20, 8700210, 8700220, 8700260, 8700260-20	8700340	8700240, 8700240-20, 8700250, 8700330, 8700270	8700300


1 – Drip chamber, which consists of: 1.1 protective cap (preserves the sterility of the medical device until it is used); 1.2 spike (provides easy piercing of various ports of infusion containers); 1.3 ventilation duct with lid and air filter - provides additional protection against contamination; 1.4 moulding ring - provides a tight connection between the upper and lower parts of the drip chamber; 1.5 the integrated dropper into the upper part of the drip chamber has standard characteristics: 20 drops = 1 ml + / - 0.1 ml (distilled water); 1.6 the lower part of the drip chamber is elastic,

so it can be used to pump the solution into the drip chamber; 1.7 liquid filter - is integrated in the lower part of the drip chamber, depending on the version of the medical device, it may have the following characteristics: - 15 µm fluid filter (Infusomat® Plus Line), 15 µm fluid filter Air-Stop (hydrophilic membrane, impermeable to air bubbles) (Infusomat® Plus SafeSet Line)). 2 – Tube - is designed for the introduction of solutions to the patient, transparent (PVC or PUR). The polymer tube with UV protection (orange transparent) is used for photosensitive drugs. 3 – Roller clamp has a connector for attaching a tube that can be used to temporarily interrupt the infusion and a for safely disposing of the tip after the infusion is completed. 4 - Peristaltic (pump) segment it is made of silicone and when exposed to the peristaltic mechanism of the pump, it provides a uniform unidirectional movement of the injected solution at a given speed. For the correct placement of the line inside the pump, it is provided: (clip (inlet/outlet) and places inside the pump, having the appropriate shape and color coding and graphic image ("asterisks") to prevent twisting of the tube during the installation of the main line in the pump.). 5 - Free Flow Protection Clamp - allows to block the flexible tube to stop the flow of the solution if necessary. 6 - Y-port - it is intended for additional bolus administration of drugs. Depending on the version, the Y-port can be needle or needle-free. 7 - Filter 0,2 µm (Sterifix) prevents the ingress of solid particles by filtering the solution. It is a non-removable component of a medical device. A filter of 0.2 microns is provided for infusion lines used for the introduction of standard infusion solutions, as well as low soluble infusion solutions and cytostatic drugs. 8 - Connector Luer Lock (type male) with protective cap - Connector Luer Lock ensures tightness of the connection between the line and the patient's catheter (due to the presence of internal thread), protective cap - depending on the version of the medical device, it can be: standard plastic or equipped with a hydrophobic membrane Prime-Stop. The presence of the PrimeStop hydrophobic membrane allows you to fill the line without removing the cap and leave the connector sterile until the line is connected to the patient.

Materials used:

PS-HI, PE, PA, SBC/SBS, PP, ABS, Borosilicate, PVC, MABS, PUR, PC, MBS, PTFE, PSU/PVP, PESU, Acrylic copolymer, SIR, SIL, IR, SB, PS, PE-LD.

Operating Instructions for infusion by compatible pump (primary method)

1. Select Infusomat® plus Line or Infusomat® plus SafeSet Line, infusion in accordance with the planned procedure.
2. Check the shelf life and integrity of the trunk packaging. Do not use if the packaging is damaged or after the expiration date.
3. Check the name, expiration date, physical properties, and dosage of the medicine.
4. Open the individual packaging, remove and expand the Infusomat® plus Line or Infusomat® plus SafeSet Line, infusion observing the rules of asepsis.
5. Close the roller clamp (Fig. 2).
6. Treat the container port, remove the protective cap from the tip and puncture the port (Fig. 1).
7. Squeezing the lower part of drip chamber, fill it with sufficient amount of the solution (Fig. 3).

8. During infusion of the polymer bottles or bags, the lid of the ventilation duct of the tip should be closed and during infusion of glass bottles the ventilation channel should be open (Fig. 4).
9. Open the roller clamp and fill the line with a solution, then close the roller clamp (Fig. 5).
10. Install the Infusomat® Plus Line according to the instructions for the peristaltic pump. Open the pump door pushing the door level. Install the green snap clamp in the appropriate pump slot (Fig. 6 pos. 2), and the orange retainer of the peristaltic segment to the corresponding place (Fig. 6 pos. 1).
11. Push the Infusion line into the air sensor on the left and into the pump based safety clamp on the right (Fig. 7).
12. After installing the line, close the pump door while simultaneously pulling the door level (Fig. 8).
13. Connect the line to the patient's catheter (Fig. 9).
14. Open the roller clamp (Fig. 10)
15. Set infusion parameters on the pump and start infusion.
16. Before removing the line from the pump, close the roller clamp (Fig. 11), disconnect the line from the patient's catheter (Fig. 12)
17. Open the pump door by pulling the door lever and remove the line from the pump.

Compatible infusion pumps according to the manufacturer's instruction for use of infusion pumps. Infusomat® Plus Lines are compatible with Infusomat compact plus (Registration Certificate № RZN 2022/16352).

Operating Instructions for infusion by gravity (secondary method)

Follow points from 1 to 9 from the chapter "Operating Instructions for infusion by compatible pump".

After that:

10. Connect the line to the patient's catheter (Fig. 9).
11. Set the required infusion rate by adjusting the roller clamp. Periodically monitor the flow rate throughout the duration of the infusion.
12. After infusion is finished, close the roller clamp (Fig. 11) and disconnect the line from the patient's catheter (Fig. 12).

Sterility

Sterilized by EO. Please refer to the information on primary packaging. Do not use if package is damaged. Visually inspect the product packaging to ensure that the sterile barrier system is intact.

Intended purpose

Single-use, sterile infusion systems for administration of infusion solutions, parenteral nutrition via one single closed line system, usable with a compatible pump (primary method) or by gravity (secondary method). Do not re-use.

Indications

8700310, 8700310-20, 8700320, 8700200, 8700200-20, 8700210, 8700220, 8700340 8700240, 8700240-20, 8700250, 8700330

Indications: for standard infusion solutions and parenteral fluids (parenteral infusion solutions and lipids)

8700300

Indications: for low soluble IV-solutions and for cytostatic solutions.

8700260, 8700260-20, 8700270

Indications: for light sensitive drugs.

Patient population

The device can be used for all patients for which infusion therapy is prescribed. Infusomat® Plus Lines can be used for adults, pediatric and neonates in hospitals and also in the outpatient market.

Selection of pump and dedicated set according to national standards and / or hospital protocols.

Intended User

Infusomat® Plus Lines are intended to be used by trained healthcare professionals, licensed and/or certified physicians, healthcare practitioners and/or specifically trained patients or caregivers according to local regulations.

Contraindications

Infusomat® plus Lines IV shall not be used in case of proven incompatibility of the drug/ solution to one of the used materials. SPCs of the drugs/solutions for contraindications should be considered.

No administration of blood products, blood components and enteral nutrition solutions.

Risks / Side effects

General risk occurring during infusion therapy including:

- Air embolism
- Local and systematic infection
- Medication error including over and under supply
- Particle contamination
- Kinking of the line may result in decrease or stop of infusion flow
- Leakage of drugs/fluids can be an additional hazard to patients, visitors or staff

Precautions

- Visual check of the Infusomat® Lines for damages before use.
- Do not re-attach the protective cap of the patient connector.
- Always refer to the manufacture's instruction for use of infusion pumps.
- If pressure infusion without infusion pump is used, remove air from container beforehand. Observe general safety pre-cautions.
- Mind the general guidelines for parallel infusions.
- Please also consider the corresponding SPCs of the drugs/ solutions regarding possible incompatibilities of drugs/ solutions and/or material of the used Infusomat® Line.
- Monitor all connections for tightness over the whole infusion procedure
- Due to the structure of the fluid filter of Infusomat® Plus SafeSet which achieves a particle retention with a filter efficiency > 80% for a particle size of 3µm, the infusion of highly viscous solutions, lipid emulsions >10% and certain suspensions can affect the flow rate. With regard to this, please also follow the information in the relevant summary of product characteristics.

Warnings

Do not use if packaging or product is damaged or contaminated or if protective caps are loose or missing.

Do not re-use: Re-use of single-use devices creates a potential risk of patient or user. It may lead to contamination and/or impairment of functional capability. Contamination and/or limited functionality of the device may lead to injury, illness or death of the patient. Administration of light sensitive drugs requires UV-protected IV administration sets.

Always check, if incompatibilities with other drugs or with respective material of different types of Infusomat® Lines is stated in the SPC (Summary of product characteristics) of the drug.

Always check, if the use of a 0.2µm filter is stated in the SPC (Summary of product characteristics) of the drug.

Duration of use

Usage time according to hospital protocols and/or guidelines. Maximum usage time of one device is up to 96h.

General technical characteristics

Article number	Length of device, mm (± 150 mm) / Length of line from peristaltic segment to Luer Lock connector, mm (± 100 mm)	Filling volume (ml) on 1 m of polymer tubing under 40 °C in accordance to EN ISO 8536-8	Storage volume (ml) on 1 m of polymer tubing under 40 °C/2 bar in accordance to EN ISO 8536-8	Tube (material colour)	Filter for liquids (15 µm)	Filter for liquids (15 µm) AirStop type	Filter 0,2 µm, (Sterifix)
8700310, 8700310-20	2400 / 1500	6	7	PVC, colorless	+	-	-
8700320	2400 / 1500	7	7	PUR, colorless.	+	-	-
8700330	3000 / 2000	6	7	PVC, colorless.	+	-	-
8700340	3000 / 2000	6	7	PVC, colorless.	+	-	-
8700200, 8700200-20	2400 / 1500	6	7	PVC, colorless.	-	+	-
8700210	3000 / 2000	6	7	PVC, colorless.	-	+	-
8700220	2400 / 1500	7	7	PUR, colorless.	-	+	-
8700240, 8700240-20	3000 / 2000	6	7	PVC, colorless.	-	+	-

8700250	3000 / 2000	7	7	PUR, colorless.	-	+	-
8700300	2400 / 1500	7	7	PUR, colorless.	-	+	+
8700260, 8700260-20	2400 / 1500	7	7	PUR, orange, UV - protect	-	+	-
8700270	3000 / 2000	7	7	PUR, orange, UV - protect	-	+	-

Internal/outer tube diameter - 4.1 x 3.0 mm ($\pm$ 0,2 mm).

Resistance to pressure up to 2.0 bar.

Luer compatibility acc. to EN ISO 80369-7/ EN ISO 80369-20

Spike compatibility acc. to EN ISO 8536-4/ EN ISO 8536-8

Scope of supply

Scope of supply of the medical device text "Infusomat® Plus Line and Infusomat® Plus Line SafeSet, infusion" is a group package that includes medical device in individual packaging – 20, 100 pcs (depending on the variants, see ch. "Name of medical device" of the instruction for use).

The group package is accompanied by instructions for use on paper - 1 pc.

Operation environment, storage and transportation conditions

Operation environment

+10°C - +40°C, humidity up to 80 %, atmospheric pressure 900 mbar - 1100 mbar

Transportation

-20°C - +40°C, humidity up to 80%, atmospheric pressure 900-1100 mbar

Storage

+15°C to +25°C, humidity up to 80%, atmospheric pressure 900-1100 mbar

Authorized organization (importer) in the Russian Federation

Authorized organization (importer) in the Russian Federation: LLC «B.Braun Medical», 191040, Saint Petersburg, Pushkinskaya 10, tel/fax (812) 320-40-04.

Disposal

Use integrated spike protection at roller clamp for disposal. Disposal of the medical device «Infusomat® Plus Line» after use is carried out in accordance with the requirements for the treatment of medical waste established by SanPiN 2.1.3684-21 for class B and/or other applicable legal acts of the Russian Federation.

Unused medical product, unused packaging should be disposed of in accordance with the procedure of disposal of solid domestic waste established by the legislation of the Russian Federation on the according to the requirements of SanPiN 2.1.3684-21 as the waste is class A and/or other applicable legal acts of the Russian Federation.

Notice to the user

If, during the use of this product or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority

If further Instructions for Use are required, they can be requested from the manufacturer or can be obtained from the B. Braun homepage: <https://eifu.bbraun.com/>.

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Symbol	Meaning of symbol
	Do not re-use
	Caution

	Consult instructions for use
	Batch number
	Sterilization by ethylene oxide
	Use-by date
	Date of manufacture
	DEHP-free
	Latex-free
	Prime-Stop
	PVC-free
	20 drops per ml
	Infusion Filter
	UV-protect
	Air-Stop
	Do not use if package is damaged
	tubing length, cm
	Filling/storage volume per metre tube length
	Non-pyrogenic
	Pressure
	Gravity
	Catalog number
	The mark of conformity to the Declaration of conformity medical devices in Russian Federation (if necessary)

	Manufacturer
	Single sterile barrier system
	Single sterile barrier system with protective packaging outside
	Country of manufacture
	Medical Device
	«Keep dry»
	«Keep away from sunlight»

CE mark

CE mark – the conformity mark of Directive 93/42/EEC and the identification number of the accredited certification body in the European Union

Manufacturer

Manufacturer: B. Braun Melsungen AG, Carl-Braun-Str. 1, 34212 Melsungen, Germany.

www.bbraun.com

Date of issue and if appl. revision date

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Chairwoman of the Supervisory Board:
Dr. Annette Belier

Executive Board:
Markus Strotmann
(Chairman)
Priv.-Doz. Dr. Stefan Ruppert
Jurgen Stahl

Corporate Office: Melsungen
Register Court:
Local Court Fritzlar
HRB 11 000
WEEE-Reg.-No. DE 42690900

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
13. Infusomat® Plus Line SafeSet for infusion by Infusomat® compact plus pump, infusion, flushing (art.№: 8700280) 25 pc./box, with instruction for use (1 pc.)

(further in the text « Infusomat® Plus Line SafeSet, infusion, flushing»)

Intended Use

Intended for continuous or intermittent IV administration by a compatible pump (primary method) and by gravity (secondary method), for flushing (emptying for underdosage prevention). Intended for IV administration of infusion solutions, liquid drugs and fat emulsions, as well as low-soluble drugs and cytotoxic drugs.

Device description


1 protective cap (preserves the sterility of the medical device until it is used) and spike (provides easy piercing of various ports of infusion containers); 2 - Snap clamp - allows you to block the flexible tube to stop the flow of the solution if necessary. 3 - Y-port needle-free (2 pcs.) - it is intended for additional administration of drugs. 4 - Drip chamber, which consists of: 4.1. moulding ring - provides a tight connection between the upper and lower parts of the drip chamber; 4.2 the integrated dropper into the upper part of the drip chamber has standard characteristics: 20 drops = 1 ml + / - 0.1 ml (distilled water); 4.3 the lower part of the drip chamber is elastic, so it can be used to pump the solution into the drip chamber; 4.4 liquid filter - 15 µm fluid filter Air-Stop (hydrophilic membrane, impermeable to air bubbles) is integrated in the lower part of the drip chamber. 5 - Tube - is designed for the introduction of solutions to the patient, transparent (PUR). 6 - Roller clamp, has a connector for attaching a tube that can be used to temporarily interrupt the infusion and a for safely disposing of the tip after the infusion is completed. 7 - Peristaltic (pump) segment it is made of silicone and when exposed to the peristaltic mechanism of the pump, it provides a uniform unidirectional movement of the injected solution at a given speed. For the correct placement of the line inside the pump, it is provided: (clip (inlet/outlet) and places inside the pump, having the appropriate shape and color coding and graphic image ("asterisks") to prevent twisting of the tube during the installation of the main line in the pump.). 8 - Free Flow Protection Clamp - allows to block the flexible tube to stop the flow of the solution if necessary. 9 - Connector Luer Lock (type male) with protective cap - Connector Luer Lock ensures tightness of the connection between the line and the patient's catheter (due to the presence of internal thread), protective cap with a hydrophobic membrane Prime-Stop - the presence of the PrimeStop hydrophobic membrane allows you to fill the line without removing the cap and leave the connector sterile until the line is connected to the patient.

Materials used:

PS-HI, PE-LD, SBC/SBS, PP, ABS, SIR, , MABS, PUR, PC, PSU/PVP, Acrylic copolymer, SIL, PE

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1. Choose an Infusomat® Plus Line SafeSet infusion, Flushing set, according to the planned procedure.
2. Check the expiry date and integrity of the package. Do not use if the package is damaged or the shelf life is expired.
3. Check the name, expiry date, physical properties and dosage of the drug.
4. Open the individual package, take out and unfold the Infusomat® Plus Line SafeSet infusion, Flushing set, observing the aseptic regulations.
5. Close the roller clamp (fig. 1).
6. Disinfect the container port, remove the protective cap from the spike and puncture the port (Fig. 2).
7. While squeezing the bottom part of the drip chamber, fill it with sufficient amount of solution (Figure 3).
8. Ventilation cap should be closed when using plastic collapsible containers or bags and opened during infusion from glass vials. (Fig. 4).
9. Open the roller clamp and prime the line with solution, then close the roller clamp (fig. 5).
10. Install the Infusomat® Plus Line SafeSet infusion, Flushing set according to the instructions for the peristaltic pump. Install the green snap clamp in the appropriate pump slot (Fig. 6 pos. 2), and the input and output retainers of the peristaltic segment to the corresponding places (Fig. 6 pos. 1).
11. Push the Infusion line into the air sensor on the left and into the pump based safety clamp on the right (Fig. 7).
12. After installing the line, close the pump door while simultaneously pulling the door level and start the administration (fig. 8).
13. Connect the line to the patient's catheter (fig. 9).
14. Open the roller clamp (fig. 10).
15. Set infusion parameters on the pump and start infusion.
16. Before removing the line from the pump, close the roller clamp (Fig. 11), disconnect the line from the patient's catheter (Fig. 12).
17. Open the pump door by pulling the door lever and remove the line from the pump (Fig. 13).

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Intended purpose

Single-use, sterile infusion systems for administration of infusion solutions, parenteral nutrition via one single closed line-system, usable with a compatible pump (primary method) or by gravity (secondary method).
Do not re-use.

Indications

For standard infusion solutions and for low soluble IV-solutions, parenteral infusion solutions (parenteral infusion solutions and lipids) and for cytostatic solutions.

Patient population

The device can be used for all patients for which infusion therapy is prescribed. Infusomat® Plus Line SafeSet, infusion, flushing can be used for adults, pediatric and neonates in hospitals and also in the outpatient market.
Selection of pump and dedicated set according to national standards and / or hospital protocols

Intended User

Infusomat® Plus Line SafeSet, infusion, flushing are intended to be used by trained healthcare professionals, licensed and/or certified physicians, healthcare practitioners and/or specifically trained patients or caregivers according to local regulations

Contraindications

Infusomat® Plus Line SafeSet, infusion, flushing shall not be used in case of proven incompatibility of the drug/IV solution to one of the used materials. SPCs of the drugs/solutions for contraindications should be considered.
No administration of blood products, blood components and enteral nutrition solutions.

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- Air embolism
- Local to systemic infection
- Medication error including over and under supply
- Particle contamination
- Kinking of the line may result in decrease or stop of infusion flow
- Leakage of drugs/fluids can be an additional hazard to patients, visitors or staff

Precautions

- Visual check of the Infusomat® Lines for damages before use.
- Do not re-attach the protective cap of the patient connector.
- Always refer to the manufacture's instruction for use of infusion pumps.
- If pressure infusion without infusion pump is used, remove air from container beforehand. Observe general safety precautions.
- Mind the general guidelines for parallel infusions.
- Please also consider the corresponding SPCs of the drugs/ solutions regarding possible incompatibilities of drugs/ solutions and/or material of the used Infusomat® Line.
- Monitor all connections for tightness over the whole infusion procedure
- Due to the structure of the fluid filter of Infusomat® SafeSet which achieves a particle retention with a filter efficiency > 80% for a particle size of 3µm, the infusion of highly viscous solutions, lipid emulsions >10% and certain suspensions can affect the flow rate. With regard to this, please also follow the information in the relevant summary of product characteristics

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Always check, if incompatibilities with other drugs or with respective material of different types of Infusomat® Lines is stated in the SPC (Summary of product characteristics) of the drug. Always check, if the use of a 0.2µm filter is stated in the SPC (Summary of product characteristics) of the drug.

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Length of device, mm (± 150 mm) / Length of line from peristaltic segment to Luer Lock connector, mm (± 100 mm) – 3000 mm / 2000 mm.

Internal/outer tube diameter - 4.1 x 3.0 mm ($\pm 0,2$ mm)

Filling volume (ml) on 1 m of polymer tubing under 40 °C in accordance to EN ISO 8536-8 / Storage volume (ml) on 1 m of polymer tubing under 40 °C/2 bar in accordance to EN ISO 8536-8: 6 ml/m 40 °C /7 ml/m 40 °C/2 bar.

Resistance to pressure up to 2.0 bar.

Luer compatibility acc. to EN ISO 80369-7/ EN ISO 80369-20

Spike compatibility acc. to EN ISO 8536-4/ EN ISO 8536-8

Scope of supply

Scope of supply of the medical device Infusomat® Plus Line SafeSet, infusion, flushing is a group package that includes: medical device in individual packaging - 25 pcs.

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Symbol	Meaning of symbol
	Do not re-use
	Caution
	Consult instructions for use

LOT	Batch number
STERILE EO	Sterilization by ethylene oxide
	Use-by date
	Date of manufacture
	DEHP-free
	Latex-free
	PVC-free
	20 drops per ml
	Infusion Filter
	Pressure
	Gravity
	Do not use if package is damaged
	tubing length, cm
VOL	Filling/storage volume per metre tube length
	Non-pyrogenic
Prime Stop	Prime-Stop
Air Stop	Air-Stop
REF	Catalog number
	The mark of conformity to the Declaration of conformity medical devices in Russian Federation
	Manufacturer
	Single sterile barrier system
	Single sterile barrier system with protective packaging outside
	Country of manufacture

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Instructions for Use (specific version of instructions for use (IFU) in English and Russian languages)

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Executive Board:
Markus Strotmann
(Chairman)
Priv.-Doz. Dr. Stefan Ruppert
Jurgen Stihl

Corporate Office: Melsungen
Register Court
Local Court Fritziar
HRB 11 000
WEEE-Reg.-No DE 42690900

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15. Infusomat® Plus Line for infusion by Infusomat® compact plus pump, transfusion (art.№: 8700360) 100 pc./box, with instruction for use (1 pc.)


(further in the text « Infusomat® Plus Line, transfusion»)

Intended Use


Intended for continuous or intermittent IV administration by a compatible pump (primary method) and by gravity (secondary method).

Intended for IV administration of blood and blood products only (transfusion).

Device description

8700350-20, 8700350, 8700350-26


8700360

1 – Drip chamber, which consists of: 1.1 protective cap (preserves the sterility of the medical device until it is used); 1.2 spike (provides easy piercing of various ports of blood and blood products containers); 1.3 ventilation inlet with lid and air filter - provides additional protection against contamination; 1.4 moulding ring - provides a tight connection between the upper and lower parts of the drip chamber; 1.5 the integrated dropper into the upper part of the drip chamber has standard characteristics: 20 drops = 1 ml + / - 0.1 ml (distilled water); 1.6 the lower part of the drip chamber is elastic, so it can be used to pump the solution into the drip chamber; 1.7 filter for blood (200 µm) - is integrated in the lower part of the drip chamber. 2 – Tube - is designed for the introduction of solutions to the patient, transparent (PVC). 3 – Roller clamp, has a connector for attaching a tube that can be used to temporarily interrupt the infusion and a for safely disposing of the tip after the infusion is completed. 4 - Peristaltic (pump) segment it is made of silicone and when exposed to the peristaltic mechanism of the pump, it provides a uniform unidirectional movement of the injected solution at a given speed. For the correct placement of the line inside the pump, it is provided: (clip (inlet/outlet) and places inside the pump, having the appropriate shape and color coding and graphic image ("asterisks") to prevent twisting of the tube during the installation of the main line in the pump.). 5 - Free Flow Protection Clamp - allows to block the flexible tube to stop the flow of the solution if necessary. 6 - Y-port needle-free - it is intended for additional bolus administration of drugs. 7 - Connector Luer Lock (type male) with protective cap - Connector Luer Lock ensures tightness of the connection between the line and the patient's catheter (due to the presence of internal thread), protective cap - equipped with a hydrophobic membrane Prime-Stop. The presence of the PrimeStop hydrophobic membrane allows you to fill the line without removing the cap and leave the connector sterile until the line is connected to the patient.

Infusomat® Plus Line, transfusion can include (depending on the variants): ABS, Acrylic copolymer/Nonwoven nylon, Borosilicate, HIPS, MABS, PSU/PVP, PA, PC, PE, PE-LD, PP, PS, PS-HI, PVC, SBC/SBS, SIL, SIR

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1. Select Infusomat® Plus Line, transfusion in accordance with the planned procedure.
2. Check the shelf life and integrity of the trunk packaging. Do not use if the packaging is damaged or after the expiration date.
3. Check the name, expiration date, physical properties, and dosage of the blood or blood product.
4. Open the individual packaging, remove and expand the Infusomat® Plus Line, transfusion observing the rules of asepsis.
5. Close the roller clamp (Fig. 1).
6. Treat the container port, remove the protective cap from the tip and puncture the port (Fig. 2).
7. Squeezing the lower part of drip chamber, fill it with sufficient amount of the blood or blood product (Fig. 3).
8. During infusion of the polymer bottles or bags, the lid of the ventilation duct of the tip should be closed and during infusion of glass bottles the ventilation channel should be open.
9. Open the roller clamp and fill the line with a blood or blood product (Fig. 4).
10. Then close the roller clamp (Fig. 5).
11. Install the Infusomat® Plus Line according to the instructions for the peristaltic pump. Install the green snap clamp in the appropriate pump slot (Fig. 6 pos. 2), and the input and output retainers of the peristaltic segment to the corresponding places (Fig. 6 pos. 1).
12. Push the Infusion line into the air sensor on the left and into the pump based safety clamp on the right (Fig. 7).
13. After installing the line, close the pump door while simultaneously pulling the door level and start the administration (Fig. 8).
14. Connect the line to the patient's catheter (Fig. 9).
15. Open the roller clamp (Fig. 10).
16. Set infusion parameters on the pump and start infusion.
17. Before removing the line from the pump, close the roller clamp (Fig. 11), disconnect the line from the patient's catheter (Fig. 12).
18. Open the pump door by pulling the door lever and remove the line from the pump.

Compatible infusion pumps according to the manufacturer's instruction for use of infusion pumps. Infusomat® plus Lines type Transfusion are compatible with Infusomat compact plus (Registration Certificate № RZN 2022/16352).

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Follow points from 1 to 10 from the chapter "Operating Instructions for infusion by compatible pump".

After that:

11. Connect the line to the patient's catheter (Fig. 9).
12. Set the required infusion rate by adjusting the roller clamp. Periodically monitor the flow rate throughout the duration of the infusion.
13. After infusion is finished, close the roller clamp (Fig. 11) and disconnect the line from the patient's catheter (Fig. 12).

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Intended purpose

Infusomat® Plus Line, transfusion is intended to be used to administer blood or blood products (whole blood, red blood cells, platelet concentrates, human albumin, granulocytes, plasma and plasma components (fresh frozen plasma, stable plasma derivatives, cryoprecipitate)), stored in flexible blood bags or rigid containers directly into a patient's circulatory system. Usable with a compatible pump (primary method) or by gravity (secondary method).

Indications

Transfusions should be administered as clinically indicated for patients with acute or chronic blood loss and those who have objective signs and symptoms of thrombocytopenia, blood coagulation disorders and / or anemia.

Transfuse blood and blood products in accordance with evidence-based indications to ensure patient safety, optimal patient outcomes, and eliminate un-necessary transfusions.

Patient population

The device can be used for all patients for which transfusion therapy is prescribed. Infusomat® Plus Lines can be used for adults, pediatric and neonates in hospitals.

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Intended User

To be used by country-specific authorized person's only (medically qualified personnel according to national regulations, especially trained in the use of the respective device, by order of a physician).

Contraindications

Infusomat® Plus Line, transfusion is not interchangeable with infusion administration sets. If specific blood sets should not be used, an additional risk to the patient, ranging from a lack of the required filtration leading to an increased risk of hemolysis and under-transfusion, is possible.

No other solutions or medications should be added to or infused through the same administration set with blood or blood products (except flushing solution, e.g. NaCl 0.9% after transfusion).

Residual Risks and Side effects

General administration set-related complications are local and systemic infection, leakage, particle contamination, air embolism and kinking that may result in decrease or stop of transfusion flow.

General risks of transfusion therapy include immuno-logic adverse reactions such as intravascular and extravascular hemolysis, febrile non-hemolytic complications, urticarial, allergic reactions or anaphylaxis as well as non-immunologic adverse reactions such as contamination, mechanical hemolysis and bacterial infections.

Many reactions become evident within 5 to 15 minutes of initiating the transfusion, but reactions can develop any time during and after a transfusion. For this reason, patients should be closely monitored during and after transfusion. Stop the transfusion immediately if signs of a reaction develop.

Precautions

- Visual check of the Infusomat® Plus Line for damages before use.
- Do not re-attach the protective cap of the patient connector.
- Always refer to the manufacture's instruction for use of infusion pumps.
- If pressure transfusion without infusion pump is used, remove air from container beforehand. Observe general safety precautions.
- Monitor all connections for tightness over the whole transfusion procedure

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Do not re-use: Re-use of single-use devices creates a potential risk of patient or user. It may lead to contamination and/or impairment of functional capability. Contamination and/or limited functionality of the device may lead to injury, illness or death of the patient.

Duration of use

In general Infusomat® Plus Line, transfusion should be changed acc. to national standards and / or hospital protocols, but at least every 6 hours.

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Article number	Length of device, mm ($\pm$ 150 mm) / Length of line from peristaltic segment to Luer Lock connector, mm ($\pm$ 100 mm)	Filling volume (ml) on 1 m of polymer tubing under 40 °C in accordance to EN ISO 8536-8	Storage volume (ml) on 1 m of polymer tubing under 40 °C/2 bar in accordance to EN ISO 8536-8	Tube (material colour)	filter for blood (200 μ m)	Y-port (needle-free)
8700350, 8700350-20, 8700350-26	2400 / 1500	6	7	PVC, colorless	+	-
8700360	3000 / 2000	6	7	PVC, colorless	+	+

Internal/outer tube diameter - 4.1 x 3.0 mm ($\pm$ 0,2 mm)

Resistance to pressure up to 2.0 bar.

Luer compatibility acc. to EN ISO 80369-7/ EN ISO 80369-20

Spike compatibility acc. to EN ISO 8536-4/ EN ISO 8536-8.

Scope of supply

Scope of supply of the medical device Infusomat® Plus Line, transfusion for using with compatible pumps and infusion by gravity" is a group package that includes medical device in individual packaging - 100 pcs.

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Operation environment, Storage and transportation conditions

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Disposal

Disposal acc. to local guidelines and or clinical protocols. Use integrated spike protection at roller clamp for disposal. Disposal of the medical device «Infusomat® Plus Line» after use is carried out in accordance with the requirements for the treatment of medical waste established by SanPiN 2.1.3684-21 for class B and/or other applicable legal acts of the Russian Federation.

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Marking

Marking on the medical device: On the top part of the drip chamber can be applied the logo of the manufacturer (in Latin «B|Braun»).

Symbols are used in accordance with GOST R ISO 15223 for marking of the individual and group packaging, for the decoding of the symbols see the Table below. Additionally, information in a foreign language can be applied on the marking, namely: «Type Transfusion - Luer Lock» (*translation from English: «Тип Трансфузионная – Люэр Лок»*); «Transfusion» (*translation from English: «Трансфузионная»*); «needle free» (*translation from English: «инъекционный порт, безыгольный»*).

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	Consult instructions for use
LOT	Batch number
STERILE EO	Sterilization by ethylene oxide
	Use-by date
	Date of manufacture
	DEHP-free
	Latex-free
Prime Stop	Prime-Stop
	20 drops per ml
	Do not use if package is damaged
	tubing length, cm
VOL	Filling/storage volume per metre tube length
	Infusion Filter
	Non-pyrogenic
	Pressure
REF	Catalog number
	The mark of conformity to the Declaration of conformity medical devices in Russian Federation
	Manufacturer
	Single sterile barrier system
	Single sterile barrier system with protective packaging outside
	Country of manufacture
MD	Medical Device
	«Keep dry»
	«Keep away from sunlight»

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