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
www.maqu

ions For Use

Flixene® ePTFE vascular grafts are offered in 2 designs; the Flixene Vascular graft. They are made entirely of polytetrafluoroethylene (ePTFE), and are available in a wide variety of sizes.

referred to as "helix") for added radial support (made from a multi-layered ePTFE filament coil).

Deployment System consisting of a pre-attached tip(s) and a flexible sheath (polyethylene).

the tapered end and the wall thickness.

The graft is a 2 layer graft employing a single layer ePTFE graft, lined with an additional layer of ePTFE for increased support. On some Advanta VXT grafts are bonded to the exterior of the graft. By following the prescribed procedure, the rings can be removed as needed.

The Advanta VXT graft is a 3 layer graft comprised of ePTFE. This graft has an additional layer of ePTFE making it Maquet's most flexible. As a result of its design, ring support is not necessary for use in the aorta or across joints.

Use

Flixene ePTFE vascular grafts are intended for use in arterial bypass, segmental bypass, and for arteriovenous vascular access.

ns

Advanta VXT, and Flixene GW grafts are contraindicated for use in the aorta, Axillofemoral bypass, or any anatomical location where a wide range of body motion. Also, tapered grafts are not recommended for use in the aorta, Axillofemoral or Axillofemoral bypass, or any anatomical location where a wide range of body motion, due to the danger of disruption during extreme body movement.

The FDA (USA) restricts this device to sale by or on the order of a qualified medical professional.

Do not use if packaging or device appear to be damaged.

Do not use any Advanta VXT or Flixene ePTFE grafts.

Do not use point suture needles as taper cut and other cutting needles on the graft material and/or host vessel.

Excessive tension to the anastomosis as elongated needle punctures or material disruption may cause needle hole bleeding, or host vessel damage.

The graft should be avoided whenever possible and limited to use with soft material.

Do not treat the Advanta VXT or Flixene ePTFE grafts by exposure to solutions such as alcohol or force aqueous solutions through the graft. The hydrophobic properties of the graft may be affected, which may result in excessive serum leakage and/or peri-graft seroma formation. Follow all handling and operative technique warnings and instructions. Do not use in extreme blood loss, loss of limb function, loss of limb, or death.

Do not peel, separate or delaminate any portion of the multi-laminate structure or the graft.

Do not use clinical data available at this time to support the use of Flixene ePTFE graft use in aortocoronary bypass or for use as a patch.

The Graft Deployment system (GDS) tip(s), Slider sheath(s), Transfer sleeve must be removed from the graft and patient after closure.

Precautions

1. Handling of the graft should be minimized and limited to clean sterile gloves.
2. Tunnel forming instruments should be sized to closely match the graft diameter to prevent formation of an oversized tunnel. An oversized tunnel may be a contributing factor in peri-graft seroma formation, and may delay or interfere with peri-graft tissue attachment and healing. When tunneling a graft which has external ring support, aggressive graft manipulation or placement within an undersized tunnel may lead to separation of the peelable external ring support.
 - a. Note: Correct sizing of the Vascular Graft Tunneler Tip to the graft is outlined in the "Bypass Technique" section of this IFU and in the Maquet Vascular Graft Tunneler Instrumentation IFU.
3. Dialysis cannulation through this graft may be instituted with caution, so long as safeguards are taken to prevent or minimize hematoma formation, pseudoaneurysm, infection and/or material disruption requiring surgical intervention.
4. Vascular access puncture sites must be adequately separated as multiple punctures in the same area may lead to disruption of the graft material and formation of a peri-graft hematoma or pseudoaneurysm requiring surgical intervention.
5. Advanta VXT and Flixene grafts should only be cut and trimmed with sharp surgical instruments to avoid reinforcement layer disruption.
6. Do not use any laser, electric, radiofrequency, or heat cautery for tailoring or cutting Advanta VXT or Flixene ePTFE graft material.
7. User should be a physician familiar with vascular graft surgical techniques and procedures before implanting a vascular graft.
8. Grafts must be fully extended prior to sizing length for implantation.
9. Upon placement of the graft, ensure the blue line on the side of the graft is straight and not twisted. A twisted line may represent a kink in how the graft is placed.

Adverse Reactions

Complications that may occur in connection with the use of any vascular graft include, but are not limited to; thrombosis, stenosis, formation of pseudoaneurysm due to excessive needle punctures, peri-graft hematoma formation, peri-graft seroma formation, excessive needle hole bleeding or weeping, infection, swelling of tissue, suture hole elongation, mechanical disruption, material separation, delamination or tearing of the graft material, suture line or host vessel which may result in extreme blood loss, loss of limb function, steal syndrome, loss of limb or possible death. The patient should be advised to contact the physician should an adverse reaction occur.

Pre-operative Device Preparation

The Advanta VXT and Flixene ePTFE vascular graft packaging consists of; a shelf box (made of fiberboard), outer tray, and inner tray which houses the graft. The inner tray and graft, housed within the outer tray, are considered sterile.

The procedure for removing the Advanta VXT or Flixene ePTFE graft for clinical use from its packaging is as follows:

1. Ensure staff member confirms the correct graft has been selected with the requesting physician.
2. Inspect both the box and the outer tray for signs of damage and that product expiration date has not been exceeded.

3. Device tracking labels are provided for hospital and patient keeping. The implant tracking card should also have a label and be returned to Maquet for processing.
4. Confirm the contents listed on the shelf box label with the tray label.
5. Upon confirming that the correct device has been selected, transfer the device into the sterile field using standard hospital technique. A non-sterile operator should peel open the outer tray exposing the inner tray so the technician within the sterile field may view the outer package.
6. The technician managing sterile product should peel the product and remove the graft once the physician is ready to implant. The grafts will be covered by a blue transfer sleeve. It is recommended the blue transfer sleeve remain on the graft until the graft is in the subcutaneous plane for implantation. The purpose of the sleeve is to minimize graft/sheath contact with foreign material.
7. It is highly recommended that graft contact with other material be minimized in an effort to minimize any potential transfer of foreign material to the operative site.
8. Do not pre-wet graft prior to placement. Inspect the graft for defects.

Standard Bypass Surgical Technique and Tunneling

1. The patient is prepped, anesthetized, and draped according to standard surgical practices. The physician should determine the site for the anastomosis sites, and dissect down to the target vessel using standard surgical technique.
2. The physician must next measure the distance of the bypass taking into account patient weight, posture, and range of motion at the placement location. Note: These factors may warrant using a longer length. It is important to make sure that the graft is long enough to create tension between the anastomosis sites.
3. Select appropriate graft diameter based upon the native vessel.
 - a. Note: If a segment of the graft is used to span a gap in the body where radial compression of the graft is required, the selected graft should have "ring" or "helix" support in the area to avoid the potential for the graft to compress when the vessel is compressed, avoiding a reduction in blood flow.
4. Once the physician has confirmed that the appropriate graft is available, he can select the appropriate Vascular Graft Tunneler Tip. To create the most appropriate size tunnel for the implantation, the following guidelines for tip selection:
 - a. For regular, non-externally supported vascular grafts and similar products, use the same size as the internal diameter of the graft.
 - b. For ringed vascular grafts or thick walled products, use the next larger size than the internal diameter of the graft.
 - c. For tapered vascular grafts, consider the larger diameter of the tapered graft for tip size selection.
5. The graft may be supplied with or without the Slider Graft Deployment System (GDS):
 - a. If the graft is not supplied with a GDS, connect the Vascular Graft Tunneler Tip by placing a suture in the Vascular Graft Tunneler Tip to the graft.
 - b. If the graft is supplied with the Slider GDS component, remove the tip from the rod of the GDS and insert the GDS into the distal end of the tunneler rod.

The Vascular Graft Tunneler Rod is 6 mm in diameter. Select the appropriate Vascular Graft Tunneler Rod model based upon anatomical location, vein required.

Insert the appropriate Vascular Graft Tunneler rod and tip, through the proximal incision site into the subcutaneous tissue. Once the tip has been passed through the tissue plane to the vein, remove the Vascular Graft Tunneler Tip by unscrewing it. Thread the GDS unit on the end of the graft into the end of the Vascular Graft Tunneler Rod, replacing the now removed Vascular Tip.

Next, pull the Vascular Graft Tunneler Rod back through the vein in a continuous and steady motion until the graft material reaches the proximal incision site. Retain the blue transfer sleeve covering the graft is pulled through, leaving it outside the body. Avoid sudden and quick pulls as this may cause the graft to be damaged. Once the graft is in position, confirm that the blue reference line on the graft is straight. A significant twist in the line represents a potential kink in the graft placement. If necessary the graft may have to be

repositioned. Next, cut the clear sheath closest to the GDS. Remove and discard the GDS. The graft may be cut using sharp surgical scissors or a scalpel. It is strongly recommended that the graft be cut with scissors to provide an even cut. Unless absolutely necessary, it is recommended not to cut until it has been tunneled and is in place, ensuring that the graft is long enough.

Next, secure the graft at the proximal incisional end using a pair of shod forceps. Make sure that when securing the graft with the forceps, the clear sheath is not secured as well. Using gentle steady pressure, pull the graft from the distal incisional site and withdraw it from the graft/patient.

Next, must now clamp off the artery/vein (physician's prerogative on the venous or arterial side) above and below the intended anastomosis. The size and type of incision will be determined by the physician based on his experience and training. An end to side anastomosis is commonly performed for arteriovenous bridge fistula

Next, using a permanent suture material (such as monofilament suture) with a non-cutting needle. Sutures should be placed 1 mm apart, and 1 mm from the edges when suturing to

the anastomosis are completed, remove the clamps and confirm the anastomosis. Should suture hole bleeding occur, apply gentle pressure with sterile gauze for 5-10 minutes in an effort to facilitate clotting. In the event that more significant bleeding occurs, additional pressure or a coagulant may be needed. The site must not be bleeding has subsided.

Next, close the tissue surgical sites using standard surgical technique. The physician should be educated by their physician on appropriate care.

Loop Graft Placement

See "Preoperative Device Preparation Steps." and "Standard Bypass Surgical Technique and Tunneling Procedure"

Once the tunnel is created, pull the graft through, making sure that the graft length exits the distal incision site.

Next, use surgical scissors carefully cut the clear sheath in the mid portion so that each end of the sheath can be removed from the proximal incision site. Do not cut any portion of the graft when cutting or separating the sheath in the mid portion of the

- Replace the appropriate Vascular Graft Tunneler Tip onto the end of the Vascular Graft Tunneler Rod and begin to tunnel from the proximal incision site to the distal incision site. Make sure that the tunnel will allow the graft to be placed in a loop without kinking.
- Remove the Vascular Graft Tunneler Tip from the rod and attach the GDS to the end of the Vascular Graft Tunneler rod. Using steady pressure, pull the graft back through the tissue plane.
- Once the graft is in position, confirm that the blue reference line on the side of the graft is straight. A significant twist in the line represents a potential kink in the graft placement. If necessary the graft may have to be repositioned.
- Next, cut the graft and clear sheath closest to the GDS. Remove and discard the GDS. The graft may be cut using sharp surgical scissors or a scalpel. It is strongly recommended that the graft be cut with scissors to provide an even cut. If a double ended GDS has been used, remove the other GDS unit at this time.
- Gently secure the graft at the middle of the loop using a pair of shod forceps. Make sure that when securing the graft with the forceps, the clear sheath is not secured as well. Using a clamp to secure either the proximal end of the clear sheath, pull on the sheath to remove it from the graft/patient. Repeat this procedure and remove the other half of the clear plastic sheath.
- Follow "Standard Bypass Surgical Technique and Tunneling Procedure" steps 11-15.

The above procedural overview is intended for informational purposes only.

Instructions on Removal of External "Rings" or "Helix" from the Advanta VXT Vascular Graft

While the rings are added to the Advanta VXT to improve kink, compression and torque resistance, a physician may choose to remove portions of it for surgical/technical reasons. The steps to do so are as follows:

- Hold the graft flat (horizontal).
- Take the very distal portion of the PTFE ring with a pair of forceps and slowly pull off the ring at a 45 degree angle, being mindful not to catch the outer soft wrap.
- If the external support rings are pulled too quickly, it is possible for a small portion of the VXT outer wrap to be removed. Should the outer wrap fray, the physician is still left with the base layer, which will be sufficient to complete the procedure.

Thrombectomy

In the event that an incision needs to be made to the Advanta VXT or Flixene graft (e.g. for the purpose of providing an access site for declotting the distal anastomosis) post implantation, the following is recommended:

- If and when a longitudinal incision is made for a thrombectomy, it is recommended to place a stay stitch at each end of the longitudinal incision prior to balloon catheter use.
- A horizontal mattress suture technique should be used for closure of the incision site of the graft.
- Do not apply excessive tension on the implanted Advanta VXT or Flixene ePTFE graft during inflated balloon catheter withdrawal. Over inflation of an embolectomy balloon or use of an inappropriate sized catheter may cause damage to the graft and anastomosis.
- Extreme care must be exercised with the use of all mechanical water jet thrombectomy and mechanical rotary brush devices so as to not disrupt or damage the Advanta VXT or Flixene ePTFE graft material, prior vascular access holes and/or suture line.

Sterilization

The Advanta VXT and Flixene ePTFE grafts are provided sterile and intended for single use only. Advanta VXT and Flixene ePTFE grafts are sterilized via Ethylene Oxide (EO) Gas Sterilization methods. Flixene ePTFE grafts are not intended for re-sterilization.

Storage Conditions

Advanta VXT and Flixene ePTFE grafts should be stored in a cool, dry place, not in direct sunlight and do not subject to excessive heat.

Symbols Used On Product Labels

REF CATALOGUE NUMBER SN SERIAL NUMBER

STERILE EO STERILIZED USING ETHYLENE OXIDE CAUTION

CONSULT INSTRUCTIONS FOR USE DO NOT RE-STERILIZE

DO NOT RESTERILIZE USE-BY DATE KEEP DRY

DO NOT USE IF PACKAGE IS DAMAGED MANUFACTURED

QTY BOX QUANTITY PEEL HERE TO OPEN

KEEP AWAY FROM SUNLIGHT

Rx Only CAUTION: FEDERAL LAW (USA) RESTRICTS THIS SALE BY OR ON THE ORDER OF A PHYSICIAN

CE MPP AUTHORIZED REPRESENTATIVE IN THE EUROPEAN

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/ (USA) restricts this device to sale by or on the order of a physician.
 ne are manufactured by Atrium Medical Corporation / 40 Continental Blvd., Merrimack, NH 03054 / +1-603-880-1433
 rk and/or Registered Trademark of MAQUET GMBH, its subsidiaries, and/or affiliates in the United States and/or other countries. Maquet is registered in United States Patent and Trademark Office. Advanta VXT and Flixene are Trademarks and/or Re-
 rp., its subsidiaries, and/or affiliates in the United States and/or other countries. Advanta VXT and Flixene are registered in United States Patent and Trademark Office.
 product names are the trademarks or registered trademarks of their respective holders. Copyright 2019 Atrium Medical Corp. or its affiliates. All rights not expressly granted are reserved.
 No. 011456

AW0114

Atrium Medical Corporation
 40 Continental Boulevard
 Merrimack, New Hampshire 03054 USA
 Phone: +1-603-880-1433

 MAQUET Critical Care AB
 Röntgenvägen 2
 SE-171 54, Solna
 Sweden

Getinge Australia Pty Ltd
 Level 2, 4 Talavera Road
 Macquarie Park NSW 2113
 Australia
 Phone: +61-(0)-2-8874-3100

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- 1.52. ADVANTA VXT, 8X70, 1GDS, FH, STR-SW;
- 1.53. ADVANTA VXT, 8X70, 1GDS, FH, STR-TW;
- 1.54. ADVANTA VXT, 8X80, 1GDS, FH, STR-TW;
- 1.55. ADVANTA VXT, 10X80, NGDS, FH, STR-SW;
- 1.56. ADVANTA VXT, 4-6X45, 1GDS, NH, STPR-SW;
- 1.57. ADVANTA VXT, 4-7X45, 1GDS, NH, STPR-SW;
- 1.58. ADVANTA VXT, 4-7X70, 1GDS, NH, STPR-SW;
- 1.59. ADVANTA VXT, 4-7X45, 2GDS, NH, STPR-SW;
- 1.60. ADVANTA VXT, 5-8X45, 1GDS, NH, STPR-SW;
- 1.61. ADVANTA VXT, 4-6X45, 1GDS, CH, STPR-SW;
- 1.62. ADVANTA VXT, 4-7X45, 1GDS, CH, STPR-SW;
- 1.63. ADVANTA VXT, 4-7X80, 1GDS, FH, LTPR-TW;
- 1.64. ADVANTA VXT, 5-8X80, 1GDS, FH, LTPR-TW;
- 1.65. FLIXENE, 6X10, NGDS, STR;
- 1.66. FLIXENE, 6X30, 1GDS, GW;
- 1.67. FLIXENE, 6X30, 1GDS, STR;
- 1.68. FLIXENE, 6X40, 1GDS, GW;
- 1.69. FLIXENE, 6X50, 1GDS, STR;
- 1.70. FLIXENE, 6X50, 1GDS, GW;
- 1.71. FLIXENE, 7X10, NGDS, STR;
- 1.72. FLIXENE, 7X30, 1GDS, GW;
- 1.73. FLIXENE, 7X40, 1GDS, GW;
- 1.74. FLIXENE, 7X50, 1GDS, STR;
- 1.75. FLIXENE, 7X50, 1GDS, GW;
- 1.76. FLIXENE, 7X80, 1GDS, GW;
- 1.77. FLIXENE, 8X50, NGDS, STR;
- 1.78. FLIXENE, 4-6X35, 2GDS, GWT-GW;
- 1.79. FLIXENE, 4-6X45, 2GDS, GWT;
- 1.80. FLIXENE, 4-6X45, 2GDS, GWT-GW;
- 1.81. FLIXENE, 4-7X30, 2GDS, GWT-GW;
- 1.82. FLIXENE, 4-7X35, 2GDS, GWT-GW;
- 1.83. FLIXENE, 4-7X45, 2GDS, GWT;
- 1.84. FLIXENE, 4-7X45, 2GDS, GWT-GW.

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