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
To be submitted to the Federal Service for Surveillance in Healthcare (Roszdravnadzor) in Russian Federation.

We, OLYMPUS MEDICAL SYSTEMS CORP., located at 2951 Ishikawa-cho, Hachioji-shi, Tokyo 192-8507, Japan, as the legal manufacturer of the Medical device "Single Use Loop Cutter FS-410U", hereby state that attached documents is true:

1. Instruction manual of Single Use Loop Cutter FS-410U
2. Addendum to Instruction manual of Single Use Loop Cutter FS-410U

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President, OLYMPUS MEDICAL SYSTEMS CORP.

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1. Country: JAPAN

This public document

2. has been signed by NAKASHIMA Yukihiro

3. acting in the capacity of Notary of the Tokyo Legal Affairs Bureau

4. bears the seal/stamp of NAKASHIMA Yukihiro, Notary

Certified .

5. at Tokyo

6. APR. 17. 2023

7. by the Ministry of Foreign Affairs

8. 23- **№ 045960**

9. Seal/stamp:

10. Signature

MAEJIMA Tadashi

For the Minister for Foreign Affairs

Registered No. 161 , 2023

NOTARIAL CERTIFICATE

This is to certify that the attached copy of document exactly corresponds with the original.

Dated this 17th day of April, 2023

NAKASHIMA Yukihiro

NAKASHIMA Yukihiro

NOTARY



7-6 AZUMA-CHO HACHIOJI-SHI TOKYO

Tokyo Legal Affairs Bureau

Single Use Loop Cutter FS-410 FS-410L, FS-410U

EndoTherapy

1 Symbols

Symbol	Description	Symbol	Description
	Refer to instructions.		Single use only
	Use by (expiration date)		Sterilized using ethylene oxide
	Sterilization lot number		Lot number
	Manufacturer		Authorized representative in the European Community
	Keep away from sunlight		Importer (into European Union)
	Do not resterilize		Keep dry
	Not Made with Natural Rubber Latex		Do not use if package is damaged
			Temperature limitation

For US Customers only

For a Symbols Glossary, visit us:

<http://www.olympus-global.com/en/common/pdf/symbolsglossary.pdf>

2 Intended Use

This instrument has been designed to be used with an Olympus endoscope. It is used to cut the residual end of the Olympus Loop used for ligating tissue and a suture within the digestive tract.

3 Instruction Manual

This instruction manual contains essential information on using this instrument safely and effectively. Before use, thoroughly review this manual and the manuals for all equipment which will be used during the procedure and use the instruments as instructed. Keep this and all related instruction manuals in a safe, accessible location. If you have any questions or comments about any information in this manual, contact Olympus.

4 User Qualifications

The operator of this instrument must be a physician or medical personnel under the supervision of a physician and must have received sufficient training in clinical endoscopic technique. This manual, therefore, does not explain or discuss clinical endoscopic procedures.

5 Signal Words

WARNING	Indicates a potentially hazardous situation which, if not avoided, could result in death or serious injury.
CAUTION	Indicates a potentially hazardous situation which, if not avoided, may result in minor or moderate injury. It may also be used to alert against unsafe practices or potential equipment damage.
NOTE	Indicates additional helpful information.


6 Precautions

Follow the warnings given below when handling the instrument. This information is to be supplemented by the warnings described in each section.

WARNING

- Do not use this instrument for any purpose other than its intended use.
- Due to the positioning relationship between the loop and loop hanger, or the suture and loop hanger, this instrument may become unable to be withdrawn from the patient body when the loop or the suture is caught by the cutter storage part. Have a full understanding of the potential of the loop or the suture caught at the distal end before using this instrument.
- If this instrument cannot be withdrawn from the patient body because the loop or the suture is caught by the distal end, refer to Chapter 12, "Emergency Treatment".

7 Nomenclature and Specifications

Nomenclature	Description
1) Slider	Push the slider to open the cutter. Pull to close.
2) Forceps cap	Prevents damage to the distal end during transportation.
3) Loop hanger	Holds the loop or the suture when cutting.
4) Cutter	Cuts the loop or the suture.
5) Cutter storage part	The cutter is stored when the slider is pulled.

8 Specifications

The compatible Olympus endoscopes are listed in the table in Section 8.2, "Specifications". New endoscopes released after the introduction of this instrument may also be compatible for use in combination with this instrument. For further details, contact Olympus.

WARNING

Use this instrument only in combination with products recommended by Olympus as shown in the Tables in Section 8.2, "Specifications". If combined with products not recommended by Olympus, patient or operator injury, malfunction, or equipment damage may result.

8.1 Environment

Operating environment	Ambient temperature	10 – 40°C (50 – 104°F)
	Relative humidity	30 – 85%
	Atmospheric pressure	700 – 1060 hPa (0.7 – 1.1 kgf/cm ²) (10.2 – 15.4 psia)
Standard storage environment (e.g. within the hospital)	Ambient temperature	5 – 40°C (41 – 104°F)
	Relative humidity	10 – 95%
	Atmospheric pressure	700 – 1060 hPa (0.7 – 1.1 kgf/cm ²) (10.2 – 15.4 psia)
Transportation environment (conditions during transportation and short-term storage)	Ambient temperature	–40 to +70°C (–40 to +158°F)
	Relative humidity	10 – 95%
	Atmospheric pressure	700 – 1060 hPa (0.7 – 1.1 kgf/cm ²) (10.2 – 15.4 psia)

8.2 Specifications

Model		FS-410L	FS-410U
Shape of the scissors			
Maximum insertion portion diameter (mm)		ø 2.4	
Working length (mm)		1650	2300
Compatible Olympus endoscopes (All of these parameters should be met.)	Model and length	Working length less than 1200 mm GIF, CF, OGF, OSF	Working length less than 1850 mm GIF, CF, SIF (SIF-10 Only), OGF, OSF, PCF
	Channel inner diameter (mm) (color code)	ø2.8, ø3.2 (yellow) ø3.7, ø6 (orange)	ø2.8, ø3.2 (yellow) ø3.7, ø4.2, ø6 (orange)
Olympus loops		MAJ-254, MAJ-339, MAJ-340, HX-400	
Designated sutures (All of these parameters should be met.)		Monofilament Barbed Absorbable 2-0, 3-0 (USP)	
Medical Devices Directive		 0197 This device complies with the requirements of Directive 93/42/EEC concerning medical devices. Classification: Class I sterile	

9 Storage

WARNING

Do not store the instrument in places where the sterile package will be damaged, wet, or improperly sealed. Otherwise, the sterile condition of the instrument may be compromised and pose an infection control risk or tissue irritation may result.

Store the instrument in the sterile package at room temperature in a clean and dry environment. Do not store in direct sunlight. Ensure that the packaged instrument is not crushed by surrounding objects during storage.

10 Preparation, Inspection and Operation

WARNING

- Do not use an instrument after the expiration date displayed on the sterile package. Doing so may pose an infection control risk or cause tissue irritation.
- Before use, prepare and inspect the instrument as instructed below. Should any irregularity be observed, do not use the instrument; use a spare instead. Damage or irregularity may compromise patient or user safety, such as punctures, hemorrhages, or mucous membrane damage, and may result in more severe equipment damage.
- When inspecting or using the instrument, always wear appropriate personal protective equipment. Otherwise, blood, mucus, and other potentially infectious material from the patient could pose an infection control risk and/or cause skin irritation.

CAUTION

- Do not coil the insertion portion with a diameter of less than 15 cm. This could damage the insertion portion.
- Never use excessive force to open or close the cutter. This could damage the instrument.

10.1 Preparation

■ **Equipment and personal protective equipment**

Prepare all equipment and personal protective equipment that will be used with the instrument in accordance with their respective instruction manuals. Appropriate personal protective equipment may include eyewear, a face mask, moisture-resistant clothing, and chemical-resistant gloves.

■ **Spare instrument**

Always have a spare instrument available.

■ **Emergency equipment**

Always have sharp pliers on hand, to cut the instrument when it cannot be withdrawn from the patient.

10.2 Inspection

Wear the personal protective equipment in accordance with their respective instruction manuals.

Before use, always inspect the instrument according to the following procedures.

If an abnormality in the instrument is detected, use a spare instrument, inspecting it thoroughly before use.

■ Inspection of the sterile package

WARNING

Do not attempt to sterilize the instrument. This could pose an infection control risk, cause tissue irritation, equipment damage, or malfunction.

Inspect the sterile package for tears, inadequate sealing, or water damage. If the sterile package shows any irregularities, the sterile condition of the instrument may have been compromised. Use a spare instead.

■ Appearance inspection

If any of following steps reveals irregularities, do not use the instrument; use a spare instead.

- 1 Unseal the sterile package and take out the instrument.
- 2 Remove the forceps cap from the distal end.
- 3 While operating the slider to open and close the cutter, confirm that the instrument is free from disconnection or looseness.
- 4 Pull the slider to make sure that the cutter closes smoothly and that it is properly aligned with the cutter storage part.
- 5 Confirm that the distal end of the instrument appears exactly as shown in the table in Section 8.2, "Specifications" and is not damaged.
- 6 Make sure that the blade is not broken or nicked.
- 7 Gently run your fingertips over the entire length of the insertion portion to check for any crushed areas, excessive bends, or other damages.
- 8 Make sure that there are no cracks in the handle.

■ Inspection of operation

If the cutter does not operate smoothly and as intended, do not use the instrument; use a spare instead.

- 1 Holding the instrument as shown in Figure 1, form a loop in the insertion portion approximately 20 cm in diameter.
- 2 Operate the slider and confirm that the cutter opens and closes smoothly.


Figure 1

10.3 Operation

The operator of the instrument must be a physician or medical personnel under the supervision of a physician and must have received sufficient training in clinical endoscopic technique.

This manual, therefore, does not explain or discuss clinical endoscopic procedures. It only describes basic operation and precautions related to the operation of this instrument.

WARNING

- Do not insert the instrument into the endoscope unless you have a clear endoscopic field of view. If you cannot see the distal end of the insertion portion in the endoscopic field of view, do not use it. This could cause patient injury, such as punctures, hemorrhages, or mucous membrane damage. It may also damage the endoscope and/or instrument.
- Do not angulate the bending section of the endoscope while the distal end of the insertion portion is extended from the distal end of the endoscope. This could cause patient injury, such as punctures, hemorrhages, or mucous membrane damage.
- Do not use the instrument to cut any objects other than the surplus of Olympus loop (MAJ-254, MAJ-339, MAJ-340, HX-400) or the designated suture. If anything other than the surplus of the loop or the designated suture is cut, the cutter may be damaged and unable to perform properly, the cut object may be caught in the cutter storage part of the instrument, and it may become difficult to safely remove the instrument from the body. Such objects other than the surplus of loop may include stent, stent wire, and loop stoppers. In this case, refer to Chapter 12, "Emergency Treatment".
- This instrument does not contain serviceable parts. Do not disassemble, modify, or attempt to repair it; patient or user injury and/or equipment damage can result.

CAUTION

When using the instrument together with a two-channel endoscope, never use electrosurgical accessories at the same time. This could result in thermal injury to the patient, operator, or assistant.

■ Inserting into the endoscope

WARNING

- When inserting the instrument into the endoscope, hold the slider firmly. Otherwise, the cutter may open and extend from the endoscope tip abruptly. This could cause patient injury, such as punctures, hemorrhages, or mucous membrane damage. It could also damage the endoscope and/or instrument.
- Do not force the instrument if resistance to insertion is encountered. Reduce the angulation until the instrument passes smoothly. This could cause patient injury, such as punctures, hemorrhages, or mucous membrane damage. It could also damage the endoscope and/or instrument.
- Do not advance or extend the instrument abruptly. This could cause patient injury, such as punctures, hemorrhages, or mucous membrane damage. It could also damage the endoscope and/or instrument.

CAUTION

- When inserting the instrument into the endoscope, hold it close to the biopsy valve and keep it as straight as possible relative to the biopsy valve. Otherwise, the insertion portion of the instrument could be damaged.
- Insert the instrument slowly. Abrupt insertion could damage the endoscope and/or instrument.

- 1 Pull the slider to close the cutter.

- 2 With the cutter closed, carefully insert the instrument into the biopsy valve. (See Figure 2)
- 3 Advance the instrument until the distal end of the insertion portion appears within the endoscopic field of view.



Figure 2

■ Cutting the loop or the suture


WARNING

- Do not force the distal end of the insertion portion against body cavity tissue. This could cause patient injury, such as punctures, hemorrhages, or mucous membrane damage.
- Position the loop or the suture on both edges of the loop hanger so as to make it plumb as possible with respect to the distal end and cut it (see Figure 3 (a)). Do not try to cut the loop or a suture that is not positioned on both edges of the loop hanger (see Figure 3 (b)). It may make cutting the loop or the suture impossible, or result in the loop or the suture getting caught in the cutter storage part of the instrument, which could make it difficult or impossible to remove from the patient. (See Figure 3 (b), (c)) If this instrument cannot be withdrawn from the patient body because the loop or the suture is caught by the distal end, refer to Chapter 12, "Emergency Treatment".


(Figures below are when cutting the loop)

(a) Correct position
(ex. loop is positioned on both edges of loop hanger.)

(b) Incorrect position
(ex. loop is positioned on only one edge of loop hanger.)

(c) The loop was caught
(ex. loop is getting caught in the cutter storage part.)

Figure 3

WARNING

- Do not cut the loop or the suture unless you have a clear endoscopic field of view. This could cause the loop or the suture getting caught in the cutter storage part of the instrument, and it may become difficult to safely remove the instrument from the body. It may also cause patient injury, such as punctures, hemorrhages, or mucous membrane damage, or the endoscope and/or instrument damage. If this instrument cannot be withdrawn from the patient body because the loop or the suture is caught by the distal end, refer to Chapter 12, "Emergency Treatment".

CAUTION

If it is difficult to simultaneously cut two loops which is in the condition other than the one shown in Figure 3 (a), cut them one at a time. Forcible cutting may damage the instrument.

- 1 To position the loop or the suture on both edges of the loop hanger angulate the endoscope or insert the instrument until it reaches the target loop or the suture.
- 2 Push the slider to open the cutter.
- 3 Position the loop or the suture on both edges of the loop hanger. (See Figure 3 (a))
- 4 Pull the slider to cut the loop or the suture.
- 5 Push the slider to release the loop or the suture being held.

■ Withdrawing the instrument from the endoscope

WARNING

Do not withdraw the instrument from the endoscope quickly. This could scatter blood, mucous, or other patient debris and pose an infection control risk.

CAUTION

- Do not withdraw the instrument from the endoscope while the cutter is open. This could damage the endoscope and/or instrument.
- If excessive resistance makes withdrawal difficult, adjust the angle of the endoscope until the instrument can be withdrawn smoothly. Forcible withdrawal could damage the instrument and/or endoscope.

- 1 Pull the slider to close the cutter.
- 2 Withdraw the instrument from the endoscope.

11 Disposal

WARNING

- After use, dispose of the instrument in an appropriate manner. If it is not properly disposed of, it could pose an infection control risk, cause user injury.
- The instrument is a single-use, disposable item. Do not reuse or attempt to sterilize it. Reusing the instrument could pose an infection control risk, cause tissue irritation, or malfunction.

After using the instrument, dispose of it in an appropriate manner.



12 Emergency Treatment

WARNING


- If this instrument cannot be withdrawn from the patient body, make a judgment from the specialist's view. Forcible withdrawal could lead to bleeding or mucous membrane damage.
- If you withdraw the endoscope by cutting the insertion portion of this instrument, withdraw it as slowly as possible. Abrupt withdrawal could lead to perforation of the loop ligation position or the suturing part, major bleeding, or mucous membrane damage.





If this instrument cannot be withdrawn from the patient body because the loop or the suture is caught by the distal end, follow the procedures description in this section.

(Figures below are an example of when cutting the loop)


If the loop or the suture cannot be disengaged from the distal end

Insert the endoscope again along the insertion portion of the instrument.

Take a proper remedial action by making a judgment from the specialist's view.

Cut the loop or the suture using a spare instrument.

Withdraw the instrument on which the loop or the suture is caught.

Manufactured by

OLYMPUS MEDICAL SYSTEMS CORP.

2951 Ishikawa-cho, Hachioji-shi, Tokyo 192-8507, Japan
Fax: (042)646-2429 Telephone: (042)642-2111

Distributed by

OLYMPUS AMERICA INC.

3500 Corporate Parkway, P.O. Box 610, Center Valley, PA 18034-0610, U.S.A.
Fax: (484)896-7128 Telephone: (484)896-5000

OLYMPUS LATIN AMERICA, INC.

5301 Blue Lagoon Drive, Suite 290 Miami, FL 33126-2097, U.S.A.
Fax: (305)261-4421 Telephone: (305)266-2332

KEYMED (MEDICAL & INDUSTRIAL EQUIPMENT) LTD.

KeyMed House, Stock Road, Southend-on-Sea, Essex SS2 5QH, United Kingdom
Fax: (01702)465677 Telephone: (01702)616333

OLYMPUS MOSCOW LIMITED LIABILITY COMPANY

Elektrozavodskaya str. 27 bld 8, 107023 Moscow, Russia
Fax: (7)495-926-7072 Telephone: (7)495-926-7077

OLYMPUS (BEIJING) SALES & SERVICE CO., LTD.

A8F, Ping An International Financial Center, No. 1-3, Xinyuan South Road,
Chaoyang District, Beijing, 100027 P.R.C.
Fax: (86)10-5976-1299 Telephone: (86)10-5819-9000

OLYMPUS KOREA CO., LTD.


(Seocho-dong, Majestar City Tower One), 2F & 3F, 12, Seocho-daero 38-gil,
Seocho-gu, Seoul, 06655, Republic of Korea
Fax: (02)6255-3210 Telephone: 1544-3200

OLYMPUS SINGAPORE PTE LTD

438B Alexandra Road #03-07/12, Alexandra Technopark Bld B, Singapore 119968
Fax: 6834-2438 Telephone: 6834-0010

OLYMPUS AUSTRALIA PTY LTD

3 Acacia Place, Notting Hill, VIC 3168, Australia
Fax: (03)9543-1350 Telephone: (03)9265-5400

OLYMPUS EUROPA SE & CO. KG

(Premises/Goods delivery) Wendenstrasse 20, 20097 Hamburg, Germany
(Letters) Postfach 10 49 08, 20034 Hamburg, Germany
Fax: (040)23773-4656 Telephone: (040)23773-0

©2018 OLYMPUS MEDICAL SYSTEMS CORP. All rights reserved.

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
Addendum to User manual for medical device: “Single Use Loop Cutter FS-410U”

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Working length, mm	2300 ± 10%
Roughness of cutters, Ra μm	Not more than 1,6
Hardness of blades, Rockwell scale, HRC	70 ± 5
Necessary force to pull of push the slider for opening and closing the scissors, N	Not more than 4,9
Device working parts connection strength, N	Not less than 19,6

NOTE: for all technical characteristics specified in the section 8.2. "Specifications" of INSTRUCTION MANUAL tolerances from declared values are +/- 5%, unless otherwise stated. The table above shows the characteristics and their values with tolerances not listed in the INSTRUCTION MANUAL

Key accessories, medical devices or non-medical devices, to be used in combination with medical device submitted for registration (annex)

Name of MD	Registration certificate
CF-LV1L	FSZ 2012/12490
CF-LV1I	FSZ 2012/12490
CF-H170L	RZN 2015/3422
CF-H170I	RZN 2015/3422
CF-H180AL	RZN 2015/3464
CF-H180AI	RZN 2015/3464
CF-H185L	RZN 2015/3464

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3. List of applied standards (Biocompatibility)

Identification	Designation
EN ISO 10993-1:2009	Biological evaluation of medical devices — Part 1: Evaluation and testing within a risk management process
EN ISO 10993-5:2009	Biological evaluation of medical devices — Part 5: Tests for in vitro cytotoxicity
EN ISO 10993-10:2010	Biological evaluation of medical devices — Part 10: Tests for irritation and skin sensitization

4. List of applied standards (Labeling)

Identification	Designation
EN ISO 15223-1:2016	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied — Part 1: General requirements
EN 1041:2008+A1	Information supplied by the manufacturer of medical devices

The warranty expiration date of medical device matches to expiration date.

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To be submitted to the Federal Service for Surveillance in Healthcare (Roszdravnadzor) in Russian Federation. We, OLYMPUS MEDICAL SYSTEMS CORP., located at 2951 Ishikawa-cho, Hachioji-shi, Tokyo 192-8507, Japan, as the legal manufacturer of the Medical device "Single Use Loop Cutter FS-410U", hereby state that attached documents is true:

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