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
М.П. / Stamp


BARDD | PERIPHERAL
VASCULAR

State of Arizona
County of Maricopa

This instrument was acknowledged before me this 9 February 2022, by Kathryn Swanstrom.
In witness whereof I herewith set my hand and official seal.

Signature of Notary Public

Notary Stamp

1. NAME OF THE MEDICAL DEVICE

Bard® Mission® disposable sterile core biopsy instruments of different types:

1. Bard® Mission® disposable sterile core biopsy instrument (14G x 6 cm);
2. Bard® Mission® disposable sterile core biopsy instrument (14G x 10 cm);
3. Bard® Mission® disposable sterile core biopsy instrument (14G x 16 cm);
4. Bard® Mission® disposable sterile core biopsy instrument (16G x 6 cm);
5. Bard® Mission® disposable sterile core biopsy instrument (16G x 10 cm);
6. Bard® Mission® disposable sterile core biopsy instrument (16G x 16 cm);
7. Bard® Mission® disposable sterile core biopsy instrument (18G x 6 cm);
8. Bard® Mission® disposable sterile core biopsy instrument (18G x 10 cm);
9. Bard® Mission® disposable sterile core biopsy instrument (18G x 16 cm);
10. Bard® Mission® disposable sterile core biopsy instrument (18G x 20 cm);
11. Bard® Mission® disposable sterile core biopsy instrument (18G x 25 cm);
12. Bard® Mission® disposable sterile core biopsy instrument (20G x 6 cm);
13. Bard® Mission® disposable sterile core biopsy instrument (20G x 10 cm);
14. Bard® Mission® disposable sterile core biopsy instrument (20G x 16 cm);
15. Bard® Mission® disposable sterile core biopsy instrument (20G x 20 cm);

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
Pic. 1 Bard® Mission® Disposable sterile Core Biopsy Instrument (Herein referred to as the "Instrument")

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Table 1. Size of the needle of the instrument Bard® Mission®.

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3.2 Warnings

1. Post-biopsy patient care may vary with the biopsy technique utilized and the individual patient's

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physiological condition. Observation of vital signs and other precautions should be taken to avoid and/or treat potential complications that may be associated with biopsy procedures.

2. When used for breast biopsy, the product is for diagnosis only.
3. The collection of multiple core biopsy samples may help to ensure the detection of any cancer tissue. A "negative" biopsy in the presence of suspicious radiographic findings does not preclude the presence of carcinoma.

4. The Instrument has been designed for single use only. Reusing this medical device bears the risk of cross-patient

contamination as medical devices – particularly those with long and small lumina, joints, and/or crevices between components – are difficult or impossible to clean once body fluids or tissues with potential pyrogenic or microbial contamination have had contact with the medical device for an indeterminable period of time. The residue of biological material can promote the contamination of the device with pyrogens or microorganisms which may lead to infectious complications.

5. Do not resterilize the Instrument. After resterilization, the sterility of the product is not guaranteed because of an indeterminable degree of potential pyrogenic or microbial contamination which may lead to infectious complications. Cleaning, reprocessing and/or resterilization of the present medical device increases the probability that the device will malfunction due to potential adverse effects on components that are influenced by thermal and/or mechanical changes.

Note: Inspect the Instrument needle components for damaged point, bent shaft or other imperfections prior to use and after each sample is collected. DO NOT USE the device if any imperfection is noted.

Note: After use, the Instrument may be a potential biohazard. Handle and dispose of in accordance with acceptable medical practice and applicable local, state, and federal laws and regulations.

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See Figure 3.

Biopsy procedure for the Bard® Mission® Disposable Core Biopsy Instrument

The biopsy procedure must be performed using appropriate aseptic techniques.

1. Prepare site as required. Adequate anesthesia should be administered prior to incision of the skin.
2. For ease of insertion, puncture the skin with a scalpel at the entry site.
3. Verify the Instrument is energized (cocked) by looking at the penetration depth indicator and confirming either a "10" or "20" is visible. **See Figure 3.**
4. With the inner stylet fully retracted, so that the sample notch is covered by the cannula, insert tip of needle to the point to be biopsied. **DO NOT** advance the inner stylet by pressing on the plunger until the Instrument is in position.
5. Press the plunger to advance the inner stylet (i.e. sample notch) into tissue. At the end of the advancement, the plunger will encounter a stop and the fire ready indicator will meet the device housing. **See Figure 4.** With imaging, verify the sample notch is in the target area to be biopsied.
6. Fire the cutting cannula by depressing the plunger past the fire ready indicator to capture the biopsy specimen in the sample notch.
7. Remove the needle from the patient and pull the plunger back to acquire the biopsy specimen. Gently advance the plunger forward to expose the biopsy specimen. **DO NOT** press the plunger past the fire ready indicator during this step.
8. If additional biopsies are required, pull back on plunger to withdraw the inner stylet and repeat steps 3 through 7.

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Figure 3: Plunger retracted to energize (cock) Instrument



Figure 4: Front of Instrument

4. TECHNICAL CHARACTERISTICS OF THE DEVICE

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4.2. Materials

Part	Material	Supplier
Bard® Mission® instrument		
Stylet	Stylet: stainless steel AISI 304 Hub: Cyclopedia C6600 PC/ABS, Colorant 8T9D071(White)	Spectra, USA
Mission Core needle cannula	Stylet: stainless steel AISI 304 Hub: Cyclopedia C6600 PC/ABS, Colorant 8T9D071(White)	Spectra, USA
Lower and Upper Case	Cyclopedia C6600 PC/ABS, Colorant 8T9D071(White) Ink: PG-NT, PMS 349C (green)	Nipro, USA
Slide Carriage	Hostaform S9363	Nipro, USA
Plunger	Acetal Celcon M-90, colorant CPOM298 Green CE5430K20 Purple CR5650K20 Pink CL 3073K20 yellow	Nipro, USA

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4.3. Packaging materials

Component	Material	Supplier
Pouch	Uncoated Tyvek: Tyvek 1073B Transparent Film: 35721K Nylon/adhesive 2635B	Bemis, USA
Needle Protection Tube	Material: Low Density Polyethylene (L.D.P.E) #640I Color: Natural	Kelpac Medical, USA
Carton	200 Lb. (or equivalent) Test EFlute Singlewall Corrugated Fiberboard, Bleached White One Side Plain Kraft (Printed) Ink: Pantone 349 Green	Flexible Packaging, USA
Shipping carton	200lb C-Flute single walled corrugated fiberboard RSC	Flexible Packaging, USA

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4.4. Requirements to the primary package materials:

- Materials shall contain acceptable level of cleanness, hard particulates and fluff;
- Materials shall not have holes, bleakness, gaps, folds or local thickening/thinning, causing the functional defects

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5.MEDICAL DEVICE LABELING AND PACKAGING

5.1 Packaging of the device. Transportation packaging.

Bard® Mission® Instruments and Kits which are packaged individually in a Tyvek pouch which is then packaged as 5 units per carton. The products are sterilized and distributed in the cartons.
Size of shipping carton depends on order.

Pouch dimensions	Long pouch (for the 20 cm and 25 cm length needles): 20.25" +/- 0.22" (L) x 4.50" +/- 0.13" (W) Short Pouch: 16.25" +/- 0.125" (L) x 4.5" +/- 0.125 (W)
Carton dimensions	Long Carton (for the 20 cm and 25 cm length needles: 4.5" +/- 0.125" (L) x 2.5" +/- 0.125" (W) x 20.25" +/- 0.125" (D) Short Carton: 4.5"/+/- 0.125" (L) x 2.5"/+/- 0.125" (W) x 16.25" +/- 0.125" (D)

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
2. Pouch and carton label.

Note: example of the label for the instrument produced in Dominican Republic

BARD® MISSION® Disposable Core Biopsy Instrument

Bard Peripheral Vascular, Inc.
C. R. Bard Inc.
1625 West 3rd Street
Tempe, AZ 85281
USA
Tel: (800) 321-4254
(480) 894-9515
Fax: (480) 966-7062

REF	Catalogue number	1406MS	
LOT	Lot number	1234567890	
	Use-by date	2022-05-28	
PD	Penetration depth	10mm and 20mm	

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	Consult instructions for use			Not made with natural rubber latex		
	Single use		Do not use if the product sterile barrier system or its packaging is compromised			
	Do not resterilize		Sterilized using Ethylene Oxide		Non-pyrogenic	
		BD, the BD Logo, Bard, and Mission are trademarks of Becton, Dickinson and Company or its affiliates.				© 2019 BD. All rights reserved.

Instrument Assembled in Dominican Republic

BARD Biopsy

Manufacturer:
Bard Peripheral Vascular, Inc.
1625 West 3rd Street
Tempe, AZ 85281 USA
Tel: 1-800-894-9515
1-800-321-4254
Fax: 1-800-966-7062
1-800-440-5376
www.crbard.com/biopsy

CE **EC REP**
2797
Authorized Representative
in the European Community
BD Switzerland Sàrl
Terre Bonne Park - A4
Route de Crassier 17
1252 Eysins - Switzerland

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
5.2. Pouch and carton label.

*Note: example of the label for the instrument
produced in Puerto Rico.*

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BARD

REF Catalogue Number **1406MS**
LOT Lot Number **ASDF1234**
Use By **2020-10-28**
PD Penetration Depth **10mm and 20mm**

(01)00801741142505

(17)201028(10)ASDF1234

Disposable Core Biopsy Instrument

Instrument pour biopsie au trocart jetable
Einmal-Stanzbiopsieinstrument
Strumento per agobiopsia monouso
Dispositivo para biopsia core desechable
Wegwerfbaar holle naaldbiopsie-instrument
Instrumento de Biopsia por Aguiha Grossa Descartável
Αναλώσιμο Εργαλείο Βιοψίας με κόπτοισα βελόνα
Instrument til grovnålsbiopsi til engangsbrug
Engångsinstrument för grovnålsbiopsi

Keräkäyttöinen paksuneulabiopsialaite
Kjernerbiopsieinstrument til engangsbrug
Jednorazowy przyrząd do biopsji gruboigłowej
Eldobható vastagű-biopsziás berendezés
Nástroj pro core biopsii k jednorázovému použití
Altlabılır Kor Biyopsi Aleti
Одноразовый инструмент для толстоигльной биопсии
拋棄式空芯切片器械
일회용 총생검 기구
Jednorazový nástroj na jadrovú biopsiu

Consult Instructions for Use

Rx Only

Single Use

Not Made With Natural Rubber Latex

Do Not Use if Package is Damaged

Do Not Resterilize

STERILE EO Sterilized Using Ethylene Oxide

Non-pyrogenic

Bard and Mission are trademarks and/or registered trademarks of C. R. Bard, Inc. or an affiliate.

Copyright ©2018 C. R. Bard, Inc. All rights reserved.

BARD Biopsy

Manufacturer:
Bard Peripheral Vascular, Inc.
1625 West 3rd Street
Tempe, AZ 85281
USA
TEL: 1-800-894-9515
1-800-321-4254
FAX: 1-800-966-7062
1-800-440-5376
www.bardbiopsy.com
0001

CE 2797
Authorized Representative
in the European Community
BD Switzerland Sàrl
Terre Bonne Park - A4
Route de Crassier 17
1262 Eysins - Switzerland
Email: EC-Rep@BD.com

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	Bard® Mission® Disposable Core Biopsy Instrument
	Contents
	Catalog number
	Use by date
	Peel to open
	Lot number
	Consult instruction for use
	Authorized Representative in the European Community
	Single use
	Do not resterilize
	Penetration depth
	Sterilized using Ethylene Oxide
	Not made with natural rubber latex
	For healthcare professionals
	Do not use if package is damaged
	Manufacturer

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During import to Russian Federation the medical device transport packages would be provided with the local Russian stickers that are prevailing when supplied to Russian Federation. In case of discrepancy of information on the sticker and on the package, it is correct to consider the information on the sticker. Local Russian sticker does not cover the original data of the manufacturer and is of the following content:

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English translation of the Russian sticker (for reference only)

Bard® Mission® disposable sterile core biopsy instruments of different types:
 Bard® Mission® disposable sterile core biopsy instrument
Manufacturer: Bard Peripheral Vascular, Inc., 1625 West 3rd Street Tempe, Arizona 85281, USA
Catalog number : XXXX
Storage conditions: temperature от 0°C до +30°C, humidity up to 75% no condensate
Use by: see on the package
Authorized representative of the manufacturer:
 Bard Rus Llc, Russian Federation, 121596, Moscow, Gorbunova street, 2, bld. 3, floor 5, premises II, room 1, office A-500.1 cabinet 5.
 tel +7499-372-50-02

RC #XXXX

Information on Russian sticker, please consider the priority for supply in the Russian Federation

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ability under this limited product warranty will be limited to repair or replacement of the defective product, in Bard Peripheral Vascular's sole discretion or refunding your net price paid. Wear and tear from normal use or defects resulting from misuse of this product are not covered by this limited warranty.

TO THE EXTENT ALLOWABLE BY APPLICABLE LAW, THIS LIMITED PRODUCT WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, WHETHER EXPRESS OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. IN NO EVENT WILL BARD PERIPHERAL VASCULAR BE LIABLE TO YOU FOR ANY INDIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES RESULTING FROM YOUR HANDLING OR USE OF THIS PRODUCT.

Some states/countries do not allow an exclusion of implied warranties, incidental or consequential damages. You may be entitled to additional remedies under the laws of your state/country.

An issue or revision date and revision number for these instructions are included for the user's information on the last page of this booklet. In the event 36 months have elapsed between this date and product use, the user should contact Bard Peripheral Vascular to see if additional product information is available.

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10.CLAIM

For questions about the circulation of the medical device in the Russian Federation, please contact Authorized Representative of the Manufacturer in the Russian Federation:
«Bard Rus Llc» Gorbunova St., 2, bld. 3, floor 5, premises II, room 1, office A-500.1 cabinet 5
Moscow, 121596, Russian Federation
Tel./Fax: +7 499 372 50 02

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11. LIST OF STANDARDS

Standard Number	Title
ASTM F1980-16	Standard Guide for Accelerated Aging of Sterile Barrier Systems for Medical Devices (Reapproved for 2011)
ASTM D4332-14	Standard Practice for Conditioning Containers, Packages, or Packaging Components for Testing
ASTM D999-08(2015)	Standard Test Methods for Vibration Testing of Shipping Containers
EN 1041:2008+A1	Information supplied by the manufacturer of medical devices
EN ISO 15223-1	Symbols for use in the labeling of medical devices
EN 556-1	Sterilization of medical devices - Requirements for medical devices to be designated "STERILE" - Part 1: Requirements for terminally sterilized medical devices
EN ISO 11737-1	Sterilization of medical devices – Microbiological methods – Part 1: Determination of a population of microorganisms on products
EN ISO 15223-1	Medical devices - Symbols to be used with medical device labels, labeling and information to be supplied - Part 1: General Requirements
EN ISO 14971	Medical devices - Application of risk management to medical devices
EN ISO 13485	Medical devices - Quality management systems - Requirements for regulatory purposes
EN ISO 11607-1	Packaging for terminally sterilized medical devices – Part 1: Requirements for materials, sterile barrier systems and packaging systems
EN ISO 11607-2	Packaging for terminally sterilized medical devices – Part 2: Validation requirements for forming, sealing and assembly processes
ISO 11135	Sterilization of health care products - Ethylene oxide - Part 1: Requirements for development, validation and routine control of a sterilization process for medical devices
EN ISO 10993-1	Biological evaluation of medical devices - Part 1: Evaluation and testing within a risk management process
ISO 14698	Cleanrooms and associated controlled environments - Biocontamination control
ISO 14644	Clean rooms and associated controlled environments.
ISTA 1G	Packaged Products 150 lb (68 kg) or Less (Random Vibration)

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