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Startup

1. Place the power switch on the back of the analyzer in the [I] position to start the analyzer.
2. Enter the username and password on the screen.
3. The startup initialization and cleaning procedures begin. After the procedures, you can perform operations on the user interface.

Daily Quality Control

Before running samples every day, run the controls in the QC interface. For details, see the section about quality control in *Operator's Manual*.

Running the Whole Blood Samples

1. Click Mode & ID in the Sample Analysis screen.
2. Select Venous Whole Blood (VWB) or Capillary Whole Blood (CWB) as required.
3. Input the sample ID for the sample to be run, then click OK.
4. Set the well-mixed sample (uncapped) below the sample probe, press the aspirate key to start the analysis. Remove the tube when you hear a beep.

NOTE: The tube should be placed vertically upward, not tilted or upside down.

Running the Predilute Samples

1. Click Mode & ID in the Sample Analysis screen.
2. Select Predilute, input the sample ID for the sample to be run, then click OK.

3. Click Add Diluent on the menu screen, present a clean centrifugal tube (uncapped) under the sample probe.

4. Press the aspirate key. When you hear a beep, remove the tube. If more portions of diluent are needed, take another clean centrifugal tube and repeat this step.
5. Click Cancel to exit the diluent dispensing.
6. Add 20µL of capillary blood to the diluent, close the tube cap and shake the tube to mix the sample.

NOTE: The centrifugal tube should be placed vertically upward, not tilted or upside down.

7. Uncap the centrifugal tube and present it to the sample probe so that the probe tip is vertically in contact with the bottom of the tube, then press the aspirate key to start the analysis. Remove the tube when you hear a beep.

Shutdown

NOTE

Shut down strictly as instructed on screen, and on the shutdown maintenance with the cleanser.

1. Click  on the menu screen, perform the shutdown procedure by following the screen instructions.
2. Turn off the power switch of the analyzer.

3. Empty the waste container and dispose of the waste properly.

Precautions

Troubleshooting

If an error occurs, the error message displays at the top of the screen. To remove the error, click the error message area, and click Remove Error. For details, see the Troubleshooting chapter in *Operator's Manual*.

When the analyzer prompts that the diluent is running out or waste is full (actually the diluent is enough and the waste is not full), please take out the corresponding sensor from the container to shake 2~3 times, then put it back. Click Remove Error to remove the error. If the error still exists, contact our customer service department.

Replacing Reagents

Once the new reagent is connected, see the section about reagent management in *Operator's Manual* for reagents replacement.


Routine Maintenance

Maintenance Scenarios	Operations
Shutdown or daily tips	Perform cleanser maintenance as instructed.
Before a short distance transport (transporting time ≤ 2h)	Perform Empty Fluidics.
Not to be used for over one week or a long distance transport (transporting time > 2h)	Perform Prepare to Ship.

Sleep

The analyzer enters sleep mode when no operation is done for a specified interval. You can touch the screen to wake up the analyzer.

For detailed instructions and precautions of the analyzer, see the *Operator's Manual*.

Startup

1. Place the power switch on the back and the left side of the analyzer in the (I) position to start the analyzer.
2. Enter the username and password on the interface.
3. The startup initialization and cleaning procedures begin. After the procedures, you can perform operations on the user interface.

Daily Quality Control

Before running samples every day, run the controls in the QC interface. For details, see the section about quality control in *Operator's Manual*.

Running the Venous Whole Blood and Capillary Whole Blood Samples

1. Click **Mode & ID** in the Sample Analysis screen.
2. Select **Venous Whole Blood (VWB)** or **Capillary Whole Blood (CWB)** as required.
3. Select the measurement mode as required and input the sample ID for the sample to be run, then click **OK**.
4. Set the well-mixed sample (uncapped) below the sample probe, press the aspirate key to run the sample. Remove the tube when you hear a beep.

NOTE: The tube shall be placed vertically upward, not tilted or upside down.

Running the Predilute Samples

1. Click **Mode & ID** in the Sample Analysis interface to set the analysis mode and select **Predilute**.
2. Select the measurement mode as required and input the sample ID for the sample to be run, then click **OK**.

3. Click **Add Diluent** on the menu interface, then present a clean centrifugal tube under the sample probe.

4. Press the aspirate key. When you hear a beep, remove the tube. If more portions of diluent are needed, take another clean centrifugal tube and repeat this step.
5. Click **Cancel** to exit the diluent dispensing.
6. Add 20µL of blood to the diluent, close the tube cap and shake the tube to mix the sample.

NOTE: The centrifugal tube shall be placed vertically upward, not tilted or upside down.

7. Uncap the centrifugal tube and present it to the sample probe so that the probe tip is vertically in contact with the bottom of the tube, then press the aspirate key to run the sample. Remove the tube when you hear a beep.

Shutdown

NOTE

To keep special protein reagent in cold temperature, do not turn off the power switch on the back of the main unit (except special occasions like transportation or long-term non-use of CRP measurement).

Shut down strictly as instructed on-screen, and do the shutdown maintenance with the cleanser.

1. Click  on the menu interface, perform the shutdown procedure by following the screen instructions.

2. Turn off the power switch at the left side of the analyzer.
3. Empty the waste container and dispose of the waste properly.

Precautions

Troubleshooting

If an error occurs, the error message displays on the top of the interface. To remove the error, double click the message area, and click Remove Error.

When the analyzer prompts that the diluent is running out or waste is full error (actually the diluent is enough and the waste is not full), please take out the corresponding sensor from the container to shake 2-3 times, then put it back. Click Remove Error to remove the error.

Replacing Reagents

Once the new reagent is connected, see the section about reagent management in *Operator's Manual* for reagents replacement.

Routine Maintenance

Maintenance Scenarios	Operations
Shutdown or daily tips	Perform cleanser maintenance procedure as instructed.
Not to be used for 1-2 weeks or before a short distance transport (transporting time ≤ 2h)	Perform Empty Fluidics.
Not to be used for more than 2 weeks or a long distance transport (transporting time > 2h)	Perform Prepare to Ship.


NOTE

For dry areas, the operation of **Empty Fluidics** should be carried out when the analyzer has not been used for two to three days, and the operation of **Prepare to Ship** should be carried out when the analyzer has not been used for more than three days.

Sleep

The analyzer enters sleep mode when no operation is done for a specified interval. You can touch the screen to wake up the analyzer.

For detailed instructions and precautions of the analyzer, see the *Operator's Manual*.

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