**Authentication of Signature**

**On Behalf of a Body Corporate, Where the  
Language of the Document is not Known to the  
Person Signing**

I, the undersigned, **Itai Izhaky, Notary** holding  
license no. 2084632, hereby certify that:

On 07/07/2020 there appeared before me at 3 Ariel  
Sharon St., Or-Yehuda, **Mrs. Rachel Dvir**, whose  
identity has been proven to me by an Israeli  
biometric ID Card No. **011035730**, issued on  
04/04/2019.

And I was convinced that the person standing before  
me fully understood the significance of the action  
and voluntarily signed the attached document,  
marked "A", on behalf of "ALLIGN  
TECHNOLOGY LTD", Private Company №  
512020934, at the address 3 Ariel Sharon St., Or-  
Yehuda; The document is drawn up in the Russian  
language, which is not known to her and is also not  
known to me.

I confirm that for the purpose of proving the right of  
the above to sign on behalf of "ALLIGN  
TECHNOLOGY LTD" I was provided with  
a Protocol of the corporation's board of directors  
meeting approved by Adv. Sarah Azulay, issued on  
06/07/2020.

And I also certify that the translator Mr. Vladimir  
Urbach (I.D. No. 307264036), from the address 5  
Einstein St. Revovot, whom I know to be fluent in  
the Russian and English languages, translated the  
said document in my presence for Mrs. Rachel Dvir  
into the English language, which is known to her and  
me, before she signed it.

The Translator Affidavit is attached and marked "1".

In witness whereof I hereby authenticate the  
signature of Mrs. Rachel Dvir by my own signature  
on this day, 07/07/2020.

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
1  
**Declaration**

I the undersigned, **Urbach Vladimir I.D No. 307264036**

After I have been warned to tell the truth and that I will be destined to be punished by Law if I will not, hereby declare in writing as follows:

1. I hereby declare that I am well acquainted with the Russian, And the English languages and that the attached document marked "B" is a correct English translation from the Russian language original document which its copy is attached hereto marked "A".

2. This is my name and this is my signature and all that is written above is true and correct.

  
Vladimir Urbach

A

## Operation manual

Intraoral Scanner iTero® Element 5D, variants:  
iTero® Element 5D, iTero® Element laptop  
configuration with accessories

### Legal Manufacturer

Align Technology, Ltd., 3 Ariel Sharon  
Boulevard, Or-Yehuda, Israel

**Rachel Ovir**

Director, Quality Assurance & Regulatory Affairs  
Align Technology Ltd.

Approve

Signature, stamp

2020

iTero element 5D

# Wheel stand configuration **operation** **manual**

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This document describes the iTero Element 5D Optical Impression Device, English language version, Updated October 2019.

## Trademarks

The following are trademarks and/or service marks of Align Technology, Inc. or one of its subsidiaries or affiliated companies and may be registered in the U.S. and/or other countries: Align Technology Inc., iTero

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
IEC 60601-1-2 Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral standard: Electromagnetic phenomena - Requirements and tests".

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IEC 60601-1 Medical electrical equipment –Part 1: General requirements for basic safety and essential performance.

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### Power supply

- Power is supplied to the system via an internal medical grade power supply.

### Battery power

- Charging – Battery will be fully charged after being plugged into a power source for 2 hours.
- With a fully charged battery, the user can scan up to 30 minutes with the iTero Element 5D scanner, without having to plug in for power.

### Electric warning

- Electric shock hazard!** Only authorized Align Technology technicians can remove external panels and covers. There are no user-serviceable parts inside.
- To avoid risk of electric shock, iTero Element 5D must only be connected to a supply mains with protective grounding.
- Only Align Technology approved Web Camera or DOK should be connected to the USB socket on the back side of the System.

### Wireless LAN

- The system comes equipped with a Wireless LAN unit.

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### Prescription health device

- The system serves as a prescription medical device and should be operated by qualified health-care providers only.

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## Cleaning & disinfection

- To avoid cross contamination, it is mandatory that after each patient session the disposable plastic sleeve be replaced and the scanning unit be disinfected.
- To avoid cross contamination, it is mandatory that after each patient session the gloves shall be removed and replaced.
- To avoid cross contamination, discard gloves when torn, contaminated, or removed for any reason.
- Dispose of scanner sleeves according to standard operating procedures or local regulations for the disposal of contaminated medical waste.

## Unpacking & installing

- The system should be unpacked and installed following Align Technology's instructions.

## Work environment

- The system should be moved between rooms with utmost care to avoid damage.
- WARNING:** Do not block the air vents on the Scanning Unit and Base Unit.
- WARNING:** System is intended for indoor use only. It should not be exposed to direct sunlight, excessive heat or humidity.

## Electro magnetic interference

- WARNING:** This device has been tested and found to comply with the requirements for medical devices according to standard IEC60601-1-2. This standard is designed to provide reasonable protection against harmful interference in a typical medical installation. However, because of the proliferation of radio-frequency transmitting equipment and other sources of electrical noise in the healthcare environments (e.g., cellular phones, mobile two-way radios, electrical appliances), it is possible that high levels of such interference due to close proximity or strength of source, may result in disruption of performance of this device. ESD may cause an interference during scanning. In case of ESD interference the system will stop scanning / displaying proper viewfinder image. If the failure occurs within scan, the user may be required to press the wand scan button to restart the scan. The viewfinder will always recover automatically.

## General

- WARNING:** No modification of this equipment is allowed.
- WARNING:** The touchscreen always needs to be in a stand while in operation!

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
**Table of contents**

|                                                                             |           |
|-----------------------------------------------------------------------------|-----------|
| <b>Chapter 1: Introduction</b>                                              | <b>2</b>  |
| About this document                                                         | 2         |
| Intended use / indication for use                                           | 2         |
| Benefits of the iTero Element 5D system                                     | 2         |
| The iTero Element 5D user interface                                         | 3         |
| <b>Chapter 2: Basic hardware features</b>                                   | <b>4</b>  |
| Front view of system                                                        | 4         |
| Back view of system                                                         | 5         |
| Scanning unit (wand)                                                        | 6         |
| <b>Chapter 3: Assembly instructions</b>                                     | <b>7</b>  |
| Step 1: Assembly                                                            | 7         |
| Assembly of iTero Element 5D counter stand                                  | 10        |
| Step 2: Make it mine process                                                | 12        |
| <b>Chapter 4: Operating instructions</b>                                    | <b>13</b> |
| End-of-day shut down                                                        | 13        |
| Moving system within the office                                             | 13        |
| <b>Chapter 5: Scanner handling, cleaning, and disinfection instructions</b> | <b>14</b> |
| Handling of the scanning unit (wand)                                        | 14        |
| Handling of the scanning unit cable                                         | 14        |
| Cleaning the scanning unit                                                  | 14        |
| Disinfecting the scanning unit                                              | 14        |
| <b>Chapter 6: Changing sleeves between patients</b>                         | <b>15</b> |
| Replacing disposable sleeves                                                | 15        |
| Scanner sleeves                                                             | 16        |
| <b>Chapter 7: Accessories (options)</b>                                     | <b>17</b> |
| Counter stand and cradle                                                    | 17        |
| Articulator                                                                 | 17        |
| <b>Chapter 8: Clinic LAN network guidelines</b>                             | <b>18</b> |
| Introduction                                                                | 18        |
| Preparations                                                                | 19        |
| Router guidelines                                                           | 19        |
| Internet connection guidelines                                              | 19        |
| Firewall                                                                    | 19        |
| Wi-Fi tips                                                                  | 20        |
| Align hostname recommendations                                              | 21        |
| <b>Appendix A: EMC declaration</b>                                          | <b>22</b> |
| <b>Appendix B: System specifications</b>                                    | <b>23</b> |

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1. The optical impression (CAD/CAM) feature of iTero Element 5D is intended/indicated for use to record the topographical images of teeth and oral tissue. Data generated from iTero may be used in conjunction with the production of dental devices (e.g., aligners, braces, appliances, etc.) and accessories.
2. iTero Element 5D software is used with the iTero scanner in capturing 3D digital impressions of teeth, oral soft tissue and structures, and bite relationship. iTero software controls the processing of the data, facilitating the integration of data, and exporting of the data for CAD/CAM fabrication of dental restorations, orthodontic devices, abutments, and accessories. In addition to scan data, various patient and case information can be imported/exported or used for simulation purposes. Other functions are available for verification and service of the system and to serve as an order management tool.
3. The Element 5D NIRI functionality, is a diagnostic aid for the detection of interproximal caries lesions above the gingiva and for monitoring the progress of such lesions.

### Benefits of the iTero Element 5D system

The iTero Element 5D system provides important advantages over existing crown-production methods, including powder-free scanning, greater crown-production accuracy, and immediate feedback during the scanning process.


Refer to our website <http://www.itero.com> to learn how the iTero Service can enhance your business by increasing patient satisfaction, improving clinical outcomes, and enhancing office efficiency.

The iTero Element 5D user interface

The iTero Element 5D system provides an intuitive user interface for performing digital scans for Restorative or Orthodontic use. The doctor is guided through the scanning sequence by means of visual and text assistance. The touchscreen and wand buttons are used to respond to screen instructions during the scanning process.

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
Chapter 2:

## Basic hardware features


Custom wheel stand hardware features: Front view of the system

Custom wheel stand hardware features: Back view of the system

## Scanning unit (wand)

## Chapter 3:

# Assembly instructions

### Step 1: Assembly

Please follow the instructions below to assemble your iTero Element 5D scanner.

- A** Wheel stand
- B** Wand with cable
- C** Wand cradle
- D** HD touch screen
- E** External battery
- F** Wand cable

1. Check the content of the box

2. Connect post to the wheel base

3. Tighten the two allen screws using the larger wrench

4. Remove cover from the back of the handle

5. Attach the wand cradle to the front of the wheel stand

6. Hold the cradle

7. Tighten the wand cradle allen screw on the back using the smaller wrench

8. Reattach the cover behind the handle


9. Remove the magnetic cover from the back of the wheel stand frame

10. Remove the battery cover

11. Slide the battery into the battery slot and tighten the thumb screws

12. Lift the HD screen to mount it

13. Turn the scanner around and tighten the thumb screw to secure the HD screen

14. Attach the power cable to the port labeled DC

15. Power cable inserted

16. Attach the magnetic back cover

17. Place the wand into the cradle


18. Attach the wand cable on the back of the HD screen

19. Attach the power cable on the bottom of the wheel stand

20. On the bottom of the wheel stand, post and secure the cable with the clip

21. Plug in the power cable

22. Position the web camera on the HD screen for remote training or support sessions

23. Plug in the webcam to the USB port at the bottom of the HD screen


24. Press button to switch on the scanner

## Assembly of iTero Element 5D counter stand


Please follow the instructions below to assemble your iTero Element 5D counter stand:


1. Insert post into frame.

2. Tighten the post using the allen wrench

3. Remove cover

4. Install new cover


5. Lift the HD screen to mount it.


6. Turn the scanner around and tighten thumb screw to secure HD screen

7. Place the cover

8. Attach the power cable to the port labelled DC

9. Attach the wand cable on the back of the HD screen

10.  
Plug in the power cable


11.  
Position the web camera on the HD screen  
for remote training or support sessions


12.  
Plug in the webcam to the USB port at the  
bottom of the HD screen

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**End-of-day shut down**

1. Close all files and applications.
2. Press and release the Power Switch on the bottom of the screen to shut down the system.

**Moving system within the office**

To ensure maximum system protection, it is recommended to have two people move the system. Follow these instructions for relocating the system:

1. Verify the scanning unit (wand) sits well inside the scanning unit cradle (holder).
2. Unplug system from the wall outlet.
3. Move the system carefully using two people.
4. Place the system at its new location and it plug into a wall outlet.

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- Soak a lint free cloth in ready to use CaviCide 1, and squeeze until the cloth is moist.
- Wipe the wand thoroughly to remove gross debris.
- The entire device needs to be visually inspected to ensure that no residual or debris remains prior to continuing to the next step.
- Soak two additional lint free cloths in CaviCide 1, and squeeze until the cloth is moist. Then wipe the device thoroughly, ensuring with special care that all the surfaces, to all edges and slots, are covered. Contact time will be at least 2 minutes.
- Then wet lint free cloths with distilled water and wipe all surfaces for at least 15-30 seconds.
- Then use dry lint free cloths to dry the surfaces.

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iTero Element 5D sleeve is intended for single patient use and shall be replaced after each patient with the purpose to avoid cross contamination.

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
Front view

## Articulator

The articulator is intended to use as a mechanical device that simulates movements of a patient's upper and lower jaws. Plaster casts of the patient's teeth and gums are placed in the device to reproduce the occlusion (bite) and articulation of the patient's jaws. The articulator is intended to fit dentures or provide orthodontic treatment.

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### Levels of Wi-Fi Internet Connectivity

Excellent  
>-50 dBm

Good  
-50 to -60 dBm

Fair  
-60 to -70 dBm

Weak  
<-70 dBm

- **IMPORTANT:** In order to achieve the best performance of your iTero Element 5D scanner, ensure that the Wi-Fi signal strength is "Excellent" or at least "Good".
- **CAUTION:** While scanning a patient, do not connect a LAN cable to iTero Element 5D – it is forbidden due to safety hazard reasons.

## Preparations

- The required Modem/Router should be configured with WPA2 Security standard, including a password.
- Ensure that your IT professional staff would be available when the scanner installation is planned to take place.
- Make sure that your Wi-Fi SSID credentials are available: Login & password.
- The minimal Wi-Fi strength signal for the system should display at least two "stripes", as shown in chapter 2, above).
- Diagnostic tool under "Settings", or the Connectivity Tool are suggested below.
- Test the local Wi-Fi connection with any PC computer, using the iTero connectivity tool (run the test as near as possible to the scanner location).
- Connectivity Tool Download at [fc1.orthocad.com/download/AlignSupport/ConnectivityCheck2.0.exe](http://fc1.orthocad.com/download/AlignSupport/ConnectivityCheck2.0.exe)
- Following are some suggestions for the office IT person, regarding what should be considered in order to prevent issues such as access or connectivity to/with the iTero scanner:
  1. Hostname recommendations related to Align services listening to ports 80 and 443.
  2. Do not prevent FTP communication since the scanner sends specific file types (.3ds and .3dc/.3dm).
  3. Disable any proxy client for data communication through TCP/IP.
  4. Do not add the scanner to any domain group.
  5. Do not run any group policy on the scanner as it may disrupt its proper functioning.

## Router guidelines

Minimum Standards: 802.11N / 802.11AC

## Internet connection guidelines

In order to achieve the best performance of your iTero Element 5D scanner, ensure that your internet connection upload speed is at least 1Mbps per scanner. Also, note that any additional devices connected to the internet in parallel to the scanner may affect the scanner's performance.

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1. Be sure to have a minimal number of walls and ceilings between the router and other network devices. Each barrier can reduce your adapter's range by 1-3 meters (3-9 feet).
2. Be sure to have a straight line, free of any partition, between network devices. Even a wall that seems rather thin can block a signal of 1 meter (3 feet) if the wall angle is shifted by only 2 degrees. To achieve the best reception, place all the devices so that the Wi-Fi signal travels straight through (90°) a wall or partition (instead of at an angle).
3. Construction materials make a difference. A solid metal door, or aluminum nails, can be very dense and may have an adverse effect on a Wi-Fi signal. Try to position access points, wireless routers, and computers so that the signal travels through drywalls or open doorways. Materials and objects such as glass, steel, metal, walls with insulation, water tanks (aquariums), mirrors, file cabinets, brick, and concrete may reduce your wireless signal.
4. Keep your iTero product away (at least 3-6 feet or 1-2 meters) from electrical devices or appliances that generate RF noise.
5. If you are using 2.4GHz cordless phones or X-10 (wireless products such as ceiling fans, remote lights, and home security systems), your wireless connection may severely degrade or entirely drop. The base of many wireless devices transmits an RF signal, even if the device is not in use. Position your other wireless devices as far as possible from your scanner and router.
6. In your area, there may be more than one active wireless network. Each network uses one or more channels. If the channel is near your system channels, the communication may gradually decline. Ask your IT department to check this, and if required, change the channel numbers used by your network.

### Align hostname recommendations

Align constantly improves its products and services. Hence, can rather commit to a Hostname, though not to a certain IP.

The following hostnames list was created to provide Align's scanners the proper operation functions, in order to be able to utilize all the advanced capabilities of the scanner performance.

Align hostnames recommendation:

| Host name                                                                                                    | Port    |
|--------------------------------------------------------------------------------------------------------------|---------|
| Mycadent.com                                                                                                 | 80, 443 |
| Myaligntech.com                                                                                              | 80, 443 |
| Export.mycadent.com                                                                                          | 80, 443 |
| Cbserver.mycadent.com                                                                                        | 80, 443 |
| Matstore.invisalign.com                                                                                      | 80, 443 |
| Matstore2.invisalign.com                                                                                     | 80, 443 |
| Matstore3.invisalign.com                                                                                     | 80, 443 |
| Matstore4.invisalign.com                                                                                     | 80, 443 |
| Matstoresg.invisalign.com                                                                                    | 80, 443 |
| Matstorechn.invisalign.com.cn                                                                                | 80, 443 |
| AWS IP range - Amazon global CDN service - IP address range varies depending on the location of the scanner. | 80, 443 |
| cloud.myitero.com                                                                                            | 443     |
| speedtest.tele2.net                                                                                          | 80      |
| alignapi.aligntech.com                                                                                       | 80, 443 |
| http://www.google.com                                                                                        | 80, 443 |
| http://www.microsoft.com                                                                                     | 80, 443 |
| http://www.yahoo.com                                                                                         | 80, 443 |
| iterosec.aligntech.com                                                                                       | 80, 443 |
| storage.cloud.aligntech.com                                                                                  | 443     |

## Appendix A:

# EMC declaration

## Summary of EMC test results for iTerо Element 5D

| Test                                                                   | Standard       | Class / Severity level                                                          | Test result |
|------------------------------------------------------------------------|----------------|---------------------------------------------------------------------------------|-------------|
| <b>Documentation (IEC 60601-1-2 sections 4 and 5)</b>                  |                |                                                                                 |             |
| General requirements for EMC                                           | section 4.1    | --                                                                              | Complies    |
| External labels                                                        | section 5.1    | --                                                                              | Complies    |
| Conformity of users' manual                                            | section 5.2.1  | --                                                                              | Complies    |
| Accuracy of technical description                                      | section 5.2.2  | --                                                                              | Complies    |
| <b>Emission (IEC 60601-1-2 / EN 60601-1-2 section 7)</b>               |                |                                                                                 |             |
| Conducted emission<br>Freq. range: 150 kHz - 30 MHz                    | CISPR 11       | Group 1 Class B 230, 120 & 100 VAC mains (50 Hz);<br>220 VAC mains (50 & 60 Hz) | Complies    |
| Radiated emission<br>Freq. range: 30 - 1000 MHz                        | CISPR 11       | Group 1 Class B                                                                 | Complies    |
| Harmonic current emission test                                         | IEC 61000-3-2  | 230 VAC mains (50 Hz & 220 V (60 Hz)                                            | Complies    |
| Voltage changes, voltage fluctuations and flicker test                 | IEC 61000-3-3  | 230 VAC mains & 220 VAC mains                                                   | Complies    |
| <b>Immunity (IEC 60601-1-2 section 8)</b>                              |                |                                                                                 |             |
| Immunity from electrostatic discharge (ESD)                            | IEC 61000-4-2  | 8 kV contact discharges & 15 kV air discharges                                  | Complies    |
| Immunity from radiated electromagnetic fields                          | IEC 61000-4-3  | 10.0 V/m; 80 MHz + 2.7 GHz, 80% AM, 1 kHz                                       | Complies    |
| Immunity from proximity field from wireless communications equipment   | IEC 61000-4-3  | List of frequencies, from 9 V/m up to 28 V/m,<br>PM (18 Hz or 217 Hz), FM 1 kHz | Complies    |
| Immunity from electrical fast transient (EFT)                          | IEC 61000-4-4  | ± 2.0 kV - on AC mains; Tr/Th - 5/50 ns, 100 kHz                                | Complies    |
| Immunity from surge                                                    | IEC 61000-4-5  | ±2.0 CM / ±1.0 kV DM on AC mains;<br>Tr/Th - 1.2/50 (8/20) µs                   | Complies    |
| Immunity from conducted disturbances induced by radio-frequency fields | IEC 61000-4-6  | 3.0, 6.0 VRMS on 230 VAC mains & wand cable;<br>0.15 + 80 MHz, 80% AM 1 kHz     | Complies    |
| Immunity from voltage dips, short interruptions and voltage variations | IEC 61000-4-11 | 0 % - 0.5 cycle & 1 cycle; 70% - 25 cycles;<br>0% - 250 cycles on AC mains      | Complies    |

## Appendix B:

## System specifications

| Item                                       | Description                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
|--------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Monitor                                    | 21.5" monitor                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| Scanner                                    | Scanner emits red laser light (680nm Class I) as well as white LED emissions, and 850nm LED emissions.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| Wireless LAN                               | LAN card provides local network communications with wireless connectivity.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| Security                                   | The iTero password is unique for each user and should not be shared with others. The iTero password policy requires high level of password strength.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| Operating power                            | 100-240VAC- 50/60 Hz – 200VA (max)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| Operating temperature                      | -8°C to 26°C / 64.4°F to 78.8°F                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| Storage/transportation temperature         | -5°C to 50°C / 23° to 122°F                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| Operating pressure & altitude              | Pressure: 520 mmHg to 771 mmHg (-69 kPa to -103 kPa)<br>Altitude: -400 feet to 10,000 feet                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| Storage/transportation pressure & altitude | Pressure: 430 mmHg to 760 mmHg (-57 kPa to -101 kPa)<br>Altitude: 0 feet to 15,000 feet                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| Relative humidity                          | Operating: 40% to 70%; Storage: 30% to 90%                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| Dimensions                                 | <div> <b>iTero HD touch monitor:</b><br/>           Height: 356 mm (-14 in)<br/>           Width: 552 mm (-21.7 in)<br/>           Depth: 65 mm (-2.5 in)         </div> <div> <b>iTero Element 5D wand:</b><br/>           Length: 346 mm (13.3 in)<br/>           Width: 50 mm (2.0 in)<br/>           Depth: 68 mm (2.7 in)         </div> <div> <b>Wheel stand:</b><br/>           Height: 1280 mm (-50 in)<br/>           Width: 645 mm (-25 in)<br/>           Depth: 625 mm (-24.5 in)         </div> <div> <b>Counter stand:</b><br/>           Height: 480 mm (-19 in)<br/>           Width: 552 mm (-21.7 in)<br/>           Depth: 220 mm (-8.7 in)         </div> |
| Net weight                                 | iTero HD touch monitor: 8.3 kg (-18.3 lbs)<br>Tero Element 5D wand: 0.47 kg (-1 lbs)<br>Wheel stand: 13.6 kg (-30 lbs)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| Shipping weight                            | -37.5 kg (-83 lbs)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |

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© 2019 Align Technology, Inc. All Rights Reserved. 206519 Rev E

# Laptop configuration **operation** **manual**

### Copyright

© 2019 Align Technology, Inc. All Rights Reserved.

The information contained in this manual is subject to change without notice.

The hardware and software described in this document are supplied under a Sales and Services Agreement and may be used only in accordance with the terms of that agreement.

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This operational manual describes the iTero® Element 5D laptop configuration Optical Impression Device. English language version. Updated May 2019.

### Trademarks

The following are trademarks and/or service marks of Align Technology, Inc. or one of its subsidiaries or affiliated companies and may be registered in the U.S. and/or other countries. Align Technology Inc., iTero

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
Align Technology B.V.  
Arlandaweg 161  
1043HS Amsterdam  
The Netherlands

### Class 1 laser compliance

This device complies with: "21 CFR 1040.10" and "EN 60825-1".

Class 1 Laser Product

### CSA compliance

This device complies with the following CSA standard for Canada and the USA:

"UL Std No. 60601-1 – Medical Electrical Equipment Part 1: General Requirements for Safety"

### FCC compliance

This device complies with Part 15 of FCC Rules and its operation is subject to the following two conditions:

1. This device may not cause harmful interference.
2. This device must accept any interference received, including interference that may cause undesired operation.

### FCC warning

Modifications to the device that are not expressly approved by the manufacturer may void your authority to operate the device under FCC Rules.

### EMC compliance

This device complies with the following EMC standard:

IEC 60601-1-2 Medical electrical equipment – Part 1-2: General requirements for basic safety and essential performance – Collateral standard: Electromagnetic phenomena – Requirements and tests".

### Safety compliance

This device complies with the following safety standard:

IEC 60601-1 Medical electrical equipment –Part 1: General requirements for basic safety and essential performance."

### CE compliance

This device complies with Council Directive 93/42/EEC for Medical Devices.

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## Symbols

The following symbols may appear on iTero Element 5D hardware components, and may also appear within this document and other iTero Element 5D literature.

Wherever this symbol appears on the device, it is recommended to refer to this document for information on proper usage of the device.

Applied part type BF. Any component on which this symbol appears is electric isolation type BF.

Separate collection of electrical waste and electronic equipment is required.

Attention: This symbol is used to highlight the fact that there are specific warnings or precautions associated with the device. Wherever this symbol appears on the device, it is mandatory to refer to safety-related information within this document.

Parts or accessories on which this symbol occurs should not be reused.

"Rx only"

CAUTION: US Federal Law restricts this device to sale by or on the order of a licensed Dentist, Orthodontist or Dental Professional. The system serves as a prescription medical device and should be operated by qualified health-care providers only.


Medical device manufacturer.

Order number.

Serial number.

Indicates the Authorized representative in the European Community.

IEC 60417-5032: Alternating current.

Indicates a medical device that needs to be protected from moisture.

Indicates the temperature limits to which the medical device can be safely exposed.

Indicates the need for the user to consult the instructions for use.

Manufacturer's batch code.

Indicates the range of atmospheric pressure to which the medical device can be safely exposed.

Indicates the range of humidity to which the medical device can be safely exposed.

Fragile, handle with care.

This side should be up.

IEC 60417-5031: Direct current.

Wand (scanning unit).

USB socket.

Electric battery.

RoHS (China).

IEC 60417-5009: STAND-BY.

CAUTION: DO NOT step on the 5D laptop configuration hub

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
### Safety instructions

Before beginning to work with the system, all users are required to read these safety instructions.

### Power supply

Power is supplied to the system via an internal medical grade power supply.

### Electric warning

- **Electric shock hazard!** Only authorized Align Technology technicians can remove external panels and covers. There are no user-serviceable parts inside.
- To avoid risk of electric shock, iTero Element 5D laptop configuration must only be connected to a supply mains with protective grounding.
- Connect the 5D laptop configuration hub device only to a laptop that is approved according to IEC60950, and to UL60950-1. Laptop and all its accessories shall be located at least 1.5m away from the patient. Do not scan a patient and touch the laptop or any of its accessories at the same time.
- Only the Align Technology scanning unit and the approved laptop should be connected to the USB sockets on the 5D laptop configuration hub.
- Only the Align Technology approved power cable should be used to connect the 5D laptop configuration hub with the AC outlet.

### Safety classifications

- Type of protection against electrical shock: Class 1.
- Degree of protection against electrical shock: Type BF.
- Degree of protection against harmful ingress of water: Ordinary.
- Equipment not suitable for use in presence of flammable anesthetic mixtures.
- Mode of operation: Continuous.

### Prescription health device

The system serves as a prescription medical device and should be operated by qualified health-care providers only.

### Scanner warnings

- The scanner emits red laser light (680nm Class 1) as well as white LED emissions, and 850nm LED emissions. Normal usage of the scanner does not present any danger to the human eye. However, doctors should refrain from shining the scanner directly into the patient's eyes.
- Avoid twisting cable, knotting cable, pulling on cable, stepping on cable.
- When the system is not in use, the scanning unit should be placed inside the cradle with the probe facing towards the cart's post and rear side of the touch screen so there will be no eye contact with the laser beam or the flickering white LED emission, and 850nm LED emissions in any case.
- The doctor should activate scanning operation only while the scanner's probe is inside the patient's mouth.
- Doctors should avoid placing the scanner in the cradle while scanning operation is still active.
- If scanner malfunction occurs, stop scanning and call service support.

### Cleaning & disinfection

- To avoid cross contamination, it is mandatory that after each patient session the disposable plastic sleeve be replaced and the scanning unit be disinfected.
- To avoid cross contamination, it is mandatory that after each patient session the gloves shall be removed and replaced.
- To avoid cross contamination discard gloves when torn, contaminated, or removed for any reason.
- Dispose of scanner sleeves according to standard operating procedures or local regulations for the disposal of contaminated medical waste.

### Unpacking & installing

The system should be unpacked and installed following Align Technology's instructions.

### Work environment

- The system should be moved between rooms with utmost care to avoid damage.
- **WARNING:** Do not block the air vents on the Scanning Unit and Base Unit.
- **WARNING:** System is intended for indoor use only. It should not be exposed to direct sunlight, excessive heat or humidity.
- The system should only be used following operating temperatures as defined in Appendix B.
- If iTero Element 5D laptop configuration has just been brought into the office from a hot or cold or humid environment, it should be set aside until it has adjusted to the room temperature to avoid internal condensation.

### Electro magnetic interference

**WARNING:** This device has been tested and found to comply with the requirements for medical devices according to standard IEC60601-1-2. This standard is designed to provide reasonable protection against harmful interference in a typical medical installation. However, because of the proliferation of radio-frequency transmitting equipment and other sources of electrical noise in the healthcare environments (e.g., cellular phones, mobile two-way radios, electrical appliances), it is possible that high levels of such interference due to close proximity or strength of source, may result in disruption of performance of this device.

ESD may cause an interference during scanning. In case of ESD interference the system will stop scanning / displaying proper viewfinder image. If the failure occurs within scan, the user may be required to press the wand scan button to restart the scan. The viewfinder will always recover automatically.

### General

**WARNING:** No modification of this equipment is allowed.

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## Table of contents

|                                                                             |           |
|-----------------------------------------------------------------------------|-----------|
| <b>Chapter 1: Introduction</b>                                              | <b>2</b>  |
| About this document                                                         | 2         |
| Intended use / indication for use                                           | 2         |
| Benefits of the iTero Element 5D laptop configuration system                | 2         |
| The iTero Element 5D laptop configuration user interface                    | 3         |
| <b>Chapter 2: Basic hardware features</b>                                   | <b>4</b>  |
| Front view of system                                                        | 4         |
| Scanning unit (wand)                                                        | 5         |
| Mobile system packaging (for transport)                                     | 6         |
| <b>Chapter 3: Assembly instructions</b>                                     | <b>7</b>  |
| Step 1: Assembly                                                            | 7         |
| Step 2: Download iTero Element 5D laptop configuration software             | 8         |
| Step 3: Make it Mine process                                                | 9         |
| <b>Chapter 4: Operating instructions</b>                                    | <b>10</b> |
| Web cam                                                                     | 10        |
| <b>Chapter 5: Scanner handling, cleaning, and disinfection instructions</b> | <b>11</b> |
| Handling of the scanning unit (wand)                                        | 11        |
| Handling of the scanning unit cable                                         | 11        |
| Cleaning the scanning unit                                                  | 11        |
| Disinfecting the scanning unit                                              | 11        |
| <b>Chapter 6: Changing sleeves between patients</b>                         | <b>12</b> |
| Replacing disposable sleeves                                                | 12        |
| Scanner sleeves                                                             | 13        |
| <b>Chapter 7: Clinic LAN network guidelines</b>                             | <b>14</b> |
| Introduction                                                                | 14        |
| Preparations                                                                | 15        |
| Router guidelines                                                           | 15        |
| Internet connection guidelines                                              | 15        |
| Firewall                                                                    | 15        |
| Wi-Fi tips                                                                  | 15        |
| <b>Appendix A: EMC declaration</b>                                          | <b>16</b> |
| <b>Appendix B: System specifications</b>                                    | <b>17</b> |

**Chapter 1:**

# Introduction

**About this operation manual**

The iTero Element 5D laptop configuration system is delivered as a proprietary, laptop based workstation for performing intra-oral scans in the doctor's office. This operational manual describes how to start and shut down the system, how to correctly handle the Scanning Unit/Wand and cable, and how to clean the Scanning Unit and replace its sleeves between patients.

**Intended use / indication for use**

The iTero Element 5D laptop configuration is an intra-oral scanner with the following features and intended/indicated for use:

1. The optical impression (CAD/CAM) feature of iTero Element 5D laptop configuration is intended/indicated for use to record the topographical images of teeth and oral tissue. Data generated from iTero may be used in conjunction with the production of dental devices (e.g., aligners, braces, appliances, etc.) and accessories.
2. iTero Element 5D laptop configuration software is used with the iTero scanner in capturing 3D digital impressions of teeth, oral soft tissue and structures, and bite relationship. iTero software controls the processing of the data, facilitating the integration of data, and exporting of the data for CAD/CAM fabrication of dental restorations, orthodontic devices, abutments, and accessories. In addition to scan data, various patient and case information can be imported/exported or used for simulation purposes. Other functions are available for verification and service of the system and to serve as an order management tool.
3. The Element 5D laptop configuration NIRI functionality, is a diagnostic aid for the detection of interproximal caries lesions above the gingiva and for monitoring the progress of such lesions.

**Benefits of the iTero Element 5D laptop configuration system**


The iTero Element 5D laptop configuration system provides important advantages over existing crown-production methods, including powder-free scanning, greater crown-production accuracy, and immediate feedback during the scanning process.

Refer to our website [www.itero.com](http://www.itero.com) to learn how the iTero Service can enhance your business by increasing patient satisfaction, improving clinical outcomes, and enhancing office efficiency.

**Software installation instructions**

Please consult the list of approved laptops and approved Antivirus applications at [www.iTero.com](http://www.iTero.com).

To ensure optimal performance use of the iTero application, only approved laptops must be used with the iTero Element 5D laptop configuration system.


## The iTero Element 5D laptop configuration user interface

The iTero Element 5D laptop configuration system provides an intuitive user interface for performing digital scans for Restorative or Orthodontic use. The doctor is guided through the scanning sequence by means of visual and text assistance. The laptop's touch screen and the wand buttons are used to respond to screen instructions during the scanning process.


One tap on the question mark will enable a transparent Help overlay that will provide a brief overview. Please note that the Headset image appears instead of the question mark while in this view. Tap anywhere to close the help screen and return to the relevant screen.


Chapter 2:

## Basic hardware features

### Front view of the system

## Scanning unit (wand)

### Mobile system packaging (for transport)

#### Transportation

Use the supplied case to move the system between offices.

#### Moving the system

To ensure maximum system protection, it is recommended to move the system with care. Follow these instructions for relocating the system:

1. Attach the blue protective sleeve onto the wand.
2. Verify the scanning unit (wand) sits well inside the case.
3. Unplug the 5D laptop configuration hub and laptop AC power cables from the wall outlet.
4. Place all items in the case (see picture above).
5. Make sure the case is kept dry to protect the system components from humidity.

### Chapter 3:

## Assembly and software installation

### Step 1: Assembly

Please follow the instructions below to assemble your iTero Element 5D laptop configuration scanner.

AC power


Wand  
(scanning unit)

USB port

- A** Hub and hub power cord
- B** Cradle
- C** Wand and wand cable
- D** USB cable to connect laptop and hub

1. Place the wand in the cradle

2. Connect hub power cord to the hub

3. Connect USB cable to the hub

4. Connect USB cable to the laptop

5. Connect wand cable to the hub

6. Plug the hub power cord into the AC power outlet

#### Two important notes:

- The hub must be connected to an AC wall outlet at all times
- The laptop should be connected to an AC wall outlet for intraoral scanning

## Step 2: Download iTero Element 5D laptop configuration software

For proper software installation and configuration of your iTero Element 5D laptop configuration scanner, please ensure your iTero wand is secure in the cradle and connected to the hub, and the hub is connected to the laptop. Please make sure your laptop is plugged into the AC wall socket during the entire software installation.

1. Laptop setup
  - A. Turn on the laptop.
  - B. Update the Windows to the latest version.
2. iTero software download
  - A. Open the registered email inbox in the browser. Look for the email "Your iTero was shipped" which includes instructions.
  - B. Click on the link to get to the page with software downloads.  
or
  - C. Alternatively, open the browser and go to (type the address into your browser) [download.itero5d.com](http://download.itero5d.com).
3. iTero software installation
  - A. Click on the "Get started" button.
  - B. Run the downloaded installation file.
  - C. Follow the instructions on the screen to complete the download of the iTero software.

It is highly recommended to regularly check for Windows updates and install them. Please enable updates.

It is also recommended to exit to Windows and shut down the computer at the end of each day to allow updates to install.

Before installing the iTero Element 5D software for the first time, please install all available Windows updates. New Windows computers should apply updates automatically.

To check for Windows Updates, open "Settings" (Windows key + I), choose "Update & Security", and then "Windows Update". Click "Check for updates" to see if there are new updates available.

**It is recommended to close other software programs while using the iTero Element 5D laptop configuration system for scanning.**


The installation process will take about 10 minutes to complete

- Step 1: Download the software  
Step 2: Connect the wand with laptop  
Step 3: Customize your account - Make it mine


**Get started**

### Step 3: Make it Mine process

1. Select language of preference and tap on the Make it Mine button to start the Wizard.

2. Follow the Wizard instructions on the screen to complete the customization of the iTero Element 5D.

## Chapter 4:

# Operating instructions

Make sure the laptop is connected to an AC wall outlet, to ensure the battery level does not reach its low battery status, at which time scanning will no longer be possible.

Power up the laptop to start iTero program.

It is recommended to keep the system in operation during office hours to allow background file transfers between the doctor's office, the doctor's partnered labs, and the Align Technology Center. It is recommended to shut down the system at the end of the day, and to reboot in the morning.

### Web cam

Connect the web cam to the laptop following the product instructions.

**Chapter 5:**

## Scanner handling, cleaning, and disinfection instructions

**Handling of the scanning unit (wand)**

The scanning unit contains delicate components and should be handled with care.

**Handling of the scanning unit cable**

- The scanner cable should be treated with care to avoid possible damage.
- Between patient sessions, it is recommended to undo any twists and knots in order to relieve all tension from the scanner cable.

The scanning unit requires proper cleaning & disinfection before first use and before each additional use, and no later than 30 minutes after the last scan.

**Cleaning and disinfecting the scanning unit**

The scanning unit (wand) should be cleaned and disinfected as follows:

- Soak a lint free cloth in ready to use CaviCide 1, and squeeze until the cloth is moist.
- Wipe the wand thoroughly to remove gross debris.
- The entire device needs to be visually inspected to ensure that no residual or debris remains prior to continuing to the next step.
- Soak two additional lint free cloths in Cavicide 1, and squeeze until the cloth is moist. Then wipe the device thoroughly, ensuring with special care that all the surfaces, to all edges and slots, are covered. Contact time will be at least 2 minutes.
- Then wet lint free cloths with distilled water and wipe all surfaces for at least 15-30 seconds.
- Then use dry lint free cloths to dry the surfaces.

**Maintenance**

The industrial camera is maintenance-free.

## Chapter 6:

Changing sleeves  
between patients

iTero Element 5D laptop configuration sleeve is intended for single patient use and shall be replaced after each patient with the purpose to avoid cross contamination.

For sleeve replacement proceed with the steps below.

**CAUTION:** Dispose of scanner sleeves according to standard operating procedures or local regulations for the disposal of contaminated medical waste.

## Replacing disposable sleeves

**Step 1**

When pulling a sleeve OFF or ON, hold the center of the sleeve.

**Step 2**

Press slightly on both sides of the disposable sleeve, pull the sleeve slowly off the scanning unit and discard.

**Step 3**


Gently slide on new sleeve onto scanning unit until it clicks into place.

**WARNING: Optical surface!**

DO NOT touch the optical surface. Contact may cause damage. If cleaning is necessary, use the wipes and anti-static cloth found inside the sleeves box. For proper use, refer to the directions found in the scanner sleeves box.

## Scanner sleeves

There are two types of sleeves intended for use with the scanner unit (wand):

### Disposable sleeve

The white sleeve is a single use sleeve for patient scanning. Always replace the white sleeve on the scanning unit between patients to avoid cross contamination. Please dispose of the white sleeve after every patient.

### Protective sleeve

The blue protective sleeve is used to protect the optical surface lens when the scanning unit is not in use. Please keep the blue sleeve in a safe place so that it does not get lost or damaged.

### Scanner sleeves packaging box

Scanner sleeves may be ordered online in boxes of 25 from the iTero store.

[www.store.itero.com](http://www.store.itero.com)

## Chapter 7:

## Clinic LAN network guidelines

### Introduction

The iTero Element 5D scanner uses the Wi-Fi internet connectivity in order to send and retrieve scans to and from the iTero cloud.

As a recommendation, it is always best to have the state-of-art available technology. Here are some helpful guidelines for the best Wi-Fi connection:

#### Levels of Wi-Fi Internet Connectivity

Excellent  
>-50 dBm

Good  
-50 to -60 dBm

Fair  
-60 to -70 dBm

Weak  
<-70 dBm

- **IMPORTANT:** In order to achieve the best performance of your iTero Element 5D laptop configuration scanner, ensure that the Wi-Fi signal strength is "Excellent" or at least "Good".
- **CAUTION:** While scanning a patient, do not connect a LAN cable to iTero Element 5D laptop configuration – it is forbidden due to safety hazard reasons.

## Preparations

- The required Modem/Router should be configured with WPA2 Security standard, including a password.
- Ensure that your IT professional staff would be available when the scanner installation is planned to take place.
- Make sure that your Wi-Fi SSID credentials are available: Login & password.
- The minimal Wi-Fi strength signal for the system should display at least two "stripes", as shown in chapter 2, above).
- Diagnostic tool under "Settings", or the Connectivity Tool are suggested below.
- Test the local Wi-Fi connection with any PC computer, using the iTero connectivity tool (run the test as near as possible to the scanner location).
- Connectivity Tool Download at [fci.orthocad.com/download/AlignSupport/ConnectivityCheck2.0.exe](http://fci.orthocad.com/download/AlignSupport/ConnectivityCheck2.0.exe)
- Following are some suggestions for the office IT person, regarding what should be considered in order to prevent issues such as access or connectivity to/with the iTero scanner:
  1. Hostname recommendations related to Align services listening to ports 80 and 443.
  2. Do not prevent FTP communication since the scanner sends specific file types (.3ds and .3dc/.3dm).
  3. Disable any proxy client for data communication through TCP/IP.
  4. Do not add the scanner to any domain group.
  5. Do not run any group policy on the scanner as it may disrupt its proper functioning.

## Router guidelines

Minimum Standards: 802.11N / 802.11AC

## Internet connection guidelines

In order to achieve the best performance of your iTero Element 5D laptop configuration scanner, ensure that your internet connection upload speed is at least 1Mbps per scanner. Also, note that any additional devices connected to the internet in parallel to the scanner may affect the scanner's performance.

## Firewall

Open Ports (in case Firewall is working):

- a. 80 - HTTP - TCP
- b. 443 - HTTPS - TCP

## Wi-Fi tips

Wi-Fi routers allow you to access your internet system using a Wi-Fi connection from essentially any place within the functional range of the wireless network. Nevertheless, the number, depth and position of walls, ceilings, or additional partitions that the wireless signals must travel through may limit the range and strength of the signal. Normal signals vary depending on the material types and background RF (radio frequency) noise in your home or business.

6. Be sure to have a minimal number of walls and ceilings between the router and other network devices. Each barrier can reduce your adapter's range by 1-3 meters (3-9 feet).
7. Be sure to have a straight line, free of any partition, between network devices. Even a wall that seems rather thin can block a signal of 1 meter (3 feet) if the wall angle is shifted by only 2 degrees. To achieve the best reception, place all the devices so that the Wi-Fi signal travels straight through (90°) a wall or partition (instead of at an angle).
8. Construction materials make a difference. A solid metal door, or aluminum nails, can be very dense and may have an adverse effect on a Wi-Fi signal. Try to position access points, wireless routers, and computers so that the signal travels through drywalls or open doorways. Materials and objects such as glass, steel, metal, walls with insulation, water tanks (aquariums), mirrors, file cabinets, brick, and concrete may reduce your wireless signal.
9. Keep your iTero product away (at least 3-6 feet or 1-2 meters) from electrical devices or appliances that generate RF noise.
10. If you are using 2.4GHz cordless phones or X-10 (wireless products such as ceiling fans, remote lights, and home security systems), your wireless connection may severely degrade or entirely drop. The base of many wireless devices transmits an RF signal, even if the device is not in use. Position your other wireless devices as far as possible from your scanner and router.
11. In your area, there may be more than one active wireless network. Each network uses one or more channels. If the channel is near your system channels, the communication may gradually decline. Ask your IT department to check this, and if required, change the channel numbers used by your network.

## Appendix A:

## EMC declaration

## Summary of EMC test results for iTero Element 5D laptop configuration

| Test                                                                   | Standard       | Class / Severity level                                                                                                                                                                             | Test result |
|------------------------------------------------------------------------|----------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------|
| <b>Documentation (IEC 60601-1-2 sections 4 and 5)</b>                  |                |                                                                                                                                                                                                    |             |
| General requirements for EMC                                           | section 4.1    | --                                                                                                                                                                                                 | Complies    |
| External labels                                                        | section 5.1    | --                                                                                                                                                                                                 | Complies    |
| Conformity of users' manual                                            | section 5.2.1  | --                                                                                                                                                                                                 | Complies    |
| Accuracy of technical description                                      | section 5.2.2  | --                                                                                                                                                                                                 | Complies    |
| <b>Emission (IEC 60601-1-2 section 7)</b>                              |                |                                                                                                                                                                                                    |             |
| Conducted emission Freq. range: 150 kHz - 30 MHz                       | CISPR 11       | Group 1 Class B<br>230, 120 & 100 VAC mains @ 50 Hz;<br>220 VAC mains @ 60 Hz                                                                                                                      | Complies    |
| Radiated emission Freq. range: 30 - 1000 MHz                           | CISPR 11       | Group 1 Class B                                                                                                                                                                                    | Complies    |
| Harmonic current emission test                                         | IEC 61000-3-2  | 230 VAC mains @ 50 Hz<br>& 220 V @ 60 Hz                                                                                                                                                           | Complies    |
| Voltage changes, voltage fluctuations and flicker test                 | IEC 61000-3-3  | 230 VAC mains @ 50 Hz<br>& 220 V @ 60 Hz                                                                                                                                                           | Complies    |
| <b>Immunity (IEC 60601-1-2 section 8)</b>                              |                |                                                                                                                                                                                                    |             |
| Immunity from electrostatic discharge (ESD)                            | IEC 61000-4-2  | 8 kV contact discharges & 15 kV air discharges                                                                                                                                                     | Complies    |
| Immunity from radiated electromagnetic fields                          | IEC 61000-4-3  | 10.0 V/m; 80 MHz, 2.7 GHz, 80% AM, 1 kHz                                                                                                                                                           | Complies    |
| Immunity from proximity field from wireless communications equipment   | IEC 61000-4-3  | List of frequencies, from 9 V/m up to 28 V/m,<br>PM (18 Hz or 217 Hz), FM 1 kHz                                                                                                                    | Complies    |
| Immunity from electrical fast transient (EFT)                          | IEC 61000-4-4  | ± 2.0 kV - on AC mains; Tr/Th - 5/50 ns, 100 kHz                                                                                                                                                   | Complies    |
| Immunity from surge                                                    | IEC 61000-4-5  | ±2.0 CM / ±1.0 kV DM<br>on 230 VAC mains @ 50 Hz<br>& 220 VAC mains @ 60 Hz;<br>Tr/Th - 1.2/50 (8/20) ms                                                                                           | Complies    |
| Immunity from conducted disturbances induced by radio-frequency fields | IEC 61000-4-6  | 3.0, 6.0 VRMS<br>on 230 VAC mains & wand cable;<br>0.15± 80 MHz, 80% AM @ 1 kHz                                                                                                                    | Complies    |
| Immunity from voltage dips, short interruptions and voltage variations | IEC 61000-4-11 | On 230 VAC & 100 VAC mains @ 50 Hz:<br>0% - 0.5 cycle & 1 cycle;<br>70% - 25 cycles; 0% - 250 cycles<br>on 220 VAC mains @ 60 Hz:<br>0% - 0.5 cycle & 1 cycle;<br>70% - 30 cycles; 0% - 300 cycles | Complies    |

## Appendix B:

## System specifications

| Item                                       | Description                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
|--------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Monitor                                    | Laptop monitor                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| Scanner                                    | Scanner emits red laser light (680nm Class II) as well as white LED emissions.                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| Operating power                            | 100-240VAC, 50/60 Hz – 40VA (max)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| Operating temperature                      | 18°C to 26°C / 64.4°F to 78.8°F                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| Storage/Transportation temperature         | -5° to 50°C 23° to 122°F                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| Operating pressure & altitude              | Pressure: 520 mmHg to 760 mmHg (-69 kPa to -101 kPa)<br>Altitude: 0 feet to 10,000 feet                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| Storage/Transportation pressure & altitude | Pressure: 430 mmHg to 760 mmHg (-57 kPa to -101 kPa)<br>Altitude: 0 feet to 15,000 feet                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| Relative humidity                          | Operating: 40% to 70%; Storage: 30% to 90%                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| Dimensions                                 | <div> <b>iTero Element 5D laptop configuration hub:</b><br/> Length: 206 mm (-8 in)<br/> Width: 94 mm (-3.7 in)<br/> Depth: 36.5 mm (-1.4 in) </div> <div> <b>iTero Element 5D wand unit:</b><br/> Length: 346 mm (13.3 in)<br/> Width: 50 mm (2.0 in)<br/> Depth: 68 mm (2.7 in) </div> <div> <b>Cradle:</b><br/> Length: 262 mm (-10 in)<br/> Width: 89 mm (-3.5 in)<br/> Depth: 52 mm (-2 in) </div> <div> <b>Case:</b><br/> Height: 326.5 mm (-13 in)<br/> Width: 455 mm (-18 in)<br/> Depth: 184 mm (-7 in) </div> |
| Net weight                                 | iTero Element 5D laptop configuration hub: ~0.5 kg (-1 lbs)<br>Scanning Unit: 0.47 kg (-1 lbs)<br>Empty Case: ~2 kg (-4.5 lbs)                                                                                                                                                                                                                                                                                                                                                                                          |
| Shipping weight                            | ~9 kg (-17.6 lbs)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |

# iTero element 5D

Align Technology, Inc.  
2820 Orchard Parkway  
San Jose, CA 95134  
USA

Align Technology B.V.  
Arlandaweg 161  
1043HS Amsterdam  
The Netherlands

Align, Invisalign, iTero, and iTero Element among others, are trademarks and/or service marks of Align Technology, Inc. or one of its subsidiaries or affiliated companies and may be registered in the U.S. and/or other countries.

© 2019 Align Technology, Inc. All Rights Reserved. 206520 Rev E

iTero element 5D

# Imaging system **instructions** **for use**

## Copyright

Copyright © by Align Technology, Inc. 2019. All rights reserved. The information contained in this manual is subject to change without notice.

The hardware and software described in this manual are supplied under a Sales and Services Agreement and may be used only in accordance with the terms of that agreement.

No part of this document may be reproduced, photocopied, stored in a retrieval system, or transmitted in any manner (electronic or mechanical) for any purpose other than the customer's normal usage, without the prior written permission of Align Technology Inc.

This document describes the iTero Element 5D Optical Impression Device. English language version. Updated in November 2019.

## Trademarks

The following are trademarks and/or service marks of Align Technology, Inc. or one of its subsidiaries or affiliated companies and may be registered in the U.S. and/or other countries:  
Align Technology Inc., iTero

Any other trademarks or registered trademarks appearing in this document are the property of their respective owners.

## USA office

Corporate Headquarters  
Align Technology, Inc.  
2820 Orchard Parkway  
San Jose, California 95134

[www.aligntech.com](http://www.aligntech.com)

Tel: +1 (408) 470-1000  
Fax: +1 (408) 470-1010

Align Technology Ltd.  
3 Ariel Sharon Boulevard  
Or-Yehuda 6037606  
Israel

Tel: +972 (3) 634-1441  
Fax: +972 (3) 634-1440

## Netherlands office

International Headquarters  
Align Technology B.V.  
Arlandaweg 161  
1043 HS Amsterdam


Tel: +31 (0) 20-586-3600  
Fax: +31 (0) 20-586-3751

## Customer support

Tel: +1 (800) 577-8767


✉ F-mail:  
[erosupport@aligntech.com](mailto:erosupport@aligntech.com)

E-mail:  
[orthocadsupport@aligntech.com](mailto:orthocadsupport@aligntech.com)

Align Technology B.V.  
Arlandaweg 161  
1043HS Amsterdam  
The Netherlands

**Table of contents**

|                                                                      |    |
|----------------------------------------------------------------------|----|
| Introduction                                                         | 2  |
| Intended use                                                         | 2  |
| Chapter 1: System registration                                       | 3  |
| Demo Mode                                                            | 3  |
| Make it Mine                                                         | 3  |
| Chapter 2: Logging in to the scanner application                     | 6  |
| Login                                                                | 6  |
| Create account                                                       | 6  |
| Forgot password                                                      | 6  |
| Chapter 3: iTero Element 5D main features                            | 7  |
| Home screen                                                          | 7  |
| Main desktop icons                                                   | 7  |
| Utility icons                                                        | 8  |
| Touch screen gestures                                                | 9  |
| Wand controls                                                        | 10 |
| Demo Mode                                                            | 10 |
| Chapter 4: Messages                                                  | 12 |
| Chapter 5: Orders                                                    | 13 |
| Chapter 6: Scan                                                      | 14 |
| Start a new scan                                                     | 14 |
| Chapter 7: Element 5D screens                                        | 17 |
| Review tool for iTero Element 5D imaging system                      | 17 |
| Using NIRI                                                           | 18 |
| Chapter 8: Patients                                                  | 21 |
| Chapter 9: iTero TimeLapse                                           | 22 |
| Chapter 10: iTero Element 5D operational manual                      | 24 |
| Chapter 11: iTero Element 5D laptop configuration operational manual | 24 |

## Introduction

The iTero Element 5D intra-oral scanner is Align Technology's latest generation integrated scanner and software system. It enables doctors to take 3D digital scans of patients' teeth and bite, make adjustments of the scan in real time and transmit the results wirelessly.

This user guide provides general information and an overview of the application.

## Intended use

The iTero Element 5D is an intra-oral scanner with the following features and intended/indicated for use:

1. The optical impression (CAD/CAM) feature of iTero Element 5D is intended/indicated for use to record the topographical images of teeth and oral tissue. Data generated from iTero may be used in conjunction with the production of dental devices (e.g., aligners, braces, appliances, etc.) and accessories.
2. iTero Element 5D software is used with the iTero Element 5D scanner in capturing 3D digital impressions of teeth, oral soft tissue and structures, and bite relationship. iTero software controls the processing of the data, facilitating the integration of data, and exporting of the data for CAD/CAM fabrication of dental restorations, orthodontic devices, abutments, and accessories. In addition to scan data, various patient and case information can be imported/exported or used for simulation purposes. Other functions are available for verification and service of the system and to serve as an order management tool.
3. The iTero Element 5D NIRI functionality, is a diagnostic aid for the detection of interproximal caries lesions above the gingiva and for monitoring the progress of such lesions.

This all-in-one system is available on a monitor with a fully interactive, touch screen display and an easy-to-use scanner. The topography of a patient's teeth can be viewed on the screen as they are being scanned, and determining the degree of occlusion of the bite can be analyzed as

the scan is completed.

In conjunction with the scanner, the application facilitates the following tasks by:

- Capturing 3D digital impressions of teeth, oral soft tissue and structures as they are being scanned.
- Enabling real-time viewing of the visual bite relationship once the scan is completed.
- Managing the processing of the data.
- Sending the data wirelessly to designated providers for the fabrication of dental restorations, orthodontic devices, abutments, and accessories, etc.
- Using the data for simulation purposes.

The application also includes administrative capabilities which enable the user to:


- Start a new scan.
- Place a new order for an existing patient or new patient.
- View orders in progress.

Review and/or track past orders.

## Chapter 1:

## System registration


**1** Demo Mode**2** Make it Mine**Demo Mode**


Press Demo Mode to familiarize yourself with the scanner's features and perform practice scans without submitting the scans. Demo Mode is designed for training new staff members and practice. Demo Mode is available anytime for dental practices to train on an iTero Element 5D device, for scanning techniques, how-to guidelines for prescription forms, case types, and to familiarize with the iTero interface. Demo Mode will feature all the aspects of the scanning process and include a wide variety of cases. Sample cases offered will include clinical cases, Invisalign cases, and restorative cases for your convenience.

When Demo Mode is in use for practice scanning, a lightly striped background and red tag in the upper left-hand corner will indicate Demo Mode. **Any scans captured in Demo Mode cannot be saved or submitted for patient treatment.** Demo Mode is available to access on the Welcome screen or at any point by pressing the iTero Element logo button on the home screen.


**Make it Mine**

Press Make It Mine to register your scanner. Please have your User Name, User Password, and Company ID ready to complete the registration process. You will receive an email from an iTero representative with login credentials and detailed information on the how to Make it Mine procedure.


1. Select the appropriate Network Connection

3. Register the system


2. Network connection verification

4. System is registered


5. Accept the License Agreement

7. System is now ready

6. Updating SW to the latest version


**Chapter 2:**

# Logging in to the scanner application

**Login**

When the scanner is powered on, the login screen for the iTero Element 5D will appear.

- 1 User information
- 2 Create Account
- 3 Forgot Password
- 4 Login

Please have your MyiTero account information ready when logging in to the iTero Scanner. You will need your name, account email, and password. Fill out all the necessary fields and click the Login button.

**Create account**

To create an account, please press the bottom button labeled Create Account. Enter an email and create an account to use on the iTero Scanner and MyiTero.com.

**Forgot password**

If you forgot your password, press the lower button labeled Forgot Password. Enter your registered account email and a temporary password will be sent to your email. Alternatively, answer your predetermined security question and gain access to the account from the scanner.


## Chapter 3:

## iTero Element 5D main features

## Home screen

The iTero Element 5D home screen provides access to all the main features of the iTero software.

Below is an image of the iTero Element 5D application's home screen.

- 1 Main home screen icons
- 2 Utility icons

## Main desktop icons

There are four main home screen icons for the iTero Element 5D software.

**New Scan:** To begin a new scan for a patient, select the New Scan icon. This home screen icon navigates to the iTero Rx form to fill out the prescription form for a patient and begin scanning.

**Patients:** To look up any patient on your system, select the Patients icon. This home screen icon navigates to the Patients page. Search for a specific patient, scroll through recently scanned patients, review scans, access iTero TimeLapse, view patient history, add new scans, or add a treatment plan.

**Orders:** To view the status of your orders, select the Orders icon. This home screen icon navigates to a list of orders in progress or past orders.

**Messages:** To read messages, select the Messages icon. This home screen icon navigates to the message inbox, where you may read on the latest notifications from the Align team regarding new updates, improved features, enhanced tools, and latest developments.

## Utility icons


**Battery:** To view the charge on the iTero external battery, press the Battery icon. A battery indicator will appear and specify how long the battery has until the iTero Element 5D may need to be charged again. When the battery icon displays a lightning bolt symbol, the battery is charging. Note: this feature only appears in supported iTero Element 5D scanners.

**Learning Center:** To access training materials and educational videos for your iTero Element 5D, press the Learning Center icon.

**Lock:** To log out of an account whenever the iTero Element 5D is not in use, press the Lock icon. This will help ensure that the dental practice is HIPAA compliant and that all medical information is secure.

**Settings:** To adjust preferences for your iTero Element 5D like wand configuration, localization, user settings, and more, press the Settings icon.

**Help:** To navigate icons for home screen and helpful hints, press the Help icon. A transparent Help overlay appears over the screen to aid in navigation of features and tools. At any time, tap the Help icon for assistance. Tap anywhere to close the Help screen and return to the relevant screen. Customer support is available on every Help overlay.

Touch screen gestures

iTero Element 5D application supports touch screen (also known as multi-touch) gestures. These gestures are predefined motions used to interact with multi-touch devices.

Examples of common touch screen gestures:

Tap

Double tap

Long press

Scroll


Rotate


Swipe

2 finger tap

Zoom in/pinch

Zoom out

Pan

## Wand controls

- 1 Side buttons
- 2 Touchpad

**Side buttons:** Press and release either of the side buttons on the wand to start or stop the scanning process. Press the buttons simultaneously to activate the touchpad.

**Touchpad:** Once activated, swipe left or right on the touchpad to navigate between scan segments. Press and hold the touchpad will enable the gyro function to rotate the model on the screen.

## Demo Mode

To access demo mode, press the iTero Element logo on the upper left side of the home screen. A menu will appear to press Demo Mode or to cancel and return to the home screen. Demo Mode will appear with all the same functionality as when taking a normal scan. No data will be permanently saved.


- 1 Demo Mode
- 2 Exit to Windows
- 3 Cancel

**Demo Mode:** Initiate Demo Mode for practice without saving any data.

**Exit to Windows:** To Exit the iTero interface and return to Windows. Note: Only available on iTero laptop configuration devices.

**Cancel:** To exit from the screen and return to normal use.

**Demo Mode tag:** The red tag indicator when using Demo Mode.


**Striped background:** A striped background will appear while scanning. It is designed to prevent scanning real patients in Demo Mode.

Practice any portion of the patient experience. **Any scans taken during Demo Mode cannot be submitted for treatment.** To exit Demo Mode, please go back to the home screen and press the iTero Element logo. The menu for Demo Mode will appear in order to exit Demo Mode.

## Chapter 4:

# Messages

1 Inbox

2 Mark as Unread

Press the Messages icon from the home screen to view notifications, updates and other messages from Align Technology. On the left side of the page, search for specific messages by subject title quickly or, press on the right side of the screen and scroll down the page to search and find a specific message. To mark any message as unread, press Mark as Unread.

Chapter 5:  
Orders

Tap on the Orders icon to view or review In Progress or Past Orders.

Orders

Dr. Sarah Smith - iTero

1

In Progress

2

Past Orders

In Progress

| ID       | Customer Name | Start Date/Time     | Scan Type   | Status     |
|----------|---------------|---------------------|-------------|------------|
| 11358871 | Jane White    | 01/30/2017 15:10:50 | Restorative | Rx Created |
|          | Anna Miller   |                     | Restorative | Scanning   |
|          | Jane Johnson  |                     | Record      | Rx Created |
|          | James Doe     | 01/27/2017 10:47:47 | Record      | Scanning   |
|          | John Smith    |                     | Record      | Rx Created |
|          | Bill Williams |                     | Record      | Rx Created |
|          | Sarah Jones   | 01/27/2017 10:30:47 | Record      | Scanning   |

Past Orders

| ID       | Customer Name | Start Date/Time | Scan Type | Status    |
|----------|---------------|-----------------|-----------|-----------|
| 11358871 | Jane Doe      | 01/27/2017      | Record    | Completed |

View Rx

Viewer

Add Rx

Invisalign Outcome Simulator

Progress Assessment

**In Progress:** Scans that have not been submitted.

**Past Orders:** Orders that have been submitted. Press on any order to check the status of the scan.

Clicking on a past order shows the following options, depending on the case type.

**View Rx:** View the prescription for this order.

**Viewer:** View the scans for this order.

**Add Rx:** Add a prescription for this order (only applicable for Orthodontic orders and available up to 21 days after a scan).

Invisalign doctors may also select Invisalign features such as Invisalign Outcome Simulator.

## Chapter 6:

# Scan

### Start a new scan

To start a new scan, press the New Scan icon from the home screen or press the new scan button from the Patient page.

New Scan

**Prescription icon:** The first step in the New Scan workflow. Fill out the prescription form with patient information, case type, and any notes you include.

**Scan icon:** The second step in the New Scan workflow. Scan the patient for a new scan to be captured.

**View and Review icon:** The third step in the New Scan workflow. Inspect the scan and examine that all adequate data has been captured.

**Send Case icon:** The final step in the New Scan workflow. Sign the confirmation form and send the case to its destination for storage or additional work needed.

- 1 Patient information
- 2 Case Type
- 3 Notes

For Restorative cases, select the case type and restoration needs for the appropriate indication. For Orthodontic cases, such as Invisalign or iRecord, select the case type and enter in patient information.

- 1 Select a lab to send the case
- 2 Select the restoration type from the dropdown

After selecting the case type for a patient, please begin scanning.

To start a scan, tap the scan icon.

- 1 Real time modeling
- 2 View finder
- 3 Color
- 4 Scan feedback
- 5 Toggle view
- 6 Enlarge view

**Scan feedback:** To view any missing tooth, soft tissue, or oral structure during the scan.

**Color:** To switch from scan to color view.

**View finder:** View the wand's current positioning.

**Toggle View:** Toggle between displaying a color image or a NIRI image in the viewfinder

**Enlarge View:** Switch to a large viewfinder in the middle of the screen

After reviewing the scan, sign the Rx form and send the scan to complete the case.

**Signature box:** Sign your signature here with a finger. Check "Save Signature" to store a signature away for future ease and use. Check "Don't show again" to eliminate this step.

**Clear:** To clear a signature from the signature box.

**Confirm & Send:** Confirm scan has been reviewed and send the scan to a lab or iTero storage.

## Chapter 7:

## iTero Element 5D screens

In addition to scanning, the iTero Element 5D scanner provides a Near Infra-Red Imaging (NIRI) mode. The NIRI and ViewFinder color images are captured during the scan, and can be reviewed, once scanning has completed.

**Review tool for iTero Element 5D imaging systems**

After scan is completed, go to View and Review icon, inspect the scan and examine that all adequate data has been captured, including NIRI and color images, on specific locations in the mouth.

**1** Review tool**2** Loupe

**Review tool:** Displays the Color and NIRI images captured during scan

**Loupe:** When loupe is dragged to the scan, the Color and NIRI images are displayed according to the specific location of the loupe on the scan

- 1 Zoom in button
- 2 Brightness and contrast toolbar

**Zoom in button:** Viewfinder window will be enlarged and only the specific image will be displayed in the viewfinder.

**Brightness and contrast toolbar:** Opens the brightness and contrast adjustment toolbar.


After review and confirmation that the scan has been captured with sufficient data, the following features can be used to view the scan in color or monochrome, edit the scan, remove any extra scanned artifacts, and erase additional scanned materials.

### Using NIRI

NIRI provides imaging of the internal tooth structure.

Structure translucency translates in the NIRI image to brightness level; the darker the object the higher is its translucency and vice versa.

Healthy enamel appears dark, translucent, while dentin or caries are less translucent and appear brighter.


Proximal lesion will appear as bright stripe/wedge-like features as in the image above, in this case proximal lesion reaching all the way to the dentin.

Changing viewing angle around such area may provide more information regarding its depth, while using the color image may help differentiate the diagnosis as well as the lesions severity.

In addition to the NIRI, the color images are also provided to give more information regarding the tooth's health and help with the diagnosis.

The NIRI above demonstrates two proximal lesions circled in red, reaching all the way to the dentin, while the color image helps determine that the lesions are proximal and not occlusal, as no white/brown spots appear in the correlating areas of the occlusal surface in the color image. The dark lines marked by the blue arrows indicate occlusal fissures.

Note: Highly opaque teeth will have low translucency of the enamel making it appear brighter, this may cause difficulty in understanding internal tooth's structures, and differentiating between enamel and dentin, thus affecting the ability to detect proximal caries. In addition, several dental restorations such as amalgam fillings are not translucent for the NIRI, meaning they may mask lesions underneath them.

NIRI along with the color images can assist in the detection of proximal caries.

## Proximal caries

Scoring provided according to ICDAS-II standard

- Sound enamel, will appear without any noticeable bright spots, meaning it will resemble a dark ring circling a bright inner circle representing the dentin.
- Enamel caries appear as a bright spot, which is not reaching the dentin.
- Caries reaching EDJ will appear as bright spot in the enamel dark ring, which connects, to the bright dentin.
- Bright spot in the enamel which is connecting significantly to the bright dentin, may indicate a dentinal/cavitated caries.

### Proximal scoring:

Clinical

Evx

0 Sound

1 Enamel

2 EDJ

3 Dentin/  
Cavitated

Chapter 8:

Patients

On the home screen, press the icon labeled Patients to locate a patient and their scans. Press on a patient to add a new scan or view previous scans.

Patients

1

2

| Chart Number | Chart Number | Date, Time             | Dr.       |
|--------------|--------------|------------------------|-----------|
| Jane Doe     | 416          | 11/10/2014 7:26:45 PM  | Dr. Jones |
| James Doe    | 4221         | 11/10/2014 9:30:11 PM  | Dr. Jones |
| Sarah Doe    | 555          | 11/10/2014 9:26:27 PM  | Dr. Jones |
| John Doe     | 5655         | 11/10/2014 9:14:30 PM  | Dr. Jones |
| Sarah Smith  | 612613       | 11/10/2014 4:49:25 PM  | Dr. Jones |
| Jane Smith   | 6561         | 11/10/2014 4:54:42 PM  | Dr. Jones |
| John Smith   | 56           | 11/24/2014 12:01:20 PM | Dr. Jones |
| James Smith  | 1            | 11/16/2014 2:11:37 PM  | Dr. Jones |
| Jack Smith   | 45           | 11/16/2014 1:52:16 PM  | Dr. Jones |
| Janice Smith | 65           | 11/11/2014 9:00:58 PM  | Dr. Jones |
| Joyce Smith  | 5            | 11/11/2014 8:37:50 PM  | Dr. Jones |

1 Patient page

2 Search bar

**Patient page:** View any patients, see historical scans associated with the patient, access iTero TimeLapse, or add a new scan.

**Search bar:** Search for specific patients using their names or chart numbers.


## Chapter 9:

# iTero TimeLapse

iTero TimeLapse compares previously captured patient's 3D scans to allow a clinician to visualize changes of a patient's teeth, tooth structures, and oral soft tissues.

Patients scanned on a regular basis can be analyzed using the iTero TimeLapse tool.

In the patient's profile, select two scans from the timeline at the bottom of the screen. Once two scans are selected, press the Compare Selected button to analyze the scans.


**Timeline:** Scans taken over time of the patient.

**Scan selection checkboxes:** Selected scans to compare and visualize in iTero TimeLapse.


**Compare selected:** Analyze scans using iTero TimeLapse.

**Color Legend:** Highlighted areas with the most changes between scans with rich color. Scale measurements are in millimeters.


**Loupe:** To view areas in detail and visualize points of interest in the Animation Window. Drag the loupe from the animation window to the area of interest in order to magnify potential treatment spots. The areas in which appear in the loupe will be magnified within the animation window.

**Animation window:** Visualize fine details of changes for the scans selected.

**Color visibility:** Toggle between monochrome and color in the Animation Window.

**Play/pause:** To play or pause the Animation Window.

**Chapter 10:**

## iTero Element 5D Operational manual

The iTero Element 5D system is delivered as a proprietary, PC-based workstation for performing intra-oral scans in the doctor's office. This document describes how to start and shut down the system, how to correctly handle the Scanning Unit/Wand and cable, and how to clean the Scanning Unit and replace its sleeves between patients.

This section is referenced from the iTero Element 5D operational manual 206519

**Chapter 11:**

## iTero Element 5D laptop configuration Operational manual

The iTero Element 5D laptop configuration system is delivered as a proprietary, laptop based workstation for performing intra-oral scans in the doctor's office. This document describes how to start and shut down the system, how to correctly handle the Scanning Unit/Wand and cable, and how to clean the Scanning Unit and replace its sleeves between patients.

This section is referenced from the iTero Element 5D laptop configuration operational manual 206520

Notes:

# iTero element 5D

Align Technology, Inc.  
2820 Orchard Parkway  
San Jose, CA 95134  
USA

Align Technology B.V.  
Arlandaweg 161  
1043HS Amsterdam  
The Netherlands

Align, Invisalign, iTero, iTero Element, and iTero Element 5D among others, are trademarks and/or service marks of Align Technology, Inc. or one of its subsidiaries or affiliated companies and may be registered in the U.S. and/or other countries.

© 2019 Align Technology, Inc. All Rights Reserved. 206721 Rev D

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**Rachel Ovir**  
Director, Quality Assurance & Regulatory Affairs  
Align Technology Ltd.

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