Authentication of Signature

On Behalf of a Body Corporate, Where the
Language of the Document is not Known to the
Person Signing

I, the undersigned, **Itai Izhaky**, Notary holding license
no. 2084632, hereby certify that:

On 20/02/2020 there appeared before me at my office,
located at 124 Eliahu Saadon Blvd. Or Yehuda, **Mrs.
Dvir Rachel**, whose identity has been proven to me by
an Israeli biometric ID Card No. **011035730**, issued on
04/04/2019.

And I was convinced that the person standing before
me fully understood the significance of the action and
voluntarily signed the attached document, marked "A",
on behalf of "Align Technology LTD", Private
Company № 51-202093-4, at the address 3 Ariel Sharon,
Israel; The document is drawn up in the Russian
language, which is not known to her and is also not
known to me.

I confirm that for the purpose of proving the right of
the above to sign on behalf of "Align Technology LTD
" I was provided with
A confirmation by Mrs. Sarah Azulay Adv., issued on
17/02/2020.

And I also certify that the translator Mr. Vladimir
Urbach (I.D. No. 307264036), from the address 5/2
Einstein, Rehovot, Israel, whom I know to be fluent in
the Russian and English languages, translated the said
document in my presence for Mrs. Dvir Rachel into the
English language, which is known to her and to me,
before she signed it.

The Translator Affidavit is attached and marked "1".

In witness whereof I hereby authenticate the signature
of Mrs. Dvir Rachel by my own signature and seal, this
day, 20/02/2020.

The document was signed in the company's offices at 3
Ariel Sharon St. Or Yehuda

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1
Declaration

I the undersigned, **Urbach Vladimir I.D No. 307264036**

After I have been warned to tell the truth and that I will be destined to be punished by Law if I will not, hereby declare in writing as follows:

1. I hereby declare that I am well acquainted with the Russian, And the English languages and that the attached document marked "B" is a correct English translation from the Russian language original document which its copy is attached hereto marked "A".

2. This is my name and this is my signature and all that is written above is true and correct.

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Rachel Dylr

Director, Quality Assurance & Regulatory Affairs

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The computer is provided with a Li-ion rechargeable battery pack. There is a danger of explosion if battery is incorrectly replaced. Replace only with same type recommended by the manufacturer. Discard used batteries according to the manufacturer's instructions.

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Power supply

- Power is supplied to the system via an internal medical grade power supply.

Battery power

- Charging - Battery will be fully charged after being plugged into a power source for 2 hours.
- With a fully charged battery, the user can scan up to 30 minutes with the iTero Element 2 scanner, without having to plug in for power.

Electric warning

- Electric shock hazard!!** Only authorized Align Technology technicians can remove external panels and covers. There are no user-serviceable parts inside.
- To avoid risk of electric shock, iTero Element 2 must only be connected to a supply mains with protective grounding.
- Only Align Technology approved Web Camera or DOK should be connected to the USB socket on the back side of the system.

Wireless LAN

- The system comes equipped with a Wireless LAN unit.

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Scanner warnings

- The scanner emits red laser light (680nm Class 1) as well as white LED emissions. Normal usage of the scanner does not present any danger to the human eye. However, doctors should refrain from shining the scanner directly into the patient's eyes.
- Avoid twisting cable, knotting cable, pulling on cable, stepping on cable.

- Usage of the scanner does not present any danger to the human eye. However, doctors should refrain from shining the scanner directly into the patient's eyes.
- When the system is not in use, the scanning unit should be placed inside the cradle with the probe facing towards the cart's post and the rear side of the touch screen so there will be no direct eye contact with the laser beam or the flickering white LED emission in any case.
- The doctor should activate scanning operation only while the scanner's probe is inside the patient's mouth.
- Doctors should avoid placing the scanner in the cradle while scanning operation is still active.

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Work environment

- The system should be moved between rooms with utmost care to avoid damage.
- Do not block the air vents on the scanning unit and base unit.
- System is intended for indoor use only. It should not be exposed to direct sunlight, excessive heat or humidity.

Electro magnetic interference

- WARNING:** This device has been tested and found to comply with the requirements for medical devices according to standard EN60601-1-2. This standard is designed to provide reasonable protection against harmful interference in a typical medical installation. However, because of the proliferation of radio-frequency transmitting equipment and other sources of electrical noise in the healthcare environments (e.g., cellular phones, mobile two-way radios, electrical appliances), it is possible that high levels of such interference due to close proximity or strength of source, may result in disruption of performance of this device.

General

- WARNING:** No modification of this equipment is allowed.
- WARNING:** The touch screen always needs to be in a stand while in operation!

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
Table of contents

Chapter 1: Introduction	2
About this operation manual	2
Intended use	2
Benefits of the iTero Element 2 system	2
The iTero Element 2 user interface	3
Chapter 2: Basic hardware features	4
Front view of the system	4
Back view of the system	5
Scanning unit (wand)	6
Chapter 3: Assembly instructions	7
Step 1: Assembly	7
Step 2: Make it Mine process	10
Chapter 4: Operating instructions	11
Initial power-up and boot	11
End-of-day shut down	11
Moving system within the office	11
Chapter 5: Scanner handling, cleaning, and disinfection instructions	12
Handling of the scanning unit (wand)	12
Handling of the scanning unit cable	12
Recommended best practices for cleaning and disinfecting	12
Cleaning and disinfectant materials for scanning unit and base unit	13
Chapter 6: Changing sleeves between patients	14
Cleaning and disinfecting the scanning unit (wand)	14
Replacing disposable sleeves	14
Scanner sleeves	15
Appendix A: EMC declaration	16
Appendix B: Hardware specifications	17

Chapter 1:

Introduction

About this operation manual

The iTero Element 2 system is delivered as a proprietary, PC-based workstation for performing intra oral scans in the doctor's office. This operation manual describes how to boot and shut down the system, how to correctly handle the scanning unit/wand and cable, and how to clean the scanning unit and replace its sleeves between patients.

Intended use

iTero Element 2 is an optical impression system (CAD/CAM) used to record the topographical images of teeth and oral tissue. Data generated from iTero may be used in conjunction with the production of dental devices (e.g., aligners, braces, appliances, etc.) and accessories.

iTero Element 2 software is used with the iTero Element 2 scanner in capturing 3D digital impressions of teeth, oral soft tissue and structures, and bite relationship. The software controls the processing of the data, facilitating the integration of data, and exporting of the data for CAD/CAM fabrication of dental restorations, Orthodontic devices, abutments, and accessories. In addition to scan data, various patient and case information can be imported/exported or used for simulation purposes. Other functions are available for verification and service of the system and to serve as an order management tool.

Benefits of the iTero Element 2 system

The iTero Element 2 system provides important advantages over existing crown-production methods, including powder-free scanning, greater crown-production accuracy, and immediate feedback during the scanning process.

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
The iTero Element 2 user interface

The iTero Element 2 system provides an intuitive user interface for performing digital scans for Restorative or Orthodontic use. The doctor is guided through the scanning sequence by means of visual and text assistance. The touch screen and wand buttons are used to respond to screen instructions during the scanning process.

 One tap on the question mark will enable a transparent Help overlay that will provide a brief overview. Please note that the Headset image appears instead of the question mark while in this view. Tap anywhere to close the help screen and return to the relevant screen.


Chapter 2:

Basic hardware features


Custom wheel stand hardware features: Front view of the system

Custom wheel stand hardware features: Back view of the system

Scanning unit (wand)

Chapter 3:

Assembly instructions

Step 1: Assembly


Please follow the instructions below to assemble your iTero Element 2 scanner:

- A** Wheel stand
- B** Wand with cable
- C** Wand cradle
- D** HD touch screen
- E** External battery
- F** Power cable

1. Check the content of the box

2. Connect post to the wheel base

3. Tighten the 2 allen screws using the larger wrench

4. Remove cover from the back of the handle

5. Attach the wand cradle to the front of the wheel stand

6. Hold the cradle

7. Tighten the wand cradle allen screw on the back using the smaller wrench

8. Reattach the cover behind the handle


9.
Remove the magnetic cover from the back of the wheel stand frame

10.
Remove the battery cover

11.
Slide the battery into the battery slot and tighten the thumb screws

12.
Lift the HD screen to mount it

13.
Turn the scanner around and tighten the thumb screw to secure the HD screen

14.
Attach the power cable to the port labeled DC


15.
Power cable inserted

16.
Attach the magnetic back cover


17.
Place the wand into the cradle

18.
Attach the wand cable on the back of the HD screen

19.
Attach the power cable on the bottom of the wheel stand

20.
On the bottom of the wheel stand, post and secure the cable with the clip

21.
Plug in the power cable


22.
Position the web camera on the HD screen for remote training or support sessions


23.
Plug in the webcam to the USB port at the bottom of the HD screen

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
2. Follow the Wizard instructions on the screen to complete the customization of the iTero Element 2.

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End-of-day shut down

1. Close all files and applications.
2. Press and release the power switch on the bottom of the screen to shut down the system.

Moving system within the office

To ensure maximum system protection, it is recommended to have two people move the system. Follow these instructions for relocating the system:

1. Verify the scanning unit (wand) sits well inside the scanning unit cradle.
2. Unplug system from the wall outlet.
3. Move the system carefully using two people.
4. Place the system at its new location and it plug into a wall outlet.

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- The scanner cable should be treated with care to avoid possible damage.

Between patient sessions, it is recommended to undo any twists and knots in order to relieve all tension from the scanner cable

Recommended best practices for cleaning and disinfecting the scanning unit, base unit, wheel stand and/or counter stand in between patients.

- Do not spray disinfectant directly on scanner system surfaces.
- Spray the disinfectant on a towel, or use disinfectant wipes for the scanning unit, and base unit.
- **Warning:** over saturation of disinfectant product on the scanner system surfaces may cause damage, including internal components.
- Follow the disinfectant manufacturers' instructions for appropriate contact time. Remove residual liquid disinfectant with a lint-free, clean cloth.
- **Note:** follow standard precautions for personal protection, as appropriate.
- **Warning:** DO NOT touch the optical surface of the scanning unit (wand).

Cleaning and disinfectant materials for scanning unit and base unit

The following cleaning and disinfectant materials are recommended for use for the scanning unit and the base unit.

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Appendix A:

EMC declaration

Summary of EMC test results for iTerо Element 2

Test	Standard	Class / Severity level	Test result
Documentation (IEC 60601-1-2 sections 4 and 5)			
General requirements for EMC	section 4.1	--	Complies
External labels	section 5.1	--	Complies
Conformity of users' manual	section 5.2.1	--	Complies
Accuracy of technical description	section 5.2.2	--	Complies
Emission (IEC 60601-1-2 / EN 60601-1-2 section 7)			
Conducted emission Freq. range: 150 kHz - 30 MHz	CISPR 11	Group 1 Class B 230, 120 & 100 VAC mains (50 Hz); 220 VAC mains (50 & 60 Hz)	Complies
Radiated emission Freq. range: 30 - 1000 MHz	CISPR 11	Group 1 Class B	Complies
Harmonic current emission test	IEC 61000-3-2	230 VAC mains (50 Hz & 220 V (60 Hz)	Complies
Voltage changes, voltage fluctuations and flicker test	IEC 61000-3-3	230 VAC mains & 220 VAC mains	Complies
Immunity (IEC 60601-1-2 section 8)			
Immunity from electrostatic discharge (ESD)	IEC 61000-4-2	8 kV contact discharges & 15 kV air discharges	Complies
Immunity from radiated electromagnetic fields	IEC 61000-4-3	10.0 V/m; 80 MHz ÷ 2.7 GHz, 80% AM, 1 kHz	Complies
Immunity from proximity field from wireless communications equipment	IEC 61000-4-3	List of frequencies, from 9 V/m up to 28 V/m, PM (18 Hz or 217 Hz), FM 1 kHz	Complies
Immunity from electrical fast transient (EFT)	IEC 61000-4-4	± 2.0 kV - on AC mains; Tr/Th - 5/50 ns, 100 kHz	Complies
Immunity from surge	IEC 61000-4-5	±2.0 CM / ±1.0 kV DM on AC mains; Tr/Th - 1.2/50 (8/20) µs	Complies
Immunity from conducted disturbances induced by radio-frequency fields	IEC 61000-4-6	3.0, 6.0 VRMS on 230 VAC mains & wand cable; 0.15÷ 80 MHz, 80% AM 1 kHz	Complies
Immunity from voltage dips, short interruptions and voltage variations	IEC 61000-4-11	0 % - 0.5 cycle & 1 cycle; 70% - 25 cycles; 0% - 250 cycles on AC mains	Complies

Appendix B:

Hardware specifications

Item	Description
Monitor	21.5" monitor
Scanner	Scanner emits red laser light (680nm Class 1) as well as white LED emissions.
Wireless LAN	LAN card provides local network communications with wireless connectivity.
Operating power	100-240VAC- 50/60 Hz – 200VA (max)
Operating temperature	18° to 26°C / 64.4° to 78.8°F
Storage/transportation temperature	-5° to 50°C / 23° to 122°F
Operating pressure & altitude	Pressure: 520 mmHg to 760 mmHg (-69 kPa to -101 kPa) Altitude: 0 feet to 10,000 feet
Storage/transportation pressure & altitude	Pressure: 430 mmHg to 760 mmHg (-57 kPa to -101 kPa) Altitude: 0 feet to 15,000 feet
Relative humidity	Operating: 40% to 70%; Storage: 30% to 90%
Dimensions	iTero HD touch monitor: Height: 356 mm (-14 in) Width: 552 mm (-21.7 in) Depth: 65 mm (-2.5 in) iTero Element wand: Length: 338.5 mm (-13 in) Width: 53.5 mm (-2 in) Depth: 69.8 mm (-3 in) Wheel stand: Height: 1280 mm (-50 in) Width: 645 mm (-25 in) Depth: 625 mm (-24.5 in)
Net weight	iTero HD touch monitor: 8.3 kg (-18.3 lbs) Tero Element wand: 0.47 kg (-1 lbs) Wheel stand: 13.6 kg (-30 lbs)
Shipping weight	-37.5 kg (-83 lbs)

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© 2018 Align Technology, Inc. All Rights Reserved. 202937 Rev D

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Operation manual

Copyright

© 2018 Align Technology, Inc. All Rights Reserved.

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The hardware and software described in this manual are supplied under a purchase agreement or purchase order, including the terms and conditions set out or referenced therein, pursuant to which the Customer has purchased the iTero Element Flex Optical Impression Device and may be used only in accordance with such terms.

No part of this manual may be reproduced, photocopied, stored in a retrieval system, or transmitted in any manner (electronic or mechanical) for any purpose other than the customer's normal usage, without the prior written permission of Align Technology Inc.

This manual describes the iTero Element Flex Optical Impression Device.

English language version. Updated January 2018.

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E-mail:
erosupport@aligntech.com

E-mail:
orthocadsupport@aligntech.com

European authorized representative

Obelis s.a.
Boulevard Général Wahis 53
B - 1030 Brussels, Belgium

www.obelis.net

Tel: +32.2.732.59.54
Fax: +32.2.732.60.03
E-mail: mail@obelis.net

Class 1 laser compliance

This device complies with: "21 CFR 1040.10" and "EN 60825-1".

Class 1 Laser Product

CSA compliance

This device complies with the following CSA standard for Canada and the USA:

"UL Std No. 60601-1 – Medical Electrical Equipment Part 1: General Requirements for Safety"

**FCC compliance**

This device complies with Part 15 of FCC Rules and its operation is subject to the following two conditions:

1. This device may not cause harmful interference.
2. This device must accept any interference received, including interference that may cause undesired operation.

**FCC warning**

Modifications to the device that are not expressly approved by the manufacturer may void your authority to operate the device under FCC Rules.

EMC compliance

This device complies with the following EMC standard:

"IEC 60601-1-2 Medical electrical equipment – Part 1-2: General requirements for basic safety and essential performance – Collateral standard: Electromagnetic phenomena – Requirements and tests".

Safety compliance

This device complies with the following safety standard:

"IEC 60601-1 Medical electrical equipment –Part 1: General requirements for basic safety and essential performance."

CE compliance

This device complies with Council Directive 93/42/EEC for Medical Devices.

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Symbols

The following symbols may appear on iTero Element Flex hardware components, and may also appear within this manual and other iTero Element Flex literature.

Wherever this symbol appears on the device, it is recommended to refer to this manual for information on proper usage of the device.

Applied part type BF. Any component on which this symbol appears is electric isolation type BF.

Parts or accessories on which this symbol occurs should not be reused.

Attention: This symbol is used to highlight the fact that there are specific warnings or precautions associated with the device. Wherever this symbol appears on the device, it is mandatory to refer to safety-related information within this manual.

"Rx only"

CAUTION: US Federal Law restricts this device to sale by or on the order of a licensed Dentist, Orthodontist or Dental Professional. The system serves as a prescription medical device and should be operated by qualified health-care providers only.

Separate collection of electrical waste and electronic equipment is required. In compliance with the European Directive on Waste Electrical and Electronic Equipment (WEEE), do not dispose this product in domestic or municipal waste. This device contains WEEE materials. Please contact EARN service. The link for the online request form <http://b2btool.earn-service.com/aligntech/select>.

Medical device manufacturer.

CAUTION: DO NOT step on the Flex hub!

CAUTION: Flex hub must be kept dry!

Order number.

Serial number.

Indicates the Authorized representative in the European Community.

Indicates the temperature limits to which the medical device can be safely exposed

Manufacturer's batch code

Indicates the need for the user to consult the instructions for use.

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Safety instructions

Before beginning to work with the system, all users are required to read these safety instructions.

The system is provided with a Flex hub containing a power supply for the scanning unit. The Flex hub must be kept dry and be guarded against breakage.

WARNING – to avoid risk of electric shock, this equipment must only be connected to a supply mains with protective earth.

Power supply

- Power to the scanning unit is supplied via an internal medical grade power supply inside the Flex hub.

Electric warning

- **Electric shock hazard!!** There are no user-serviceable parts inside the Flex hub.
- To avoid risk of electric shock, iTero Element Flex must only be connected to a supply mains with protective grounding.
- Connect the Flex hub device only to a laptop that is approved according to IEC60950, and to UL60950-1. Laptop and all its accessories shall be located at least 1.5m away from the patient. Do not scan a patient and touch the laptop or any of its accessories at the same time.
- Only the Align Technology scanning unit and the approved laptop should be connected to the USB sockets on the Flex hub.
- Only the Align Technology approved power cable should be used to connect the Flex hub with the AC outlet.

Safety classifications

- Type of protection against electrical shock: Class 1.
- Degree of protection against electrical shock: Type BF.
- Degree of protection against harmful ingress of water: Ordinary.
- Equipment not suitable for use in presence of flammable anesthetic mixtures.
- Mode of operation: Continuous.

Prescription health device

- The system serves as a prescription medical device and should be operated by qualified health-care providers only.

Scanner warnings

- The scanner emits red laser light (680nm Class 1) as well as white LED emissions. Normal usage of the scanner does not present any danger to the human eye. However, doctors should refrain from shining the scanner directly into the patient's eyes.
- Avoid twisting cable, knotting cable, pulling on cable, stepping on cable.
- When the system is not in use, the scanning unit should be placed inside the cradle with the probe facing down so there will be no eye contact with the laser beam in any case.

- Doctors should activate the scanning operation only while scanner's probe is inside the patient's mouth.
- Doctors should avoid placing the scanner in the cradle while scanning operation is still active.

Cleaning & disinfection

- To avoid cross contamination, it is mandatory that after each patient session the disposable plastic sleeve be replaced and the scanning unit be disinfected.
- Dispose of scanner sleeves according to standard operating procedures or local regulations for the disposal of contaminated medical waste.

Unpacking & installing

- The system should be unpacked and installed following Align Technology's instructions.

Work environment

- The system should be moved between offices inside its designated case to avoid damage.
- Do not block the air vents on the scanning unit.
- System is intended for indoor use only. It should not be exposed to direct sunlight, excessive heat or humidity.
- The system should only be used following operating temperatures as defined in Appendix B.
- If iTero Element Flex has just been brought into the office from a hot or cold or humid environment, it should be set aside until it has adjusted to the room temperature to avoid internal condensation.

Electro magnetic interference

- **WARNING:** This device has been tested and found to comply with the requirements for medical devices according to standard IEC60601-1-2. This standard is designed to provide reasonable protection against harmful interference in a typical medical installation. However, because of the proliferation of radio-frequency transmitting equipment and other sources of electrical noise in the healthcare environments (e.g., cellular phones, mobile two-way radios, electrical appliances), it is possible that high levels of such interference due to close proximity or strength of source, may result in disruption of performance of this device.
- In the event of a sharp drop in the input mains voltage, the system will not operate but will remain safe for the user. The system will return to its working state after the voltage returns to its nominal values.

General

- **WARNING:** No modification of this equipment is allowed.

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
Table of contents

Chapter 1: Introduction	2
About this operation manual	2
Intended use	2
Benefits of the iTero Element Flex system	2
The iTero Element Flex user interface	3
Chapter 2: Basic hardware features	4
Front view of the system	4
Scanning unit (wand)	5
Mobile system packaging (for transport)	6
Moving the system	6
Chapter 3: Assembly and software installation	7
Step 1: Assembly	7
Step 2: Download iTero Element Flex software	8
Step 3: Make it Mine process	9
Chapter 4: Operating instructions	10
Initial power-up and boot	10
Web cam	10
Chapter 5: Scanner handling, cleaning, and disinfection instructions	11
Handling of the scanning unit (wand)	11
Handling of the scanning unit cable	11
Recommended best practices for cleaning and disinfecting	11
Cleaning and disinfectant materials for scanning unit and base unit	12
Chapter 6: Changing sleeves between patients	13
Cleaning and disinfecting the scanning unit (wand)	13
Replacing disposable sleeves	13
Scanner sleeves	14
Appendix A: EMC declaration	15
Appendix B: Hardware specifications	16

Chapter 1:

Introduction

About this operation manual

The iTero Element Flex system is delivered as a proprietary, laptop-based workstation for performing intra oral scans in the doctor's office. This operation manual describes how to boot and shut down the system, how to correctly handle the scanning unit/wand and cable, and how to clean the scanning unit and replace its sleeves between patients.

Intended use

iTero Element Flex is a mobile optical impression system (CAD/CAM) used to record the topographical images of teeth and oral tissue. The system can be connected to laptops approved by Align Technology to reduce office space requirements and increase system mobility.

Data generated from iTero may be used in conjunction with the production of dental devices (e.g., aligners, braces, appliances, etc.) and accessories.

iTero Element Flex software is used with the iTero Element scanner in capturing 3D digital impressions of teeth, oral soft tissue and structures, and bite relationship. The software controls the processing of the data, facilitating the integration of data, and exporting of the data for CAD/CAM fabrication of dental restorations, Orthodontic devices, abutments, and accessories. In addition to scanning data, various patient and case information can be imported/exported or used for simulation purposes. Other functions are available for verification and service of the system and to serve as an order management tool.

Benefits of the iTero Element Flex system

The iTero Element Flex system provides important advantages over existing crown-production methods, including powder-free scanning, greater crown-production accuracy, and immediate feedback during the scanning process.

Refer to our website <http://www.itero.com> to learn how the iTero Service can enhance your business by increasing patient satisfaction, improving clinical outcomes, and enhancing office efficiency.

Software installation instructions

Please consult the list of approved laptops and approved Antivirus applications at www.iTero.com.


To ensure optimal performance use of the iTero application, only approved laptops must be used with the iTero Element Flex system.

The iTero Element Flex user interface

The iTero Element Flex system provides an intuitive user interface for performing digital scans for Restorative or Orthodontic use. The doctor is guided through the scanning sequence by means of visual and text assistance. The laptop's touch screen and the wand buttons are used to respond to screen instructions during the scanning process.

One tap on the question mark will enable a transparent Help overlay that will provide a brief overview. Please note that the Headset image appears instead of the question mark while in this view. Tap anywhere to close the help screen and return to the relevant screen.

Chapter 2:


Basic hardware features

Front view of the system


Laptop with touch screen

Flex hub

Wand and cradle

Scanning unit (wand)

Mobile system packaging (for transport)

Transportation

Use the supplied case to move the system between offices.

Moving the system

To ensure maximum system protection, it is recommended to move the system with care. Follow these instructions for relocating the system:

1. Attach the blue protective sleeve onto the wand.
2. Verify the scanning unit (wand) sits well inside the case.
3. Unplug the Flex hub and laptop AC power cables from the wall outlet.
4. Place all items in the case (see picture above).
5. Make sure the case is kept dry to protect the system components from humidity.

Chapter 3:

Assembly and software installation

Step 1: Assembly

Please follow the instructions below to assemble your iTero Element Flex scanner:

- A** Hub and hub power cord
- B** Cradle
- C** Wand and wand cable
- D** USB cable to connect laptop and hub

1. Place the wand in the cradle

2. Connect hub power cord to the hub

3. Connect USB cable to the hub

4. Connect USB cable to the laptop

5. Connect wand cable to the hub

6. Plug the hub power cord into the AC power outlet

Step 2: Download iTero Element Flex software

For proper software installation and configuration of your iTero Element Flex scanner, please ensure your iTero wand is secure in the cradle and connected to the hub, and the hub is connected to the laptop. Please make sure your laptop is plugged into the AC wall socket during the entire software installation.

1. Laptop setup

- Turn on the laptop.
- Update the Windows to the latest version.

2. iTero software download

- Open the registered email inbox in the browser. Look for the email "Your iTero was shipped" which includes instructions.
- Click on the link to get to the page with software downloads.
or
- Alternatively, open the browser and go to (type the address into your browser) download.itero.com.

3. iTero software installation


- Click on the "Get started" button.
- Run the downloaded installation file.
- Follow the instructions on the screen to complete the download of the iTero software.

It is highly recommended to regularly check for Windows updates and install them. Please enable updates.


It is also recommended to exit to Windows and shut down the computer at the end of each day to allow updates to install.

Before installing the iTero Element Flex software for first time, please install all available Windows updates. New Windows computers should apply updates automatically.


To check for Windows Updates, open "Settings" (Windows key + I), choose "Update & Security", and then "Windows Update". Click "Check for updates" to see if there are new updates available.

Step 3: Make it Mine process

1. Select language of preference and tap on the Make it Mine button to start the Wizard.

2. Follow the Wizard instructions on the screen to complete the customization and installation of the iTero Element Flex.

Chapter 4:

Operating instructions

Make sure the laptop is connected to an AC wall outlet, to ensure the battery level does not reach its low battery status, at which time scanning will no longer be possible.

Power up the laptop to start iTero program.

It is recommended to keep the system in operation during office hours to allow background file transfers between the doctor's office, the doctor's partnered labs, and the Align Technology Center. It is recommended to shut down the system at the end of the day, and to reboot in the morning.

Web cam

 Connect the web cam to the laptop following the product instructions.

Chapter 5:

Scanner handling, cleaning, and disinfection instructions

Handling of the scanning unit (wand)

- The scanning unit contains delicate components and should be handled with care.

Handling of the scanning unit cable

- The scanner cable should be treated with care to avoid possible damage.
- Avoid twisting, knotting, pulling or stepping on the cable, etc.
- Between patient sessions, it is recommended to undo any twists and knots in order to relieve all tension from the scanner cable.

Recommended best practices for cleaning and disinfecting the scanning unit in between patients.

- Do not spray disinfectant directly on scanner system surfaces.
- Spray the disinfectant on a towel, or use disinfectant wipes for the scanning unit.
- **Warning:** over saturation of disinfectant product on the scanner system surfaces may cause damage, including internal components.
- Follow the disinfectant manufacturers' instructions for appropriate contact time. Remove residual liquid disinfectant with a lint-free, clean cloth.
- **Note:** follow standard precautions for personal protection, as appropriate.
- **Warning:** DO NOT touch the optical surface of the scanning unit (wand).

Cleaning and disinfectant materials for scanning unit and cradle

The following cleaning and disinfectant materials are recommended for use for the scanning unit and the cradle.

Description	pH	Manufacturer P/N	Manufacturer
Birex® Quat disinfectant wipes	7.6	BI 240	Biotrol Intl.
CaviCide AF	12.7	13-800	Metrex
CaviCide CaviWipe	12.5	13-1000 13-1100	Metrex
CaviCide 1 CaviWipe 1	12.5	13-5000 13-5100	Metrex
Clorox Healthcare® hydrogen peroxide cleaner disinfecting liquid	2-3	30828, 30829	Clorox® Healthcare™
Clorox Healthcare® hydrogen peroxide cleaner disinfectant wipes	2-3	30824, 30825	Clorox® Healthcare™
Opti-Cide 3® liquid	7.6	DOCS12-024, DOCP04-128	Biotrol Intl.
Opti-Cide 3® wipes	7.6	DOCW06-100	Biotrol Intl.
OPTIM 33TB liquid	2.5-3.5	OPT33-1GAL, OPT-33-1QT	SciCan Inc.
OPTIM 33TB wipes	2.5-3.5	OPT33-W10X10, OPT-33-W12	SciCan Inc.
ProSpray ProSpray wipes	10	PSC240 PSW-1	Certol
Webcol® alcohol prep pads	7	5110	Medtronic

Chapter 6:

Changing sleeves between patients

Cleaning and disinfecting the scanning unit (wand)

To avoid cross contamination, it is essential that after each patient you fully clean and disinfect the scanning unit and the disposable sleeve. First spray disinfectant material on towel or use disinfectant wipes to clean the scanning unit and scanning unit cradle. Then proceed with the steps below to remove the used sleeve and attach a new disposable sleeve.

CAUTION: Dispose of scanner sleeves according to standard operating procedures or local regulations for the disposal of contaminated medical waste.

Replacing disposable sleeves

Step 1

When pulling a sleeve OFF or ON, hold the center of the sleeve.

Step 2

Press slightly on both sides of the disposable sleeve, pull the sleeve slowly off the scanning unit and discard.

Step 3

Gently slide on new sleeve onto scanning unit until it clicks into place.

**WARNING: Optical surface!**

DO NOT touch the optical surface. Contact may cause damage. If cleaning is necessary, use the wipes and anti-static cloth found inside the sleeves box. For proper use, refer to the directions found in the scanner sleeves box.

Scanner sleeves

There are two types of sleeves intended for use with the scanner unit (wand):

Disposable sleeve

The white sleeve is a single use sleeve for patient scanning. Always replace the white sleeve on the scanning unit between patients to avoid cross contamination. Please dispose of the white sleeve after every patient.

Protective sleeve

The blue protective sleeve is used to protect the optical surface lens when the scanning unit is not in use. Please keep the blue sleeve in a safe place so that it does not get lost or damaged.

Scanner sleeves packaging box

Scanner sleeves may be ordered online in boxes of 25 from the iTero store www.store.itero.com, where available.

Appendix A:

EMC declaration


Summary of EMC test results for iTero Element Flex

Test	Standard	Class / Severity level	Test result
Documentation (IEC 60601-1-2 sections 4 and 5)			
General requirements for EMC	section 4.1	--	Complies
External labels	section 5.1	--	Complies
Conformity of users' manual	section 5.2.1	--	Complies
Accuracy of technical description	section 5.2.2	--	Complies
Emission (IEC 60601-1-2 section 7)			
Conducted emission Freq. range: 150 kHz - 30 MHz	CISPR 11	Group 1 Class B 230, 120 & 100 VAC mains @ 50 Hz; 220 VAC mains @ 60 Hz	Complies
Radiated emission Freq. range: 30 - 1000 MHz	CISPR 11	Group 1 Class B	Complies
Harmonic current emission test	IEC 61000-3-2	230 VAC mains @ 50 Hz & 220 V @ 60 Hz	Complies
Voltage changes, voltage fluctuations and flicker test	IEC 61000-3-3	230 VAC mains @ 50 Hz & 220 V @ 60 Hz	Complies
Immunity (IEC 60601-1-2 section 8)			
Immunity from electrostatic discharge (ESD)	IEC 61000-4-2	8 kV contact discharges & 15 kV air discharges	Complies
Immunity from radiated electromagnetic fields	IEC 61000-4-3	10.0 V/m; 80 MHz , 2.7 GHz, 80% AM, 1 kHz	Complies
Immunity from proximity field from wireless communications equipment	IEC 61000-4-3	List of frequencies, from 9 V/m up to 28 V/m, PM (18 Hz or 217 Hz), FM 1 kHz	Complies
Immunity from electrical fast transient (EFT)	IEC 61000-4-4	± 2.0 kV - on AC mains; Tr/Th - 5/50 ns, 100 kHz	Complies
Immunity from surge	IEC 61000-4-5	± 2.0 CM / ± 1.0 kV DM on 230 VAC mains @ 50 Hz & 220 VAC mains @ 60 Hz; Tr/Th - 1.2/50 (8/20) ms	Complies
Immunity from conducted disturbances induced by radio-frequency fields	IEC 61000-4-6	3.0, 6.0 VRMS on 230 VAC mains & wand cable; 0.15 - 80 MHz, 80% AM @ 1 kHz	Complies
Immunity from voltage dips, short interruptions and voltage variations	IEC 61000-4-11	On 230 VAC & 100 VAC mains @ 50 Hz: 0% - 0.5 cycle & 1 cycle; 70% - 25 cycles; 0% - 250 cycles on 220 VAC mains @ 60 Hz: 0% - 0.5 cycle & 1 cycle; 70% - 30 cycles; 0% - 300 cycles	Complies

Appendix B:

Hardware specifications

Item	Description
Monitor	Laptop monitor
Scanner	Scanner emits red laser light (680nm Class 1) as well as white LED emissions.
Operating power	100-240VAC- 50/60 Hz – 40VA (max)
Operating temperature	18° to 26°C / 64.4° to 78.8°F
Storage/Transportation temperature	-5° to 50°C 23° to 122°F
Operating pressure & altitude	Pressure: 520 mmHg to 760 mmHg (-69 kPa to -101 kPa) Altitude: 0 feet to 10,000 feet
Storage/Transportation pressure & altitude	Pressure: 430 mmHg to 760 mmHg (-57 kPa to -101 kPa) Altitude: 0 feet to 15,000 feet
Relative humidity	Operating: 40% to 70%; Storage: 30% to 90%
Dimensions	iTero Element Flex hub: Length: 206 mm (-8 in) Width: 94 mm (-3.7 in) Depth: 36.5 mm (-1.4 in) iTero Element wand: Length: 338.5 mm (-13 in) Width: 53.5 mm (-2 in) Depth: 69.8 mm (-3 in) Cradle: Length: 262 mm (-10 in) Width: 89 mm (-3.5 in) Depth: 52 mm (-2 in) Case: Height: 326.5 mm (-13 in) Width: 455 mm (-18 in) Depth: 184 mm (-7 in)
Net weight	iTero Element Flex hub: ~0.5 kg (~1 lbs) iTero Element wand: 0.47 kg (~1 lbs) Empty Case: ~2 kg (~4.5 lbs)
Shipping weight	~8 kg (~17.6 lbs)

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© 2018 Align Technology, Inc. All Rights Reserved. 204164 Rev D

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Align Technology, Inc.
2820 Orchard Parkway
San Jose, CA 95134
USA

Align Technology B.V.
Ariandaweg 161
1043HS Amsterdam
The Netherlands

Align, Invisalign, iTero, and iTero Element among others, are trademarks and/or service marks of Align Technology, Inc. or one of its subsidiaries or affiliated companies and may be registered in the U.S. and/or other countries.

© 2019 Align Technology, Inc. All Rights Reserved. 203410 Rev C

Copyright

© 2019 Align Technology, Inc. All Rights Reserved

The information contained in this manual is subject to change without notice.

The hardware and software described in this manual are supplied under a Sales and Services Agreement and may be used only in accordance with the terms of that agreement.

No part of this manual may be reproduced, photocopied, stored in a retrieval system, or transmitted in any manner (electronic or mechanical) for any purpose other than the customer's normal usage, without the prior written permission of Align Technology Ltd.

English language version. Updated in April 2019.

Trademarks

Align, Invisalign, iTero, and iTero Element among others, are trademarks and/or service marks of Align Technology, Inc. or one of its subsidiaries or affiliated companies and may be registered in the U.S. and/or other countries.

Any other trademarks or registered trademarks appearing in this manual are the property of their respective owners.

USA office

Corporate Headquarters
Align Technology, Inc.
2820 Orchard Parkway
San Jose, California 95134

www.aligntech.com

Tel: +1 (408) 470-1000
Fax: +1 (408) 470-1010

Align Technology Ltd.
3 Ariel Sharon Boulevard
Or-Yehuda 6037606
Israel

Tel: +972 (3) 634-1441
Fax: +972 (3) 634-1440

Netherlands office

International Headquarters
Align Technology B.V.
Arlandaweg 161
1043 HS Amsterdam

Tel: +31 (0) 20-586-3600
Fax: +31 (0) 20-586-3751

Customer support

Tel: +1 (800) 577-8767

E-mail:
iterosupport@aligntech.com

E-mail:
orthocadsupport@aligntech.com

Align Technology B.V.
Arlandaweg 161
1043HS Amsterdam
The Netherlands

Table of contents

Introduction	2
Intended use	2
1.0 System registration	3
Demo Mode	3
Make it Mine	4
Connecting to hidden networks	4
2.0 Logging in to the scanner application	6
Login	6
Create account	6
Forgot password	6
3.0 iTero Element main features	7
Home screen	7
Main desktop icons	7
Utility icons	8
Touch screen gestures	9
Wand controls	10
Demo Mode	11
4.0 Messages	12
5.0 Orders	13
6.0 Scan	14
7.0 Patients	19
8.0 iTero TimeLapse	20

Introduction

The iTero Element intraoral scanner is Align Technology's latest generation integrated scanner and software system. It enables doctors to take 3D digital scans of patients' teeth and bite, make adjustments of the scan in real time and transmit the results wirelessly.

This user guide provides general information and an overview of the application.

Intended use

The iTero Element is an optical impression system (CAD/CAM) used to record the topographical images of teeth and oral tissue. Data generated from iTero may be used in conjunction with the production of dental devices (e.g., aligners, braces, appliances, etc.) and accessories.

iTero Element software is used with the iTero Element scanner capturing 3D digital impressions of teeth, oral soft tissue and structures, and bite relationship. The software controls the processing of the data, facilitating the integration of data, and exporting of the data for CAD/CAM fabrication of dental restorations, orthodontic devices, abutments, and accessories. In addition to scan data, various patient and case information can be imported/ exported or used for simulation purposes. Other functions are available for verification and service of the system and to serve as an order management tool.

The iTero Element is an optical impression system built using computer-aided design and manufacturing (CAD/CAM) methods. The system is used to scan and record the topographical images of teeth and oral tissue. Upon the completion of a scan, data generated from the iTero Element can be sent wirelessly for the production of dental devices and accessories.

This all-in-one system is available on a monitor with a fully interactive, touch screen display and an easy-to-use scanner. The topography of a patient's teeth can be viewed on the screen as they are being scanned, and determining the degree of occlusion of the bite can be analyzed as the scan is completed.


In conjunction with the scanner, the application facilitates the following tasks by:

- Capturing 3D digital impressions of teeth, oral soft tissue and structures as they are being scanned.
- Enabling real-time viewing of the visual bite relationship once the scan is completed.
- Managing the processing of the data.
- Sending the data wirelessly to designated providers for the fabrication of dental restorations, orthodontic devices, abutments, and accessories, etc.
- Using the data for simulation purposes.

The application also includes administrative capabilities which enable the user to:

- Start a new scan.
- Place a new order for an existing patient or new patient.
- View orders in progress.
- Review and/or track past orders.

1.0 System registration

1 Demo Mode

2 Make it Mine

Demo Mode

Press Demo Mode to familiarize yourself with the scanner's features and perform practice scans without submitting the scans. Demo Mode is designed for training new staff members and practice. Demo Mode is available anytime for dental practices to train on an iTero Element device, for scanning techniques, how-to guidelines for prescription forms, case types, and to familiarize with the iTero interface. Demo Mode will feature all the aspects of the scanning process and include a wide variety of cases. Sample cases offered will include clinical cases, Invisalign cases, and restorative cases for your convenience.

When Demo Mode is in use for practice scanning, a lightly striped background and red tag in the upper left-hand corner will indicate Demo Mode. Any scans captured in Demo Mode cannot be saved or submitted for patient treatment. Demo Mode is available to access on the Welcome screen or at any point by pressing the iTero logo button on the home screen.

Make it Mine


Press Make It Mine to register your scanner. Please have your User Name, User Password, and Company ID ready to complete the registration process. You will receive an email from an iTero representative with login credentials and detailed information on the how to Make it Mine procedure.

Connecting to hidden networks

When setting up the iTero scanner, dental practices may connect to private internet connections on their Wi-Fi. These private or hidden networks are regular networks that do not present their SSID or network name.


All private networks will be labeled "Hidden Networks" as options under listed networks. When connecting to a hidden network, clinicians will need to input the network name as well as the security key or password.

2.0 Logging in to the scanner application

Login

When the scanner is powered on, the login screen for the iTero Element will appear.

- 1 User information
- 2 Create account
- 3 Forgot password
- 4 Login

Please have your MyiTero account information ready when logging in to the iTero scanner. You will need your name, account email, and password. Fill out all the necessary fields and click the Login button.

Create Account

To create an account, please press the bottom button labeled Create Account. Enter an email and create an account to use on the iTero scanner and MyiTero.com.

Forgot Password


If you forgot your password, press the lower button labeled Forgot Password. Enter your registered account email and a temporary password will be sent to your email. Alternatively, answer your predetermined security question and gain access to the account from the scanner.

3.0 iTero Element main features

Home screen

The iTero Element home screen provides access to all the main features of the iTero software.

Below is an image of the iTero Element application's home screen.

- 1 Main home screen icons
- 2 Utility icons

Main desktop icons

There are four main home screen icons for the iTero Element software.

New Scan - To begin a new scan for a patient, select the New Scan icon. This home screen icon navigates to the iTero Rx form to fill out the prescription form for a patient and begin scanning.

Patients - To look up any patient on your system, select the Patients icon. This home screen icon navigates to the Patients page. Search for a specific patient, scroll through recently scanned patients, review scans, access iTero TimeLapse, view patient history, add new scans, or add a treatment plan.

Orders - To view the status of your orders, select the Orders icon. This home screen icon navigates to a list of orders in progress or past orders.

Messages - To read messages, select the Messages icon. This home screen icon navigates to the message inbox, where you may read on the latest notifications from the Align team regarding new updates, improved features, enhanced tools, and latest developments.

Utility icons

Battery - To view the charge on the iTero external battery, press the Battery icon. A battery indicator will appear and specify how long the battery has until the iTero Element may need to be charged again. When the battery icon displays a lightning bolt symbol, the battery is charging. Note: this feature only appears in supported iTero Element scanners.

Learning Center - To access training materials and educational videos for your iTero Element, press the Learning Center icon.

Lock - To log out of an account whenever the iTero Element is not in use, press the Lock icon. This will help ensure that the dental practice is HIPAA compliant and that all medical information is secure.

Settings - To adjust preferences for your iTero Element like wand configuration, localization, user settings, and more, press the Settings icon.

Help - To navigate icons for home screen and helpful hints, press the Help icon. A transparent Help overlay appears over the screen to aid in navigation of features and tools. At any time, tap the Help icon for assistance. Tap anywhere to close the Help screen and return to the relevant screen. Customer support is available on every Help overlay.

Touch screen gestures

iTero Element application supports touch screen (also known as multi-touch) gestures. These gestures are predefined motions used to interact with multi-touch devices.

Examples of common touch screen gestures:

Tap

Double tap

Long press

Scroll

Rotate

Swipe

2 finger tap

Zoom in/pinch

Zoom out

Pan

Wand controls

- 1 Side buttons
- 2 Touchpad

Side buttons - Press and release either of the side buttons on the wand to start or stop the scanning process. Press the buttons simultaneously to activate the touchpad.

Touchpad - Once activated, swipe left or right on the touchpad to navigate between scan segments. Press and hold the touchpad will enable the gyro function to rotate the model on the screen.

Demo Mode

To access demo mode, press the iTero Element logo on the upper left side of the home screen. A menu will appear to press Demo Mode or to cancel and return to the home screen. Demo Mode will appear with all the same functionality as when taking a normal scan. No data will be permanently saved.


- 1 Demo mode
- 2 Exit to Windows
- 3 Cancel

Demo Mode - Initiate Demo Mode for practice without saving any data.

Exit to Windows - To Exit the iTero interface and return to Windows.
Note: Only available on iTero Flex devices.

Cancel - To exit from the screen and return to normal use.


Demo Mode tag - The red tag indicator when using Demo Mode.

Striped background - A striped background will appear while scanning. It is designed to prevent scanning real patients in Demo Mode.

Practice any portion of the patient experience. Any scans taken during Demo Mode cannot be submitted for treatment. To exit Demo Mode, please go back to the home screen and press the iTero Element logo. The menu for Demo Mode will appear in order to exit Demo Mode.


4.0 Messages

Press the Messages icon from the home screen to view notifications, updates and other messages from Align Technology. On the left side of the page, search for specific messages by subject title quickly or, press on the right side of the screen and scroll down the page to search and find a specific message. To mark any message as unread, press Mark as Unread.

5.0 Orders

Tap on the Orders icon to view or review In Progress or Past Orders.

In Progress – Scans that have not been submitted.

Past Orders – Orders that have been submitted. Press on any order to check the status of the scan.

Clicking on a past order shows the following options, depending on the case type.

View Rx – View the prescription for this order.

Viewer – View the scans for this order.

Add Rx – Add a prescription for this order (only applicable for Orthodontic orders and available up to 21 days after a scan).

Invisalign doctors may also select Invisalign features such as Invisalign Outcome Simulator.

6.0 Scan

To start a new scan, press the New Scan icon from the home screen or press the new scan button from the Patient page.


Prescription icon - The first step in the New Scan workflow. Fill out the prescription form with patient information, case type, and any notes to include.

Scan icon - The second step in the New Scan workflow. Scan the patient for a new scan to be captured.

View and Review icon - The third step in the New Scan workflow. Inspect the scan and examine that all adequate data has been captured.

Send Case icon - The final step in the New Scan workflow. Sign the confirmation form and send the case to its destination for storage or additional work needed.

- 1 Patient information
- 2 Case type
- 3 Notes
- 4 Treatment stage

For Restorative cases, select the case type and restoration needs for the appropriate indication.

Note: For implant cases (scan body indication), availability of scan bodies may differ from country to country. Even if a brand or model of scan bodies are available on the iTero Rx, it is the responsibility of the dental professionals to use parts available in their regions.


For Orthodontic cases, such as Invisalign or iRecord, select the case type and enter in patient information.

- 1 Select a lab to send the case
- 2 Select the restoration type from the dropdown

After selecting the case type for a patient, please begin scanning.

To start a scan, tap the scan icon.


- 1 Real time modeling
- 2 View finder
- 3 Color
- 4 Scan feedback

Scan feedback - To view any missing tooth, soft tissue, or oral structure during the scan.

Color - To switch from scan to color view.


View finder - View the wand's current positioning.

During scan processing, the system will inform clinicians if re-scanning is needed for additional data capture and/or bite imaging.

After review and confirmation that the scan has been captured with sufficient data, the following features can be used to view the scan in color or monochrome, edit the scan, remove any extra scanned artifacts, and erase additional scanned materials.

Occlusal map – View a patient's occlusal map with color indications.

Eraser - To erase any selected parts of the scan in order to re-scan for better data capture.

Edge trim – Remove any excess or unnecessary artifact from the scan.

After reviewing the scan, sign the Rx form and send the scan to complete the case.


Signature box – Sign your signature here with a finger. Check “Save Signature” to store a signature away for future ease and use. Check “Don’t show again” to eliminate this step.

Clear – To clear a signature from the signature box.

Confirm & Send – Confirm scan has been reviewed and send the scan to a lab or iTerо storage.

7.0 Patients

On the home screen, press the icon labeled Patients to locate a patient and their scans. Press on a patient to add a new scan or view previous scans

Patient page – View any patients, see historical scans associated with the patient, access iTerо TimeLapse, or add a new scan.


Search bar – Search for specific patients using their names or chart numbers.

8.0 iTero TimeLapse

iTero TimeLapse compares previously captured patient's 3D scans to allow a clinician to visualize changes of a patient's teeth, tooth structures, and oral soft tissues.

Patients scanned on a regular basis can be analyzed using the iTero TimeLapse tool. In the

Patient's profile, select two scans from the timeline at the bottom of the screen. Once two scans are selected, press the Compare Selected button to analyze the scans.

Timeline – Scans taken over time of the patient.


Scan selection checkboxes – Selected scans to compare and visualize in iTero TimeLapse.

Compare Selected – Analyze scans using iTero TimeLapse.

Color Legend – Highlighted areas with the most changes between scans with rich color. Scale measurements are in millimeters.

Loupe – To view areas in detail and visualize points of interest in the Animation Window. Drag the loupe from the animation window to the area of interest in order to magnify potential treatment spots. The areas in which appear in the loupe will be magnified within the animation window.

Animation window – Visualize fine details of changes for the scans selected.

Color visibility – Toggle between monochrome and color in the Animation Window.

Play/pause – To play or pause the Animation Window.

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