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
INSTRUCTIONS FOR USE

INCRAFT AAA Stent Graft System, in modifications

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
DOLPH McGRATH,
1, Old Waterford Road,
Cionmel,
CO. TIPPERARY, IRELAND
Notary Public for the
County of Tipperary.
Commissioned for Life.

2016

CORDIS CASHEL, CAHIR ROAD, CASHEL, CO. TIPPERARY

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INCRAFT AAA Stent Graft System, in modifications:

I. INCRAFT AAA Stent Graft System, in a set:

- 1) INCRAFT AAA Delivery System 4,7 mm x 54 cm with Aortic Bifurcate Prosthesis of following sizes 22 mm; 26 mm; 30 mm – 1 pc.
- 2) Stent implant card – 1 pc.
- 3) Instruction for use – 1 pc.

II. INCRAFT AAA Stent Graft System, in a set:

- 1) INCRAFT AAA Delivery System 5,3 mm x 54 cm with Aortic Bifurcate Prosthesis of following sizes 34 mm – 1 pc.
- 2) Stent implant card – 1 pc.
- 3) Instruction for use – 1 pc.

III. INCRAFT AAA Stent Graft System, in a set:

- 1) INCRAFT AAA Delivery System 4,0 mm x 77 cm with Aortic Bifurcate Prosthesis of following sizes: 10 mm x 8,2 cm; 10 mm x 10,1 cm; 10 mm x 12 cm; 10 mm x 13,8 cm; 13 mm x 8,2 cm; 13 mm x 10,1 cm; 13 mm x 12,0 cm; 13 mm x 13,8 cm; 16 mm x 8,2 cm; 16 mm x 10,1 cm; 16 mm x 12,0 cm; 16 mm x 13,8 cm; 20 mm x 8,2 cm; 20 mm x 10,1 cm; 20 mm x 12,0 cm; 20 mm x 13,8 cm – 1 pc.
- 2) Stent implant card – 1 pc.
- 3) Instruction for use – 1 pc.

IV. INCRAFT AAA Stent Graft System, in a set:

- 1) INCRAFT AAA Delivery System 4,3 mm x 77 cm with Aortic Bifurcate Prosthesis of following sizes 24 mm x 10,1 cm; 24 mm x 12,0 cm; 24 mm x 13,8 cm – 1 pc.
- 2) Stent implant card – 1 pc.
- 3) Instruction for use – 1 pc.

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Address of Manufacturing Facilities:

Cordis de Mexico, S.A. DE C.V.
Circuito Interior Norte No. 1820, Parque
Industrial Salvarcar, Cd. Juarez, Chihuahua,
Mexico 32574

**Russian Federation Authorized
Representative:**

Johnson & Johnson LLC
Russian Federation, 17, Krylatskaya street
building 3, Moscow, 121614, Russia
Tel.: + 7 (495) 580 77 77;
Fax: + 7 (495) 580 78 78

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The INCRAFT AAA Stent Graft System in modifications (hereinafter INCRAFT AAA Stent-Graft System) is intended for the endovascular treatment of patients with infrarenal abdominal aortic aneurysms

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
Each prosthesis is constructed of a seamless, low porosity, woven polyester graft supported by a series of short, electropolished, laser-cut, self-expanding nitinol stent-rings throughout the entire length. The nitinol stent-rings are sutured to the inner surface of the graft material. In addition to the stents being visible under fluoroscopy, radiopaque markers are sewn onto each component to aid visualization and to facilitate accurate placement.

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The Aortic Bifurcate Prosthesis is deployed first into the cranial portion of the infra-renal aorta. It has a flared bare trans-renal stent with 8 or 10 laser-cut barbs depending on the cranial diameter. The barbs help keep the prosthesis in place.

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
The Aortic Bifurcate Prosthesis has one main trunk with two sealing stent and a taper stent that divides into the ipsilateral and contralateral legs, supported by a series of z-stents. While the diameter of the trunk varies by code, the lengths of the trunk and legs, as well as the diameters of the legs are constant. The Aortic Bifurcate Prosthesis is manufactured in four trunk diameter sizes (22, 26, 30 and 34mm). Refer to table 1.

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The Iliac Limb Prostheses are deployed into the legs of the Aortic Bifurcate Prosthesis and into the ipsilateral and contralateral iliac vessels. The overlap between the Aortic Bifurcate Prosthesis and the Iliac Limb Prosthesis can vary between 2cm and 5cm on the ipsilateral side, and between 2cm and 4cm on the contralateral side.

The Iliac Limb Prostheses could also be used as iliac extensions by placing them into a previously deployed Iliac Limb Prosthesis to gain additional exclusion

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length.

The Iliac Limb prosthesis has a series of z-stents cranially, one or more taper stents (if other than a straight configuration), and a diamond sealing stent caudally. The cranial diameter is always constant while the length and the caudal diameter of the Iliac Limb Prosthesis could vary by code.


The Iliac Limb Prostheses are available in 5 different distal diameters (10, 13, 16, 20 and 24mm) and in four different lengths (8, 10, 12, and 14cm). Referto Table 2.

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Each prosthesis is loaded into a delivery system which facilitates controlled deployment of the prosthesis into the intended locations under fluoroscopic guidance (see Figure 4). There are two variations of the delivery system: one for the Aortic Bifurcate Prosthesis, and one for the Iliac Limb Prosthesis.

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The Aortic Bifurcate Delivery System has an integrated sheath introducer along with a hemostatic valve to facilitate component exchanges during the procedure.

The size of the integrated sheath introducer varies depending on the diameter of the prosthesis it contains. For prosthesis diameters of 22, 26, and 30mm, the inner diameter of the integrated sheath introducer is 13F (outer diameter of 14F). For the prostheses diameter of 34mm, the inner diameter of the integrated sheath introducer is 15F (outer diameter of 16F).

The outer surface of the integrated sheath introducer has a lubricious (hydrophilic) coating at the distal end to facilitate introduction into the vasculature.

Each Aortic Bifurcate Delivery System handle is labeled with "AB" to indicate that it contains an Aortic Bifurcate Prosthesis and a number that indicates the trunk diameter of the Aortic Bifurcate Prosthesis.

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The Delivery System of the Iliac Limb Prosthesis is similar to that of the Aortic Bifurcate Prosthesis except for its size and that it does not have an integrated sheath introducer. The Iliac Limb Delivery System can be delivered through the integrated sheath introducer of the Aortic Bifurcate System.

The Iliac Limb Delivery System has a 12F outer diameter for prosthesis diameters between 10mm and 20mm, and a 13F outer diameter for the 24mm diameter prosthesis (refer to Table 2).

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The outer surface of each Iliac Limb Delivery System has a lubricious (hydrophilic) coating at the distal end to facilitate introduction into the vasculature.

Each Iliac Limb Delivery System handle is labeled with "IL" to indicate that it contains an Iliac Limb Prosthesis and two numbers in the format AAxBB that indicate the size of Iliac Limb Prosthesis where AA equals the Limb Diameter in mm and BB equals the Limb Length in cm.

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The INCRAFT AAA Stent Graft System is intended for the endovascular treatment of patients with infrarenal abdominal aortic aneurysms with the following characteristics:

- Femoral access vessels should be adequate to fit the selected delivery system
- Proximal neck length $\geq 10\text{mm}$
- Aortic neck diameters $\geq 17\text{mm}$ and $\leq 31\text{mm}$
- Aortic neck suitable for suprarenal fixation
- Infrarenal and suprarenal neck angulation $\leq 60^\circ$
- Iliac fixation length $\geq 15\text{mm}$
- Iliac diameters $\geq 7\text{mm}$ and $\leq 22\text{mm}$
- Minimum overall AAA treatment length (proximal landing location to distal landing location) $\geq 128\text{mm}$
- Morphology suitable for aneurysm repair

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The INCRAFT AAA Stent Graft System is contraindicated for:

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- Patients with a known contraindication to undergoing angiography or anticoagulation.

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Carefully observe all warnings and cautions noted throughout these instructions as failure to do so may result in injury to the patient.

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• The use of the INCRAFT AAA Stent Graft System requires that physicians be specially trained in endovascular abdominal aortic aneurysm repair techniques, including experience with high resolution fluoroscopy and radiation safety.

• A vascular surgical team should be available while the implant procedure is in progress in case a conversion to an open surgical repair is required.

• An increase in aneurysm size and/or persistent endoleak may lead to aneurysm rupture.

• For single use only. Do not resterilize or re-use. Re-use, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failures, which in turn may lead to injury, illness or death of the patient.

• Note product "Use By" date.

• After use, all components used and packaging materials may be a potential biohazard. Handle and dispose of in

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accordance with the accepted medical practice and with applicable local, state and federal laws and regulations.

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- Irregular calcification and/or plaque may compromise the fixation and sealing of the implant, especially at the cranial and caudal sealing zones.

- Inadequate seal zone length may result in increased risk of leakage into the aneurysm or migration of the prosthesis.

- All patients should be advised that endovascular treatment of infra-renal Abdominal Aortic Aneurysms requires lifelong, regular follow-up to assess their health and the performance of the implanted endovascular prosthesis.

- When selecting a bifurcate prosthesis, attention should be given to the abdominal treatment length from the lowest renal artery to the aortic bifurcation. If this length is less than the length of the contralateral length of the aortic bifurcate prosthesis (8.6 cm) then it could result in increased difficulty with cannulating the contralateral gate.

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Non-clinical testing has demonstrated that the INCRAFT AAA Stent Graft is MR Conditional. It can be scanned safely in both 1.5Tesla (T) and 3T magnetic resonance (MR) systems under the following conditions:

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- Spatial gradient field ≤ 2500 Gauss/cm (25T/m)
- Maximum whole-body-averaged specific absorption rate (SAR) of 4 W/kg for 15 minutes of scanning
- Patients with an INCRAFT AAA Stent Graft implanted in the abdominal aorta may safely undergo MRI for Normal Mode and First Level Controlled Operating Mode of the MR System, as defined in IEC Standard 60601-2-33.

In non-clinical testing, the INCRAFT AAA Stent Graft produced a maximum temperature rise of less than 6°C during 15 minutes of continuous MR scanning in the first level controlled mode at a whole-body averaged specific absorption rate (SAR) of 4W/kg using a static field of 1.5T and 3T and a GE Signa HDx MR scanning system. The calculations do not take into consideration the cooling effects of blood flow, and therefore, actual in-vivo rises are expected to be lower. The effect of heating in the MRI environment for fractured struts is not known.

The INCRAFT AAA Stent Graft should not migrate in the 1.5T and 3T MRI environments. Non-clinical testing of an individual stent, two overlapping stents and prosthesis assemblies in the 1.5T and 3T MR systems demonstrated magnetically induced torque and translational attraction on the stents at levels less than those imposed by a patient's daily activity. MRI at 1.5T or 3T may be performed immediately following the implantation of the prosthesis. Non-clinical testing has not been performed to rule out the possibility of prosthesis migration at magnetic field strengths above 3T.

MR image quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the

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stent. Optimization of MR imaging parameters is recommended. The image artifact extends approximately 2mm and 5mm from the device, both inside and outside the device lumen when scanned in nonclinical testing using the sequence spin echo and gradient echo respectively, in a 3 Tesla GE Signa HDx MR system with a whole body coil.

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Potential adverse events and potential device risks include, but are not limited to following:

1. Amputation
2. Anesthesia complications
3. Aneurysmal enlargement
4. Aneurysm sac rupture
5. Aortic damage (perforation, dissection, bleeding, rupture)
6. Arterial or venous thrombosis
7. Bleeding, hematoma or coagulopathy
8. Bowel complications (e.g., ileus, transient ischemia, infarction, necrosis)
9. Cardiac arrhythmia
10. Cardiac complications
11. Cardiac failure or infarction
12. Claudication
13. Contrast toxicity
14. Death
15. Edema
16. Embolism
17. Endoleak
18. Fever
19. Gastrointestinal complications
20. Genitourinary complications (e.g., ischemia, erosion, fistula, incontinence, hematuria)
21. Hematoma (surgical)
22. Hepatic failure
23. Impotence
24. Infection
25. Lymphatic complications
26. Neurological complications (e.g. CVA, TIA)
27. Paralysis or paraparesis
28. Post-implant syndrome

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Warning:

Failure to comply with sizing requirements outlined in the table above may result in prostheses leaks, compromised flow, prostheses migration, compromised long term device durability or other complication or adverse event.

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INCRAFT AAA Stent-Graft System consists of:

1. Stent-graft pre-installed into the delivery system – 1 pc.
2. Stent implant card – 1 pc.
3. Instruction for use – 1 pc.

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contraindicated, an alternative anticoagulant should be considered.

2. Following standard practices, gain vascular access on both ipsilateral and contralateral sides.

3. Following standard practices, perform angiography with a diagnostic catheter placed from the contralateral side (anteroposterior, oblique and lateral views as necessary) to confirm patient eligibility, correct device component sizing and deployment locations. Leave the diagnostic catheter in place at the level of the renal arteries.

4. Always use fluoroscopy while advancing and removing guidewires, catheters, and the INCRAFT AAA Stent Graft System.

5. From the ipsilateral side advance a 0.035" (0.89mm) stiff guidewire to the level of the renal arteries.

Cautions:

- Failure to use an appropriate guidewire may result in vessel trauma and compromise deliverability and/or performance of the delivery system.

- When advancing the guidewires, catheters, and the INCRAFT AAA Stent Graft System into the abdominal aorta, do not disturb the thrombus mass within the aneurysm. Doing so may dislodge emboli, which can cause distal embolization. If distal embolization should occur, use conventional treatment methods.

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Prepare each Delivery System as needed. Use Figure 4 (below) as a reference.

Remove the appropriately sized Aortic Bifurcate and two Iliac Limb Delivery

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
The Delivery System components are identified in Figure 4.

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Figure 4. Delivery System Component Identification Guide

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
Radiopaque markers provide a reference for proper alignment when deploying the prosthesis components (Figures 5). When describing the orientation of this product, cranial refers to the portion of the prosthesis that is closer to the head of the patient. Caudal refers to the portion of the prostheses that is closer to the foot of the patient.

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The following Warnings and Precautions apply throughout the implant procedure:

- Ensure that the Delivery System Handle and Delivery System Sheath are parallel with the patient's leg. Excessive angulation where the White Handle Component meets the Delivery System Sheath may prevent Delivery System Sheath retraction.

- Do not handle the delivery system by the Gold Handle Component since rotation

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of the Gold Handle Component during positioning may cause premature deployment of the prosthesis.

- Before deployment ensure that the Delivery System is straight and without any slack.

- Do not torque the Delivery System more than 90° without confirming rotational response at the distal end of the device (Contralateral Side Marker- Figure 5).

- Do not start deployment until the delivery system is accurately placed within the vasculature and ready for deployment.

- When positioning the loaded Delivery System, hold only the White Handle Component.

- Do not rotate the Gold Handle Component in a counter clockwise direction, as this may result in an inability to deploy the prosthesis.

- When deploying the prosthesis, be sure to hold the White Handle Component of the Delivery System firmly against a stationary object.

- Prosthesis components cannot be resheathed or drawn back into the Delivery System without compromising the system, even if the prosthesis component is only partially deployed.

- If the outer sheath is accidentally withdrawn exposing the prosthesis, the device will prematurely deploy and may be incorrectly positioned.

- Failure to position the Bifurcate Cranial Edge Markers within the healthy infrarenal aortic neck may result in prosthesis leaks or require a bail-out procedure.

- Failure to position the prosthesis below the lowest renal ostium may result in occlusion of the renal arteries.

- Failure to position the Limb Caudal Edge Markers cranial to the internal iliac artery origin may result in occlusion of the internal iliac artery.

- Always use fluoroscopy to verify the prostheses is completely released from the Delivery System. Incomplete retraction of the Delivery System Sheath or

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incomplete displacement (pull) of the Fixation Release Wire could lead to dislodgement of the prosthesis when you remove the Delivery System subcomponent from the patient.

- Use fluoroscopic guidance to advance the Delivery System and to detect kinking or alignment problems with the stent graft system. Do not use excessive force to advance or withdraw the Delivery System when resistance is encountered. If the Delivery System kinks during insertion, do not attempt to deploy the stent graft component. Remove the device and insert a new Delivery System.

- Prosthesis migration or incorrect prosthesis deployment may require surgical intervention.

- Exercise particular care in areas that are difficult to navigate, such as areas of stenosis, intravascular thrombus, calcification or tortuosity, or where excessive resistance is experienced, as vessel or catheter damage could occur. Consider performing balloon angioplasty at the site of a narrowed or stenotic vessel, and then attempt to gently reintroduce the catheter delivery system. Also exercise care with device selection and correct placement/ positioning of the device in the presence of anatomically challenging situations such as areas of significant stenosis, intravascular thrombus, calcification, tortuosity and/or angulation which can affect successful initial treatment of the aneurysm.

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1. Using the existing diagnostic catheter previously introduced in the contralateral side, locate the renal arteries, aortic bifurcation, and the iliac bifurcations.

2. If not already present on the ipsilateral side, insert a 0.035"(0.89mm) stiff guidewire. Position the tip of the guidewire in the descending thoracic aorta. Wet the Aortic Bifurcate Delivery System (tip and approximately 25cm of the sheath) with saline to activate the hydrophilic coating.

3. Immediately prior to use, re-check that the Sheath Hemostasis Valve screw cap is secure.

Caution:

Failure to ensure the Sheath Hemostasis Valve screw cap is secure may result in connection separation and subsequent incomplete deployment of the prosthesis.

4. Insert the Aortic Bifurcate Delivery System over the ipsilateral guidewire and, under fluoroscopy, advance the Aortic Bifurcate Delivery System beyond the renal arteries and pull back such that the bottom end of the Aortic Bifurcate Cranial Edge Markers are just below the lowest renal artery.

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- Once cranial position has been identified, do not move the patient or imaging equipment, as it may compromise accuracy of prosthesis placement.

- The diagnostic catheter can be removed prior to deployment. However, if it is not removed until after deployment, ensure that the tip is straightened (for example, if it is pigtail catheter) with a guidewire before removal so that the prosthesis is not subject to migration.

- When aligning the position of prosthesis, be sure the fluoroscope is angled perpendicularly to the center line of the infrarenal aorta to avoid parallax or other source of visualization error that could impact proper positioning. Some cranial-caudal angulation of the image intensifier tube may be necessary to achieve this, especially if there is anterior angulation of the aneurysm neck.

1. Using the Contralateral Side Marker as a reference, orient the Aortic Bifurcate Prosthesis so that the Contralateral Side Marker is aligned with either the Contralateral Iliac Artery or the desired position (refer to Figure 5).

2. Ensure that the prosthesis position has been maintained with respect to the lowest renal artery. Confirm with injection through angiographic catheter as necessary. If needed, correct position of the Bifurcate Cranial Edge Markers under fluoroscopic guidance.

Caution:

Prior to deployment ensure that the Sheath Hemostasis Valve screw cap is secure.

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Prosthesis, proceed as follows: Place your hand firmly on the White Handle Component (refer to Figure 4) of the Aortic Bifurcate Delivery System to maintain the Delivery System position with respect to the renal arteries (use fluoroscopy and / or angiography to confirm position). While holding the White Handle Component with one hand, slowly turn the Gold Handle Component (refer to Figure 4) of the Aortic Bifurcate Delivery System in the clockwise direction as indicated by the directional arrow on the Gold Handle Component until the sheath begins to retract (use fluoroscopy to confirm that retraction has begun).

Warning:

Ensure that the position of the Aortic Bifurcate Cranial Edge Markers does not change while retracting the sheath. Once the trans-renal stent is unsheathed, the barbs on the trans-renal stent are exposed and may be engaged with the aortic wall. Repositioning or rotating the Aortic Bifurcate Delivery System with the barbs exposed may damage the trans-renal stent or the Aortic Bifurcate Delivery System, or cause injury to the artery or cause partial/complete coverage of aortic side-branches. Repositioning of the prosthesis may also result in peripheral and aortic side-branch vessel embolization. Failure to position the bottom edge of the Aortic Bifurcate Cranial Edge Markers below the lowest renal ostium may result in occlusion of the renal arteries. Rotating the delivery system after the sheath has commenced retracting may compromise the operator's ability to pull the secondary release wire in step 9.

Caution:

Ensure that the Delivery System Handle

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and Delivery System Sheath are parallel with the patient's leg. Excessive angulation where the White Handle Component meets the Delivery System Sheath may prevent Delivery System Sheath retraction.

4. Under fluoroscopy, keep turning the Gold Handle Component until the transrenal stent begins to expand. Assess device positioning relative to the renal arteries (refer to Figure 6) and make minor adjustments if needed. Continue retraction of the sheath until at least the Contralateral Leg of the Aortic Bifurcate Prosthesis is unsheathed (refer to Figure 5) as confirmed by fluoroscopy where the sheath tip is just caudal to the four Contralateral Leg-Gate Markers.

Warning:

Do not retract the guidewire beyond the contralateral leg-gate marker, as guidewire access can be lost or result in a kink of the delivery system.

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
1. While holding the White Handle Component stationary, locate the Fixation Release Wire on the handle of the Aortic Bifurcate Delivery System (refer to Figure 4). Unclip the release wire retainer and pull the fixation release wire until the cranial end of the trans-renal stent is released as confirmed by fluoroscopy.

Caution:

While activating the fixation release wire ensure that the white handle is held against a steady object to avoid compromising placement of the stent graft.

2. Under fluoroscopy, If the prosthesis has not been completely unsheathed, continue turning the Gold Handle Component (refer to Figure 4) until it is completely unsheathed. Use fluoroscopy to confirm complete deployment of the Aortic Bifurcate Prostheses.

3. Using standard techniques, slowly reposition the diagnostic catheter from its current position to a location directly above the aortic bifurcation.

4. Upon full deployment of Aortic Bifurcate Prosthesis, hold the integrated sheath introducer firmly to prevent movement while disconnecting the Sheath Hemostasis Valve to separate the sheath into two parts. While holding the sheath firmly, remove the Aortic Bifurcate Delivery System by the handle under fluoroscopic guidance to ensure safe removal of the device, leaving the integrated sheath introducer in place for use when deploying the ipsilateral limb and maintaining guidewire access.

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Caution:

If upon removal of the aortic bifurcate prosthesis delivery system excessive bleeding is seen through the valve, reposition the tip of the delivery system into the valve until the ipsilateral iliac limb prosthesis delivery system is inserted.

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The sequence of implanting the Ipsilateral and Contralateral Iliac Limb Prostheses could vary based on local practice and clinical situation.

The description below describes the situation where the Ipsilateral Iliac Limb Prosthesis is implanted first.

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1. Visualize the ipsilateral internal iliac artery origin and the caudal landing zone in the ipsilateral iliac artery by standard techniques.

2. Wet the outer sheath of the Iliac Limb Delivery System (tip to approximately 25cm of the sheath) with saline to assure activation of the hydrophilic coating.

3. Under fluoroscopic guidance, introduce the Iliac Limb Delivery System over the 0.035" (0.89mm) stiff guidewire through the integrated sheath introducer left in place from the Aortic Bifurcate Prosthesis deployment so that the Limb Cranial Edge Marker of the Iliac Limb Prostheses is aligned with the Maximum Overlap Marker of the Aortic Bifurcate Prostheses while ensuring that no cranial migration of the Aortic Bifurcate Prosthesis occurs.

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Warning:

Cranial migration of the Aortic Bifurcate Prosthesis may result in renal occlusion.

4. Under fluoroscopic guidance pull the sheath Hemostasis Valve back until the tip of the integrated Sheath Introducer is at least 1cm below the Limb Caudal Edge Marker to ensure that it does not interfere with deployment of the Iliac Limb Prostheses while also ensuring that the Delivery System Sheath Tip Marker remains intravascular, particularly if in the proximity of the vascular access site (refer to Figure 7).

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
Warning:

Failure to pull the Delivery System sheath back may result in incomplete deployment of the Iliac Limb Prosthesis.

5. Under fluoroscopic guidance, make final caudal adjustment of the Iliac Limb Delivery System by retracting the Iliac Limb Delivery System to the desired deployment position in the iliac artery, ensuring that the Limb Caudal Edge Marker is aligned with the iliac landing zone and the Limb Cranial Edge Marker is located between the Ipsilateral Minimum Overlap Marker and the Maximum Overlap Marker (refer to Figure 5).

Note:

The Limb Caudal Graft Edge is 0.5-1.5mm below the Limb Caudal Edge Marker.

Warning:

The Limb Cranial Edge Marker of the Ipsilateral Iliac Limb Prosthesis must be between the Minimum Overlap Marker and Maximum Overlap Marker of the Aortic Bifurcate Prosthesis. Failure to do so may result in prosthesis leaks or flow restrictions through one or both Iliac Limb Prostheses.

Caution:

Before deployment, ensure that the Delivery System Handle and Delivery System Sheath are parallel with the patient's leg. Excessive angulation where the White Handle Component meets the Delivery System Sheath may prevent Delivery System Sheath retraction.

6. To deploy the Ipsilateral Iliac Limb Prosthesis, proceed as follows: Place your hand firmly on the White Handle Component (refer to Figure 4) of the Iliac

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Limb Delivery System to maintain the Iliac Limb Delivery System position with respect to the distal landing zone (use fluoroscopy to confirm position). While holding the White Handle Component with one hand, slowly turn the Gold Handle Component (refer to Figure 4) in the clockwise direction as indicated by the directional arrow on the Gold Handle Component to retract the sheath.

Caution:

Ensure that the position of the Limb Cranial Edge Marker does not change while retracting the sheath.

Warning:

Do not advance the Iliac Limb Delivery System once the distal tip of the sheath is retracted. Doing so could displace or damage the Aortic Bifurcate Prosthesis.

7. Keep turning the Gold Handle Component until the Iliac Limb Prosthesis is fully deployed. Continue to turn the Gold Handle Component until the Sheath Tip Marker is at least 1cm beyond the Limb Caudal Marker, as confirmed by fluoroscopy.

8. Unclip the Release Wire Retainer and pull it (refer to Figure 4) until it comes completely out of the handle, thereby releasing and securing the cranial end of the Iliac Limb Prosthesis.

Warning:

Do not advance or withdraw the Iliac Limb Delivery System after the Iliac Limb Prosthesis has been deployed until the Fixation Release Wire has been retracted. Doing so could result in adverse events, including damage to the vessel wall or prosthesis.

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Note:

Unlike with Aortic Bifurcate Prosthesis, the fixation release wire for the Iliac Limb Delivery System must be completely removed from the Delivery System.

Caution:

Failure to remove the fixation release wire completely may result in caudal migration of the Iliac Limb during withdrawal of the delivery system.

9. Under fluoroscopy, carefully remove the Iliac Limb Delivery System from the patient, ensuring that the Iliac Limb Delivery System does not dislodge the prosthesis and guide wire access is maintained. Leave the Delivery System Sheath in place to use for balloon tacking after system deployment is completed.

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1. Using standard practice, cannulate the contralateral leg of the Bifurcated Aortic Prosthesis guided by the Contralateral Leg-Gate Markers (refer to Figure 5).

Notes:

- Consider rotating the C-arm to a different angle to ensure that the Contralateral Leg of the Aortic Bifurcate can be uniquely visualized.

- Ensure the Contralateral Leg has been cannulated. Oblique imaging may be helpful in determining actual guidewire position relative to the leg opening.

- In the event that the contralateral leg is cannulated prior to placing the ipsilateral limb, ensure that the appropriate

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limb is cannulated prior to deployment.

1. Exchange the guidewire to a 0.035" (0.89mm) stiff guidewire on the contralateral side.

2. Using standard techniques, visualize the contralateral internal iliac artery origin and the caudal landing zone in the contralateral iliac artery.

3. Remove the contralateral diagnostic catheter, if used, under fluoroscopic guidance. Unless already present, introduce an ipsilateral diagnostic catheter.

4. Wet the outer sheath of the Iliac Limb Delivery System Sheath (tip to approximately 25cm of the sheath) with saline to ensure activation of the hydrophilic coating.

5. Under fluoroscopic guidance, introduce the Iliac Limb Delivery System over a 0.035" (0.89mm) stiff guidewire so that the Cranial Marker of the Iliac Limb Prosthesis is aligned with the Maximum Overlap Marker of the Bifurcate Prostheses while ensuring that no migration of the Aortic Bifurcate Prosthesis occurs.

Warning:

Cranial migration of the Aortic Bifurcate Prosthesis may result in renal occlusion.

6. Under fluoroscopic guidance, make final caudal adjustment of the Iliac Limb Delivery System by retracting the Iliac Limb Delivery System to the desired deployment position in the iliac artery, ensuring that the Limb Caudal Edge Marker is aligned with the iliac landing zone and the Limb Cranial Edge Marker is located between the contralateral leg Minimum Overlap Marker and the Maximum Overlap

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Marker (refer to Figure 5).

Warning:

The Limb Cranial Edge Marker of the Contralateral Iliac Limb Prosthesis must be between the Minimum Overlap Marker and Maximum Overlap Marker (overlap zone) of the Aortic Bifurcate Prosthesis. Failure to do so may result in prosthesis leaks or flow restrictions through one or both Iliac Limb Prostheses.

Note:

If the Contralateral Iliac Limb Delivery System will not freely advance into the Aortic Bifurcate Prosthesis, remove the Contralateral Iliac Limb Delivery System and re-introduce the guidewire into the contralateral leg. The guidewire may have been advanced between a stent and the graft material, or advanced into the ipsilateral leg.

Caution:

Ensure that the Delivery System Handle and Delivery System Sheath are parallel with the patient's leg. Excessive angulation where the White Handle Component meets the Delivery System Sheath may prevent Delivery System Sheath retraction.

7. To deploy the Contralateral Iliac Limb Prosthesis, proceed as follows: Place your hand firmly on the White Handle Component (refer to Figure 4) of the Iliac Limb Delivery System handle to maintain the Delivery System position with respect to the distal landing zone (use fluoroscopy to confirm position). While holding the White Handle Component with one hand, slowly turn the Gold Handle Component (refer to Figure 4) in the clockwise direction as indicated by the directional arrow on the Gold Handle Component to retract the sheath.

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Delivery System must be completely removed from the Delivery System.

Warning:

Failure to remove the fixation release wire completely may result in caudal migration of the Iliac Limb during withdrawal of the delivery system.

3. Under fluoroscopy, carefully remove the Iliac Limb Delivery System from the patient, ensuring that the Delivery System does not dislodge the prosthesis and guide wire access is maintained. If a sheath introducer was used, leave the sheath in place to use for balloon tacking.

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Perform the following steps using appropriate endovascular and surgical techniques.

1. Under fluoroscopy, use an appropriately sized and compatible compliant balloon catheter to tack the cranial and caudal seal zones and the overlap regions of the modular components as well as any areas of potential restricted blood flow.

Caution:

Be careful not to displace the prostheses upon introducing and retracting the balloon catheter.

Note:

Care should be taken when inflating the balloon, especially with calcified, tortuous, stenotic, or otherwise diseased vessels. Ensure to not inflate the balloon outside or the graft material. Inflate balloon slowly. It is recommended that a backup balloon be available.

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Warnings:

- Over inflation of balloon can cause graft tears and/or vessel dissection or rupture.

- When expanding the prostheses, there is an increased risk of vessel injury and/or rupture, and possible patient death, if the balloon's proximal and distal radiopaque markers are not completely within the covered (graft fabric) portion of the prosthesis.

- Do not tack the trans-renal stent of the Aortic Bifurcate Prosthesis as it may cause vessel injury and vessel rupture and could snag onto the trans-renal stent of the Aortic Bifurcate Prosthesis.

- Cranial migration of the Aortic Bifurcate Prosthesis may result in renal artery occlusion. Caudal migration of the Iliac Limb Prosthesis may result in occlusion of the internal iliac arteries.

At the completion of the procedure, perform angiography to assess the prosthesis for proximal, modular overlap and distal endoleaks and to verify position of the implanted prosthesis in relation to the aneurysm, the renal as well as internal iliac arteries.

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Note:

Similar to iliac limbs, the amount of overlap is adjustable, however, minimum overlap is recommended as to prevent the potential multiple overlap of the bifurcate leg, iliac limb, and iliac extension, which may result in flow restrictions.

Note:

When using iliac leg extension, the minimum overlap varies by the iliac limb size being extended as per the Table 3.

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Note:

It is not possible to extend 10mm iliac limb.

1. Select the appropriate device as per Table 3.

Caution:

After removal from the package take care not to pull or snag the fixation release wire (refer to Figure 4) when handling and preparing the Delivery System.

2. Prior to use, flush the guidewire lumen with heparinized saline through the guidewire lumen flush connector of the Iliac Limb Delivery System (refer to Figure 4) until saline is observed exiting through the distal end of the device.

Note:

It is acceptable for small amounts of saline to be observed coming from locations other than the distal end of the device.

3. Visualize the caudal landing zone in the iliac artery by standard techniques.

4. Wet the outer sheath of the Iliac Limb Delivery System (tip and approximately 25cm of the sheath) with saline to assure activation of the hydrophilic coating.

5. Under fluoroscopic guidance, introduce the Iliac Limb Delivery System over the 0.035" (0.89mm) stiff guidewire. If using the integrated Sheath introducer left in place from the Aortic Bifurcate Prosthesis deployment, refer to steps 3-4 of the "Implant the Ipsilateral Iliac Limb Prosthesis" section.

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Note:

If deploying the Iliac Limb Prosthesis near the vascular access site, take extra care to ensure that sheath remains intravascular.

6. Under fluoroscopy, advance the Iliac Limb Delivery System beyond the targeted deployment site but not beyond the Maximum overlap marker of the aortic bifurcate and pull back such that the Extension (Limb) Cranial Edge


Marker is more cranial than three short non-tapered z-stents (refer to Figure 8).

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Caution:

There is no minimum marker band on the iliac limb being extended.

1. Under fluoroscopic guidance, ensure proper overlap as Figure 8.

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2. If a catheter sheath introducer is used, ensure it is pulled back beyond the caudal edge marker of the iliac extension.

3. To deploy the Iliac Limb Prosthesis, proceed as follows: Place your hand firmly on the White Handle Component (refer to Figure 4) of the Iliac Limb Delivery System to maintain the Iliac Limb Delivery System position with respect to the distal landing zone (use fluoroscopy to confirm position). While holding the White Handle Component with one hand, slowly turn the Gold Handle Component (refer to Figure 4) in the clockwise direction as indicated by the directional arrow on the Gold Handle Component to retract the sheath.

Caution:

Ensure that the position of the Limb Extension Cranial Edge Marker does not change while retracting the sheath.

Warning:

Do not advance the Iliac Limb Delivery System once the distal tip of the sheath is retracted. Doing so could displace or damage the Aortic Bifurcate Prosthesis or Limb Prosthesis previously deployed.

4. Keep turning the Gold Handle Component until the Iliac Limb Prosthesis

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is fully deployed. Continue to turn the Gold Handle Component until the Sheath Tip Marker is at least 1 cm beyond the Limb Caudal Marker, as confirmed by fluoroscopy.

5. Unclip the Release Wire Retainer and pull it (refer to Figure 4) until it comes completely out of the handle, thereby releasing and securing the cranial end of the Iliac Limb Prosthesis.

Warning:

Do not advance or withdraw the Iliac Limb Delivery System after the prosthesis has been deployed until the Fixation Release Wire has been retracted. Doing so could result in adverse events including damage to the vessel walls or prosthesis.

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Caution:

Failure to remove the fixation release wire completely may result in caudal migration of the Iliac Limb during withdrawal of the Delivery System.

6. Under fluoroscopy, carefully remove the Iliac Limb Delivery System from the patient, ensuring that the Iliac Limb Delivery System does not dislodge the prosthesis and guide wire access is maintained. Leave the Delivery System Sheath in place to use for balloon tacking.

7. Introduce appropriately sized and compatible molding balloon and tack the overlap and seal areas of the Iliac limb extension.

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Caution:

Be careful not to displace the prosthesis upon introducing and retracting the balloon catheter.

Note:

Care should be taken when inflating the balloon, especially with calcified, tortuous, stenotic, or otherwise diseased vessels. Inflate slowly. Ensure to not inflate the balloon outside the graft material. It is recommended that a backup balloon be available.

Warnings:

- Over inflation of balloon can cause graft tears and/or vessel dissection or rupture.
- When expanding the prosthesis, there is an increased risk of vessel injury and/or rupture, and possible patient death, if the balloon's proximal and distal radiopaque markers are not completely within the covered (graft fabric) portion of the prosthesis.

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It is recommended that physicians conduct regular examinations and imaging for the patient's lifetime. Follow-up imaging should be decided based upon the physician's clinical assessment of the patient pre- and post-implantation of the stent graft. After endovascular graft placement, patients should be regularly monitored for perigraft flow, aneurysm growth or changes in the structure or position of the endovascular graft. Annual imaging is recommended, including

- abdominal radiographs to examine device integrity (stent fracture, separation between bifurcated device and proximal cuffs or limb extensions, if applicable); and

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• contrast and non-contrast CT to examine aneurysm changes, perigraft flow, patency, tortuosity and progressive disease. If renal complications or other factors preclude the use of image contrast media, abdominal radiographs and duplex ultrasound may provide similar information.

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There is no express or implied warranty, including without limitation any implied warranty of merchantability or fitness for a particular purpose, on the product(s) described in this publication. Under no circumstances shall Manufacturer be liable for any direct, incidental, or consequential damages other than as expressly provided by specific law. No person has the authority to bind Manufacturer to any representation or warranty except as specifically set forth herein. Descriptions or specifications in Manufacturer printed matter, including this publication, are meant solely to generally describe the product at the time of manufacture and do not constitute any express warranties. Manufacturer will not be responsible for any direct, incidental, or consequential damages resulting from reuse of the product.

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INCRAFT AAA Stent-Graft System is transported by all modes of transport in covered vehicles in accordance with the goods transportation rules operating on each type of transport.

Medical devices are transported and stored in a cool, dry, dark environment at a temperature from - 40 °C to + 54 °C and a relative humidity of 20% to 80% (non-condensing).

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Shelf life of the medical device is 2 years. It is not allowed to use the device after expiration date indicated on the device package.

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Dispose of all used devices in accordance with hospital policy for bio hazardous materials.

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