071-00001 Rev 21

Instructions for Use
Cordis S.M.A.R.T.® Flex Vascular Stent System

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
Upute za upotrebu
Sustav za primjenu vaskularnog stenta Cordis S.M.A.R.T.® Flex

Navodila za uporabo
Sistem vaskularnega stenta Cordis S.M.A.R.T.® Flex

DOLPH McGRATH,
1, Old Waterford Road,
Cionmel,
CO. TIPPERARY, IRELAND
Notary Public for the
County of Tipperary.
Commissioned for Life.

DOLPH McGRATH
15th July 2016


Signed on behalf of Cordis Cashel
Nicola Parry
QA/RA Lead

Cordis Cashel,
Cahir Road, Cashel,
Co. Tipperary, Ireland.
Registration No. 424180

Content / Table des matières / Inhalt / Contenido / Contenido / Conteudo / Inhoud / Innhald / Sisallus / Innehåll / Innhold / Περιεχόμενα / Obsah / Spis treści / Obsah / Съдържание / Cuprins / Sadržaj / Tarkib / Saturs / Turmyn / İçindekiler / Содержание / 目录 / 目錄 / 目録 / 目次 / Isi / Sadržaj / Vsebina	
English	Page
Français	Page
Deutsch	Seite
Italiano	Pagina
Español	Página
Português	Página
Nederlands	Pagina
Dansk	Side
Suomi	Sivu
Svenska	sid
Norsk	Side
Ελληνικά	Σελίδα
Cestina	Strana
Magyar	Oldal
Polski	Strona
Slovensky	Strana
Български	Страница
România	Pagina
Eesti	Lk
Latviešu	lpp
Lietuvių k	P
Türkçe	Sayfa
Pyjama	ᲞᲠᲣᲚᲘᲗ
ភាសាខ្មែរ	ទំព័រ
ភាសាខ្មែរ	ទំព័រ
ភាសាខ្មែរ	ទំព័រ
Bahasa Indonesia	Halaman
Hrvatski	Stranica
Slovensčina	Strana

[illegible]

Figure 1 - S.M.A.R.T.® Flex Vascular Stent System / Figure 1 - Système d'endoprothèse vasculaire S.M.A.R.T.® Flex / Abbildung 1 - S.M.A.R.T.® Flex-Stentssystem für Gefäße / Figura 1 - Sistema stent vascolare S.M.A.R.T.® Flex / Figura 1: Sistema de stent vascular S.M.A.R.T.® Flex / Figure 1: Sistema de stent vascular S.M.A.R.T.® Flex / Abeefteling 1 - S.M.A.R.T.® Flex vasculair systeem / Figure 1 - S.M.A.R.T.® Flex vaskulär stentsystem / Kuva 1 - S.M.A.R.T.® Flex-verisuontenstijärjestelmä / Fig. 1 - S.M.A.R.T.® Flex vaskulær stentsystem / Figure 1 - S.M.A.R.T.® Flex vaskulært stentsystem / Ekvivo 1 - Συστήμα αγγειοκάθετην S.M.A.R.T.® Flex / Ölsbräck 1 - System cívlného stentu S.M.A.R.T.® Flex / 1. ábra - S.M.A.R.T.® Flex vaskuláris stentszűrrendszere / Rysunek 1 - System naczyńnowej stentu S.M.A.R.T.® Flex / Obrazok 1 - Vaskulyarny stentovyy sistem S.M.A.R.T.® Flex / Dnyrnyia 1 - Cúdovalno stentno sústema S.M.A.R.T.® Flex / Figura 1 - Sistemul de stent vascular S.M.A.R.T.® Flex / Joons1 - S.M.A.R.T.® Flex vaskulaare stendisüsteem / Silika 1 - Sistem vaskularnega stenta S.M.A.R.T.® Flex / 1. attels - S.M.A.R.T.® Flex vaskulara stenta silika / 1 pav. - Kraujagimiai stento sistema S.M.A.R.T.® Flex / Şekil 1 - S.M.A.R.T.® Flex Vasküler Stent Sistemi / Pvc. 1 - Набор для установки сосудистых стентов S.M.A.R.T.® Flex / figo 1 - S.M.A.R.T.® Flex 血管支架系统 / 圖 1 - S.M.A.R.T.® Flex 血管支架系統 / 그림 1 - S.M.A.R.T.® 플렉스 혈관 스텐트 시스템 / yslr1 - S.M.A.R.T.® Flex Vaskulær Stent System / Gambar 1 - S.M.A.R.T.® Flex Vascular Stent System / Слика 1 - Васкуларни систем за стент S.M.A.R.T.® Flex

- [illegible]

English

Stent: Manufactured by ethylene oxide gas. Do not use if package is opened or damaged. For single use only. Re-use may increase the risk of injury, inadequate performance, and infection. Non-pyrogenic. Radiopaque. MR Conditional.

These recommendations are designed to serve only as a general guideline. They are not intended to supersede institutional protocols or professional clinical judgment concerning patient care.

Stent Implant Patient Information Card

Designed for the patient to carry in their wallet for reference along with their healthcare cards. This patient information card provides information pertaining to the stent including the model and lot number of the stent, the date of the procedure and the location of the stent in their leg. The card also provides manufacturing information and MR Conditions.

INTENDED USE

The stent is intended as a treatment for atherosclerotic peripheral femoral artery lesions and proximal popliteal lesions.

CONTRAINDICATIONS

- Inability to tolerate antithrombotic or antiplatelet therapies
- Pregnancy
- Hypercoagulable state
- 100% occluded superficial femoral artery that can not be successfully crossed
- Not intended for use in lesions that cannot be pre-dilated
- Not intended for use in a perforated vessel

WARNINGS

- Persons allergic to nickel (nickel titanium) may suffer an allergic reaction to this implant.
- Appropriate antiplatelet/anticoagulant therapy should be administered pre and post procedure.
- The long term outcome of repeated dilation of endovascularized stents is unknown at present.
- When multiple stents are used resulting in stent-to-stent contact, they should be of similar composition.

PRECAUTIONS

- This product should only be used by physicians trained and experienced in diagnostic and interventional vascular techniques. Standard techniques for interventional vascular procedures should be employed.
- Use of the S.M.A.R.T.® Flex Vascular Stent System requires fluoroscopic control.
- Do not use power injection systems with the delivery system.
- Use in patients with a history of contrast sensitivity is not recommended unless the patient can be adequately pre-medicated.
- The effect of heating in the MRI environment for overlapping stents or stents with fractured struts is not known.
- Safety and effectiveness has not been demonstrated in:
 - Lesions that are either totally or densely calcified
 - Pediatric patients

System Handling - Precautions

- Do not attempt to remove the stent from the delivery system before use.
- Store at ambient room conditions out of direct sunlight.
- Do not use beyond the "Use By" date.
- Do not use if the system appears to be damaged.
- Do not expose the delivery system to organic solvents.
- Do not use if the outer sheath of the device cannot be flushed.
- Do not use if the inner lumen of the guidewire tube cannot be flushed.

Stent Placement - Precautions

- Use caution when crossing a deployed stent with any adjunct device.
- If resistance is encountered at any time during the insertion procedure, do not force passage. Resistance may cause damage to stent or system.
- If resistance occurs during movement through the sheath, carefully withdraw the stent system.
- Once stent deployment has begun, the stent must be fully deployed. The system is not designed for stent repositioning or recapturing.
- If resistance is felt when initially retracting the outer deployment sheath, do not force deployment. Carefully withdraw the stent system without deploying the stent.

POTENTIAL ADVERSE EVENTS - Potential hazards and side effects include, but are not limited to:

- Abrupt stent closure
- Allergic reaction to nitinol
- Anginal/coronary ischemia
- Aneurysm or Pseudoaneurysm in vessel or vessel access site
- Arrhythmia
- Atheroembolization (Blue Toe Syndrome)
- AV fistula formation
- Death
- Embolization (air, plaque, thrombus, device or other)
- Hematomas/Hemorrhage
- Hypotension/hypertension
- Infection and/or pain at the access site
- Infection secondary to contamination of the stent may lead to sepsis and hemodynamic instability
- Intimal injury/dissection
- Ischemia requiring intervention (bypass or amputation of toe, foot or leg)
- Limb loss
- Myocardial infarction (MI)
- Overstretching the artery may result in rupture
- Renal insufficiency or failure
- Restenosis of the stented segment
- Septicemia/bacteremia
- Stent malapposition/migration
- Stent structure fracture
- Stroke
- Thrombosis/thrombus
- Tissue necrosis
- Vasospasm
- Vessel perforation or rupture
- Worsened claudication/rest pain

DEVICE DESCRIPTION

The S.M.A.R.T.® Flex Vascular Stent System is a self-expanding stent (S.M.A.R.T.® Flex Stent) pre-mounted and delivered with a simple retractable sheath over the wire delivery system. The S.M.A.R.T.® Flex Stent is a nearly fully connected stent made from superelastic nitinol tubing and is constructed through the integration of helically wound struts with helical flexible coils. Both the strut elements and the helical coils provide radial stiffness and an expansion mechanism. The S.M.A.R.T.® Flex Stent is very bendable and will conform with minimal to no fish scaling. Upon deployment, the stent achieves its predetermined diameter and exerts a constant, gentle outward force to maintain an unobstructed lumen.

The S.M.A.R.T.® Flex Vascular Stent System, as shown in Figure 1, is composed of a 6F outer sheath (1) and a coaxial inner shaft. The inner shaft is composed of a pusher tube (2) and a guidewire tube (7). The pusher tube (2) axially abuts the proximal end of the stent (9) and the guidewire tube (7) is coaxial to the stent (9) and pusher tube (2). The guidewire tube (7) terminates distally in a catheter tip (8) and originates proximally in a Luer hub (6) designed to accept a 0.035" guidewire. The lumen created by the inside diameter of the guidewire tube (7) is flushed prior to the procedure by injecting fluid via the proximal Luer hub (6). The 6F outer sheath (1) connects proximally to a handle (3) that has a Y connection with female luer (4) and a Tuohy Borsst valve (5) at the proximal portion of the pusher tube (2). The self-expanding stent (9) is constrained within the space between the outer sheath (1) and the guidewire tube (7). This space is flushed prior to the procedure by injecting fluid via the Y connection with female luer (4).

Stent positioning about the target lesion is achieved prior to deployment utilizing the distal stent markers (10) and the proximal placement marker (12). These sets of markers indicate the approximate position of the stent after deployment. The proximal stent markers (11) are proximal to the proximal placement marker (12) and indicate the amount the stent will shorten during deployment. The stent system is designed to shorten from proximal to distal during deployment, such that once the distal portion of the stent has been placed at the target location, it will not move during the remainder of the stent deployment and most stent shortening occurs from the proximal end. The distal and proximal radiopaque stent markers (10 and 11) are attached to each end of the stent. The proximal placement marker (12) is a radiopaque band attached to the guidewire tube (7).

For stent deployment, the Tuohy Borsst valve (5) must be loosened to allow free movement of the pusher tube (2). To retract the outer sheath (1) and deploy the stent (9), remove any slack in the system, grip the outer sheath (1) and handle (3) with one hand and the pusher tube (2) with the other. Move the outer sheath (1) proximally relative to the pusher tube (2), while holding the pusher tube (2) in a fixed position. The delivery system has a certain amount of compliance and there may be some lag between the proximal movement of the outer sheath (1) / handle (3) and the retraction of the outer sheath (1) at the distal end of the stent (9). The position of the delivery system may need to be adjusted to keep the distal stent markers (10) at the target site during this lag period. Once the distal end of the stent (9) starts to deploy continue to retract the proximal outer sheath (1) and handle (3) while holding the pusher tube (2) still until the stent (9) is fully deployed and the proximal stent markers (11) have expanded. Note: The proximal placement marker (12) may move during the stent deployment and may not coincide with the location of the proximal stent markers (11).

PROCEDURE

- The percutaneous placement of the stent should be done in an angiography procedure room under fluoroscopic guidance. Diagnostic angiography should be performed to map out the extent of the lesion and the collateral flow and measure the length of the target lesion and diameters of the reference vessel (proximal and distal to the lesion). Patient preparation and sterile precautions should be the same as for any angioplasty procedure.
- Select a stent size diameter according to the reference vessel diameter measured and the stent size selection table below.

Stent Diameter (mm)	Reference Vessel Diameter (mm)
3	3.5-4.5
4	4.5-5.5
5	5.5-6.5
6	6.5-7.5

- Select Stent Length - Measure the length of the target lesion to determine the length of stent required. Select a stent length that allows for the proximal and distal aspects of the stent to cover the entire target lesion.
 - The Maximum Constrained Length is the length of the stent in at the smallest indicated reference vessel diameter. This length is indicated on the system with the proximal placement marker (see (12) in Fig. 1).
 - The Minimum Constrained Length is the length of the stent in at the largest indicated reference vessel diameter.

Labeled Stent Length (mm)	Maximum Constrained Length (mm) Per Proximal placement marker	Minimum Constrained Length (mm)
30	33	32
40	45	42
60	66	63
80	88	84
100	110	105
120	132	126
150	165	157
200	210	200

- Access the treatment site utilizing the appropriate accessory equipment. The stent delivery system is compatible with a 6F (or larger) introducer sheath.

- Introduce an appropriate 035" guidewire through the access catheter or introducer across the distal segment of the target lesion. It is recommended to use a sheath to cross over the aortic bifurcation. It is recommended that the patient will receive a therapeutic dose of heparin such that the ACT is greater than 250.
- Pre-Dilate the Lesion - Using standard balloon angioplasty techniques dilate the lesion to a size equivalent to the reference vessel diameter. Remove the balloon catheter from the patient maintaining lesion access with the guidewire.
- Prepare Stent Delivery System
 - Open the outer box to reveal the pouch containing the stent and delivery catheter.
 - After careful inspection of the pouch, looking for damage to the sterile barrier, carefully peel open the outer pouch and extract the inner pouch. Carefully peel open the inner pouch and remove the backboard with carrier tube which holds the stent delivery system. **Warning: DO NOT use if pouch is opened or damaged.**
 - Set the backboard with carrier tube on a flat surface. Remove the stent/delivery system from the carrier tube. Examine the device for damage. If it is suspected that the sterility has been compromised or the device is damaged do not use the device.
 - If the distal tip (8) is not seated in the outer sheath (1), loosen the Tuohy Borsst valve (5) on the handle (3) and retract the pusher tube (2) such that the distal tip (8) will seat in the outer sheath (1). Tighten the Tuohy Borsst valve by rotating the valve hub clockwise.
 - Check to ensure that the Tuohy Borsst valve (5) on the handle is tightened on the pusher tube (2).
 - Use a 3-10 cc syringe to flush the outer sheath (1) with sterile heparinized saline through the female luer (4) on the handle. Flush until only a few drops of saline exit the distal end of the outer sheath. Complete system flushing may require 2-3 flushings with a 1cc syringe. **Warning: If the outer sheath (1) cannot be flushed do not use the device.**
 - Use a 3-10 cc syringe to flush the inner lumen of the guidewire tube with sterile heparinized saline through the proximal Luer hub (6) attached to the pusher tube (2), until saline flows out of the guidewire lumen at the distal tip. **Warning: If the inner lumen of the guidewire tube (7) cannot be flushed do not use the device.**
- Introduce the Stent Delivery System
 - Advance the device over the guidewire and through the introducer to the target site. If resistance is met during delivery system introduction, the system should be withdrawn and another system used.
 - Under fluoroscopic guidance, position the distal stent markers (10) distal to the target lesion and proximal placement marker (12) proximal to the target lesion.
 - Straighten the proximal part of the delivery system as much as possible and keep the handle in a stable position.
 - The Tuohy Borsst valve (5) on the handle (3) must be loosened to allow free movement of the pusher tube (2).
- Deploy the Stent
 - Stent deployment must be performed under fluoroscopic guidance.
 - Grip the outer sheath (1) and handle (3) with one hand and the pusher tube (2) with the other.
 - Move the outer sheath (1) proximally relative to the pusher tube (2), while holding the pusher tube (2) in a fixed position.
 - The position of the delivery system may need to be adjusted to keep the distal stent markers (10) at the appropriate location.
 - Once the distal end of the stent (9) starts to deploy continue to retract the proximal outer sheath (1) and handle (3) while holding the pusher tube (2) still until the stent (9) is fully deployed and the proximal stent markers (11) have expanded.
- Note: The proximal placement marker (12) may move during the stent deployment and may not coincide with the location of the proximal stent markers (11) after deployment.
- If resistance is felt during retraction of the outer sheath (1), do not force deployment. Carefully withdraw the stent system without deploying the stent.
- Under fluoroscopic guidance, withdraw the entire delivery system as one unit, over the guidewire while keeping the guidewire in position, and out of the body. Remove the delivery device from the guidewire. **Warning: Do not forcefully remove delivery system from introducer/ sheath - If resistance is met, remove sheath and delivery system as a unit.**
- Perform an arterial angiogram to verify full deployment. If incomplete expansion exists within the stent at any point along the lesion, post deployment balloon dilation can be performed at the discretion of the physician. Remove the balloon from the patient.
- When clinically appropriate, remove the guidewire, sheath, and any other accessory equipment from the body and achieve hemostasis of the access site.

MR COMPATIBILITY INFORMATION
The S.M.A.R.T.® Flex Stent (IMPLANT) was determined to be MR-Conditional according to the terminology specified in the American Society for Testing and Materials (ASTM) International, Designation: F2503-05, Standard Practice for Marking Medical Devices and Other Items for Safety in the Magnetic Resonance Environment. ASTM International, 100 Bar Harbor Drive, PO Box C700, West Conshohocken, Pennsylvania, 19428

Non-clinical testing demonstrated that the IMPLANT is MR Conditional. A patient with this implant can be scanned safely immediately after placement under the following conditions:

Static Magnetic Field

- Static magnetic field of 3-Tesla or less
- Highest spatial gradient magnetic field of 720-Gauss/cm

MR-Related Heating

In non-clinical testing of 8x150 and 8x200mm samples, the S.M.A.R.T.® Flex Stent produced the following maximum temperature rises and specific absorption rate (SAR) measurement using calorimetry during MRI performed for 15-min of scanning (i.e. per pulse sequence) in 1.5-Tesla/64-MHz (Magnetom, Siemens Medical Solutions, Malvern, PA. Software: Magnets4, Version Syngo MR 2002B DHMS Active-shielded, horizontal field scanner) and 3-Tesla/128-MHz (Excite, HDx, Software: 14X.M5, General Electric Healthcare, Milwaukee, WI) MR systems:

Maximum Measurement	1.5-T	3-T
MR system reported, whole body averaged SAR	3.8-W/kg	1.0-W/kg
Calorimetry measured values, whole body averaged SAR	3.1-W/kg	2.8-W/kg
Highest temperature change	2.7°C	4.4°C

Artifact Information

MR image quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the S.M.A.R.T.® Flex Stent. Therefore, optimization of MR imaging parameters to compensate for the presence of this device may be necessary. The maximum artifact size (i.e. as seen on the gradient echo pulse sequence) extends approximately 10-mm relative to the size and shape of this implant.

Pulse Sequence	T1 SE	T1 SE	GRE	GRE
Signal Void Size	2792	102	4254	302
Plane Orientation	Parallel	Perpendicular	Parallel	Perpendicular

WARRANTY DISCLAIMER

Although this product has been manufactured under carefully controlled conditions, Cordis Catheter has no control over the conditions under which this product is used. Cordis Catheter therefore disclaims all warranties, both express and implied, with respect to the product including, but not limited to, any implied warranty of merchantability or fitness for a particular purpose. Cordis Catheter shall not be liable to any person or entity for any medical expenses or any direct, indirect, or consequential damages caused by any use, defect, failure or malfunction of the product, whether a claim for such damages is based upon warranty, contract, tort or otherwise. No person has any authority to bind Cordis Catheter to any representation or warranty with respect to the product. The exclusions and limitations set out above are not intended to and should not be construed so as to contravene mandatory provisions of applicable law. If any part or term of this Disclaimer of Warranty is held to be illegal, unenforceable or in conflict with applicable law by a court of competent jurisdiction, the validity of the remaining portions of this Disclaimer of Warranty shall not be affected, and all rights and obligations shall be construed and enforced as if this Disclaimer of Warranty did not contain the particular part or term held to be invalid.

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Dolph McGrath
25th July 2016

Nicola Parry
Signed on behalf of Cordis Cashel
Nicola Parry
QA/RA Lead

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
- [illegible]

English

Sterile: Sterilized by ethylene oxide.

Do not resterilize - single use only. Re-use may increase the risk of injury, inadequate performance, and infection.

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DEVICE DESCRIPTION

The **S.M.A.R.T.® Flex** Vascular Stent System is a self-expanding stent (**S.M.A.R.T.® Flex Stent**) pre-mounted and delivered with a simple retractable sheath over-the-wire delivery system. The **S.M.A.R.T.® Flex Stent** is a fully connected stent made from superelastic nitinol tubing and is constructed through the integration of helically wound struts with helical flexible coils. Both the strut elements and the helical coils provide radial stiffness and an expansion mechanism. The **S.M.A.R.T.® Flex Stent** is very flexible and will conform with minimal to no fish scaling. Upon deployment, the stent achieves its predetermined diameter and exerts a constant, gentle outward force to maintain an unobstructed lumen.

INDICATIONS FOR USE

The **S.M.A.R.T.® Flex Stent** is intended for use in the common and external iliac arteries to improve luminal diameters in patients with symptomatic vascular stenotic and/or occlusive diseases.

CONTRAINDICATIONS

- Inability to tolerate antithrombotic or antiplatelet therapies
- Pregnancy
- Hypercoagulable state
- Not intended for use in lesions that cannot be pre-dilated
- Not intended for use in a perforated vessel

WARNINGS

- Persons allergic to nitinol (nickel titanium) may suffer an allergic reaction to this implant.
- Appropriate antiplatelet/anticoagulant therapy should be administered pre and post procedure.
- The long term outcome following repeat dilation of endothelialized stents is unknown at present.
- When multiple stents are used resulting in stent-to-stent contact, they should be of similar composition

PRECAUTIONS

- This product should only be used by physicians trained and experienced in diagnostic and interventional vascular techniques. Standard techniques for interventional vascular procedures should be employed.
- Use of the **S.M.A.R.T.® Flex** Vascular Stent System requires fluoroscopic control.
- Do not use power injection systems with the delivery system.
- Use in patients with a history of contrast sensitivity is not recommended unless the patient can be adequately pre-medicated.
- The effect of heating in the MRI environment for overlapping stents or stents with fractured struts is not known.
- Safety and effectiveness has not been demonstrated in:
 - Lesions that are either totally or densely calcified
 - Pediatric patients

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POTENTIAL ADVERSE EVENTS Potential hazards and side effects include, but are not limited to:

- Abrupt stent closure
- Allergic reaction to nitinol
- Angina/coronary ischemia
- Aneurysm or Pseudoaneurysm in vessel or vessel access site
- Arrhythmia
- Atheroembolization (Blue Toe Syndrome)
- AV Fistula formation
- Death
- Embolization (air, plaque, thrombus, device or other)
- Fever
- Hematoma/Hemorrhage
- Hypotension/hypertension
- Infection and/or pain at the access site
- Infection secondary to contamination of the stent may lead to sepsis and hemodynamic instability
- Intimal injury/dissection
- Ischemia requiring intervention (bypass or amputation of toe, foot or leg)
- Limb loss
- Myocardial infarction (MI)
- Overstretching the artery may result in rupture
- Renal insufficiency or failure
- Restenosis of the stented segment
- Septicemia/bacteremia
- Stent malapposition/migration
- Stent structure fracture
- Stroke
- Thrombosis/thrombus
- Tissue necrosis
- Vasospasm
- Vessel perforation or rupture
- Worsened claudication/rest pain

MRI COMPATIBILITY INFORMATION

Non-clinical testing demonstrated that the IMPLANT is **MR Conditional**. A patient with this implant can be scanned safely immediately after placement under the following conditions:

- Normal operating mode of the MR system and use of whole body transmit coil.
- Maximum time for WB-SAR scan 15 min.
- 1.5-Tesla/64-MHz MR System reported whole body average SAR value of 3.8-W/kg was associated with a calculated whole body average SAR value of 3.1-W/kg
- 3-Tesla/128-MHz MR System reported whole body average SAR value of 3.0-W/kg was associated with a calculated whole body average SAR value of 2.8-W/kg
- Static magnetic field of 3-Tesla or 1.5-Tesla
- Spatial gradient magnetic field of 720-Gauss/cm or less

MRI-Related Heating

In non-clinical testing, the IMPLANT produced the following temperature rises during MRI performed for 15-min in 1.5-Tesla (1.5-Tesla/64-MHz, Magnetom; Siemens Medical Solutions, Malvern, PA, Software Numaris4, Version Syngo MR 2002B DHHS) and 3-Tesla (3-Tesla/128-MHz, Excite, Software G3.0-0528, General Electric Healthcare, Milwaukee, WI) MR systems, as follows:

Highest Temperature Change	MRI Condition	Stent Configuration
+2.4°C	1.5-T/64-MHz	Single 10mm X 100mm
+3.5°C	3-T/128-MHz	Single 10mm X 100mm

Therefore, the MRI-related heating experiments for the IMPLANT at 1.5- and 3-Tesla using a transmit/receive RF body coil of MR system reported whole body averaged SARs of 3.8-W/kg (i.e., associated with a calorimetry value of 3.1-W/kg) and 3.0-W/kg (i.e., associated with a calorimetry value of 2.8-W/kg), respectively, indicated that the greatest amount of heating that occurred in association with these specific conditions was equal to or less than 3.5°C in the maximum size 10mm x 100mm single stent configuration. The effect of heating in the MRI environment for overlapping stents or stents with fractures is not known.

Artifact Information

MR image quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the IMPLANT. Therefore, optimization of MR imaging parameters to compensate for the presence of this implant may be necessary.

Dimensions: IMPLANT, 8mm diameter, 150mm length

Pulse Sequence	T1-SE	T1-SE	GRE	GRE
Signal Void Size	1,585 mm ²	76 mm ²	1,870 mm ²	95 mm ²
Plane Orientation	Parallel	Perpendicular	Parallel	Perpendicular

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- Measure the length of the target lesion to determine the length of stent required. Select a stent length that allows for the proximal and distal aspects of the stent to cover the entire target lesion.
 - The **Maximum Constrained Length** is the length of the stent in the smallest indicated reference vessel diameter.
 - The **Minimum Constrained Length** is the length of the stent in the largest indicated reference vessel diameter.

Labeled Stent Length (mm)	Maximum Constrained Length (mm)	Minimum Constrained Length (mm)
20	22	21
30	33	32
40	45	42
60	66	63
80	88	84
100	110	105

- Access the treatment site utilizing the appropriate accessory equipment. The stent delivery system is compatible with a 6Fr (or larger) introducer sheath.
- Using a 6Fr (or larger) introducer sheath, introduce an appropriate .035" guidewire through the access catheter or introducer across the distal segment of the target lesion. It is recommended to use a sheath to cross over the aortic bifurcation. It is recommended that the patient will receive a therapeutic dose of heparin such that the ACT is greater than 250.
- Pre-dilate the lesion:** Using standard balloon angioplasty techniques dilate the lesion minimally to a size equivalent to the reference vessel diameter up to the labeled stent size. Remove the balloon catheter from the patient maintaining lesion access with the guidewire.
- Prepare the stent delivery system:**
 - Open the outer box to reveal the pouch containing the stent and delivery catheter.
 - After careful inspection of the pouch, looking for damage to the sterile barrier, carefully peel open the outer pouch and extract the inner pouch. Carefully peel open the inner pouch and remove the backboard with carrier tube which holds the stent delivery system. **Warning: DO NOT use if either the inner or outer pouch is opened or damaged.**

- Set the backboard with carrier tube on a flat surface.

Remove the stent/delivery system from the carrier tube. Examine the device for damage. If it is suspected that the stent has been compromised or the device is damaged, do not use the device.

- If the distal tip (8) is not seated in the outer sheath (1), loosen the Tuohy Bort valve (5) on the handle (3) and retract the pusher tube (2) such that the distal tip (8) will seat in the outer sheath (1). Tighten the Tuohy Bort valve by rotating the valve hub clockwise.
- Check to ensure that the Tuohy Bort valve (5) on the handle is tightened on the pusher tube (2) to prevent premature stent deployment and leakage during flushing.
- Use a 1-3 cc syringe to flush the outer sheath (1) with sterile heparinized saline through the female luer (4) on the handle. Flush until only a few drops of saline exit the distal end of the outer sheath. Complete system flushing may require 2-3 flushings with a 1cc syringe. **Warning: If the outer sheath (1) cannot be flushed do not use the device.**
- Use a 3-10 cc syringe to flush the inner lumen of the guidewire tube with sterile heparinized saline through the proximal luer hub (6) attached to the pusher tube (2), until saline flows out of the guidewire lumen at the distal tip. **Warning: If the inner lumen of the guidewire tube (7) cannot be flushed do not use the device.**

- Introduce the stent delivery system:** Advance the device over the guidewire and through the introducer to the target site. If resistance is met during delivery system introduction, the system should be withdrawn and another system used.

- Under fluoroscopic guidance, position the distal stent markers (10) distal to the target lesion and proximal placement marker (12) proximal to the target lesion.

- Straighten the proximal part of the delivery system to remove as much slack as possible and keep the handle in a stable position.

- The Tuohy Bort valve (5) on the handle (3) must be loosened to allow free movement of the pusher tube (2).

- Deploy the stent:**

- Stent deployment must be performed under fluoroscopic guidance.
- Grip the outer sheath (1) and handle (3) with one hand and the pusher tube (2) with the other.
- Move the outer sheath (1) proximally relative to the pusher tube (2), while holding the pusher tube (2) in a fixed position. **Warning: Do not move pusher tube (2) in distal direction to facilitate deployment.**
- The position of the delivery system may need to be adjusted to keep the distal stent markers (10) at the appropriate location.
- Once the distal end of the stent (9) starts to deploy continue to retract the proximal outer sheath (1) and handle (3) while holding the pusher tube (2) still until the stent (9) is fully deployed and the proximal stent markers (11) have expanded.

Note: The proximal placement marker (12) may move during the stent deployment and may not coincide with the location of the proximal stent markers (11) after deployment.

- If resistance is felt during retraction of the outer sheath (1), do not force deployment. Carefully withdraw the stent system without deploying the stent.

- Under fluoroscopic guidance, withdraw the entire delivery system over the guidewire as one unit, while keeping the guidewire in position and out of the body. Remove the delivery device from the guidewire.

Warning: Do not forcefully remove delivery system from introducer/sheath - if resistance is met, remove sheath and delivery system as a unit.

- Perform an arterial angiogram to verify full deployment. If incomplete expansion exists within the stent at any point along the lesion, post deployment balloon dilation can be performed at the discretion of the physician. Remove the balloon from the patient.

Warning: Do not attempt to reposition stent after start of deployment.

- When clinically appropriate, remove the guidewire, sheath, and any other accessory equipment from the body and achieve hemostasis at the access site.

Patient Materials

Designed for the patient to carry in their wallet for reference along with their healthcare cards. This patient information card provides information pertaining to the stent including the model and lot number of the stent, the date of the procedure, and the location of the stent in their leg. The card also provides manufacturing information and MR Conditions.

WARRANTY DISCLAIMER

Although this product has been manufactured under carefully controlled conditions, **Cordis Cashel** has no control over the conditions under which this product is used. **Cordis Cashel** therefore disclaims all warranties, both express and implied, with respect to the product including, but not limited to, any implied warranty of merchantability or fitness for a particular purpose. **Cordis Cashel** shall not be liable to any person or entity for any medical expenses or any direct, incidental or consequential damages caused by any use, defect, failure or malfunction of the product, whether a claim for such damages is based upon warranty, contract, tort or otherwise. No person has any authority to bind **Cordis Cashel** to any representation or warranty with respect to the product. The exclusions and limitations set out above are not intended to, and should not be construed so as to contravene mandatory provisions of applicable law. If any part or term of this Disclaimer of Warranty is held to be illegal, unenforceable or in conflict with applicable law by a court of competent jurisdiction, the validity of the remaining portions of this Disclaimer of Warranty shall not be affected, and all rights and obligations shall be construed and enforced as if this Disclaimer of Warranty did not contain the particular part or term held to be invalid.

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Russian Federation Authorized Repre- sentative:

Johnson & Johnson LLC
Russian Federation, 17, Krylatskaya street
building 2, Moscow, 121614
Tel.: + 7 (495) 580 77 77;
Fax: + 7 (495) 580 78 78

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*S.M.A.R.T. Flex. Vascular Stent System
for ileal artery*

- Preinfarction angina in patients com-
plained for frequent and intensive attacks
of stenocardia.
- Severe transmural heart attack posing
threat for life.
- Maintenance of artificial system of blood
circulation with tendency for narrowing
- Atherosclerosis

CORDIS CASHEL, CAHIR ROAD, CASHEL, CO. TIPPERARY

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This medical device is used only in the conditions of medical and therapeutic institutions by qualified medical personnel.

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The medical S.M.A.R.T. Flex not negatively affect humans and the environment during use, transportation, and storage through following the instructions for use. Special personnel protection and the patient in use of device it is not required.

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The medical device is transported using all available modes of transportation in the covered vehicles according to the transportation rules forced on each transportation mode.

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Medical device is stored at a room temperature, at a place protected from direct sunlight.

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The medical device is transported in cool, dry, dark place, at temperature from -40°C to $+54^{\circ}\text{C}$ and air relative humidity from 20% to 80% (without condensation).

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S.M.A.R.T. Flex system in operation are stable against influence of climatic factors for temperature condition from $+32^{\circ}\text{C}$ to $+42^{\circ}\text{C}$ and to influence of biological liquids.

CORDIS CASHEL, CAHIR ROAD, CASHEL, CO. TIPPERARY

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After use all used components and packing materials can be potentially biologically dangerous. Disposal is carried out according to the accepted medical practice, and also according to local, state and federal laws and provisions.

CORDIS CASHEL, CAHIR ROAD, CASHEL, CO. TIPPERARY

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071-00001 Rev 21

of medical device

S.M.A.R.T. Flex vascular stent system for femoral
and popliteal arterial segments

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

2017

DOLPH McGRATH,
1, Old Waterford Road,
Cionnail.

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CO. TIPPERARY, IRELAND
Notary Public for the
County of Tipperary.
Commissioned for Life.

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S.M.A.R.T. Flex. Vascular Stent System,
in the design styles:

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I. S.M.A.R.T. Flex. Vascular Stent System, for femoral and popliteal arterial segment comprised of

1. 80 cm long stent delivery system with S.M.A.R.T. Flex stent of the following design styles:
5 mm x 30 mm, 5 mm x 40 mm, 5 mm x 60 mm, 5 mm x 80 mm, 5 mm x 100 mm, 5 mm x 120 mm, 5 mm x 150 mm, 6 mm x 30 mm, 6 mm x 40 mm, 6 mm x 60 mm, 6 mm x 80 mm, 6 mm x 100 mm, 6 mm x 120 mm, 6 mm x 150 mm, 7 mm x 30 mm, 7 mm x 40 mm, 7 mm x 60 mm, 7 mm x 80 mm, 7 mm x 100 mm, 7 mm x 120 mm, 7 mm x 150 mm, 7 mm x 200 mm, 8 mm x 30 mm, 8 mm x 40 mm, 8 mm x 60 mm, 8 mm x 80 mm, 8 mm x 100 mm, 8 mm x 120 mm, 8 mm x 150 mm, 8 mm x 200 mm – 1 pc.
2. Personal card of the patient with implanted stent - 1 pc.
3. Instructions for Use - 1 pc.

II. S.M.A.R.T. Flex. Vascular Stent System, for femoral and popliteal arterial segment comprised of

1. 120 cm long stent delivery system with S.M.A.R.T. Flex stent of the following design styles:
5 mm x 30 mm, 5 mm x 40 mm, 5 mm x 60 mm, 5 mm x 80 mm, 5 mm x 100 mm, 5 mm x 120 mm, 5 mm x 150 mm, 6 mm x 30 mm, 6 mm x 40 mm, 6 mm x 60 mm, 6 mm x 80 mm, 6 mm x 100 mm, 6 mm x 120 mm, 6 mm x 150 mm, 7 mm x 30 mm, 7 mm x 40 mm, 7 mm x 60 mm, 7 mm x 80 mm, 7 mm x 100 mm, 7 mm x 120 mm, 7 mm x 150 mm, 7 mm x 200 mm, 8 mm x 30 mm, 8 mm x 40 mm, 8 mm x 60 mm, 8 mm x 80 mm, 8 mm x 100 mm, 8 mm x 120 mm, 8 mm x 150 mm, 8 mm x 200 mm – 1 pc.

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2. Personal card of the patient with
implanted stent - 1 pc.
3. Instructions for Use - 1 pc.

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Burpee Materials Technology, LLC
15 Christopher Way, Eatontown, New
Jersey 07724, USA (США)

Address of Manufacturing Facilities:

Burpee Materials Technology, LLC
15 Christopher Way, Eatontown, New
Jersey 07724, USA

Cordis Cashel

Cahir Road, Cashel, Co. Tipperary, Ire-
land (Ирландия)

Cordis Cashel

Cahir Road, Cashel, Co. Tipperary, Ire-
land

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Authorized Representative of the manufacturer in the Russian Feder- ation:

Johnson & Johnson LLC
Russian Federation, Krylatskaya street 17,
building 2, Moscow, 121614
Tel.: +7 (495) 580 77 77;
Fax: +7 (495) 580 78 78

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S.M.A.R.T. Flex. Vascular Stent System for femoral and popliteal arterial segment:

- Acute and subacute thrombosis
- Thrombosis after surgeries
- Thrombosis after endovascular intervention
- Thrombosis of dialysis fistulas

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In case of adverse events that threaten the patient's life, the physician may decide to remove the product.

Deimplantation of the product is carried out using standard surgical procedures.

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The manufacturer warrants conformance of this medical device with the declared characteristics within expiry term provided that all requirements of the manufacturer's use instructions are met.

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I. S.M.A.R.T. Flex. Vascular Stent System, for iliac artery comprised of

1. 80 cm long stent delivery system with S.M.A.R.T. Flex stent of the following design styles:
9 mm x 20 mm, 9 mm x 30 mm,
9 mm x 40 mm, 9 mm x 60 mm,
9 mm x 80 mm, 9 mm x 100 mm, 10 mm x 20 mm,
10 mm x 30 mm, 10 mm x 40 mm, 10 mm x 60 mm,
10 mm x 80 mm, 10 mm x 100 mm – 1 pc.
2. Personal card of the patient with implanted stent - 1 pc.
3. Instructions for Use - 1 pc.

II. S.M.A.R.T. Flex. Vascular Stent System, for iliac artery comprised of

1. 120 cm long stent delivery system with S.M.A.R.T. Flex stent of the following design styles:
9 mm x 20 mm, 9 mm x 30 mm,
9 mm x 40 mm, 9 mm x 60 mm,
9 mm x 80 mm, 9 mm x 100 mm, 10 mm x 20 mm,
10 mm x 30 mm, 10 mm x 40 mm, 10 mm x 60 mm,
10 mm x 80 mm, 10 mm x 100 mm – 1 pc.
2. Personal card of the patient with implanted stent - 1 pc.
3. Instructions for Use - 1 pc.

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Burpee Materials Technology, LLC
15 Christopher Way, Eatontown, New Jersey
07724, USA (США)

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Burpee Materials Technology, LLC
15 Christopher Way, Eatontown, New Jersey
07724, USA

Cordis Cashel

Cahir Road, Cashel, Co. Tipperary, Ireland
(Ирландия)

Cordis Cashel

Cahir Road, Cashel, Co. Tipperary, Ireland

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Johnson & Johnson LLC
Russian Federation, Krylatskaya street 17,
building 2, Moscow, 121614
Tel.: + 7 (495) 580 77 77;
Fax: + 7 (495) 580 78 78

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© Cordis Corporation, 2014
S.M.A.R.T.® is a registered trademark of Cordis