JP

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CERTIFICATE

中国国际贸易促进委员会暨中国国际商会
China Council for the Promotion of International Trade is China Chamber of International Commerce

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China Council for the Promotion of International Trade
China Chamber of International Commerce

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THIS IS TO CERTIFY THAT: the seal of LYNMED MEDICAL TECHNOLOGY (BEIJING) CO., LTD. on the annexed DOCUMENT is genuine.

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证明地址 Website for verifying the certificate: <http://www.rzccpit.com/validate.html>

Date: 1st of June, 2022

Approved by

王鹏 (signature)

Wang, General manager

Lyncmed Medical Technology (Beijing) Co.Ltd.

INSTRUCTION FOR USE
FOR THE MEDICAL DEVICE

CLEAN+SAFE® Non-Sterile Nitrile Rubber Medical Examination Gloves

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1. Medical device name and versions

CLEAN+SAFE® non-sterile nitrile rubber medical examination gloves

Versions:

1. CLEAN+SAFE® non-sterile nitrile rubber medical examination gloves: powder-free with textured fingertips (blue). Sizes: XS, S, M, L, XL;
2. CLEAN+SAFE® non-sterile nitrile rubber medical examination gloves: powder-free with textured fingertips (purple). Sizes: XS, S, M, L, XL;

Hereinafter referred to as gloves, devices, CLEAN+SAFE gloves, CLEAN+SAFE® nitrile rubber gloves.

2. Information about the developer, manufacturer and authorized representative of the medical device manufacturer

Developer

Lyncmed Medical Technology (Beijing) Co. Ltd.

Address: Room 1601, Building No. 2, Zhubang 2000 Business Building, Balizhuangxili 99, Chaoyang District, 100022 Beijing, P.R. China
Tel. +86-010-86468058, E-mail: Contactus@LyncMed.com

Manufacturer

Lyncmed Medical Technology (Beijing) Co. Ltd.

Address: Room 1601, Building No. 2, Zhubang 2000 Business Building, Balizhuangxili 99, Chaoyang District, 100022 Beijing, P.R. China
Tel. +86-010-86468058, E-mail: Contactus@LyncMed.com

Authorized representative of manufacturer (producer)

PROTECO, JSC

Legal address: 192102, St. Petersburg, Salova st. 57, build. 1, lit. I.

Business address: 196128, St. Petersburg, Varshavskaya st. 5, bldg. 2, office 401, Business Center "Varshavsky"

Tel./fax: +7 (812) 635-88-90; E-mail: info@protecodent.ru

3. Intended use of the medical device with indication of an intended user

3.1. Intended use of the medical device

Clean+Safe® gloves are intended to provide two-sided barrier to protect both patients and healthcare personnel from cross-contamination during medical research, diagnostic and therapeutic procedures, as well as for use in cosmetology. Gloves are single use non-sterile devices.

3.2. Intended users

INSTRUCTION FOR USE
FOR THE MEDICAL DEVICE

CLEAN+SAFE® Non-Sterile Nitrile Rubber Medical Examination Gloves

Clean+Safe® nitrile rubber gloves are recommended for

- healthcare personnel in cosmetology clinics and healthcare facilities, such as hospitals, maternity hospitals, clinics and doctors' offices (including private) and other healthcare facilities.
- individuals using the device for its intended purpose specified by the manufacturer.

4. Device description and design

Clean + Safe® nitrile gloves, ambidextrous, are made of nitrile rubber, have a beaded wrist cuff that snugs the wrist. Bead is rolled inside the glove and eliminates roll-ups and accumulation of infectious agents between the bead and the external surface. The device has a powder-free surface, is a fingertip-textured. The device is single chlorinated to prevent the sides of gloves from sticking together and to ensure easy removal from the hand. The gloves are presented in two colors- blue and purple.

5. Indications, contraindications and possible side effects associated with intended use of medical device

5.1. Indications

The gloves are intended to protect hands and maintain hygiene standards in examination of patients, diagnostic and other manipulations where gloves must be used, in order to prevent infection of the patient and the physician.

5.2. Contraindications

Contraindication for use of Clean+Safe® nitrile rubber gloves is individual intolerance to raw materials components. If irritation occurs, stop using the device and seek medical attention.

5.3 Side effects

Side-effects that may occur during normal conditions of intended use: non-allergic contact dermatitis due to the individual intolerance to raw materials, including additives used in gloves production.

6. Precautions

Read instructions before use. Do not use after expiration date.

Do not wear the same pair of gloves for the care of more than one patient. When wearing gloves, change them during care of the same patient if moving from one body site to another one with different bacterial flora.

Hand hygiene must be performed before putting on and after removing gloves. After handrubbing or handwashing, let your hands dry completely before putting on gloves.

If gloves are contaminated with secretions, blood, etc. in order to avoid contamination of the hands during their removal, use a swab (cloth) moistened with a solution of a disinfectant (or antiseptic) to remove visible contaminants. Remove gloves; immerse them in a disinfectant solution and dispose of them. Treat hands with an antiseptic.

Gloves should be removed immediately, if they are damaged (or there is a suspicion that their integrity is compromised) or there is indication for hand hygiene.

INSTRUCTION FOR USE
FOR THE MEDICAL DEVICE

CLEAN+SAFE® Non-Sterile Nitrile Rubber Medical Examination Gloves

2. Usage procedure

Take out a glove from its original box.


2. Touch only restricted surface of the glove corresponding to the wrist (at the top of the cuff).

3. Don the first glove.

4. Take the second glove with bare hand and touch only restricted surface of glove corresponding to wrist.

5. To avoid touching the skin of the forearm with gloved hand, turn external surface of the glove to be donned on the folded fingers of the gloved hand, thus permitting to glove the second hand.

INSTRUCTION FOR USE
FOR THE MEDICAL DEVICE

CLEAN+SAFE® Non-Sterile Nitrile Rubber Medical Examination Gloves

Once gloved, hands should not touch surrounding objects the contact with which is an indication for glove change.

7. After use pinch one glove at the wrist level to remove it, without touching the skin of the forearm, and peel away from the hand, thus allowing the glove to turn inside out.

8. Hold the removed glove in the gloved hand and slide the fingers of the other hand between the glove and the wrist. Remove the second glove by rolling it down the hand and fold into the first glove.

9. Discard the removed gloves.

10. Then perform hand hygiene by rubbing them with hygienic antiseptic agents or by washing them with soap and water.

8. Technical and functional features

Physical properties:

1. Elongation at break – not less than 500 % (before accelerated aging), not less than 400 % after accelerated aging; force at break – not less than 6,0 H (before and after accelerated aging).
2. Hermeticity, inspection level /AQL – G1/2.5.
3. Residual powder content: not more than 2 mg/glove.

Dimensions:

Dimension	Specification				
Size	XS	S	M	L	XL
Width, mm	70-79	80-89	90-99	100-109	110-119
Length, mm	Not less than 240				
Thickness, mm	Palm – not less than 0,05				

INSTRUCTION FOR USE
FOR THE MEDICAL DEVICE
CLEAN+SAFE® Non-Sterile Nitrile Rubber Medical Examination Gloves

Finger tip – not less than 0,05

Packaging

Medical devices "Clean+Safe® non-sterile nitrile rubber medical examination gloves" are packed in a cardboard box in the amount – 100 pcs. (50 pairs).

Description of symbols on a package:

Table 1. Symbols and descriptions

Symbol	Description	Symbol	Description
	Manufacturer		Manufacturing date
	Lot number		Keep dry
	Non-sterile		Product must not be - disposed of as household waste
	Keep away from sunlight		Do not reuse
	Consult instructions for use		Reference number (article)
	CE compliance mark		Product does not contain natural rubber latex
	Do not use if package is damaged		

10. Requirements for maintenance and repair of medical device

The devices are disposable and easy to use, do not require special skills and training. They do not need service and repair.

11. Sterilization and disinfection

Clean+Safe® non-sterile nitrile rubber medical examination gloves are single use non-sterile devices. Sterilization and disinfection are not applicable.

INSTRUCTION FOR USE
FOR THE MEDICAL DEVICE
CLEAN+SAFE® Non-Sterile Nitrile Rubber Medical Examination Gloves

12. Conditions of storage, operation and transportation

Storage. Devices should be stored in well-ventilated areas away from heat, direct sunlight, direct exposure of lighting systems and electrical equipment. Packages must be protected from contamination, dust, moisture and mechanical damage.

Storage conditions:

- temperature: 10 to +30 °C;
- relative humidity: up to 80% at 25 °C;
- atmospheric pressure: no specific requirements.

Transportation. Gloves packed in shipping box may be transported by all types of covered vehicles according to applicable regulations for the carriage of goods by this type of transport, observing the precautions indicated on the shipping box.

Transportation conditions:

- temperature: -30 to 50 °C;
- relative humidity: up to 98% at 35 °C; (non-condensing);
- atmospheric pressure: no specific requirements.

Operation. Devices are used in healthcare facilities, ambulance vehicles, in field and home conditions, to protect the hands of medical staff when using the product for its intended purpose.

Operation conditions:

- temperature: 10 to 35 °C;
- relative humidity: up to 80% at 25 °C;
- atmospheric pressure: no specific requirements.

13. Disposal

The medical device "CLEAN + SAFE® Non-sterile Nitrile Rubber Medical Examination Gloves" after their intended use are classified as epidemiologically hazardous waste (products contaminated with blood and/or other body fluids of patients).

In this case, a medical organization performs only collection and temporary storage of the above waste. Waste is collected to single-use bags or containers with color marking adopted for hazardous medical waste in a room or in places for temporary storage of medical waste before further transportation to a destruction or disposal site.

Transportation of medical waste from medical facilities, and their disposal is provided by specialized organizations licensed to treat epidemiologically hazardous medical waste.

Consumer packing, rejected and expired unused products are classified as epidemiologically safe waste (non-hazardous waste not contacting body fluids of patients), which is close by its composition to municipal solid wastes (MSW).

Non-hazardous waste is collected in reusable containers or disposable bags (excluding bags for hazardous waste). The procedure for collection and storage of non-hazardous waste is determined in accordance with the waste management scheme in each specific organization.

Transportation and disposal are performed considering requirements adopted for a given

INSTRUCTION FOR USE
FOR THE MEDICAL DEVICE

CLEAN+SAFE® Non-Sterile Nitrile Rubber Medical Examination Gloves

tion, in accordance with the sanitary legislation of a particular country for the maintenance of populated areas and treatment of production and consumption waste.

14. Shelf life

The manufacturer guarantees compliance of all parameters and characteristics of the gloves throughout the shelf life: 5 years, subject to the packaging integrity and compliance with storage, transportation and operation conditions, in accordance with the instructions.

15. List of international standards that the medical device complies with

93/42/EEC, EN 455-1, EN 455-2, EN 455-3, EN 455-4, EN 1041, ISO 9001, EN ISO 10993-1, EN ISO 10993-5, EN ISO 10993-10, ISO 11193-1 (Except "Thickness" and "Force at break" parameters), EN ISO 14971, EN ISO 13485, EN ISO 15223-1, EN 16523-1, EN ISO 21171, ASTM D412-16, ASTM D573-04, ASTM D3577-09e1, ASTM D3767-03, ASTM D5151-06, ASTM D6124-06, ASTM D6319-10 (In the part of "Thickness" parameter).

Reference standards are used in the latest edition: the organization maintains and updates the list of reference standards as a record of its own management system in accordance with ISO 13485.

INSTRUCTION FOR USE
FOR THE MEDICAL DEVICE
CLEAN+SAFE® Non-Sterile Nitrile Rubber Medical Examination Gloves

This document contains 8 pages
Peng Wang 王鹏 (signature)
General Manager,
Lynmed Medical Technology (Beijing) Co., Ltd.

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