“APPROVED”

Ontex BVBA

Group Sustainability & Scientific Affairs Manager

(position)

hoot

ame)

ture)

2017

rals)

B 00141573


B 00141573

APOSTILLE (Convention de La Haye du 5 octobre 1961)	
1. Land/Pays/Land :	BELGIË – BELGIQUE – BELGIEN.
2. Deze openbare akte is ondertekend door : Le présent acte public a été signé par : Diese öffentliche Urkunde ist unterschrieben von:	Van Den Bossche, Henri
3. Handelend in hoedanigheid van : Agissant en qualité de : In seiner/ihrer Eigenschaft als:	Notaris/Notaire/Notar
4. Is voorzien van het zegel van Est revêtu du sceau de Sie ist versehen mit dem Siegel des/der:	Notaris/Notaire/Notar Buggenhout
Voor echt verklaard / Attesté / Bestätigt	
5. Te Brussel/A Bruxelles/In Brüssel	6. Op/Le/Am : 21/11/2017
7. Door FOD Buitenlandse Zaken, Buitenlandse Handel en Ontwikkelingssamenwerking Par le SPF Affaires étrangères, Commerce extérieur et Coopération au Développement Durch FOD Auswärtige Angelegenheiten, Außenhandel und Entwicklungszusammenarbeit	
8. Onder Nr. / Sous le n°/ Unter Nr. :	9805803921891744
9. Stempel/Sceau/Stempel :	10. Ondertekening/Signature/ Unterschrift : Jan Van de Velde
Prijs/Prix/ Preis : 20 EUR	
Deze Apostille waarborgt de authenticiteit van de inhoud van het document: niet. Cette Apostille ne garantit pas l'authenticité du contenu du document. Diese Apostille dient nicht dem Beweis der Authentizität des Inhalts des Dokuments. Deze Apostille controleren? - Vérifier cette Apostille? - Diese Apostille überprüfen?: http://legalweb.diplomatie.be	

I, the undersigned, Henri Van den Bossche, notary holding residence at Buggenhout, confirm that the signature on this document is the signature of Mr. Bart Waterschoot and is well known by me.

This confirmation of the signature only concerns its authenticity and does not imply a certification by the notary of the capacity or the powers of the aforementioned signatory.

"APPROVED"

Ontex BVBA

Group Sustainability & Scientific Affairs Manager

(position)

Bart Waterschoot


(name)

(signature)

« 30 » October 2017

«day» month (numerals)

Exploitation documentation

**DIRECTIONS FOR USE
OF THE MEDICAL PRODUCT**

**Urological pads «Optio» in the following size
(see Paragraph 1 on page 2)**

logical pads with


**Ontex BVBA», Genthof 5, 9255 Buggenhout, Belgium
2017**

an den Bossche, notary holding residence at Buggenhout, confirm that the
is the signature of Mr. Bart Waterschoot and is well known by me.

ature only concerns its authenticity and does not imply a certification by the
powers of the aforementioned signatory.

2017 at Buggenhout.

1. Name of the medical product (hereinafter "MP")

Urological pads «Optio» in the following sizes (hereinafter "Pads, Products"):

- Urological pads «Optio» Ultra for women;
- Urological pads «Optio» Mini for women;
- Urological pads «Optio» Normal for women;
- Urological pads «Optio» Extra for women;
- Urological pads «Optio» Super for women;
- Urological pads «Optio» Extra Plus for women;
- Urological pads «Optio» Maxi for women;
- Urological pads «Optio» Super plus for women;
- Urological pads «Optio» Extra for men;
- Urological pads «Optio» Super for men;
- Urological pads «Optio» Level 1 for men;
- Urological pads «Optio» Level 2 for men.

2. Information on the medical product

Name of the operating company: Ontex BVBA

Registered address of the company: «Ontex BVBA», Genthof 5, 9255 Buggenhout, Belgium.

Name of the manufacturing companies:

1. «Ontex BVBA», Belgium.
2. «Ontex BVBA Eeklo», Belgium.
3. «Ontex Sante France SAS», France.
4. «Ontex CZ s.r.o.», Czech Republic.

Registered address of the company:

1. Genthof 5, 9255 Buggenhout, Belgium.
2. Korte Moeie 53, 9900 Eeklo, Belgium.
3. Quai du Rivage, 62119 Dourges, France.
4. Turnov, Vesecko cp. 491, PSC 51101, Czech Republic.

Company profile of the manufacturer:

Production, sale and distribution of hygienic and medical products both under own brands and retail brands.

Authorised representative of the manufacturer: LLC «Ontex RU», 105064, Moscow, Zemlyanoy val St, 9, www.ontex.ru, Russia.

3. Description of the product

Urological pads «Optio» in the following size (see Paragraph 1 on page 2) - single use. Anatomic shape of the pads makes them easy to use and avoids the feeling of their constant presence.

3.1 Description of the product structure


Picture urological pads with elastic rubber band

Picture urological pads



Picture urological pads for men


Products contain the following layers:

Schematic image – a sample of urological pads.

Schematic image – a sample of urological pads with elastic rubber band

Schematic image – a sample of urological pads for men

"Urological pads "Optio" for women, size Ultra, "Urological pads "Optio" for women, size Ultra, "Urological pads "Optio" for men, size Extra, Urological pads "Optio" for men, the size of the pads have in its composition elastic gum.

Products, which include elastic bands:

3.2 Designation of the product given by the producer

Products protect clothes and bedding from wetting and contamination by biological secretions of patients.

3.4 Contra indications for use of the products

3.5 Side effects of the products

3.6 Scope of product application

3.7 Use instructions of the product

- wash hands before and after the change of the pads;
- change the pad every 3-4 hours or more often, if the pad soaks quicker;
- do not tear the used pad, don't wash it, don't throw to the sewage line;
- don't cut or change the form of the product because it may influence its effectiveness;
- wash genitals or clean the intimate zone with the help of special napkins for intimate hygiene during each pads change.

3.8 Methods of product application

3.9 Completeness

Products are supplied completed in accordance with the engineering and technical documents and order requirements. Completeness of the product delivery should comply with the information specified in the Table No 1.

Table No 1. Completeness of the delivery

No	Name	Number (pcs.)
1.	Urological pads «Optio» in the following size (see Paragraph 1 on page 2)	from 8 to 30

4. Technical characteristics of the product

4.1 Size of the products

The size of the products should comply with the sizes specified in the Table No 2, 2.1.

Table No 2. Size of the products

Name of the sizes of the products	Sizes	
	Main size of the product*, mm	Absorbing layer*, mm
Ultra for women	190x75	160x55
Mini for women	227x92	205x73
Normal for women	242x92	220x73
Extra for women	257x92	235x73
Super for women	280x100	255x80
Extra Plus for women	350x125	320x100
Maxi for women	380x125	350x100
Super Plus for women	405x155	370x125
Extra for men	240x80	190x80
Super for men	410x80	340x80
Level 1 for men	305x185	275x145
Level 2 for men	305x185	275x145

*acceptable variations on target indicators specified in the table are of +-10%

Table No 2.1 Size of the products

Nane MP	Elastic headband		Protective tape	
	Length *, mm	Width *, mm	Length *, mm	Width *, mm
Urological pads «Optio» Ultra for women	-	-	160	48
Urological pads «Optio» Mini for women	-	-	175	48
Urological pads «Optio» Normal for women	130	10	219	48
Urological pads «Optio» Extra for women	160	10	267	48
Urological pads «Optio» Super for women	200	10	318	48
Urological pads «Optio» Extra Plus for women	190	10	324	48

Urological pads «Optio» Maxi for women	215	10	343	
Urological pads «Optio» Super plus for women	230	10	368	
Urological pads «Optio» Extra for men	-	-	126	
Urological pads «Optio» Super for men	-	-	168	42
Urological pads «Optio» Level 1 for men	190	6	190	48
Urological pads «Optio» Level 2 for men	185	6	190	48

*acceptable variations on target indicators specified in the table are of +-5%

4.2 Materials of the products

Materials specified in the Table No 3 are used for manufacturing of the products.

Table No 3. Materials of the product

Layer of the products	Materials of the products
Upper covering layer	Hydrophilic fiber-bonded material of polypropylene under trademark «NW 21113» produced by «INNOWO GmbH» Germany
Distributing layer	Fiber-bonded material of polypropylene under trademark «LOFT SOFT» produced by «LIBELTEX», Belgium
Protective layer	Polyethylene film under trademark «WHI 90» produced by «TRIOPLANEX FRANCE SA», France Hydrophobic fiber-bonded material of polypropylene under trademark «LA-PEE012W-SPU013W» produced by «EXTEN SA», Switzerland.
Absorbing layer	Cellulose bleached without chlorine under trademark «GP Cellulose» produced by «GP Cellulose GmbH», Switzerland. Polymer superabsorbent under trademark «AQUAKEEP D60B» produced by «SUMITOMO SEIKA EUROPA SA», France.
Elastic rubber band	Elastomeric polyurethane thread under trademark "Spandex" produced by «HYOSUNG SPANDEX (JIAXING) CO. LTD», the Republic of Korea.
Glue	Hot-melt glue «Bostik» on the base of styrole copolymer under trademark «H20028» produced by «Bostik Nederland B.V.», The Netherlands.

4.3 Weight and functional performance of the product

Weight and functional performance of the products are specified in the Table No 4.

Table No 4. Weight and functional performance of the products

Name of the sizes of the products	Specification					
	Weight of the products*, g	Minimal quantity of the absorber, g	Absorptive capacity, not less than, g	Water holding capacity, not less, g	Rewetting, no more, g	Absorption rate, sec
Ultra for women	5	1	100	24	5	<70

Mini for women	9.5	2	170	48	5	<70
Normal for women	10.5	2.5	200	55	5	<90
Extra for women	11.5	2.5	270	75	5	<100
Super for women	15	3	330	97	5	<100
Extra Plus for women	22.5	5	530	143	5.5	<100
Maxi for women	34.5	6	800	210	6	<110
Super Plus for women	34.5	6	559	152	6	<100
Extra for men	19	4	560	153	6	<100
Super for men	31	4	700	205	7	<105
Level 1 for men	19	4	560	154	6	<70
Level 2 for men	31	4	700	206	7	<75

*acceptable variations on target indicators specified in the table are of $\pm 10\%$

5. Classification of the medical product

"Urological pads «Optio» in the following size» (see Paragraph 1 on page 2)" belong to the 1 class according to Annex 2 of the Order No 4n "On approval of nomenclature classification of the medical products" as of June 06, 2012 of the Ministry of Health of the Russian Federation with subsequent amendments.

Type of contact with human body: brief contact with intact skin and mucous membranes of the human body.

6. Marking of the products

Ready-made products are marked on the package or on the labels; place of marking application is specified in the design documents.

Marking of the separate tare contains:

- The name of the manufacturing company and (or) its trademark;
- The address of the manufacturing company, including the country of origin;
- The name of the products;
- The size of the product;
- The number of products in the package;
- The materials used in the manufacture;
- The date of manufacturing of the products specified in the batch number on the side of the packaging in the format: LOT XX: hour/day/year/min/month;
- The expiration date;
- The expiry date (see batch number: EXP month/year);
- Method of application (drawings and/or words);
- The usage instruction;
- The designation of the product;
- Contra indications for use of the products;
- Storage conditions;

- The information on disposal of the products;
- The bar code of the product (if available);
- The information about registration products;
- The information on the single use of the products;
- The name of the importing company on the territory of the Russian Federation;
- The address of the importing company on the territory of the Russian Federation;

The of signs and pictograms on the label:

- "Do not reuse"

- "Protect from moisture"

- Disposal: do not dispose into the sewers

- The Loop

PAP

- The Loop

- the level of absorbance

It is allowed to mark other data including information about consumer properties and the use of the products.

Data are printed or embossed that ensure their readability and preservation (stability) during the storage, transportation, handling and intended usage of the products.

Marking of the transportation tare with the following information:

- the name of the products;
- the article of the products;
- the number of the products in the transportation tare;
- the name and the address of the manufacturing company;
- the name and the address of the company authorized by the manufacturer, of the importer;
- the bar code of the transportation tare;
- the date of manufacturing of the products (month, year);
- the expiry date;
- conditions of storage, transportation, recycling possibility;
- icons and symbols according to the legislation requirements of the consignee country.

7. Packaging

Package should provide the protection of the products from the contamination, mechanical damage, weather impact in transportation and storage as well as the most complete usage of vehicle capacity (tonnage) and easy performance of loading-unloading operations.

Boxes of corrugated cardboard and other packing materials ensuring the protection of the products may be used as transportation tare.

It is allowed to use the other packing materials including those that are manufactured by the manufacturing company which provide the necessary durability and protection of the products in transportation and storage.

8. Transportation and storage

Products are transported by all kinds of vehicles in accordance with the cargo transportation rules operating on the transport of given type at a temperature of from -50 °C to +50 °C a relative humidity of 80%. Material is not hazardous at any kind of transportation.

Ready-made products should be stored in sheltered storage facilities. Storage conditions should exclude the impact of moisture and aggressive environment on the products.

To preserve the quality product should be preserved in original package at the room temperature from +5 °C to +25 °C and relative humidity no more than 70 % in childproof and moisture protected sites. Don't expose the products to direct sunlight.

The temperature of using from +10 °C up to +45 °C.

Shelf life from date of manufacture - 5 years, the size of the Ultra shelf life - 3 years from date of manufacture.

9. Safety requirements

Work with products does not require special safety measures.

Structure of the products does not contain materials which are hazardous for the health of the person during the operation.

Products are explosive safe, inflammable when exposed to the open fire without explosions.

In case of fire materials of the products and package should be extinguished with foam, extinguishing powder, carbon dioxide, water, sand.

10. Requirements for environment protection

Requirements for environment protection during the application of the product

Products do not release toxic substances to the environment and do not impact (at direct contact) the health of the person and the environment during the transportation, storage and operation. Products are explosive safe, inflammable when exposed to the open fire without explosions. Special environment protective measures from harmful influence are not required during the operation.

Requirements for environment protection and disposal of the products

The main possible harmful impact on the environment is the contamination of the atmosphere, soils and waters of the inhabited areas as a result of:

- unorganized waste landfilling on the territory of manufacturing company or outside its territory;
- random landfilling in unspecified places.

Products and materials used for their manufacture do not impose risk for the health of people and environment both during the operation and after its completion.

11. Maintenance and repair

Not applicable. Products are not objects of maintenance and repair.

12. Disposal

It is allowed to dispose waste and materials during the manufacture by contract with company which has the required license. (To dispose with the help of the method approved by the local legislation).
The product can be disposed as household waste after its usage. It is not allowed to throw the used products into the toilet.
Products are classified under the waste class A – epidemiological waste safe according SanPiN 2.1.7290-10.

13. Warranty liabilities

Manufacturer assures the compliance of the products with the requirements of the product specifications and guarantees the replacement of the products in case of nonconformity.
Warranty period for the product is 5 years from the moment of manufacture/packing. Date of manufacture is specified with news ticker on the side of the package.
After the expiration date this medical product should be disposed.

14. Information about reclamations

If defects of the products appear during warranty period as well as at primary acceptance of the products, the owner of the product should send a reclamation to the address of the selling organization – request for the replacement according to the form No 1.

Form No 1.

The address where the representative of the company carrying out warranty liabilities should arrive, phone number;		
Date	Short description of the defect	Date of the sending of the reclamation

Send the reclamations to the address:
LLC «Ontex RU», 105064, Moscow, Zemlyanoy val St, 9, Russia.
Postal address: LLC «Ontex RU», 105064, Moscow, Zemlyanoy val St, 9, Russia.
Tel.: 8 495 646 92 21; website: www.ontex.ru

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