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CERTIFICATE

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China Council for the Promotion of International Trade
China Chamber of International Commerce

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THIS IS TO CERTIFY THAT: the seal of ASSURE
TECH. (HANGZHOU) CO., LTD. on the annexed DOCUMENT is
genuine.

China Council for the Promotion
of International Trade

授权签字: 
Authorized
Signature: Guo Xiaofu

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证明书查询网址 Website for verifying the certificate: <http://www.rzccpit.com/validate.html>

INSTRUCTIONS FOR USE ON A MEDICAL DEVICE
FOR IN VITRO DIAGNOSTICS

**Rapid test for detecting drugs and their metabolites in urine by
immunochromatographic method ("Test system OLLTEST® 1")**

Instructions for use VER 05.02.2023

2023

Content

1. Description of medical device	4
1.1. Name of medical device	4
1.2. Intended use	10
1.3. Scope and potential users.....	12
1.4. Indications and contraindications of use.....	13
1.5. Diagnostic value	15
1.6. Description of analytes	15
1.6.1. Amphetamine (AMP)	15
1.6.2. Barbiturates (BAR)	16
1.6.4. Cocaine (COC).....	17
1.6.5. Methylenedioxymethamphetamine (MDMA)	17
1.6.7. Methamphetamine (MET).....	18
1.6.8. Methadone (MTD).....	18
1.6.9. Opiates (OPI)	19
1.6.10. Phencyclidine (PCP)	19
1.6.11. Marijuana (THC)	19
2. Principle of the method	20
3. Medical device specifications.....	21
3.1. Technical characteristics of medical device	21
3.4. Analytical specificity.....	23
3.5 Cross reactivity	23
3.6 Reproducibility	23
4. Precautions	23
5 Equipment and materials that are not supplied in sets of test systems, but are necessary for analysis.....	24
6 Samples	24
7 Method of analysis	25
8 Interpretation of results	27
8.1 Test results for drug detection.....	27
8.2 Test results to establish the fact of falsification	27
9 Quality control	28

10 Limitations of the use of the test system	28
11 Operating conditions, storage and transportation	29
11.1 Operating conditions	29
11.2 Storage conditions.....	29
11.3 Transport conditions	29
12 Procedure for disposal and destruction of the test system	29
13. Information on the developer and manufacturer of the medical device	30
14. Authorized representative of the manufacturer in the Russian Federation.....	30
15 Warranty	30
16 Bibliography	30
Appendix A. Test composition.....	32
Appendix B. Sample correlation	40
Appendix C. The report on Analytical sensitivity	43
Appendix D. The report on Analytical specificity	45
Appendix E. Non cross-reactive compounds.....	51

1. Description of medical device

1.1. Name of medical device

Express test for detecting drugs and their metabolites in urine by immunochromatographic method ("Test system OLLTEST[®] 1")

Execution options:

1. Rapid test for the detection of 11 drugs and their metabolites in the urine by the immunochromatographic method with a strip to control the falsification of urine. Contents:

- A disposable test panel for the determination of 11 types of drugs and their metabolites: amphetamine, barbiturates, benzodiazepines, cocaine, synthetic cathinones, marijuana, methamphetamine, methadone, methylenedioxymethamphetamine, phencyclidine, opiates with a strip to control urine falsification according to pH, SG, creatinine - 25 pcs.

- Color table for establishing the fact of falsification - 1 pc.

- Instructions for use - 1 pc.

2. Express test for the detection of 11 drugs and their metabolites in the urine by the immunochromatographic method

Contents:

- Disposable test panel for the determination of 11 types of drugs and their metabolites: amphetamine, barbiturates, benzodiazepines, cocaine, synthetic cathinones, marijuana, methamphetamine, methadone, methylenedioxymethamphetamine, phencyclidine, opiates - 25 pcs.

- Instructions for use - 1 pc.

3. Express test for the detection of 6 drugs and their metabolites in the urine by the immunochromatographic method with a strip to control the falsification of urine. Contents:

- A disposable test panel for the determination of 6 types of drugs and their metabolites: amphetamine, cocaine, methamphetamine, opiates, marijuana,

methadone with a strip to control urine falsification according to the parameters creatinine, pH, specific gravity - 25 pcs.

- Color table for establishing the fact of falsification - 1 pc.
- Instructions for use - 1 pc.

4. Express test for detection of 6 drugs and their metabolites in urine by immunochromatographic method.

Contents:

- Disposable test panel for the determination of 6 types of drugs and their metabolites: amphetamine, cocaine, methamphetamine, opiates, marijuana, methadone - 25 pcs.

- Instructions for use - 1 pc.

5. Express test for the detection of 6 drugs and their metabolites in the urine by the immunochromatographic method with a strip to control the falsification of urine. Contents:

- A disposable test panel for the determination of 6 types of drugs and their metabolites: amphetamine, barbiturates, benzodiazepines, opiates, marijuana, methadone with a strip to control urine falsification according to the parameters creatinine, pH, specific gravity - 25 pcs.

- Color table for establishing the fact of falsification - 1 pc.
- Instructions for use - 1 pc.

6. Express test for detecting 6 drugs and their metabolites in urine by immunochromatographic method

Contents:

- Disposable test panel for the determination of 6 types of drugs and their metabolites: amphetamine, barbiturates, benzodiazepines, opiates, marijuana, methadone - 25 pcs.

- Instructions for use - 1 pc.

7. Express test for the detection of 2 drugs and their metabolites in the urine by the immunochromatographic method with a strip to control the falsification of urine. Contents:

- A disposable test panel for the determination of 2 types of drugs substances and their metabolites: barbiturates, benzodiazepines with a strip to control urine falsification by the parameters creatinine, pH, specific gravity - 25 pcs.

- Color table for establishing the fact of falsification - 1 pc.

- Instructions for use - 1 pc.

8. Express test for the detection of 2 drugs and their metabolites in the urine by the immunochromatographic method

Contents:

- Disposable test panel for the determination of 2 types of drugs and their metabolites: barbiturates, benzodiazepines - 25 pcs.

- Instructions for use - 1 pc.

9. Express test for the detection of 3 drugs and their metabolites in the urine by the immunochromatographic method with a strip to control the falsification of urine. Contents:

- A disposable test panel for the determination of 3 types of drugs and their metabolites: barbiturates, synthetic cathinones, phencyclidine with a strip to control urine falsification according to the parameters creatinine, pH, specific gravity - 25 pcs.

- Color table for establishing the fact of falsification - 1 pc.

- Instructions for use - 1 pc.

10. Express test for detection of 3 drugs and their metabolites in urine by immunochromatographic method

Contents:

- Disposable test panel for the determination of 3 types of drugs and their metabolites: barbiturates, synthetic cathinones, phencyclidine - 25 pcs.

- Instructions for use - 1 pc.

11. Express test for the detection of 2 drugs and their metabolites in the urine by the immunochromatographic method with a strip to control the falsification of urine.

Contents:

- A disposable test panel for the determination of 2 types of drugs and their metabolites: synthetic cathinones, phencyclidine with a strip to control urine falsification in terms of creatinine, pH, specific gravity - 25 pcs.

- Color table for establishing the fact of falsification - 1 pc.

- Instructions for use - 1 pc

12. Express test for detection of 2 drugs and their metabolites in urine by immunochromatographic method.

Contents:

- Disposable test panel for the determination of 2 types of drugs and their metabolites: synthetic cathinones, phencyclidine - 25 pcs.

- Instructions for use - 1 pc.

13. Express test for the detection of 2 drugs and their metabolites in the urine by the immunochromatographic method with a strip to control the falsification of urine.

Contents:

- A disposable test panel for the determination of 2 types of drugs and their metabolites: cocaine, phencyclidine with a strip to control the falsification of urine in terms of creatinine, pH, specific gravity - 25 pcs.

- Color table for establishing the fact of falsification - 1 pc.

- Instructions for use - 1 pc.

14. Express test for detection of 2 drugs and their metabolites in urine by immunochromatographic method

Contents:

- Disposable test panel for the determination of 2 types of drugs and their metabolites: cocaine, phencyclidine - 25 pcs.

- Instructions for use - 1 pc.

15. Express test for the detection of 2 drugs and their metabolites in the urine by the immunochromatographic method with a strip to control the falsification of urine.

Contents:

- A disposable test panel for the determination of 2 types of drugs and their metabolites: barbiturates, phencyclidine with a strip to control urine falsification according to the parameters creatinine, pH, specific gravity - 25 pcs.
- Color table for establishing the fact of falsification - 1 pc.
- Instructions for use - 1pc.

16. Rapid test for the detection of 2 drugs and their metabolites in urine by immunochromatographic method

Contents:

- Disposable test panel for the determination of 2 types of drugs and their metabolites: barbiturates, phencyclidine - 25 pcs.
- Instructions for use - 1pc.

17. Rapid test for the detection of 3 drugs and their metabolites in the urine by the immunochromatographic method with a strip to control the falsification of urine.

Contents:

- A disposable test panel for the determination of 3 types of drugs and their metabolites: methamphetamine, cocaine, methylenedioxymethamphetamine with a strip to control the falsification of urine according to the parameters creatinine, pH, specific gravity - 25 pcs.
- Color table for establishing the fact of falsification - 1 pc.
- Instructions for use - 1pc.

18. Rapid test for detection of 3 drugs and their metabolites in urine by immunochromatographic method.

Contents:

- Disposable test panel for the determination of 3 types of drugs and their metabolites: methamphetamine, cocaine, methylenedioxymethamphetamine - 25 pcs.

- Instructions for use - 1 pc.

19. Rapid test for the detection of 3 drugs and their metabolites in the urine by the immunochromatographic method with a strip to control the falsification of urine.

Contents:

- A disposable test panel for the determination of 3 types of drugs and their metabolites: barbiturates, methadone, methylenedioxymethamphetamine with a strip to control urine falsification by the parameters creatinine, pH, specific gravity - 25 pcs.

- Color table for establishing the fact of falsification - 1 pc.

- Instructions for use - 1 pc.

20. Rapid test for the detection of 3 drugs and their metabolites in the urine by the immunochromatographic method.

Contents:

- Disposable test panel for the determination of 3 types of drugs and their metabolites: barbiturates, methadone, methylenedioxymethamphetamine - 25 pcs.

- Instructions for use - 1 pc.

21. Rapid test for the detection of 3 drugs and their metabolites in the urine by the immunochromatographic method with a strip to control the falsification of urine.

Contents:

- Disposable test panel for the determination of 3 types of drugs and their metabolites: amphetamine, opiates, marijuana with a strip to control the falsification of urine by the parameters of creatinine, pH, specific gravity - 25 pcs.
- Color table for establishing the fact of falsification - 1 pc.
- Instructions for use - 1 pc.

22. Rapid test for the detection of 3 drugs and their metabolites in the urine by immunochromatographic method.

Contents:

- Disposable test panel for the determination of 3 types of drugs and their metabolites: amphetamine, opiates, marijuana - 25 pcs.
- Instructions for use - 1 pc.

1.2. Intended use

Rapid test for detecting drugs and their metabolites in urine by the immunochromatographic method ("Test-system OLLTEST® 1") (hereinafter referred to as the test system) is intended for the qualitative simultaneous determination of 2 to 11 drugs and their metabolites.

Functional purpose - drug screening. Screening analyzes of various drugs in urine range from simple immunochromatographic tests to complex analytical procedures. The speed and sensitivity of immunochromatographic tests ensure their wide use in the analysis of urine for the presence of various drugs. The test system is a device for carrying out immunochromatographic analysis for the qualitative

detection of substances, without the need for additional devices. The calculated cut-off concentrations and the calibrator for the determined drugs are shown in Table 1.

Table 1 - Calculated cut-off concentrations and calibrator

Parameter	Name	Calibrator	Cut-off (ng/ml)
AMP	Amphetamine	D- Amphetamine	300
BAR	Barbiturates	Secobarbital	2000
BZO	Benzodiazepines	Oxazepam	100
COC	Cocaine	Benzoyllecgonine	100
MDMA	Methylenedioxymethamphetamine	3,4-methylenedioxymethamphetamine	500
α-PVP	Synthetic cathinones	Alpha-pirrolidinovalerophenon	500
MET	Methamphetamine	D-Methamphetamine	500
MTD	Methadone	Methadone	200
OPI	Opiates	Morphine	2000
PCP	Phencyclidine	Phencyclidine	25
THC	Marijuana	11 nor-Δ ⁹ -THC-9-COOH	25

The test system is an immunochromatographic express analysis for the qualitative preliminary detection of any combinations of drugs in human urine sample. Each test uses monoclonal antibodies to selectively detect elevated urinary drug levels. The use of a test system for the detection of drugs and their metabolites in urine gives good results if the content of drugs in urine reaches a certain cut-off concentration. This test is recommended by the United States Substance Abuse and Mental Health Services Administration (SAMHSA).

Specific pathology, condition or risk factor:

The test system is used for medical examination of persons with suspicion of drug use, pre-trip and post-trip examinations, when persons are admitted to the hospital and other cases associated with increased danger.

Analyzed sample type:

- freshly collected urine stored at a temperature of 18 to 25°C for up to 4 hours;
- urine stored in a refrigerator at a temperature of 2 to 8°C for up to 48 hours;
- frozen urine, stored at temperatures from minus 20 to minus 30 ° C for up to 30 days.

The results of the analysis are not affected by population and demographic aspects of the use of the test system.

Options 1, 3, 5, 7, 9, 11, 13, 15, 17, 19, 21 are equipped with a test strip for detecting sample falsification.

Falsification is tampering with a urine sample to alter test results. The use of impurities can lead to false negative drug test results, either interfering with the screening test and / or destroying drugs present in the urine. Dilution can also be used to produce false negative drug test results.

One of the best ways to check for tampering or dilution is to measure some of the characteristics of the urine, such as creatinine, pH, and urine specific gravity.

Creatinine (CRE): Sample dilution tests. Creatinine is a byproduct of creatine and is an amino acid found in muscle tissue and found in urine. A person may try to interfere with a drug test by drinking excessive amounts of water or diuretics such as herbal teas to flush out the body. Creatinine and specific gravity are two methods of testing for dilution and washing, which are the most common mechanisms used to bypass drug testing. Low creatinine and specific gravity levels may indicate diluted urine. A lack of creatinine (<5 mg / dL) indicates that the sample does not match human urine.

pH: tests for the presence of acidic or alkaline impurities in the urine. Normal pH levels should be between 4.0 and 9.0. Values outside this range may indicate that the sample has been falsified.

Specific gravity: tests for sample dilution. The normal range is 1.003 to 1.030. Values outside this range may be the result of sample dilution or tampering.

1.3. Scope and potential users

Application: clinical laboratory diagnostics. The test system is intended for use in a laboratory of a medical institution for carrying out an express analysis to identify drugs and their metabolites in human urine.

Potential user: a laboratory specialist with a higher medical education, as well as a laboratory assistant with a secondary specialized medical education.

For in vitro diagnostic use only.

The test system is intended for single use.

1.4. Indications and contraindications of use

Indications: Identification of drugs and their metabolites in human urine.

Option 1 is designed to detect 11 types of drugs and their metabolites: amphetamine, barbiturates, benzodiazepines, cocaine, synthetic cathinones, marijuana, methamphetamine, methadone, methylenedioxymethamphetamine, phencyclidine, opiates, as well as to control the falsification of urine by the parameters of creatinine, pH, SG.

Option 2 is designed to detect 11 types of drugs and their metabolites: amphetamine, barbiturates, benzodiazepines, cocaine, synthetic cathinones, marijuana, methamphetamine, methadone, methylenedioxymethamphetamine, phencyclidine, opiates.

Option 3 is intended for the determination of 6 types of drugs and their metabolites: amphetamine, cocaine, methamphetamine, opiates, marijuana, methadone, as well as for the control of urine falsification by the parameters creatinine, pH, specific gravity

Option 4 is designed to determine 6 types of drugs and their metabolites: amphetamine, cocaine, methamphetamine, opiates, marijuana, methadone.

Option 5 is designed to determine 6 types of drugs and their metabolites: amphetamine, barbiturates, benzodiazepines, opiates, marijuana, methadone, as well as to control the falsification of urine by the parameters creatinine, pH, specific gravity.

Option 6 is designed to determine 6 types of drugs and their metabolites: amphetamine, barbiturates, benzodiazepines, opiates, marijuana, methadone.

Option 7 is intended for the determination of 2 types of drugs and their metabolites: barbiturates, benzodiazepines, as well as for the control of urine falsification by the parameters creatinine, pH, specific gravity.

Option 8 is intended for the determination of 2 types of drugs and their metabolites: barbiturates, benzodiazepines.

Option 9 is intended for the determination of 3 types of drugs and their metabolites: barbiturates, synthetic cathinones, phencyclidine, as well as for the control of urine falsification in terms of creatinine, pH, specific gravity.

Option 10 is intended for the determination of 3 types of drugs and their metabolites: barbiturates, synthetic cathinones, phencyclidine.

Option 11 is designed to determine 2 types of drugs and their metabolites: synthetic cathinones, phencyclidine, as well as to control the falsification of urine by the parameters creatinine, pH, specific gravity.

Option 12 is designed to determine 2 types of drugs and their metabolites: synthetic cathinones, phencyclidine.

Option 13 is designed to determine 2 types of drugs and their metabolites: cocaine, phencyclidine, as well as to control the falsification of urine by the parameters creatinine, pH, specific gravity.

Option 14 is designed to determine 2 types of drugs and their metabolites: cocaine, phencyclidine.

Option 15 is intended for the determination of 2 types of drugs and their metabolites: barbiturates, phencyclidine, as well as for the control of urine falsification according to the parameters creatinine, pH, specific gravity.

Option 16 is designed to determine 2 types of drugs and their metabolites: barbiturates, phencyclidine.

Option 17 is designed to determine 3 types of drugs and their metabolites: methamphetamine, cocaine, methylenedioxymethamphetamine, as well as to control the falsification of urine according to the parameters creatinine, pH, specific gravity.

Option 18 is designed to determine 3 types of drugs and their metabolites: methamphetamine, cocaine, methylenedioxymethamphetamine.

Option 19 is designed to determine 3 types of drugs and their metabolites: barbiturates, methadone, methylenedioxymethamphetamine, as well as to control the falsification of urine by the parameters of creatinine, pH, specific gravity.

Option 20 is intended for the determination of 3 types of drugs and their metabolites: barbiturates, methadone, methylenedioxymethamphetamine.

Option 21 is designed to determine 3 types of drugs and their metabolites: amphetamine, opiates, marijuana, as well as to control the falsification of urine in terms of creatinine, pH, specific gravity.

Option 22 is designed to determine 3 types of drugs and their metabolites: amphetamine, opiates, marijuana.

Contraindications: It has no contraindications within the scope of the intended use established by the manufacturer.

1.5. Diagnostic value

The test system is used to detect narcotic substances and their metabolites in urine and can be used in studies of schoolchildren and students (order of the Ministry of Health of the Russian Federation 581n dated 06.10.2014), employees employed in public service (decree of the government of the Russian Federation No. 565 dated 04.07.2013, Order of the Ministry of Health of the Russian Federation 988n of December 22, 2016, order of the Ministry of Health of the Russian Federation No. order of the Ministry of Health of the Russian Federation No. 266 of 04/30/2019), employees of hazardous and hazardous industries, when conducting pre-trip and pre-shift inspections on transport (order of the Ministry of Health of the Russian Federation No. 835n of 12/15/2014, order Ministry of Health of the Russian Federation 39n of 01/29/2016), during medical examinations of athletes, during medical examination (order of the Ministry of Health of the Russian Federation No. 933n of 12/18/2015, order of the Ministry of Health of the Russian Federation No. 344 of 06/15/2015, order of the Ministry of Health of the Russian Federation No. 441 of 06/30/2016) and in other cases where the result of this test is required. This analysis provides only a preliminary analytical result. A more accurate alternative chemical method should be used to confirm the analytical result. Gas chromatography with mass spectrometer (GC / MS) or liquid chromatography with mass spectrometer (LC / MS) are the preferred confirmatory methods. Clinical considerations and professional judgment should apply to any drug abuse test result, especially if preliminary positive results are indicated. Test strips for detecting falsified urine is a semi-quantitative test based on color comparison to determine creatinine, pH, specific gravity of human urine. This test provides only preliminary data. A more accurate alternative chemical method must be used to obtain a confirmed analytical result. Results must be sent to the laboratory for confirmation.

1.6. Description of analytes

1.6.1. Amphetamine (AMP)

Amphetamines belong to the category of potent adrenomimetic drugs used for therapeutic purposes. By the chemical composition, amphetamines belong to the human natural catecholamines: adrenaline and norepinephrine. High doses of amphetamine lead to increased excitation of the central nervous system and cause a state of euphoria and increased alertness, reduce appetite and cause a feeling of increased energy and strength. The cardiovascular response to amphetamines includes increased blood pressure and cardiac arrhythmias. More acute reactions are

the appearance of anxiety disorders, paranoia, hallucinations and psychotic behavior. The effect of amphetamines lasts 2-4 hours after use, the half-life of the drug from the body is 4-24 hours. About 30% of amphetamines are excreted in the urine in an unchanged form with the remains of hydroxylation products and deamidated derivatives.

1.6.2. Barbiturates (BAR)

Barbiturates weaken the activity of the central nervous system. For therapeutic purposes, they are used as sedatives, hypnotics and anticonvulsants. Barbiturates are usually taken orally in the form of capsules or tablets. The effect resembles the state of alcohol intoxication. Long-term use of barbiturates leads to immunological tolerance and physical dependence. Quick-acting barbiturates taken at 400 mg / day for 2–3 months form a clinically significant degree of physical dependence. Symptoms of withdrawal can be quite severe and can be fatal. A very small amount (less than 5%) of barbiturates is excreted in the urine in unchanged form. The detection period for barbiturates in urine is 4-7 days.

1.6.3. Benzodiazepines (BZO)

Benzodiazepines are medications often prescribed for the symptomatic treatment of anxiety and sleep disorders. They act through specific receptors using gamma-aminobutyric acid (GABA). Being safer and more effective, benzodiazepines have replaced barbiturates in the treatment of anxiety disorders and sleep disorders. Benzodiazepines are also used as sedatives before surgical and medical procedures and for the treatment of epilepsy and alcohol withdrawal.

With regular use of benzodiazepines (for example, daily) for several months, the risk of physical dependence increases, especially when taking large doses. Abrupt cessation of admission can lead to symptoms such as trouble sleeping, indigestion, poor health, loss of appetite, sweating, trembling, weakness, anxiety disorders, and altered perceptions. Only a negligible amount of most benzodiazepines (less than

1%) are excreted unchanged in the urine; concentration in the urine from paired drugs. The detection period of benzodiazepines in urine is 3-7 days.

1.6.4. Cocaine (COC)

Cocaine is a potent central nervous system and local anesthetic. Initially, it causes excessive energy and excitement, gradually leading to tremors, increased sensitivity and cramping. In large quantities, cocaine causes fever, immunity, shortness of breath and loss of consciousness.

Often, cocaine is taken on its own through inhalation through the nose, intravenous injection and smoking pure cocaine. It is rapidly excreted in the urine as benzoylegonin (1,2). Benzoylegonin, the main metabolite of cocaine, has a longer biological half-life (5-8 hours) than cocaine (0.5-1.5 hours), and can be detected within 24-48 hours after exposure to cocaine (2).

1.6.5. Methylenedioxymethamphetamine (MDMA)

Methylenedioxymethamphetamine (ecstasy) is a designer drug first synthesized in 1914 by the German pharmaceutical company for the treatment of obesity. Those who took the drug often reported side effects, such as increased muscle tension and sweating. Although methylenedioxymethamphetamine is clearly not a stimulant, like amphetamine it can increase blood pressure and heart rate. Methylenedioxymethamphetamine causes changes in perception, which is manifested in an increase in sensitivity to light, difficulty in focusing and blurred vision in some people who consume the substance. It is believed that the mechanism of action of methylenedioxymethamphetamine is manifested through the release of the neurotransmitter Serotonin. Methylenedioxymethamphetamine can also release dopamine, but it is generally accepted that this is a secondary effect of the drug (Nichols and Oberlender, 1990). The most common effect of methylenedioxymethamphetamine, which occurs in almost all people who took a moderate dose of the drug, was to compress the jaws.

1.6.6. Synthetic cathinones (a-PVP)

A-Pyrrolidinovalerophenone is a recreational psychoactive substance, the main effect of which is psychostimulation. When taking the substance, a person experiences the following effects: euphoria, talkativeness, fussiness, increased concentration, as well as a general pleasant state of health. The intense sensation of pleasure lasts 30-45 minutes, reaching a peak 1.5 hours after ingestion. The main effects disappear after about 4 hours, the general "desirable" experience lasts more than 6-8 hours

1.6.7. Methamphetamine (MET)

Methamphetamine is an addictive stimulant that intensely stimulates certain brain systems. The chemical composition of methamphetamine is very similar to amphetamine, but the effect of methamphetamine on the central nervous system is much stronger. Methamphetamine produced in illegal laboratories is very addictive. The drug is taken orally by injection or inhalation. High doses lead to increased excitation of the central nervous system and cause a state of euphoria and increased alertness, reduce appetite and cause a feeling of increased energy and strength. The cardiovascular response to methamphetamine includes increased blood pressure and cardiac arrhythmias. More acute reactions include the appearance of anxiety disorders, paranoia, hallucinations, psychotic behavior, depression and exhaustion.

The effect of methamphetamine lasts 2-4 hours after use, the half-life of the drug from the body is 9-24 hours. Methamphetamine is excreted in the urine as amphetamine, an oxidized and deaminated derivative. However, 10-20% of methamphetamine is excreted unchanged in the urine. Thus, the presence of the parent compound in the urine indicates the use of methamphetamine. Typically, methamphetamine can be detected in urine within 3-5 days, depending on the pH level of the urine.

1.6.8. Methadone (MTD)

Methadone is a long-acting painkiller, the effect of which lasts from twelve to forty-eight hours. Methadone relieves consumers of the stress associated with finding illegal heroin, the dangers of injections and mood swings caused by most opiates. Long-term use of methadone in large doses can lead to very long withdrawal

symptoms. Compared to other opiates, when they take a withdrawal period from a week to ten days, the withdrawal period for methadone users can be 5-6 weeks.

1.6.9. Opiates (OPI)

Opiates include any drug derived from opium poppy, including natural substances, morphine and codeine, and semi-synthetic drugs such as heroin. An opioid is a more generic name that refers to any drug that acts on opioid receptors.

Opioid analgesics constitute a large group of substances that control pain through inhibition of the central nervous system. Large doses of morphine can cause high levels of tolerance and physical dependence and lead to drug addiction and abuse. Morphine is excreted in a non-metabolized form and is the main metabolic product of codeine and heroin. Morphine can be found in urine within a few days after taking a dose of opiates.

1.6.10. Phencyclidine (PCP)

Phencyclidine, also known as PCP, is a hallucinogen that was first used in the 1950s as an anesthetic for surgical interventions. PCP was withdrawn from use because patients using the substance were raving and hallucinating. Phencyclidine is used in the form of powders, capsules and tablets. The powder can be consumed by inhalation or smoking after mixing with marijuana or herbal impurities. Phencyclidine is most commonly used by inhalation, but can also be used intravenously, intranasally, and orally. After taking small doses, mental and abilities are accelerated, the receiver thinks and acts quickly, and also experiences mood swings from euphoria to depression. Self-mutilation is one of the most dangerous effects of phencyclidine. PCP is excreted in the urine unchanged (from 4% to 19%) and in the form of conjugates (from 25% to 30%), with a half-life of 12 hours.

1.6.11. Marijuana (THC)

THC (Δ^9 - tetrahydrocannabinol) is the main active component of cannabinoids (marijuana).

If you smoke or take orally, it causes a state of euphoria. Consumers have a weakened memory and slow comprehension. They may also experience short-term confusion and fear. Long-term use of large doses can lead to behavioral disorders. The peak effect when smoking marijuana is manifested in 20-30 minutes and lasts 90-120 minutes after smoking one cigarette. An increased level of metabolism in the urine is observed during the exposure time and can be detected within 3-10 days after smoking. The main metabolite excreted in urine is 11-nor- Δ^9 - tetrahydrocannabinol - 9-carboxylic acid (Δ^9 - THC-COOH).

2. Principle of the method

A rapid test for the detection of drugs and their metabolites in urine by the immunochromatographic method works on the principle of a one-step immunochemical assay, in which chemically labeled narcotic substances (drug-protein conjugates) compete with the molecules of the narcotic substance in the urine for binding to a limited number of antibody binding sites. Drug-protein conjugates are pre-applied to the test strip membranes in the test strip area. At one end of each test strip is a pad containing the antibody-colloidal gold conjugate. In the absence of a drug in the urine, a solution of the colored antibody-colloidal gold conjugate rises up the membrane chromatographically using capillary forces and reaches the immobilized drug-protein conjugates in the test strip zone. The stained antibody-colloidal gold conjugate then binds to the drug-protein conjugates. A colored line appears due to the formation of the antibody / drug-protein conjugate complex. Thus, the formation of a colored line in the area of the test strip indicates a negative result. In the presence of a drug in urine, the drug / metabolite antigen competes with the drug-protein conjugate for binding to antibodies in the test strip area. With a sufficient concentration of the narcotic substance, saturation of all antibody binding sites occurs. This prevents the binding of the antibody-colloidal gold conjugate with the drug-protein conjugate and the formation of a colored band in the test area. Thus, the absence of a colored band in the test area indicates a positive result.

A control strip with a different antigen / antibody reaction is added to the immunochromatographic membrane strip in the control area (C) to indicate that the test is performed correctly. This control line should always appear regardless of the presence of the drug or metabolite. If the control line does not appear, discard the test strip.

3. Medical device specifications

3.1. Technical characteristics of medical device

Overall dimensions of the tablet: length (74.0 ± 1.0) mm, width (51.0 ± 0.5) mm, thickness (5.0 ± 0.3) mm. Tablet mass (11.5 ± 0.5) g.

Test-strips for detecting drugs are represented in three different versions:

- For joint detection of three substances;
- For joint detection of two substances;
- For detection of one substance.

Types of test-strips for detecting drugs:

- For the joint detection of COC/MET/THC
- For the joint detection of BAR/MTD
- For the joint detection of AMP/BZO
- For detection of AMP
- For detection of BAR
- For detection of BZO
- For detection of COC
- For detection of OPI
- For detection of THC
- For detection of MDMA
- For detection of MET
- For detection of MTD
- For detection of a-PVP
- For detection of PCP

The structure of the test strip for the detection of narcotic substances and their metabolites is shown in Figure 1:

A - Substrate of the sample zone (fiberglass, brand: Glass Fiber 8964, manufacturer - Assure Tech, China);

B - Substrate of the conjugate application zone (polyester fiber mesh membrane, brand: Polyester 6613, manufacturer - Assure Tech, China);

C - Nitrocellulose membrane (nitrocellulose, brand: CN140, manufacturer - Assure Tech, China);

D - Absorbent plate (filter paper, cellulose, brand: AN3 Filter Paper, manufacturer - Assure Tech, China)

E - Lower substrate pad (polyvinyl chloride, grade PVC-01, manufacturer Anhui beichi, China, brand).

- Width of test strip - (3 ± 1) mm

The reagents used in the test system are listed in Appendix A.


Figure 1 - The structure of the test strip for the detection of narcotic substances and their metabolites

3.2 Comparison (correlation) of samples

Comparison of samples was carried out on urine samples. All samples were confirmed by gas chromatography-mass spectrometry (GC / MS) or high performance liquid chromatography (HPLC). Samples were randomly selected and analyzed using a test system. Samples were classified as positive or negative after 5 minutes. The results of the analysis are given in Appendix B.

3.3 Analytical sensitivity

Analytical sensitivity was determined by adding a control sample of the drug to a negative urine sample at $\pm 50\%$ of the cut-off value, $\pm 25\%$ of the cut-off value and

the cut-off value. Samples of the drug were randomly selected and encoded. The results were confirmed by GC / MS. The analysis of each sample of the drug was performed 10 times in accordance with the instructions for use of the drug. The results of the analysis were considered after 5 minutes after adding the sample.

Most samples containing the drug at a concentration below the cut-off concentration show a negative result, samples containing the drug at a concentration equal to the cut-off concentration are defined as positive or negative. All samples containing the drug in concentrations exceeding the cut-off concentration show a positive result. The results are presented in Appendix C.

3.4. Analytical specificity

Data of analytical specificity is represented in Appendix D.

3.5 Cross reactivity

The aim of the study was to determine the cross-reactivity of a test with substances in urine that does not contain analyte. The following substances did not show cross-reactivity when tested using a test system at a concentration of 100 µg / ml. The results of cross-reactivity study are given in Appendix E.

3.6 Reproducibility

This study was conducted in 2 series per day and lasted 25 days for each group of three parties. It was carried out by three specialists who were not aware of the sample identification system. Each of the three operators analyzed 2 aliquots of urine samples with different concentrations of the narcotic substance (2 series per day). In total, each specialist conducted 50 tests for each concentration. The results are shown in Appendix F.

4. Precautions

Test system is intended only for in vitro diagnostics.

The package, containing test system must be sealed. Do not use the device if the package is spoiled.

An analysis presents a possible risk of infection. When working with biological materials (or materials of biological origin), personal protective equipment, such as gloves and glasses, should be used.

To avoid contamination, use a separate urine collection container for each sample.

In the manufacture of the test system, materials of animal origin are used, namely, fractions of bovine serum albumin (BSA), goat antibodies, monoclonal mouse antibodies. The materials are non-toxic, non-toxic and not in direct contact with the human body. The materials are not a pharmacological substance or drug. The processes of processing, conservation, testing and processing of tissues, cells and substances of animal origin are designed in such a way as to ensure the maximum level of biological safety for users, patients and other people. Tissues, cells and biological substances of animal origin used in the manufacture of the test system were obtained from animals that passed veterinary control and were monitored, taking into account the intended use of these tissues.

5 Equipment and materials that are not supplied in sets of test systems, but are necessary for analysis

- Stopwatch or timer;
- Urine collection container;
- Medical gloves;
- Laboratory centrifuge;
- Refrigerator with temperature from 2 to 4 ° C;
- Freezer with temperature from minus 20 ° C to minus 30 ° C;
- Aseptic wipes (if necessary).

6 Samples

The test system is designed for use with urine samples.

A urine sample should be collected in a clean and dry container. Urine collected at any time of the day can be used.

Urine samples with visible sediments are centrifuged for 10 minutes. at an acceleration of 1500 g (revolutions per minute), filter or settle until a clear supernatant is obtained.

Fresh urine can be stored at room temperature from 18 to 25°C for 4 hours or stored in a refrigerator from 2 to 8°C for 48 hours until the test. For long-term storage, samples can be frozen and stored at temperatures from minus 20 to minus 30 ° C for up to 30 days.

Before analysis, the urine samples that were stored at a temperature of 2 to 8 ° C should be kept at room temperature from 18 to 25 ° C for at least 30 minutes.

Frozen test urine samples should be kept at room temperature from 18 to 25 ° C for at least 60 minutes.

7 Method of analysis

To obtain reliable results, strict adherence to the Instructions for the use of the test system is necessary.

Before use, verify the integrity of the packaging of the test system by visual inspection.

Do not use the product if the packaging is damaged or poorly sealed.

The test system, if stored at a temperature of 2 to 8 ° C, must be maintained at room temperature from 18 to 25 ° C for at least 10 minutes before analysis. The test system, in the case of storage at a temperature of from 8 to 18 ° C, before analysis should be maintained at room temperature from 18 to 25 ° C for at least 5 minutes. If the test system is stored at a temperature of 25 to 30 ° C, analysis can begin immediately.

Open the individual packaging of the test system immediately before the study.

Remove the test device from the sealed packaging and use it as soon as possible.

After opening the individual sealed packaging of the test system, the analysis should be carried out no later than 5 minutes provided that the room temperature is from 18 to 25 ° C and humidity is not more than 80%.

Dip the substrate area to apply the sample on the test strips into the urine sample (Figure 2).

Remove the test system and place it on a flat, horizontal, clean and dry surface.

For falsification studies, visually compare the color of the reaction substrate with the color chart on the color table (Figure 3). Results should be read after 2 minutes. Do not interpret the results later than 5 minutes.


Fig. 2 – Analysis

Parameter	Normal				Abnormal				
Specific gravity	 1,003	 1,005	 1,015	 1,025	 1,000	 1,035			
Creatinine	 100 Mg/dl		 20 Mg/dl		 0 Mg/dl		 10 Mg/dl		
pH	 4,0	 7,0		 9,0	 2,0	 3,0	 10,0	 11,0	 12,0

Figure 3 - Color table for interpreting the results of falsification tests

Interpret test results after 5 minutes from the start of the analysis. Do not interpret the results later than 10 minutes after adding the test sample to the test system.

8 Interpretation of results

8.1 Test results for drug detection

Positive result (Figure 4): in the control area (C), only one color bar appears. There is no visible color bar in the test area (T).

Negative result (Figure 5): two colored stripes appear on the membrane. One lane appears in the control region (C), and the other lane appears in the test region (T).

Invalid result (Figure 6): the control strip is not displayed. The results of any test that did not give a control strip at the indicated reading time should be discarded. Please familiarize yourself with the procedure and repeat the analysis with the new test system. If the problem persists, immediately stop using the test system and contact the authorized representative of the manufacturer in the Russian Federation.

Fig. 4 – positive result

Fig. 5 – negative result

Fig. 6 – invalid

NOTE:

The color intensity in the test region (T) may vary depending on the concentration of analytes present in the sample. Therefore, any shade of color in the test area should be considered negative. Please note that this is only a qualitative test and cannot be used for quantification.

Insufficient sample volume, non-compliance with the methodology or expired test system expiration date are the most probable reasons for the lack of a control line.

8.2 Test results to establish the fact of falsification

To interpret the results, use the color table (Figure 3), which is part of the test system.

NOTE. Test strips for detecting fraud (urine) are designed to identify samples with deviations from the norm by the main parameters: creatinine, pH, specific gravity

Creatinine: Normal creatinine levels are between 20 and 350 mg / dl. In rare cases, with some kidney diseases, its amount may decrease.

Specific Gravity: Elevated levels of protein in the urine can lead to abnormally high specific gravities.

Normal pH levels should range from 4.0 to 9.0.

9 Quality control

The test system provides internal control of the analysis procedure with different antigen / antibody responses in the control region (C) on each strip. This control line should always appear regardless of the presence or absence of a drug or metabolite in the sample. If the control line does not appear, discard the test strip. The presence of this control band in the control region serves to confirm a sufficient volume of the test sample, the correct technique for performing analysis procedures, and to verify that optimal chromatographic conditions are achieved.

It is possible to use external control materials to ensure proper operation of the test system. Samples of control materials are commercially available. When conducting quality control using external control materials, it is necessary to be guided by the analysis methodology of clause 7 of this Instructions for use. If the control values do not fall within the established limits, then the analysis results are invalid.

10 Limitations of the use of the test system

The test system is intended for in vitro laboratory diagnostics and should only be used for the qualitative determination of drugs.

The test system is intended for analysis only using human urine.

The test system is intended for single use only.

A positive result from any of the tests indicates only the presence of the drug / metabolite and does not establish the degree of intoxication.

There is a possibility that a technical or procedural error, accompanying substances, as well as factors not listed, may interfere with the test and lead to false results.

If a drug / metabolite is detected in a urine sample, the analysis will not show the frequency of drug use and the effect of food and drug administration on the result of the analysis.

11 Operating conditions, storage and transportation

11.1 Operating conditions

After opening the package, the test system should be used for 1 hour at room temperature from 18 to 25 ° C and relative humidity not more than 80%.

The test system, if stored at a temperature of 2 to 8 ° C, must be maintained at room temperature from 18 to 25 ° C for at least 10 minutes before analysis.

To obtain reliable results, strict adherence to the Instructions for Use is required.

11.2 Storage conditions

The test system should be stored in the manufacturer's packaging at a temperature of 2 to 30°C and humidity not more than 60% until the expiration date indicated on the package.

Do not freeze. The product is sensitive to moisture and should be used immediately after opening. Any test in improperly sealed packaging should be discarded.

11.3 Transport conditions

Transport the test system by all types of covered transport in accordance with the transportation rules applicable to a particular type of transport. Transport at a temperature of 2 ° C to 30 ° C.

Handle carefully and avoid strong pressure on the packaging during transport and storage.

Products transported in violation of the temperature regime are not subject to use.

12 Procedure for disposal and destruction of the test system

Destruction of the test system is carried out by organizations that have the appropriate license, on specially equipped sites, landfills and premises in accordance with the requirements stipulated by existing Federal laws, and in compliance with the mandatory requirements for environmental protection.

Samples may contain pathogens. When handling biological fluids, follow generally accepted safety rules and follow the rules of your local laboratory.

Test systems are epidemiologically hazardous waste.

Packaging refers to epidemiologically safe waste that is similar in composition to municipal solid waste.

13. Information on the developer and manufacturer of the medical device

Assure Tech. (Hangzhou) Co., Ltd.

Address: Building 4, No.1418-50, Moganshan Road, Gongshu District, Hangzhou, 310011 Zhejiang, P.R. China

phone: +86-571-81022698, fax: +86571-8886920, email: contact@diareagent.com

14. Authorized representative of the manufacturer in the Russian Federation

LLC «Skylab», Address: 125212, Moscow, Leningradskoye sh, № 37, building 1, floor 1, tel: +7(926)176-91-36, website: www.sky-lab.su, TIN 7743301018.

15 Warranty

The manufacturer guarantees that the test system meets the established requirements of the standards and documentation of the manufacturer. Subject to the conditions of transportation, storage and operation, the test system is stable throughout the expiration date.

The warranty shelf life of the test system is 24 months from the date of production.

Claims, wishes regarding the quality of the test system should be sent to the LLC «Skylab», Address: 125212, Moscow, Leningradskoye sh, № 37, building 1, floor 1, tel: +7(926)176-91-36, website: www.sky-lab.su, TIN 7743301018.

16 Bibliography

1. Baselt RC. The distribution of toxic substances in the human body. 2nd Edition Davis: Biomedical Publications; 1982.

2. Hawks RL, Chiang CN, eds. Examination of urine for the presence of narcotic substances. Rockville: Department of Health and Human Services, National Institute on Drug Abuse; 1986.
3. Substance Abuse and Mental Health Services Administration. Mandatory Guidelines for Federal Workplace Drug Testing Programs. 53 Federal Register; 1988.
4. McBay AJ. Drug-analysis technology - pitfalls and problems of drug testing. Clin Chem. 1987 Oct; 33 (11 Suppl): 33B-40B.
5. Gilman AG, Goodman LS, Gilman A, eds. Goodman and Gilman's The Pharmacological Basis of Therapeutics. 6th ed. New York: Macmillan; 1980.

Appendix A. Test composition

Component	Reagent	Amount	Quantity	Manufacturer
Test strip (BAR2000/MTD200)				
Control zone (Internal control components)	Goat anti rabbit pAb	$\cong 0.5 \mu\text{g}$	Buffer solution: Phosphate buffer pH 5,5-9,3	Surko Bio Tech
	Goat anti mouse pAb	$\cong 0.5 \mu\text{g}$	Preservative: Sodium azide below 0,1%	
Test zone (Main reaction components)	Antigen: BAR-BSA MTD-BSA	$\cong 0.3 \mu\text{g}$ $\cong 0.3 \mu\text{g}$	Buffer solution: Phosphate buffer pH 5,5-9,3 Preservative: Sodium azide below 0,1%	Surko Bio Tech
Conjugate Pad	Mouse Anti BAR	$\cong 0.4 \text{ мкг}$	Buffer solution: Phosphate buffer pH 5,5-9,3	Surko Bio Tech
	Mouse Anti MTD	$\cong 0.4 \text{ мкг}$	Preservative: Sodium azide below 0,1%	
Test strip (COC100/MET500/ THC25)				
Control zone (Internal control components)	Goat anti rabbit pAb	$\cong 0.5 \mu\text{g}$	Buffer solution: Phosphate buffer pH 5,5-9,3	Surko Bio Tech
	Goat anti mouse pAb	$\cong 0.5 \mu\text{g}$	Preservative: Sodium azide below 0,1%	
Test zone (Main reaction components)	Antigen: COC-BSA MET-BSA THC-BSA	$\cong 0.3 \mu\text{g}$ $\cong 0.3 \mu\text{g}$ $\cong 0.3 \mu\text{g}$	Buffer solution: Phosphate buffer pH 5,5-9,3 Preservative: Sodium azide below 0,1%	Surko Bio Tech

Conjugate Pad	Mouse Anti COC	$\cong 0.4 \text{ мкг}$	Buffer solution: Phosphate buffer pH 5,5-9,3	Surko Bio Tech
	Mouse Anti MET	$\cong 0.4 \text{ мкг}$	Preservative: Sodium azide below 0,1%	
	Mouse Anti THC	$\cong 0.4 \text{ мкг}$		
Test strip (AMP300/BZO100)				
Control zone (Internal control components)	Goat anti rabbit pAb	$\cong 0.5 \text{ }\mu\text{g}$	Buffer solution: Phosphate buffer pH 5,5-9,3	Surko Bio Tech
	Goat anti mouse pAb	$\cong 0.5 \text{ }\mu\text{g}$	Preservative: Sodium azide below 0,1%	
Test zone (Main reaction components)	Antigen: AMP-BSA BZO-BSA	$\cong 0.3 \text{ }\mu\text{g}$ $\cong 0.3 \text{ }\mu\text{g}$	Buffer solution: Phosphate buffer pH 5,5-9,3 Preservative: Sodium azide below 0,1%	Surko Bio Tech
Conjugate Pad	Mouse Anti AMP	$\cong 0.4 \text{ }\mu\text{g}$	Buffer solution: Phosphate buffer pH 5,5-9,3	Surko Bio Tech
	Mouse Anti BZO	$\cong 0.4 \text{ }\mu\text{g}$	Preservative: Sodium azide below 0,1%	
Test strip AMP300				
Control zone (Internal control components)	Goat anti rabbit pAb	$\cong 0.5 \text{ }\mu\text{g}$	Buffer solution: Phosphate buffer pH 5,5-9,3	Surko Bio Tech
	Goat anti mouse pAb	$\cong 0.5 \text{ }\mu\text{g}$	Preservative: Sodium azide below 0,1%	
Test zone (Main reaction components)	Antigen: AMP-BSA	$\cong 0.3 \text{ }\mu\text{g}$	Buffer solution: Phosphate buffer pH 5,5-9,3	Surko Bio Tech

			Preservative: Sodium azide below 0,1%	
Conjugate Pad	Mouse Anti AMP	$\cong 0.4 \text{ мкг}$	Buffer solution: Phosphate buffer pH 5,5-9,3 Preservative: Sodium azide below 0,1%	Surko Bio Tech
Test strip (BAR2000)				
Control zone (Internal control components)	Goat anti rabbit pAb	$\cong 0.5 \text{ μg}$	Buffer solution: Phosphate buffer pH 5,5-9,3 Preservative: Sodium azide below 0,1%	Surko Bio Tech
	Goat anti mouse pAb	$\cong 0.5 \text{ μg}$		
Test zone (Main reaction components)	Antigen: BAR-BSA	$\cong 0.3 \text{ μg}$	Buffer solution: Phosphate buffer pH 5,5-9,3 Preservative: Sodium azide below 0,1%	Surko Bio Tech
Conjugate Pad	Mouse Anti BAR	$\cong 0.4 \text{ μg}$	Buffer solution: Phosphate buffer pH 5,5-9,3 Preservative: Sodium azide below 0,1%	Surko Bio Tech
Test strip (BZO100)				
Control zone (Internal control components)	Goat anti rabbit pAb	$\cong 0.5 \text{ μg}$	Buffer solution: Phosphate buffer pH 5,5-9,3 Preservative: Sodium azide below 0,1%	Surko Bio Tech
	Goat anti mouse pAb	$\cong 0.5 \text{ μg}$		

Test zone (Main reaction components)	Antigen: BZO-BSA	$\cong 0.3 \mu\text{g}$	Buffer solution: Phosphate buffer pH 5,5-9,3 Preservative: Sodium azide below 0,1%	Surko Bio Tech
Conjugate Pad	Mouse Anti BZO	$\cong 0.4 \mu\text{g}$	Buffer solution: Phosphate buffer pH 5,5-9,3 Preservative: Sodium azide below 0,1%	Surko Bio Tech
Test strip (COC100)				
Control zone (Internal control components)	Goat anti rabbit pAb	$\cong 0.5 \mu\text{g}$	Buffer solution: Phosphate buffer pH 5,5-9,3 Preservative: Sodium azide below 0,1%	Surko Bio Tech
	Goat anti mouse pAb	$\cong 0.5 \mu\text{g}$	Buffer solution: Phosphate buffer pH 5,5-9,3 Preservative: Sodium azide below 0,1%	
Test zone (Main reaction components)	Antigen: COC-BSA	$\cong 0.3 \mu\text{g}$	Buffer solution: Phosphate buffer pH 5,5-9,3 Preservative: Sodium azide below 0,1%	Surko Bio Tech
Conjugate Pad	Mouse Anti COC	$\cong 0.4 \mu\text{g}$	Buffer solution: Phosphate buffer pH 5,5-9,3 Preservative: Sodium azide below 0,1%	Surko Bio Tech
Test strip (OPI2000)				
Control zone (Internal control components)	Goat anti rabbit pAb	$\cong 0.5 \mu\text{g}$	Buffer solution: Phosphate buffer pH 5,5-9,3	Surko Bio Tech

	Goat anti mouse pAb	$\cong 0.5 \mu\text{g}$	Preservative: Sodium azide below 0,1%	
Test zone (Main reaction components)	Antigen: OPI-BSA	$\cong 0.3 \mu\text{g}$	Buffer solution: Phosphate buffer pH 5,5-9,3 Preservative: Sodium azide below 0,1%	Surko Bio Tech
Conjugate Pad	Mouse Anti OPI	$\cong 0.4 \mu\text{g}$	Buffer solution: Phosphate buffer pH 5,5-9,3 Preservative: Sodium azide below 0,1%	Surko Bio Tech
Test strip THC25				
Control zone (Internal control components)	Goat anti rabbit pAb	$\cong 0.5 \mu\text{g}$	Buffer solution: Phosphate buffer pH 5,5-9,3 Preservative: Sodium azide below 0,1%	Surko Bio Tech
	Goat anti mouse pAb	$\cong 0.5 \mu\text{g}$		
Test zone (Main reaction components)	Antigen: THC-BSA	$\cong 0.3 \mu\text{g}$	Buffer solution: Phosphate buffer pH 5,5-9,3 Preservative: Sodium azide below 0,1%	Surko Bio Tech
Conjugate Pad	Mouse Anti THC	$\cong 0.4 \mu\text{g}$	Buffer solution: Phosphate buffer pH 5,5-9,3 Preservative: Sodium azide below 0,1%	Surko Bio Tech
Test strip MET500				

Control zone (Internal control components)	Goat anti rabbit pAb	$\cong 0.5 \mu\text{g}$	Buffer solution: Phosphate buffer pH 5,5-9,3	Surko Bio Tech
	Goat anti mouse pAb	$\cong 0.5 \mu\text{g}$	Preservative: Sodium azide below 0,1%	
Test zone (Main reaction components)	Antigen: MET-BSA	$\cong 0.3 \mu\text{g}$	Buffer solution: Phosphate buffer pH 5,5-9,3 Preservative: Sodium azide below 0,1%	Surko Bio Tech
Conjugate Pad	Mouse Anti MET	$\cong 0.4 \mu\text{g}$	Buffer solution: Phosphate buffer pH 5,5-9,3 Preservative: Sodium azide below 0,1%	Surko Bio Tech
Test strip (MTD200)				
Control zone (Internal control components)	Goat anti rabbit pAb	$\cong 0.5 \mu\text{g}$	Buffer solution: Phosphate buffer pH 5,5-9,3	Surko Bio Tech
	Goat anti mouse pAb	$\cong 0.5 \mu\text{g}$	Preservative: Sodium azide below 0,1%	
Test zone (Main reaction components)	Antigen: MTD-BSA	$\cong 0.3 \mu\text{g}$	Buffer solution: Phosphate buffer pH 5,5-9,3 Preservative: Sodium azide below 0,1%	Surko Bio Tech
Conjugate Pad	Mouse Anti MTD	$\cong 0.4 \mu\text{g}$	Buffer solution: Phosphate buffer pH 5,5-9,3	Surko Bio Tech

			Preservative: Sodium azide below 0,1%	
Test strip (a-PVP500)				
Test zone (Main reaction components)	Antigen: a-PVP -BSA	$\cong 0.3 \mu\text{g}$	Buffer solution: Phosphate buffer pH 5,5-9,3 Preservative: Sodium azide below 0,1%	Surko Bio Tech
Control zone (Internal control components)	Goat anti rabbit pAb	$\cong 0.5 \mu\text{g}$	Buffer solution: Phosphate buffer pH 5,5-9,3 Preservative: Sodium azide below 0,1%	Surko Bio Tech
	Goat anti mouse pAb	$\cong 0.5 \mu\text{g}$		
Conjugate Pad	Mouse Anti a-PVP	$\cong 0.4 \mu\text{g}$	Buffer solution: Phosphate buffer pH 5,5-9,3 Preservative: Sodium azide below 0,1%	Surko Bio Tech
Test strip (MDMA500)				
Control zone (Internal control components)	Goat anti rabbit pAb	$\cong 0.5 \mu\text{g}$	Buffer solution: Phosphate buffer pH 5,5-9,3 Preservative: Sodium azide below 0,1%	Surko Bio Tech
	Goat anti mouse pAb	$\cong 0.5 \mu\text{g}$		
Test zone (Main reaction components)	Antigen: MDMA -BSA	$\cong 0.3 \mu\text{g}$	Buffer solution: Phosphate buffer pH 5,5-9,3 Preservative: Sodium azide below 0,1%	Surko Bio Tech

Conjugate Pad	Mouse Anti MDMA	$\cong 0.4 \mu\text{g}$	Buffer solution: Phosphate buffer pH 5,5-9,3 Preservative: Sodium azide below 0,1%	Surko Bio Tech
Test strip (PCP25)				
Control zone (Internal control components)	Goat anti rabbit pAb	$\cong 0.5 \mu\text{g}$	Buffer solution: Phosphate buffer pH 5,5-9,3 Preservative: Sodium azide below 0,1%	Surko Bio Tech
	Goat anti mouse pAb	$\cong 0.5 \mu\text{g}$		
Test zone (Main reaction components)	Antigen: PCP -BSA	$\cong 0.3 \mu\text{g}$	Buffer solution: Phosphate buffer pH 5,5-9,3 Preservative: Sodium azide below 0,1%	Surko Bio Tech
Conjugate Pad	Mouse Anti PCP	$\cong 0.4 \mu\text{g}$	Buffer solution: Phosphate buffer pH 5,5-9,3 Preservative: Sodium azide below 0,1%	Surko Bio Tech

Appendix B. Sample correlation

AMP:

GC/MS

“OLLTEST® 1 Test System”		+	-
	+	125	0
	-	5	130

Positive agreement with GC/MS: $125/(125+5) = 96.1\%$

Negative agreement with GC/MS: $130/(130+0) = 100\%$

Total agreement with GC/MS: $(125+130)/(125+5+0+130) = 98.1\%$

BAR:

GC/MS

“OLLTEST® 1 Test System”		+	-
	+	105	3
	-	2	113

Positive agreement with GC/MS: $105/(105+2) = 98.1\%$

Negative agreement with GC/MS: $113/(113+3) = 97.4\%$

Total agreement with GC/MS: $(105+113)/(105+2+3+113) = 97.8\%$

BZO:

GC/MS

“OLLTEST® 1 Test System”		+	-
	+	47	1
	-	2	50

Positive agreement with GC/MS: $47/(47+2) = 95.9\%$

Negative agreement with GC/MS: $50/(50+1) = 98\%$

Total agreement with GC/MS: $(47+50)/(47+2+1+50) = 97\%$

COC:

GC/MS

“OLLTEST® 1 Test System”		+	-
	+	165	5
	-	3	260

Positive agreement with GC/MS: $165/(165+3) = 98.2\%$

Negative agreement with GC/MS: $260/(260+5) = 98.1\%$

Total agreement with GC/MS: $(165+260)/(165+3+5+260) = 98.2\%$

MDMA:

GC/MS

“OLLTEST [®] 1 Test System”		+	-
	+	85	0
	-	0	151

Positive agreement with GC/MS: $85/(85+0) = 100\%$

Negative agreement with GC/MS: $151/(151+0) = 100\%$

Total agreement with GC/MS: $(85+151)/(85+0+0+151) = 100\%$

a-PVP:

GC/MS

OLLTEST [®] 1 Test System”		+	-
	+	52	0
	-	2	66

Positive agreement with GC/MS: $52/(52+2) = 96.3\%$

Negative agreement with GC/MS: $66/(66+0) = 100\%$

Total agreement with GC/MS: $(52+66)/(52+66+2+0) = 98.3\%$

MET:

GC/MS

OLLTEST [®] 1 Test System”		+	-
	+	63	0
	-	2	54

Positive agreement with GC/MS: $63/(63+2) = 96.9\%$

Negative agreement with GC/MS: $54/(54+0) = 100\%$

Total agreement with GC/MS: $(63+54)/(63+2+0+54) = 98.3\%$

MTD:

GC/MS

OLLTEST [®] 1 Test System”		+	-
	+	58	0

	-	0	20
--	---	---	----

Positive agreement with GC/MS: $58/(58+1) = 97.3\%$

Negative agreement with GC/MS: $20/(20+0) = 100\%$

Total agreement with GC/MS: $(58+20)/(58+1+0+20) = 98.7\%$

OPI:

GC/MS

OLLTEST® 1 Test System"		+	-
	+	205	4
	-	5	250

Positive agreement with GC/MS: $205/(205+5) = 97.6\%$

Negative agreement with GC/MS: $250/(250+4) = 98.4\%$

Total agreement with GC/MS: $(205+250)/(205+5+4+250) = 98.1\%$

PCP:

GC/MS

OLLTEST® 1 Test System"		+	-
	+	220	0
	-	5	250

Positive agreement with GC/MS: $220/(220+5) = 97.8\%$

Negative agreement with GC/MS: $250/(250+0) = 100\%$

Total agreement with GC/MS: $(220+250)/(220+5+0+250) = 98.9\%$

THC:

GC/MS

OLLTEST® 1 Test System"		+	-
	+	60	1
	-	2	59

Positive agreement with GC/MS: $60/(60+2) = 96.8\%$

Negative agreement with GC/MS: $59/(59+1) = 98.3\%$

Total agreement with GC/MS: $(60+59)/(60+2+1+59) = 97.5\%$

Appendix C. The report on Analytical sensitivity

Concentration	n	AMP		BAR		BZO		COC	
Result		-	+	-	+	-	+	-	+
Negative sample	50	50	0	50	0	50	0	50	0
50% cut-off	50	50	0	50	0	50	0	50	0
75% cut-off	50	50	0	50	0	50	0	50	0
cut-off	50	20	30	11	39	12	38	23	27
125% cut-off	50	0	50	0	50	0	50	0	50
150% cut-off	50	0	50	0	50	0	50	0	50
3X cut-off	50	0	50	0	50	0	50	0	50

Concentration	n	MDMA		a-PVP		MET	
Result		-	+	-	+	-	+
Negative sample	50	50	0	50	0	50	0
50% cut-off	50	50	0	50	0	50	0
75% cut-off	50	50	0	50	0	50	0
cut-off	50	13	37	16	34	10	40
125% cut-off	50	0	50	0	50	0	50
150% cut-off	50	0	50	0	50	0	50
3X cut-off	50	0	50	0	50	0	50

Concentration	n	MTD		OPI		PCP		THC	
Result		-	+	-	+	-	+	-	+

Negative sample	50	50	0	50	0	50	0	50	0
50% cut-off	50	50	0	50	0	50	0	50	0
75% cut-off	50	50	0	50	0	50	0	50	0
cut-off	50	6	44	23	27	9	41	11	39
125% cut-off	50	0	50	0	50	0	50	0	50
150% cut-off	50	0	50	0	50	0	50	0	50
3X cut-off	50	0	50	0	50	0	50	0	50

Appendix D. The report on Analytical specificity

<u>Amphetamine 300 related compounds</u>	<u>Concentration ng/ml</u>
d- Amphetamine	300
l- Amphetamine	50,000
Mephentermine hemisulfate salt	>100,000
3,4-Methylenedioxyamphetamine (MDA)	625
Phentermine	625
Paramethoxyamphetamine (PMA)	625
Paramethoxymethamphetamine (PMMA)	>100,000
Tyramine	>100,000
<u>Barbiturates 2000 related compounds</u>	
Secobarbital	2,000
Allobarbital	15,000
Alphenal	450
Amobarbital	10,000
Aprobarbital	1,000
Butabarbital	400
Butalbital	10,000
Butethal	1,000
Cyclopentobarbital	10,000

Pentobarbital	1,300
Phenobarbital	500
<u>Benzodiazepines 100 related compounds</u>	
Oxazepam	100
Alprazolam	42
Bromazepam	208
Chlordiazepoxide	833
Clobazam	21
Clonazepam	833
Clorazepate	1,110
Desalkflurazepam	83
Diazepam	83
Estazolam	1,667
Fentanyl	>100,000
Flunitrazepam	125
Flurazepam	>100,000
Lorazepam	417
Lormetazepam	417
Medazepam	>100,000
Midazolam	>100,000

Nitrazepam	8,333
Norchlordiazepoxide	83
Nordiazepam	167
Prazepam	>100,000
Temazepam	21
Triazolam	1,667
<u>Cocaine 100 related compounds</u>	
Benzoylecgonine	100
<u>MDMA 500 related compounds</u>	
3,4-Methylenedioxyamphetamine	500
d-Amphetamine	>100,000
l- Amphetamine	>100,000
d-Methamphetamine	>100,000
l- Methamphetamine	>100,000
3,4- Methylenedioxyamphetamine	2,500
3,4- Methylenedioxyethylamphetamine	156
Paramethoxyamphetamine	50,000
Paramethoxymethamphetamine	>100,000
<u>a-PVP 500 related compounds</u>	
a-PVP	500

3,4-Methylenedioxypropiovalerone HCl (MDPV)	40
PVP	>100,000
<u>Methamphetamine 500 related compounds</u>	
d- Methamphetamine	500
Chloroquine	12,500
Fenfluramine	12,500
l- Methamphetamine	3,125
Mephentermine hemisulfate salt	25,000
3,4- Methylenedioxyethylamphetamine	12,500
3,4-Methylenedioxymethamphetamine	1,875
Paramethoxymethamphetamine	625
(-)-Ephedrine	2,000
<u>Methadone 200 related compounds</u>	
Methadone	200
(-)-alpha-methadol	1,500
Doxylamine	3,500
LAAM HCl	6500
<u>Opiates 2000 related compounds</u>	
Morphine	2,000
Acetylcodeine	1,563

Buprenorphine	25,000
Codeine	500
Diacetylmorphine (Heroin)	1,250
Dihydrocodeine	1,563
Dihydrocodeine	800
Hydromorphone	25,000
Hydrocodone	50,000
Merperidine	>100,000
6- Monoacetylmorphine (6-MAM)	1,250
Morphine-3-β-d- glucuronide	12,500
Nalorphine Hydrochloride	>100,000
Oxycodone	>100,000
Oxymorphone	>100,000
Rifampicine	>100,000
Thebaine	50,000
<u>Phencyclidine 25 related compounds</u>	
Phencyclidine	25
Hydrocodone	12,500
Hydromorphone	6,250
4- hydroxyphencyclidine	75

Marijuana 25 related compounds

11-ноρ-Δ9-THC-9-COOH	25
11-ноρ-Δ8-THC-9-COOH	15
Δ8- Tetrahydrocannabinol	7,500
Δ9- Tetrahydrocannabinol	7,500
Cannabinol	10,000

Appendix E. Non cross-reactive compounds.

№	Name	№	Name	№	Name
1	(-)-Ephedrine (Except MET)	17	Chlorpheniramine	32	Oxalic Acid
2	(+)-Naproxen	18	Creatine	33	Penicillin – G
3	(+/-)-Ephedrine (Except MET)	19	Dextromethorphan	34	Pheniramine
4	4-Dimethylaminoantipyrene	20	Dextrophan tartrate	35	Phenothiazine
5	Acetaminophen	21	Dopamine	36	Procaine
6	Acetone	22	Erythromycin	37	Protonix
7	Albumin	23	Ethanol	38	Pseudoephedrine
8	Amitriptyline	24	Furosemide	39	Quinidine
9	Ampicillin	25	Glucose	40	Ranitidine
10	Aspartame	26	Guaiacol Glyceryl Ether	41	Sertraline
11	Aspirin	27	Hemoglobin	42	Tyramine
12	Benzocaine	28	Ibuprofen	43	Vitamin C (Ascorbic acid)
13	Bilirubin	29	Imipramine	44	Trimeprazine
14	b- Phenylethyl-amine	30	Isoproterenol	45	Venfaxine
15	Caffeine	31	Lidocaine		
16	Chloroquine				

Conclusion: The components listed in the above table show no cross-reactivity after 5 minutes in the assay at 100 µg / ml

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3. Substance Abuse and Mental Health Services Administration. Mandatory Guidelines for Federal Workplace Drug Testing Programs. 53 Federal Register; 1988.

4. McBay AJ. Drug-analysis technology--pitfalls and problems of drug testing. Clin Chem. 1987 Oct; 33 (11 Suppl): 33B-40B.

5. Gilman AG, Goodman LS, Gilman A, eds. Goodman and Gilman's The Pharmacological Basis of Therapeutics. 6th ed. New York: Macmillan; 1980.

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