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Operating Requirement:

Initial use of the product, before using the product pls check the packaging are complete:

- Product : 1 set
- User manual : 1 pc

Steps:

Sensor selection for SpO₂ measurement depends on the patient's age. The LEDs are placed in one part and the photodetector is placed in another part.

Please follow instructions below to use the adult finger SpO₂ sensor:

- Insert the sensor to the SpO₂ socket.
- Turn on the monitor.

·Attach the sensor to an appropriate site on the patient's finger

·The measurements will be displayed on the monitor's screen.

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Cable exists above finger

NOTE

Ensure the SpO₂ sensor is placed on the finger correctly. The LED part of the sensor should be to the back of the patient hand and the photodetector part to the inside. Make sure you insert the finger to a suitable depth into the sensor so that the fingernail is just opposite to the light emitted from the sensor.

To acquire accurate results, please read data until the sensor is steadily placed. Readings may not be accurate when either the SpO₂ sensor or the patient is moving.

If the accuracy of any measurement does not seem correct, first check the patient's vital signs by an alternate method. Then check the instrument for correct functioning.

Inaccurate measurements may be caused by :

- Incorrect sensor application or use;
- High-frequency electrical noise, such as noise from electrosurgical apparatus connected to the system;
- Significant levels of dysfunctional haemoglobins (e.g. carboxyhaemoglobin or methemoglobin);
- Significant concentrations of dysfunctional haemoglobin, such as carboxyhaemoglobin and methemoglobin;
- Intravascular dyes such as indocyanine green or methylene blue;
- Exposure to excessive illumination, such as surgical lamps (especially ones with a xenon light source), bilirubin lamps, fluorescent lights, infrared heating lamps, or direct sunlight (exposure to excessive illumination can be corrected by covering the sensor with a dark material);
- Excessive patient movement;
- Venous pulsations;
- SpO₂ is too low;
- improper sensor installation or incorrect contact position with the patient;
- Placement of a sensor on the same extremity with a blood pressure cuff, arterial

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catheter, or intravascular line.

Loss of pulse signal can occur in the following situation:

- The sensor is too tight;
- There is excessive illumination from light sources such as surgical lamp, a bilirubin lamp, or sunlight;
- A blood pressure cuff is inflated on the same extremity as the one with a SpO₂ sensor attached;
- The patient has hypotension, severe vasoconstriction, severe anaemia, or hypothermia;
- There is arterial occlusion close to the sensor;
- The patient is in cardiac arrest or in shock.

Contraindications:

Contraindications have not been discovered

Precautions:

NOTE

not perform SpO₂ monitoring and NIBP measurements on the same arm simultaneously, Obstruction of blood flow during NIBP measurements may adversely affect the reading of the SpO₂ value.

WARNING:

- Check if the sensor cable is in normal condition before monitoring. Do not use the SpO₂ sensor once the package or the sensor is found damaged.
- Remove the SpO₂ sensor from the patient after measurement.
- As with all medical equipment, carefully position cabling to reduce the possibility of patient injury. Cables of electrical surgical equipment should not be wound around the SpO₂ sensor.
- Do not put the sensor on extremities with the arterial catheter or venous syringes.
- If no pulse is found or the reading is unreasonable, first check the patient's condition, and then check the sensor installation and connection with the monitor, finally ask the qualified engineer to check the device and the SpO₂ sensor to ensure proper working order.
- Do not use the monitor to measure patients whose pulse rate is lower than 30bpm, this may cause incorrect results.
- Prolonged and continuous monitoring may increase a change in dermal condition such as abnormal sensitivity, rubescence, vesicle, repressive putrescence, and so on. It is ...

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important to check the skin more frequently with a child and patient with poor perfusion or immature dermatographia by light collimation and positioning correctly according to changes of the skin. Check the sensor placement 2 ~ 3 hours and move it when the skin deteriorates.

- Make sure no contamination or scar exists in the site where the sensor is placed. Otherwise, the measured result may be incorrect because the signal received by the sensor is affected.
- When used on different patients, the monitor is prone to cross contamination, which should be prevented and controlled by the user. Disinfection is recommended before using the SpO₂ sensor on other patients.

• Please use the SpO₂ sensor supplied by the monitor or compatible to the monitor.

CAUTION:

- SpO₂ sensors are fragile. Avoid from pressure and knocks
- The product should be used together with the compatible devices, regular calibration and maintenance should be carried out according to the compatible devices' requirements; Persons who implements the preventive inspection and maintenance for the devices and accessories should received special training program for the related products, regular maintenance cycle of half year recommended.
- Do not use the sensor during MRI operating, otherwise, the induced currents produced will burn the patients. In addition, the sensor will also affect the MRI imaging, and MRI equipment will also affect sensor accuracy.
- Do not transform or modify the sensor, such actions would damage the the performance or accuracies of the sensor.
- Do not disassemble or repair the sensor by yourself, otherwise, damage of the product or injury of the operator may be caused, and a direct result of the failure of the product warranty.
- Product with disposable label, please make sure that the package is intact before opening. When the package is damaged, pls do not use it! Non-repetitive use!
- When using the sensor, please make sure that the sensor cable is placed properly in order to avoid entanglement or suffocation of the patient.
- The maximum contact temperature between sensor and patient should not exceed 41° C.
- The possible impact on the application of a medical device for people with implanted into the human body medical devices, pregnant women, women during breastfeeding, children, adults with chronic diseases have not been identified.

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• Information on the possible impact of the use of a medical device on the ability to drive vehicles, mechanisms have not been identified.

Cleaning

1) Finger Sensor

- Clean the finger sensor surface with a cloth dampened with an isopropyl alcohol 70% solution or bleach 10% solution. After cleaning wipe the sensor down thoroughly with a clean water dampened cloth to rinse.
- Ensure the finger sensor windows are clean and dry before reuse.
- Do not immerse the finger sensor.

• Do not attempt to sterilize the finger sensor.

• After cleaning the finger sensor, verify the sensor is physically intact, with no broken/frayed wires or damaged parts. Ensure the connectors are clean and dry, with no signs of contamination or corrosion. Do not use a broken or damaged sensor or one with wet, contaminated or corroded connectors.

2) Y or C-Sensor

- The Y or C-sensor may be immersed-up to, but not including, the connector, in a gluteraldehyde solution, 10% bleach solution.
- Rinse thoroughly with water and dry before use.
- Do not attempt to sterilize Y or C-sensor.
- Do not immerse the connector on the Y or C-sensor.

• After cleaning the Y or C-sensor, verify the sensor is physically intact, with on broken/frayed wirs or damaged parts. Ensure the connectors are clean and dry, with no signs of contamination or corrosion. Do not use a broken or damaged sensor or one with wet, contaminated or corroded connectors.

Transport and storage environmental requirements:

The storage environment for the pacakged SpO₂ sensor, cables for SpO₂ senor should be in a clean, no corrosive gas and well-ventilated room, with temperature at - 20°C +55 °C and relative humidity 93%. Keep out of reach of children. In transport condition, proper stacking should follow the lable instruction on the outer box, and collision, severe vibration and direct rain and snow shower hit should be prevented.

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Product should be repaired by an authorized representative of the manufacturer in Russia.

Disposal:

All hospital waste must be disposed of in accordance with applicable rules and regulations.

Valid period of the product:

one year.

On the territory of the Russian Federation on issues relating to the quality of products manufactured by Shenzhen EMS Life Technology, LTD, China.

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