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Product name: NIBP cuff series

Product Model:

Manufacturer: Shenzhen EMS LIFE Technology Co.,LTD


Address: Shi Dai Jing Yuan, Fu Yong Town, Baoan District, Shenzhen P.R China. Postcode: 518103

Application Scope:

To measure the non-invasive blood pressure.


Product Working Principles of Function and Structure:

The monitor measures non-invasive blood pressure by the oscillometric method. It detects the change of blood volume in the artery and change in pressure when inflating the cuff. When the cuff is deflating, the measurements will be detected by the change of pressure. The product consists of bandage, hose and/or bladder three parts. Product structure as below:

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Key technical parameter:

1. The use life for bladder is more than 30,000 times.

Joint strength between hose and bladder: 3~5KN

Temperature: 0°C ~ 45°C

Relative Humidity: ≤95%:

Operating Requirement:

Steps:

It is important to select the proper size blood pressure cuff and place the cuff correctly.

1. Initial use of the product, before using the product pls check the packaging are complete:

Product : 1 set

User manual : 1 pc

2. Cuff selection:

1) Identify patient limb circumference.

2) Select appropriate cuff (limb circumference is identified on each cuff).

3. Cuff Placement:

1) Place cuff firmly around extremity being used.

2) Marking on the cuff should match the artery location.

3) Cuff should be 2 to 5 cm (1 to 2 inches) above the elbow if using brachial artery.

4) The air pad should be over the brachial artery. Tubing is immediately to the right or left of the brachial artery to prevent kinking when the elbow is bent. The cuff is usually placed around the patient arm, but if the patient is too fat, select a thigh cuff and place around his thigh.

Ensure the tube between the cuff and the monitor is not compressed or restricted.

The arm should be placed on a surface that level with the heart. The palm of the hand should face up. Shown as Figure 2 and Figure 3.

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
Figure 2 Placement of the Cuff


Figure 3 Placement of the Cuff

To measure blood pressure: (Figure 2 and 3)

- 1) Insert the sensor's connector to the monitor's NIBP socket;
 - 2) Turn on the monitor to enter the parameter monitoring interface;
- Place the cuff correctly around the arm;
Press blood pressure measurement button to measure the NIBP;

The readings will be displayed on the screen.

Contraindications:

Contraindications have not been discovered

Precautions

Cuff width should be forty percent of the limb's girth or two thirds of the upper arm length. The length of the cuffs inflated part should be fifty to eighty percent of the limb's girth. The wrong cuff will lead to an inaccurate reading.

The arm should be placed on a surface that is level with the heart. If this cannot be done amend the measured results as follows:

1. If the cuff is higher than the level of the patient's heart, you should add 0.9mmHg (0.10Kpa) per cm or 2.2mmHg (0.25Kpa) per inch.
2. If the cuff is lower than the level of the patient heart, you should reduce 0.9mmHg (0.10Kpa) per cm or 2.2mmHg (0.25Kpa) per inch.

If you doubt at the veracity of the readings, please look over the patient's life character before checking the patient monitor function.

For use on behalf of
SHENZHEN EMS LIFE TECHNOLOGY CO., LIMITED
深圳市易美生科技有限公司

WARNING:

- The measurement is not recommended for patients experiencing seizure or tremors, or patients with hypotension, hypertension, arrhythmias, or extremely high or low heart rate.

Do not place the cuff on a position where skin or tissue is damaged.

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Do not place the cuff on a limb being used for A-V Fistulas, intravenous infusion or any area where circulation is compromised or has the potential to be compromised

The measure patient selected in the setup menu and the cuff size used must be correct to obtain reliable NIBP data and to prevent overpressure in pediatric use.

When continues mode used long time, the skin under the cuff may occur nerve dandified or blood lost. Periodically check patient limb circulation distal to the cuff. Check frequently when using Auto NIBP in 1 and 2 minute intervals and continues mode.

The possible impact on the application of a medical device for people with implanted into the human body medical devices, pregnant women, women during breastfeeding, children, adults with chronic diseases have not been identified.

Information on the possible impact of the use of a medical device on the ability to drive vehicles, mechanisms have not been identified.

Disinfection and Cleaning:

- 1) Prepared enzyme detergent and 10% of the bleach , were put into different spray bottle.
- 2) Evenly spray the cleaner on the surface of cuff and hose, if the dirty is dry, the cleaner should be stayed on the surface for one minute.
- 3) Wipe the surface with a clean soft cloth. For the visible stains or irregular surface, use a soft hair brush to clean it.
- 4) Rinse by plenty of distilled water.
- 5) When disinfection, spray 10% of the bleach on the cuff surface, and place the product for 5 minutes.
- 6) Wipe off excess solvent night and rinse with distilled water, then let it dry.

Icon Interpretation:

Product serial number label
Lifetime mark
Production date marked
Latex free

For and on behalf of
SHENZHEN FMS LIFE TECHNOLOGY CO., LTD.
深圳市易美生科技有限公司

Wang Jiming

Transport and storage environmental requirements:

The storage environment for the NIBP cuff should be in a clean, dry and well-ventilated room, with temperature at -5°C + 40°C and relative humidity < 85%, toxic, harmful and corrosive materials should not be stroed together The valid storage period is two years. In transport condition, the transport vehicles should be clean and hygiene, proper stacking

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should follow the label instruction on the outer box, and collision, severe vibration and direct rain and snow shower hit should be prevented. Toxic, harmful and corrosive materials should not be mixed together for transportation

Disposal:

All hospital waste must be disposed of in accordance with applicable rules and regulations.

Valid period of the product:

Six months.

On the territory of the Russian Federation on issues relating to the quality of products manufactured by Shenzhen EMS Life Technology, LTD, China, please contact:

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