Medtronic

Medtronic, Inc.
710 Medtronic Parkway
Minneapolis, MN 55432-5604
USA
www.medtronic.com


Date: 17 August 2023

Declaration Letter

We, **Medtronic, Inc.**, with a registered business address at address: 710 Medtronic Pkwy., Minneapolis, MN 55432, USA, hereby affirm that the attached documentation, Technical Manual and Appendix Technical Manual for Attain Performa S MRI SureScan 4298 are true and unaltered copies of the original documents are provided.

Sincerely,
Medtronic Inc.

By: *Jennifer Anderson* 17 Aug 2023
Jennifer Anderson
Regulatory Affairs Director
Medtronic Inc.

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The following list includes trademarks or registered trademarks of Medtronic in the United States and possibly in other countries. All other trademarks are the property of their respective owners.

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1 Description

The Medtronic Attain Performa MRI SureScan 4298 steroid-eluting, quadripolar electrode, transvenous, over-the-wire LV lead is designed for pacing via a cardiac vein. The lead has been tested for use in the Magnetic Resonance Imaging (MRI) environment. All lead lengths for this lead model are MR Conditional. This lead contains 4 electrodes designed to function as cathodes or anodes, depending on how the device LV pacing polarity is programmed:

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The electrodes are positioned relative to the 2 distal curves of the lead to facilitate contact with the cardiac veins.

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The Attain Performa MRI SureScan 4298 LV lead can be positioned with the aid of a guide wire, a stylet, an inner catheter, or an inner catheter plus a guide wire or a stylet.

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1.1 Medtronic SureScan system

The Model 4298 lead is part of the Medtronic SureScan system. The SureScan system includes a Medtronic SureScan device connected to Medtronic SureScan leads. Labeling for SureScan system components displays the SureScan logo and the MR Conditional symbol.

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3 Contraindications

Before performing an MRI scan, refer to the SureScan MRI technical manual for MRI-specific contraindications. The Model 4298 lead is contraindicated as follows:

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- Rust stylets are not recommended with this lead due to the risk of conductor coil/insulation perforation.
- Verify lead length on the lead label on the connector to choose an appropriate stylet kit (downsized knob) length when repositioning. Always choose a stylet kit (downsized knob) 3 cm shorter than the lead length. For example, choose a stylet kit (downsized knob) with stylets 75 cm long for a lead 78 cm long.
- Lead removal may result in evulsion of the endocardium, valve, or vein.
- Lead junctions may separate, leaving the lead tip and bare wire in the heart or vein.
- Chronic repositioning may adversely affect the low-threshold performance of a steroid-eluting lead.
- Cap abandoned leads to avoid transmitting electrical signals.
- For leads that have been severed, seal the remaining lead end and suture the lead to adjacent tissue.
- If a lead is removed and repositioned, inspect it carefully for insulator or conductor coil damage before repositioning.

Magnetic resonance imaging (MRI) – An MRI is a type of medical imaging that uses magnetic fields to create an internal view of the body. If certain criteria are met and the warnings and precautions provided by Medtronic are followed, patients with an MRI Conditional device and lead system are able to undergo an MRI scan; for details, refer to the SureScan MRI technical manual that Medtronic provides for an MRI Conditional device.

Diathermy treatment (Including therapeutic ultrasound) – Diathermy is a treatment that involves the therapeutic heating of body tissues. Diathermy treatments include high frequency, short wave, microwave, and therapeutic ultrasound. Except for therapeutic ultrasound, do not use diathermy treatments on cardiac device patients. Diathermy treatments may result in serious injury or damage to an implanted device and lead system. Therapeutic ultrasound (including physiotherapy, high intensity therapeutic ultrasound, and high intensity focused ultrasound), is the use of ultrasound at higher energies than diagnostic ultrasound to bring heat or agitation into the body. Therapeutic ultrasound is acceptable if treatment is performed with a minimum separation distance of 15 cm between the applicator and the implanted device and lead system, as long as the ultrasonic beam is pointing away from the device and lead system.

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To implant the Attain Performa MRI SureScan 4298 LV lead in a selected cardiac vein, a compatible delivery system is required, such as a Medtronic delivery system. A compatible delivery system includes a guide catheter and either a hemostasis valve or an introducer valve that can be removed or that allows passage over the lead connector. Contact a Medtronic representative for further information regarding compatible delivery systems.

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
1. Insert a straight stylet or a guide wire into the lead to vary the shape of the distal end of the lead (Figure 2).
- Note:** When a stylet is fully inserted, the distal tip of the stylet does not reach the distal tip of the lead.

Figure 2.

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6.6 Lead placement

The Attain Performa MRI SureScan 4298 LV lead can be placed with the aid of a guide wire, a stylet, an inner catheter, or an inner catheter plus a guide wire or a stylet.

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
1. Select a guide wire.
If the patient has tortuous anatomy, a more flexible guide wire is recommended. If additional support is needed, use a firmer guide wire.
If you are using an Attain Hybrid guide wire, the steps involving the insertion tool (Step 2 through Step 4) do not apply.
2. Insert the guide wire into the lead by placing the distal (flexible) end of the guide wire into the lead connector pin using the guide wire insertion tool included in the package (Figure 6). To prevent the lead from dislodging from the guide wire insertion tool, grasp the lead and guide wire insertion tool firmly between thumb and forefinger.
Caution: To minimize the risk of damaging the guide wire, be sure that the flexible section of the guide wire is fully inserted into the lead before removing the guide wire insertion tool from the lead.
Note: Be sure to remove the guide wire insertion tool before lead implant.


Figure 6.

- 1 Guide wire
- 2 Lead connector pin

3. Disengage the guide wire insertion tool from the lead connector pin.
4. Remove the guide wire insertion tool by sliding the tool off the end of the guide wire.
5. Position the guide wire steering handle (Figure 7):
For an Attain Hybrid guide wire, preload the guide wire steering handle onto the guide wire before loading the guide wire into the connector pin on the proximal end of the lead.
For other guide wires:
 - a. Advance the guide wire steering handle over the proximal (rigid) end of the guide wire.
 - b. Tighten the guide wire steering handle onto the guide wire near the lead connector pin.

Figure 7.

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
Using fluoroscopy for guidance, fixate the lead by wedging both curves of the lead into the cardiac vein (Figure 9).

Figure 9.

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
6.10 Taking electrical measurements

The Attain Performa MRI SureScan 4298 LV lead was designed to provide pacing via selectable electrodes. The 16 available pacing polarities are shown in Figure 10.

Figure 10.

- 1 Extended bipolar pacing polarities
- 2 Bipolar (reversible) pacing polarities

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Note: If you are monitoring the LV EGM signal when the LV2 to LV3 or LV3 to LV2 pacing polarity is selected, the signal amplitude may be attenuated compared to other pacing configurations such as LV1 to RVcoil. This signal attenuation is an expected characteristic of the Attain Performa MRI SureScan 4298 LV lead with the short spacing between electrodes LV2 and LV3.

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
Figure 12.

1 Removing the AccuRead 2.0 tool from the connector pin

2 Removing the AccuRead 2.0 tool from the stylet using the slit on the side of the tool

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1. If used, remove the guide wire and guide wire insertion tool. Replace the guide wire with a straight stylet (downsized knob). Insert the straight stylet into the lead to the mid-coronary sinus.
2. Remove the guide catheter from the lead. See the delivery system product documentation for details.
Note: For Medtronic slittable delivery systems, use a slitter compatible with a 1.75 mm (5.3 French) lead body.

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8 CE mark of conformity

The following CE mark of conformity applies to the Model 4298 lead.

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**Authorized Representative in the European
Community**
Medtronic B.V.
Earl Bakkenstraat 10
6422 RJ Heerlen
The Netherlands
+31 45 566 8000

Europe/Middle East/Africa
Medtronic International Trading Sàrl
Route du Molliat 31
Case Postale 84
CH-1131 Tolochenaz
Switzerland
+41 21 802 7000

Australia

Medtronic Australasia Pty Ltd
97 Waterloo Road
North Ryde, NSW 2113
Australia

Canada

Medtronic of Canada Ltd
99 Hereford Street
Brampton, Ontario L6Y 0R3
Canada
+1 905 460 3900

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1. Name of medical product	Cardiac vein pacing lead Attain Performa. Cardiac vein pacing lead ATTAIN PERFORMA of various version: Cardiac vein pacing lead ATTAIN PERFORMA MRI SURESCAN 4298 size: 78 cm, 88 cm, composed: • lead with anchoring sleeve; • guide wire insertion tool; • guide wire clip; • guide wire steering handle; • stylets - no more than 4 pcs; • AccuRead 2.0 analyzer cable interface tools - no more than 2 pcs; • product documentation: Technical Manual; Limited Warranty; Stickers.
2. Information about manufacturer of medical device (name, location address, mailing address, phone number, address of place of production)	Design Site: Medtronic Inc., USA 8200 Coral Sea St., Mounds View, Minnesota 55112 Tel.: +1 800-633-8766; www.medtronic.com Legal Manufacturer: Medtronic Inc., USA 710 Medtronic Parkway, Minneapolis, Minnesota 55432, USA Tel.: +1 800-633-8766 www.medtronic.com Authorized Representative of Manufacturer: Medtronic, LLC, 123112, Moscow, in.ter.c. Presnensky Municipal District, Presnenskaya nab., 10, floor 9, office III, room 41 Tel.: (495) 580-73-77, Fax: (495) 580-73-78, www.medtronic.ru Holder of Registration Certificate in Russia (name, address, tel, fax, e-mail, website): Medtronic, LLC,

	123112, Moscow, in.ter.c. Presnensky Municipal District, Presnenskaya nab., 10, floor 9, office III, room 41 Tel.: (495) 580-73-77, Fax: (495) 580-73-78, www.medtronic.ru Address of Manufacturing Site for medical device: Medtronic Puerto Rico Operations Co., Villalba Rd. 149, KM. 56.3, Call Box 6001, Villalba, PR 00766, USA
3. Prescription of a medical device	Intended for left ventricular stimulation through the vein of the heart in the left ventricle, as part of the biventricular pacing system.
4. Information on potential consumers of medical device	This medical device is intended for use in clinical practice by qualified physicians and medical specialists.
5. Sterilization method, Sterilization Validation, Packaging Validation for Sterile devices	Sterilization – Medtronic has sterilized the package contents with ethylene oxide before shipment. This lead is for single use only and is not intended to be resterilized. Device is sterile and non-pyrogenic. - Sterilized by EO. - Sterility Assurance Level (SAL)- 10^{-6}.
6. Procedure and conditions for the disposal or destruction of a medical device	Dispose of all single-use accessories according to local environmental requirements. After use, this product may become biohazardous. Unused products with expired shelf life are classified as Class A medical waste. Products extracted from the patient's body are classified as Class B medical waste.
7. Class risk depending on the potential risk of a medical product, in accordance with the nomenclature	Classification of medical device: class risk 3; B - products whose failure reduces the effectiveness or delays the treatment and diagnostic process in non-critical situations, or increases the burden on medical or maintenance personnel; After use, this product may be a potential biohazard; Unused, expired devices are classified as medical waste class «A»; Leads extracted from the patient's body are medical waste Class B; Class of electrical safety, laser safety, and other types of hazards: not applicable; Type of working part: not applicable; The device is sterile, not re-sterilized.

8. Methods and apparatus disinfection and sterilization cleaning	Not applicable. Sterile. For single use only.	
9. Shelf life	Shelf life no more 24 months from the date of sterile packaging.	
10. Operating conditions	In a medical and preventative treatment facility. The device should be used only by or under the supervision of a trained physician.	
		Service conditions
	Operating temperature	from + 32°C to + 42°C
	Humidity rate	80%-100%
	Atmospheric pressure	from 700 hPa to 1060 hPa
11. Conditions of transportation and storage	Storage temperature – Store at 25°C (77°F). Excursions from this storage temperature are permitted in the range of 15 to 30°C (59 to 86°F). Transient spikes up to 40°C (104°F) are permitted as long as they do not exceed 24 hours.	
12. Warranty	Disclaimer of warranty and general warning	
	Although the Medtronic lead, hereafter referred to as “Product”, has been carefully designed, manufactured and tested prior to sale, the Product may fail to perform its intended function satisfactorily for a variety of reasons. Warnings contained in the Product labeling provide more detailed information and are considered an integral part of this disclaimer of warranty. Medtronic, therefore, disclaims all warranties, both express and implied, with respect to the Product.	
	Medtronic shall not be liable for any incidental or consequential damages caused by any use, defect, or failure of the Product, whether the claim is based on warranty, contract, tort, or otherwise.	
	The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this disclaimer of warranty is held by any court of competent jurisdiction to be illegal, unenforceable, or in any conflict with applicable law, the validity of the remaining portion of the disclaimer of warranty shall not be affected, and all rights and obligations shall be construed and enforced as if this disclaimer of warranty did not contain the particular part or term held to be invalid.	
	General warning	
	Medtronic implantable leads are implanted in the extremely hostile environment of the human body. Leads are necessarily very small in diameter and must still be very flexible, which unavoidably reduces their potential performance or longevity.	
	Leads may fail to function for a variety of causes, including, but not limited to:	

**Medtronic**

Medtronic Inc.
710 Medtronic Parkway N.E.
Minneapolis, MN 55432
USA

	Medical complications, body rejection phenomena, allergic reaction, fibrotic tissue, or failure of leads by breakage or by breach of their insulation covering. In addition, despite the exercise of all due care in design, component selection, manufacture, and testing prior to sale, leads may be easily damaged before, during, or after insertion by improper handling or other intervening acts. Consequently, no representation or warranty is made that failure or cessation of function of leads will not occur; or that the body will not react adversely to the implantation of leads; or that medical complications (including perforation of the heart) will not follow implantation of leads; or that the lead will, in all cases, restore adequate cardiac function. Warranty shelf life: The warranty covers the entire shelf life of the medical device.
13. The list of national and international regulatory documents/standards, which matches medical device	EN 45502-1, EN 45502-2-1, EN ISO 15223-1, EN ISO 14971, EN ISO 11607-1, EN 1041, EN ISO 11135, EN 980, EN 62366, EN ISO 10993-1, ISO 10993-3, ISO 10993-4, ISO 10993-5, ISO 10993-6, ISO 10993-7, ISO 10993-10, ISO 10993-11, ISO 27186.
14. Customer Service	For more information, contact your local Medtronic representative, or contact Medtronic Inc. as listed below: Medtronic Inc., USA Information for patients: Toll-free +1 800 551 5544 (7 a.m. – 6 p.m., Monday–Friday, central time) Fax +1 763 514 1855 Representative in Russia: Medtronic, LLC, 123112, Moscow, in.ter.c. Presnensky Municipal District, Presnenskaya nab., 10, floor 9, office III, room 41 Tel.: (495) 580-73-77, Fax: (495) 580-73-78, www.medtronic.ru

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

Date: 17 August 2023


Declaration Letter

We, **Medtronic, Inc.**, with a registered business address at address: 710 Medtronic Pkwy., Minneapolis, MN 55432, USA, hereby affirm that the attached documentation, Technical Manual and Appendix Technical Manual for Attain Performa S MRI SureScan 4398 are true and unaltered copies of the original documents are provided.

Sincerely,
Medtronic Inc.

By:  17 Aug 2023
Jennifer Anderson
Regulatory Affairs Director
Medtronic Inc.

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Contents

1	Description	3
2	Indications	4
3	Contraindications	4
4	Warnings and precautions	4
5	Potential adverse events	6
6	Implant procedure	7
7	Specifications	13
8	CE mark of conformity	15
9	Explanation of symbols on package labeling	15
10	Medtronic disclaimer of warranty	16
11	Service	16

1 Description

The Medtronic Attain Performa Straight MRI SureScan 4398 steroid-eluting, quadripolar electrode, transvenous, over-the-wire LV lead is designed for pacing via a cardiac vein. The lead has been tested for use in the Magnetic Resonance Imaging (MRI) environment. All lead lengths for this lead model are MR Conditional. This lead contains 4 electrodes designed to function as cathodes or anodes, depending on how the device LV pacing polarity is programmed:

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¹ IS4-LLLL refers to an International Connector Standard (ISO 27186:2010) whereby pulse generators and leads so designated are assured of a basic mechanical fit.

LLLL defines the lead connector contacts as low voltage (L).

² Technology developed by NASA.

See Section 7.2, "Specifications drawing (nominal)", page 14 for an illustration of the lead electrodes and connector contacts.

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Each electrode contains a Monolithic controlled release device (MCRD) for elution of steroid to reduce inflammatory response within the cardiac vein. The MCRDs contain a combined total dosage of <1.0 mg of dexamethasone acetate steroid. Upon exposure to body fluids, the steroid elutes from the MCRDs. The steroid suppresses the inflammatory response that is believed to cause threshold rise typically associated with implanted pacing electrodes.

The lead features 2 polyurethane lines located between electrode LV1 and electrode LV2 to facilitate contact with the cardiac veins.

The outer insulation of the lead is polyurethane and the inner insulation is Si-polyimide (SI-PI)². The SI-PI is applied as a coating to the conductor wire before coiling.

The Attain Performa Straight MRI SureScan 4398 LV lead can be positioned with the aid of a guide wire, a stylet, an inner catheter, or an inner catheter plus a guide wire or a stylet.

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1.1 Medtronic SureScan system

The Model 4398 lead is part of the Medtronic SureScan system. The SureScan system includes a Medtronic SureScan device connected to Medtronic SureScan leads. Labeling for SureScan system components displays the SureScan logo and the MR Conditional symbol.

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3 Contraindications

Before performing an MRI scan, refer to the SureScan MRI technical manual for MRI-specific contraindications. The Model 4396 lead is contraindicated as follows:

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Magnetic resonance imaging (MRI) – An MRI is a type of medical imaging that uses magnetic fields to create an internal view of the body. If certain criteria are met and the warnings and precautions provided by Medtronic are followed, patients with an MR Conditional device and lead system are able to undergo an MRI scan; for details, refer to the SureScan MRI technical manual that Medtronic provides for an MR Conditional device.

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Potential adverse event	Indicator of potential adverse event	Corrective actions to consider
Lead dislodgement ^a	Intermittent or continuous loss of capture or LV EGM signal integrity (including sensing) ^a	Reprogram the LV pacing polarity. Reposition the lead.
Lead dislodgement ^a	Intermittent or continuous oversensing	Reprogram the LV pacing polarity. Reposition the lead.
Lead conductor fracture	Intermittent or continuous loss of capture or LV EGM signal integrity (including sensing) ^a	Replace the lead. Reprogram the LV pacing polarity.

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To implant the Attain Performa Straight MRI SureScan 4398 LV lead in a selected cardiac vein, a compatible delivery system is required, such as a Medtronic delivery system. A compatible delivery system includes a guide catheter and either a hemostasis valve or an introducer valve that can be removed or that allows passage over the lead connector. Contact a Medtronic representative for further information regarding compatible delivery systems.

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
- To access the subclavian vein, use a preferred method based on professional experience.

Caution: Certain anatomical abnormalities, such as thoracic outlet syndrome, may precipitate pinching and subsequent fracture of the lead.

Caution: Insertion should be done as far lateral as possible to avoid clamping the lead body between the clavicle and the first rib (Figure 1).

Figure 1.

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
1. Insert a stylet or a guide wire into the lead to vary the shape of the distal end of the lead.
Note: When a stylet is fully inserted, the distal tip of the stylet does not reach the distal tip of the lead.
2. Insert the lead into the delivery system. See the delivery system product documentation for additional information.

6.6 Lead placement

The Attain Performa Straight MRI SureScan 4398 LV lead can be placed with the aid of a guide wire, a stylet, an inner catheter, or an inner catheter plus a guide wire or a stylet.

Stylet delivery – If the patient's anatomy features gentle vein angulation off the coronary sinus and the cardiac vein branch is not tortuous (Figure 2), a stylet may be used for lead delivery.


Figure 2.

- 1 Coronary sinus
- 2 Cardiac vein

Guide wire delivery – If the patient's anatomy features acute vein angulation off the coronary sinus and the cardiac vein branch is tortuous (Figure 3), a guide wire may be used for lead delivery.

Figure 3.

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Caution: Do not use a sharp object to impart a curve to the distal end of the stylet. Imparting a curve to the stylet can be accomplished with a smooth-surface, sterile instrument (Figure 4).

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6.8 Placing the lead using a guide wire

Warning: Do not insert the proximal end of the guide wire through the lead tip seal without using the guide wire insertion tool. Inserting the guide wire without the guide wire insertion tool could damage the lead tip seal or the conductor coil or insulation.

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
- Select a guide wire.
If the patient has tortuous anatomy, a more flexible guide wire is recommended. If additional support is needed, use a firmer guide wire.
If you are using an Attain Hybrid guide wire, the steps involving the insertion tool (Step 2 through Step 4) do not apply.
- Insert the guide wire into the lead by placing the distal (flexible) end of the guide wire into the lead connector pin using the guide wire insertion tool included in the package (Figure 5). To prevent the lead from dislodging from the guide

wire insertion tool, grasp the lead and guide wire insertion tool firmly between thumb and forefinger.

Caution: To minimize the risk of damaging the guide wire, be sure that the flexible section of the guide wire is fully inserted into the lead before removing the guide wire insertion tool from the lead.

Note: Be sure to remove the guide wire insertion tool before lead implant.


Figure 5.

- Guide wire
- Lead connector pin

- Disengage the guide wire insertion tool from the lead connector pin.
- Remove the guide wire insertion tool by sliding the tool off the end of the guide wire.
- Position the guide wire steering handle (Figure 6):
For an Attain Hybrid guide wire, preload the guide wire steering handle onto the guide wire before loading the guide wire into the connector pin on the proximal end of the lead.
For other guide wires:
 - Advance the guide wire steering handle over the proximal (rigid) end of the guide wire.
 - Tighten the guide wire steering handle onto the guide wire near the lead connector pin.

Figure 6.

- Guide wire steering handle
- Lead connector pin

- Attach the guide wire clip to the guide wire and secure it within the sterile field. Medtronic recommends securing the guide wire clip to the patient's sterile surgical drape.

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
Insert the guide wire into the lead by placing the proximal (rigid) end of the guide wire into the distal lead tip using the guide wire insertion tool included in the package (Figure 7).

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Caution: If the distal end of the guide wire becomes severely kinked or twisted, it may be difficult to withdraw it back through the lead. Therefore, if there is an indication that the distal end of the guide wire has become damaged, or if there is significant resistance in guide wire passage, remove the lead and guide wire together as a unit. Remove the guide wire from the lead and insert a new guide wire into the lead. Do not use excessive force to retract the guide wire from the lead.

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
Using fluoroscopy for guidance, fixate the lead into the cardiac vein (Figure 8).

Figure 8.

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
6.10 Taking electrical measurements

The Attain Performa Straight MRI SureScan 4396 LV lead was designed to provide pacing via selectable electrodes. The 16 available pacing polarities are shown in Figure 9.

Figure 9.

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Note: If you are monitoring the LV EGM signal when the LV2 to LV3 or LV3 to LV2 pacing polarity is selected, the signal amplitude may be attenuated compared to other pacing configurations such as LV1 to RVcoil. This signal attenuation is an expected characteristic of the Attain Performa Straight MRI SureScan 4398 LV lead with the short spacing between electrodes LV2 and LV3.

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Note: Before taking electrical measurements, it is recommended that you retract the stylet or guide wire, inside the lead lumen, to a point proximal of all electrodes. This action allows the lead tip to resume its normal shape, enabling appropriate electrode-tissue contact to occur.

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
1. Ensure that the lead connector is inserted into the AccuRead 2.0 tool completely. If the AccuRead 2.0 tool is properly attached, the connector pin is accessible (see Figure 10).

Figure 10.

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Measurement recommended	Left ventricle
Maximum acute stimulation thresholds ^b	3.0 V
Minimum acute LV EGM signal amplitude	4.0 mV ^c

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
4. If electrical measurements do not stabilize to acceptable levels, it may be necessary to reposition the lead and repeat the testing procedure.

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5. Check for phrenic nerve stimulation by pacing at 10 V and a pulse width setting greater than 0.5 ms. Then observe for diaphragmatic contracting either by fluoroscopy or direct abdominal palpitation. Further testing may include patient positional changes to simulate upright chronic conditions. If phrenic nerve stimulation occurs, reduce the voltage until a phrenic nerve stimulation threshold is determined. Phrenic nerve stimulation may necessitate changing the pacing polarity or repositioning of the lead.
6. The AccuRead 2.0 tool may be removed from a guide wire or stylet using the slit on the side of the tool (see Figure 11).

Figure 11.

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1. If used, remove the guide wire and guide wire insertion tool. Replace the guide wire with a straight stylet (downsized knob). Insert the straight stylet into the lead to the mid-coronary sinus.
2. Remove the guide catheter from the lead. See the delivery system product documentation for details.
Note: For Medtronic slittable delivery systems, use a slitter compatible with a 1.75 mm (5.3 French) lead body.
3. Carefully and completely remove the stylet or Attain Hybrid guide wire. When removing the stylet, grip the lead firmly close to the distal end of the connector pin; gripping the lead at this location helps prevent possible lead tip dislodgement.
4. Ensure that the lead connector is inserted into the AccuRead 2.0 analyzer cable interface tool and then repeat the electrical measurements. See Section 6.10, "Taking electrical measurements", page 10.

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- Use only nonabsorbable sutures to anchor the lead.
- Do not attempt to remove or cut the anchoring sleeve.
- Do not use the anchoring sleeve tabs for suturing.
- During lead anchoring, take care to avoid dislodging the lead tip.
- Do not secure sutures so tightly that they damage the vein, lead, or anchoring sleeve (Figure 12).
- Do not tie a suture directly to the lead body (Figure 12).

Figure 12.

Anchor the lead using all 3 grooves:

1. Position the anchoring sleeve against or near the vein.
2. Secure the anchoring sleeve to the lead body by tying a suture firmly in each of the 3 grooves (Figure 13).

Figure 13.

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
- Ensure that the leads do not leave the device at an acute angle.
- Do not grip the lead or device with surgical instruments.
- Do not coil the lead. Coiling the lead can twist the lead body and may result in lead dislodgement (Figure 14).

Figure 14.

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1. To prevent undesirable twisting of the lead body or bending of the lead body at an acute angle, rotate the device to wrap the excess lead length loosely (Figure 15).

Figure 15.

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7 Specifications

7.1 Specifications (nominal)


Parameter	4396
Serial number prefix	QUB
Type	Quadripolar electrode
Chamber paced	Left ventricle
Length	78 and 88 cm
Connector	IS4-LLLL
Material	Conductor: 25% Ag-con-MP35N
	Insulators: Polyurethane (outer) SI-polyimide (SI-PI) (inner) ^a
	Electrodes: Platinum iridium alloy with titanium nitride coating
	Connector pin: MP35N
	Connector rings: MP35N
	Molded tip seal: Silicone rubber
Electrode configuration	Radiused, titanium-nitride-coated, steroid-eluting
Diameter	Lead body: 1.75 mm (5.3 French)
	Electrodes: 1.70 mm (5.1 French)
	Spacer between electrodes LV1 (distal) and LV2: 1.30 mm (3.9 French)
	Spacer between electrodes LV2 and LV3: 1.57 mm (4.7 French)
	Spacer between electrodes LV3 and LV4 (proximal): 1.30 mm (3.9 French)
Medtronic delivery system (recommended inner diameter)	1.90 mm (5.7 French)
Diagnostic guide wire (recommended diameter)	0.36 mm to 0.46 mm
Electrode surface area	5.8 mm ²

Parameter		4398
Distance between electrodes	Electrode LV1 (distal) to LV2:	21 mm
	Electrode LV2 to LV3:	1.3 mm
	Electrode LV3 to LV4 (proximal):	21 mm
Conductor resistance	LV1	22 $\pm$ 5 Ω (78 cm) 24 $\pm$ 6 Ω (88 cm)
	LV2	19 $\pm$ 4 Ω (78 cm) 21 $\pm$ 5 Ω (88 cm)
	LV3	18 $\pm$ 4 Ω (78 cm) 21 $\pm$ 4 Ω (88 cm)
	LV4	17 $\pm$ 4 Ω (78 cm) 20 $\pm$ 4 Ω (88 cm)
Steroid		Dexamethasone acetate
Combined total of steroid		<1.0 mg
	Steroid binder	Silicone rubber

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8 CE mark of conformity

The following CE mark of conformity applies to the Model 4398 lead.

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Symbol	Explanation
CE0123	Conformité Européenne (European Conformity). This symbol means that the device fully complies with applicable European Union acts.
	Do not use if package is damaged
	Do not reuse
	Upper limit of temperature
	Sterilized using ethylene oxide
	Consult instructions for use
	Consult instructions for use at this website
	Date of manufacture
	Manufacturer
	Use by

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


64 Srpski

Printing instructions: doc: #163256
Category: Leads single pig.
Size: 7.125 x 8 in (181 x 203 mm)

502111-005

M958923A026 A Medtronic Confidential

Composed: 2015-05-01 07:12:29
XSL-Stylesheet C 18-MAR-2015

Medtronic

Medtronic, Inc.
710 Medtronic Parkway
Minneapolis, MN 55432
USA
www.medtronic.com
+1 763 514 4000

**Authorized Representative in the European
Community**

Medtronic B.V.
Earl Bakkenstraat 10
6422 RJ Heerlen
The Netherlands
+31 45 586 8000

Europe/Middle East/Africa
Medtronic International Trading Sàrl
Route du Molliau 31
Case Postale 84
CH-1131 Tolochenaz
Switzerland
+41 21 802 7000

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1. Name of medical product	Cardiac vein pacing lead Attain Performa. Cardiac vein pacing lead ATTAIN PERFORMA of various version: Cardiac vein pacing lead ATTAIN PERFORMA STRAIGHT MRI SURESCAN 4398 size: 78 cm, 88 cm, composed: • lead with anchoring sleeve; • guide wire insertion tool; • guide wire clip; • guide wire steering handle; • stylets - no more than 4 pcs; • AccuRead 2.0 analyzer cable interface tools - no more than 2 pcs; • product documentation: Technical Manual; Limited Warranty; Stickers.
2. Information about manufacturer of medical device (name, location address, mailing address, phone number, address of place of production)	Design Site: Medtronic Inc., USA 8200 Coral Sea St., Mounds View, Minnesota 55112 Tel.: +1 800-633-8766; www.medtronic.com Legal Manufacturer: Medtronic Inc., USA 710 Medtronic Parkway, Minneapolis, Minnesota 55432, USA Tel.: +1 800-633-8766 www.medtronic.com Authorized Representative of Manufacturer: Medtronic, LLC, 123112, Moscow, in.ter.c. Presnensky Municipal District, Presnenskaya nab., 10, floor 9, office III, room 41 Tel.: (495) 580-73-77, Fax: (495) 580-73-78, www.medtronic.ru Holder of Registration Certificate in Russia (name, address, tel, fax, e-mail, website):

	Medtronic, LLC, 123112, Moscow, in.ter.c. Presnensky Municipal District, Presnenskaya nab., 10, floor 9, office III, room 41 Tel.: (495) 580-73-77, Fax: (495) 580-73-78, www.medtronic.ru Address of Manufacturing Site for medical device: Medtronic Puerto Rico Operations Co., Villalba Rd. 149, KM. 56.3, Call Box 6001, Villalba, PR 00766, USA
3. Prescription of a medical device	Intended for left ventricular stimulation through the vein of the heart in the left ventricle, as part of the biventricular pacing system.
4. Information on potential consumers of medical device	This medical device is intended for use in clinical practice by qualified physicians and medical specialists.
5. Sterilization method, Sterilization Validation, Packaging Validation for Sterile devices	Sterilization – Medtronic has sterilized the package contents with ethylene oxide before shipment. This lead is for single use only and is not intended to be resterilized. Device is sterile and non-pyrogenic. - Sterilized by EO. - Sterility Assurance Level (SAL)- 10^{-6}.
6. Procedure and conditions for the disposal or destruction of a medical device	Dispose of all single-use accessories according to local environmental requirements. After use, this product may become biohazardous. Unused products with expired shelf life are classified as Class A medical waste. Products extracted from the patient's body are classified as Class B medical waste.
7. Class risk depending on the potential risk of a medical product, in accordance with the nomenclature	Classification of medical device: class risk 3; B - products whose failure reduces the effectiveness or delays the treatment and diagnostic process in non-critical situations, or increases the burden on medical or maintenance personnel; After use, this product may be a potential biohazard; Unused, expired devices are classified as medical waste class «A»; Leads extracted from the patient's body are medical waste Class B; Class of electrical safety, laser safety, and other types of hazards: not applicable; Type of working part: not applicable; The device is sterile, not re-sterilized.

8. Methods and apparatus for disinfection and sterilization cleaning	Not applicable. Sterile. For single use only.	
9. Shelf life	Shelf life no more 24 months from the date of sterile packaging.	
10. Operating conditions	In a medical and preventative treatment facility. The device should be used only by or under the supervision of a trained physician.	
		Service conditions
	Operating temperature	from + 32°C to + 42°C
	Humidity rate	80%-100%
11. Conditions of transportation and storage	Atmospheric pressure	from 700 hPa to 1060 hPa
	Storage temperature – Store at 25°C (77°F). Excursions from this storage temperature are permitted in the range of 15 to 30°C (59 to 86°F). Transient spikes up to 40°C (104°F) are permitted as long as they do not exceed 24 hours.	
12. Warranty	Disclaimer of warranty and general warning Although the Medtronic lead, hereafter referred to as “Product”, has been carefully designed, manufactured and tested prior to sale, the Product may fail to perform its intended function satisfactorily for a variety of reasons. Warnings contained in the Product labeling provide more detailed information and are considered an integral part of this disclaimer of warranty. Medtronic, therefore, disclaims all warranties, both express and implied, with respect to the Product. Medtronic shall not be liable for any incidental or consequential damages caused by any use, defect, or failure of the Product, whether the claim is based on warranty, contract, tort, or otherwise. The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this disclaimer of warranty is held by any court of competent jurisdiction to be illegal, unenforceable, or in any conflict with applicable law, the validity of the remaining portion of the disclaimer of warranty shall not be affected, and all rights and obligations shall be construed and enforced as if this disclaimer of warranty did not contain the particular part or term held to be invalid. General warning Medtronic implantable leads are implanted in the extremely hostile environment of the human body. Leads are necessarily very small in diameter and must still be very flexible, which unavoidably reduces their potential performance or longevity. Leads may fail to function for a variety of causes, including, but not limited to:	

	Medical complications, body rejection phenomena, allergic reaction, fibrotic tissue, or failure of leads by breakage or by breach of their insulation covering. In addition, despite the exercise of all due care in design, component selection, manufacture, and testing prior to sale, leads may be easily damaged before, during, or after insertion by improper handling or other intervening acts. Consequently, no representation or warranty is made that failure or cessation of function of leads will not occur; or that the body will not react adversely to the implantation of leads; or that medical complications (including perforation of the heart) will not follow implantation of leads; or that the lead will, in all cases, restore adequate cardiac function. Warranty shelf life: The warranty covers the entire shelf life of the medical device.
13. The list of national and international regulatory documents/standards, which matches medical device	EN 45502-1, EN 45502-2-1, EN ISO 15223-1, EN ISO 14971, EN ISO 11607-1, EN 1041, EN ISO 11135, EN 980, EN 62366, EN ISO 10993-1, ISO 10993-3, ISO 10993-4, ISO 10993-5, ISO 10993-6, ISO 10993-7, ISO 10993-10, ISO 10993-11, ISO 27186.
14. Customer Service	For more information, contact your local Medtronic representative, or contact Medtronic Inc. as listed below: Medtronic Inc., USA Information for patients: Toll-free +1 800 551 5544 (7 a.m.– 6 p.m., Monday–Friday, central time) Fax +1 763 514 1855 Representative in Russia: Medtronic, LLC, 123112, Moscow, in.ter.c. Presnensky Municipal District, Presnenskaya nab., 10, floor 9, office III, room 41 Tel.: (495) 580-73-77, Fax: (495) 580-73-78, www.medtronic.ru

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Date: 17 August 2023

Declaration Letter


We, **Medtronic, Inc.**, with a registered business address at address: 710 Medtronic Pkwy., Minneapolis, MN 55432, USA, hereby affirm that the attached documentation, Technical Manual and Appendix Technical Manual for Attain Performa S MRI SureScan 4598 are true and unaltered copies of the original documents are provided.


Sincerely,
Medtronic Inc.

By:

 17 Aug 2023

Jennifer Anderson
Regulatory Affairs Director
Medtronic Inc.

17 Aug 2023

Medtronic

ATTAIN PERFORMA™ S MRI SURESCAN™ 4598

Steroid-eluting, quadripolar electrode, transvenous, over-the-wire, cardiac vein pacing lead

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AccuRead, Attain, Attain Hybrid, Medtronic, Performa, SureScan

Contents

1	Description	3
2	Indications	4
3	Contraindications	4
4	Warnings and precautions	4
5	Potential adverse events	6
6	Implant procedure	7
7	Specifications	14
8	CE mark of conformity	15
9	Explanation of symbols on package labeling	15
10	Medtronic disclaimer of warranty	16
11	Service	16

1 Description

The Medtronic Attain Performa S MRI SureScan 4598 steroid-eluting, quadripolar electrode, transvenous, over-the-wire LV lead is designed for pacing via a cardiac vein. The lead has been tested for use in the Magnetic Resonance Imaging (MRI) environment. All lead lengths for this lead model are MR Conditional. This lead contains 4 electrodes designed to function as cathodes or anodes, depending on how the device LV pacing polarity is programmed:

- electrode LV1, the distal electrode, positioned near the distal tip of the lead
- electrode LV2, positioned 21 mm proximal to electrode LV1
- electrode LV3, positioned 1.3 mm proximal to electrode LV2
- electrode LV4, the proximal electrode, positioned 21 mm proximal to electrode LV3

See Section 6.10, "Taking electrical measurements", page 10 for information about LV pacing polarity selections.

The Medtronic IS4-LLLL¹ four-pole inline connector on the lead facilitates device connection during implant. The tip of the connector pin has a white band indicator that can be used for visual confirmation of proper lead connection to the device.

The connector contacts align with the lead electrodes LV1 to LV1, LV2 to LV2, and so on:

- connector contact LV1, the connector pin, positioned at the proximal tip of the lead (aligns to electrode LV1)
- connector contact LV2, positioned distal to connector pin LV1 (aligns to electrode LV2)
- connector contact LV3, positioned distal to connector contact LV2 (aligns to electrode LV3)
- connector contact LV4, positioned distal to connector contact LV3 (aligns to electrode LV4)

¹ IS4-LLLL refers to an International Connector Standard (ISO 27186:2010) whereby pulse generators and leads so designated are assured of a basic mechanical fit. LLLL defines the lead connector contacts as low voltage (L).

² Technology developed by NASA.

See Section 7.2, "Specifications drawing (nominal)", page 15 for an illustration of the lead electrodes and connector contacts.

The distal tip of the lead allows a guide wire to pass through to aid in cardiac vein selection. The tip contains a silicone rubber membrane, which seals the lead inner lumen to reduce blood ingress.

Each electrode contains a Monolithic controlled release device (MCRD) for elution of steroid to reduce inflammatory responses within the cardiac vein. The MCRDs contain a combined total dosage of <1.0 mg of dexamethasone acetate steroid. Upon exposure to body fluids, the steroid elutes from the MCRDs. The steroid suppresses the inflammatory response that is believed to cause threshold rise typically associated with implanted pacing electrodes.

The electrodes are positioned relative to the 3 distal curves of the lead to facilitate contact with the cardiac veins.

The outer insulation of the lead is polyurethane and the inner insulation is SI-polyimide (SI-PI)². The SI-PI is applied as a coating to the conductor wire before coiling.

The Attain Performa S MRI SureScan 4598 LV lead can be positioned with the aid of a guide wire, a stylet, an inner catheter, or an inner catheter plus a guide wire or a stylet.

To implant the lead in a selected cardiac vein, a compatible delivery system is required. A compatible delivery system includes a guide catheter and either a hemostasis valve or an introducer valve that can be removed or that allows passage over the lead connector. Contact a Medtronic representative for further information regarding compatible delivery systems.

1.1 Medtronic SureScan system

The Model 4598 lead is part of the Medtronic SureScan system. The SureScan system includes a Medtronic SureScan device connected to Medtronic SureScan leads. Labeling for SureScan system components displays the SureScan logo and the MR Conditional symbol.

SureScan logo

MR Conditional symbol. The Medtronic SureScan system is MR Conditional and is designed to allow implanted patients to undergo an MRI scan under the specified MRI conditions for use.

The MRI SureScan feature permits a mode of operation that allows a patient with a SureScan device to be safely scanned by an MRI machine while the device continues to provide appropriate pacing. When programmed to On, MRI SureScan operation

disables arrhythmia detection, magnet mode, and all user-defined diagnostics. Before performing an MRI scan, refer to the SureScan MRI technical manual for important information about procedures and MRI-specific warnings and precautions.

1.2 Contents of package

Leads and accessories are supplied sterile. Each package contains the following items:

- 1 lead with anchoring sleeve
- 1 guide wire insertion tool
- 1 guide wire clip
- 1 guide wire steering handle
- 4 stylets
- 2 AccuRead 2.0 analyzer cable interface tools
- product documentation

1.3 Accessory descriptions

Dispose of all single-use accessories according to local environmental requirements.

AccuRead 2.0 tool – The AccuRead 2.0 analyzer cable interface tool guides the connection of analyzer cable terminals to the lead connector contacts. This connection facilitates accurate electrical measurements during implant and prevents possible connector damage.

Anchoring sleeve – An anchoring sleeve secures the lead to prevent it from moving and protects the lead insulation and conductors from damage caused by tight sutures.

Guide wire clip – A guide wire clip secures the excess guide wire and helps to protect and maintain the sterility of the guide wire.

Guide wire insertion tool – A guide wire insertion tool provides additional control when inserting a guide wire into the lead connector pin or the lead tip.

Guide wire steering handle – A guide wire steering handle is used only with guide wires 0.46 mm or less in diameter. The steering handle provides additional guide wire steering and rotation control.

Stylet – A stylet provides additional stiffness and controlled flexibility for maneuvering the lead into position. Each stylet knob is labeled with the stylet diameter and corresponding lead length.

2 Indications

The lead has application as part of a Medtronic biventricular pacing system.

3 Contraindications

Before performing an MRI scan, refer to the SureScan MRI technical manual for MRI-specific contraindications. The Model 4598 lead is contraindicated as follows:

Coronary vasculature – This lead is contraindicated for patients with coronary venous vasculature that is inadequate for lead placement, as indicated by venogram.

Steroid use – The lead is contraindicated in patients for whom a single dose of 1.0 mg of dexamethasone acetate may be contraindicated.

4 Warnings and precautions

A complete SureScan system is required for use in the MRI environment. Before performing an MRI scan, refer to the SureScan MRI technical manual for MRI-specific warnings and precautions.

Note: Medical procedure warnings and precautions that pertain to the Medtronic implanted system are provided in the manual that is packaged with the device or on the Medtronic Manual Library website (www.medtronic.com/manuals).

Inspecting the sterile package – Inspect the sterile package with care before opening it.

- If the seal or package is damaged, contact a Medtronic representative.
- Do not store this product above 40 °C.
- Do not use the product after its expiration date.

Sterilization – Medtronic has sterilized the package contents with ethylene oxide before shipment. This lead is for single use only and is not intended to be resterilized.

Single use – The lead and accessories are for single use only.

Necessary hospital equipment – Keep external defibrillation equipment nearby for immediate use during acute lead system testing, the implant procedure, or whenever arrhythmias are possible or intentionally induced during post-implant testing.

External defibrillation and cardioversion – External defibrillation and cardioversion are therapies that deliver an electrical shock to the heart to convert an abnormal heart rhythm to a normal rhythm.

Medtronic cardiac devices are designed to withstand exposure to external defibrillation and cardioversion. While damage to an implanted system from an external shock is rare, the probability increases with increased energy levels. These procedures may also temporarily or permanently elevate pacing thresholds or temporarily or permanently damage the myocardium. If external defibrillation or cardioversion are required, consider the following precautions:

- Use the lowest clinically appropriate energy.
- Position the patches or paddles a minimum of 15 cm away from the device.
- Position the patches or paddles perpendicular to the device and lead system.
- If an external defibrillation or cardioversion is delivered within 15 cm of the device, use a Medtronic programmer to evaluate the device and lead system.

Line-powered and battery-powered equipment – An implanted lead forms a direct current path to the myocardium.

4 English

During lead implant and testing, use only battery-powered equipment or line-powered equipment specifically designed for this purpose to protect against fibrillation that may be caused by alternating currents. Line-powered equipment used in the vicinity of the patient must be properly grounded. Lead connector pins must be insulated from any leakage currents that may arise from line-powered equipment.

Concurrent devices – Output pulses, especially from unipolar devices, may adversely affect device sensing capabilities. If a patient requires a separate stimulation device, either permanent or temporary, allow enough space between the leads of the separate systems to avoid interference in the sensing capabilities of the devices. Previously implanted pulse generators and implantable cardioverter defibrillators should generally be explanted.

Steroid use – It has not been determined whether the warnings, precautions, or complications usually associated with injectable dexamethasone acetate apply to the use of this highly localized, controlled-release device.

Pregnancy – Dexamethasone acetate has been shown to be teratogenic in many species when given in doses equivalent to the human dose. There are no adequate and well-controlled studies in pregnant women. Dexamethasone acetate should be used during pregnancy only if the potential benefit justifies the potential risk to the fetus. Studies in mice, rats, and rabbits have shown that adrenocorticoids increase the incidence of cleft palate, placental insufficiency, and spontaneous abortions, and can decrease the intrauterine growth rate.

Nursing mothers – Systemically administered corticosteroids appear in human milk and could suppress growth, interfere with endogenous corticosteroid production, or cause other untoward effects in nursing infants. Because of the potential for serious adverse reactions in nursing infants from corticosteroids, a decision should be made whether to discontinue nursing or to use a non-steroidal lead, taking into account the importance of the lead and the drug to the mother.

Handling steroid monolithic controlled release devices (MCRDs) – Avoid reducing the amount of steroid available before implanting the lead. Reducing the available amount of steroid may adversely affect low-threshold performance.

- Do not allow the electrode surfaces to come in contact with surface contaminants.
- Do not wipe or immerse the electrodes in fluid, except blood, at the time of implant.

Bipolar pacing – If the surface area of the selected anode electrode is equal to or less than the surface area of the selected cathode, higher pacing thresholds or anodal stimulation may result. This lead uses 4 equal-sized electrodes; therefore, bipolar pacing configurations may result in higher pacing thresholds or anodal stimulation.

If therapy cannot be delivered via the LV1 electrode, the LV2, LV3, and LV4 electrodes are available for use with specific Medtronic devices. Refer to the appropriate Medtronic Cardiac Resynchronization Therapy (CRT-D) system manual for use of available LV lead pacing polarity options.

Handling the stylet – Handle the stylet with care at all times.

- To minimize the likelihood of trauma to the vein and to maintain lead flexibility while advancing the lead through the vein, keep the stylet withdrawn 1 to 2 cm or select a more flexible stylet.
- Do not use excessive force or surgical instruments when inserting a stylet.
- Avoid overbending, kinking, or blood contact on stylets.
- Use a new stylet when blood or other fluids accumulate on the stylet. Accumulated fluids may cause damage to the lead or difficulty in passing the stylet through the lead.
- Curving the distal end of the stylet prior to insertion into the lead will achieve a curvature at the distal end of the lead. Do not use a sharp object to impart a curve to the distal end of the stylet.

Handling the guide wire – Handle the guide wire with care at all times.

- Do not insert the proximal end of the guide wire through the lead tip seal without using the guide wire insertion tool. Inserting the guide wire without the guide wire insertion tool may damage the lead.
- Damage to a guide wire may prevent the guide wire from performing with accurate torque response and control and may cause vessel damage. For additional information about vessel damage and other potential adverse events, refer to the technical manual packaged with the appropriate guide wire.
- If the distal end of the guide wire becomes severely kinked or twisted, it may be difficult to withdraw it back through the lead. Therefore, if there is an indication that the distal end of the guide wire has become damaged, or if there is significant resistance in guide wire passage, remove the lead and guide wire together as a unit. Remove the guide wire from the lead and insert a new guide wire into the lead. Do not use excessive force to retract the guide wire from the lead. Refer to the product documentation packaged with the guide wire for additional information.

Handling the lead – Handle the lead with care at all times.

- Rust stylets are not recommended with this lead due to the risk of conductor coil or insulation perforation.
- If a stylet is used for lead positioning, use only the stylets packaged with the lead or in a stylet kit (downsized knob). Other stylets may extend beyond the lead tip causing patient injury.
- If the lead is damaged, do not implant it. Return the lead to a Medtronic representative.
- Protect the lead from materials that shed particles such as lint and dust. Lead insulators attract these particles.
- Handle the lead with sterile surgical gloves that have been rinsed in sterile water or a comparable substance.
- Do not severely bend, kink, or stretch the lead.
- Do not use surgical instruments to grasp the lead or connector pin.
- Do not immerse leads in mineral oil, silicone oil, or any other liquid, except blood, at the time of implant.

English 5

- Use an anchoring sleeve with all leads. Ensure that the anchoring sleeve is positioned close to the lead connector pin, to prevent inadvertent passage of the sleeve into the vein. If it is necessary to wipe the lead before insertion, ensure that the anchoring sleeve remains in position.
- Do not force the guide catheter or leads if significant resistance is encountered. Use of guide catheters or leads may cause trauma to the heart.

Chronic repositioning or removal – Chronic repositioning or removal of leads may be difficult because of fibrotic tissue development. The clinical study was not designed to evaluate the removal of left ventricular leads from the coronary venous vasculature. If a lead must be removed or repositioned, proceed with extreme caution. Return all removed leads to Medtronic.

- Rust stylets are not recommended with this lead due to the risk of conductor coil/insulation perforation.
- Verify lead length on the lead label on the connector to choose an appropriate stylet kit (downsized knob) length when repositioning. Always choose a stylet kit (downsized knob) 3 cm shorter than the lead length. For example, choose a stylet kit (downsized knob) with stylets 75 cm long for a lead 78 cm long.
- Lead removal may result in avulsion of the endocardium, valve, or vein.
- Lead junctions may separate, leaving the lead tip and bare wire in the heart or vein.
- Chronic repositioning may adversely affect the low-threshold performance of a steroid-eluting lead.
- Cap abandoned leads to avoid transmitting electrical signals.
- For leads that have been severed, seal the remaining lead end and suture the lead to adjacent tissue.
- If a lead is removed and repositioned, inspect it carefully for insulator or conductor coil damage before repositioning.

Magnetic resonance imaging (MRI) – An MRI is a type of medical imaging that uses magnetic fields to create an internal view of the body. If certain criteria are met and the warnings and precautions provided by Medtronic are followed, patients with an MR Conditional device and lead system are able to undergo an MRI scan; for details, refer to the SureScan MRI technical manual that Medtronic provides for an MRI Conditional device.

Diathermy treatment (including therapeutic ultrasound) – Diathermy is a treatment that involves the therapeutic heating of body tissues. Diathermy treatments include high frequency, short wave, microwave, and therapeutic ultrasound. Except for therapeutic ultrasound, do not use diathermy treatments on cardiac device patients. Diathermy treatments may result in serious injury or damage to an implanted device and lead system. Therapeutic ultrasound (including physiotherapy, high intensity therapeutic ultrasound, and high intensity focused ultrasound), is the use of ultrasound at higher energies than diagnostic ultrasound to bring heat or agitation into the body. Therapeutic ultrasound is acceptable if treatment is performed with a minimum separation distance of 15 cm between the applicator and the implanted device and lead system, as long as the ultrasonic beam is pointing away from the device and lead system.

5 Potential adverse events

The potential adverse events related to the use of transvenous leads include, but are not limited to, the following patient-related conditions:

- air embolism
- cardiac dissection
- cardiac perforation
- cardiac tamponade
- coronary sinus dissection
- death
- endocarditis and pericarditis
- erosion through the skin
- extracardiac muscle or nerve stimulation
- fibrillation or other arrhythmias
- heart block
- heart wall or vein wall rupture
- hematoma/seroma
- infection
- myocardial irritability
- myopotential sensing
- pericardial effusion
- pericardial rub
- pneumothorax
- body rejection phenomena including local tissue reaction
- threshold elevation
- thrombosis
- thrombotic embolism
- valve damage (particularly in fragile hearts)

Additional potential adverse events related to the lead and the programmed parameters include, but are not limited to, the following:

Potential adverse event	Indicator of potential adverse event	Corrective actions to consider
Lead dislodgement ^a	Intermittent or continuous loss of capture or LV EGM signal integrity (including sensing) ^b	Reprogram the LV pacing polarity. Reposition the lead.
Lead dislodgement ^a	Intermittent or continuous oversensing	Reprogram the LV pacing polarity. Reposition the lead.
Lead conductor fracture	Intermittent or continuous loss of capture or LV EGM signal integrity (including sensing) ^a	Replace the lead. Reprogram the LV pacing polarity.

6 English

Potential adverse event	Indicator of potential adverse event	Corrective actions to consider
Lead conductor insulation failure	Intermittent or continuous loss of capture or LV EGM signal integrity (including sensing) ^a	Replace the lead. Reprogram the LV pacing polarity.
Threshold elevation or exit block	Loss of capture ^a	Adjust the implantable device output. Reprogram the LV pacing polarity. Replace or reposition the lead.

^a Transient loss of capture or LV EGM signal integrity (including sensing) may occur following surgery until lead stabilization takes place. If stabilization does not occur, lead dislodgement may be suspected.

Implant techniques that may damage the lead include, but are not limited to, the following techniques:

Implant techniques that may damage the lead	Possible effects on the lead	Corrective action to consider
Forcing the lead through the introducer/delivery system	Electrode, conductor coil, or insulation damage	Replace the lead.
Use of too medial of an approach with venous introducer resulting in clavicle and first rib binding	Conductor coil fracture, insulation damage	Replace the lead.
Using too stiff a stylet	Conductor coil/insulation perforation	Replace the lead.
Puncturing the pericardium or tendon when using subclavian introducer approach resulting in binding	Conductor coil fracture, insulation damage	Replace the lead.
Advancing the lead through the non-coronary central access veins without the stylet or guide wire fully inserted	Tip distortion or insulation perforation	Replace the lead.
Inserting the proximal end of the guide wire through the lead tip seal without using the guide wire insertion tool	Lead tip seal damage or conductor coil/insulation damage	Replace the lead.

6 Implant procedure

Warning: Before implanting a SureScan system, consider the risks associated with removing previously implanted leads. Abandoned leads or previously implanted leads not tested for MRI compatibility compromise the ability to safely scan the SureScan system during MRI scans.

Warning: Do not force the guide catheter or lead if significant resistance is encountered. Use of guide catheters or leads may cause trauma to the heart.

To implant the Attain Performa S MRI SureScan 4598 LV lead in a selected cardiac vein, a compatible delivery system is required, such as a Medtronic delivery system. A compatible delivery system includes a guide catheter and either a hemostasis valve or an introducer valve that can be removed or that allows passage over the lead connector. Contact a Medtronic representative for further information regarding compatible delivery systems.

Proper surgical procedures and sterile techniques are the responsibility of the medical professional. The implant procedures described in this manual are provided for information only. Each physician must apply the information in these instructions according to professional medical training and experience.

6.1 Placing the right ventricular lead

When deciding which ventricular lead to place first, consider the ease of coronary sinus cannulation and the need for backup pacing.

6.2 Preparing the delivery system


Prepare the delivery system for lead implantation according to the instructions in the product documentation packaged with the delivery system.

6.3 Venous access

Warning: Backup pacing should be readily available during implant. Use of the delivery system or leads may cause heart block.

- To access the subclavian vein, use a preferred method based on professional experience.
Caution: Certain anatomical abnormalities, such as thoracic outlet syndrome, may precipitate pinching and subsequent fracture of the lead.
Caution: Insertion should be done as far lateral as possible to avoid clamping the lead body between the clavicle and the first rib (Figure 1).

Figure 1.

1 Suggested entry site

English 7

2. Introduce a J-shaped introducer guide wire and percutaneous introducer sheath.
3. Introduce the guide catheter assembly to access the coronary sinus.

See the delivery system product documentation for additional information.

6.4 Obtaining venograms

Before placing a lead in the coronary sinus, obtain venograms. Venograms are recommended to assess a route for passage and a site for final placement based on the size, shape, location, and tortuosity of the veins. Also, venograms may be useful in identifying suspected coronary sinus trauma. For information on obtaining a venogram by using a venogram balloon catheter, see the product documentation packaged with an appropriate venogram balloon catheter.

6.5 Inserting the lead into the delivery system

Warning: If a stylet is used for lead positioning, use only the stylets packaged with the lead or in a stylet kit (downsized knob). Always use a stylet that is 3 cm shorter than the lead length listed on the IS4 connector label. Other stylets may extend beyond the lead tip causing injury or perforation of the cardiac vein or heart.

Warning: Do not force the lead if significant resistance is encountered during lead passage. The use of guide catheters or leads may cause trauma to the heart.


Caution: Use care when handling the lead during insertion.

- Do not severely bend, kink, or stretch the lead.
- Do not use surgical instruments to grasp the lead or connector pin.

The following steps may be used to insert the lead:

1. Insert a straight stylet or a guide wire into the lead to vary the shape of the distal end of the lead (Figure 2).
Note: When a stylet is fully inserted, the distal tip of the stylet does not reach the distal tip of the lead.

Figure 2.

- 1 Stylet fully withdrawn
- 2 Stylet fully inserted


2. Insert the lead into the delivery system. See the delivery system product documentation for additional information.

6.6 Lead placement

The Attain Performa S MRI SureScan 4598 LV lead can be placed with the aid of a guide wire, a stylet, an inner catheter, or an inner catheter plus a guide wire or a stylet.

Stylet delivery – If the patient's anatomy features gentle vein angulation off the coronary sinus and the cardiac vein branch is not tortuous (Figure 3), a stylet may be used for lead delivery.


Figure 3.

- 1 Coronary sinus
- 2 Cardiac vein

Guide wire delivery – If the patient's anatomy features acute vein angulation off the coronary sinus and the cardiac vein branch is tortuous (Figure 4), a guide wire may be used for lead delivery.

Figure 4.

- 1 Tortuous cardiac vein branch with gentle angulation from the coronary sinus
- 2 Tortuous cardiac vein branch with acute angulation from the coronary sinus

6.7 Placing the lead using a stylet

Warning: Do not force the lead if significant resistance is encountered during lead passage.

Warning: To minimize the likelihood of trauma to the vein and to maintain lead flexibility while advancing the lead through the vein, keep the stylet withdrawn 1 to 2 cm or select a more flexible stylet.

Caution: Do not use a sharp object to impart a curve to the distal end of the stylet. Imparting a curve to the stylet can be accomplished with a smooth-surface, sterile instrument (Figure 5).

Figure 5.

Note: If it is difficult to advance the stylet around a bend, consider using a different stylet. More flexible stylets are recommended for tortuous anatomies. Firmer stylets are recommended where additional support is needed.

There are many techniques that may be used to advance the lead into a cardiac vein using a stylet. The choice of technique is left to the discretion of the physician.

6.8 Placing the lead using a guide wire

Warning: Do not insert the proximal end of the guide wire through the lead tip seal without using the guide wire insertion tool.

Inserting the guide wire without the guide wire insertion tool could damage the lead tip seal or the conductor coil or insulation.

Warning: Damage to a guide wire may prevent the guide wire from performing with accurate torque response and control and may cause vessel damage. For additional information about vessel damage and other potential adverse events, refer to the technical manual packaged with the appropriate guide wire.

Caution: Use care when positioning the guide wire. Refer to the product documentation packaged with the guide wire for additional information.


Notes:

- Medtronic recommends using guide wires 0.36 mm to 0.46 mm in diameter. Contact a Medtronic representative for further information about recommended guide wires.
- Consider soaking the guide wire in a heparin solution before insertion to minimize the risk of thrombus formation during use.
- If you are using an Attain Hybrid guide wire, the procedure for preparing the guide wire for use differs from other guide wires because of the attached proximal knob.

Use the following steps to prepare the guide wire for use:

1. Select a guide wire.
 If the patient has tortuous anatomy, a more flexible guide wire is recommended. If additional support is needed, use a firmer guide wire.
 If you are using an Attain Hybrid guide wire, the steps involving the insertion tool (Step 2 through Step 4) do not apply.
2. Insert the guide wire into the lead by placing the distal (flexible) end of the guide wire into the lead connector pin using the guide wire insertion tool included in the package (Figure 6). To prevent the lead from dislodging from the guide wire insertion tool, grasp the lead and guide wire insertion tool firmly between thumb and forefinger.
Caution: To minimize the risk of damaging the guide wire, be sure that the flexible section of the guide wire is fully inserted into the lead before removing the guide wire insertion tool from the lead.
Note: Be sure to remove the guide wire insertion tool before lead implant.

Figure 6.


- 1 Guide wire
- 2 Lead connector pin

3. Disengage the guide wire insertion tool from the lead connector pin.

English 9

4. Remove the guide wire insertion tool by sliding the tool off the end of the guide wire.
5. Position the guide wire steering handle (Figure 7):
 For an Attain Hybrid guide wire, preload the guide wire steering handle onto the guide wire before loading the guide wire into the connector pin on the proximal end of the lead.
 For other guide wires:
 - a. Advance the guide wire steering handle over the proximal (rigid) end of the guide wire.
 - b. Tighten the guide wire steering handle onto the guide wire near the lead connector pin.

Figure 7.

- 1 Guide wire steering handle
- 2 Lead connector pin

6. Attach the guide wire clip to the guide wire and secure it within the sterile field. Medtronic recommends securing the guide wire clip to the patient's sterile surgical drape.


As an alternate approach in situations where the guide wire is already in place, the lead can be loaded over the guide wire using the guide wire insertion tool.

Insert the guide wire into the lead by placing the proximal (rigid) end of the guide wire into the distal lead tip using the guide wire insertion tool included in the package (Figure 8).

Notes:

- There may be slight resistance as the guide wire passes through the lead tip seal.
- Grasp the lead and guide wire insertion tool between your thumb and forefinger to prevent the lead from dislodging from the guide wire insertion tool while you are inserting the guide wire.
- Be sure to remove the guide wire insertion tool, through the slit in the guide wire insertion tool, before lead implant.

Figure 8.

- 1 Guide wire
- 2 Lead distal tip

Warning: Do not force the lead if significant resistance is encountered during lead passage.

Caution: If the distal end of the guide wire becomes severely kinked or twisted, it may be difficult to withdraw it back through the

lead. Therefore, if there is an indication that the distal end of the guide wire has become damaged, or if there is significant resistance in guide wire passage, remove the lead and guide wire together as a unit. Remove the guide wire from the lead and insert a new guide wire into the lead. Do not use excessive force to retract the guide wire from the lead.

Notes:

- If the lead is not advancing, or if the lead and guide wire seem to be sticking together, there may be thrombus on the guide wire at the lead tip. Remove and inspect the lead and guide wire. Consider using a new guide wire. Reinsert the lead and the guide wire as described in previous steps.
- If it is difficult to advance the guide wire around a bend, consider using a different guide wire. More flexible guide wires are recommended for tortuous anatomies. Firmer guide wires are recommended where additional support is needed.


There are many techniques that may be used to advance the lead into a cardiac vein using a guide wire. The choice of technique is left to the discretion of the physician.

6.9 Fixating the lead

Fixation is achieved the same way regardless of lead placement.

Using fluoroscopy for guidance, fixate the lead by wedging the 3 curves of the lead into the cardiac vein (Figure 9).

Figure 9.

For optimal electrical performance and fixation, place the proximal electrode in the selected vein, not in the coronary sinus.


Note: If the selected vein is large, it may be necessary to position the lead in a smaller cardiac vein in order to achieve lead fixation. For optimal electrical performance, Medtronic recommends electrode-tissue contact.

6.10 Taking electrical measurements

The Attain Performa S MRI SureScan 4598 LV lead was designed to provide pacing via selectable electrodes. The 16 available pacing polarities are shown in Figure 10.

10 English

Figure 10.

- 1 Extended bipolar pacing polarities
- 2 Bipolar (reversible) pacing polarities

Extended bipolar pacing is available using this lead in combination with a compatible CRT-D system and RV defibrillation lead. The 4 available pacing polarities include:

- LV1 to RVcoil
- LV2 to RVcoil
- LV3 to RVcoil
- LV4 to RVcoil

Bipolar pacing is available using this lead in combination with a compatible CRT-D system. The polarity of each of the 6 electrode pairs can be reversed to yield a total of 12 bipolar pacing polarities:

- LV1 to LV2, LV2 to LV1
- LV1 to LV3, LV3 to LV1
- LV1 to LV4, LV4 to LV1
- LV2 to LV3, LV3 to LV2
- LV2 to LV4, LV4 to LV2
- LV3 to LV4, LV4 to LV3

Note: If you are monitoring the LV EGM signal when the LV2 to LV3 or LV3 to LV2 pacing polarity is selected, the signal amplitude may be attenuated compared to other pacing configurations such as LV1 to RVcoil. This signal attenuation is an expected characteristic of the Attain Performa S MRI SureScan 4598 LV lead with the short spacing between electrodes LV2 and LV3.

Caution: Before taking electrical or defibrillation efficacy measurements, move objects made of conductive materials, such as guide wires or stylets, away from all electrodes.

Note: Before taking electrical measurements, it is recommended that you retract the stylet or guide wire, inside the lead lumen, to a point proximal of all electrodes. This action allows the lead tip to

resume its normal shape, enabling appropriate electrode-tissue contact to occur.


The AccuRead 2.0 analyzer cable interface tool is used to facilitate accurate electrical measurements during implant.

Caution: Use the AccuRead 2.0 cable interface tool when connecting analyzer cable terminals to the lead connector. This tool enables you to take accurate electrical measurements at implant while reducing the risk of connector damage, electrical bridging, or electrical shorting. The potential for connector damage, bridging, or shorting is due to variations in analyzer cable terminals and due to the width and proximity of the contacts (rings and pin) on the IS4 connector.

Use the following steps to take electrical measurements:

1. Ensure that the lead connector is inserted into the AccuRead 2.0 tool completely. If the AccuRead 2.0 tool is properly attached, the connector pin is accessible (see Figure 11).

Figure 11.

2. Attach a surgical cable to the lead connector. Use the AccuRead 2.0 cable interface tool to guide the alignment of the cable clips with the contacts on the lead connector. The tool helps to ensure that accurate readings are obtained. See the lead drawing in Chapter 7 for information about the alignment of the lead connector contacts with the lead electrodes.
3. Use an implant support instrument, such as a pacing system analyzer, for obtaining electrical measurements. For information on the use of the implant support instrument, consult the product documentation for that device. Satisfactory lead placement is indicated by low stimulation thresholds and adequate sensing of intracardiac signal amplitudes.

Notes:

- A low stimulation threshold provides a desirable safety margin, allowing for a possible rise in thresholds that may occur within 2 months following implant.
- Adequate acute LV EGM signal amplitudes ensure that the lead is properly sensing intrinsic cardiac signals. Minimum signal requirements depend on the sensitivity capabilities of the device. Acceptable acute signal amplitudes for the lead must be greater than the minimum device sensing capabilities. Be sure to include an adequate safety margin to account for lead maturity.³

Table 1. Recommended measurements at implant^a

Measurement recommended	Left ventricle
Maximum acute stimulation threshold ^b	3.0 V
Minimum acute LV EGM signal amplitude	4.0 mV ^c


^a Assuming 500 Ω resistance.

^b At a pulse width setting of 0.5 ms.

^c The LV EGM signal amplitude for the LV2 to LV3 and LV3 to LV2 pacing polarities may be attenuated compared to other pacing polarities such as LV1 to RVcoil.

4. If electrical measurements do not stabilize to acceptable levels, it may be necessary to reposition the lead and repeat the testing procedure.
5. Check for phrenic nerve stimulation by pacing at 10 V and a pulse width setting greater than 0.5 ms. Then observe for diaphragmatic contracting either by fluoroscopy or direct abdominal palpitation. Further testing may include patient positional changes to simulate upright chronic conditions. If phrenic nerve stimulation occurs, reduce the voltage until a phrenic nerve stimulation threshold is determined. Phrenic nerve stimulation may necessitate changing the pacing polarity or repositioning of the lead.
6. The AccuRead 2.0 tool may be removed from a guide wire or stylet using the slit on the side of the tool (see Figure 12).

Figure 12.

- 1 Removing the AccuRead 2.0 tool from the connector pin
- 2 Removing the AccuRead 2.0 tool from the stylet using the slit on the side of the tool

6.11 Removing the guide catheter from the lead

Note: If an Attain Hybrid guide wire was used for placement of the lead, go to Step 2. It is not necessary to replace the Attain Hybrid guide wire with a stylet for catheter removal.

Once the lead is in the final position, remove the guide catheter from the lead:

1. If used, remove the guide wire and guide wire insertion tool. Replace the guide wire with a straight stylet (downsized knob). Insert the straight stylet into the lead to the mid-coronary sinus.
2. Remove the guide catheter from the lead. See the delivery system product documentation for details.
Note: For Medtronic slidable delivery systems, use a sifter compatible with a 1.75 mm (5.3 French) lead body.

³ The LV EGM signal amplitude for the LV2 to LV3 and LV3 to LV2 pacing polarities may be attenuated compared to other pacing polarities such as LV1 to RVcoil.

3. Carefully and completely remove the stylet or Attain Hybrid guide wire. When removing the stylet, grip the lead firmly close to the distal end of the connector pin; gripping the lead at this location helps prevent possible lead tip dislodgement.
4. Ensure that the lead connector is inserted into the AccuRead 2.0 analyzer cable interface tool and then repeat the electrical measurements. See Section 6.10, "Taking electrical measurements", page 10.

6.12 Anchoring the lead

Caution: Use care when anchoring the lead.

- Use only nonabsorbable sutures to anchor the lead.
- Do not attempt to remove or cut the anchoring sleeve.
- Do not use the anchoring sleeve tabs for suturing.
- During lead anchoring, take care to avoid dislodging the lead tip.
- Do not secure sutures so tightly that they damage the vein, lead, or anchoring sleeve (Figure 13).
- Do not tie a suture directly to the lead body (Figure 13).

Figure 13.

Anchor the lead using all 3 grooves:

1. Position the anchoring sleeve against or near the vein.
2. Secure the anchoring sleeve to the lead body by tying a suture firmly in each of the 3 grooves (Figure 14).

Figure 14.

3. Use at least 1 additional suture in 1 of the grooves to secure the anchoring sleeve and lead body to the fascia.

6.13 Connecting the lead

Caution: Before connecting the lead to the device, always remove implant tools such as a stylet, guide wire, or AccuRead 2.0 cable interface tool. Failure to remove implant tools may result in lead failure.

Connect the lead to an implantable device.


1. Insert the lead connector into the device connector block. Consult the product documentation packaged with the implantable device for instructions on proper lead connections.
2. Take electrical measurements through the device.

6.14 Placing the device and lead into the pocket

Caution: Use care when placing the device and leads into the pocket.

- Ensure that the leads do not leave the device at an acute angle.
- Do not grip the lead or device with surgical instruments.
- Do not coil the lead. Coiling the lead can twist the lead body and may result in lead dislodgement (Figure 15).

Figure 15.

Use the following steps to place the device and leads into the pocket:

1. To prevent undesirable twisting of the lead body or bending of the lead body at an acute angle, rotate the device to wrap the excess lead length loosely (Figure 16).

Figure 16.

2. Insert the device and leads into the pocket.
3. Before closing the pocket, verify sensing, pacing, cardioversion, and defibrillation efficacy.

6.15 Post-implant evaluation

After implant, monitor the patient's electrocardiogram until the patient is discharged. If a lead dislodges, it usually occurs during the immediate postoperative period.

Recommendations for verifying that the lead is properly positioned include x-rays and pacing thresholds taken at pre-hospital discharge, 3 months after implant, and every 6 months thereafter.

In the event of a patient death, explant all implanted leads and devices and return them to Medtronic with a completed Product Information Report form.

7 Specifications

7.1 Specifications (nominal)


Parameter	4898
Serial number prefix	QUC
Type	Quadrupole electrode
Chamber paced	Left ventricle
Length	78 and 88 cm
Connector	IS4-LLLL
Material	Conductor: 25% Ag-core-MP35N Insulators: Polyurethane (outer) SI-polyimide (SI-Pi) (inner) ^a Electrodes: Platinum iridium alloy with titanium nitride coating Connector pin: MP35N Connector rings: MP35N Molded tip seal: Silicone rubber
Electrode configuration	Radial, titanium-nitride-coated, steroid-eluting
Diameter	Lead body: 1.75 mm (5.3 French) Electrodes: 1.70 mm (5.1 French) Spacer between electrodes LV1 (distal) and LV2: 1.30 mm (3.9 French) Spacer between electrodes LV2 and LV3: 1.57 mm (4.7 French) Spacer between electrodes LV3 and LV4 (proximal): 1.30 mm (3.9 French) Medtronic delivery system (recommended inner diameter): 1.90 mm (5.7 French) Diagnostic guide wire (recommended diameter): 0.36 mm to 0.46 mm Electrode surface area: 5.8 mm ² Distance between electrodes: Electrode LV1 (distal) to LV2: 21 mm Electrode LV2 to LV3: 1.3 mm

Parameter	4898
	Electrode LV3 to LV4 (proximal): 21 mm
Conductor resistance	LV1: 22 ±5 Ω (78 cm) 24 ±6 Ω (88 cm) LV2: 19 ±4 Ω (78 cm) 21 ±5 Ω (88 cm) LV3: 18 ±4 Ω (78 cm) 21 ±4 Ω (88 cm) LV4: 17 ±4 Ω (78 cm) 20 ±4 Ω (88 cm)
Steroid	Dexamethasone acetate
Combined total of steroid	<1.0 mg
	Steroid binder: Silicone rubber

^aTechnology developed by NASA.

7.2 Specifications drawing (nominal)

Figure 17.

- 1 Electrode LV1 (the distal electrode): nominally 5.8 mm² geometric pacing surface area
- 2 Electrode LV2: nominally 5.8 mm² geometric pacing surface area
- 3 Electrode LV3: nominally 5.8 mm² geometric pacing surface area
- 4 Electrode LV4 (the proximal electrode): nominally 5.8 mm² geometric pacing surface area
- 5 Anchoring sleeve

- 6 IS4 connector
- 7 Contact LV4
- 8 Contact LV3
- 9 Contact LV2
- 10 Contact LV1

8 CE mark of conformity

The following CE mark of conformity applies to the Model 4598 lead.

CE0123
2013

9 Explanation of symbols on package labeling

Refer to the package labels to see which symbols apply to this product.

Table 2. Explanation of symbols on package labeling

Symbol	Explanation
CE0123	Conformité Européenne (European Conformity). This symbol means that the device fully complies with applicable European Union acts.
	Do not use if package is damaged
	Do not reuse
	Upper limit of temperature
STERILE EO	Sterilized using ethylene oxide
	Consult instructions for use
	Consult instructions for use at this website
	Date of manufacture
	Manufacturer
	Use by

English 15

Table 2. Explanation of symbols on package labeling (continued)

Symbol	Explanation
	Reorder number
	Authorized representative in the European Community
	Serial number
	Open here
	Package contents
	Product documentation
	Caution
	Accessories
	Inner diameter
	Lead length
	Catheter
	Catheter delivered OTW
	Lead
	Transvenous cardiac vein pacing lead
	Transvenous ventricular lead
	Transvenous atrial lead
	Transvenous lead with one defibrillation electrode

Table 2. Explanation of symbols on package labeling (continued)

Symbol	Explanation
	Transvenous lead with two defibrillation electrodes
	Pace
	Sense
	Defibrillation
	Extendable and retractable screw-in
	Tined
	Steroid-eluting
	Lead introducer
	Lead introducer with guide wire
	J-shaped
	SureScan logo
	MR Conditional

10 Medtronic disclaimer of warranty

For complete disclaimer of warranty information, see the accompanying disclaimer of warranty document.

11 Service

Medtronic employs highly trained representatives and engineers located throughout the world to serve you and, upon request, to provide training to qualified hospital personnel in the use of Medtronic products. Medtronic also maintains a professional staff to provide technical consultation to product users. For more information, contact your local Medtronic representative, or call or write Medtronic at the appropriate telephone number or address listed on the back cover.

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64 Srpski

Printing instructions: doc #163256
Category: Leads single pkg.
Size: 7.125 x 8 in (181 x 203 mm)

502115-008

M958924A028 A Medtronic Confidential

Composed: 2015-05-01 07:18:06
XSL-Stylesheet C 18-MAR-2015

Medtronic

Medtronic, Inc.
710 Medtronic Parkway
Minneapolis, MN 55432
USA
www.medtronic.com
+1 763 514 4000

**Authorized Representative in the European
Community**
Medtronic B.V.
Earl Bakkenstraat 10
6422 PJ Heerlen
The Netherlands
+31 45 566 8000

Europe/Middle East/Africa
Medtronic International Trading Sàrl
Route du Molliat 31
Case Postale 84
CH-1131 Tolochenaz
Switzerland
+41 21 802 7000

Australia
Medtronic Australasia Pty Ltd
97 Waterloo Road
North Ryde, NSW 2113
Australia

Canada
Medtronic of Canada Ltd
99 Hereford Street
Brampton, Ontario L6Y 0R3
Canada
+1 905 460 3800

Technical manuals
www.medtronic.com/manuals

© 2015 Medtronic, Inc.
M958924A028 A
2015-05-01

M958924A028

Appendix technical manual

When the device is delivered to the Russian Federation, the information contained in this appendix to the technical manual is a priority

1. Name of medical product	Cardiac vein pacing lead Attain Performa. Cardiac vein pacing lead ATTAIN PERFORMA of various version: Cardiac vein pacing lead ATTAIN PERFORMA S MRI SURESCAN 4598 size: 78 cm, 88 cm, composed: • lead with anchoring sleeve; • guide wire insertion tool; • guide wire clip; • guide wire steering handle; • stylets - no more than 4 pcs; • AccuRead 2.0 analyzer cable interface tools - no more than 2 pcs; • product documentation: Technical Manual; Limited Warranty; Stickers.
2. Information about manufacturer of medical device (name, location address, mailing address, phone number, address of place of production)	Design Site: Medtronic Inc., USA 8200 Coral Sea St., Mounds View, Minnesota 55112 Tel.: +1 800-633-8766; www.medtronic.com Legal Manufacturer: Medtronic Inc., USA 710 Medtronic Parkway, Minneapolis, Minnesota 55432, USA Tel.: +1 800-633-8766 www.medtronic.com Authorized Representative of Manufacturer: Medtronic, LLC, 123112, Moscow, in.ter.c. Presnensky Municipal District, Presnenskaya nab., 10, floor 9, office III, room 41 Tel.: (495) 580-73-77, Fax: (495) 580-73-78, www.medtronic.ru Holder of Registration Certificate in Russia (name, address, tel, fax, e-mail, website): Medtronic, LLC,

	123112, Moscow, in.ter.c. Presnensky Municipal District, Presnenskaya nab., 10, floor 9, office III, room 41 Tel.: (495) 580-73-77, Fax: (495) 580-73-78, www.medtronic.ru Address of Manufacturing Site for medical device: Medtronic Puerto Rico Operations Co., Villalba Rd. 149, KM. 56.3, Call Box 6001, Villalba, PR 00766, USA
3. Prescription of a medical device	Intended for left ventricular stimulation through the vein of the heart in the left ventricle, as part of the biventricular pacing system.
4. Information on potential consumers of medical device	This medical device is intended for use in clinical practice by qualified physicians and medical specialists.
5. Sterilization method, Sterilization Validation, Packaging Validation for Sterile devices	Sterilization – Medtronic has sterilized the package contents with ethylene oxide before shipment. This lead is for single use only and is not intended to be resterilized. Device is sterile and non-pyrogenic. - Sterilized by EO. - Sterility Assurance Level (SAL)- 10^{-6}.
6. Procedure and conditions for the disposal or destruction of a medical device	Dispose of all single-use accessories according to local environmental requirements. After use, this product may become biohazardous. Unused products with expired shelf life are classified as Class A medical waste. Products extracted from the patient's body are classified as Class B medical waste.
7. Class risk depending on the potential risk of a medical product, in accordance with the nomenclature	Classification of medical device: class risk 3; B - products whose failure reduces the effectiveness or delays the treatment and diagnostic process in non-critical situations, or increases the burden on medical or maintenance personnel; After use, this product may be a potential biohazard; Unused, expired devices are classified as medical waste class «A»; Leads extracted from the patient's body are medical waste Class B; Class of electrical safety, laser safety, and other types of hazards: not applicable; Type of working part: not applicable; The device is sterile, not re-sterilized.
8. Methods and apparatus for	Not applicable. Sterile. For single use only.

disinfection and sterilization cleaning									
9. Shelf life	Shelf life no more 24 months from the date of sterile packaging.								
10. Operating conditions	In a medical and preventative treatment facility. The device should be used only by or under the supervision of a trained physician. <table> <tr> <td></td><td>Service conditions</td></tr> <tr> <td>Operating temperature</td><td>from + 32°C to + 42°C</td></tr> <tr> <td>Humidity rate</td><td>80%-100%</td></tr> <tr> <td>Atmospheric pressure</td><td>from 700 hPa to 1060 hPa</td></tr> </table>		Service conditions	Operating temperature	from + 32°C to + 42°C	Humidity rate	80%-100%	Atmospheric pressure	from 700 hPa to 1060 hPa
	Service conditions								
Operating temperature	from + 32°C to + 42°C								
Humidity rate	80%-100%								
Atmospheric pressure	from 700 hPa to 1060 hPa								
11. Conditions of transportation and storage	Storage temperature – Store at 25°C (77°F). Excursions from this storage temperature are permitted in the range of 15 to 30°C (59 to 86°F). Transient spikes up to 40°C (104°F) are permitted as long as they do not exceed 24 hours.								
12. Warranty	Disclaimer of warranty and general warning Although the Medtronic lead, hereafter referred to as “Product”, has been carefully designed, manufactured and tested prior to sale, the Product may fail to perform its intended function satisfactorily for a variety of reasons. Warnings contained in the Product labeling provide more detailed information and are considered an integral part of this disclaimer of warranty. Medtronic, therefore, disclaims all warranties, both express and implied, with respect to the Product. Medtronic shall not be liable for any incidental or consequential damages caused by any use, defect, or failure of the Product, whether the claim is based on warranty, contract, tort, or otherwise. The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this disclaimer of warranty is held by any court of competent jurisdiction to be illegal, unenforceable, or in any conflict with applicable law, the validity of the remaining portion of the disclaimer of warranty shall not be affected, and all rights and obligations shall be construed and enforced as if this disclaimer of warranty did not contain the particular part or term held to be invalid. General warning Medtronic implantable leads are implanted in the extremely hostile environment of the human body. Leads are necessarily very small in diameter and must still be very flexible, which unavoidably reduces their potential performance or longevity. Leads may fail to function for a variety of causes, including, but not limited to: Medical complications, body rejection phenomena, allergic reaction, fibrotic tissue, or failure of leads by breakage or by breach of their insulation covering. In addition, despite the exercise of all due care in design, component								

**Medtronic**

Medtronic Inc.
710 Medtronic Parkway N.E.
Minneapolis, MN 55432
USA

	selection, manufacture, and testing prior to sale, leads may be easily damaged before, during, or after insertion by improper handling or other intervening acts. Consequently, no representation or warranty is made that failure or cessation of function of leads will not occur; or that the body will not react adversely to the implantation of leads; or that medical complications (including perforation of the heart) will not follow implantation of leads; or that the lead will, in all cases, restore adequate cardiac function. Warranty shelf life: The warranty covers the entire shelf life of the medical device.
13. The list of national and international regulatory documents/standards, which matches medical device	EN 45502-1, EN 45502-2-1, EN ISO 15223-1, EN ISO 14971, EN ISO 11607-1, EN 1041, EN ISO 11135, EN 980, EN 62366, EN ISO 10993-1, ISO 10993-3, ISO 10993-4, ISO 10993-5, ISO 10993-6, ISO 10993-7, ISO 10993-10, ISO 10993-11, ISO 27186.
14. Customer Service	For more information, contact your local Medtronic representative, or contact Medtronic Inc. as listed below: Medtronic Inc., USA Information for patients: Toll-free +1 800 551 5544 (7 a.m.– 6 p.m., Monday–Friday, central time) Fax +1 763 514 1855 Representative in Russia: Medtronic, LLC, 123112, Moscow, in.ter.c. Presnensky Municipal District, Presnenskaya nab., 10, floor 9, office III, room 41 Tel.: (495) 580-73-77, Fax: (495) 580-73-78, www.medtronic.ru

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