SEELER & DR.FELLMETH
NOTARE

NOTARE SEELER & DR. FELLMETH · Q7, 23 · 68161 MANNHEIM

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
Date: 26 January 2023

i.V.

Andrea Weber
Manager Global Regulatory Affairs

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English

Intended use

The S1 Rinse Solution is intended to be used with every cobas b 221 system instrument configuration in a near-patient testing environment for the washing of all tubing and the system. The empty S1 Rinse Solution bottle can be used for the collection of liquid waste instead of the W Waste Container.

Not for self-testing.

Summary

S1 Rinse Solution is used in critical professional near-patient environments such as e.g. intensive care areas, emergency departments and operating rooms as well as the hospital laboratory. For more information on the handling of the reagent refer to the latest cobas b 221 system Instructions for Use.

Test principle

Exchanging an S1 Rinse Solution

⚠ WARNING! Perform QC measurements on 3 levels after each S1 Rinse Solution exchange on the cobas b 221 system. Failure to follow QC protocols or ignoring QC results may lead to incorrect patient results.

⚠ CAUTION! Do not insert damaged S1 Rinse Solution in the cobas b 221 system. The cobas b 221 S1 Rinse Solution must be installed according to instructions given in the cobas b 221 system Instructions for Use.

Reagents

Specimen collection is an aqueous solution of detergents and preservatives.

User training

Do not carry out operation tasks or maintenance actions unless you have received training from Roche Diagnostics. Leave tasks that are not described in the user documentation to trained Roche Service representatives.

Precautions and warnings

For in vitro diagnostic use by health care professionals. Exercise the normal precautions required for handling all laboratory reagents. Dispose of waste according to accepted laboratory instructions and procedures. Apply all relevant local disposal regulations to determine the safe disposal. Safety data sheet available for professional user on request.

The cobas b 221 system components classified as follows in accordance with the Regulation (EC) No. 753/2008.

EUH 208 May produce an allergic reaction.

Handling and disposal

Empty containers should be taken to an approved waste handling site for disposal.

- Wear appropriate personal protective equipment, including, but not limited to, the following items: Eye protection, lab coat and lab gloves.
- When handling S1 Rinse Solution, keep it upright to avoid spillage.
- If you spill contents from the S1 Rinse Solution, wipe it up immediately.
- If contents from S1 Rinse Solution come into contact with your skin or eye, wash the affected area immediately with plenty of water.

The empty S1 Rinse Solution container can be used as a W Waste Container. Remove the label from the empty S1 Rinse Solution bottle, and immediately discard it.

Information regarding the label removal process can be found in the cobas b 221 system Instructions for Use.

Storage and stability

Unopened in the closed box:

- up to the stated expiration date at +2°C to +30°C (+35°F to +86°F).

In-use:

- up to the stated expiration date, for a maximum of 42 days when installed in the cobas b 221 system analyzer.

The operating conditions of the analyzer are further described in the cobas b 221 system Instructions for Use.

Specimen collection and preparation

Specimen collection and preparation is described in the Instructions for Use of the cobas b 221 system.

Materials required (but not provided)

The minimum required components are listed below, but not limited to:

- One of the following instruments:
- REF 0333171001 cobas b 221 <2> system
- REF 0333173601 cobas b 221 <3> system
- REF 0333175401 cobas b 221 <4> system

All remaining supplies contain suitable electrodes and a cassette or suitable dummy electrodes and a cassette. Refer to the cobas b 221 system Instructions for Use for further details. The following supplies must be present:

- REF 032603184 S1 Flush Pack
- REF 0314405401 W Waste Container or Empty S1 Flush Pack

In case of cobas b 221 <3> system, the following solution must be used:

- REF 032603184 S3 Flush Pack A

Any level of the following Quality Control material:

- REF 0333192001 COMBITROL PLUS B, Level 1
- REF 0333192101 COMBITROL PLUS B, Level 2
- REF 0333191901 COMBITROL PLUS B, Level 3
- REF 0333191960 COMBITROL PLUS B, Level 1 US only
- REF 0333191960 COMBITROL PLUS B, Level 2 US only
- REF 0333191960 COMBITROL PLUS B, Level 3 US only

In case of using an AU680001 COOKMSS Verification Material:

- REF 0333192001 AUTOQC
- The following QC materials are to be used with the AUTO QC module:
- REF 0333191901 AUTO-CONTROL PLUS B, Level 1
- REF 0333192101 AUTO-CONTROL PLUS B, Level 2
- REF 0333191901 AUTO-CONTROL PLUS B, Level 3

Refer to the cobas b 221 system Instructions for Use for further details and for details of other accessories.

For further information, please refer to the cobas b 221 system Instructions for Use, and the Method Sheets of all necessary components.

Any serious incident that has occurred in relation to the device should be notified to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

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