October 29, 2018

CERTIFICATION OF TRUTHFUL AND ACCURATE STATEMENT

I, Nivedita S Rao, in my capacity as Sr. Regulatory Affairs Specialist – Abbott Vascular certify that the attached **XIENCE Alpine Everolimus Eluting Coronary Stent System (EECSS), Instructions For Use** is a true, correct and complete copy of the original document in my control and that to the best of my knowledge.

Nivedita S Rao

Sr. Regulatory Affairs Specialist

State of California
County of Santa Clara


A notary public or other officer completing this certificate verifies only the identity of the individual who signed the document, to which this certificate is attached, and not the truthfulness, accuracy, or validity of that document.

On October 29, 2018, before me, Joanna Marie Lacuna Manuel, Notary Public, personally appeared Nivedita S Rao, who proved to me on the basis of satisfactory evidence to be the person whose name is subscribed to the within instrument and acknowledged to me that she executed the same in her authorized capacity, and that by her signature on the instrument the person, or the entity upon behalf of which the person acted, executed the instrument.

I certify under PENALTY OF PERJURY under the laws of the State of California that the foregoing paragraph is true and correct.

WITNESS my hand and official seal.

Joanna Marie Lacuna Manuel, Notary Public

INSTRUCTION FOR USE
XIENCE Alpine Everolimus Eluting Coronary Stent System**Table of Contents**

1.0	MEDICAL DEVICE DESCRIPTION
1.1.	Medical Device Name
1.2.	Details of Manufacturers
1.3.	Technical Specification
1.4.	How Supplied
2.0	INDICATIONS
3.0	CONTRAINDICATIONS
4.0	WARNINGS
5.0	PRECAUTIONS
5.1	Stent Handling
5.2	Stent Placement
5.3	Use in Conjunction with Other Procedures
5.4	Stent / System Removal
5.5	Post Implant
5.6	Use in Special Populations
5.6.1	Pregnancy
5.6.2	Lactation
5.6.3	Pediatric use
5.7	Magnetic Resonance Imaging (MRI) Safety Information
5.8	Drug Interactions
5.9	Immune Suppression Potential
5.10	Lipid Elevation Potential
6.0	DRUG INFORMATION
6.1	Interactions with Drugs or Other Substances
6.2	Pregnancy
6.3	Lactation
7.0	POTENTIAL ADVERSE EVENTS
8.0	PATIENT SELECTION AND TREATMENT
8.1	Individualization of Treatment
9.0	CLINICIAN USE INFORMATION
9.1	Inspection Prior to Use
9.2	Materials Required
9.3	Preparation
9.3.1	Packaging Removal
9.3.2	Guide Wire Lumen Flush
9.3.3	Delivery System Preparation
9.4	Delivery Procedure

- 9.5 Deployment Procedure
- 9.6 Removal Procedure
- 9.7 Post-deployment Dilation of Stent Segments
- 10.0 SPIRIT AND XIENCE FAMILY OF CLINICAL TRIALS
 - 10.1 Pre-Market Clinical Trials
 - Table 10.1-1: SPIRIT Family of XIENCE V Clinical Trial Designs (Pre-market)
 - Table 10.1-2: SPIRIT Family of Clinical Trials Angiographic Results (Pre-market)
 - Table 10.1-3: SPIRIT Family of Clinical Trials Principal 1-Year Clinical Outcomes (Pre-market)
 - Table 10.1-4: SPIRIT Family of Clinical Trials Principal Clinical Outcomes from Latest Follow-up (Pre-market)
 - 10.2 Post-Market Clinical Trials
 - Table 10.2-1: SPIRIT and XIENCE V Family Clinical Trial Designs (Post-market)
 - Table 10.2-2: SPIRIT and XIENCE V Family of Trials Principal Clinical Outcomes (Post-market)
 - 10.3 Pooled Clinical Trials Analysis for Stent Thrombosis Subsequent to DAPT Interruption / Discontinuation
 - Table 10.3-1: Stent Thrombosis (ARC Definite / Probable) Subsequent to DAPT Interruption / Discontinuation through 730 days Versus No Interruption; Pooled Data from 7 Pre- and Post-market Trials
- 11.0 STERILIZATION CONDITIONS
- 12.0 PACKAGING, USAGE
- 13.0 TRANSPORTATION AND STORAGE CONDITIONS
- 14.0 ENVIRONMENT PROTECTION REQUIREMENTS
- 15.0 DISPOSAL
- 16.0 DISCLAIMER WARRANTY AND LIMITATION OF REMEDY
- 17.0 PATENTS AND TRADEMARKS
- 18.0 GRAPHICAL SYMBOLS FOR MEDICAL DEVICE LABELING

1.0 MEDICAL DEVICE DESCRIPTION

1.1 Medical Device Name

The XIENCE Alpine Everolimus Eluting Coronary Stent System (hereinafter "XIENCE Alpine Stent System", "XIENCE Alpine Stent", "EECSS").

1.2 Details of Manufacturers

Legal Manufacturer: ABBOTT VASCULAR

3200 Lakeside Drive, Santa Clara, California, 95054, USA

Tel.: 408 845 3000, Fax: 408 845 3333

Manufacturing Sites:

- 26531 Ynez Road, Temecula, CA 92591, USA;
- Cashel Road, Clonmel, County Tipperary, Ireland;

Authorized Representative in Russia

Technoproject Ltd.

191 014 St. Petersburg, Mayakovskogo str., 22-24 lit. A, room 16N, Russian Federation
tel.: +7 (812) 272-97-87, 275-82-33

1.3 Technical Specification


Figure 1. XIENCE Alpine Everolimus Eluting Coronary Stent System



Figure 2. Flushing tool

The XIENCE Alpine Stent System includes the following sizes: 2.0, 2.25, 2.5, 2.75, 3.0, 3.25, 3.5, and 4.0 mm diameters; 8, 12, 15, 18, 23, 28, 33, and 38 mm stent lengths. The 33 and 38 mm lengths are omitted for the 2.0 and 2.25 mm diameters.

XIENCE Alpine Everolimus Eluting Coronary Stents are balloon-expandable stents composed of Co-Cr-W-Ni (ASTM F90) with a PBMA and PVDF-HFP/Everolimus coating. The Small and Medium stent designs have the same 3-link, short/long crest design. In order to allow for different expansion ranges, each design varies the number of short and long crests per ring. The Small Stents have 3 short and 3 long crests per ring. The Medium Stents have 3 short and 6 long crests per ring.

The XIENCE Alpine Delivery Catheter has a rapid-exchange (RX) co-axial design with the balloon and stent at the distal end of the catheter. The proximal lumen provides for inflation of the balloon with contrast medium and the central distal lumen permits a guide wire to facilitate advancement of the catheter. The distal and intermediate portions of the delivery catheter, the tip, and tapers of the balloon are coated with HYDROCOAT™ Hydrophilic Coating.

Radiopaque markers are positioned on the inner member underneath the balloon to provide accurate positioning of the stent/balloon in the artery. The balloon is designed to deliver an expandable stent of known diameter and length at specified pressures.

Markers located on the proximal outer shafts (95 cm and 105 cm proximal to the distal tip) help the physician gauge the delivery catheter position relative to the guiding catheter tip. A shaft color change denotes the guide wire exit notch. Working catheter length is 145 cm.

An adaption arm on the proximal end of the delivery catheter provides access to the inflation lumen. It is designed with a Luer-lock fitting to facilitate connection to an inflation device.

The available dose of Everolimus on the stent varies by size as follows:

Table 1.

Stent Diameter (mm)	Stent Length (mm)	Drug Dose (µg)
2.0, 2.25, 2.5, 2.75, 3.0, 3.25	8	40
2.0, 2.25, 2.5, 2.75, 3.0, 3.25	12	60
2.0, 2.25, 2.5, 2.75, 3.0, 3.25	15	74
2.0, 2.25, 2.5, 2.75, 3.0, 3.25	18	88
2.0, 2.25, 2.5, 2.75, 3.0, 3.25	23	109
2.0, 2.25, 2.5, 2.75, 3.0, 3.25	28	137
2.5, 2.75, 3.0, 3.25	33	157
2.5, 2.75, 3.0, 3.25	38	185
3.5, 4.0	8	50
3.5, 4.0	12	75
3.5, 4.0	15	91
3.5, 4.0	18	116
3.5, 4.0	23	141
3.5, 4.0	28	174
3.5, 4.0	33	199
3.5, 4.0	38	232

In vitro Device Specifications are presented in a table below:

Table 2.

Stent Diameter (mm)	Stent Length (mm)	Minimum Guiding Catheter Compatibility (ID)	* <i>in vitro</i> Stent Nominal Pressure		Rated Burst Pressure (RBP)		Stent Free Area (%)
			(atm)	kPa	(atm)	kPa	
2.0	8, 12, 15, 18, 23, 28	5F (0.056" / 1.42 mm)	10	1013	18	1824	81
2.25	8, 12, 15, 18, 23, 28	5F (0.056" / 1.42 mm)	10	1013	18	1824	83
2.5	8, 12, 15, 18, 23, 28, 33, 38	5F (0.056" / 1.42 mm)	10	1013	18	1824	84
2.75	8, 12, 15, 18, 23, 28, 33, 38	5F (0.056" / 1.42 mm)	10	1013	18	1824	86
3.0	8, 12, 15, 18, 23, 28, 33, 38	5F (0.056" / 1.42 mm)	10	1013	18	1824	87
3.25	8, 12, 15, 18, 23, 28, 33, 38	5F (0.056" / 1.42 mm)	10	1013	18	1824	88
3.5	8, 12, 15, 18, 23, 28, 33, 38	5F (0.056" / 1.42 mm)	10	1013	18	1824	85
4.0	8, 12, 15, 18, 23, 28, 33	5F (0.056" / 1.42 mm)	10	1013	18	1824	87
4.0	38	6F (0.066" / 1.68 mm)	10	1013	18	1824	87

* Assure full deployment of the stent. Deployment pressures should be based on lesion characteristics.

The technical characteristics of the flushing tool are presented in the table 3.

Table 3.

Needle outer diameter	0.0122" ± 0.0004"
Blunt Point	Square cut

1.4 How supplied

Sterile – This device is sterilized with ethylene oxide gas. Non-pyrogenic. Do not use if the package is open or damaged.

This single-use device cannot be reused on another patient, as it is not designed to perform as intended after the first usage. Changes in mechanical, physical, and / or chemical characteristics introduced under conditions of repeated use, cleaning, and / or resterilization may compromise the integrity of the design and / or materials, leading to contamination due to narrow gaps and / or spaces and diminished safety and / or performance of the device. Absence of original labeling may lead to misuse and eliminate traceability. Absence of original packaging may lead to device damage, loss of sterility, and risk of injury to the patient and / or user.

Contents – One (1) XIENCE Alpine Everolimus Eluting Coronary Stent System; one (1) flushing tool.

2.0 INDICATIONS

The XIENCE Alpine Everolimus Eluting Coronary Stent System is indicated for improving coronary luminal diameter in the following:

- Patients with symptomatic ischemic heart disease due to discrete *de novo* native coronary artery lesions.
- For restoring coronary flow in patients experiencing acute myocardial infarction who present within 12 hours of symptom onset.
- For the treatment of patients with concomitant diabetes, acute coronary syndrome, dual vessel lesions (two lesions in two different epicardial vessels), lesions residing within small coronary vessels; lesions where treatment results in the jailing of side branches (lesions with a side branch < 2 mm in diameter or an ostial stenosis < 50%); for the treatment of elderly patients (age ≥ 65), and for treatment of both men and women.
- For the treatment of patients presenting with in-stent restenosis in coronary artery lesions; chronic total occluded coronary artery lesions (defined as coronary artery lesions with TIMI flow 0 and lasting longer than 3 months); and coronary artery bifurcation lesions.

In all cases, the treated lesion length should be less than the nominal stent length (8 mm, 12 mm, 15 mm, 18 mm, 23 mm, 28 mm, 33 mm, or 38 mm) with a reference vessel diameter of ≥ 2.00 mm and ≤ 4.25 mm.

3.0 CONTRAINDICATIONS

The XIENCE Alpine Everolimus Eluting Coronary Stent System is contraindicated for use in:

- Patients who cannot tolerate, including allergy or hypersensitivity to, procedural anticoagulation or the post procedural antiplatelet regimen.
- Patients with hypersensitivity or contraindication to everolimus or structurally related compounds, or known hypersensitivity to stent components (cobalt, chromium, nickel, tungsten, acrylic, fluoropolymers), or with contrast sensitivity.

4.0 WARNINGS

- For single use only. Do not resterilize or reuse. Note the product "Use by" date on the package.
- It is not recommended to treat patients having a lesion that prevents complete inflation of an angioplasty balloon.
- Antiplatelet therapy should be administered post-procedure (See *Section 8.1, Individualization of Treatment*).
- This product should not be used in patients who are not likely to comply with the recommended antiplatelet therapy
- Judicious selection of patients is necessary, since the use of this device carries the associated risk of stent thrombosis, vascular complications, and / or bleeding events.

5.0 PRECAUTIONS

5.1 Stent Handling

- Implantation of the stent should be performed only by physicians who have received appropriate training.
- Stent placement should only be performed at centers where emergency coronary artery bypass graft (CABG) is available.
- **The foil pouch is not a sterile barrier.** The inner header bag (pouch) within the foil pouch is the sterile barrier. **Only the contents of the inner pouch should be considered sterile.** The outside surface of the inner pouch is NOT sterile.

- To confirm sterility has been maintained, ensure that the inner package sterile barrier has not been opened or damaged prior to use.
- Care should be taken to control the guiding catheter tip during stent delivery, deployment, and balloon withdrawal. Before withdrawing the stent delivery system, visually confirm complete balloon deflation by fluoroscopy to avoid guiding catheter movement into the vessel and subsequent arterial damage.
- Special care must be taken not to handle or in any way disrupt the stent on the balloon. This is most important during catheter removal from packaging, placement over the guide wire, and advancement through the rotating hemostatic valve adapter and guiding catheter hub.
- **Do not manipulate, touch, or handle the stent** with your fingers, as this may cause coating damage, contamination, or dislodgement of the stent from the delivery balloon.
- Use only the appropriate balloon inflation media. Do not use air or any gaseous medium to inflate the balloon, as this may cause uneven expansion and difficulty in deployment of the stent. If gaseous medium is used and balloon rupture occurs, there is the potential of causing air embolism and / or vessel injury.

5.2 Stent Placement

- Use guiding catheters which have lumen sizes that are suitable to accommodate the stent delivery system.
- Do not prepare or pre-inflate the delivery system prior to stent deployment other than as directed. Use balloon purging technique described under *Section 9.3.3, Delivery System Preparation*.
- When pre-dilatation is performed, an appropriate balloon size should be used. Failure to do so may increase the difficulty of stent placement and cause procedural complications.
- The decision to pre-dilate the lesion with an appropriate sized balloon should be based on patient and lesion characteristics. Direct stenting in less complex coronary lesions has been shown to be as effective and safe as stenting with pre-dilation for device lengths up to 28 mm in real-world settings. **If pre-dilation is performed**, limit the length of pre-dilation by the PTCA balloon to avoid creating a region of vessel injury that is outside the boundaries of the implanted stent.
- When introducing the delivery system into the vessel, do not induce negative pressure on the delivery system. This may cause dislodgement of the stent from the balloon.
- Do not torque the catheter more than one (1) full turn.
- Implanting a stent may lead to dissection of the vessel distal and / or proximal to the stent, and may cause abrupt closure of the vessel requiring additional

intervention (CABG, surgery, further dilation, placement of additional stents, or other).

- An unexpanded stent may be retracted into the guiding catheter one time only. An unexpanded stent should not be reintroduced into the artery once it has been pulled back into the guiding catheter. Subsequent movement in and out through the distal end of the guiding catheter should not be performed, as the stent may be damaged or dislodged during retraction back into the guiding catheter.
- Should resistance be felt at any time during removal of the undeployed coronary stent system, please refer to the steps provided in *Section 6.4 Stent / System Removal*.
- Do not expand the stent if it is not properly positioned in the vessel. (See *Section 6.4 Stent / System Removal – Precautions*.)
- The inflated balloon diameter of the system used to deploy the stent should approximate the diameter of the vessel. Oversizing of the stent can result in a ruptured vessel. To ensure full expansion of the stent, the balloon should be inflated to a minimum of nominal pressure.
- Do not exceed RBP as indicated on product label. Monitor balloon pressures during inflation. Use of pressures higher than specified on product label may result in a ruptured balloon with possible intimal damage and dissection.
- When performed, post-dilatation should be performed at high pressure with a noncompliant balloon.
- Under-expansion of the stent may result in stent movement. Care must be taken to properly size the stent to ensure that the stent is in full contact with the arterial wall upon deflation of the balloon. All efforts should be made to assure that the stent is not under dilated. Refer to *Section 9.0 Clinician Use Information*.
- Placement of a stent has the potential to compromise side branch patency.
- Stent retrieval methods (use of additional wires, snares, and / or forceps) may result in additional trauma to the coronary vasculature and / or the vascular access site. Complications may include bleeding, hematoma, or pseudoaneurysm.
- When treating multiple lesions within the same vessel, stent the distal lesion prior to stenting the proximal lesion. Stenting in this order obviates the need to cross the proximal stent during placement of the distal stent, and reduces the chance of damaging or dislodging the proximal stent.
- When multiple DES stents are required, only stent material with similar composition (e.g., XIENCE Everolimus Eluting Coronary family Stents with the identical cobalt-chromium stent substrate and identical drug-eluting polymer coating) should be used. Potential interaction with other drug-eluting stents or

coated stents has not been evaluated, and should be avoided. Placing multiple stents of different metals in contact with each other may increase the potential for corrosion *in vivo*, although *in vitro* corrosion tests using an L-605 CoCr alloy stent in combination with a 316L stainless steel alloy stent did not appear to increase corrosion.

- The extent of the patient's exposure to drug and polymer is directly related to the number of stents implanted. A patient can receive up to four XIENCE Xpedition, Alpine Everolimus Eluting Coronary Stents or other everolimus eluting coronary stents from the XIENCE family (i.e., XIENCE V, XIENCE PRIME, XIENCE Xpedition, XIENCE Alpine) depending on the number of vessels treated and the lesion length. Those patients receiving bailout stenting will receive additional XIENCE family stents. The use of multiple XIENCE family stents will result in the patient receiving larger amounts of drug and polymer.
- The safety and effectiveness of the XIENCE Alpine Everolimus Eluting Coronary Stent in patients with prior brachytherapy of the target lesion or the use of brachytherapy for treated site restenosis in a XIENCE Everolimus Eluting Coronary family stent have not been established. Both vascular brachytherapy and the XIENCE Everolimus Eluting Coronary family stents alter arterial remodeling. The potential combined effect on arterial remodeling by these two treatments is not known.

5.3 Use in Conjunction with Other Procedures

- While vessel preparation in complex lesions may include the use of various mechanical atherectomy devices, the safety and effectiveness of the XIENCE Alpine stents have not been established in clinical trials with use of either mechanical atherectomy devices (directional atherectomy catheters, rotational atherectomy catheters) or laser angioplasty catheters.

5.4 Stent / System Removal

- **Stent delivery system removal prior to stent deployment:**

If removal of a stent system is required prior to deployment, ensure that the guide catheter is coaxially positioned relative to the stent delivery system, and cautiously withdraw the stent delivery system into the guiding catheter. Should unusual resistance be felt at any time when withdrawing the stent into the guide catheter, the stent delivery system and the guide catheter should be removed as a single unit. This should be done under direct visualization with fluoroscopy.

- **Withdrawal of the stent delivery system / post-dilation balloon from the deployed stent:**

1. Deflate the balloon by pulling negative on the inflation device. Larger and longer balloons will take more time (up to 30 seconds) to deflate than smaller and shorter balloons. Confirm balloon deflation under fluoroscopy and wait 10 – 15 seconds longer.
2. Position the inflation device to "negative" or "neutral" pressure.

3. Stabilize guide catheter position just outside coronary ostium and anchor in place. Maintain guide wire placement across stent segment.
4. Gently remove the stent delivery system / post-dilation balloon with slow and steady pressure.
5. Tighten the rotating hemostatic valve.

Notes:

1. **If during withdrawal of the catheter from the deployed stent, resistance is encountered, use the following steps to improve balloon rewrap:**
 - Re-inflate the balloon up to nominal pressure.
 - Repeat steps 1 through 5 above.
2. After successful withdrawal of the balloon from the deployed stent, should any resistance be felt at any time when withdrawing the stent delivery system or post-dilatation balloon into the guide catheter, remove the entire system as a single unit.
 - Failure to follow these steps and / or applying excessive force to the delivery system can potentially result in loss or damage to the stent and / or delivery system components.
 - If it is necessary to retain guide wire position for subsequent artery / lesion access, leave the guide wire in place and remove all other system components.

5.5 Post Implant

- If necessary to cross a newly deployed stent with a guide wire, balloon, delivery system, or imaging catheters, exercise care to avoid disrupting the stent geometry.
- Subsequent restenosis may require repeat dilatation of the arterial segment containing the stent. The long-term outcome following repeat dilatation of stents is unknown at present.
- If the patient requires imaging, see *Section 5.7 Magnetic Resonance Imaging (MRI) Safety Information*.

5.6 Use in Special Populations**5.6.1 Pregnancy**

Pregnancy Category C: see *Section 6.2 Pregnancy*. This product has not been tested in pregnant women or in men intending to father children. Effects on the developing foetus have not been studied. Effects of a XIENCE V on pre-natal and post-natal rat development were no different than the controls. When administered at oral doses of 0.1 mg/kg or above to animals, everolimus has shown reproductive toxicity effects including embryo toxicity and foetotoxicity¹. Effective contraception is recommended to be initiated

before implanting and continued for one year after implantation. While there is no contraindication, the risks and reproductive effects are unknown at this time¹.

5.6.2 Lactation

See *Section 6.3 Lactation*. It is unknown whether everolimus is distributed in human milk. A decision should be made whether or not to discontinue nursing prior to stent implantation, considering the importance of the stent to the mother.

5.6.3 Pediatric use

The safety and effectiveness of the XIENCE Alpine stent in pediatric subjects have not been established.

5.7. Magnetic Resonance Imaging (MRI) Safety Information

Non-clinical testing has demonstrated that the XIENCE Alpine stent, in single and in overlapped configurations up to 71 mm in length, is MR Conditional. A patient with this device can be safely scanned in a MR system meeting the following conditions:

- Static magnetic field of 1.5 or 3 T
- Maximum Spatial gradient field of 2500 gauss/cm
- Maximum MR System reported whole-body-averaged specific absorption rate (SAR) of 2.0 W/kg (normal operating mode)

Under the scan conditions defined above, the XIENCE Alpine stent is expected to produce a maximum temperature rise of less than 3.3°C after 15 minutes of continuous scanning.

In non-clinical testing, the image artifact caused by the device extends approximately 13 mm from the XIENCE Alpine stent when imaged with a gradient pulse sequence and a 3T MRI system.

5.8 Drug Interactions

See *Section 6.1 Interactions with Drugs or Other Substances*. Several drugs are known to affect everolimus metabolism, and other drug interactions may also occur. Everolimus is known to be a substrate for both cytochrome P4503A4 (CYP3A4) and P-glycoprotein (PgP). Everolimus absorption and subsequent elimination may be influenced by drugs that affect these pathways. Everolimus has also been shown to reduce the clearance of some prescription medications when administered orally along with cyclosporine (CsA). Formal drug interaction studies have not been performed with the XIENCE Alpine stent because of limited exposure to everolimus eluted from the stent. Therefore, due consideration should be given to the potential for both systemic and local drug interactions in the vessel wall, when deciding to place the XIENCE Alpine stent in a

¹ Certican® UK label Mar 2015, Afinitor® EU authorization SPC Dec 2014, Votubia® EU SPC Sept 2014, Afinitor® US label Jan 2015, and Zortress® US label Sept 2015. Refer to www.MHRA.gov.uk, www.ema.europa.eu, and www.fda.gov for the most recent versions of these SPC/labels.

patient taking a drug with known interaction with everolimus, or when deciding to initiate therapy with such a drug in a patient who has recently received the XIENCE Alpine stent.

5.9 Immune Suppression Potential

Everolimus, the XIENCE Alpine stent active ingredient, is an immunosuppressive agent. Immune suppression was not observed in the SPIRIT and XIENCE family of clinical trials. However, for patients who receive several XIENCE Alpine devices simultaneously, it may be possible for everolimus systemic concentrations to approach immunosuppressive levels temporarily, especially in patients who also have hepatic insufficiency or who are taking drugs that inhibit CYP3A4 or P-glycoprotein. Therefore, consideration should be given to patients taking other immunosuppressive agents or who are at risk for immune suppression.

5.10 Lipid Elevation Potential

Oral everolimus use in renal transplant and advanced renal cell carcinoma patients was associated with increased serum cholesterol and triglyceride levels, which in some cases required treatment. The effect was seen with both low and high dose prolonged oral therapy in a dose related manner. When used according to the indications for use, exposure to systemic everolimus concentrations from the XIENCE Alpine stent is expected to be significantly lower than concentration exposure usually obtained in transplant patients. Increased serum cholesterol and triglyceride levels were not observed in the SPIRIT and XIENCE family of clinical trials. Oral administration of everolimus in combination with cyclosporine has been associated with increased serum cholesterol and triglyceride levels.

6.0 DRUG INFORMATION

6.1 Interactions with Drugs or Other Substances

Everolimus is extensively metabolized by the cytochrome P4503A4 (CYP3A4) in the gut wall and liver, and is a substrate for the countertransporter P-glycoprotein. Everolimus has also been shown to reduce the clearance of some prescription medications when it was administered orally along with cyclosporine (CsA). Hence, everolimus, when prescribed as an oral medication, may interact with other medications that include (but are not restricted to) inhibitors and inducers of CYP3A4 isozymes; absorption and subsequent elimination of everolimus may be influenced by drugs that affect these pathways. Formal drug interaction studies have not been performed with the XIENCE Alpine or XIENCE V stents because of limited systemic exposure to everolimus eluted from XIENCE V. However, consideration should be given to the potential for both systemic and local drug interactions in the vessel wall when deciding to place the XIENCE Alpine stent in a subject taking a drug with known interaction with everolimus.

Everolimus, when prescribed as an oral medication, may interact with the following including, but not limited to:

- CYP3A4 / P-glycoprotein isozyme inhibitors
 - Antifungal agents (e.g., fluconazole, ketoconazole, itraconazole, posaconazole, voriconazole)
 - Macrolide antibiotics (e.g., erythromycin, clarithromycin, telithromycin)
 - Calcium channel blockers (e.g., verapamil, nifedipine, diltiazem)
 - Protease inhibitors (e.g., ritonavir, atazanavir, saquinavir, darunavir, indinavir, nelfinavir, amprenavir, fosamprenavir)
 - Other (e.g., cyclosporine, nefazodone, cisapride, metoclopramide, bromocriptine, cimetidine, danazol, sildenafil, terfenadine, astemizole, grapefruit / grapefruit juice, digoxin)
- CYP3A4 / P-glycoprotein isozyme inducers
 - Antibiotics (e.g., rifampin, rifabutin, ciprofloxacin, ofloxacin)
 - Anticonvulsants (e.g., carbamazepine, phenobarbital, phenytoin)
 - Non-nucleoside reverse transcriptase inhibitors (e.g., efavirenz, nevirapine)
 - Glucocorticoids (e.g., dexamethasone, prednisone, prednisolone)
 - HMGCoA reductase inhibitors (simvastatin, lovastatin)
 - Other (e.g., St. John's Wort)

For more detailed drug interaction information, please reference the most recent everolimus drug label¹.

6.2 Pregnancy

Pregnancy Category C: There are no adequate everolimus or XIENCE Alpine stent-related studies in pregnant women. Effects of a similar stent, XIENCE V, on pre-natal and post-natal rat development were no different than the controls. When administered at oral doses of 0.1 mg/kg or above to animals, everolimus has shown reproductive toxicity effects including embryotoxicity and fetotoxicity². Effective contraception is recommended to be initiated before implanting a XIENCE Alpine stent and continued for one year post-implantation. The XIENCE Alpine stent should be used in pregnant women only if potential benefits of the stent outweigh potential risks.

The safety of the XIENCE Alpine stent has not been evaluated in males intending to father children.

6.3 Lactation

It is unknown whether everolimus is distributed in human milk. Also, everolimus pharmacokinetic and safety profiles have not been determined in infants. Consequently, mothers should be advised of potential serious adverse reactions to everolimus in nursing infants. Prior to XIENCE Alpine stent implantation, decisions should be made regarding whether to discontinue nursing or conduct an alternate percutaneous coronary intervention procedure.

² Certican® UK label Mar 2015, Afinitor® EU authorization SPC Dec 2014, Votubia® EU SPC Sept 2014, Afinitor® US label Jan 2015, and Zortress® US label Sept 2015. Refer to www.MHRA.gov.uk, www.ema.europa.eu, and www.fda.gov for the most recent versions of these SPC/labels.

7.0 POTENTIAL ADVERSE EVENTS

Adverse events that may be associated with PCI treatment procedures and the use of a stent in native coronary include, but are not limited to, the following:

- Allergic reaction or hypersensitivity to latex, contrast agent, anesthesia, device materials (cobalt, chromium, nickel, tungsten, acrylic, and fluoropolymers), and drug reactions to everolimus, anticoagulation, or antiplatelet drugs
- Vascular access complications which may require transfusion or vessel repair, including:
 - Catheter site reactions
 - Bleeding (ecchymosis, oozing, hematoma, hemorrhage, retroperitoneal hemorrhage)
 - Arteriovenous fistula, pseudoaneurysm, aneurysm, dissection, perforation / rupture
 - Embolism (air, tissue, plaque, thrombotic material or device)
 - Peripheral nerve injury
 - Peripheral ischemia
- Coronary artery complications which may require additional intervention, including:
 - Total occlusion or abrupt closure
 - Arteriovenous fistula, pseudoaneurysm, aneurysm, dissection, perforation / rupture
 - Tissue prolapse / plaque shift
 - Embolism (air, tissue, plaque, thrombotic material, or device)
 - Coronary or stent thrombosis (acute, subacute, late, very late)
 - Stenosis or restenosis
- Pericardial complications which may require additional intervention, including:
 - Cardiac tamponade
 - Pericardial effusion
 - Pericarditis
- Cardiac arrhythmias (including conduction disorders, aspecific, atrial and ventricular arrhythmias)
- Cardiac ischemic conditions (including myocardial ischemia, myocardial infarction (including acute), coronary artery spasm, and unstable or stable angina pectoris)
- Stroke / cerebrovascular accident (CVA) and Transient Ischemic Attack (TIA)
- System organ failures:
 - Cardio-respiratory arrest
 - Cardiac failure
 - Cardiopulmonary failure (including pulmonary edema)
 - Renal Insufficiency / failure
 - Shock
- Blood cell disorders (including Heparin Induced Thrombocytopenia (HIT))
- Hypotension / hypertension
- Infection
- Nausea and vomiting
- Palpitations, dizziness, and syncope
- Chest pain
- Fever

- Pain
- Death

Adverse events associated with daily oral administration of everolimus in doses varying from 1.5 mg to 10 mg daily can be found in the Summary of Product Characteristics (SPC) and labels for the drug³. The risks described below include the anticipated adverse events relevant for the cardiac population referenced in the contraindications, warnings and precaution sections of the everolimus labels / SPCs and / or observed at incidences $\geq 10\%$ in clinical trials with oral everolimus for different indications. Please refer to the drug SPCs and labels¹ for more detailed information and less frequent adverse events:

- Abdominal pain
- Anemia
- Angioedema (increased risk with concomitant ACE inhibitor use)
- Arterial thrombotic events
- Bleeding and coagulopathy (including Hemolytic Uremic Syndrome [HUS], Thrombotic Thrombocytopenic Purpura [TTP], and thrombotic microangiopathy – increased risk with concomitant cyclosporine use)
- Constipation
- Cough
- Diabetes mellitus
- Diarrhea
- Dyspnea
- Embryo-fetal toxicity
- Erythema
- Erythroderma
- Headache
- Hepatic Artery Thrombosis (HAT)
- Hepatic disorders (including hepatitis and jaundice)
- Hypersensitivity to everolimus active substance, or to other rapamycin derivatives
- Hypertension
- Infection (bacterial, fungal, viral or protozoan infections, including infections with opportunistic pathogens). Polyoma virus-associated nephropathy (PVAN) JC virus associated progressive multiple leukoencephalopathy (PML), fatal infections and sepsis have been reported in patients treated with oral everolimus.
- Kidney arterial and venous thrombosis
- Laboratory test alterations (elevations of serum creatinine, proteinuria, hypokalemia; hyperglycemia, dyslipidemia including hypercholesterolemia and hypertriglyceridemia; abnormal liver function tests; decreases in hemoglobin, lymphocytes, neutrophils, and platelets)
- Lymphoma and skin cancer
- Nausea
- Nephrotoxicity (in combination with cyclosporine)
- Non-infectious pneumonitis (including interstitial lung disease)

³ Certican® UK label Mar 2015, Afinitor® EU authorization SPC Dec 2014, Votubia® EU SPC Sept 2014, Afinitor® US label Jan 2015, and Zortress® US label Sept 2015. Refer to www.MHRA.gov.uk, www.ema.europa.eu, and www.fda.gov for the most recent versions of these SPC/labels.

- Oral ulcerations
- Pain
- Pancreatitis
- Pericardial effusion
- Peripheral edema
- Pleural effusion
- Pneumonia
- Pyrexia
- Rash
- Renal failure
- Upper respiratory tract infection
- Urinary tract infection
- Venous thromboembolism
- Vomiting
- Wound healing complications (including wound infections and lymphocele)

8.0 PATIENT SELECTION AND TREATMENT

8.1 Individualization of Treatment

The risks and benefits described above should be considered for each patient before using a device from the family of XIENCE EECSS. Patient selection factors to be assessed should include judgment regarding risk of antiplatelet therapy. Special consideration should be given to those patients with recently active gastritis or peptic ulcer disease.

Antiplatelet drugs should be used in combination with the XIENCE EECSS per the guidelines from the American College of Cardiology, American Heart Association, and Society for Cardiovascular Angiography and Interventions (ACC/AHA/SCAI) and ESC guidelines. Physicians should use the information from the large body of clinical evidence for XIENCE stents, coupled with the current literature on drug-eluting stents, current guidelines and the specific needs of individual patients, to determine the specific antiplatelet / anticoagulation regimen to be used for their patients in general practice.

Current guidelines for the DAPT discontinuation should be followed and are recommended. The decision to interrupt or discontinue DAPT is the responsibility of the treating physician, taking into consideration the individual patient's condition. In case an unanticipated interruption or discontinuation of DAPT is required any time after one month following XIENCE coronary stent implantation, two-year data from the XIENCE coronary clinical trials show low stent thrombosis rates and no observed increased risk for stent thrombosis. Refer to Table 10.3-1 for rate information.

It is very important that the patient comply with the post-procedural antiplatelet recommendations. Premature discontinuation of prescribed antiplatelet medication could result in a higher risk of thrombosis, myocardial infarction, or death. Prior to percutaneous coronary intervention (PCI), if a surgical or dental procedure is anticipated that requires early discontinuation of antiplatelet therapy, the interventionalist and patient should carefully consider whether a drug-eluting stent and its associated recommended

antiplatelet therapy is the appropriate PCI choice. Following PCI, should a surgical or dental procedure be recommended, the risks and benefits of the procedure should be weighed against the possible risk associated with premature discontinuation of antiplatelet therapy.

Patients who require premature discontinuation of antiplatelet therapy secondary to significant active bleeding, should be monitored carefully for cardiac events and, once stabilized, have their antiplatelet therapy restarted as soon as possible per the discretion of their treating physicians.

9.0 CLINICIAN USE INFORMATION

9.1 Inspection Prior to Use

1. Carefully inspect the sterile package before opening and check for damage to the sterile barrier. Do not use if the integrity of the sterile package has been compromised.
2. Do not use after the "Use by" date.
3. Tear open the foil pouch and remove the inner pouch.

Note: The outside of the inner pouch is NOT sterile. Open the inner pouch and pass or drop the product into the sterile field using an aseptic technique.

4. Prior to using the XIENCE Alpine EECSS, carefully remove the system from the package and inspect for bends, kinks, and other damage. Verify that the stent does not extend beyond the radiopaque balloon markers. Do not use if any defects are noted. However, do not manipulate, touch, or handle the stent, which may cause coating damage, contamination, or stent dislodgement from the delivery balloon.

Note: At any time during use of the XIENCE Alpine EECSS, if the stainless steel proximal shaft has been bent or kinked, do not continue to use the catheter.

9.2 Materials Required

- Appropriate arterial sheath
- Appropriate guiding catheter(s)
- 2 – 3 syringes (10 – 20 ml)
- 1,000 u/500 ml heparinized normal saline (HepNS)
- Rotating hemostatic valve with appropriate minimum inner diameter (0.096 inch [2.44 mm])
- 0.014 inch (0.36 mm) x 175 cm (minimum length) guide wire
- Torque device
- Guide wire introducer
- Contrast diluted 1:1 with heparinized normal saline
- Inflation device
- Appropriate sized pre-dilatation angioplasty balloon

- Appropriate sized post-dilatation noncompliant angioplasty balloon
- Three-way stopcock
- Appropriate anticoagulation and antiplatelet drugs

9.3 Preparation

9.3.1 Packaging Removal

Note: The foil pouch is not a sterile barrier. The inner header bag (pouch) within the foil pouch is the sterile barrier. Only the contents of the inner pouch should be considered sterile. The outside surface of the inner pouch is NOT sterile.

5. Carefully remove the delivery system from its protective tubing for preparation of the delivery system. When using a Rapid Exchange (RX) system, do not bend or kink the hypotube during removal.
6. Remove the product mandrel and protective stent sheath by grasping the catheter just proximal to the stent (at the proximal balloon bond site), and with the other hand, grasp the stent protector and gently remove distally. If unusual resistance is felt during product mandrel and stent sheath removal, do not use this product and replace with another. Follow product returns procedure for the unused device.

9.3.2 Guide Wire Lumen Flush

1. Flush the guide wire lumen with HepNS using the flushing tool supplied with the product. Insert the flushing tool into the tip of the catheter and flush until fluid exits the guide wire exit notch.

Note: Avoid manipulation of the stent while flushing the guide wire lumen, as this may disrupt the placement of the stent on the balloon.

9.3.3 Delivery System Preparation

1. Prepare an inflation device / syringe with diluted contrast medium.
2. Attach an inflation device / syringe to the stopcock; attach it to the inflation port of the product. Do not bend the product hypotube when connecting to the inflation device / syringe.
3. With the tip down, orient the delivery system vertically.
4. Open the stopcock to delivery system; pull negative for 30 seconds; release to neutral for contrast fill.
5. Close the stopcock to the delivery system; purge the inflation device / syringe of all air.
6. Repeat steps 3 through 5 until all air is expelled. If bubbles persist, do not use the product.
7. If a syringe was used, attach a prepared inflation device to stopcock.
8. Open the stopcock to the delivery system.
9. Leave on neutral.

Note: While introducing the delivery system into the vessel, do not induce negative pressure on the delivery system. This may cause dislodgement of the stent from the balloon.

Note: If air is seen in the shaft, repeat *Section 9.3.3 Delivery System Preparation*, steps 3 through 5, to prevent uneven stent expansion.

9.4 Delivery Procedure

1. Prepare the vascular access site according to standard practice.
2. The decision to pre-dilate the lesion with an appropriate sized balloon should be based on patient and lesion characteristics. **If pre-dilatation is performed**, limit the length of pre-dilation by the PTCA balloon to avoid creating a region of vessel injury that is outside the boundaries of the XIENCE Alpine stent.
3. For long lesions, size the stent to the diameter of the most distal portion of the vessel.

Note: If choosing between two stent diameters for tight lesions, choose the smaller diameter stent and inflate. See product label for compliance information.

4. Maintain neutral pressure on the inflation device attached to the delivery system. Open the rotating hemostatic valve as wide as possible.
5. Backload the delivery system onto the proximal portion of the guide wire while maintaining guide wire position across the target lesion.
6. Carefully advance the delivery system into the guiding catheter and over the guide wire to the target lesion. When using a Rapid Exchange (RX) system, be sure to keep the hypotube straight. Ensure guiding catheter stability before advancing the stent system into the coronary artery.

Note: If unusual resistance is felt before the stent exits the guiding catheter, do not force passage. Resistance may indicate a problem and the use of excessive force may result in stent damage or dislodgement. Maintain guide wire placement across the lesion and remove the delivery system and guiding catheter as a single unit.

7. Advance the delivery system over the guide wire to the target lesion under direct fluoroscopic visualization. Utilize the radiopaque balloon markers to position the stent across the lesion. Perform angiography to confirm stent position. If the position of the stent is not optimal, it should be carefully repositioned or removed (see *Section 9.6 Removal Procedure*). The balloon markers indicate both the stent edges and the balloon shoulders. Expansion of the stent should not be undertaken if the stent is not properly positioned in the target lesion.

Note: If removal of a stent system is required prior to deployment, ensure that the guiding catheter is coaxially positioned relative to the stent delivery system, and cautiously withdraw the stent delivery system into the guiding catheter. Should unusual resistance be felt at any time when withdrawing the stent towards the guide catheter, the stent delivery system and the guiding catheter should be removed as a single unit. This should be done under direct visualization with fluoroscopy.

8. Tighten the rotating hemostatic valve. The stent is now ready to be deployed.

9.5 Deployment Procedure

CAUTION: Refer to product label for *in vitro* stent inner diameter, nominal pressure, and RBP.

1. Prior to deployment, reconfirm the correct position of the stent relative to the target lesion using the radiopaque balloon markers.
2. Deploy the stent slowly by pressurizing the delivery system in 2 atm increments, every 5 seconds, until stent is completely expanded. Fully expand the stent by inflating to nominal pressure at a minimum. Accepted practice generally targets an initial deployment pressure that would achieve a stent inner diameter ratio of about 1.1 times the reference vessel diameter (refer to product label for *in vitro* stent inner diameter, nominal pressure, and RBP).
3. For long lesions, size the stent to the diameter of the most distal portion of the vessel and expand stent to nominal pressure at minimum. Maintain pressure for 30 seconds. If necessary, the delivery system can be repressurized or further pressurized to assure complete apposition of the stent to the artery wall.
4. Maintain pressure for 30 seconds for full expansion of the stent. Fluoroscopic visualization during stent expansion should be used in order to properly judge the optimum stent diameter as compared to the proximal and distal native coronary artery diameters (reference vessel diameters). Optimal stent expansion and proper apposition requires that the stent be in full contact with the arterial wall.

Note: See *Section 9.6 Removal Procedure* for instruction on withdrawal of stent delivery system.

5. If necessary, the delivery system can be repressurized or further pressurized to assure complete apposition of the stent to the artery wall.

Note: Do not exceed the labeled rated burst pressure (RBP) of 18 atm (1824 kPa).

6. Fully cover the entire lesion and balloon-treated area (including dissections) with the XIENCE Alpine stent, allowing for adequate stent coverage into healthy tissue proximal and distal to the lesion.
7. Deflate the balloon by pulling negative on the inflation device for 30 seconds. Confirm complete balloon deflation before attempting to move the delivery system. If unusual resistance is felt during stent delivery system withdrawal, pay particular attention to guiding catheter position.

Note: See *Section 9.6 Removal Procedure* for instruction on withdrawal of stent delivery system.

8. Confirm stent position and deployment using standard angiographic techniques. For optimal results, the entire stenosed arterial segment should be covered by the stent. Fluoroscopic visualization during stent expansion should be used in order to properly judge the optimum expanded stent diameter as compared to the proximal and distal coronary artery diameter(s). Optimal expansion requires that the stent be in full contact with the artery wall. Stent wall contact should be verified through routine angiography or intravascular ultrasound (IVUS).
9. If the deployed stent size is still inadequate with respect to reference vessel diameter, a larger balloon may be used to further expand the stent. If the initial angiographic appearance is suboptimal, the stent may be further expanded using

a low profile, high pressure, noncompliant balloon dilation catheter. If this is required, the stented segment should be carefully recrossed with a prolapsed guide wire to avoid disrupting the stent geometry. Deployed stents should not be left underdilated.

CAUTION: Do not dilate the stent beyond the following limits.

<u>Nominal Stent Diameter</u>	<u>Dilation Limit</u>
2.0 – 2.5 mm	3.25 mm
2.75 – 3.25 mm	3.75 mm
3.5 – 4.0 mm	4.50 mm

10. If more than one XIENCE Alpine stent is needed to cover the lesion and balloon-treated area, it is suggested that, to avoid the potential for gap restenosis, the stents be adequately overlapped. To ensure that there are no gaps between stents, the balloon marker bands of the second XIENCE Alpine stent should be positioned inside the deployed stent prior to expansion.
11. Reconfirm stent position and angiographic results. Repeat inflations until optimal stent deployment is achieved.

9.6 Removal Procedure

Withdrawal of the stent delivery catheter from the deployed stent:

1. Deflate the balloon by pulling negative pressure on the inflation device. Larger and longer balloons will take more time (up to 30 seconds) to deflate than smaller and shorter balloons. Confirm balloon deflation under fluoroscopy and wait 10 to 15 seconds longer.
2. Position inflation device on "negative" or "neutral" pressure.
3. Stabilize guide catheter position just outside coronary ostium and anchor in place. Maintain guide wire placement across stent segment.
4. Gently remove the stent delivery system with slow and steady pressure.
5. Tighten the rotating hemostatic valve.

If during withdrawal of the stent delivery catheter resistance is encountered, use the following steps to improve balloon rewrap:

- Re-inflate the balloon up to nominal pressure.
- Repeat steps 1 through 5 above.

Post stent delivery system withdrawal; Stent deployment confirmation

1. Confirm stent position and deployment using standard angiographic techniques. For optimal results, the entire stenosed arterial segment should be covered by the stent. Fluoroscopic visualization during stent expansion should be used in order to properly judge the optimum expanded stent diameter as compared to the

proximal and distal coronary artery diameter(s). Optimal expansion requires that the stent be in full contact with the artery wall. Stent wall contact should be verified through routine angiography or intravascular ultrasound (IVUS).

2. If more than one XIENCE Alpine stent is needed to cover the lesion and balloon-treated area, it is suggested that, to avoid the potential for gap restenosis, the stents be adequately overlapped.
3. To ensure that there are no gaps between stents, the balloon marker bands of the second XIENCE Alpine stent should be positioned inside the deployed stent prior to expansion.
4. Reconfirm stent position and angiographic results to assess stented area. Repeat inflations until optimal stent deployment is achieved. If post-dilation is necessary, ensure that the final stent diameter matches the reference vessel diameter. Ensure that the stent wall is in contact with the artery wall.

9.7 Post-deployment Dilation of Stent Segments

1. All efforts should be taken to ensure that the stent is not under-dilated.
2. If the deployed stent size is still inadequate with respect to the vessel diameter, or if full contact with the vessel wall is not achieved, a larger balloon may be used to expand the stent further. The stent may be further expanded using a low profile, high pressure, and noncompliant balloon catheter. If this is required, the stented segment should be recrossed carefully with a prolapsed guide wire to avoid dislodging the stent. The balloon should be centered within the stent and should not extend outside of the stented region.

CAUTION: Do not dilate the stent beyond the following limits.

<u>Nominal Stent Diameter</u>	<u>Dilation Limit</u>
2.0 – 2.5 mm	3.25 mm
2.75 – 3.25 mm	3.75 mm
3.5 – 4.0 mm	4.50 mm

10.0 SPIRIT AND XIENCE FAMILY OF CLINICAL TRIALS

The XIENCE Alpine EECSS is based on the predicate devices XIENCE V EECSS, XIENCE PRIME EECSS and XIENCE Xpedition EECSS.

The XIENCE Alpine EECSS uses the identical stent platform, identical drug coating formulation, identical drug primer, identical nominal total drug content, and the identical stent contacting balloon materials as the XIENCE PRIME EECSS and XIENCE Xpedition EECSS.

The XIENCE Alpine EECSS differs from the XIENCE Xpedition EECSS only in the stent delivery system. The XIENCE Alpine stent delivery system utilizes the same principle of operation and materials as other Abbott Vascular RX stent systems and coronary dilatation catheters.

Compared to the XIENCE V EECSS, the XIENCE Alpine EECSS has the same stent security specification, same stent placement on the balloon between the balloon markers, a similar tip entry profile, and a similar taper length for XIENCE Alpine stent sizes up to 28 mm in length. Based on the identical nature of the XIENCE Alpine stent to the XIENCE PRIME stent, performance of the XIENCE Alpine EECSS can be predicted to be similar to the performance of XIENCE V and XIENCE PRIME. Therefore, clinical trial data for XIENCE V and XIENCE PRIME are summarized in this section.

10.1 Pre-market Clinical Trials

Principal safety and effectiveness of the XIENCE V stent has been established from a series of pre-market clinical trials. SPIRIT III RCT was the pivotal randomized clinical trial (RCT) that established the non-inferiority of the XIENCE V stent to the TAXUS® Express® stent (TAXUS stent). The SPIRIT IV trial is a prospective, randomized, active-controlled, single-blinded, multi-center evaluation of the XIENCE V stent compared to the TAXUS Express stent (TAXUS stent) in the treatment of up to three *de novo* lesions ≤ 28 mm in length in native coronary arteries with RVD ≥ 2.5 mm to ≤ 4.25 mm. The SPIRIT IV trial established the superiority of the XIENCE V stent compared to the TAXUS Express stent. The SPIRIT Small Vessel (SV) Registry is a prospective, single-arm, open-label, US multi-center registry study that established the safety and effectiveness of the XIENCE V 2.25 mm diameter stent. SPIRIT PRIME is a prospective, open-label, multi-center nonrandomized clinical trial with two study arms using the core size XIENCE PRIME and XIENCE PRIME LL stent system that established the safety and effectiveness of the XIENCE PRIME and XIENCE PRIME LL stents. Tables 10.1-1 through 10.1-4 present the trial designs, angiographic results (for studies that required angiographic follow-up), and principal clinical outcomes at 1-year and from latest follow-up, respectively.

Table 10.1-1: SPIRIT Family of XIENCE V Clinical Trial Designs (Pre-market)

	SPIRIT III RCT	SPIRIT IV	SPIRIT Small Vessel Registry	SPIRIT PRIME Clinical Trial	
				Core Size Registry	Long Lesion Registry
Study Type / Design	- Multi-center - Randomized - Single-blinded - Active-control	- Multi-center - Randomized - Single-blinded - Active-control	- Multi-center - Open-label - Single-arm	- Multi-center - Open-label - Single-arm	- Multi-center - Open-label - Single-arm
Number of Subjects Enrolled	Total: 1,002 XIENCE V: 668 TAXUS Express Control: 334	Total: 3,690 XIENCE V: 2,460 TAXUS Express Control: 1,230 ¹	Total: 150 2.25 mm XIENCE V	Total: 400 XIENCE PRIME	Total: 100 XIENCE PRIME
Treatment	Up to two <i>de novo</i> lesions in different epicardial vessels	Up to three <i>de novo</i> lesions, maximum of two lesions per epicardial vessel	Up to two <i>de novo</i> lesions in different epicardial vessels	Up to two <i>de novo</i> lesions in different epicardial vessels	Up to two <i>de novo</i> lesions in different epicardial vessels
Lesion Size	RVD: ≥ 2.5 ≤ 3.75 mm Length: ≤ 28 mm	RVD: ≥ 2.5 ≤ 4.25 mm ² Length: ≤ 28 mm	RVD: ≥ 2.25 < 2.50 mm Length: ≤ 28 mm	RVD: ≥ 2.25 ≤ 4.25 mm Length: ≤ 22 mm	XIENCE PRIME CS: RVD: ≥ 2.25 ≤ 4.25 mm Length: ≤ 22 mm XIENCE PRIME LL: RVD: ≥ 2.5 ≤ 4.25 mm Length: > 22 mm and ≤ 32 mm
Primary Endpoint	In-segment late loss at 240 days	Ischemia-driven target lesion failure at 1 year (composite of cardiac death, target vessel MI or ischemia-driven TLR)	TLF (target lesion failure) at 1 year	TLF (target lesion failure) at 1 year	TLF (target lesion failure) at 1 year
Co-Primary Endpoint	TVF at 270 days	None	None	None	None
Clinical Follow-up	30, 180, 240, 270 days, 1 to 5 years	30, 180, 270 days, 1 to 3 years	30 days, 240 days, 1 to 3 years	30, 180 days, 1 to 3 years	30, 180 days, 1 to 3 years
Angiographic Follow-up	240 days (N = 564)	None	240 days (N = 69)	None	None

¹ In the TAXUS stent arm, there was 1 subject who received 1 TAXUS® Liberté® stent.

² RVD ≥ 2.5 mm to ≤ 3.75 mm and stent sizes up to 3.5 mm until 4.0 mm TAXUS is commercially available.

**Table 10.1-2: SPIRIT Family of Clinical Trials
 Angiographic Results (Pre-market)**

Angiographic Results	SPIRIT III RCT 240 Days		SPIRIT Small Vessel 240 Days
	XIENCE V (N = 376) (M = 427)	TAXUS (N = 188) (M = 220)	2.25 mm XIENCE V (N = 69) (M = 69)
In-Stent Late Loss (mm)	0.16 ± 0.41 (342)	0.30 ± 0.53 (158)	0.20 ± 0.40 (52)
In-Segment Late Loss (mm)	0.14 ± 0.39 (343)	0.26 ± 0.46 (158)	0.16 ± 0.41 (52)
In-Stent Binary Restenosis	2.3% (8/343)	5.7% (9/158)	3.8% (2/52)
In-Segment Binary Restenosis	4.7% (16/344)	8.9% (14/158)	9.6% (5/52)

Notes:

- Data are mean (mm) ± SD or % (n/N).
- N is total number of patients. M is total number of lesions.
- SPIRIT III and SV 240-day include follow-up window (240 + 28 days).

**Table 10.1-3: SPIRIT Family of Clinical Trials
 Principal 1-Year Clinical Outcomes (Pre-market)**

	SPIRIT IV		SPIRIT III RCT		SPIRIT Small Vessel	SPIRIT PRIME Clinical Trial	
	XIENCE V (N = 2458)	TAXUS (N = 1229)	XIENCE V (N = 669)	TAXUS (N = 333)	2.25 mm XIENCE V (N = 144)	Core Size Registry (N = 401)	Long Lesion Registry (N = 104)
TLF	4.0% (97/2416)	6.8% (81/1195)	5.3% (35/655)	9.7% (31/319)	8.1% (11/136)	4.5% (18/399)	7.7% (8/104)
TVF	5.5% (134/2416)	7.7% (92/1195)	8.5% (56/655)	11.6% (37/319)	11.0% (15/136)	N/A	N/A
MACE	4.1% (98/2416)	6.9% (82/1195)	6.0% (39/655)	10.3% (33/319)	8.1% (11/136)	4.5% (18/399)	7.7% (8/104)
All Death	1.0% (25/2416)	1.3% (15/1195)	1.2% (8/657)	1.3% (4/320)	1.5% (2/136)	0.8% (3/399)	1.0% (1/104)
Cardiac Death	0.4% (10/2416)	0.4% (5/1195)	0.8% (5/657)	0.9% (3/320)	1.5% (2/136)	0.3% (1/399)	0.0% (0/104)
MI	1.9% (45/2416)	3.1% (37/1195)	2.7% (18/655)	4.1% (13/319)	1.5% (2/136)	1.8% (7/399)	4.8% (5/104)
Cardiac Death or MI	2.2% (54/2416)	3.3% (39/1195)	3.4% (22/655)	4.7% (15/319)	2.9% (4/136)	2.0% (8/399)	4.8% (5/104)
Ischemia-Driven TLR	2.3% (56/2416)	4.6% (55/1195)	3.4% (22/655)	5.6% (18/319)	5.1% (7/136)	2.5% (10/399)	2.9% (3/104)
Ischemia-Driven TVR, Non TL	2.2% (54/2416)	2.4% (29/1195)	3.2% (21/655)	4.7% (15/319)	5.9% (8/136)	2.8% (11/399)	2.9% (3/104)
Stent Thrombosis							
ARC (Definite / Probable)	0.29% (7/2391)	1.10% (13/1181)	0.9% (6/650)	0.6% (2/316)	1.5% (2/136)	0.5% (2/399)	0.0% (0/104)
ARC (Definite)	0.3% (6/2385)	0.8% (10/1183)	0.8% (5/650)	0.3% (1/317)	0.7% (1/138)	0.5% (2/399)	0.0% (0/104)

Notes:

- All counts presented in this table are subject counts. Subjects are counted only once for each event for each time period.
- 1-year includes the follow-up window (365 + 28 days) for all trials.
- TLF includes cardiac death, MI attributed to target vessel and ischemia-driven TLR. SPIRIT SV and PRIME used clinically indicated TLR definition rather than ischemia-driven TLR.

- TVF includes cardiac death, MI, ischemia-driven TLR and TVR, non-target lesion. SPIRIT SV and PRIME used clinically indicated TLR and TVR definition rather than ischemia-driven TLR and TVR definition, which was used for SPIRIT II, SPIRIT III, and SPIRIT IV.
- MACE includes cardiac death, MI and ischemia-driven TLR.

Table 10.1-4: SPIRIT Family of Clinical Trials
Principal Clinical Outcomes from Latest Follow-up (Pre-market)

	SPIRIT IV 3 Years		SPIRIT III RCT 5 Years		SPIRIT Small Vessel 2 Years	SPIRIT PRIME Clinical Trial 1 Year	
	XIENCE V (N = 2458)	TAXUS (N = 1229)	XIENCE V (N = 669)	TAXUS (N = 333)	2.25 mm XIENCE V (N = 144)	Core Size Registry (N = 401)	Long Lesion Registry (N = 104)
TLF	9.5% (223/2348)	11.9% (138/1158)	13.4% (81/605)	20.6% (59/286)	8.3% (11/133)	4.5% (18/399)	7.7% (8/104)
TVF	13.3% (312/2348)	14.5% (168/1158)	20.3% (123/605)	26.6% (76/286)	12.0% (16/133)	N/A	N/A
MACE	9.8% (231/2348)	12.3% (142/1158)	14.4% (87/605)	22.0% (63/286)	8.3% (11/133)	4.5% (18/399)	7.7% (8/104)
All Death	3.4% (81/2348)	5.2% (60/1158)	6.0% (37/621)	10.3% (31/300)	1.5% (2/133)	0.8% (3/399)	1.0% (1/104)
Cardiac Death	1.4% (34/2348)	1.9% (22/1158)	2.7% (17/621)	4.3% (13/300)	1.5% (2/133)	0.3% (1/399)	0.0% (0/104)
MI	3.1% (73/2348)	4.7% (55/1158)	4.6% (28/605)	7.0% (20/286)	1.5% (2/133)	1.8% (7/399)	4.8% (5/104)
Cardiac Death or MI	4.5% (105/2348)	6.0% (70/1158)	7.1% (43/605)	11.2% (32/286)	3.0% (4/133)	2.0% (8/399)	4.8% (5/104)
Ischemia-Driven TLR	6.3% (148/2348)	7.9% (92/1158)	8.9% (54/605)	12.9% (37/286)	5.3% (7/133)	2.5% (10/399)	2.9% (3/104)
Ischemia-Driven TVR, Non TL	5.6% (132/2348)	5.4% (63/1158)	8.8% (53/605)	11.9% (34/286)	6.8% (9/133)	2.8% (11/399)	2.9% (3/104)
Stent Thrombosis							
ARC (Definite / Probable)	0.62% (14/2263)	1.73% (19/1098)	1.5% (9/582)	1.9% (5/268)	1.5% (2/132)	0.5% (2/399)	0.0% (0/104)
ARC (Definite)	0.49% (11/2263)	1.28% (14/1098)	1.2% (7/582)	0.7% (2/268)	0.8% (1/132)	0.5% (2/399)	0.0% (0/104)

Notes:

- All counts presented in this table are subject counts. Subjects are counted only once for each event for each time period.
- Data includes the follow-up window of + 28 days for all trials.
- TLF includes cardiac death, MI attributed to target vessel and ischemia-driven TLR. SPIRIT SV and PRIME used clinically indicated TLR definition rather than ischemia-driven TLR.
- TVF includes cardiac death, MI, ischemia-driven TLR and TVR, non-target lesion.
- MACE includes cardiac death, MI, and ischemia-driven TLR.

10.2 Post-Market Clinical Trials

The XIENCE V USA study is a prospective, multi-center, FDA-mandated post-market study to evaluate the continued safety and effectiveness of the XIENCE V EECSS in real-world settings after it was commercialized in the U.S., and also to support the FDA DAPT initiative. The objective of the SPIRIT V Single Arm Study (SAS) was to continue the assessment of the performance of XIENCE V EECSS in the treatment of patients with *de novo* coronary artery lesions. XIENCE V India is a prospective, open-label, multi-center, observational, single-arm registry to evaluate XIENCE V EECSS's continued safety and effectiveness during commercial use in real-world settings in India. Tables 10.2-1 through 10.2-2 present clinical trial designs, principal clinical outcomes at 1 year and from latest follow-up, respectively.

The results from these post-marketing clinical trials demonstrate safety and effectiveness of XIENCE V in real-world settings. In addition, XIENCE V improved patient-reported outcomes (including better quality of life, reduced angina frequency, improved angina stability, and reduced physical limitation) at 6 months and the improvements were sustained for 1 year in coronary artery disease patients.

Table 10.2-1: SPIRIT and XIENCE V Family Clinical Trial Designs (Post-market)

	XIENCE V USA Phase I Cohort	XIENCE V USA Long- Term Follow- up Cohort	XIENCE V USA AV-DAPT Cohort	SPIRIT V (SAS)	XIENCE V India
Study Type / Design	- Multi-center - Prospective - Single-arm	- Multi-center - Prospective - Single-arm	- Multi-center - Randomized - Double-blinded - Placebo control	- Multi-center - Prospective - Single-arm	- Multi-center - Prospective - Single-arm
Number of Subjects Enrolled	8040	4663	868	2663	977
Treatment	Per site standard care	Per site standard care	Patients were randomized to receive either thienopyridine or placebo treatment for additional 18 months along with aspirin	Maximum of one <i>de novo</i> , native target lesion per major epicardial vessel or side branch (no prior stent implant, no prior brachytherapy), maximum of 4 planned EES	Per site standard care
Lesion Size	No angiographic restrictions			RVD = ≥ 2.25 ≤ 4.0 mm Length ≤ 28 mm by visual estimation	No angiographic restrictions
Primary Endpoint	ARC definite and probable stent thrombosis up to 1 year	ARC definite and probable stent thrombosis from year 1 to 5	MACE (composite of all death, MI, and stroke) 12–33 months	Composite rate of all death, MI, TVR at 30-day	ARC ST (Definite / probable) 1 year and yearly after through 3 years
Co-Primary Endpoint	Cardiac death or any MI at 1 year	Cardiac death or any MI from year 1 to 5	ARC definite and probable ST 12– 33 months	None	Cardiac death and any MI at 1 year
Clinical Follow-up	14, 30, 180 days, and 1 year	2, 3 and 4 years	15, 24, 30, and 33 months	30 days and 1 and 2 years	14, 30, 180 days and 1, 2, and 3 years
Angiographic Follow-up	None	None	None	None	None

Table 10.2-2: SPIRIT and XIENCE V Family of Trials
Principal Clinical Outcomes (Post-market)

	XIENCE V USA Phase I 1 Year	SPIRIT V (SAS)		XIENCE V India	
		1 Year	2 Years	1 Year	2 Years
	XIENCE V (N = 8040)	XIENCE V (N = 2663)	XIENCE V (N = 2663)	XIENCE V (N = 990)	XIENCE V (N = 990)
TLF (ARC)	9.4% (707/7522)	5.25% (138/2627)	7.49% (192/2562)	2.4% (24/986)	3.4% (32/942)
TLF	6.8% (513/7505)	N/A	N/A	N/A	N/A
All Death, MI (ARC) and TVR	N/A	7.04% (185/2627)	10.34% (265/2562)	N/A	N/A
Cardiac Death or MI (ARC)	7.2% (545/7522)	4.23% (111/2627)	5.74% (147/2562)	1.9% (19/986)	2.8% (26/942)
Cardiac Death or MI	3.3% (249/7505)	N/A	N/A	N/A	N/A
TLR	4.6% (349/7522)	1.90% (50/2627)	3.04% (78/2562)	1.2% (12/986)	1.5% (14/942)
TVR, non TLR	2.3% (176/7522)	1.45% (38/2627)	2.26% (58/2562)	0.1% (1/986)	0.1% (1/942)
All Death	2.6% (194/7522)	1.71% (45/2627)	2.97% (76/2562)	0.9% (9/986)	1.7% (16/942)
Cardiac Death	1.4% (108/7522)	1.10% (29/2627)	1.87% (48/2562)	0.9% (9/986)	1.7% (16/942)
All MI (ARC)	6.3% (475/7522)	3.54% (93/2627)	4.45% (114/2562)	1.3% (13/986)	1.6% (15/942)
All MI	2.2% (162/7505)	N/A	N/A	N/A	N/A
Stent Thrombosis					
ARC (Definite / Probable)	0.81% (60/7380)	0.65% (17/2607)	0.79% (20/2523)	0.51% (5/986)	0.53% (5/939)
ARC (Definite)	0.54% (40/7380)	N/A	N/A	0.41% (4/986)	0.43% (4/939)

Notes:

- All counts presented in this table are subject counts. Subjects are counted only once for each event for each time period.
- XIENCE V USA 1-year data include 42-day window or through randomization date if occurred earlier than 407 days for the second enrollment phase. XIENCE V India 1-year data include 47-day window. SPIRIT V 1-year data include 28-day window.
- SPIRIT V 2-year data include 28-day window. XIENCE V India 2-year data include 47-day window.
- TLF (ARC) includes cardiac death, MI attributed to target vessel (per ARC definition), clinically indicated TLR. TLF includes cardiac death, MI attributed to target vessel (per protocol definition), clinically indicated TLR.

10.3 Pooled Clinical Trials Analysis for Stent Thrombosis Subsequent to DAPT Interruption / Discontinuation

The XIENCE family of stents have demonstrated safety and effectiveness in multiple trials in standard / low-risk and extended-risk real-world patients and demonstrated low stent thrombosis (ST) rates.

This includes patients who discontinued or interrupted dual antiplatelet therapy (DAPT) after 1 month post coronary stent implantation. Clinical trials described below and / or in the published literature have been analyzed for stent thrombosis correlated with DAPT. The data (Table 10.3-1) are from a pooled patient level analysis of the following pre- and post-market trials:

Pre-market: SPIRIT II, SPIRIT III, SPIRIT IV

Post-market: SPIRIT V, SPIRIT Women, XIENCE V USA, XIENCE V India

Table 10.3-1: Stent Thrombosis (ARC Definite / Probable) Subsequent to DAPT Interruption / Discontinuation through 730 days Versus No Interruption; Pooled Data from 7 Pre- and Post-market Trials

	Pooled 7 Trials (2 year results) N = 11,219		
	DAPT Interrupted ≤30 days	DAPT Interrupted >30–730 days	No Interruption of DAPT
Subsequent ST (Def / Prob)	2.51% (12/478)	0.38% (13/3445)	0.66% (44/6648)

11.0 STERILIZATION CONDITIONS

The XIENCE Alpine Everolimus Eluting Coronary Stent System is supplied sterile.

The XIENCE Alpine Everolimus Eluting Coronary Stent System is sterilized with ethylene oxide gas.

Disposable medical product. Do not re-sterilize. Do not use if the package is open or damaged. Medical device is non-pyrogenic.

Each lot of products produced by ABBOTT VASCULAR is inspected to be in compliance with product specifications including sterilization process validation. A Declaration of Conformity with recognized global standards is issued by ABBOTT VASCULAR. Products which are not compliant with product specifications and the required global standards cannot be distributed.

12.0 PACKAGING, USAGE

Each product is individually packed.

The XIENCE Alpine Stent System has a protective sheath that covers the stent/balloon area. The protective sheath (manufactured in-house by Abbott Vascular) is a polymer-

based material (composed of Linear Low Density Polyethylene (LLDPE) used to cover and protect the drug coated stent. The finished product stent system is inserted into a high-density polyethylene (HDPE) dispenser coil. The dispenser coil is placed into a Polyethylene (LLDPE)/Nylon/Polyethylene (LLDPE)/Tyvek® inner header bag, which is heat sealed and given a product label. The inner header bag with a Tyvek portion allows for steam, ethylene oxide (EO), and oxygen gas exchange during sterilization while maintaining the sterile barrier. Header bags are placed into individual cartons and then a corrugated box for shipping and EO sterilization. An Instruction for Use booklet (in Russian) is placed in each product carton post sterilization, during addendum labeling prior to shipping product to the region.

The product packaging ensures sterility and functionality until the expiry date stated on the product labels. Labels are placed on the primary and secondary packaging for accurate identification of the product, determination of its size, the lot number and expiration date. The labels also include symbols regarding product warnings and precautions.

13.0 TRANSPORTATION AND STORAGE CONDITIONS

Devices packaged by manufacturer may be transported by all modes of transport in covered vehicles in accordance with the rules of transportation of goods for each type of transport. The temperature at transportation should be between -20°C to +50°C with a range of $\pm 2^\circ\text{C}$ and Relative Humidity from 30% to 75% with a range of $\pm 5\%$.

Products should be stored at 15°C – 30°C in an enclosed space in the original manufacturer's packaging for the labeled shelf life.

Shelf life of the XIENCE Alpine Everolimus Eluting Coronary Stent System is 3 years (36 months).

14.0 ENVIRONMENT PROTECTION REQUIREMENTS

The XIENCE Alpine Everolimus Eluting Coronary Stent System does not adversely impact the environment

15.0 DISPOSAL

Products, which had contact with blood and/or other biological liquids of patients should be disposed and eliminated in accordance with the requirements for handling medical waste, established by SanPiN 2.1.7.2790-2010 for class B and other normative legal acts of the Russian Federation.

Unused products, including the ones which expired (which had no contact with the biological liquids of patients) should be disposed in accordance with the procedure of disposal of solid domestic waste established by the legislation of the Russian Federation for class A by SanPiN 2.1.7.2790.

16.0 DISCLAIMER OF WARRANTY AND LIMITATION OF REMEDY

Abbott Vascular has exercised reasonable care in the manufacture of this device; however, numerous factors beyond Abbott Vascular's control can affect this device and the results obtained from its use. **ABBOTT VASCULAR, THEREFORE, MAKES NO WARRANTIES REGARDING THIS DEVICE AND EXPRESSLY DISCLAIMS ALL WARRANTIES, WHETHER EXPRESS OR IMPLIED BY OPERATION OF LAW OR OTHERWISE, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. ABBOTT VASCULAR SHALL NOT BE LIABLE TO ANY PERSON OR ENTITY FOR ANY INCIDENTAL OR CONSEQUENTIAL LOSS, DAMAGE, OR EXPENSE, DIRECTLY OR INDIRECTLY ARISING FROM THE USE OF THIS DEVICE. NO PERSON HAS ANY AUTHORITY TO ASSUME ANY LIABILITY ON BEHALF OF ABBOTT VASCULAR WITH RESPECT TO THE DEVICE OR TO MAKE ANY GUARANTEE, REPRESENTATION, OR WARRANTY WITH RESPECT TO THE DEVICE.** The exclusions and limitations set forth above are not intended to and should not be construed so as to contravene mandatory provisions of applicable law. If any term of this Disclaimer of Warranty is held to be illegal, unenforceable or in conflict with applicable law by a court of competent jurisdiction, the validity of the remaining portions of this Disclaimer of Warranty shall not be affected.

Descriptions or specifications in any Abbott Vascular printed material are intended solely to generally describe the device at the time of manufacture and do not constitute any express warranties.

ABBOTT VASCULAR WILL NOT BE RESPONSIBLE FOR ANY DIRECT, INCIDENTAL OR CONSEQUENTIAL LOSS, DAMAGE, OR EXPENSE RESULTING FROM REUSE, REPROCESSING, OR RE-STERILIZATION OF THIS DEVICE.

17.0 PATENTS AND TRADEMARKS

This product and / or its use may be covered by one or more of the following United States Patents: 6,179,810;; 6,440,990; 6,629,994; 6,656,220; 6,746,423; 6,887,219; 6,887,510; 6,890,318; 6,929,657; 6,939,373; 6,957,152; 7,549,975; 7,828,766; 7,833,193; 7,906,066; 7,947,207; 8,043,553; 8,052,638; 8,075,583; 8,173,062; 8,221,112; 8,221,444; 8,308,711; 8,382,738; 8,388,575; 8,388,602; 8,394,055; 8,444,608; 8,444,802; 8,535,596; 8,540,927; 8,613,722. Other U.S. patents pending. Other patents issued and pending outside of the U.S.

CERTICAN, AFINITOR, VOTUBIA and ZORTRESS are trademarks of the NOVARTIS group of companies.

AMERICAN COLLEGE OF CARDIOLOGY is a trademark of the American College of Cardiology Foundation.

AMERICAN HEART ASSOCIATION is a trademark of the American Heart Association, Inc.

TAXUS and TAXUS LIBERTÉ are trademarks of the Boston Scientific group of companies.

XIENCE, XIENCE V, XIENCE PRIME, XIENCE PRO, XIENCE Xpedition, and XIENCE Alpine are trademarks of the Abbott group of companies.

Abbott Vascular
3200 Lakeside Drive
Santa Clara, CA 95054 USA

CUSTOMER SERVICE

TEL: (800) 227-9902
FAX: (800) 601-8874
Outside USA TEL: (951) 914-4669
Outside USA FAX: (951) 914-2531

18.0 GRAPHICAL SYMBOLS FOR MEDICAL DEVICE LABELING

 Manufacturer	 Catalogue number	 French size	 EXCURSIONS Excursions permitted to temperature range
	 Sterilized using ethylene oxide	 Consult instructions for use	
 Use by	 Batch code	 Do not reuse	 Temperature limitation
 Date of manufacture	 Inner diameter	 Outer diameter	
 Non-pyrogenic	 Contents (numeral represents quantity of units inside)	 Guiding catheter	 Flushing tool
 MR Conditional	 Do not resterilize	 Stent length	 Do not use if package is damaged

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©2018 Abbott