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OPERATOR MANUAL for Medical Device: Arctic Sun 5000 system and accessories for management of normothermia and hypothermia

Manufacturer: Medivance Inc., USA

State of Colorado
County of Boulder

Signed before me on March 14, 2022
by Scott Light

Leanne Sloan
(Notary's official signature)

Notary Public
(Title of office)

October 04, 2022
(Commission Expiration)

LEANNE SLOAN
NOTARY PUBLIC
STATE OF COLORADO
NOTARY ID 20184039166
MY COMMISSION EXPIRES OCTOBER 04, 2022

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М.П. / Stamp
« March 14 » 2022г.

Medivance

a wholly owned subsidiary of
Bard Medical Division
C.R. Bard, Inc.

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The Arctic Sun 5000 Temperature Management System is a thermal regulating system, indicated for monitoring and controlling temperature patients of all ages.

3.3. Classification of the device

Class of potential risk of use of the medical device	26
Product classification by type and duration of contact	System and accessories –no direct contact with patient skin, pads -long-term up to 30 days.
Frequency of use	Multiple use medical device, pads – single use
Software classification (ISO 62304)	Class B
Protection against electric shock	System – class 1, accessories – type BF
protection against ingress of water and particulate matter IPX0	IPX0
Mode of operation	Continuous
ME EQUIPMENT and ME SYSTEMS intended for use in an OXYGEN RICH ENVIRONMENT	Not intended for use in an-oxygen rich environment
Software	version 3.0.2, March 2018

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3.4 WARNINGS

- Do not use the Arctic Sun system in the presence of flammable agents because an explosion and/or fire may result.
- Do not use high frequency surgical instruments or endocardial catheters while the Arctic Sun 5000 Temperature Management System is in use.
- There is a risk of electrical shock and hazardous moving parts. There are no user serviceable parts inside. Do not remove covers. Refer servicing to qualified personnel.

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Power cord has a hospital grade plug. Grounding reliability can only be achieved when connected to an equivalent receptacle marked "hospital use" or "hospital grade".

When using the Arctic Sun 5000 system, note that all other thermal conductive systems, such as water blankets and water gels, in use while warming or cooling with the Arctic Sun 5000 Temperature Management System may actually alter or interfere with patient temperature control. Do not place ArcticGel Pads over transdermal medication patches as warming can increase drug delivery, resulting in possible harm to the patient.

The Arctic Sun 5000 system is not intended for use in the operating room environment.

5 CAUTIONS

This product is to be used by or under the supervision of trained, qualified medical personnel.

Use only sterile water. The use of other fluids will damage the Arctic Sun 5000 Temperature Management System.

When moving the Arctic Sun 5000 Temperature Management System always use the handle to lift the controller over an obstacle to avoid over balancing.

The patient's bed surface should be located between 75 cm and 150 cm above the floor to ensure proper flow and minimize risk of leaks.

The clinician is responsible to determine the appropriateness of custom parameters. When the system is powered off, all changes to parameters will revert to the default unless the new settings have been saved as new defaults in the Advanced Setup screen. For small patients (≤ 30 kg) it is recommended to use the following settings: Water Temperature High Limit $\leq 40^{\circ}\text{C}$; Water Temperature Low Limit $\geq 10^{\circ}\text{C}$; Control Strategy = 1. It is recommended to use the Patient Temperature High and Patient Temperature Low alert settings.

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- The operator must continuously monitor patient temperature when using Manual Control and adjust the temperature of the water flowing through the pads accordingly. Patient temperature will not be controlled by the Arctic Sun 5000 Temperature Management System in Manual Control.

- Due to the system's high efficiency, Manual Control is not recommended for long duration use. The operator is advised to use the automatic therapy modes (e.g. Control Patient, Cool Patient, Rewarm Patient) for automatic patient temperature monitoring and control.

- The Arctic Sun 5000 Temperature Management System will monitor and control patient core temperature based on the temperature probe attached to the system. The clinician is responsible for correctly placing the temperature probe and verifying the accuracy and placement of the patient probe at the start of the procedure.

- Manufacturer supplies temperature simulators (fixed value resistors) for testing, training and demonstration purposes. Never use this device, or other method, to circumvent the normal patient temperature feedback control when the system is connected to the patient. Doing so exposes the patient to the hazards associated with severe hypo- or hyper-thermia.

- Manufacturer recommends measuring patient temperature from a second site to verify patient temperature. Medivance recommends the use of a second patient temperature probe connected to the Arctic Sun 5000 Temperature Management System Temperature 2 input as it provides continuous monitoring and safety alarm features. Alternatively, patient temperature may be verified periodically with independent instrumentation.

- The displayed temperature graph is for general information purposes only and is not intended to replace standard medical record documentation for use in therapy decisions.

- Patient temperature will not be controlled, and alarms are not enabled in Stop Mode. Patient temperature may increase or decrease with the Arctic Sun 5000 Temperature Management System in Stop Mode.

- Carefully observe the system for air leaks before and during use. If the pads fail to prime or a significant continuous air leak is observed in the pad return line, check connections. If needed, replace the leaking pad. Air leakage may result in lower flow rates and potentially decrease the performance of the system.

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The ArcticGel Pads are non-sterile for single patient use. Do not reprocess or sterilize. If used in a sterile environment, pads should be placed according to the physician's request, either prior to the sterile preparation or sterile draping. ArcticGel Pads should not be placed on a sterile field.

• Use pads immediately after opening. Do not store pads once the kit has been opened.

• Do not place ArcticGel Pads on skin that has signs of ulceration, burns, hives, or rash.

• While there are no known allergies to hydrogel materials, caution should be exercised with any patient who has a history of skin allergies or sensitivities. Examine the packaging before use. Do not use in case of damage.

• Do not allow circulating water to contaminate the sterile field when patient lines are disconnected.

• The water content of the hydrogel affects the pad's adhesion to the skin and conductivity, and therefore, the efficiency of controlling patient temperature. Periodically check that pads remain moist and adherent. Replace pads when the hydrogel no longer uniformly adheres to the skin. Replacing pads at least every 5 days is recommended.

• Do not puncture the ArcticGel Pads with sharp objects. Punctures will result in air entering the fluid pathway and may reduce performance.

• If accessible, examine the patient's skin under the ArcticGel Pads often, especially those at higher risk of skin injury. Skin injury may occur as a cumulative result of pressure, time and temperature. Possible skin injuries include bruising, tearing, skin ulcerations, blistering, and necrosis. Do not place bean bag or other firm positioning devices under the ArcticGel Pads.

Do not place positioning devices under the pad manifolds or patient lines.

• The rate of temperature change and potentially the final achievable patient temperature is affected by many factors.

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Treatment application, monitoring and results are the responsibility of the attending physician.

If the patient does not reach target temperature in a reasonable time or the patient is not able to be maintained at the target temperature, the skin may be exposed to low or high water temperatures for an extended period of time which may increase the risk for skin injury. Ensure that pad sizing/coverage and custom parameter settings are correct for the patient and treatment goals, water flow is greater than or equal to 2.3 liters per minute and the patient temperature probe is in the correct place. For patient cooling, ensure environmental factors such as excessively hot rooms, heat lamps, and heated nebulizers are eliminated and patient shivering is controlled. Otherwise, consider increasing minimum water temperature, modifying target temperature to an attainable setting or discontinuing treatment. For patient warming, consider decreasing maximum water temperature, modifying target temperature to an attainable setting or discontinuing treatment.

• Due to underlying medical or physiological conditions, some patients are more susceptible to skin damage from pressure and heat or cold. Patients at risk include those with poor tissue perfusion or poor skin integrity due to diabetes, peripheral vascular disease, poor nutritional status, steroid use or high dose vasopressor therapy. If warranted, use pressure relieving or pressure reducing devices under the patient to protect from skin injury.

• Do not allow urine, antibacterial solutions or other agents to pool underneath the ArcticGel Pads. Urine and antibacterial agents can absorb into the pad hydrogel and cause chemical injury and loss of pad adhesion. Replace pads immediately if these fluids come into contact with the hydrogel.

• Do not place ArcticGel Pads over an electrosurgical grounding pad. The combination of heat sources may result in skin burns.

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Power On

- 1) Turn the power ON by activating the Power Switch.
- 2) The control module will automatically go through a brief self-test of the independent safety alarm.
- 3) A New User Training module option is available from the start up screen.
- 4) When the self-test is complete, the **Patient Therapy Selection** screen will appear on the control panel.

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Fill Reservoir

- 1) Fill the reservoir with sterile water only.
- 2) Four liters of water will be required to fill the reservoir at initial installation.
- 3) Add one vial of Arctic Sun 5000 Temperature Management System Cleaning Solution to the sterile water.
- 4) From the **Patient Therapy Selection** screen, press either the **Normothermia** or **Hypothermia** button, under the New Patient heading.

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6) The **Fill Reservoir** screen will appear. Follow the directions on the screen.

Functional Verification

Perform the following functional verification procedure after initial setup and installation of the control module.

1) Power **On** the control module

2) From the **Patient Therapy Selection** screen, press the **Hypothermia** button to display the **Hypothermia** therapy screen.

3) From the **Hypothermia** therapy screen, press the **Manual Control** button to open the **Manual Control** window.

4) Use the Up and Down arrows to set the **Manual Control** water target temperature to 40°C and the duration to 30 minutes.

5) Press the **Start** button to initiate **Manual Control**. Allow at least 3 minutes for the system to stabilize.

6) Monitor the flow rate and water temperature in the **System** status area on the **Hypothermia** therapy screen.

7) Verify that the flow rate reaches at least 1.5 liters/minute.

8) Verify that the water temperature increases to 30°C.

9) Press the **Stop** button.

10) Set the Manual Control water target temperature to 4°C and the duration to 30 minutes.

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5. Attach the pad's line connectors to the patient line manifolds. Begin circulating water through the pads using either Patient Temperature Control Mode (automatic) or Water Temperature Control Mode (manual). If the pads fail to prime or a significant continuous air leak is observed in the pad return line, check connections, then if needed replace the leaking pad.

6. Once the pads are primed, assure the flow rate displayed on the control panel is greater than 1.7 liters per minute, which is the minimum flow rate for a full pad kit (4).

7. When finished, empty water from pads. Cold temperature increases the adhesiveness of the hydrogel. For ease of removal, leave pads on the patient for approximately 15 minutes to allow the hydrogel to warm. Slowly remove pads from the patient and discard.

Directions for ArcticGel Small Universal Pad

• The ArcticGel Small Universal Pads are only for use with the Arctic Sun Temperature Management System. See Operators Manual for detailed instructions on system use.

• Select the proper number of pads for the patient size and clinical indication. However, the rate of temperature change and potentially the final achievable temperature is affected by pad surface area, patient size, pad placement and water temperature range. Best system performance will be achieved by using the maximum number of pads.

• The following sizing chart is provided as guidance. Modify the number of pads and locations as necessary to meet specific clinical needs.

Patient Weight	Number of Pads
2.5 - 5 kg (5.5 - 10 lbs)	1
5 - 10 kg (10 - 22 lbs)	2
10 - 16 kg (22 - 35 lbs)	3

• Place the pads on healthy, clean skin only. Locate pad edges away from articulating areas of the body to avoid irritation. Place pads to allow for full respiratory excursion.

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• The pad surface must be contacting the skin for optimal energy transfer efficiency. The Small Universal Pads are provided with a cloth liner over the hydrogel. The pads may be used with the cloth liner in place or removed to expose the hydrogel adhesive. If the cloth liner is used, secure the pads with the Velcro straps to ensure good contact with the skin. The pads may be removed and reapplied if necessary.

• Due to the small patient size and the potential for rapid patient temperature change

• It is recommended to use the following settings:

Water Temperature High Limit: $\leq 40^{\circ}\text{C}$ (104°F)

Water Temperature Low Limit: $\geq 10^{\circ}\text{C}$ (50°F)

Control Strategy: 2

• It is recommended to use the Patient Temperature High and Patient Temperature Low alert settings.

• Place a core patient temperature probe and connect to the Arctic Sun Temperature Management System patient temperature 1 connector for continuous patient temperature feedback. A rectal or esophageal temperature probe is recommended.

• Attach the pad's line connectors to the fluid delivery line manifolds. Begin treating the patient.

• If the pads fail to prime or a significant continuous air leak is observed in the pad return line, check connections, then if needed replace the leaking pad.

• Once all the pads are primed, assure the steady state flow rate displayed on the control panel is appropriate for the number of pads connected.

Number Pads 1 2 3

Minimum Flow Rate l/m 0.9 1.7 2.4

• When finished, purge water from pads. Slowly remove pads from the patient and discard.

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Temperature Probe Placement

Patient temperature control with the Arctic Sun 5000 Temperature Management System requires patient temperature feedback provided by an indwelling patient temperature probe connected to the Patient Temperature I connector on the back of the Control Module.

Patient Therapy Selection

Use the **Patient Therapy Selection** screen to initiate a **New Patient**, **Continue a Current Patient**, or access the **Advanced Setup** screen.

New Patient - Normothermia

Select **Normothermia** if the therapy goal is to maintain a patient temperature at a pre-defined target temperature for an indefinite period of time. Press the **Normothermia** button to display the **Normothermia** therapy screen.

New Patient - Hypothermia

Select **Hypothermia** to reduce and maintain a patient temperature at a set target temperature for a defined period of time, then slowly re-warm the patient at a controlled re-warming rate. Press the **Hypothermia** button to display the **Hypothermia** therapy screen.

Additional Protocol Option

Two additional protocols (**Hypothermia** or **Normothermia**) may be visible on the **Patient Therapy Selection** screen.

Current Patient

The **Continue Current Patient** button and the date and time that the current therapy was paused will display on the **Patient Therapy Selection** screen if a patient therapy was paused within the past 6 hours.

Press the **Continue Current Patient** button to resume a paused patient therapy.

Initiate Normothermia (Control Patient & Rewarm Patient)

Normothermia therapy is initiated and managed, and patient temperature is automatically controlled to a set target temperature from the **Control Patient** window in the **Normothermia** therapy screen. The **Control Patient** window displays the patient target temperature and the duration since the initiation of normothermia therapy.

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To initiate Normothermia therapy:

- 1) From the **Patient Therapy Selection** screen, press the **Normothermia** button to display the **Normothermia** therapy screen.
- 2) The default patient target temperature will display in the **Control Patient** window.
- 3) To modify the patient target temperature, press the **Adjust** button to display the **Control Patient-Adjust** window.
- 4) **Control Patient** to: Use the Up and Down arrows to set the desired patient target temperature to control the patient.
- 5) **Rewarm at a Rate** of: Use the Up and Down arrows on the right of the screen to set the rewarming rate.
- 6) Press the **Save** button to save the new settings and close the **Control Patient-Adjust** window
- 7) Press **Start**, in the **Control Patient** window to initiate therapy. You will hear a tone and then a voice stating "Therapy Started". Additionally, the **Control Patient** window and the Arctic Sun 5000 Temperature Management System icon will blink, indicating that therapy is in progress.

Initiate Hypothermia (Cool Patient and Rewarm Patient)

Hypothermia therapy is initiated and managed, and patient temperature is automatically controlled to a set target temperature from the **Cool Patient** and **Rewarm Patient** windows in the **Hypothermia** therapy screen.

The **Cool Patient** window displays the cooling phase patient target temperature and the length of time remaining in the cooling phase of the **Hypothermia** therapy.

The **Rewarm Patient** window displays the rewarming phase patient target temperature and the length of time remaining in the rewarming phase of the **Hypothermia** therapy.

To initiate hypothermia therapy:

From the **Patient Therapy Selection** screen, press the **Hypothermia** button to display the **Hypothermia** therapy screen.

1. Cool Patient Settings

- The default patient target temperature and duration will display in the **Cool Patient** window.
- To modify the patient target temperature and duration, press the **Adjust** button to display the **Cool Patient-Adjust** window.

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• **Cool Patient To:** Use the Up and Down arrows on the left side to set the desired patient target temperature to cool the patient

• **Cool Patient For:** Use the Up and Down arrows on the right side to set the cooling duration to cool the patient before rewarming begins.

• Press the **Save** button to save the new settings and close the **Cool Patient-Adjust** window

2. Rewarm Patient Settings

• The default patient target temperature and duration will display in the **Rewarm Patient** window.

• To change the rewarming phase patient target temperature and rewarming rate, press the **Adjust** button in the **Rewarm Patient** window

to display the **Rewarm Patient-Adjust** screen. Use the Up and Down arrows on the left side to set the desired final patient target temperature.

• **Rewarm Patient To:** Use the Up and Down arrows on the right side to set the desired final patient target temperature.

• **Rewarm at a Rate of:** Use the Up and Down arrows in the center of the screen to set the rewarming rate.

• **Rewarm Patient From:** When cooling a patient, adjustment of the **Rewarm Patient From** setting on the left side of the screen is disabled and defaults to the **Cool Patient** target temperature.

• When rewarming a patient, the **Rewarm Patient From** adjustment is enabled and the value can be modified. The **Rewarm Patient From** setting is the temperature to which the system is currently controlling the patient. The **Rewarm Patient From** temperature will automatically increase as the rewarming process continues. This feature allows the rewarming procedure to be optimized by allowing complete control of the rewarming ramp.

• Using the **Rewarm Patient From** temperature, the **Rewarm Patient To** temperature and the rewarming rate settings, the system will calculate and display the rewarming duration and the date/time at which the patient will reach the final rewarming target temperature.

• Press the **Save** button to save the new settings and close the **Rewarm Patient-Adjust** window.

3. Initiate Patient Cooling

• Press **Start**, in the **Cool Patient** window to initiate therapy. You will hear a tone and then a voice stating "Therapy Started". Additionally, the **Cool Patient** window and the Arctic Sun 5000 Temperature Management System icon will blink, indicating that therapy is in progress.

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4. Initiate Patient Rewarming

• Upon completion of the cooling phase, there are two options for initiation of patient rewarming, either Automatically or Manually, depending on the **Rewarming Begins** setting in **Hypothermia Settings**.

• If **Rewarming Begins** is set to **Automatically**, the rewarming process starts automatically when the **Cool Patient** therapy is complete and the duration reaches zero.

• If **Rewarming Begins** is set to **Manually**, the rewarming process starts when the **Start** button is pressed in the **Rewarm Patient** window. The

cooling process will continue until the **Rewarm Patient Start** button is pressed. An **Alert** will occur when the **Cool Patient** duration reaches zero. When the **Rewarm Patient** duration timer reaches zero, the system will continue to control the patient to the target temperature until the **Stop** button is pressed. Once in **Normothermia**, the timer will reset and begin tracking the duration of **Normothermia** therapy.

End Therapy

• From the **Normothermia** therapy or **Hypothermia** therapy screen, press the **Stop** button to terminate water circulation to the pads.

• Press the **Empty Pads** button, and follow the instructions to purge the pads of water.

• Disconnect the pads from the Fluid Delivery Line.

• Slowly and carefully remove pads from the patient skin.

• Discard the used pads in accordance with hospital procedures for medical waste.

• Press the power switch **Off**.

If power is lost while the power switch is in the On position, an audible alert will be issued until it is switched Off. This alerts the user that the treatment may have been accidentally stopped.

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
Normothermia Settings

Use the **Normothermia Settings** screen to view the current settings and modify the settings for the following parameters. To modify any parameter setting, press the **Adjust** button to the right of the parameter.

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To access the Normothermia Settings screen:

- 1) Press **Adjust** on the **Control Patient** window.
- 2) Press the **More** button on the **Control Patient Adjust** window.
- 3) The **Normothermia Settings** screen will be displayed.
- 4) To save the new settings as the current patient therapy settings, press the **Close** button. For instructions on saving the settings as the system defaults, see **Advanced Setup**.

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
NOTE: Verify that the Patient Temperature 1 probe is correctly positioned and correctly connected to the system. If the patient temperature change is less than 0.15°C after the first hour of therapy, the system will generate *Alert 116 – Patient Temperature 1 Change Not Detected*. If Alert 116 is not acknowledged after 5 minutes, the system will generate *Alarm 117 – Patient Temperature 1 Change Not Detected*. Alarm 117 will stop therapy and an audible alarm will sound. After Alarm 117 is acknowledged, therapy will require a restart.

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Hypothermia Settings

Use the **Hypothermia Settings** screen to view the current settings and modify the settings for the following parameters. To modify any parameter setting, press the **Adjust** button to the right of the parameter.

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- Cooling Begins
- Rewarming Begins

Water Temperature Settings

- Pre-Condition Water
- Manual Control
- High Water Limit
- Low Water Limit

Patient Temperature Settings

- High Patient Alert
- Low Patient Alert
- Control Strategy

Display Settings

- Temperature Units
- Temperature Units Adjust
- Patient Temp 2

To access the Hypothermia Settings screen:

- 1) Press **Adjust** on the **Cool Patient** window or the **Rewarm Patient** window.
- 2) Press the **More** button on the **Cool Patient Adjust** window or **Rewarm Patient Adjust** window.
- 3) The **Hypothermia Settings** screen will be displayed.

4) To save the new settings as the current patient therapy settings, press the **Close** button. For instructions on saving the settings as the system defaults, see **Advanced Setup**.

NOTE: Verify that the Patient Temperature 1 probe is correctly positioned and correctly connected to the system. If the patient temperature change is less than 0.15°C after the first hour of therapy, the system will generate **Alert 116 – Patient Temperature 1 Change Not Detected**. If Alert 116 is not acknowledged after 5 minutes, the system will generate **Alarm 117 – Patient Temperature 1 Change Not Detected**. Alarm 117 will stop therapy and an audible alarm will sound. After Alarm 117 is acknowledged, therapy will require a restart.

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4.4. Advanced Setup

Use the **Advanced Setup** screen to view the current settings and modify the settings for the following parameters. To modify any parameter setting, press the **Adjust** button to the right of the parameter.

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Alarms

An Alarm is denoted by an audio signal that repeats every 10 seconds until the Alarm is cleared. The Alarm screen will appear that displays the alarm number, alarm title, a description of the problem or conditions that triggered the alarm, and solutions and instructions for troubleshooting and resolving the alarm condition. If certain Alarm conditions are not acknowledged by the operator within 2 minutes, a Reminder tone will sound. All Alarm settings are maintained in the event of a power interruption.

Main Safety Alarms

While there are multiple alarms and safety features in the Arctic Sun 5000 Temperature Management System, there are five main safety alarms that will place the device into Stop mode until the condition is addressed.

Alarm Specification

High Patient Temperature 39.5°C

Low Patient Temperature 31.0°C

High Water Temperature 42.5°C / 44°C

Low Water Temperature 3.0°C / 3.5°C

System Self-Test Failure At device power ON

Each time the Arctic Sun 5000 Temperature Management System is powered On, a system self test for the independent safety alarm is automatically run.

This test simulates a “water high temperature” fault situation on both the primary and secondary water temperature sensors. Both the primary and secondary safety systems must respond to the fault and be verified by the opposing safety system. If either safety systems do not respond appropriately either an alarm 80 or 81 will be issued. Contact Customer Support.

Non-Recoverable Alarms

If an Alarm condition occurs that prevents proper use of the device or proper patient treatment (such as the five main safety alarms discussed above), the system is placed into Stop mode and will not allow therapy to continue.

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This type of Alarm is known as Non-Recoverable. If this situation occurs, cycle the device power (turn device Off then On). If the alarm recurs contact Customer Support.

Recoverable Alarms

Other Alarms that temporarily Stop the device until the user is able to correct the cause and clear the Alarm are classified as Recoverable. If the condition that initiated the alarm is not addressed and problem persists, the Alarm will recur.

If a Recoverable Alarm occurs:

- 1) When an alarm is issued the device is placed into **Stop** mode.
- 2) Read the displayed instructions.
- 3) Note the **Alarm** number.
- 4) Press the **Close** button to clear the alarm.
- 5) Follow the instructions to correct the alarm condition. Perform the actions in the order listed until the alarm condition is resolved.
- 6) Once you have cleared the alarm, press the Start button in the therapy window to restart therapy. You will hear a tone and a voice stating "Therapy Started". Additionally, the active therapy window and the Arctic Sun 5000 Temperature Management System icon will blink.
- 7) If the condition does not resolve, contact Customer Support.

Alerts

Alerts are denoted by an audio signal that repeats every 25 seconds. The Alert screen will appear that displays the alert number, alert title, a description of the problem that triggered the alert, and solutions and instructions for troubleshooting and resolving the alert condition.

If an Alert occurs:

- 1) Read the displayed instructions.
- 2) Note the Alert number.
- 3) Press the **Close** button to clear the alert.
- 4) Follow the instructions to correct the alert condition. Perform the actions in the order listed until the alarm condition is resolved. If the condition does not resolve, contact Customer Support.

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5) Refer to the Arctic Sun 5000 Temperature Management System Help screens for additional information regarding alarms and alerts.

5. MAINTENANCE AND SERVICE

Maintenance Schedule

Procedure	Interval	Who performs
Clean external surfaces	As required	Medical personnel
Inspect connectors and cables	6 month	Authorized service support
Clean the condenser	6 month	Authorized service support
Replenish Cleaning solution	6 month	Authorized service support
Inspect screen protector	6 month	Authorized service support
Calibration	Every 2000 hours or 250 uses, whichever occurs first, as indicated by system display	Authorized service support
Inspect fluid delivery line	6 month	Authorized service support

Cleaning and Maintenance

Routine cleaning and preventive maintenance should be performed on the Arctic Sun 5000 Temperature Management System control module every 6 months at a minimum. This consists of cleaning the external surfaces, accessories and chiller condenser, inspecting the device, and replenishing the internal cleaning solution that suppresses microorganism growth in the water reservoir and hydraulic circuit.

External Surfaces

Clean the exterior body of the control module, fluid delivery lines, power cords and temperature cables using a

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soft cloth and mild detergent or disinfectant according to hospital protocol.

Medivance has qualified and approves the use of the following types of disinfectants for exterior surfaces: sodium hypochlorite, isopropyl alcohol, and quaternary ammonium

Condenser

- A dirty chiller condenser will significantly reduce the cooling capacity of the control module.

- To clean the condenser, wipe the dust from the exterior grill using a soft cloth. Depending on the quality of your institution's air, periodically remove the back cover and vacuum or brush the condenser fins. At a minimum the condenser fins should be cleaned annually. Maintenance activities should be performed by qualified personnel.

Device Inspection

- Periodically inspect the external areas of the device for damaged, loose or missing parts, and frayed or twisted power cords and cables.

- Discontinue using the device displaying one or more of the above conditions until the problem is corrected and has been verified to be operating correctly.

Replenish Internal Cleaning Solution

Contact Customer Service to order internal cleaning solution.

To replenish the internal cleaning solution:

- 1) Drain the reservoir.

- Turn control module power Off.

- Attach the drain line to the two drain ports on the back of the control module. Place the end of the drain line into a container. The water will passively drain into the container.

- 2) Refill the reservoir.

- From the Hypothermia therapy screen or the Normothermia therapy screen, press the Fill Reservoir button.

- The Fill Reservoir screen will appear. Follow the directions on the screen.

- Add one vial of Arctic Sun 5000 Temperature Management System cleaning solution to the first bottle of sterile water.

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The filling process will automatically stop when the reservoir is full. Continue to replace the bottles of sterile water until the filling process stops.

When the Fill Reservoir process is complete, the screen will close.

Inspect Screen Protector

The Control Panel's touchscreen is supplied with a disposable screen protector. If it becomes damaged, it can be removed by lifting the edge and carefully peeling it from the screen. To ensure dust and particulates are removed, clean the touchscreen using isopropyl alcohol. Remove the blue liner from the screen protector. Then carefully apply the protector to the screen with the liner side down against the screen.

Inspect Fluid Delivery Line

1. Power On the system
2. From the Patient Therapy Selection screen press the Hypothermia button to display the Hypothermia therapy screen.
3. From the Hypothermia therapy screen, press the Manual Control button to open the Manual Control window.
4. Set the Manual Control water target temperature to 28°C and the duration to 30 minutes.
5. Connect a shunt to a set of fluid delivery line ports.
6. Press the Help button and then press the Help Index button. Select the topic Maintenance and Service and subtopic System Diagnostics then press the Display button. Verify that inlet pressure is -7 ± 0.2 .
7. Repeat on all valves. If inlet pressure is out of range, replace the two valves that the shunt is connected to.
8. Ensure that the shunt is removed before device is put back in service.

Preventative Maintenance

Use of the Arctic Sun 5000 Temperature Management System in excess of 2,000 hours, without conducting preventative maintenance, may result in failure of certain system components and failure of the system to function as intended. To maintain system performance, the Arctic Sun 5000 Temperature Management System requires periodic service of key components.

Preventive maintenance should be performed by Authorized representatives of the service support. Service includes

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regulating of electronic components, cleaning solution replenish, change of old components, software update (if necessary).

Contact AR "Bard Rus" LLC to confirm contacts of Authorized service provider.

Software Update

Software updates will be provided on a flash drive. Installation of software updates will be performed via the USB port on the front of the control module. The software update feature will automatically initiate if the control module detects the proper files on a flash drive inserted in the USB port at power On.

To install software update:

- 1) Insert flash drive into the USB port.
- 2) An image of a timer will display while the software update is being installed, and will disappear when the software installation process is completed.
- 3) After installation, the new software version will display in the Software Version field in Advanced Setup.

Service

Contact AR "Bard Rus" LLC for technical support and customer service instructions to enable appropriately qualified technical personnel to repair those parts of the equipment that are considered repairable.

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5.1.1. Calibration Test Unit

A separate device, the Calibration Test Unit (CTU) is required to perform periodic calibration on the Arctic Sun 5000 Temperature Management System.

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When to Perform a Calibration or Calibration Check

1. Calibration is recommended after 2000 hours of operation or 250 uses, whichever occurs first. The status of Calibration is available in the Advanced Settings screen.

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2. In addition, Calibration may be required after replacing certain components (Manifold Harness, Tank Temperature Sensor Harness, Processor Circuit Card, Isolation Circuit Card, Input/Output Circuit Card).

3. A Calibration Check confirms the device flow, ability to heat and cool, and the temperature sensing systems are all within specification. During the calibration check errors may be displayed with diagnostic information assisting with performance or calibration issues. After successful completion of a calibration check, a report is displayed showing a pass or fail status of all parameters checked.

5.1.2. Performing a Calibration

To perform a calibration on the Arctic Sun 5000 Temperature Management System, press the Advanced Setup button on the Therapy Selection screen. Press the Start button next to Calibration and follow the on-screen instructions.

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6. ELECTROMAGNETIC COMPATABILITY

Medical electrical equipment needs special precautions regarding electromagnetic compatibility. Ensure that the Arctic Sun 5000 Temperature Management System is installed and used according to the electromagnetic compatibility information provided. The following are guidance and manufacturer's declarations regarding electromagnetic compatibility for the Arctic Sun 5000 Temperature Management System.

- The use of accessories or cables other than those specified or sold by Medivance (shown below) is not recommended. Use of unapproved accessories or cables may result in increased emissions or in decreased immunity of the Arctic Sun 5000 Temperature Management System.

- If the Arctic Sun 5000 Temperature Management System is used directly adjacent to or stacked with other equipment, the user should periodically observe the Arctic Sun 5000 Temperature Management System device to verify it operates normally in that environment.

- Portable and mobile RF communications equipment can affect Medical Electrical Equipment.

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Accessories and cables approved by Medivance for use with the Arctic Sun® Temperature Management System	Part #
Temperature In-Cable - Nellcor	735-02
Temperature In-Cable - Bard®	735-03
Temperature In-Cable - Rusch	735-04
Temperature In-Cable - GE	735-05
Temperature In-Cable - Philips	735-06
Temperature Out-Cable - Nellcor	735-52
Temperature Out-Cable - Bard®	735-53
Temperature Out-Cable - Rusch	735-54
Temperature Out-Cable - GE	735-55
Temperature Out-Cable - Philips	735-56
Power Cord, US, Canada, Mexico	733-00
Power Cord, Continental Europe	733-01
Power Cord, UK, Ireland	733-02
Power Cord, Australia, New Zealand	733-03
Power Cord, Mainland China	733-04
Power Cord, Brazil	733-05
Power Cord, Switzerland	733-07
Power Cord, South Africa	733-08
Transmission Interface Module (TIM) Kit	760-00
Transmission Interface Module (TIM)	761-00
RS232 Cord	762-00

1.1 EN/IEC 60601-1-2 Table 1

Guidance and Manufacturer's Declaration - Electromagnetic Emissions

The Arctic Sun® Temperature Management System is intended for use in the electromagnetic environment specified below. The customer or the end user of the Arctic Sun® Temperature Management System should assure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The Arctic Sun® Temperature Management System uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.

1.1 EN/IEC 60601-1-2 Table 1 (cont.)

Guidance and Manufacturer's Declaration - Electromagnetic Emissions

Emissions Test	Compliance	Electromagnetic environment guidance
RF emissions CISPR 11	Class A	The Arctic Sun® Temperature Management System unit is suitable for use in all establishments other than domestic, establishments and those directly connected to the public low-voltage power supply network that supplies buildings for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ Flicker emissions IEC 61000-3-3	Complies	

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12 EN/IEC 60601-1-2 Table 2			
Guidance and Manufacturer's Declaration – Electromagnetic Immunity			
The Arctic Sun® Temperature Management System unit is intended for use in the electromagnetic environment specified below. The customer or the end user of the Arctic Sun® Temperature Management System unit should assure it is used only in such an environment.			
Immunity Test	IEC60601 test level	Compliance Level	Intended Electromagnetic Environment
Electromagnetic discharge (ESD) IEC 61000-4-2	±8kV contact ±8kV air	±8kV contact ±8kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	±2kV for power supply lines ±1kV for input/output lines	±2kV for power supply lines ±1kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1kV differential mode (line-line) ±2kV common mode (line-earth)	±1kV differential mode (line-line) ±2kV common mode (line-earth)	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% UT (>95% dip in UT) for 0.5 cycle 40% UT (60% dip in UT) for 5 cycles 70% UT (30% dip in UT) for 25 cycles <5% UT (>95% dip in UT) for 5 seconds	<5% UT (>95% dip in UT) for 0.5 cycle 40% UT (60% dip in UT) for 5 cycles 70% UT (30% dip in UT) for 25 cycles <5% UT (>95% dip in UT) for 5 seconds	Mains power quality should be that of a typical commercial or hospital environment. If the user of the Arctic Sun® Temperature Management System unit requires continued operation during power mains interruptions, it is recommended that the Arctic Sun® Temperature Management System unit be powered from an uninterruptible power supply with sufficient capacity to run the unit for the maximum required time of interruption.
Power frequency (50/60Hz) magnetic field IEC 61000-4-8	3A/m	3A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Note: UT is the a.c. mains voltage prior to application of the test level.

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1.3 EN/IEC 60601-1-2:2007 Sub-clause 5.2.2.2 Table 4:

Guidance and Manufacturer's Declaration – Electromagnetic Immunity

The Arctic Sun® Temperature Management System unit is intended for use in the electromagnetic environment specified below. The customer or the end user of the Arctic Sun® Temperature Management System unit should assure it is used in such an environment.

Immunity test	IEC60601 test level	Compliance Level	Intended Electromagnetic Environment
Conducted RF EC 61000-4-8 Radiated RF EC 61000-4-3	3Vrms 150kHz to 80MHz 3V/m 80MHz to 2.5GHz	3Vrms 150kHz to 80MHz 3V/m 80MHz to 2.5GHz	Portable and mobile RF communications equipment should be used no closer to any part of the Arctic Sun® Temperature Management System unit, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = 1.2 \sqrt{P}$ $d = 1.2 \sqrt{P}$ 80MHz to 800 MHz $d = 2.3 \sqrt{P}$ 800MHz to 2.5GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended minimum separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ^a , should be less than the compliance level in each frequency range. ^b Interference may occur in the vicinity of equipment marked with the following symbol: 

NOTE 1 At 80MHz and 800MHz, the higher frequency range applies

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from objects, structures and people.

Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Arctic Sun® Temperature Management System unit is used exceeds the applicable RF compliance level above, the Arctic Sun® Temperature Management System unit should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the Arctic Sun® Temperature Management System unit.

^a Over the frequency range 150kHz to 80MHz, field strengths should be less than 3V/m.

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
7. MEDICAL DEVICE PACKAGING

The Arctic Sun 5000 System is packaged in a custom-designed configuration that has been validated to meet packaging and shipping requirements. The following diagrams show the construction of the device package.

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Package is made from corrugated carton 200 # C Flute Brown.

Box is shipped on the wooden pallet.

Accessories are packed in the individual pouches made of polyethylene and then box of whited carton.

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Parameter	Specification
Therapy Modes	Normothermia: Control Patient, Rewarm Patient Hypothermia: Cool Patient, Rewarm Patient
Heater Capacity	2500 BTU/hr / 750 Watts
Circulating Fluid	Sterile Water
Reservoir Capacity	3.5 liters
Water Flow Rate	5 liters per minute
Patient Probe Type	YSI 400 Series compatible
Patient Temperature Inputs	Patient Temp 1: control, monitor, alarm Patient Temp 2: monitor, alarm
Patient Temperature Display Range	10°C to 44°C 50°F to 111.2°F 0.1°C /°F increments
Patient Temperature Measurement Accuracy	±0.4°C (10°C to 32°C) ±0.2°C (32°C to 38°C) ±0.4°C (38°C to 44°C) Includes ±0.1°C external probe
Responses of the PCLCS (Physiologic Closed-Loop Control System)	Settling Time: ~4.5 hrs Relative Overshoot: <0.5°C Command Overshoot: <0.5°C Response Time: Warming (max) 33°C to 37°C: ~6 hrs Cooling 37°C to 33°C: ~2 hrs Steady State Deviation: 0 Tracking Error: 0 Note: All values derived from testing in simulated use.
Patient Temperature Control Range	32°C to 38.5°C 89.6°F to 101.3°F 0.1°C /°F increments
Water Temperature Display Range	3°C to 45°C / 37.4°F to 113.0°F 0.1°C /°F increments
Water Temperature Control Range (Manual)	4°C to 42°C / 39.2°F to 107.6°F 1°C /°F increments
High Water Temperature Limit	36°C to 42°C / 96.8°F to 107.6°F 1°C /°F increments
Low Water Temperature Limit	4°C to 25°C / 39.2°F to 77°F 1°C /°F increments
Time to heat water from 20°C to 37°C	8 minutes (approximate)
Sound Pressure	Alarm Tone: 70dB to 80dB at 1 meter, repeats every 10 seconds Alert Tone: 63dB to 71dB at 1 meter, repeats every 25 seconds Reminder Tone: 65dB at 3 meters, 0.5 seconds on/20 seconds off
Mains Input	100-120VAC, 50-60Hz, 11A 220-240VAC, 50-60Hz, 5.5A
Leakage Current	<300 µA
Operating Relative Humidity Range	5% to 70% non-condensing
Storage Relative Humidity Range	5% to 95% non-condensing
Operating Temperature Range	10°C to 27°C / 50°F to 80°F
Storage Temperature Range	-30°C to 50°C / -20°F to 120°F
Atmospheric Pressure Range	60 kPa to 110 kPa
Dimensions	Height: 35 inches (89 cm) Width: 14 inches (36 cm) Depth: 18.5 inches (47 cm)
Weight	Empty: 43 kg / 95 lbs ; Filled: 47 kg / 103 lbs

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Arctic Sun 5000 temperature management system is indendent for use with ArcticGel pads only.

- Arctic Sun 5000 system can be used with temperature probes YSI 400 Series compatible and probes of other manufacturers in accordance with the technical specifications and instruction for use. Temperaure probes are not supplied with the sytem.

Table to of temperature and resistance correlation for comparatable sensors:

Temperature, °C	Resistance, Ω
32	1,667
33	1,599
34	1,533
35	1,471
36	1,412
37	1,355
38	1,301
39	1,249

Cable for sensor YSI 400 series 1 pc. is included into the kit of the device and is supplied with the ArcticSun 5000 system. For connection of the sensors cabels with the connectors GE, Philips, Rusch are used and ordered separately and are not supplied with the system.

List of the temperature sensors compatable with the ArcticSun 5000 system consist of, but not limited to the following manufacturers:

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
- Temperature probe sensor for general indication 9F/12F (not more than 10 pcs) – accessory to the medical device “Patient monitor Carescape B450 with accessories”. Registration certificate RZN 2016/3956 date 15.04.2016, p/78, Manufacturer: GE Healthcare Finland Oy, Kuortaneenkatu 2, 00510 Helsinki, Finland
- Sensor for general indication, series 400, pediatric/adult (not more than 10 pc.) accessory to the medical device “Patient monitor Carescape B450 with accessories”. Registration certificate RZN 2016/3956 date 15.04.2016, p/78, Manufacturer: GE Healthcare Finland Oy, Kuortaneenkatu 2, 00510 Helsinki, Finland
- Temperature sensors multiple use (not more than 5 pc.), Accessory to the medical device “Module Patient monitor IntelliVue Mx with accessories”, p.46, Registration certificate FSZ 2012/13060 date 13.05.2021, manufacturer: Philips Medizin systeme Boblingen GmbH, Hewlett-Packard-Strasse 2, 71034 Boblingen, Germany
- Temperature sensors single use (not more than 5 pc.), Accessory to the medical device “Module Patient monitor IntelliVue Mx with accessories”, p.47, Registration certificate FSZ 2012/13060 date 13.05.2021, manufacturer: Philips Medizin systeme Boblingen GmbH, Hewlett-Packard-Strasse 2, 71034 Boblingen, Germany
- Temperature sensors YSI 400 series, 5 types, manufacturer “YSI Incorporated”, USA, accessory to the medical device “Resuscitation and anesthesiology Monitor for control of the of the physiology parameters MITAR-01-“R-D” per TU 9441-002-24149103-2002”, registration certificate FSR 2012/13343 date 19.02.2021, Manufacturer “NPP Monitor” LLC, 344068, Russia, Rostov-on-Don, Krasnokursantskaya street, 104a

Ой”, Финляндия, GE Healthcare Finland Oy, Kuortaneenkatu 2, 00510 Helsinki, Finland

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


FORM 007-00014

1.1 Symbols present on the system, pads and accessories labeling:

	For the safe and effective use of this device, the operator must consult the accompanying documents prior to use.
	Indicates the Authorized Representative in the European Community.
	This symbol adjacent to the patient connections means that the thermal probe connection is a "Defibrillator-Proof, Type BF Applied Part", per standard IEC 60601-1 and affords the degree of patient protection defined in that standard for this type of applied part.
	Per ETL Intertek, models of the Arctic Sun® Temperature Management System that bear the ETL Monogram conform to AAMI ES60601-1, IEC 60601-1-8, IEC 60601-10, IEC 60601-2-35 and are certified to CSA C22.2 No. 60601-1.
	Indicates high temperature part or component. The maximum temperature of this internal component allowed by the protective system is also listed.
	Indicates that only sterile water should be used when filling the Arctic Sun® Temperature Management System Control Module.
	Identifies Patient Temperature 1, the patient temperature probe input for monitoring and control.
	Identifies Patient Temperature 2, the patient temperature probe input for monitoring
	Identifies Patient Temperature Out, the patient temperature output to an external hospital monitor.
	Identifies the drain port.
	Indicates electrical hazard
	Identifies the storage temperature range.
	Identifies the storage relative humidity range.
	Manufacturer
	Date of Manufacture
	Do not re-use.
	Risk of overbalance due to pushing, leaning, resting, etc.
	Arctic Sun® Temperature Management System must be disposed of properly. DO NOT dispose of into the garbage.
	Identifies mechanical hazard

RxOnly	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Unit quantity
REF	Catalog number
SN	Serial number
	CE mark
	Fragile
	Upward
	Keep dry
	Use by date
	Time of application
	Do not reuse
	Do not puncture or cut
	MR-safe
	Lot
	Non-sterile
	Not made with natural latex
	Consult instruction for use
	Medical device

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
BD Switzerland Sàrl
Terre Bonne Park - A4
Route de Crassier 17
1262 Eysins, Switzerland
Made in Mexico
PK 7650756 03/2020

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BD Switzerland Sàrl
Terre Bonne Park - A4
Route de Crassier 17
1262 Eysins, Switzerland
Made in Mexico
PK 7650758 03/2020

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BD Switzerland Sàrl
Terre Bonne Park - A4
Route de Cressier 17
1282 Eysins, Switzerland
Made in Mexico
PK 7640760 03/2020

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BD Switzerland Sàrl
Terre Bonne Park - A4
Route de Cressier 17
1282 Eysins, Switzerland
Made in Mexico
PK 7640767 03/2020

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
EC REP 80 Switzerland Sàrl
Terre Bonne Park - A4
Route de Crassier 17
1262 Eysins, Switzerland
Made in Mexico
PK 7650794 03/2020

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

EC REP 80 Switzerland Sàrl
Terre Bonne Park - A4
Route de Crassier 17
1262 Eysins, Switzerland
Made in Mexico
PK 7650766 03/2020

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Small Universal

ARCTICGEL™ Pad

ARCTIC SUN®
TEMPERATURE MANAGEMENT SYSTEM

Small Universal	Lille universel	通用小	Žgħir Universal
Universelle petite	Small universal	소형 공용	Universală mică
Klein (Universalgröße)	Pieni universaali	通用小	Malý univerzálny
Small universale	Liten universal	Малък, универсален	Majhno univerzalno
Pequeño universal	Malý univerzálny	Malo univerzalno	
Klein universeel	Extra univerzális	Väike kõigile sobiv	
Universal pequeno	Malý univerzální	Mazs universāls	
Μικρό γενικής χρήσης	Универсальный малый	Mažas Universalus	

2.5-16 kg
(5.5-35 lbs.)

Medivance®

REF 318-01

LOT NGWW0000

2020-12-28

Medivance, Inc.
321 South Taylor Avenue, Suite 200
Louisville, CO 80027 USA
Phone: +1-303-926-1917
Toll-Free: +1-800-526-4455

BD Switzerland Sàrl
Terre Bonne Park – A4
Route de Crassier 17
1262 Eysins, Switzerland
Made in Mexico
PK 1650768 03/2020

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OPERATOR MANUAL for Medical Device:

Arctic Sun 5000 system and accessories for management of normothermia and hypothermia

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Signed before me on 14th, 2022
by Scott Light

[Signature]
(Notary's official signature)

Notary Public
(Title of office)

October 04, 2024
(Commission Expiration)

LEAHNE SLOW
NOTARY PUBLIC
STATE OF COLORADO
MY COMMISSION EXPIRES OCTOBER 04, 2024

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[Signature]
(Indonesian/signature)

M. II. / Stamp
« March 14 » 2022r.

Medivance
 a wholly owned subsidiary of
 Bard Medical Devices
 C.A. Bard Inc.
 321 South Taylor Ave., Suite 200
 Louisville, CO 80027
 (303) 826-1917

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