Approved by CERTEST BIOTEC S.L.
Pol. Industrial Río Gállego II, Calle J, №1
50840, San Mateo de Gállego, Zaragoza (Spain)
The place of signature Zaragoza
CEO Nelson Gonçalves Fernandes
Company stamp

11/12/22

certest

Certest Biotec, S.L.
Pol. Ind. Río Gállego II, 1
50840 San Mateo de Gállego
Zaragoza (Spain)
Tel. (+34) 976 520 354
www.certest.es

CIF: B-50946854

Nelson Gonçalves Fernandes

Package insert

a latex turbidimetric assay for the quantitative *in vitro* detection of haemoglobin (faecal occult blood) in human stool samples (FOB Turbilatex), in versions, manufactured by CERTEST BIOTEC S.L., Spain

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PAPEL EXCLUSIVO PARA DOCUMENTOS NOTIFICADOS

HF6784593

Package insert of a latex turbidimetric assay for the quantitative *in vitro* detection of haemoglobin (faecal occult blood) in human stool samples (FOB Turbilatex), in versions.

Medical product name and related implementations

FOB Turbilatex

FOB Turbilatex is produced in next options:

FOB Turbilatex, version I, consists of:

1. FOB Turbilatex Reagent 1 – 1x68 mL;
2. FOB Turbilatex Reagent 2 – 1x17 mL;
3. FOB Turbilatex Calibrator 0 – 1x0,5 mL;
4. FOB Turbilatex Calibrator 1 – 1x0,5 mL;
5. FOB Turbilatex Calibrator 2 – 1x0,5 mL;
6. FOB Turbilatex Calibrator 3 – 1x0,5 mL;
7. FOB Turbilatex Calibrator 4 – 1x0,5 mL;
8. FOB Turbilatex Calibrator 5 – 1x0,5 mL;
9. FOB Turbilatex Control 1 – 9x0,5 mL;
10. FOB Turbilatex Control 2 – 9x0,5 mL;
11. Instruction for use – 1 pc.;
12. Value sheet – 1 pc.;
13. Bags – 3 pcs

14. Universal Turbilatex Sample Collection Vial for *in vitro* diagnosis - if necessary, in versions:

I. Universal Turbilatex Sample Collection Vial- no more 1000 pcs. , consists of:

- 1) Universal Turbilatex Sample Collection Vial- no more 1000 pcs.
- 2) Leaflet– 1 pc.

II. Universal Turbilatex Sample Collection Vial- no more 1000 pcs., consists of:

- 1) Universal Turbilatex Sample Collection Vial- no more 1000 pcs.;
- 2) Caps- no more 1000 pcs.;
- 3) Leaflet - 1 pc.

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IV. FOB Turbilatex, version IV, consists of:

1. FOB Turbilatex Reagent 1 – 1x20 mL;
2. FOB Turbilatex Reagent 2 – 1x5 mL;
3. Instruction for use – 1 pc.
4. Universal Turbilatex Sample Collection Vial for in vitro diagnosis - if

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Medical product purpose

FOB Turbilatex is intended for the quantitative determination of hemoglobin in human fecal samples in biochemical analyzers by the latex immunoturbidimetric method in clinical diagnostic laboratories and research practice.

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Scope

Clinical *in vitro* laboratory diagnostics, research practice.

Biochemical analysers used

FOB Turbilatex fob turbilatex is intended for use on automatic photometers (biochemical analyzers) registered and approved for circulation in the territory of the Russian Federation, providing immunoturbidimetric analysis based on measuring the intensity of light with a wavelength of 505 nm - recommended, secondary -800 nm, passed through a cuvette containing a colloidal solution (suspension formed by particles). The FOB Turbilatex kit is intended for use on the biochemical analyzers presented in the table.

Name of the medical device	Biochemical Analyzer
FOB Turbilatex	JEOL BioMajesty
	Furuno CA
	Dirui CS
	iMagic-S7
	DDS-240
	URIT-8280

Functional use

FOB Turbilatex is an auxiliary tool to make a diagnosis of lost gastrointestinal blood (gastrointestinal bleeding), diagnosis of precancer intestinal losses and

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of early-stage colorectal cancer by detecting the loss of intact or nearly intact haemoglobin from the intestine.

Medical Product Potential consumers Information

FOB Turbilatex is intended to be used only for *in vitro* diagnostics in clinical diagnostic laboratories and research practice. The product should be used only by instructed and trained personnel (specially qualified specialists: clinical laboratory diagnostics doctor, medical laboratory technician (paramedic-laboratory assistant), clinician, researcher).

General description and medical product structure

FOB Turbilatex is produced in various options. Complete set of FOB Turbilatex includes the following components:

Reagents (Reagent 1 and Reagent 2);

FOB Turbilatex Calibrator (calibrator 0 to 5) – liquid calibrator containing protein (human hemoglobin) at different concentrations, concentration is indicated on the instruction for use of the product, value sheet and on the label of each vial;

FOB Turbilatex Control 1 and 2 – liquid control materials containing protein (human hemoglobin) at different concentrations. Concentration is indicated on the label of each vial and value sheet.

Universal Turbilatex Sample Collection Vial for *in vitro* diagnostics

FOB Turbilatex components and concentrations:

Reagent 1	Buffer	<300 mM, pH 8,4
	Polyester	<100 г/л
	Sodium Azide	<0,1%
Reagent 2	Buffer	<100 mM, pH 8
	Protein	<1%

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Haemoglobin Concentration in calibrations:

Calibrator	Concentration
FOB Turbilatex Calibrator 0	0 ng/mL
FOB Turbilatex Calibrator 1	50 ng /mL
FOB Turbilatex Calibrator 2	100 ng /mL
FOB Turbilatex Calibrator 3	250 ng /mL
FOB Turbilatex Calibrator 4	500 ng /mL
FOB Turbilatex Calibrator 5	1000 ng /mL

Target analyte description and evaluations of detection

There is normally no blood in the feces. Its appearance may indicate a pathological process in organism, including a malignant neoplasm in the

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digestive system. In the early stages, pathology often develops without obvious symptoms, and it can be difficult to notice occult blood without tests. Faecal analysis of occult blood allows to detect the smallest blood impurities, thus, diagnose diseases in the early stages, that is extremely important if cancer is suspected (in particular, colorectal cancer — CRC). CRC affects both men and women of all racial and ethnic groups, most often occurs in people aged 50 and older. Men CRC is the most cancer spread type after prostate and lung cancer, women CRC — after breast and lung cancer. Screening with faecal occult blood tests is based on the concept that important target colonic neoplasm, such as early-stage cancer and large adenomatous polyps. Colorectal cancer is also associated with local acute inflammatory reaction being visualized, in some cases, by white cell neutrophil scanning.

In some studies it has been found 5.1 micrograms/g of faecal hemoglobin are equal or exceed cut-off value allowing the identification of gastrointestinal bleeding in adults, for its diagnosis extra diagnostic invasive procedures are required.

When gastrointestinal blood is lost, the stool will contain a combination of intact or nearly intact haemoglobin, intact heme and heme-derived porphyrins in amounts that depend on the site and amount of bleeding and the transit time through the gut. Haemoglobin is the iron-containing oxygen-transport protein in the red blood cells of all vertebrates that may be leaked into gastrointestinal tract and then discharged with the faeces in gastrointestinal bleeding diseases. FOB Turbilatex tests detect intact or nearly intact human haemoglobin, being a very specific technique for detecting loss of blood from the lower intestine.

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Indications and contraindications

FOB Turbilatex is used as an auxiliary tool to make a diagnosis of lost gastrointestinal blood (gastrointestinal bleeding), diagnosis of precancer intestinal losses and of early-stage colorectal cancer. For diagnostic purposes, the results should always be associated with the patient's medical background, anamnesis, lab results and other available data. The test result is not a diagnosis and should be interpreted only by a doctor.

There are no contraindications to medical kit usage.

Assay Sample Type

FOB Turbilatex is assayed in human stool sample.

Specimen collection and stability

1. Collect enough quantity of faeces in clean and dry containers (no preservatives or transport media).
2. Fecal samples should not be collected during menstrual period, bleeding haemorrhoids, blood in urine or they have strained during bowel movement.
3. If not immediately tested, freeze the stored samples at 2-8°C for 3 days or at -20°C maximum 6 months.
4. Before testing frozen samples are totally thawed, and brought to room temperature (15-30°C). Homogenise stool samples as thoroughly as possible prior to preparation.
5. Do not use stool specimens beyond the expiry allowed. Avoid storage and usage of contaminated samples. Repeated samples freezing is not allowed.

Equipment and materials required for FOB Turbilatex using

Materials and instruments that could be required:

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- Lab centrifugation;
- Vortex;
- Automated micropipettes with removable tips;
- Disposable gloves;
- Freezer +2—+8 °C ;
- Container for disposal of used materials;
- Universal Turbilatex Sample Collection Vial
- Specimen collection container;
- «Eppendorf» Plastic micro-vials, capacity up to 2,0 ml (if necessary).

All medical kits used for joint usage should have valid registration certificates for legal circulation on the Russian Federation territory. There are no restrictions on the joint FOB Turbilatex usage with the above-mentioned products.

Others of Universal Turbilatex Sample Collection Vial besides recommended by FOB Turbilatex manufacturer is not allowed.

Analysis

Frequency of vial use

The product is intended for single use.

Components preparing for use

All components of FOB Turbilatex kit (Reagent 1, Reagent 2, calibrators and control materials, vials) are ready for use, no preliminary preparation is required.

Samples preparation for use

Use the special Sample Collection Vial from FOB Turbilatex list to collect the sample. Perform sample preparation following the procedure below.

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
Filter

Cap 1: Stool sample
collection zone.

Universal Turbilatex Sample
Collection Vial

1. Put the sample collection vial with
cap 1 up and take out it.

Cap 1 with the stick to take the sample.

Stool sample

2. Introduce the stick into 4 different zones
of the sample

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
3. Close cap 1.

4. Shake the sample collection vial with diluent+sample in order to get good sample dispersion. Use a vortex (1 minute).


5. Open cap 2 and add enough volume to follow with the analytical procedure to assay the sample, and then close cap 2.

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Collect enough quantity of feces (1-2 g or mL for liquid sample) in clean and dry container (no preservatives or transport media). Homogenize stool samples as thoroughly as possible prior to preparation.

Take out the cap 1 of Universal Turbilatex Sample Collection Vial and use the stick to pick up enough sample quantity for procedure. Introduce the stick only one time into 4 different zones of the stool sample and add it to the sample dilution vial. If the sample is liquid (type 7, Bristol stool form scale) add 20 μ L of sample into the sample collection vial by a micropipette.

Close the vial (cap 1) with the diluent and stool sample. Shake the vial with the cap 1 down to assure good sample dispersion, using a vortex (1 minute). Let diluted sample to sit for some minutes. The sample collection vial with diluted sample can be stored for 7 days in the temperature 2-8°C.

After given time expires, open the cap 2 and add enough volume to follow with analytical procedure to assay the sample, and then close cap 2. The

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required volume of a biological sample in a secondary test tube consists of the dead volume of the sample in an «Eppendorf» Plastic micro-vials, the target volume of the sample directly dosed into the reaction cuvette and the extra volume of the sample. Thus, the volume of the sample sufficient to perform the analysis is determined by the analyzer used and the type of micro-sample used.

Universal Turbilatex Sample Collection Vial use for stool samples (option II):

← Cap 1: Stool sample collection zone.

1. Universal Turbilatex Sample Collection Vial.

2. Pull out the shaped stick by turning the screw-cap (white) counter clockwise.

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
3. Introduce the stick into 4 different zones of the sample

4. Shake the sample dilution vial with diluent+sample in order to get good sample dispersion. Use a vortex (1 minute). Centrifuge.

5. Remove both the cap and the blue separator by turning them to the opposite side.

6. Load vials in the analyser. Close with the additional cap

Collect enough quantity of faeces (1-2 g or mL for liquid sample) in clean and dry container (no preservatives or transport media). Homogenize stool samples as thoroughly as possible prior to preparation.


Pull out the shaped stick by turning the white screw-cap counter-clockwise. Do not remove the blue plastic. Replace the stick into 4 different zones of stool

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
sample and repeatedly rotate the stick in to make the sample as homogeneous as possible. Make sure all grooves are filled with faecal material.

Before placing the stick into the tube, remove the stool excess by rotating the stick on the internal wall of the stool container.

Fully insert the stick with faecal material into the tube containing the extraction solution through the separator (blue). While turning the screw cap clockwise until completely closed, keep the separator (blue) by the fingers to avoid its opening. The Universal Turbilatex Sample Collection Vial with diluted sample can be stored for 7 days in the refrigerator (2- 8°C).

Samples from fridge must be kept at room temperature before using.

Shake each vial by vortex for 30-60 seconds to properly homogenate the content. Make sure if the grooves of the devices are visually empty of any stool material, otherwise repeat the vortexing for other 30-60 seconds. Proceed to next step regardless if there is some residual stool sample in the grooves. Vortex time is 120 seconds maximumly.

Centrifuge the vial for 10 minutes at 1000–3000 g or allow it to stand upright for at least 15 minutes to remove the possible residuals of faecal material, which might damage the aspirating and dispensing systems of the automated clinical chemistry instruments. After centrifugation, remove the vial from the centrifuge.

After the completion of the previous steps, remove the white screw cap and the blue separator by rotating the separator clockwise.

The faecal extract (supernatant) is now ready to be tested. Insert the opened vial directly into the the sample rotor of the clinical chemistry analyser or transfer the supernatant to «Eppendorf» Plastic micro-vial..

Close the collection device by the provided additional screw cap and store or discard.

Sample performing for haemoglobin level measuring with the FOB Turbilatex kit can be only using special Turbilatex test tubes manufactured by CERTEST BIOTEC S.L., Spain.

Vials produced by other manufacturers are not suitable to use with the FOB Turbilatex kit.

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Analytical method description and test principle

Antibodies against human haemoglobin covalently bound to latex particles in Reagent 2 and interact with sample haemoglobin. As a result, antigen-antibody immune complexes and cloudy solution are formed. If no antigen (haemoglobin) presents in the sample the mixture will keep its appearance as a smooth suspension. Suspension turbidity intensity is directly proportional to haemoglobin concentration in the sample. Such intensity of turbidity is measured by analyzer.

Sample haemoglobin concentration is determined by a calibration curve done by haemoglobin calibrators of various levels. FOB Turbilatex Calibrator 0 is used till the zero point.

Quality control is implemented control materials - FOB Turbilatex Control 1 and FOB Turbilatex Control 2, assigned for the analyte being determined.

Measurement Procedure

Determination method: particle enhancement immunoturbidimetric test.

Haemoglobin concentration determination is a photometric measurement of the antigen-antibody reaction between human haemoglobin antibodies bound to latex particles and haemoglobin in sample.

Basic settings for the biochemical automatic analyzer CA -800 by Furuno Electric Co., Ltd.:

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Wavelength	510 nm
Temperature	37 °C
Measurement	End point
Sample/calibrator volume	6 mkL
Reagent 1 volume	60 mkL
Reagent 2 volume	15 mkL
Absorbance of blank	35–37 Cycle
Incubation	
Absorbance	64–66 Cycle

Result is automatically calculated by analyzer in ng/mL.

Note: During FOB Turbilatex using on other biochemical analyzers, request adaptations from an authorized representative of the manufacturer in the Russian Federation.

Calibration

For photometric systems calibration use only the following materials:
FOB Turbilatex Calibrator 0–5 (options I and III).

For calibration the assigned haemoglobin values for each level are used that indicated in FOB Turbilatex instructions and values sheet and on the label of each vial.

Haemoglobin concentrations in samples are obtained from a calibration curve done with six calibrators of different levels, including a zero-concentration calibrator.

Calibration stability is 1 month. Calibrate the system at least once a month is extremely recommended by manufacturer.

Recalibrate the photometric systems in case of below:

- when the controls are out of the assigned range for the specified lot of control material used;
- new lot of reagent is used;
- other cases provided by lab internal rules.

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The calibrator is traced to cyanomethemoglobin, that is used as an international standard for immunodiagnostic assays through an internal standard derived from this international standard and stabilized accordingly.

Quality Control

Quality control is to verify the analyte determination results in each analytical series, done directly in the lab and allowing to identify random and systematic errors of analytical stage. It is recommended to assay FOB Turbilatex Control 1 and FOB Turbilatex Control 2 measurements in each samples batch.

Haemoglobin concentration in controls is in the range of low and high values. The assigned haemoglobin values indicated on vial label and in relevant values sheet refer only to the specified lot of control materials.

Results of control materials analysis are intended to use for definition of the acceptability of analytical series results, which can be used for the diagnosis and disease or treatment planning prediction after. The results reliability issue for most tests can be solved via the regular use of control materials and appropriate statistical methods. Each lab should have corrective actions in case of control range deviations.

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Haemoglobin concentration in samples is calculated by a calibration curve using an appropriate mathematical model that is recommended in adaptations to a specific biochemical analyzer.

Limitations

1. FOB Turbilatex must be used in human stool samples only.
2. FOB Turbilatex is intended for professional *in vitro* diagnostic use only.
3. This test may be less sensitive for detecting upper gastrointestinal bleeding.
4. Technical errors and violations during the analysis procedure, as well as the possible presence of various substances in the stool samples, can affect the analysis (see the corresponding document section), this may lead to erroneous results.
5. Do not mix reagents, calibrators and control materials of CERTEST BIOTEC S.L., Spain manufactured that have different lots, and analogue components of another supplier.
6. Do not use sample collection vials that are not recommended by manufacturer.
7. During a medical device using follow the instructions strictly. Failure to follow the manufacturer's recommendations and product handling may lead to erroneous results and its misinterpretation.

Analytical characteristics

Results obtained in definite labs can be differ from below analytical characteristics.

Measuring range (Linearity)

Haemoglobin measuring range is 10-1000 ng hHb/mL.

If haemoglobin concentration in sample exceeds upper values of measuring range, diluent and extracted sample should be diluted by new diluent from

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another Sample collection vial manufactured by CERTEST BIOTEC S.L., Spain, in the ratio of 1: 3 and multiply the result by dilution.

Prozone effect

Prozone effect has not been observed to a concentration of human haemoglobin till 10 µg of hHb/mL. Studies using higher concentrations have not been carried out.

Sensitivity (Limit of detection)

Limit of detection (LOD): 8 ng hHb/mL.

Limit of quantification (LOQ): 10 ng hHb/mL.

Clinical sensitivity and specificity

Haemoglobin in feces measuring results were evaluated by comparing of an immunoturbidimetric assay using FOB Turbilatex kit manufactured by CERTEST BIOTEC S.L., Spain, and an immunological method using available competitive kits for detection of haemoglobin (faecal occult blood) in human stool samples.

The results obtained:

	FOB Turbilatex / Competitor 1	FOB Turbilatex / Competitor 2	Total
True positive results	22	48	70
False positive results	0	6	6
False negative results	1	2	3
True negative results	33	44	77
Amount of samples researched	56	100	156
Sensitivity, %	95,7 (78,1–99,9)	96,0 (86,3–99,5)	95,9 (88,5–99,1)
Specificity, %	100 (89,4–100,0)	88,0 (75,7–95,5)	92,8 (84,9–97,3)
PPV, %	100 (84,6–100,0)	88,9 (77,4–95,8)	92,1 (83,6–97,0)
NPV, %	97,1 (84,7–99,9)	95,7 (85,2–99,5)	96,3 (89,4–99,2)
LR+	63,75 (4,062–1000)	8,00 (3,769–16,98)	13,27 (6,129–28,71)
LR-	0,063 (0,013–0,299)	0,045 (0,012–0,177)	0,044 (0,015–0,134)

Comments:

- Sensitivity – probability of test result to be positive in case of disease.
- Specificity – probability of test result to be negative in case of disease absence.

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- PPV – predictive value of positive result (probability of disease presence in case of positive test).
- NPV – predictive value of negative result (no disease probability in case of negative test).
- LR+ – Positive likelihood ratio – ratio of positive test result probability in case of disease presence.
- LR- – Negative likelihood ratio – ratio of negative test result probability in case of disease and negative test result probability in case of disease absence.

The results reflected high sensitivity and specificity of human hemoglobin determination using the FOB Turbilatex by CERTEST BIOTEC S.L. kit, Spain.

Measurement accuracy

Inter-series accuracy (reproducibility):

Haemoglobin concentration, ng/mL	BS200E			Biolis 24i			ChemWellT		
	Average value	CO, ng/mL	CV, %	Average value	CO, ng/mL	CV, %	Average value	CO, ng/mL	CV, %
50	50,3	7,7	15	50,8	5,9	12	50,1	6,2	12
250	245,5	8,9	4	251,2	8,4	3	275,4	20,6	7
500	505,2	18,2	4	15,0	15,0	3	521,1	40,8	8

In-series accuracy:

Haemoglobin concentration, ng/mL	BS200E			Biolis 24i			ChemWellT		
	Average value, ng/ml	CO, ng/ml	CV, %	Average value, ng/ml	CO, ng/ml	CV, %	Average value, ng/ml	CO, ng/ml	CV, %

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Interferences

Substances are time to time in human stool samples do not affect the measurement. No interferences were founded against other substances tested: BSA, , Ascorbic acid, Hemin (concentration equal or lower to 0.25 mg/mL).

Cross reactivity

No cross reactivity was founded against other faecal markers occasionally present in faeces, such as human calprotectin, bovine and pig haemoglobin, human lactoferrin, bovine and human transferrin.

The cross reactions with bovine, pig haemoglobins and bovine transferrin are not present at concentrations below 0.2 µg/mL, 0.9 µg/mL and 250 µg/mL respectively. Take note that these concentrations are higher that the concentrations that could be found on the human faecal samples.

Expected values (reference interval) and results interpretation

In some studies it has been found haemoglobin (faecal occult blood) in human stool samples level ≥ 5.1 mcg/g used as cut-off, allowing to identify adults with gastrointestinal bleeding, who requires extra diagnostic invasive procedures for diagnosis. This value is equivalent to 50ng hHb/ml diluent from vial of FOB Turbilatex.

Positive result detecting abnormal presence of human haemoglobin in stool samples	Haemoglobin in human stool samples level is higher or equal to threshold value by clinical lab.
Recommended levels of hemoglobin in feces	Diagnostics – 5 mkg /g; Screening – 20 mkg/g.

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Received results interpretation by FOB Turbilatex:

below 50 ng/ml	Normal values, no gastrointestinal bleeding.
50–200 ng/ml	Slight bleeding in the gastrointestinal tract. Monitoring and patient follow-up is recommended.
above 200 ng/ml	Severe bleeding in gastrointestinal tract.

FOB Turbilatex cut-off value: 50-200 ng hHb/mL (x100 dilution factor, 5-20 µg hHb/g in stool samples).

Positive results determine the presence of human Haemoglobin (hHb) in stool samples. Positive result should be accompanied by additional diagnostic invasive procedures to determine precise cause and blood in the stool source.

If symptoms or the situation persists, haemoglobin detection should be conducted by invasive methods. Negative results do not exclude bleeding, as some polyps and colorectal cancers may bleed intermittently, or they do not bleed during certain disease stages. Moreover, blood may be unevenly distributed in stool samples.

Pre-packing and Packing

Reagent 1 and Reagent 2 are packaged in plastic vials (Reagent 1: 68 or 20 mL; Reagent 2: 17 or 5 mL), closed by plastic screw caps. FOB Turbilatex Calibrator 0-5 are packaged in six plastic vials of 0.5 mL each, closed by plastic screw caps. FOB Turbilatex Control 1 and FOB Turbilatex Control 2 are packaged in plastic vials (2 or 9 vials for each level) of 0.5 mL each, closed

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by plastic screw caps. The test tubes are packed in bags, if necessary, the bags are packed in a cardboard box (optional).

Tests number per kit is determined by reagents consumption per measurement that is directly related to biochemical analyzer used and its technical features. Bottles configuration is compatible with the reagent rotor of some biochemical analyzers series, including BioMajesty, Furuno CA, Dirui CS, URIT, OLYMPUS AU, Beckman Coulter AU, Siemens Advia, Biolyzer 300/600, Sapphire 800, this allows them to be put into the rotor without transfusion into other vials or through special adapters, provided with the analyzer. Reagent bottles have a barcode (optional) containing reagent information, lot number and expiry date. The barcode is read by device (analyzer) automatically.

Each production item is completed with instructions for use done by the manufacturer for the Russian market in accordance with valid standards. It contains all the necessary information for the correct measurement and results interpretation.

FOB Turbilatex (options I, III) is completed with values sheet, that includes the concentration of analyte in calibrators, range of control materials for quality control during FOB Turbilatex using.

FOB Turbilatex components are packed in a cardboard box to protect from mechanical and climatic factors during transportation/storage and it guarantees the convenience during loading and unloading operations.

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Storage and Stability

FOB Turbilatex should be stored in the transport or individual packaging of manufacturer in unopened vials indoors, in a refrigerator (industrial, laboratory), in a place protected from light. Permitted temperature of FOB Turbilatex storage and Universal Turbilatex Sample Collection Vial is below. Freezing of product and its components of Kit and accessories is not allowed. Contaminated FOB Turbilatex components are not for storage and usage. In case of storage conditions product can be stored and used during the expiration date indicated on the label. Products stored in violation of specified storage conditions are not subject to use. Stop to use it after expiry indicated on the product label.

Storage temperature: Reagent 1, Reagent 2, calibrators, controls – from +2 to +8 °C; vials for collecting stool samples – from +2 to +30 °C.

To get reliable results, strict compliance with instructions is necessary.

General precautions

1. For *in vitro* diagnostic use only.
2. Do not use the product if pack is damaged.
3. Do not use open vials or damaged pack.
4. FOB Turbilatex contains biological material. Follow universal safety regulations and Good Lab Practice during product applying as potential infectious specimens.
5. All the specimen should be considered potentially hazardous. Should be handled in the same manner as an infectious agent.
6. Sample contains preservatives such as sodium azide (<0.1%). Don't swallow! Avoid any contact with skin or mucous membrane.
7. Follow Good Laboratory Practices, wear protective clothing, use disposal gloves, goggles and mask to avoid such components as vials, test

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tubes, caps, probes and test tubes and vials contents. Do not eat, drink or smoke in the working area.

8. Store and transport the product only in accordance with relevant document section conditions. Stop using the product if these terms are violated.

9. Stop using product after expiry indicated on the label.

10. The diagnosis of analysis results can only be made by a doctor and only by complete information, including other patient tests and clinical condition results. For diagnostic purposes, the results should be always associated with the patient's medical background, anamnesis, laboratory results and other available data.

11. For disposal of unused reagents, reagents with expired shelf life, as well as used reagents, including test tubes containing samples of the analyzed material, reaction mixtures, vials with reagents, tools and objects contaminated with biological material (store and transport) in accordance with the requirements of SanPiN 2.1.3684 -21.

12. Biological material, including tools and objects contaminated with the material, belong to the hazard class of medical waste B.

Information about medicinal product presence in a medical device

The product is not a medicinal product and it does not contain substances used in the pharmaceutical industry in its compound.

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Environmental protection requirements and medical device disposal procedure

1. Product components do not contain toxic and explosive substances to pose environmental threat and medical personnel health.
2. After test completing remove the used vials of reagents, calibrators, controls and samples from the analyzer.
3. FOB Turbilatex components contain biological material. Product Disposal as well as tools and objects contaminated by biological material, should be in protective clothing (gloves, medical gown).
4. Disposal and destruction of the FOB Turbilatex Kit (used, unused, unusable, including with damaged packaging, due to expiration date and/or non-compliance with the regulated regime of storage, transportation and use) should be carried out in accordance with the requirements of SanPiN 2.1.3684-21 for the disposal of medical waste of class "B".
5. Biological material including tools and objects contaminated by material, belong to medical waste B hazard class.
6. To dispose of the kit components, use suitable containers to avoid environmental pollution.

Installation, adjustment, calibration and other actions necessary for putting a medical device into operation

Not applicable.

Requirement for premises to install (assembling) medical device is supposed to be

Not applicable.

Requirement for persons to be qualified and trained to install (assembling) medical device is supposed to be

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Requirements for checking of medical device correct installation and its readiness for safe operation

Not applicable. Product is not subject to installation and maintenance.

Information about sterile medical device condition

FOB Turbilatex is a non-sterile medical device. FOB Turbilatex usage is in non-sterile *in vitro* conditions in analyzer. All precautions are aimed to prevent product components from human body and environment entering, as they include biological material that is considered potentially dangerous. Precautions for product handling are described in the relevant section of this document.

Components sterilization and disinfection - not applicable.

Information of human and animal origin materials

CERTEST BIOTEC S.L. duly registered in accordance with Spanish law, specializing in the research, development and production of "in vitro" products for clinical diagnostics. The company confirms that all reagents and/or components of the FOB Turbilatex do not contain animal or human cells, tissues and/or their derivatives.

Manufacturer informs compound of this kit contains components of biological origin. Protein (human hemoglobin) in calibrators and controls was obtained from healthy, properly tested donors and does not bring any danger. Animal origin materials (antibodies) used in production were accompanied by an official veterinary certificate of animal health state.

Animal origin components are from healthy animals under veterinarian supervision, they are not susceptible to infectious, contagious and parasitic diseases. Nevertheless, product should be handled the same way as a potentially dangerous vector of infectious diseases, using personal protection, as well as thoroughly washing hands after product usage.

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Information about nature, type, intensity and distribution of medical device radiation

Not applicable.

Medical device Marking

Symbols list of FOB Turbilatex on the label below.

Symbol	Symbol description	Symbol	Symbol description
	Medical product for <i>in vitro</i> diagnostics		Caution! Refer to the instructions
	Consult instructions for use		Use by
	Keep away from sunlight		Catalogue number
	Keep dry		Lot number
	Temperature limitation		Compound
	Manufacturer		Barcode
	Biological risk		QR-code

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	Logo		CE mark
	Quantity in the package		Prohibition of re-use

Label also includes "For professional use".

Transportation

Product transportation should be realized by all types of covered transport in accordance with the current rules of transportation. Transportation can be in the temperature condition: FOB Turbilatex : from +2 to +8 °C; Universal Turbilatex Sample Collection Vial: from +2 to +30 °C.

Transportation is possible only in compliance with storage conditions. Product and accessories freezing is not allowed. Products transported in temperature violation are not subject to use.

Products in transport packaging are resistant to the mechanical factors effects that occur during transportation.

Guarantee

Manufacturer confirms it has not had any indication to date that this product caused or contributed to adverse events, undesirable phenomena or corrective actions related to safety.

Product has not warranty service.

Manufacturer guarantees the quality, effectiveness and safety, unacceptable risk absence of harm to life, human health and the environment during FOB

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Turbilatex using for its intended purpose under the conditions provided by instruction.

Manufacturer's warranty obligations apply from the date of manufacture and for medical device expiry, subject to proper transportation, storage and operation.

If it is necessary to get advisory support or in case of products quality questions, contact the authorized representative of the manufacturer in Russian Federation territory.

List of the Russian Federation national standards for products confirming the quality of medical device

Standard	Name
ISO P 51352-2013	In vitro diagnostic medical devices. Test methods
ISO P 51088-2013	Medical devices for <i>in vitro</i> diagnostics. Reagents, reagent kits, test systems, control materials, nutrient media. Requirements for products and supporting documentation
EN ISO 23640-2015	<i>In vitro</i> diagnostic medical devices - Evaluation of stability of <i>in vitro</i> diagnostic reagents
ISO 18113-1-2015	<i>In vitro</i> diagnostic medical devices - Information supplied by the manufacturer (labelling) - Part 1: Terms, definitions and general requirements
ISO 18113-2-2015	<i>In vitro</i> diagnostic medical devices — Information supplied by the manufacturer (labelling) - Part 2: <i>In vitro</i> diagnostic reagents for professional use.
ISO 15223-1-2020	Medical devices — Symbols to be used with information to be supplied by the manufacturer — Part 1: General requirements
ISO 14971-2021	Medical devices. Application of risk management to medical devices.

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ISO 17511-2011	In vitro diagnostic medical devices. Measurement of quantities in biological samples. Metrological traceability of values assigned to calibrators and control materials.
SanPiN 2.1.3684-2021	Sanitary and epidemiological requirements for the maintenance of urban and rural settlements territories, for water objects, drinking water and drinking water supply to the population, atmospheric air, soils, residential premises, operation of industrial public premises, organization and conduct of sanitary and anti-epidemic (preventive) measures

References

1. Towler BP, Irwig L, Glasziou P, Weller D, Kewenter J. Screening for colorectal cancer using the faecal occult blood test, Hemoccult. Cochrane Database Syst Rev. 2000;(2)CD001216.
2. Ransohoff DF and Lang CA. Screening for colorectal cancer with the Fecal Occult Blood Test: a background paper. Ann Intern Med. 1997; 126: 811-822.
3. Ransohoff DF and Lang CA. Suggested technique for Faecal Occult Blood testing and interpretation in colorectal cancer screening. Ann Intern Med. 1997; 126: 808-810.

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