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Document No. ST-K-0176

Date of Issue: MAY-22-2023

*«I confirm the accuracy, correctness and
reliability of the content of the texts»*

President

Sentech Korea Corp.

SENTECH KOREA CORP.

D. J. Yoo

D. J. YOO / PRESIDENT

Do Joon Yoo

« 12 » May 2023

Design file taking into account the operational suitability
for a medical device:

«Ethanol Vapor Concentration Analyzer in Exhaled Air Dingo with Accessories»

Productions

Sentech Korea Corp., Korea

2023

Table of contents

1. Characteristics and purpose of a medical device	3
2. Technical description of the medical device	6
3. Identification of the characteristics of a medical device that affect the safety of its use	14
4. Operational Suitability Specification	19
5. Risk assessment	22
6. Residual risk	27
7. Conclusions	27
8. Conclusion	28

Introduction

This document is designed for a medical device "Ethanol vapor concentration analyzer in exhaled air Dingo with accessories" (hereinafter referred to as the analyzer) to achieve sufficient operational suitability of the medical device(s).

The purpose of the risk analysis:

- a) to determine the risks associated with the operational suitability of the product;
- b) assess and limit these risks;
- c) check the measures taken to limit risks.

This file is made according to the requirements of GOST R IEC 62366-1-2021 "Medical devices. Part 1. Design of medical devices taking into account operational suitability".

In order to achieve a sufficient (safe) level of operational suitability of the analyzer for the patient/user, an action plan was developed according to the methodology given in Annexes E, G and H to the standard, so as to have the most complete list of possible risks available.

This document includes the specification of operational suitability, the main operational functions, use cases, user interface requirements for the main operational functions, risk analysis, risk assessment and the acceptability of residual risks.

During the risk analysis, the following research methods were used: the study of literature, documents and results of activities in the field under study; monitoring, study and generalization of experience;

systematic and comparative analysis; formation of expert assessments.

1. Characteristics and purpose of a medical device

1.1 General description, including intended use of the medical device

The Dingo Exhaled Air Ethanol Vapor Concentration Analyzer with Accessories is designed for express measurement of the mass concentration of ethanol vapors in a selected sample of exhaled air.

Class depending on the degree of potential risk of use for medical purposes in accordance with the nomenclature classifier of medical devices approved by the Order of the Ministry of Health of the Russian Federation dated June 06, 2012 No. 4n: 1.

Scope of application: The analyzer is used in professional practice by medical workers in the implementation of pre-trip and post-trip inspection of drivers, pre-shift and post-shift inspection of employees of organizations, as well as ambulance workers. The analyzer can be used for medical examination for intoxication.


Recommendations on areas of application in the field of state regulation of ensuring the uniformity of measurements:

- implementation of activities in the field of health;
- implementation of measures of state control and supervision in the field of road safety (if there is a portable printer for printing measurement results).

Method of application: according to the instruction manual.

1.2 General description of the main functional elements of the main type representatives of the medical device

A) Structurally, the analyzer is made in the form of a monoblock. The appearance of the analyzer and the designations of the elements are shown in picture 1.

Pic. 1 The appearance of the analyzer

1. LED	6. On/off/confirmation button	11. Batteries
2. IR port	7. Battery compartment button	12. USB port
3. Display	8. Battery cover with marking sticker	13. Strap hole
4. Up button	9. Sensor unit	14. The hole of the mouthpiece retainers
5. Down/menu button	10. Sensor block latches	15. Mouthpiece

On the front side of the breathalyzer there are on/off buttons, menu navigation and a backlit graphic display.

The battery compartment cover is located on the back of the device. Under it there are batteries and a removable sensor unit.

In the upper end of the analyzer there is a recess for the mouthpiece, which can be located both left-hand and right-hand orientation for the convenience of the operator.

The breathalyzer is used for rapid non-invasive testing, exhaled air is used as the analyzed medium. The device, in a state of readiness for testing, is brought to the mouth of the subject, from whom it is

required to exhale into a disposable mouthpiece. The pressure from the exhaled air activates the sampling device, after the analysis of which the result is displayed digitally. Based on the result of the analysis, the operator decides on further actions in relation to the subject.

The analyzer is simple and easy to operate, its work is based on modern achievements of microelectronics. The analyzer is controlled using three buttons on the front panel and a menu system.

The analyzer power button is located on the front panel. When the analyzer is turned on, voltage is applied to the analyzer circuit, and the autotesting procedure begins.

The graphical display shows the measurement results, messages about the operating modes of the analyzer, the state of charge of the batteries and auxiliary information.

The analyzer uses an electrochemical sensor to measure the mass concentration of ethanol vapors in the analyzed air, manufactured by Sentech Korea Corp., Korea.

The analyzer's microprocessor controls all operating modes and converts the output signals of the measuring sensor into readings. The microprocessor fully controls all stages of the measurement, and any violation of the procedure is displayed on the display as a corresponding warning. The built-in software of the analyzers was developed by the manufacturer specifically to solve the problem of measuring the mass concentration of ethanol vapors in exhaled air, as well as to display the measurement results on the display. Identification of the embedded software is performed by displaying the version on the analyzer display in the administrator menu in the "Information" section in the "Device info" subsection.

All stages of the analyzer are accompanied by sound signals.

When performing measurements, disposable plastic mouthpieces are used. The mouthpiece is delivered in an individual package, which is opened immediately before taking measurements.

The measurement results are printed out on a portable printer. For communication between the analyzer and the printer, a communication cable from the supplied package is used. For analyzers in the Dingo E-200(V) version, it is possible to connect to the printer via the Bluetooth interface.


2. Technical description of the medical device

2.1 Main technical characteristics

2.1.1 The analyzers comply with the requirements of GOST R 50444-2020, GOST R IEC 60601-1-2010, GOST R IEC 62304-2013, GOST R ISO/IEC 9126-93, GOST 28195-89, GOST R ISO/IEC 12119-2000, GOST R 51188-98, GOST R IEC 60601-1-6-2014 and a set of design documentation.

Name	Photo	Description
Analyzer case		Overall dimensions (L×W×C), mm not more than: 150×85×10 Weight, g not more than: 20
AA battery		Alkaline (alkaline) batteries or accumulators, rated voltage — 1.5 V Weight, g not more than: 24
Disposable mouthpiece		The mouthpieces are made of food-grade polystyrene, when tested, they are attached directly to the body of the breathalyzer Overall dimensions (L×W×C), mm not more than: 101×15×17 Weight, g not more than: 5

Case for transportation and storage		The case is made of durable plastic, a separate cell is provided inside for each element, which prevents their collision and, accordingly, accidental damage. Overall dimensions (L×W×C), mm not more than: 230×55×195 Weight, g not more than: 340
A case for transporting and storing the analyzer with a printer		The case is made of durable plastic, a separate cell is provided inside for each element, which prevents their collision and, accordingly, accidental damage. Overall dimensions (L×W×C), mm not more than: 330×120×300 Weight, g not more than: 1460
Arm strap		The strap is used to prevent the device from falling in order to avoid unpleasant consequences Strap length, mm max: 240 Weight, g not more than: 2
Funnel mouthpiece		Special mouthpieces for repeated use, designed to conduct an exhalation test to obtain an estimated result, are made of polymer material, attached directly to the body of the breathalyzer Overall dimensions (L×Ø×C), mm not more than: 44×45×33 Weight, g not more than: 5

Printer		The portable printer is used to print out the test results, powered by 4 AA alkaline batteries, or from an external power source via an adapter with a stabilized output voltage of 7.5 V / 2.0A Overall dimensions (L×W×C), mm not more than: 131×76×53 Weight, g not more than: 255
Cable for connecting the analyzer to a PC		The cable is used for data exchange between the analyzer and the PC through a special communication program. (item 7) Cable length, mm max: 1530 Weight, g not more than: 61 Input connector interface: USB Output connector interface: HDMI
Cable for connecting the analyzer to the printer		The cable is used to transmit the test results to print. Cable length, mm max: 1520 Weight, g not more than: 60 Input connector interface: 6P6C Output connector interface: HDMI
Adapter to the printer from 220 V		The adapter is used to connect the printer to a household AC power supply. Cable length, mm max: 1890 Weight, g no more: 121 Alternating current: 110-220V, 50/60 Hz, 0.5A Constant current: 7.5 V, 2.0A

AA battery for the printer		Alkaline (alkaline) batteries or accumulators are used for autonomous power supply of the printer; rated voltage — 1.5 V Weight, g not more than: 24
Paper to the printer		Standard thermal paper with a width of 57 ± 0.1 mm, is used to print out the test results

Name of the parameter	Specifications	
	Dingo E-200	Dingo E-200(B)
General information:		
Weight, g, no more	200	
Overall dimensions (H×W×D), mm no more	133×65×33	
Sensor Type	electrochemical removable	
Purge mode	1. Measuring mode - through a disposable mouthpiece 2. Screening mode - through a funnel mouthpiece, or without the use of mouthpieces	
Calibration interval	1 year	
Displaying results	4 – bit digital display	
Range of indications, mg/l	От 0,000 до 2,000	
The price of the lower division of the scale	0,005	

discharge, mg/l		
Parameters of the analyzed gas mixture when the sample is fed to the input of the analyzers (automatic sampling mode):		
- consumption of the analyzed gas mixture, l/min, not less	9	
- the volume of the sample of the analyzed gas mixture, l, not less	1,2	
Time of preparation for work after switching on, s., no more	20	
Measurement time after sampling, sec., no more	10	
Preparation time for operation after analysis of a gas mixture with a mass concentration of ethanol 0.50 mg/l, sec., no more	120	
The time interval of the analyzers without adjusting the readings, months, no more	12	
The number of measurements on the analyzers without replacing the batteries, at least	3000	
Communication with the printer	1. Via USB port via a wired cable	1. Via USB port via a wired cable; 2. Via Bluetooth protocol
Power supply:		
Analyzer power supply	AA type battery	
Number of batteries, pcs.	2	
Operating voltage, V	3	

Printer power supply	AA type battery
Number of batteries, pcs.	4
Operating voltage, V	6
Operating conditions:	
Operating temperature range	From +10°C to +40°C
Relative humidity	From 10% to 90%
Atmospheric pressure	From 84,0 kPa to 106,7 kPa
Transportation and storage conditions:	
Ambient temperature	From -10°C to + 50°C
Relative humidity	From 10% to 90% (without condensation)
Atmospheric pressure	60 kPa – 140 kPa
Special requirements for transportation and storage conditions	Keep the appliance away from alcohol-containing liquids

Measurement range of ethanol mass concentration, mg/l	Limits of permissible error	
	Absolute	Relative
From 0,000 to 0,500 include.	±0,050 mg/l	-
More than 0,500 to 1,500	-	±10%

Notes:

- 1) In the analyzers, the minimum interval of readings is programmatically set, which are displayed on the analyzer display and paper media in the form of zero readings: from 0.000 to 0.030 mg /l.
- 2) When checking the readings of analyzers using ethanol/ nitrogen gas mixtures in pressure cylinders, a conversion factor of the readings K^{GS} equal to 1.1 is used
- 3) There is no additional error due to the presence of unmeasurable components.

2.2 The scope of delivery of the medical device

I. Ethanol Vapor Concentration Analyzer in exhaled Dingo air with accessories consisting of:

1. Dingo E-200/ Dingo E-200 breathalyzer(V) – 1 pc.
2. AA battery - 2 pcs.
3. Disposable mouthpiece – 5 pcs.
4. Case for transportation and storage – 1 pc.
5. Hand strap – 1 pc.
6. Passport – 1 pc.
7. Instruction manual – 1 pc.

Accessories:

1. Analyzer case – 1 pc. (if necessary)
2. Funnel mouthpiece - 1 pc. (if necessary)
3. Printer – 1 pc. (if necessary)
4. A case for transporting and storing the analyzer with a printer – 1 pc. (if necessary)
5. Cable for connecting the analyzer to a PC – 1 pc. (if necessary)
6. Cable for connecting the analyzer to the printer – 1 pc. (if necessary)
7. Adapter to the printer from 220V – 1 pc. (if necessary)
8. AA battery to the printer – 4 pcs. (if necessary)
9. Paper to the printer – 1 roll (if necessary)

2.3 Marking

2.3.1 The marking of the medical device is made in accordance with the requirements of GOST R 50444, GOST R IEC 60601-1, GOST R ISO 15223-1, in Russian, in printed form and is resistant to aggressive environments in acceptable operating and maintenance conditions of the medical device.

2.3.2 On the back cover of the analyzer there is an individual marking of the product, which contains

the following information:

- name of the product model;
- trademark;
- serial number;
- date of manufacture in the format month/year;
- Refer to the operating instructions;
- the sign of the measuring instrument;
- rated voltage;
- degree of protection according to IP classification;
- the name of the manufacturer and the country of manufacture.

There is an individual marking on the back of the printer, which contains the following information:

- name of the printer model;
- serial number;
- DC value;
- The name and address of the manufacturer and the country of manufacture.

2.3.3 Transport marking according to GOST 14192. A marking is applied to the transport container (packaging), which contains the following information:

- product name;
- trademark;
- "special disposal" sign;
- production date in month/year format;
- overall dimensions of the package;
- gross weight;
- registration certificate number;
- barcode;
- storage conditions;
- transportation conditions;
- name of the manufacturer;
- address of the manufacturer;
- authorized representative of the manufacturer;
- manipulative signs.

The marking is applied to a paper label.

The handling signs corresponding to the values are applied to the transport package: "Fragile. Caution", "Top", "Keep away from moisture".

The designation of storage conditions and other additional inscriptions are applied to the container or

label in places free from transport markings.

2.4 Packaging

The operational documentation is packed in a plastic bag.

Disposable mouthpiece, funnel mouthpiece are packed in a HDPE plastic bag with a density of 925 kg/m³, 0.1 mm thick, manufactured by Shenzhen Xing Sheng Li Imports and Exports Co., Ltd. (China).

The analyzer with accessories is packed in a case and transport packaging – a corrugated box made of five-layer corrugated cardboard with a thickness of 6 mm, a density of 150 g/m³, white color, overall dimensions (L × W×C) not more than 400×250×305 mm, gross weight not more than 2.5 kg, manufactured by Shenzhen Xing Sheng Li Imports and Exports Co., Ltd." (China).

3. Identification of the characteristics of a medical device that affect the safety of its use

To identify the characteristics of the analyzer related to operational suitability, which may have an impact on safety, a list of questions has been compiled in accordance with GOST R IEC 62366 Appendix E. The answers to the compiled questions helped to conduct a comprehensive and most complete identification of the dangers associated with the use of a medical device.

E.2.1 Is the medical device sterile or is it intended for sterilization by the user, or are other types of microbiological treatment necessary?	
Factors considered:	
- is the medical device intended for single or multiple use;	<i>Designed for repeated use</i>
- any restrictions on the number of cycles of repeated applications;	<i>Does not require consideration</i>
- the method of the sterilization process used.	<i>Does not require consideration</i>
Is the medical device labeled appropriately to inform the user whether it is intended for single or multiple use?	<i>Does not require consideration</i>
Are any handling or service life restrictions clearly indicated on the packaging?	<i>Specified</i>
The operational documentation should clearly indicate the appropriate methods and means for cleaning or sterilization, and provide information on the frequency of cleaning.	Disinfection of the analyzer must be carried out after each use by wiping the outer surfaces twice with a cloth made of calico or gauze moistened with a 3% solution of hydrogen peroxide, in accordance with MU-287-113. The napkins should be wrung out. It is not recommended to use products containing aldehydes.

	Mouthpieces are disposable, are not subject to reuse and processing, they should be disposed of in accordance with SanPiN 2.1.3684-21. A disposable mouthpiece and a disposable mouthpiece with a check valve must be disposed of in accordance with Class B
E.2.2 Are measurements being performed?	<u>YES</u>
Factors considered:	
- measured variables;	<i>Requires consideration</i>
- reliability;	<i>Requires consideration</i>
- accuracy of measurement results;	<i>Requires consideration</i>
- awareness of the frequency with which measurements need to be carried out;	<i>Does not require consideration</i>
- parameter affecting the result;	<i>Does not require consideration</i>
- necessary consumables;	<i>Requires consideration</i>
- how to handle the results and interpret them.	<i>Requires consideration</i>
E.2.3 Is the medical device intended for use in conjunction with medicines or other medical technologies?	<u>NO</u>
Factors considered:	
- identification of any medications or other medical technologies that may be involved in the interaction;	<i>Does not require consideration</i>
- possible problems associated with such interaction;	<i>Does not require consideration</i>
- as well as the patient's reaction to the treatment.	<i>Does not require consideration</i>
E.2.4 Is there an unwanted release of energy or substances?	<u>NO</u>
Factors considered:	
- noise and vibration;	<i>Does not require consideration</i>
- heat;	<i>Does not require consideration</i>
- radiation (including ionizing, non-ionizing and ultraviolet/visible/infrared radiation);	<i>Does not require consideration</i>
- leakage currents;	<i>Requires consideration</i>
- electric and/or magnetic fields;	<i>Requires consideration</i>
- the negative impact of noise, vibration, heat and waste (for example, exhaust gases) on the physiology and psychology	<i>Does not require consideration</i>

of patients and third parties.	
E.2.5 Is the medical device sensitive to environmental influences?	<u>YES</u>
Factors considered:	
- shock, vibration;	<i>Requires consideration</i>
- pressure, temperature and humidity;	<i>Requires consideration</i>
- illumination;	<i>Does not require consideration</i>
- sensitivity to changes in power supply and cooling;	<i>Does not require consideration</i>
- electromagnetic interference.	<i>Does not require consideration</i>
E.2.6 Are there any special consumables or accessories related to the medical device?	<u>YES</u>
Factors considered:	
- shock, vibration;	<i>Requires consideration</i>
- pressure, temperature and humidity;	<i>Does not require consideration</i>
E.2.7 Is maintenance and/or calibration necessary?	<u>YES</u>
Factors considered:	
- is preventive maintenance and/or calibration carried out by the user or the responsible organization;	<i>Requires consideration</i>
- whether special substances or equipment are necessary for proper maintenance and/or calibration.	<i>Requires consideration</i>
E.2.8 Does the medical device have a limited shelf life?	<u>YES</u>
Factors considered:	
- marking;	<i>Requires consideration</i>
- indicators;	<i>Does not require consideration</i>
- methods of disposal of a medical device.	<i>Requires consideration</i>
E.2.9 Are there any delayed and/or long-term consequences of using?	<u>NO</u>
Factors considered:	
- ergonomic and cumulative effects;	<i>Does not require consideration</i>
- prolonged exposure to vibration, noise, heat, gases, as well as insufficient ergonomics (wear of joints, muscles and nerves, etc.).	<i>Does not require consideration</i>
E.2.10 What mechanical forces can affect a medical device?	-
Factors considered:	
- information about whether the responsible organization, the user, or controls the mechanical forces affecting the	<i>Does not require consideration</i>

medical device;	
- it is controlled by the user in interaction with other persons.	<i>Requires consideration</i>
E.2.11 Is the medical device intended for single use?	<u>NO</u>
Factors considered:	
- is the medical device self-destructible after use;	<i>Does not require consideration</i>
- is it obvious that a medical device has been used?	<i>Does not require consideration</i>
E.2.12 Does the installation or use of a medical device require special training or special skills?	<u>YES</u>
Factors considered:	
- novelty of the medical device;	<i>Does not require consideration</i>
- probable skills and training of the person installing the medical device.	<i>Requires consideration</i> Qualified personnel familiar with the analyzer operation process
E.2.13 What information on safe use will be provided to the user?	-
Factors considered:	
- will the information be provided to the user directly by the manufacturer or will third parties be involved, such as installation/installation specialists, caregivers, healthcare professionals or pharmacists, and whether training is implied;	<i>Requires consideration</i>
- commissioning and transfer to the user and the feasibility / possibility of installation by persons who do not have special skills.	<i>Requires consideration</i>
E.2.14 Can the design features of the user interface lead to operational/application errors?	-
Factors considered:	
- design features of the user interface that can lead to errors in operation/application.	<i>Does not require consideration</i>
E.2.15 Is a medical device used in an environment with a lot of distractions?	<u>YES</u>
E.12.16 Does the medical device have connecting parts or accessories?	<u>YES</u>
Factors considered:	
- the probability of incorrect connection;	<i>Requires consideration</i> <i>The product has a disposable mouthpiece, which must be attached.</i>

- differentiation;	<i>Does not require consideration</i>
- similarity with compounds in other products;	<i>Does not require consideration</i>
- the strength of the connection.	<i>Requires consideration</i>
E.12.17 Does the medical device have a control interface?	<u>YES</u>
Factors considered:	
- spatial distribution of controls, their coding, grouping, mutual arrangement;	<i>Requires consideration</i>
- types of responses;	<i>Requires consideration</i>
- gross mistakes;	<i>Does not require consideration</i>
- differentiation of management;	<i>Does not require consideration</i>
- image quality;	<i>Does not require consideration</i>
- discreteness or continuity of operation of controls;	<i>Requires consideration</i>
- reversibility of settings or actions.	<i>Requires consideration</i>
E.12.18 How information is displayed on a medical device?	-
Factors considered:	
- conditions, orientation, visual capabilities of the user;	<i>Does not require consideration</i>
- populations, clarity of the information provided;	<i>Requires consideration</i>
- blocks, color coding and availability of critical information.	<i>Requires consideration</i>
E.12.19 Is the medical device controlled using the menu?	<u>YES</u>
Factors considered:	
- difficulty and number of menu levels;	<i>Requires consideration</i>
- sufficiency of information about the status of the menu;	<i>Does not require consideration</i>
- location of settings;	<i>Requires consideration</i>
- navigation method, the number of steps in one action;	<i>Does not require consideration</i>
- - clarity consistency and memory problems;	<i>Does not require consideration</i>
- the ratio of importance and availability of management functions.	<i>Requires consideration</i>
E.12.20 Will the medical device be used by persons with special needs?	<u>YES</u>
Factors considered:	
- mental and physical capabilities;	<i>Requires consideration</i>
- user skills and training;	<i>Requires consideration</i>
- ergonomics issues;	<i>Requires consideration</i>
- the environment in which the medical device will be used;	<i>Requires consideration</i>

- who will be installed and whether the patient can control or influence the use of a medical device.	<i>Requires consideration</i>
E.12.21 How can a medical device be intentionally used incorrectly?	-
Factors considered:	
- incorrect use of connectors	<i>Requires consideration</i>
- disabling security elements or alarm systems;	<i>Requires consideration</i>
- disregard of the MANUFACTURER's recommendations for maintenance.	<i>Requires consideration</i>
E.12.22 Is the medical device intended for use as a mobile or portable?	<u>YES</u>
Factors considered:	
- clamps, handles, wheels, brake;	<i>Does not require consideration</i>
- mechanical stability and durability.	<i>Requires consideration</i>

4. Operational Suitability Specification

The Serviceability specification is the documentation defining the user interface requirements regarding serviceability.

1. DESCRIPTION
The Dingo Exhaled Air Ethanol Vapor Concentration Analyzer with Accessories is designed for express measurement of the mass concentration of ethanol vapors in a selected sample of exhaled air. Analyzers can be used both for self-monitoring and for checking third parties. Analyzers are portable automatic devices of cyclic action. The work of the analyzers is fully automated, all stages of preparation and testing are accompanied by an audible alarm and text messages.
2. SUMMARY OF THE MEDICAL DEVICE APPLICATION SPECIFICATION
2.1. Medical appointment
The Dingo Exhaled Air Ethanol Vapor Concentration Analyzer with Accessories is designed for express measurement of the mass concentration of ethanol vapors in a selected sample of exhaled air.
2.2. Patient population
Age: from 12 years
Weight: no restrictions.

There are no contraindications to use.
Note: Diabetes, a low-calorie diet and other factors that lead to an increase in the concentration of ketones in the body may manifest themselves during testing in the form of small concentrations of ethanol that do not exceed the permissible threshold of 0.16 mg/l.
Nationality: any.
Patient's condition: - - condition is not important.
2.3. Body part or type of tissue to which it is applied or interacts with
Area of measurement/exposure: surface contact with human skin in the process of taking an air sample to analyze the concentration of ethanol vapors in the exhaled air.
Conditions: used as an analyzer of the concentration of ethanol vapors in the exhaled air of the patient.
2.4. Intended User
Education: - minimum: unlimited. - maximum: unlimited.
Knowledge: - minimum: reading and understanding Arabic numerals when writing them. - maximum: unlimited.
Understanding of the language: Russian.
Acceptable disadvantages: - moderate decrease in visual acuity or vision corrected; - - color blindness.
2.5. Application
General provisions: - - Use both indoors and outdoors.
Operating conditions of a medical device Operating temperature range: from +10°C to +40°C Relative humidity: from 10% to 90% Atmospheric pressure: From 84.0 kPa to 106.7 kPa
During operation, the analyzer must be serviced by qualified personnel – specialists of the buyer (or third parties hired by the buyer and agreed with the manufacturer). The frequency of maintenance of the analyzer is carried out at least once every 12 months.
Frequency of use: short-term regular use.
Mobility: mobile, portable product.
3. MAIN WORKING FUNCTIONS


Frequently used functions: - Screening mode; - Printing of measurement results; - Active test; - View memory; - Cleaning and disinfection.
SECURITY-related functions: - Protection from moisture and dust; - Electrical safety class II; - - cleaning and disinfection.
4. USAGE SCENARIOS
(See paragraph 1.2)
5. USER INTERFACE REQUIREMENTS FOR BASIC OPERATIONAL FUNCTIONS
Screening mode: - Easy button detection, menu navigation; - the location eliminates the random selection of an inappropriate function. - easy reading of control icons, test results.
Printing of measurement results: - availability of a portable/built-in printer; - easy to connect to a portable printer; - - fast printing immediately after the test for the built-in printer. - the ability to re-print.
Active test: - Easy button detection, menu navigation; - location eliminates the random selection of an inappropriate function; - speed of analysis; - easy reading of control icons, test results.
Memory view: - - availability of memory for a large number of measurement results; - easy reading of the results on the product screen; - sort by date and time.
Cleaning, disinfection and sterilization: - no surrounding protrusions or crevices for easy cleaning; - protection of the control unit, dust/moisture protection of the housing.
6. REQUIREMENTS FOR DETERMINING WHETHER THE MAIN OPERATIONAL FUNCTIONS ARE EASY FOR THE USER TO UNDERSTAND
Screening mode: - the color scheme, the location of the buttons, their designation do not affect the user's understanding of the functions; - easy attachment of a disposable mouthpiece.
Printing measurement results:

- the type of printer, the location of the buttons do not affect the user's understanding of the functions.
Active Test: - the color scheme, the location of the buttons, their designation do not affect the user's understanding of the functions; - easy attachment of a disposable mouthpiece.
Viewing memory: - easy sorting of measurement results by date and time; - large font, convenient placement of information on the screen.
Cleaning, disinfection and sterilization: - control unit, housing meet the requirements of moisture protection; - the case has no protrusions, cracks, hard-to-reach areas

5. Risk assessment

For a complete risk analysis of the operational suitability of a medical device, the entire range of application scenarios and related factors that can lead to harm were considered. For the evaluated medical devices, many hazards will never lead to harm and can be excluded from further consideration. The definition and relationship of the terms are shown in the figure and in Table F1. In accordance with GOST R IEC 62366 Annex F.

Term	Value from ISO 14971
Harm	Physical injury or damage to human health, or property, or the environment
Danger	Potential source of harm
DANGEROUS SITUATION	Circumstances in which people, property, or the environment are exposed to one or more hazards
Risk	The combination of the probability of causing harm and the severity of this harm
SEVERITY	Measure of possible consequences of danger

Note 1 is the probability of a dangerous situation.

Note 2 is the probability that a dangerous situation will lead to harm.

Figure F.1 - Graphical representation of the relationship between danger, sequence of events, dangerous situation and harm

Evaluation criteria:

To analyze the risk of operational suitability of a medical device, a qualitative method described in Appendix D of GOST ISO 14971 was used.

Determine 3 levels of probability (Table B.1) and 3 levels of severity (Table B.2). Fill in a 3×3 matrix (Table B.3) to describe the probability and severity of risk for each dangerous situation. Each element of the matrix represents a subgroup of the full set of possible risks.

The elements of the matrix are formed by the intersection of the range of possible probabilities and the range of possible severity (possible consequences).

Table D.1 - Levels of severity of harm by qualitative characteristics

Determination of severity	Description
Significant (3)	Death, or loss of body functions, or a change in the anatomical structure of the body
Moderate (2)	Reversible or minor injury (lesion)

Negligible (1)	No injury (lesion) or very minor injury (lesion)
----------------	--

Table D.2 - Levels of probability of harm by qualitative characteristics

Determination of probability	Description
High (3)	A very likely, frequent risk
Average (2)	Possible, but not frequent risk
Low (1)	Usually an impossible, rare, long-term risk

When using probability designations in the rows of the matrix, and severity designations in the columns, a 3x3 risk matrix can be created (Table D.3). The risks that are defined (R1, R2, R3 ...) are entered into the corresponding elements of the matrix.

Table D.3 Risk assessment matrix

Levels of severity of harm

Probability levels of harm	Negligible (1)	Moderate (2)	Significant (3)
High (3)			
Average (2)			
Low (1)			

Symbols:

- gray color is an unacceptable risk;
- white color - acceptable risk


thus, the risks falling into the area of the gray matrix will be unacceptable.

To assess the risk, Table 3 "Risk assessment and risk management measures introduced at the MI production stage" considers possible dangerous situations with an indication of constructive, protective and informational measures to prevent them. This table also provides an assessment of the probability of a dangerous situation and possible harm, as well as a conclusion about its admissibility.

Table 3. "Risk assessment and risk management measures introduced at the MI production stage"

№	Dang	Descriptio	Risk assessment	Harm	User	Residual risk
---	------	------------	-----------------	------	------	---------------

	er	n of the scenario*	Proba bility of harm	Seve rity of har m	Risk		Interface Risk Manage ment tools	Proba bility of harm	Seve rity of har m	Risk
1.	Analyzer	<u>The USER has not properly read the operating manual.</u> <u>The USER subjected the medical device to unskilled repair and maintenance.</u>	2	2	Invalid	Electric shock, user injury, product failure	1. Operation manual.	1	2	Acceptable
2.	Analyzer	<u>The USER has not properly read the operating manual.</u> <u>The USER has not familiarized himself with the assembly and operation requirements.</u>	2	2	Invalid	Electric shock, operator injury, product failure	1. Operation manual.	1	2	Acceptable
3.	Analyzer	<u>The USER has not properly read the operating manual.</u> <u>The USER</u>	2	1	Acceptable	Product failure	1. Operation manual.	2	1	Acceptable

		<u>did not ensure the use of the medical device in the specified electromagnetic environment.</u>								
4.	Analyzer	<u>The USER ignores the operating instructions of the medical device.</u> <u>The USER does not use the medical device for its intended purpose</u>	1	2	Acceptable		1. Operating room staff trained to work with a medical device	1	1	Acceptable
5.	Analyzer	<u>THE USER is not aware of the need for maintenance.</u> <u>The USER operates the medical device for a long time without maintenance.</u>	1	2	Acceptable		 The operating manual contains data on the frequency of maintenance and a list of necessary checks.	1	1	Acceptable

6.	Analyzer surfaces	<u>The cleaning disinfection and sterilization procedures were performed incorrectly</u> <u>Contact with surfaces; reuse of a disposable mouthpiece</u>	2	2	Invalid	Bacterial infection	1. The operating manual contains all the necessary requirements for processing methods and conditions.	1	2	Acceptable
* The described scenario includes a sequence of events, at least one of which is an OPERATION/APPLICATION ERROR, and influencing factors that, starting from DANGER, through a dangerous situation, lead to HARM. The underlined part of the description indicates an OPERATION/APPLICATION ERROR. The italicized part of the description indicates A DANGEROUS SITUATION.										

6. Residual risk

Table D.4 Residual risk matrix

Levels of severity of harm Probability levels of harm	Negligible (1)	Moderate (2)	Significant (3)
High (3)			
Average(2)	P ₃		
Low(1)	P ₅ , P ₆ , P ₈	P ₁ , P ₂ , P ₄ , P ₇	

P₁-P₈ is the risk number according to Table 3.

Symbols:

- gray color - unacceptable risk; - white color - acceptable risk

As you can see, all the predicted risks are in the range of acceptable values and are acceptable. Therefore, the residual risk can be considered insignificant.

7. Conclusions

Using Table 3 "Assessment of harm due to hazards related to operational suitability", a list of known and predicted hazards has been compiled, the sequence of events of at least one of which is an operation/application error and through a dangerous situation leads to harm, the predicted

consequences associated with the use of a medical device, the probability and severity of harm have been determined.

The residual risks of all hazards in a certain scenario were considered and evaluated, none of the hazards remained outside the risk analysis. The completeness of risk management is confirmed by a commission consisting of qualified technical personnel.

As a result of the risk analysis, 8 predicted dangerous situations were identified, the residual cumulative risk of MI is acceptable.

Nevertheless, the residual risk requires attention from the manufacturer and user of the medical device. It is necessary to provide for additional collection of information from consumers.

8. Conclusion

The design process has been carried out taking into account the operational suitability of a medical device: An analyzer for the concentration of ethanol vapors in the exhaled air of a Dingo with accessories. The residual cumulative risk of MI is acceptable.

Appendix A

List of documents referenced

The designation of the document to which the reference is given	Name of the document
GOST 8273-75	Wrapping paper. Technical conditions.
GOST 8828-89	Base paper and double-layer waterproof packaging paper. Technical conditions.
GOST 9569-2006	Waxed paper. Technical conditions.
GOST 10354-82	Polyethylene film. Technical conditions.
GOST 14192-96	Marking of goods.
GOST R ISO 15223-1-2014	Medical products. Symbols used for labeling on medical devices, labels and accompanying documentation. Part 1. Basic requirements.
GOST 31508-2012	Medical products. Classification depending on the potential risk of use. General requirements.
GOST R 50444-2020	Medical devices, apparatuses and equipment. General technical conditions.
GOST R IEC 60601-1-2010	Medical electrical products. Part 1. General safety requirements considering the main functional characteristics.
GOST R IEC 60601-1-2-2014	Medical electrical products. Part 1-2. General safety requirements, considering the main functional characteristics. Parallel standard. Electromagnetic compatibility. Requirements and tests

GOST R IEC 62366-1-2021	Medical products. Part 1. Design of medical devices considering operational suitability
GOST ISO 14971-2021	Medical products. Application of risk management to medical devices
MU-287-113 dated 30.12.1998	Guidelines for disinfection, pre-sterilization cleaning and sterilization of medical devices.
GOST R IEC 60601-1-6	Medical electrical products. Part 1-6. General safety requirements, taking into account the main functional characteristics. Additional standard. Operational suitability
GOST 30893.1	Basic rules of interchangeability. General tolerances. Limit deviations of linear and angular dimensions with unspecified tolerances

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PARK DASOM
attorney-in-fact of
YOO DO JOON president of
SENTECH KOREA CO.LTD

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appeared before me and admitted
said principal's subscription to
the attached Ethanol Vapor Concentration
Analyzer in Exhaled Air Dingo with
Accessories Productions

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This is hereby attested on this
23rd day of May 2023 at this office.

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Name of the office
Dae Seung Law and Notary Office Inc

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
Belong to Uijeongbu
District Prosecutor's Office

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Address of the office
SamyeongB/D 1F
9, Sicheong-ro, Paju-si Gyeonggi-Do, Korea

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Signature of the Notary Public
BAN SUNG CHAN

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This office has been authorized by the
Minister of Justice, the Republic of
Korea, to act as Notary Public Since
7, Feb. 2020 Under Law No.118.

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Do Joon Yoo

« 12 » May 2023

Ethanol Vapor Concentration Analyzer in Exhaled Air Dingo with accessories

Passport

**Registration certificate No. RZN 2014/1689
Registered in the Federal Information Fund
for Ensuring the Uniformity of Measurements No. 80132-20
2023**

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Authorized representative of the manufacturer in the territory of the Russian Federation:

Limited Liability Company «SIMS-2» (SIMS-2 LLC),
125430, Moscow, Mitinskaya street, 16, fl. 10, room 1012B, com. from 15 to 18
Phone: +7(495)792-31-90
Web-site: www.sims2.ru
E-mail address: info@sims2.ru

Testing Center:

Federal State Unitary Enterprise "All-Russian Research Institute of Metrology named after D.I.

Mendelev" (FSUE "VNIIM named after D.I. Mendelev»)

Address: 19, Moskovsky Ave., Saint Petersburg, 190005, Russia.

Phone: +7 812 251-7601

Web-site: <https://www.vniim.ru>

E-mail address: info@vniim.ru

Certificate of accreditation of FSUE "VNIIM named after D.I. Mendelev" for testing measuring instruments for type approval No. RA.RU.311541 dated 28.01.2016.

2. Operating conditions

The analyzer can be used both indoors and outdoors if the environmental parameters correspond to the operating conditions specified in Table 1.

Table 1 Operating conditions of the analyzer

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3. Metrological and technical characteristics

The metrological and technical characteristics of the ethanol vapor concentration analyzer in the exhaled air of the Dingo are given in Table 2.

Table 2 Metrological and technical characteristics of the analyzer

Name of the parameter	Specifications	
	Dingo E-200	Dingo E-200(B)

General information:		
Weight, g, no more	200	
Overall dimensions (H×W×D), mm no more	133×65×33	
Sensor Type	electrochemical removable	
Purge mode	1. Measuring mode - through a disposable mouthpiece 2. Screening mode - through a funnel mouthpiece, or without the use of mouthpieces	
Calibration interval	1 year	
Displaying results	4 – bit digital display	
Range of indications, mg/l	От 0,000 до 2,000	
The price of the lower division of the scale discharge, mg/l	0,005	
Parameters of the analyzed gas mixture when the sample is fed to the input of the analyzers (automatic sampling mode):	 	
- consumption of the analyzed gas mixture, l/min, not less		
- the volume of the sample of the analyzed gas mixture, l, not less	9	1,2
Time of preparation for work after switching on, s., no more	20	
Measurement time after sampling, sec., no more	10	
Preparation time for operation after analysis of a gas mixture with a mass concentration of ethanol 0.50 mg/l, sec., no more	120	
The time interval of the analyzers without adjusting the readings, months, no more	12	
The number of measurements on the analyzers without replacing the batteries, at least	3000	
Communication with the printer	1. Via USB port via a wired cable	1. Via USB port via a wired cable; 2. Via Bluetooth protocol
Power supply:		
Analyzer power supply	AA type battery	
Number of batteries, pcs.	2	
Operating voltage, V	3	

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4. Software

The analyzer has built-in software (software) designed to process measurement information. The firmware provides the following basic functions:

- processing of measurement information;
- diagnostics of the analyzer hardware;
- setting up the analyzer;
- display of measurement results on the display.

The influence of the embedded software on the metrological characteristics of the analyzers is considered when they are normalized. The level of software protection is "average" in accordance with P 50.2.077—2014.

The identification data of the software is shown in Table 3.

Table 3 Software identification data

Identification data (signs)	Meaning	
	Dingo E-200	Dingo E-200 (B)
Software identification name	ALP-1_v7.06	ALP-1_v7.05
Version number (identification number) software, not lower	Version 7.06	Version 7.05
Digital ID of the software	BAD1B438	927BF844
Algorithm for obtaining a digital identifier	CRC32	CRC32
Note – The values of the digital identifiers of the software specified in the table refer only to the firmware file of the specified version.		

5. Completeness of the measuring instrument

The scope of delivery of the analyzer is shown in Table 4.

Table 4 Scope of delivery of the analyzer

No	Name	Quantity
1	Breathalyzer Dingo E-200/ Dingo E-200(B)	1 pc.
2	AA type battery	2 pc.
3	Disposable mouthpiece	5 pc.
4	Arm strap	1 pc.
5	Passport	1 pc.
6	User Manual	1 pc.
Accessories:		
1	Analyzer case	1 pc. (if necessary)
2	Funnel mouthpiece	1 pc. (if necessary)
3	Printer	1 pc. (if necessary)
4	Case for transportation and storage	1 pc. (if necessary)
5	A case for transporting and storing the analyzer with a printer	1 pc. (if necessary)
6	Cable for connecting the analyzer to a PC	1 pc. (if necessary)
7	Cable for connecting the analyzer to the printer	1 pc. (if necessary)
8	Adapter to the printer from 220V	1 pc. (if necessary)
9	AA type battery for the printer	4 pc. (if necessary)
10	Paper for the printer	1 roll (if necessary)

6. Certificate of acceptance

Ethanol Vapor Concentration Analyzer in Exhaled Air Dingo with Accessories,

serial number _____

release date _____

complies with technical documentation and is recognized as serviceable.

 signature of the person responsible for acceptance

7. Initial verification

Verification date _____

The Verifier:

Signature

Name

Verification mark _____ Verification is completed

Verification is carried out according to the MP document-242-2327-2019 "GSI. Analyzers of ethanol vapor concentration in exhaled air Dingo E-200, Dingo E-200(B). Verification methodology", approved by FSUE "VNIIM named after D.I. Mendeleev" on August 22, 2019. The verification methodology is provided to organizations accredited for verification upon request.

Basic means of verification: working standards of 1 or 2 categories in accordance with

The state verification scheme for measuring the content of components in gas and gas condensate media, approved by Rosstandart Order No. 2664 dated 12/14/2018 – generators of gas mixtures of ethanol vapors in the air complete with standard samples of the composition of aqueous ethanol solutions VRE-2 GSO 8789-2006. Limits of permissible relative error $\pm 5\%$.

It is allowed to use similar means of verification that ensure the determination of metrological characteristics of the verified SI with the required accuracy.

The verification mark is applied to the verification certificate and/or to the passport (during the initial verification before commissioning).

8. Warranty conditions

SIMS-2 LLC, being the official representative of the manufacturer, guarantees that the analyzer you purchased has no manufacturing defects in terms of materials and components at the time of sale and undertakes to repair the failed elements free of charge during the entire warranty period, with the exception of routine maintenance caused by natural wear of the sensor sensor.

The warranty period of operation for a medical device is 12 months.

The warranty period of storage for a medical device is 3 months.

Attention! The warranty does not apply to chemical power supplies.

The analyzer is accepted for warranty service only if there is a passport (operating manual) indicating the serial (factory) number, date of sale and with the seal of the trading organization.

SIMS-2 LLC reserves the right to refuse free warranty service in the following cases:

1. Documents confirming the date of sale of the analyzer are lost or incorrectly filled in.
2. The analyzer was subjected to unauthorized opening.

3. The analyzer was used in violation of the rules of operation.
4. The analyzer has traces of mechanical damage caused by blows, falling, or attempts to open.
5. The analyzer has damage caused by foreign objects, liquids or insects entering the device.
6. The analyzer has damage caused by the use of non-standard components and accessories.

In case of refusal of warranty service, the buyer is issued a certificate of technical expertise with justification of the reason for the refusal.

During the warranty period, if factory defects are detected, the user can contact the authorized representative of the manufacturer in the Russian Federation - SIMS-2 LLC at the address: 125430, Moscow, Mitinskaya str., 16, fl. 10, room 1012B, com. from 15 to 18, phone: (495) 792-31-90 (multi-channel), email address: info@sims2.ru.

WARRANTY CARD

For repair (replacement) during the warranty period

Ethanol Vapor Concentration Analyzer in Exhaled Air Dingo

(serial number)

Purchased _____
(date, signature, stamp of the trading organization)

Manufacturer:

Sentech Korea Corp.,
Address: 21-6, Jimok-ro, 75 beon-gil, Paju-si, Gyeonggi-do, Korea.
Phone: (82 31) 80714400
Web-site: www.sentechkorea.com
E-mail address: sentech@paran.co.kr

Place of production:

Sentech Korea Corp.,
Address: 21-6, Jimok-ro, 75 beon-gil, Paju-si, Gyeonggi-do, Korea
Phone: (82 31) 80714400
Web-site: www.sentechkorea.com
E-mail address: sentech@paran.co.kr

Limited liability company «ARIDES»

(ARIDES LLC), Republic of Armenia

Address: 111 Raffi str., Malatia-Sebastia, Yerevan, 0064

Phone: +37411 26 99 50

Web-site: arides.am

E-mail address: info@arides.am

9. Maintenance

Maintenance and repair of the analyzer, including warranty repairs, is carried out on the territory of the Russian Federation by the authorized distributor of the manufacturer - SIMS-2 LLC at the address: 125430, Moscow, Mitinskaya Street, 16, fl. 10, room 1012B, com. from 15 to 18.

tel: (495) 792-31-90 (multichannel), e-mail: info@sims2.ru .

Maintenance of the analyzer is carried out in order to ensure constant serviceability and readiness for operation.

Daily maintenance of the analyzer includes an external inspection.

Periodic maintenance of the analyzer during the entire period of operation includes:

- correction of the analyzer readings – when checking if necessary, or when the inscription "Calibration expired" appears on the analyzer display;
- analyzer verification – 1 time per year;

The analyzer readings are adjusted during verification as necessary or when the inscription "Calibration expired" appears on the analyzer display. The manufacturer has programmed the analyzer to issue this message every 365 days.

After the correction of the readings, the analyzer is necessarily subject to verification.

Note – The message on the analyzer display "Calibration expired" is informational. At the time of measurement, the analyzer must be verified, the date of verification is indicated in the certificate of verification or the passport of the analyzer (during initial verification before commissioning).

Instructions for correcting readings are provided by the official representative of the manufacturer in Russia - SIMS-2 LLC, upon a separate request to organizations accredited for verification or service centers.

It is recommended to note the work on correcting the readings in the analyzer passport (in the maintenance section).

The analyzers are verified according to the MP document-242-2327-2019 "GSI. Analyzers of ethanol vapor concentration in exhaled air Dingo E-200, Dingo E-200 (V). Verification method". The interval between verifications is 1 year.

The analyzers are verified by accredited legal entities or individual entrepreneurs.

Persons who are familiar with GOST R 54794-2011, GOST R 50444-2020, Rosstandart Order No.

2664 dated 12/14/2018, the operating manual of verifiable analyzers and reference measuring instruments, who are qualified at least as an engineer and have passed safety instructions, are allowed to work with analyzers and carry out verification.

Basic means of verification:

– Working ethanols of the 1st or 2nd category in accordance with the State verification scheme for measuring the content of components in gas and gas condensate media, approved by Rosstandart Order No. 2664 dated 12/14/2018 – generators of gas mixtures of ethanol vapors in the air complete with standard samples of the composition of aqueous ethanol solutions VRE-2 GSO 8789-2006.

Limits of permissible relative error $\pm 5\%$.

or

– Working ethanols of the 1st category in accordance with the State verification scheme for measuring the content of components in gas and gas condensate media, approved by Rosstandart Order No. 2664 dated 12/14/2018 – standard samples of the composition of ethanol/nitrogen gas mixtures in pressure cylinders GSO 10338-2013. The limits of the relative error at $P = 0.95 \pm (2-4)\%$.

It is allowed to use similar means of verification that ensure the determination of the metrological characteristics of the analyzers being verified with the required accuracy.

The verification mark is applied to the verification certificate or to the passport (during the initial verification before commissioning).

Attention! The responsibility for the timely verification is borne by the owner of the analyzer.

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appeared before me and admitted
said principal's subscription to
the attached Ethanol Vapor Concentration
Analyzer in Exhaled Air Dingo with
accessories Passport

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