I hereby certify accuracy, correctness and
reliability of this document
translated into Russian

Directions for medical device use

LifeStream balloon expandable endovascular covered stent- graft with catheter

clearStream
MEDICAL TECHNOLOGIES LTD.

30 MAR 2016

MOYNE UPPER
ENNISCORTHY, CO. WEXFORD.

«APPROVED BY»
“ClearStream Technologies, Ltd.” Ireland

Fiona Ní Mhulláin,
Senior RA/QA Manager
ClearStream Technologies Ltd.

Signed:

F Ní Mhulláin

Date:

30 March 16

Subscribed and sworn before me
by the said *Fiona Ní Mhulláin*
at Enniscorthy in the County of
Wexford this *30* day of *March* 2016

Anthony H. Ensor
Anthony H. Ensor, Solicitor and
Notary Public for the County of
Wexford

CONTENTS

1. Name of the medical device	3
2. Medical device manufacturer information	3
3. Purpose of the medical device assigned by the manufacturer	3
4. Description of the medical device	4
5. Technical description of the medical device	5
6. Directions for use	11
7. Magnetic resonance imaging (MRI) information.....	13
8. Biocompatibility verification	14
9. Safety requirements	14
10. Environmental protection requirements while applying the medical device ..	15
11. Potential adverse effects of the medical device use	17
12. Methods of sterilization	17
13. Labeling, operation conditions, transportation and storage	17
14. Packaging	19
15. Warranty	19
16. Shelf life	20
17. Procedure for the medical device disposal or destruction	20

1. Name of the medical device

LifeStream balloon expandable endovascular covered stent-graft with a catheter, design variants

1. LifeStream Balloon Expandable Endovascular Covered Stent-graft with a 80mm catheter, sizes: 5 mm x 26 mm, 5 mm x 37 mm, 6 mm x 16 mm, 6 mm x 26 mm, 6 mm x 37 mm, 6 mm x 58 mm, 7 mm x 16 mm, 7 mm x 26 mm, 7 mm x 37 mm, 7 mm x 58 mm, 8 mm x 16 mm, 8 mm x 26 mm, 8 mm x 37 mm, 8 mm x 58 mm, 9 mm x 38 mm, 9 mm x 58 mm, 10 mm x 38 mm, 10 mm x 58 mm, 12 mm x 38 mm, 12 mm x 58 mm.

2. LifeStream Balloon Expandable endovascular Covered Stent-graft with a 135mm catheter, sizes: 5 mm x 26 mm, 5 mm x 37 mm, 6 mm x 16 mm, 6 mm x 26 mm, 6 mm x 37 mm, 6 mm x 58 mm, 7 mm x 16 mm, 7 mm x 26 mm, 7 mm x 37 mm, 7 mm x 58 mm, 8 mm x 16 mm, 8 mm x 26 mm, 8 mm x 37 mm, 8 mm x 58 mm, 9 mm x 38 mm, 9 mm x 58 mm, 10 mm x 38 mm, 10 mm x 58 mm, 12 mm x 38 mm, 12 mm x 58 mm.

2. Medical device manufacturer information

Manufacture /Designer: "ClearStream Technologies, Ltd.", Moyne Upper, Enniscorthy, Co. Wexford, Ireland, Tel: +353 0 53923 7111

Medical device production facility location: Moyne Upper, Enniscorthy, Co. Wexford, Ireland

Authorized representative of the manufacturer in the Russian Federation:

Limited liability company «Bard Rus» (Bard Rus, LLC), Building 204, 2 Ul. Gorbunova, Moscow, Russia, 121596. Phone: +7 (499) 372-5002, fax: +7 (499) 372-5003

3. Purpose of the medical device assigned by the manufacturer

LifeStream balloon expandable endovascular covered stent-graft with a catheter (hereinafter, LifeStream stent-graft, stent-graft, device, implant) is intended for treatment of atherosclerotic lesions in common and external iliac arteries areas.

Indications for use:

- Atherosclerotic lesions in common and external iliac arteries.

Contraindications:

- Allergy or hypersensitivity to stainless steel, nickel, chromium, or ePTFE;
- Uncorrected bleeding disorders;
- Contraindications to anticoagulation and (or) antiplatelet agents;
- Damage that prevents complete extension of the stent-graft;
- Lesions in which the lumen diameter post balloon angioplasty is insufficient for the passage of the catheter;
- Lesion locations subject to external compression.

Conditions for application:


This medical device should be used by specially personnel trained in specialized medical institutions only..

4. Description of the medical device

LifeStream balloon expandable endovascular covered stent-graft with a catheter, design variants, is comprised of an electropolished balloon-expandable stent-graft made from stainless steel and encapsulated between two layers of ePTFE. (See Fig. 1)

Fig. 1 LifeStream balloon expandable endovascular covered stent-graft

The photographic image of the stent-graft, expanded with the balloon to a nominal size

The covered stent-graft is supplied pre-mounted on a delivery system consisting of a 80 cm or 135 cm long catheter which uses the balloon to expand the stent-graft in the area of the iliac artery to be manipulated (See Fig. 2)

Fig. 2 The catheter of the LifeStream balloon expandable endovascular covered stent-graft

Two radiopaque markers (A and B) located on the balloon catheter indicate the length of the balloon and the ends of the covered stent-graft and aid in accurate stent-graft deployment. The proximal end of the catheter includes an inflation female luer lock hub (C) and a guide wire female luer lock hub (D).

A luer lock hub meets the requirements of ISO 594: the inner smaller diameter (for information, as this is the specification for the syringe hub with exterior thread): 3,925 -3,990 mm, the inner diameter of the larger orifice of the taper: 4,270 -4,315 mm, the minimum length of the tapered hub is 7.500 mm. The maximum diameter of the hub outer part should not exceed 6.73 mm. Maximum bend radius near the orifice should not exceed R 0.5 mm.

The catheter is compatible with a guide wire with a diameter of 0.035 inch (0.89 mm). Pre-mounted stent-grafts of different diameters and lengths are commercially available.

Contents

The box of the LifeStream balloon expandable endovascular covered stent-graft with a catheter includes:

- Stent-graft pre-mounted on a catheter - 1.
- Instructions for use - 1.
- Patient Implant Card - 1.

Medical devices required for the placement of LifeStream balloon expandable endovascular covered stent-graft with a catheter:

- Heparinized saline;
- Sterile syringes;
- 0.035" (0.89 mm) guide wire with a length at least twice as long as the endovascular system;
- Introducer sheath with appropriate inner diameter;
- Balloon angioplasty catheter for the pre-and/or post-dilatation, if required;
- Inflation device with pressure monitoring capability;
- Diagnostic catheters and accessories;
- Contrast medium.

5. Technical description of the medical device

Stent-graft protective cover

Stent-graft

Balloon

Catheter coupling hub

Catheter

Distal radiopaque mark

Proximal radiopaque mark

Y-shaped connector

Materials:

- Stent-graft - stainless steel 316L, ePTFE coated - high molecular weight F-107/FX-9015;
- Balloon - Vestamid L2140 polymer;
- Catheter - 83.5% PEBAX 7233, 11.5% barium sulfate, 1.0% erucamide Aldrich, purple dye: 4.0% Pantone PMS 267C;
- Catheter coupling hub - Wacker LR 3003/80A US, dye 2% PMS267C;
- Y-shaped connector - polycarbonate Makrolon RX 1805-451118;
- Radiopaque markers (distal and proximal radiopaque marks) - 90% platinum and 10% iridium;
- Stent-graft protective cover - teflon, Translucent green dye.

Rated burst pressure for all sizes - 1216 kPa

Nominal pressure for all sizes - 811 kPa.

The scale of pressure depending on the diameter of the balloon

kPa	5.0 mm	6.0 mm	7.0 mm	8.0 mm	9.0 mm	10.0 mm	12.0 mm
811	5.1	6.0	7.1	8.2	9.1	10.1	12.0
912	5.1	6.0	7.2	8.3	9.2	10.2	12.0
1013	5.1	6.1	7.2	8.4	9.2	10.2	12.1
1115	5.2	6.1	7.3	8.4	9.3	10.3	12.2
1216	5.2	6.2	7.3	8.5	9.4	10.4	12.3

Balloon compliance chart and stent-graft dimensions

Expanded stent-graft diameter, mm	Stent-graft length, mm	Nominal sizes (after inflation, 8 atm. (811 kPa)), mm		Rated sizes (after inflation, 12 atm. (1216 kPa)), mm	
		Mean diameter, mm	Average length, mm	Mean diameter, mm	Average length, mm
5	26	4.4	26.4	5.0	26.2
5	37	5.0	36.8	5.3	36.4
6	16	5.5	16.1	5.9	15.9
6	26	5.7	25.9	6.0	25.7
6	37	5.7	36.5	6.0	36.1
6	58	5.9	56.8	6.2	56.4
7	16	6.9	15.5	7.3	15.4
7	26	6.8	25.2	7.1	24.9
7	37	6.6	35.6	7.1	35.3
7	58	6.7	55.7	7.1	55.1
8	16	7.7	15.3	8.2	15.0
8	26	7.9	24.6	8.3	24.3
8	37	7.9	34.7	8.3	33.9
8	58	7.8	55.5	8.4	54.4
9	38	8.6	36.5	9.0	36.2
9	58	8.6	56.0	9.0	55.9
10	38	9.7	35.7	10.2	34.8
10	58	9.6	55.4	10.1	54.5
12	38	11.5	33.6	12.2	33.5
12	58	11.4	53.9	12.0	52.6

*Notes: approximate values

*Permissible deviations $\pm 5\%$, unless otherwise indicated.

Balloon and stent-graft compliance chart (nominal pressure and post rated burst pressure inflation) are indicated on the device label.

The diameters of the distal and proximal ends of the catheters

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External diameter of the introducer compatible with the stent-graft

Diameter of the expanded stent-graft, mm	Stent-graft length, mm	Introducer sheath external diameter
5	26	6 F (1.98 mm)
5	37	6 F (1.98 mm)
6	16	6 F (1.98 mm)
6	26	6 F (1.98 mm)
6	37	6 F (1.98 mm)
6	58	6 F (1.98 mm)
7	16	6 F (1.98 mm)
7	26	6 F (1.98 mm)
7	37	6 F (1.98 mm)
7	58	7 F (2.31 mm)
8	16	7 F (2.31 mm)
8	26	7 F (2.31 mm)
8	37	7 F (2.31 mm)
8	58	7 F (2.31 mm)
9	38	7 F (2.31 mm)
9	58	7 F (2.31 mm)
10	38	8 F (2.64 mm)
10	58	8 F (2.64 mm)
12	38	8 F (2.64 mm)
12	58	8 F (2.64 mm)

Force transmitted by the stent-graft to surrounding tissues: 0.10-0.25 N/mm

Expanded stent-graft diameter, mm	Stent-graft length, mm	Forces transmitted by the stent-graft to surrounding tissues
5	26	0.10 N/mm
5	37	
6	16	0.12 N/mm
6	26	
6	37	
6	58	
7	16	0.15 N/mm
7	26	
7	37	
7	58	
8	16	0.18 N/mm
8	26	
8	37	
8	58	
9	38	0.20 N/mm
9	58	
10	38	0.22 N/mm
10	58	
12	38	0.25 N/mm
12	58	

The forces required to deform the stent-graft: 0.5-1.25 N

Expanded stent-graft diameter, mm	Stent-graft length, mm	The forces required to deform the stent-graft
5	26	0.5 N
5	37	
6	16	0.8 N
6	26	
6	37	
6	58	
7	16	1.05 N
7	26	
7	37	
7	58	
8	16	1.10 N
8	26	
8	37	
8	58	
9	38	1.15 N
9	58	
10	38	1.20 N
10	58	
12	38	1.25 N
12	58	

Stent-graft specifications

Range of expected lumen diameters, mm	Percentage of the stent-graft shortening after implantation	Minimum effective length, mm (± 2 mm)	Maximum effective length (stent-graft working length), mm	Stent-graft length in the delivery system, mm
5.1 - 5.2	$\leq 15\%$	22.0 ± 2 mm (vessel diameter - 5.2 mm)	26.0 ± 2 mm (vessel diameter - 5.1 mm)	26.0 ± 3 mm
	$\leq 15\%$	31.0 ± 2 mm (vessel diameter - 5.2 mm)	37.0 ± 2 mm (vessel diameter - 5.1 mm)	37.0 ± 3 mm
6.0 - 6.2	$\leq 17\%$	13.0 ± 2 mm (vessel diameter - 6.2 mm)	16.0 ± 2 mm (vessel diameter - 6.0 mm)	16.0 ± 3 mm
	$\leq 17\%$	22.0 ± 2 mm (vessel diameter - 6.2 mm)	26.0 ± 2 mm (vessel diameter - 6.0 mm)	26.0 ± 3 mm
	$\leq 17\%$	31.0 ± 2 mm (vessel diameter - 6.2 mm)	37.0 ± 2 mm (vessel diameter - 6.0 mm)	37.0 ± 3 mm
	$\leq 17\%$	48.0 ± 2 mm (vessel diameter - 6.2 mm)	58.0 ± 2 mm (vessel diameter - 6.0 mm)	58.0 ± 3 mm
7.1 - 7.3	$\leq 19\%$	13.0 ± 2 mm (vessel diameter - 7.3 mm)	16.0 ± 2 mm (vessel diameter - 7.1 mm)	16.0 ± 3 mm

Range of expected lumen diameters, mm	Percentage of the stent-graft shortening after implantation	Minimum effective length, mm (± 2 mm)	Maximum effective length (stent-graft working length), mm	Stent-graft length in the delivery system, mm
	$<19\%$	21.0 $\pm$ 2 mm (vessel diameter - 7.3 mm)	26.0 $\pm$ 2 mm (vessel diameter - 7.1 mm)	26.0 $\pm$ 3 mm
	$<19\%$	30.0 $\pm$ 2 mm (vessel diameter - 7.3 mm)	37.0 $\pm$ 2 mm (vessel diameter - 7.1 mm)	37.0 $\pm$ 3 mm
	$\leq 19\%$	47.0 $\pm$ 2 mm (vessel diameter - 7.3 mm)	58.0 $\pm$ 2 mm (vessel diameter - 7.1 mm)	58.0 $\pm$ 3 mm
8.2 - 8.5	$<21\%$	13.0 $\pm$ 2 mm (vessel diameter - 8.5 mm)	16.0 $\pm$ 2 mm (vessel diameter - 8.2 mm)	16.0 $\pm$ 3 mm
	$\leq 21\%$	21.0 $\pm$ 2 mm (vessel diameter - 8.5 mm)	26.0 $\pm$ 2 mm (vessel diameter - 8.2 mm)	26.0 $\pm$ 3 mm
	$<21\%$	29.0 $\pm$ 2 mm (vessel diameter - 8.5 mm)	37.0 $\pm$ 2 mm (vessel diameter - 8.2 mm)	37.0 $\pm$ 3 mm
	$<21\%$	46.0 $\pm$ 2 mm (vessel diameter - 8.5 mm)	58.0 $\pm$ 2 mm (vessel diameter - 8.2 mm)	58.0 $\pm$ 3 mm
9.1 - 9.4	$\leq 23\%$	29.0 $\pm$ 2 mm (vessel diameter - 9.4 mm)	38.0 $\pm$ 2 mm (vessel diameter - 9.1 mm)	38.0 $\pm$ 3 mm
	$<23\%$	45.0 $\pm$ 2 mm (vessel diameter - 9.4 mm)	58.0 $\pm$ 2 mm (vessel diameter - 9.1 mm)	58.0 $\pm$ 3 mm
10.1 - 10.4	$\leq 24\%$	29.0 $\pm$ 2 mm (vessel diameter - 10.4 mm)	38.0 $\pm$ 2 mm (vessel diameter - 10.1 mm)	38.0 $\pm$ 3 mm
	$\leq 24\%$	44.0 $\pm$ 2 mm (vessel diameter - 10.4 mm)	58.0 $\pm$ 2 mm (vessel diameter - 10.1 mm)	58.0 $\pm$ 3 mm
12.0 - 12.3	$\leq 25\%$	29.0 $\pm$ 2 mm (vessel diameter - 12.3 mm)	38.0 $\pm$ 2 mm (vessel diameter - 12.0 mm)	38.0 $\pm$ 3 mm
	$\leq 25\%$	44.0 $\pm$ 2 mm (vessel diameter - 12.3 mm)	58.0 $\pm$ 2 mm (vessel diameter - 12.0 mm)	58.0 $\pm$ 3 mm


Stent-graft wall thickness, mm	Stent-graft profile (d, external diameter)	Profile form
0.07 mm	5.1 $\pm$ 0.3 mm	
	6.0 $\pm$ 0.3 mm	
	7.1 $\pm$ 0.3 mm	
	8.2 $\pm$ 0.3 mm	
	9.1 $\pm$ 0.3 mm	
	10.1 $\pm$ 0.3 mm	
	12.0 $\pm$ 0.3 mm	

Fig.1

Expanded stent-graft diameter, mm	Stent-graft length, mm	Radial strength of the stent-graft	Minimum bend radius of the stent-graft	Relation between metal surface area and intima in the stented segment
5	26	Minimum strength – 0.150 N, Maximum strength – 0.340 N	15.4 mm	≤30%
5	37		17.0 mm	
6	16			
6	26			
6	37		19.4 mm	
6	58			
7	16			
7	26		20.0 mm	
7	37			
7	58			
8	16		22.0 mm	
8	26			
8	37			
8	58			
9	38	Minimum strength – 0.300 N, Maximum strength – 0.340 N	22.0 mm	
9	58		24.2 mm	
10	38			
10	58		25.0 mm	
12	38			
12	58			

Stent-graft functioning factors:

- Probability of contact, galvanic, spot and fretting corrosion, probability of oxidation: probability is low and guaranteed by the application of bioinert materials;
- Probability of relative stent-graft and tissue dislodgement: where the stent-graft diameter is selected depending on the vessel diameter, there is no plane-parallel compression, probability of relative stent-graft and tissue dislodgement is negligible;
- Expected period for the stent-graft ingrowth, penetration into the tissue, tissue perforation, stent-graft dislodgement and migration - about 14 weeks;
- Information on biodegradation: none, due to the use of bioinert materials.

6. Directions for use

Site access and preparation

1. Using standard techniques access the artery and place the introducer with a corresponding outer diameter and the guide wire with a diameter of 0.035 inch (0.89 mm) to the target lesion site.
2. Perform diagnostic angiography to confirm the site of implantation and measure the reference vessel diameter and lesion length.

Covered stent-graft size selection

3. Select a covered stent-graft diameter that is approximately 5-20% larger than the largest reference vessel diameter at the proximal or distal target site.
Referred to the sizing table on the packaging label for appropriate selection of the colored stent-graft diameter and length.

Endovascular system preparation

4. Carefully remove the selected device from the package.
5. Inspect the covered stent-graft for adherence to the balloon and stent-graft placement between the balloon markers. Do not use the device, if the stent-graft is not centered and (or) does not firmly adhere to the balloon.
6. Flush delivery system guide wire lumen with sterile saline mixture until saline drops from the distal end of the endovascular system.

Air evacuation

7. Use the recommended 20 cc or smaller luer-lock syringe with a minimum of 5 cc's sterile saline mixture for aspirating this device.
8. With the distal balloon tip pointing down and positioned below the level of the syringe, pull negative pressure until all air is expelled.
9. Induce a negative pressure to remove all air from the balloon and inflation lumen. Repeat until all the air is expelled.
10. Carefully release to neutral. Allow the inflation lumen to fill with the diluted contrast medium and maintain a neutral pressure.

Important: Do not apply positive pressure to the balloon.

11. Attach the prefilled inflation device to the inflation lumen of the catheter hub intended to inflate the balloon, ensuring no air bubbles remain at the catheter connection.
12. Make sure that the covered stent-graft is still centered between the two radiopaque markers on the balloon catheter.

The introduction of the device and placement of the covered stent-graft

13. Advance the device over the guide wire into the introducer sheath.
14. Further advance the device to the target treatment site placement with the introducer sheath and position the stent-graft across the lesion. Verify that the covered stent-graft is still centered between the two markers of the balloon. Slowly retract the introducer sheath while maintaining the position of the covered stent-graft. Ensure the introducer sheath is retracted far enough to not compromise the balloon expansion and the stent-graft release.

15. Slowly inflate the device to the nominal pressure, expanding the covered stent-graft. Confirm complete extension via fluoroscopic visualization. A 15-30 second inflation time is recommended.

Important: Do not exceed the Rated Burst Pressure of the delivery system.

16. After the covered stent-graft deployment, apply negative pressure to the balloon until it is fully deflated. By using fluoroscopic visualization and maintaining negative pressure, withdrawal the delivery system from the body with the guide wire remaining across the lesion.

Important: After deployment of the covered stent-graft, if resistance is felt during the device withdrawal, carefully remove it and the guide wire as a single unit.

If, while removing the device and guide wire as a single unit, resistance is still there, remove the guide wire, device and introducer sheath as a single unit.

17. Confirm optimal covered stent-graft adherence to the vessel walls using standard angiographic techniques. If the covered stent-graft has not achieved full wall position, a post dilation with an appropriately sized balloon should be performed. 5-8 mm devices may be post-dilated with balloons up to 10 mm in diameter. 9-12 mm devices may be post-dilated with balloons up to 12 mm in diameter.

Note: Blood flow impedance remaining after the covered stent-graft deployment (for example, residual proximal or distal stenosis, or splitting, or weak peripheral outflow) may increase the risk of the stent-graft thrombosis. After the procedure was completed, inflow and outflow of blood should be evaluated, and, if necessary, additional steps should be taken (e.g., additional post-dilatation, additional stenting or distal bypass grafting) to support the sufficient inflow and outflow.

18. Remove the guide wire and introducer sheath from the patient.

19. Close the incision as appropriate.
20. Dispose of the device in accordance with the instructions for use.

Placement of two overlapped covered stent-grafts

21. In the event that two stent-grafts must be placed overlapped, place the distal stent-graft first. Ensure approximately 10 mm overlap zone.

Patient implant information card

A Patient Implant Information Card is provided with this device. The patient data, implant data and hospital data should be accurately recorded on the card and given to the patient.

Apply one of the peel-off stickers found on the product label on the pouch to the indicated area on the Patient implant information card. This peel-off sticker contains important information about the patient's covered stent-graft implant. The patient should carry the Scott with them and provide to any medical personnel caring for the patient in the future.

7. Magnetic Resonance Imaging (MRI) information

LifeStream stent-graft was determined to be MR conditional according to the terminology specified in the American Society for Testing and Materials (ASTM) International, Designation: F2503-08 - Standard Practice for Marking Medical Devices and Other Items for Safety in the Magnetic Resonance Environment. ASTM International, 100 Barr Harbor Drive, PO Box C700. West Conshohocken, Pennsylvania, in 2008.

Non-clinical testing demonstrated that the LifeStream stent-graft is MR Conditional. A patient with this implant can be scanned safely immediately after the implant placement under the following conditions:

- Static magnetic field of 3 Tesla or 1.5 Tesla;
- Spatial gradient of the magnetic field not exceeding 3000 gauss/cm;
- Maximum MR system reported whole-body-averaged Specific Adsorption Rate (SAR) of 2 W/kg in the normal operating mode.

In non-clinical testing, the LifeStream Stent-graft produced a temperature rise of 3.5 °C at a maximum MR system-reported whole body averaged specific absorption rate (SAR) of 2-W/kg for 15 minutes of continuous MR scanning in a 3-Tesla MR system using a transmit/receive body coil (Signa HDx, General Electric Healthcare, Milwaukee, WI). The actual in-vivo temperature rise is expected to be less as this calculation did not include the cooling due to blood flow in the lumen of the stent-graft and blood perfusion in the tissue outside the stent-graft. The effect of heating in the MRI environment for stent-grafts with fractured struts is not known.

In non-clinical testing, the LifeStream stent-graft produced a temperature rise of 3.4 °C at a maximum MR system-reported whole body averaged specific absorption rate (SAR) of 2-W/kg for 15 minutes of continuous MR scanning in a 1.5-Tesla MR test coil (Model 46-306600G2, General Electric Healthcare, Milwaukee, WI). The actual in-vivo rise is expected to be less as this calculation did not include the cooling due to blood flow in the lumen of the stent-graft and blood perfusion in the tissue outside the stent-graft. The effect of heating in the MRI environment for stent-grafts with fractured struts is not known.

MR image quality may be compromised if the area of interest is located in the exact same area or relatively close to the position of the LifeStream stent-graft. Therefore, optimization of MR imaging parameters to compensate for the presence of this implant may be necessary.

Artifact information:

Image artifact of the LifeStream stent-graft was tested according to ASTM F2119-07 in a GE Signa HDx 3-Tesla cylindrical bore scanner. Maximum artifact beyond the stent-graft was 13 mm for the spin echo sequence and 16 mm for the gradient echo sequence, both inside and outside the device lumen. The lumen was obscured. Imaging parameters may need to be adjusted for artifact optimization.

It is recommended that patients with a stent-graft register the MR conditions with the MedicAlert Foundation (www.medicalert.org).

8. Biocompatibility verification

ISO 10993-1 (2009) defines a LifeStream stent-graft as an implantable device permanently in contact with the circulating blood. The standard recommends conducting the following stent-graft tests: cytotoxicity, sensitization, excitability or intradermal reactivity, acute systemic toxicity, hemocompatibility, genotoxicity, and testing by applying muscle implantation. ISO 10993-1 (2009) defines a delivery system (catheter) as an externally attachable product and contacting circulating blood during a limited period of time. The standard recommends conducting the following catheter tests: cytotoxicity, sensitization, excitability or intradermal reactivity, acute systemic toxicity, and hemocompatibility.

Biocompatibility test results confirmed absence of toxic or other adverse reactions. Moreover, the stent-graft and catheter material platform confirms their clinical safety for many years. For these reasons, we conclude that the device meets the requirements of ISO 10993-1.

9. Safety requirements

Warnings

- The LifeStream stent-graft is supplied sterile and is intended for single use only. Do not re-sterilize and/or reuse the device. Reuse, re-sterilization, reprocessing and/or repackaging may create a risk of patient or user infection, compromise the structural integrity and/or essential material and design characteristics of the device, which may lead to device failure, and/or injury, illness or death of the patient.
- Patients with allergies or hypersensitivity to stainless steel, nickel, chromium or ePTFE may suffer an allergic response to this implant.
- Do not use if packaging is damaged.
- Use the device prior to the USE BY date specified on the package.
- Stenting across a vessel side branch may impede blood flow and hinder or prevent future procedures.
- Stenting into a bifurcation may compromise future diagnostic or therapeutic procedures.
- Should excessive resistance be felt at any time during the insertion process, do not force passage.
- Attempts to retract the stent-graft into the sheath may result in dislodgement and embolization of the stent-graft. Maintain guide wire placement across the lesion and withdraw the device only until the proximal end of the stent-graft is aligned with the tip of the sheath. Do not attempt to remove an unexpanded stent-graft through the sheath. Attempting to remove an unexpanded covered stent-graft by pulling it back into the sheath may result in stent-graft dislodgement.
- Do not exceed the maximum rated burst pressure since this increases the potential for balloon rupture and vessel damage. The use of a pressure monitoring device is recommended to prevent over pressurization.
- Use only diluted contrast medium for balloon inflation. Do not use air or any gaseous medium to inflate the balloon as this may cause uneven expansion and difficulty in deployment of the stent-graft. For the contrast/saline solution, a ratio of 50/50 is recommended.

- The stent-graft cannot be repositioned after it is deployed.
- Do not retract the balloon until the balloon is fully deflated under vacuum.
- Do not expose the stent-graft to temperatures higher than 260 °C (500 °F). ePTFE decomposes at elevated temperatures, producing highly toxic decomposition products.
- Use of a laser on or around the surface of the stent-graft may result in damage to the stent-graft and could create toxic fumes, which may harm the patient or operator.
- After use, the catheter is a potential biohazard. Handle and dispose of it in accordance with accepted medical practice and applicable local, state and federal laws and regulations.

Precautions

- The device should only be used by physicians who are trained in endovascular procedures and are familiar with the complications, side effects and hazards of peripheral vascular interventions.
- Anatomical variances may complicate the procedure; use caution when advancing the endovascular system through tortuous or difficult anatomy.
- Prior to device use, refer to the Stent-graft sizing table on the label and read the Instructions for Use.
- Do not use if the delivery system cannot be properly flushed.
- Crossing the implant with a catheters or other adjunct devices can result in stent-graft dislodgement or damage.
- This device has not been tested for use in overlapped conditions with stent-grafts or stent-grafts from other manufacturers.
- Store in a cool and dry place. Keep away from sunlight.

10. Environmental protection requirements while applying the medical device

The medical device does not have a negative impact on people and the environment during its use, transportation and storage.

11. Potential adverse effects of the medical device use

Potential patient/device adverse effects:

- Abscess
- Allergic/anaphylactoid reaction
- Amputation
- Aneurysm/pseudoaneurysm
- Angina/coronary ischemia
- Arterial occlusion/thrombus, near the puncture site
- Arterial occlusion/thrombus, remote from puncture site
- Arterial occlusion/restenosis of the treated vessel
- Arteriovenous fistula
- Arrhythmia
- Balloon rupture
- Blockage of major collateral artery or arterial branch
- Bypass surgery
- Stent-graft dislodgement from balloon during tracking procedure
- Stent-graft misplacement during placement procedure
- Stent-graft migration post placement procedure
- Stent-graft insufficient wall apposition
- Stent-graft deformation/kink/fracture
- Death
- Distal embolization





- Drug reaction or allergic reaction to medication, substances or materials used for the procedure (e.g. anticoagulation or antiplatelet agent, contrast medium, stent-graft or catheter materials)
- Edema
- Fever
- Hemorrhage/bleeding
- Hematoma and/or bleeding at puncture (access) site
- Hypotension/hypertension
- Inability to introduce/withdraw endovascular system
- Inability to track endovascular system to the target lesion
- Inability to inflate the balloon/deploy stent-graft
- Infection at access site
- Infection at or around implant
- Inflammation
- Ischemia/infarction of tissue/organ
- Myocardial infarction
- Pain
- Radiation injuries
- Renal insufficiency/failure/toxicity
- Respiratory arrest
- Restenosis in the treatment area/stent-graft edge
- Sepsis
- Shock
- Stroke/Transient Ischemic Attack (TIA)
- Thromboembolic event/thrombosis
- Vasospasm
- Vessel wall trauma, perforation/dissection/rupture

12. Method of sterilization

This medical device is supplied sterile, using proven methods (sterilization using ethylene oxide) that ensures the required level of sterility (SAL) 10^{-6} . This medical device is intended for single use only, should be considered sterile if the packaging is closed and intact. Do not use the device if there is any doubt about its sterility.

13. Labeling, operation conditions, transportation and storage

	Stent-graft length in uninflated condition	NP	Nominal pressure
	Stent-graft diameter in expanded condition		See Instructions for use
	Recommended guide wire		Do not use if packaging is damaged
	Catheter length	RBP	Rated burst pressure

	Diameters of the distal and proximal ends of the catheters		Contents
	Use by ...		Do not reuse
	Length of the stent-graft in expanded condition		Keep away from sunlight
	Recommended introducer sheath		Keep dry
	Reference number		MRI compatible
	Lot code		Does not contain natural latex
	Manufacturer		Non- pyrogenous
	For physicians only		Ethylene oxide sterilization
	Compliant with European Directive 93/42/EEC.		Do not re-sterilize

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(01) 0 000000 00000 0 (17) 150724 (10) XXXXXXXX

Operation

Stent-grafts are resistant to climatic factors under temperatures within the range of + 32 °C to + 42 °C and exposure to body fluids.

Transportation

Stent-grafts can be transported by all modes of transport, in covered vehicles, in accordance with the requirements and regulations for transportation of goods applicable to each mode of transport. In transit, stent-grafts are resistant to climatic factors under temperatures from -50 °C to + 50 °C, relative humidity from 60% to 95% (non-condensing humidity) and mechanical vibration loads within the frequency range from 10 to 55 Hz, displacement amplitude of 0.35 mm, and shock loads with peak acceleration of 10g at shock load lasting for not more than 16 ms.

Storage

In storage, stent-grafts are resistant to climatic factors under temperatures ranging from -50 to + 50 °C and relative humidity of 60% to 95% (non-condensing humidity). Stent-grafts should be stored in places protected from direct sunlight, heat and frost.

14. Packaging

All devices are supplied in a Tyvek protective sealed primary packaging.

Then, in the device in the primary packaging is placed onto a tray made of high density polyethylene and sealed with polythene film.

After that, the tray with finished products is placed into the bleached corrugated cardboard box, 1 in each, with Instructions for Use and the Patient Implant Card.

For transport purposes, boxes are placed inside large boxes made of thick corrugated board, 20 in each.

Overall package dimensions: (L x H x W) 320 x 268 x 14 mm ($\pm 10\%$)

Weight: 400 g ($\pm 10\%$)

15. Warranty

ClearStream Technologies, Ltd during one year from the date of the original purchase warrants to the primary purchaser of this product that this product free from material and workmanship defects. Liability under this limited product warranty is limited to repair or replacement of a defective product or full refunding of the purchase price at ClearStream Technologies, Ltd.'s discretion. This limited warranty does not cover normal wear and tear, as well as faults arising from use of the product for other purposes.

To the extent permitted by law, this limited product warranty excludes all other warranties, express or implied, in particular, any implied warranties of merchantability or fitness for a particular purpose. Under no circumstances ClearStream Technologies, Ltd. shall not be liable for any indirect, incidental or consequential damages arising out of this product handling or use.

In some countries, exclusion of implied warranties, incidental or consequential damages is not applicable. Pursuant to such countries laws, the user can be provided with additional remedies.

16. Shelf life

Shelf life of the medical device is 1 year. This medical device should not be used after the USE BY date indicated on the package.

17. Procedure for the medical device disposal or destruction

After use, the catheter sheath is a potential biological hazard. Handle and dispose of the catheter sheath in accordance with accepted medical practice and applicable local, state and federal laws and regulations

**I hereby certify accuracy, correctness and
reliability of this document
translated into Russian**

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