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The device design ensures parallel reduction of the jaws, which guarantees the high-quality formation of staples of the same height along the entire staple line length.

Product Codes Table

WARNING: The blue (standard) and green (thick) Echelon Contour Curved Cutter Reloads may only be used with The Echelon Contour Curved Cutter but can be used interchangeably.

Device codes:

Instrument Code	Reload Code	Tissue Thickness	Reload Color	Number of staples	Closed Staple Height
GCS40B	GCR40B	Standard	Blue	46	1.5 mm
GCS40G	GCR40G	Thick	Green	46	2.0 mm

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INSTRUCTIONS FOR USE

Verify compatibility of all instruments and accessories prior to using the instrument.

1. Determine the proper reload. Tissue thickness should be carefully evaluated prior to using the instrument and reload. Refer to the Product Codes Table for proper reload selection.

WARNING: Use of an inappropriate reload may result in improper staple formation.

2. Using sterile technique, remove the instrument with reload from the package. To avoid damage, do not flip the instrument or reload into the sterile field.

3. Examine the reload for the presence of a staple retainer (Illustration 1-12). The staple retainer ensures proper reload alignment and protects the staple leg points during shipping and transportation. If the retainer is not in place, discard the reload.

WARNING: The staple retainer ensures proper reload alignment and protects the staple leg points during shipping and transportation. If the reload is removed from the device without the staple retainer in place, the reload cannot be used and should be discarded. Failure to do so may result in malformed staples, incomplete cut line, bleeding, leakage from the staple line, and/or difficulty removing the device. Instrument damage or failure may result.

4. Remove the staple retainer by pushing on the pad labeled “Push to Remove” (Illustration 2). The loaded device cannot be used until the staple retainer is removed.

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Caution: After removing the staple retainer, observe the surface of the reload. The reload must be replaced with another reload if any colored drivers are visible.

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If colored drivers are visible, the reload may not contain staples. Attempting to close the device when the colored drivers are visible will initiate the lockout feature and prevent closure.

5. Check that the reload is held securely within the jaws. To ensure the reload is fully seated, the user may push on the blue or green cartridge (Illustration 1-10) where indicated by the arrow in Illustration 10. The user may notice audible and/or tactile feedback when the reload is correctly seated.

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6. Once the device is ready for use, grasp the instrument handle (Illustration 1-6) in the palm of the hand with the fingers around the closure trigger (Illustration 1-7). Position tissue to be cut and stapled between the reload anvil in the yellow washer (Illustration 1-11) and the green or blue cartridge (Illustration 1-10).

7. To capture tissue, manually advance the retaining pin (Illustration 1-13) using the retaining pin actuator (Illustration 1-4) on top of the handle (Illustration 3a). If the retaining pin was not previously manually advanced, it will automatically advance upon closure (Illustration 3b).

WARNING: Tissue that is too wide or thick to fit comfortably in the jaw aperture may indicate there is too much tissue to transect in a single firing. Attempting to force the trigger to complete the closing stroke with too much tissue or thickened tissue may result in poor staple line integrity with possible leakage or disruption. In addition, instrument damage or failure may result. It may be appropriate to complete the transection with more than one firing. See below for Instructions for Multiple Firings.

WARNING: Make sure tissue to be stapled is properly positioned in the jaws before stapling. Bunching, stretching, or uneven loading of tissue could result in leakage, lack of hemostasis, or disruption of staple line.

WARNING: If the retaining pin inadvertently pierces an unintended area of tissue during closure, this may result in a leak or trauma at the site of piercing. Be sure to carefully reposition the tissue into the appropriate area.

8. Squeeze the closure trigger and handle together until the closure trigger is latched. When the closure trigger is latched, the firing trigger moves to the ready-

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to-fire position. The reload has clamped onto the tissue which is ready to be cut and stapled. Ensure that the closing stroke is completed. Listen for an audible click (Illustration 4). All fingers should be completely removed from the closing trigger to ensure that the closing system holds.

9. An audible click halfway through the closure stroke indicates that the device is in the intermediate position. In the intermediate position, the pin is fully seated in the reload anvil in the yellow washer capturing the tissue, and the jaws are partially open. Reposition tissue within the instrument if desired. The closure trigger can be released at this point, and the instrument will maintain its intermediate jaw opening to allow for final positioning of the tissue to be cut and stapled. The retaining pin is fully engaged.

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WARNING: Tissue that is difficult to close in the jaws of the device may indicate there is too much tissue to transect in a single firing. Attempting to force the trigger to complete the closing stroke with too much tissue or thickened tissue may result in poor staple line integrity with possible leakage or disruption. In addition, instrument damage or failure may result. It may be appropriate to complete the transection with more than one firing. See below for "Instructions for Multiple Firings."

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If the tissue needs to be repositioned within the instrument before stapling, open the jaws by squeezing the closure trigger and pushing the release button with the thumb. The closure trigger will return to the fully open position and the jaws will release the tissue, but the retaining pin will remain extended. Tissue can now be repositioned.

10. After closing, wait 15 seconds to allow for adequate tissue compression.

11. With the handle and closure trigger grasped in the palm of the hand, place all fingers around the firing trigger. Before firing, check to ensure the retaining pin is seated in the reload anvil in the yellow washer. Always ensure complete closing prior to firing.

WARNING: Before firing, verify that the intended tissue is within the closed jaws of the instrument. If the pin is not properly positioned, staples may not form properly, which may result in leakage or disruption of the staple line.

WARNING: When positioning the stapler on the application site, ensure that no obstructions such as clips, stents, or guide wires are within the instrument jaws. Firing over an obstruction may result in incomplete cutting action and/or improperly formed staples.

12. Completely fire the instrument by squeezing the firing trigger until it touches the closure trigger (plastic to plastic), signifying that the instrument has fired the staples and knife blade, and the tissue has been transected (Illustration 5). The firing trigger automatically returns to the intermediate position as soon as it is released, leaving the closure trigger in the fully closed position.

WARNING: Make sure that the release button is not pressed during firing as proper formation of staples may be compromised.

WARNING: Care should be taken not to fire this instrument when the resulting staple line will be in tension as this may result in leakage at the staple line.

WARNING: The firing stroke must be completed. Do not partially fire the instrument. Incomplete firing can result in malformed staples, incomplete cut line, bleeding, and leakage from the staple line, and/or difficulty removing the device.

13. After firing the instrument, squeeze the closure trigger and push the release button with the thumb. To open the jaws of the instrument, release the closure trigger while simultaneously holding the release button (Illustration 6). The

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reload drivers are exposed to indicate that the reload has been used. If the instrument becomes locked onto tissue and will not open by pressing the release button, the device can be opened by depressing the release button and pulling the closure trigger simultaneously.

Caution: Examine the staple line for hemostasis and proper staple formation. If hemostasis is not present, appropriate techniques should be used to achieve hemostasis.

14. Retract the retaining pin manually (Illustration 7). Ensure that the staple line and the cut line are complete. Remove the instrument.

Instructions for Multiple Firings

1. If it is determined two firings will be needed to transect the tissue, for the first firing, ensure the device reaches more than half way across the tissue that is to be transected.

2. To capture tissue, manually advance the retaining pin (Illustration 1-13) using the retaining pin actuator (Illustration 1-4) on top of the handle (Illustration 3a) until it pierces the tissue. If the retaining pin was not previously manually advanced, it will automatically advance upon closure (Illustration 3b).

3. Squeeze the closure trigger and handle together until the closure trigger is latched. When the closure trigger is latched, the firing trigger moves to the ready-to-fire position. The reload has clamped onto the tissue which is ready to be cut and stapled. Ensure the closing stroke is completed. Listen for an audible click (Illustration 4). All fingers should be completely removed from the closing trigger to ensure the closing system holds.

WARNING: Do not grasp the firing trigger (Illustration 1-8) before the instrument is to be fired. Doing so may result in malformed staples, incomplete cut line, bleeding, and leakage from the staple line.

Caution: Do not grasp the firing trigger before the instrument is to be fired, doing so could initiate the lockout feature; in this case the device will need to be released from the tissue. The instrument may be partially fired and will need to be reloaded before using on tissue.

WARNING: Attempting to force the trigger to complete the closing stroke with too much tissue or thickened tissue may result in poor staple line integrity with

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possible leakage or disruption. In addition, instrument damage or failure may result.

4. If the tissue needs to be repositioned within the instrument before stapling, open the jaws by squeezing the handle and closure trigger together slightly, push the release button (Illustration 1-5), and slowly release grasp of the closure trigger (Illustration 6). The closure trigger will return to the fully open position and the jaws will release the tissue, but the retaining pin will remain extended. Tissue can now be repositioned. If device is repositioned, ensure the pierced tissue will be removed with the specimen.

Caution: The Echelon Contour Curved Cutter is designed with a feature which prevents closure if a used reload, no reload, or an incorrectly inserted reload is in the instrument. The used or misloaded reload module should then be replaced with a new reload or, if no reload was present, a reload should be loaded in the instrument. After following the above steps, any device which does not close either partially or completely should not be used.

5. After closing, wait 15 seconds to allow for adequate tissue compression.

6. With the handle and closure trigger grasped in the palm of the hand, place all fingers around the firing trigger. Before firing, check to ensure the retaining pin is seated in the reload anvil in the yellow washer. Always ensure complete closing prior to firing.

WARNING: Before firing, verify that the intended tissue is in the closed jaws of the instrument. If the pin is not properly positioned, staples may not form properly, which may result in leakage or disruption of the staple line.

WARNING: When positioning the stapler on the application site, ensure that no hard objects such as clips, stents, or guide wires are within the instrument jaws. Firing over an obstruction may result in incomplete cutting action and/or improperly formed staples.

7. Completely fire the instrument by squeezing the firing trigger until it touches the closure trigger (plastic to plastic), signifying that the instrument has fired the staples and knife blade, and the tissue has been transected (Illustration 5). The firing trigger automatically returns to the intermediate position as soon as it is released, leaving the closure trigger in the fully closed position. Make sure that the release button is not pressed during firing as proper formation of staples may be compromised.

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WARNING: The firing stroke must be completed. Do not partially fire the instrument. Incomplete firing can result in malformed staples, incomplete cut line, bleeding, and leakage

from the staple line, and/or difficulty removing the device.

8. After firing the instrument, squeeze the closure trigger and push the release button with the thumb. To open the jaws of the instrument, release the closure trigger while simultaneously holding the release button (Illustration 6). The reload drivers are exposed to indicate that the reload has been used. If the instrument becomes locked onto tissue and will not open by pressing the release button, the device can be opened by depressing the release button and simultaneously pulling the closure trigger open.

9. Examine the staple line for proper staple formation.

10. Following the first firing, the tissue is transected up to the retaining pin and the staple line extends beyond the pin. The staple line extension temporarily seals against the potential of lumen content until the transection can be completed with a second firing. At this time, the device should be reloaded for the second firing.

11. See "Instructions for Reload" below for reload removal, device cleaning, reload selection, and installation.

12. When positioning the instrument for the second firing, the second staple line may be positioned where the first line ended (Illustration 11a), or it may be offset to control the location of the staple line crossing (Illustration 11b).

13. Advance the retaining pin to capture the remaining tissue within the jaw opening. See "Instructions for Use" above for device positioning, closure, and firing. After clamping, to ensure adequate tissue compression, wait 15 seconds prior to firing.

14. Release the device and examine the tissue for hemostasis and proper staple formation. If hemostasis is not present, appropriate techniques should be used to achieve hemostasis.

Retract the retaining pin manually (Illustration 7). Ensure that the staple line and the cut line are complete. Remove the instrument.

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Instructions for Reload

1. Remove the used reload from the instrument. Grasp the top of the reload near the grasping feature and pull out in a slight arcing motion, unsnapping the reload from the jaws (Illustration 8). Properly discard the used reload.

2. Thoroughly clean the instrument by completely submerging the curved head (Illustration 1-1)

and reload housing (Illustration 1-2) in a sterile solution and swish vigorously. After swishing, visually inspect the curved head and reload housing and remove all residual staples and/or foreign matter from instrument prior to reloading.

WARNING: Failure to properly clean the instrument prior to reloading could damage the instrument as well as compromise cut line and staple formation.

3. Determine the proper reload. Tissue thickness should be carefully evaluated prior to using the instrument and reload. Refer to the Product Codes Table for proper reload selection.

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4. Using sterile technique, remove the reload from the package. To avoid damage, do not flip the reload into the sterile field.

5. Examine the new reload for the presence of a staple retainer. If the retainer is not in place, discard the reload. Do not remove the staple retainer until the reload is inserted into the device. The staple retainer cannot be removed until the reload is correctly loaded into the device. The loaded device cannot be used until the staple retainer is removed.

6. Depress the release button to ensure the instrument is in the open position. Ensure that the retaining pin actuator (Illustration 1-4) has been retracted. Load the device with the appropriate reload by aligning the reload with the device. The tracks on each side of the reload should be used as guides to align the reload within the jaws of the instrument (Illustration 9). When the reload (staple retainer still in place) is properly aligned, push the reload into the instrument until it is fully seated. Check that the reload is held securely within the jaws. The user may notice audible and/or tactile feedback when the Echelon Contour Curved Cutter Reload is correctly seated (Illustration 10).

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WARNING: Do not remove the yellow washer (Illustration 1-11) as it is critical to the structural integrity of the reload. Removing the washer may result in malformed staples, incomplete cut line, bleeding, and leakage from the staple line.

WARNING: The staple retainer ensures proper reload alignment and protects the staple leg points during shipping and transportation. If the reload is removed from the device without the staple retainer in place, the reload cannot be used and should be discarded. Failure to do so may result in malformed staples, incomplete cut line, bleeding, and leakage from the staple line, and/or difficulty removing the device. In addition, instrument damage or failure may result.

7. Remove the staple retainer by pushing on the pad labeled “Push to Remove” (Illustration 2) and discard. The loaded device cannot be used until the staple retainer is removed.

WARNING: Do not remove the yellow washer (Illustration 1-11) as it is critical to the structural integrity of the reload. Removing the washer may result in malformed staples, incomplete cut line, bleeding, and leakage from the staple line.

Caution: After removing the staple retainer, observe the surface of the new reload. The reload must be replaced with another reload if any colored drivers are visible. If colored drivers are visible, the reload may not contain staples. Attempting to close the device when the colored drivers are visible will initiate the lockout feature and prevent closure.

8. Check that the reload is held securely within the jaws. To ensure the reload is fully seated, the user may push on the blue or green cartridge (Illustration 1-10) where indicated by the arrow in Illustration 10. The user may notice audible and/or tactile feedback when the reload is correctly seated.

Caution: Attempting to close the device with an improperly seated reload will initiate the lockout feature and prevent closure.

9. The instrument is now reloaded and ready to use.

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WARNINGS AND PRECAUTIONS

- Minimally invasive procedures should be performed only by persons having adequate training and familiarity with minimally invasive techniques. Consult medical literature relative to techniques, complications, and hazards prior to performance of any minimally invasive procedure.
- Minimally invasive instruments may vary in diameter from manufacturer to manufacturer. When minimally invasive instruments and accessories from different manufacturers are employed together in a procedure, verify compatibility prior to initiation of the procedure.
- A thorough understanding of the principles and techniques involved in laser, electrosurgical, and ultrasonic procedures is essential to avoid shock and burn hazards to both patient and medical personnel and damage to the device or other medical instruments. Ensure that electrical insulation or grounding is not compromised. Do not immerse electrosurgical instruments in liquid unless the instruments are designed and labeled to be immersed.
- Pre-operative radiotherapy may result in changes to tissue characteristics. These changes may, for example, cause the tissue thickness to exceed the indicated range for the selected staple. Careful consideration should be given to any pre-surgical treatment the patient may have undergone, which may require alterations to surgical technique or alternative surgical procedures.
- The blue (standard) and green (thick) Echelon Contour™ Curved Cutter Reloads may only be used with The Echelon Contour™ Curved Cutter but can be used interchangeably.
- Do not fire device without a loaded reload.
- Use of an inappropriate reload may result in improper staple formation.
- Attempting to load the device without a staple retainer could result in misalignment, creating malformed staples, incomplete cut line, bleeding, and leakage from the staple line, and/or difficulty removing the device. In addition, instrument damage or failure may result.
- The staple retainer ensures proper reload alignment and protects the staple leg points during shipping and transportation. If the reload is removed from the device without the staple retainer in place, the reload cannot be used and should be discarded. Failure to do so may result in malformed staples, incomplete cut line, bleeding, and leakage from the staple line, and/or difficulty removing the device. In addition, instrument damage or failure may result.

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- Do not remove the yellow washer as it is critical to the structural integrity of the reload. Removing the washer may result in malformed staples, incomplete cut line, bleeding, and leakage from the staple line.

- After removing the staple retainer, examine the surface of the new reload. The reload must be replaced with another reload if any colored drivers are visible. If colored drivers are visible, the reload may not contain staples. Attempting to close the device when the colored drivers are visible will initiate the lockout feature and prevent closure.

- Make sure tissue to be stapled is properly positioned in the jaws before stapling. Bunching, stretching, or uneven loading of tissue could result in leakage, lack of hemostasis, or disruption of staple line.

- Ensure that unintentional tissue is not cut.

- Care should be taken not to fire this instrument when the resulting staple line will be in tension.

- Tissue that is too wide or thick to fit comfortably in the jaw aperture may indicate there is too much tissue to transect in a single firing. Attempting to force the trigger to complete the closing stroke with too much tissue or thickened tissue may result in poor staple line integrity with possible leakage or disruption. In addition, instrument damage or failure may result. It may be appropriate to complete the transection with more than one firing.

- Tissue that is difficult to close in the jaws of the device may indicate there is too much tissue to transect in a single firing. Attempting to force the trigger to complete the closing stroke with too much tissue or thickened tissue may result in poor staple line integrity with possible leakage or disruption. In addition, instrument damage or failure may result. It may be appropriate to complete the transection with more than one firing.

- Before firing, verify that the intended tissue is in the closed jaws of the instrument. If the pin is not properly positioned, staples may not form properly, which may result in leakage or disruption of the staple line.

- When positioning the stapler on the application site, ensure no hard objects such as clips, stents, or guide wires are within the instrument jaws. Firing over an obstruction may result in incomplete cutting action and/or improperly formed staples.

- The firing stroke must be completed. Do not partially fire the instrument. Incomplete firing can result in malformed staples, incomplete cut line, bleeding, and leakage from the staple line, and/or difficulty removing the device.

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- Examine the staple line for hemostasis and proper staple formation. If hemostasis is not present, appropriate techniques should be used to achieve hemostasis.

- Do not grasp the firing trigger (Illustration 1-8) before the instrument is to be fired. Doing so may result in malformed staples, incomplete cut line, bleeding, and leakage from the staple line.

- If the retaining pin inadvertently pierces an unintended area of tissue during closure, this may result in a leak or trauma at the site of piercing. Be sure to carefully reposition the tissue into the appropriate area.

- Do not grasp the firing trigger before the instrument is to be fired. Doing so could initiate the lockout feature. In this case the device will need to be released from the tissue. The instrument may be partially fired and will need to be reloaded before using on tissue.

- Attempting to force the trigger to complete the closing stroke with too much tissue or thickened tissue may result in poor staple line integrity with possible leakage or disruption. In addition, instrument damage or failure may result.

- The Echelon Contour Curved Cutter is designed with a feature which prevents closure if a used reload, no reload, or an incorrectly inserted reload is in the instrument. The used or misloaded reload should then be replaced with a new reload or, if no reload was present, a reload should be loaded in the instrument. After following the above steps, any device which does not close either partially or completely should not be used.

- Failure to properly clean the instrument prior to reloading could damage the instrument and compromise cut line and staple formation of subsequent firings.

- When using other technologies in the procedure, observe the precautions suggested by the original equipment manufacturer to avoid the hazards associated with their use.

- Do not use the Echelon Contour Curved Cutter or Reloads after the expiration date or if the package is damaged. This may render the equipment inoperable or non-sterile.

- Dispose of all opened instruments whether used or unused. This device is packaged and sterilized for single use only.

- Reuse and improper reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure that in turn may result in patient injury, illness, or death.

- Reuse and improper reprocessing or resterilization of single-use devices may create a risk of contamination and/or cause infection or cross-infection, including,

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but not limited to, the transmission of infectious diseases. Contamination may lead to injury, illness, or death.

- Instruments or devices which come into contact with bodily fluids may require special disposal handling to prevent biological contamination. Disposal of all medical device instruments and accessories should be performed according to local requirements and regulations.
- Attempting to close the device with an improperly seated reload will initiate the lockout feature and prevent closure.
- Discard used items containing knives or blades in "sharps" containers.

MRI Safety Information

Non-clinical testing has demonstrated the staples made of titanium (Ti3Al2.5V) alloy are MR Conditional. A patient with the implanted staples can be safely scanned in an MR system meeting the following conditions:

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- Maximum spatial field gradient of 3,000 gauss/cm (30 T/m)
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2 W/kg (Normal Operating Mode) Under the scan conditions defined above, the staples are expected to produce a maximum temperature rise of less than 2° C after 15 minutes of continuous scanning.

In non-clinical testing, the image artifact caused by the device extends approximately 4 mm from the staple line when imaged with a gradient echo pulse sequence and a 3.0 T MRI system.

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TRANSPORT AND STORAGE

Transportation conditions:

Designation	Range
Temperature	Maximum + 60°C
Relative humidity	Less than 60 %
Barometric pressure	No specific conditions

Storage conditions:

Designation	Range
Temperature	Maximum + 25°C
Relative humidity	Less than 60 %
Barometric pressure	No specific conditions

Keep Dry.

Keep Away from Heat.

OPERATING ENVIRONMENT CONDITIONS

Hospital, outpatient surgical center.

Operation conditions	
Humidity	20 % to 60 %
Temperature	20 °C to 24 °C
Barometric pressure	No specific conditions

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LIST OF NATIONAL STANDARDS APPLIED BY THE MANUFACTURER OF A MEDICAL DEVICE

STANDARDS	DESCRIPTION
EN ISO 13485	Medical devices - Quality management systems - Requirements for regulatory purposes
EN ISO 14971	Medical devices - Application of risk management to medical devices
EN ISO 15223-1	Medical devices - Symbols to be used with medical device labels, labeling and information to be supplied -Part 1: General requirements
EN 1041	Information provided with medical devices
EN ISO 20417	Medical devices Information to be supplied by the manufacturer
EN ISO 16061	Instrumentation for use in association with non-active surgical implants - General requirements
EN ISO 14630	Non-active surgical implants - General requirements
EN 62570	Standard practice for marking medical devices and other items for safety in the magnetic resonance environment
EN 62366	Medical Devices- Application of usability engineering to medical devices
EN ISO 10993-1	Biological evaluation of medical devices -Part 1: Evaluation and testing within a risk management process
EN ISO 10993-2	Biological evaluation of medical devices - Part 2: Animal welfare requirements

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