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CERTIFICATE

中国国际贸易促进委员会暨中国国际商会
China Council for the Promotion of International Trade is China Chamber of International Commerce

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China Council for the Promotion of International Trade
China Chamber of International Commerce

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THIS IS TO CERTIFY THAT: the seal of Shanghai MicroPort EP MedTech Co., Ltd. on the annexed INSTRUCTIONS FOR USE is genuine.

China Council for the Promotion
of International Trade CCPIT

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证明书查询网址 Website for verifying the certificate: <http://www.rzccpit.com/validate.html>

Shanghai MicroPort EP MedTech Co., Ltd.
CEO

Sun, Yi Yong

(Stamp, Signature)

Date:

Feb 6, 2023

INSTRUCTIONS FOR USE

for Medical Device

«PathBuilder transseptal access kit, in versions»

Official manufacturer: Shanghai MicroPort EP MedTech Co. Ltd., Building 23&28, Lane 588,
Tianxiong Rd., 201318 Shanghai, People's Republic of China

Instructions for use

PathBuilder transseptal access kit with Transseptal Guiding Introducer

Shanghai MicroPort EP MedTech Co., Ltd.

【WARNING】

- Do not reuse this device. The biological materials and foreign bodies attached cannot be removed thoroughly after use, so any reuse of this product may lead to adverse effect of patients.
- The "Use By" date on the device package and the sterile packaging should be inspected prior to use. Do not use the device if past the "Use By" date or the packaging appears damaged.
- Do not alter this device in any way.
- Do not use excessive force to advance or withdraw the device when resistance is encountered.
- Procedures present the potential for significant x-ray exposure, which can result in acute radiation injury as well as increased risk for somatic and genetic effects, to both patients and physician due to the x-ray beam intensity and duration of the fluoroscopic imaging. Procedures should only be performed after adequate attention has been given to the potential radiation exposure associated with the procedure, and steps taken to minimize this exposure.
- When used, device shall be disposed of as medical waste according local laws and regulations.
- Only those physicians who are trained in transseptal procedures and catheter delivery systems should use this device.
- If the patient with a history of heparin induced thrombocytopenia (HIT), bivalirudin is a preferred alternative to heparin for anticoagulation during procedures.

【Packing List】

PathBuilder transseptal access kit, in versions:

1. PathBuilder transseptal access kit with Transseptal Guiding Introducer, composed of:

1. Transseptal Guiding Introducer PathBuilder, in versions: diameter 8F(2.67mm); 8.5F(2.83mm), length: 630 mm; curve configuration: L0; L1; L2; L3; L4; R0; R1; R2; R3; R4 – 1pc.;
2. Dilator in versions: length: 670mm – 1pc.;
3. Guidewire, length 180cm, diameter 0.81mm – 1pc.;
4. Transseptal needle PathBuilder 18G(1.21mm), in versions: M, L, length 710mm – 1pc. (if necessary);
5. Stylet – 1pc. (if necessary);
6. Instructions for use

2. PathBuilder transseptal access kit with Transseptal Guiding Introducer, composed of:

1. Transseptal Guiding Introducer PathBuilder, in versions: diameter 8F(2.67mm); 8.5F(2.83mm), length: 810mm; curve configuration: L0; L1; L2; L3; L4; R0; R1; R2; R3; R4 – 1pc.;
2. Dilator in versions: length: 850mm – 1pc.;
3. Guidewire, length 180cm, diameter 0.81mm – 1pc.;
4. Transseptal needle PathBuilder 18G(1.21mm) in variants: M, L, length 890mm – 1pc. (if necessary);
5. Stylet – 1pc. (if necessary);
6. Instructions for use

【Overview】

Pathbuilder Transseptal Guiding Introducer manufactured by Shanghai MicroPort EP MedTech Co., Ltd. is designed to introduce various cardiovascular catheters into each heart chamber for further electrophysiology operation. Different pathways can be built according to various procedure needs to introduce the electrophysiology catheters into different part of heart chambers.

Pathbuilder Transseptal Guiding Introducer is comprised of a sheath, a dilator and a guidewire. On the tip of sheath is a marker band to observe the position of catheter at any time; on the sheath is provided with valves to provide hemostasis during dilator/catheter exchange dilator; on its sideport is provided with three-way stopcock that can be used for aspiration, fluid infusion, blood sampling and pressure monitoring. The dilator inner lumen is also designed with a special geometry as its distal end to limit the exposure of needle.

【PRODUCT SPECIFICATION】

Pathbuilder Transseptal Guiding Introducer can be classified into different specifications based on the curve, diameter as well as length of the body. Different color codes are attached near the sideport of sheath to distinguish tube diameter.

Pathbuilder Transseptal Guiding Introducer has two specifications of 8F(2.67mm) and 8.5F(2.83mm) according to its diameter, which are identified by the white or blue color code attached near the sideport of sheath (See Fig.1,2.); each specification has ten models of L0, L1, L2, L3, L4, R0, R1, R2, R3 and R4 according to its curve shape; in each model, there are two specifications of length that are 63cm and 81cm. Physicians can make choices based on actual needs.

Table 1: Specification table of Pathbuilder Transseptal Guiding Introducer

Model	Curve	Radius of curve (mm)	Color code	Sheath		Dilator		Guidewire	
				Length (cm/mm)	Diameter (F/mm)	Length (cm/mm)	Distal end	Length (cm)	Outer diameter (inch/mm)
TS L0 63 8	L0	25	White	63/630	8/2.67	67/670	Curved	180	0.032/0.81
TS L1 63 8	L1	25		63/630	8/2.67	67/670	Curved	180	0.032/0.81
TS L2 63 8	L2	25		63/630	8/2.67	67/670	Curved	180	0.032/0.81
TS L3 63 8	L3	25		63/630	8/2.67	67/670	Curved	180	0.032/0.81
TS L4 63 8	L4	25		63/630	8/2.67	67/670	Curved	180	0.032/0.81
TS L0 63 85	L0	25	Blue	63/630	8.5/2.83	67/670	Curved	180	0.032/0.81
TS L1 63 85	L1	25		63/630	8.5/2.83	67/670	Curved	180	0.032/0.81

TS L2 63 85	L2	25		63/630	8.5/2.83	67/670	Curved	180	
TS L3 63 85	L3	25		63/630	8.5/2.83	67/670	Curved	180	0.032/0.81
TS L4 63 85	L4	25		63/630	8.5/2.83	67/670	Curved	180	0.032/0.81
TS L0 81 8	L0	25	White	81/810	8/2.67	85/850	Curved	180	0.032/0.81
TS L1 81 8	L1	25		81/810	8/2.67	85/850	Curved	180	0.032/0.81
TS L2 81 8	L2	25		81/810	8/2.67	85/850	Curved	180	0.032/0.81
TS L3 81 8	L3	25		81/810	8/2.67	85/850	Curved	180	0.032/0.81
TS L4 81 8	L4	25		81/810	8/2.67	85/850	Curved	180	0.032/0.81
TS L0 81 85	L0	25		81/810	8/2.67	85/850	Curved	180	0.032/0.81
TS L1 81 85	L1	25	Blue	81/810	8.5/2.83	85/850	Curved	180	0.032/0.81
TS L2 81 85	L2	25		81/810	8.5/2.83	85/850	Curved	180	0.032/0.81
TS L3 81 85	L3	25		81/810	8.5/2.83	85/850	Curved	180	0.032/0.81
TS L4 81 85	L4	25		81/810	8.5/2.83	85/850	Curved	180	0.032/0.81
TS R0 63 8	R0	25		81/810	8.5/2.83	85/850	Curved	180	0.032/0.81
TS R0 63 8	R0	25	White	63/630	8/2.67	67/670	Straight	180	0.032/0.81
TS R1 63 8	R1	25		63/630	8/2.67	67/670	Straight	180	0.032/0.81
TS R2 63 8	R2	25		63/630	8/2.67	67/670	Straight	180	0.032/0.81
TS R3 63 8	R3	25		63/630	8/2.67	67/670	Straight	180	0.032/0.81
TS R4 63 8	R4	25		63/630	8/2.67	67/670	Straight	180	0.032/0.81
TS R0 63 85	R0	25	Blue	63/630	8.5/2.83	67/670	Straight	180	0.032/0.81
TS R1 63 85	R1	25		63/630	8.5/2.83	67/670	Straight	180	0.032/0.81
TS R2 63 85	R2	25		63/630	8.5/2.83	67/670	Straight	180	0.032/0.81
TS R3 63 85	R3	25		63/630	8.5/2.83	67/670	Straight	180	0.032/0.81
TS R4 63 85	R4	25		63/630	8.5/2.83	67/670	Straight	180	0.032/0.81
TS R0 81 8	R0	25	White	81/810	8/2.67	85/850	Straight	180	0.032/0.81
TS R1 81 8	R1	25		81/810	8/2.67	85/850	Straight	180	0.032/0.81
TS R2 81 8	R2	25		81/810	8/2.67	85/850	Straight	180	0.032/0.81
TS R3 81 8	R3	25		81/810	8/2.67	85/850	Straight	180	0.032/0.81
TS R4 81 8	R4	25		81/810	8/2.67	85/850	Straight	180	0.032/0.81
TS R0 81 85	R0	25	Blue	81/810	8.5/2.83	85/850	Straight	180	0.032/0.81
TS R1 81 85	R1	25		81/810	8.5/2.83	85/850	Straight	180	0.032/0.81
TS R2 81 85	R2	25		81/810	8.5/2.83	85/850	Straight	180	0.032/0.81
TS R3 81 85	R3	25		81/810	8.5/2.83	85/850	Straight	180	0.032/0.81
TS R4 81 85	R4	25		81/810	8.5/2.83	85/850	Straight	180	0.032/0.81


Fig. 1: 8F(2.67mm) white color code attached near the sideport of sheath


Fig. 2: 8.5F(2.83mm) blue color code attached near the sideport of sheath

【INDICATIONS FOR USE】

PathBuilder Transseptal Guiding Introducer is indicated for introducing various cardiovascular catheters into the heart via femoral vein, including the left side of the heart through the interatrial septum, and when used with cardiovascular catheters, for assisting catheters in reaching the target site.

【INDICATIONS】

Atrial fibrillation

【CONTRAINDICATIONS】

- 1) Intra-atrial septal patch;
- 2) Atrial myxoma;
- 3) Myocardial infarction within the last two weeks;
- 4) Unstable angina;
- 5) Cerebral vascular accident (CVA) within the last two weeks;
- 6) Patients who do not tolerate anticoagulation therapy;
- 7) Patients with an active infection;
- 8) Left atrial thrombus which is presence after more than three weeks of anticoagulation

【TARGET GROUP】

18-75 years old adults, men or non-pregnant women

【Risks and Side Effects】

Complications that may occur during use of this device include scratches of endocardium and vessel wall cause reactions like inflammation or abscess; toxic reaction; complications caused by vascular puncture (such as pseudoaneurysm, arteriovenous fistula, local hematoma, or ecchymosis and infection); radiation injury; thrombosis / air embolism; injury to coronary artery; serious coagulation reactions; pericardial tamponade; cardiac / artery perforation; pericardial effusion; cerebral embolism / cerebral apoplexy / stroke; thoracotomy to remove foreign body.

【Accessories Required for Using Product】

- Transseptal needle (supplied in kit if necessary or to be purchased as PathBuilder transseptal access kit with Transseptal Needle)

PathBuilder Transseptal Guiding Introducer	PathBuilder Transseptal Needle		
	Model	Diameter	Length
TS L0 63 8	TNM7118	18Ga (1.21mm)	710mm
TS L1 63 8	TNL7118	18Ga (1.21mm)	710mm
TS L2 63 8		18Ga (1.21mm)	710mm
TS L3 63 8		18Ga (1.21mm)	710mm
TS L4 63 8		18Ga (1.21mm)	710mm
TS L0 63 85		18Ga (1.21mm)	710mm
TS L1 63 85		18Ga (1.21mm)	710mm
TS L2 63 85		18Ga (1.21mm)	710mm
TS L3 63 85		18Ga (1.21mm)	710mm

TS L4 63 85		18Ga (1.21mm)	710mm
TS L0 81 8	TNM8918	18Ga (1.21mm)	890mm
TS L1 81 8	TNL8918	18Ga (1.21mm)	890mm
TS L2 81 8		18Ga (1.21mm)	890mm
TS L3 81 8		18Ga (1.21mm)	890mm
TS L4 81 8		18Ga (1.21mm)	890mm
TS L0 81 85		18Ga (1.21mm)	890mm
TS L1 81 85		18Ga (1.21mm)	890mm
TS L2 81 85		18Ga (1.21mm)	890mm
TS L3 81 85		18Ga (1.21mm)	890mm
TS L4 81 85		18Ga (1.21mm)	890mm
TS R0 63 8	TNM7118	18Ga (1.21mm)	710mm
TS R1 63 8	TNL7118	18Ga (1.21mm)	710mm
TS R2 63 8		18Ga (1.21mm)	710mm
TS R3 63 8		18Ga (1.21mm)	710mm
TS R4 63 8		18Ga (1.21mm)	710mm
TS R0 63 85		18Ga (1.21mm)	710mm
TS R1 63 85		18Ga (1.21mm)	710mm
TS R2 63 85		18Ga (1.21mm)	710mm
TS R3 63 85		18Ga (1.21mm)	710mm
TS R4 63 85		18Ga (1.21mm)	710mm
TS R0 81 8	TNM8918	18Ga (1.21mm)	890mm
TS R1 81 8	TNL8918	18Ga (1.21mm)	890mm
TS R2 81 8		18Ga (1.21mm)	890mm
TS R3 81 8		18Ga (1.21mm)	890mm
TS R4 81 8		18Ga (1.21mm)	890mm
TS R0 81 85		18Ga (1.21mm)	890mm
TS R1 81 85		18Ga (1.21mm)	890mm
TS R2 81 85		18Ga (1.21mm)	890mm
TS R3 81 85		18Ga (1.21mm)	890mm
TS R4 81 85		18Ga (1.21mm)	890mm

- Cardiovascular catheter (to be purchased separately)
Ablation or diagnostic catheters (OD: 7F/2.33mm and 8F/2.67mm)
- Syringe (to be purchased separately)
Syringes with 6% (Luer) taper lock fitting

【DIRECTIONS FOR USE】

1. Prepare and assemble equipment

Preparing the equipment requires the following items:

- One transseptal sheath and dilator,
- One matching transseptal needle and a stainless steel stylet,
- One 0.032"(0.81mm), 180cm guidewire with a 3mm"J"tip,
- Syringes for aspiration and flushing,
- Sterile heparinized saline
- Flush the transseptal sheath and dilator with sterile heparinized saline and insert the dilator fully into the sheath.

- Remove the stylet from the needle and flush the needle with sterile heparinized saline. Re-insert the stylet into the needle and lock it onto the hub.
- Insert the needle assembly into the sheath/dilator. Due to the stop feature of the dilator, there will be a gap between the dilator hub and the needle pointer flange. (See Fig.3. Magnification is 2.9).


Fig. 3

- Withdraw the needle assembly until the tip of the stylet is just within the tip of the dilator.
2. **Advance sheath/dilator assembly into superior vena cava.**
 - Introduce the prepared guidewire from femoral vein (right fem oral preferred) to superior vena cava. Note: 0.032 inch(0.81mm) is the maximum guidewire diameter that can be used with the transseptal dilator.
 - Insert the transseptal sheath/dilator assembly over the guidewire and advance the assembly into the superior vena cava.
 3. **Position the needle/stylet assembly inside the sheath/dilator assembly.**
 - Remove the guidewire from the dilator
 - Aspirate and then flush the dilator with clean heparinized saline, ensuring that no air enters the bloodstream.
 - Separate the sheath and dilator by withdrawing the dilator by a distance sufficient to accommodate the needle curve. This will facilitate passage of the needle curve through the dilator and sheath hubs.
 - Insert the needle/stylet assembly into the dilator, letting it rotate freely as it advances.
 - After the needle curve is advanced beyond the hemostasis valve hub of the sheath, reconnect the sheath and dilator by sliding the sheath back over the dilator slowly.(do not advance the dilator)
 - Advance the needle/stylet assembly and confirm that the needle tip is inside the dilator by fluoroscopy.

- After removing the stylet, attach a syringe to the needle hub and aspirate until blood return is observed. Then discard the syringe.
- Fill the syringe and flush the needle with clean heparinized saline, ensuring that the syringe is attached to the needle hub and no air enters the bloodstream.

4. Engage the fossa ovalis

- Visualize and identify anatomic landmarks. Set the fluoroscopy unit to an appropriate angle, parallel to the plane of the mitral valve and orthogonal to the plane of the septum. This will typically be left anterior oblique (LAO).
- Adjust the needle pointer so that the needle is perpendicular to the fossa ovalis (typically between 3:00 and 5:00 o'clock, as viewed from the foot end of the patient). (See Fig.4. Magnification is 2.9).


Fig. 4

- Confirm that the needle tip is inside the dilator by fluoroscopy. Then drag the entire assembly slowly. Prevent any movement of the assembly parts relative to each other. Maintain the previous orientation of the needle pointer.
- Observe the tip of the dilator during the drag for medial (or rightward) movement, indicating the tip has engaged the fossa ovalis.

5. Puncture the fossa ovalis

- Once the correct location is confirmed, advance the needle across the interatrial septum.
- Under pressure monitoring, entry into the left atrium is confirmed when the pressure tracing shows a left atrial pressure waveform.

If there is any resistance to needle advancement, re-evaluate the anatomic landmarks.

Caution: if pericardial or aortic entry occurs, do not advance the dilator over the needle. If the needle has penetrated the pericardium or aorta, it must be withdrawn. Monitor vital signs closely.

6. Advance the sheath/dilator assembly

- While maintaining a fixed needle position, advance the sheath/dilator assembly over the needle.

aspirate until blood
ensuring that the

7. Advance the sheath over the fixed dilator and needle into the left atrium.

- Withdraw the needle into the dilator until it is just inside the tip. Maintain the position of the dilator across the septum.
- With the dilator in a fixed location, advance the sheath over the dilator.

8. Withdraw the needle and the dilator

- Disconnect any attachments to the needle hub
- Immediately aspirate blood through the side arm of the sheath after withdrawing the needle and dilator slowly at the same time. The blood should be arterial blood.
- Under pressure monitoring, the pressure tracing still shows a left atrial pressure waveform. The sheath is now in place in the left atrium.

9. Advance the catheter

- Insert electrode catheter into the sheath slowly and observe extension length of catheter by fluoroscopy.

10. Remove the sheath

- After removing the catheter, re-insert the guidewire over the sheath. Then re-insert the dilator over the guidewire and keep the sheath straight. At last remove the sheath, the dilator and the guidewire as a whole.

【STERILIZATION】

Pathbuilder™ Transseptal Guiding Introducer has been sterilized by epoxy ethane (ethylene oxide) before delivery.

【STORAGE & TRANSPORTATION REQUIREMENTS】

The product shall be stored in a shady and cool, dry, clean, and well-ventilated warehouse which is in natural air circulation environment. In transit, the product shall be protected from heavy load, direct sunlight and rain or as specified in the ordering contract.

Situation	Ambient temperature	Relative air humidity	Atmospheric pressure
Transportation	-40°C to 55°C	10% to 90%	50kPa to 106kPa
Storage	18°C to 28°C	10% to 90%	50kPa to 106kPa
Working	5°C to 40°C	10% to 80%	70kPa to 106kPa

【PERIOD OF VALIDITY】

The sterilization period of validity is three years under specified storage conditions.

【MANUFACTURING DATA】

See label statement for service life.

【DISPOSAL】

Dispose the devices according to local environmental requirements.

Environmental Monitoring

Class of hazard for medical wastes by reference to specification of morphological composition
(in accordance with Russian Sanitary Regulations and Standards 2.1.3684-21) Class B

【SYMBOLS DESCRIPTION】

	SYMBOL	DESCRIPTION
1.		DO NOT REUSE
2.		CONSULT INSTRUCTIONS FOR USE
3.		PROTECT FROM HEAT SOURCE AND RADIATION SOURCE
4.		KEEP DRY
5.		USE BY
6.		BATCH CODE
7.		CATALOGUE NUMBER
8.		STERILIZED USING ETHYLENE OXIDE
9.		DO NOT USE IF PACKAGE IS DAMAGED
10.		QUANTITY OF PRODUCT CONTAINED 1
11.		DATE OF MANUFACTURE
12.		MANUFACTURER
13.		AUTHORISED REPRESENTATIVE IN THE EUROPEAN COMMUNITY
14.		TEMPERATURE LIMITATION

【AFTER-SALES SERVICE】

With “providing the medical sector with high quality and efficacy medical products” as its top operational objective, Shanghai MicroPort EP MedTech Co., Ltd. (hereinafter referred to as MicroPort EP) guarantees that its products are free of defects in materials or manufacturing when the clients receive them. For other questions relating to the products, please directly consult the company.

【SOLEMN STATEMENT】

Shanghai MicroPort EP MedTech Co., Ltd. (hereinafter referred to as MicroPort EP Co.) expressly states herein that its Pathbuilder Transseptal Guiding Introducer is disposable and cannot be reused. MicroPort EP Co. will not recommend, indicate and imply in any manner the reusability of the system, and will not assume the responsibility for any accident or product damage resulting from reuse.

Pathbuilder Transseptal Guiding Introducer can be connected and used only with the compatible devices specified herein, and MicroPort EP Co. will not assume the responsibility for damage to product device, procedure failure and the like resulting from operating mistakes or any other human factors.

【GENERAL STANDARDS】

Shanghai MicroPort EP MedTech Co., Ltd. declares that PathBuilder transseptal access kit with Transseptal Guiding Introducer meets the requirements ISO 11070, EN ISO 13485, EN ISO 14971. A complete list of requirements is available upon request.

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Shanghai MicroPort EP MedTech Co., Ltd.

Address: Building 23&28, Lane 588, Tianxiong Rd., 201318 Shanghai, PEOPLE'S REPUBLIC OF CHINA

Postcode: 201318

Tel: +86 (21) 60969600

Fax: +86 (21) 20903925

E-mail: customerservice@everpace.com

Website: www.everpace.com

【Authorized representative IN the European community】

MicroPort Medical B.V.

Address: Paasheuvelweg 25, 1105BP Amsterdam, The Netherlands

Contact person: HH Chang

E-mail: cs@microport-int.com

Tel: +31 (0)20 545 0100-8

Fax: +31(0)205450109

【Authorized representative IN Russian Federation】

Synergy Consulting LLC

Address: Russian Federation, 109548, Moscow, Shosseynaya, 1, bld. 2, suite I, floor 4, office 26

E-mail: Info@synergyc.ru

The copyright of this document is owned by Shanghai MicroPort EP MedTech Co., Ltd.
Copying is not allowed without its consent.

A-C6G01-019, Rev.: A, Revision date: 2023-01

Instructions for use

**PathBuilder transseptal access kit
with Transseptal Needle**

Shanghai MicroPort EP MedTech Co., Ltd.

【WARNING】

- Do not reuse this device. The biological materials and foreign bodies attached cannot be removed thoroughly after use, so any reuse of this product may lead to adverse effect of patients.
- The "Use By" date on the device package and the sterile packaging should be inspected prior to use. Do not use the device if past the "Use By" date or the packaging appears damaged.
- Do not alter this device in any way.
- Do not use excessive force to advance or withdraw the device when resistance is encountered.
- Procedures present the potential for significant x-ray exposure, which can result in acute radiation injury as well as increased risk for somatic and genetic effects, to both patients and physician due to the x-ray beam intensity and duration of the fluoroscopic imaging. Procedures should only be performed after adequate attention has been given to the potential radiation exposure associated with the procedure, and steps taken to minimize this exposure.
- When used, device shall be disposed of as medical waste according local laws and regulations.
- Only those physicians who are trained in transseptal procedures and catheter delivery systems should use this device.
- If the patient with a history of heparin induced thrombocytopenia (HIT), bivalirudin is a preferred alternative to heparin for anticoagulation during procedures.

【Packing List】

PathBuilder transseptal access kit, in variants:

1. PathBuilder transseptal access kit with Transseptal needle, composed of:

- Transseptal needle PathBuilder 19G(1.10mm) S, length 560 mm – 1pc;
- Stylet – 1pc;
- Instructions for use

2. PathBuilder transseptal access kit with Transseptal needle, composed of:

1. Transseptal needle PathBuilder 19G(1.10mm) M, length 560 mm – 1pc;
2. Stylet – 1pc;
3. Instructions for use

TN M 71
TN M 89
TN M 98
TN L 71
TN L 98

3. PathBuilder transseptal access kit with Transseptal needle, composed of:


1. Transseptal needle PathBuilder 18G(1.21mm) M, length 710 mm; 890 mm; 980 mm—1pc;
2. Stylet – 1pc;
3. Instructions for use

4. PathBuilder transseptal access kit with Transseptal needle, composed of:

1. Transseptal needle PathBuilder 18G(1.21mm) L, length 710 mm; 890 mm; 980 mm—1pc;
2. Stylet – 1pc;
3. Instructions for use

【Overview】

Pathbuilder Transseptal Needle manufactured by Shanghai MicroPort EP MedTech Co., Ltd. is indicated for atrial septal puncture when accessing from femoral vein and to build a pathway from the right atrium to the left atrium together with the transseptal guiding introducer. The product structure diagram is as shown below:

Pathbuilder Transseptal Needle is consisted of thin section, curved section, main body, pointer flange, hub and stylet. The pointer flange indicates the curve orientation of the Pathbuilder Transseptal Needle. The thin section of the Pathbuilder Transseptal Needle fits well with that of the dilator which fixes the maximum length of the needle into the dilator.

【PRODUCT SPECIFICATION】

Pathbuilder Transseptal Needle can be classified into different specifications based on curve shape, diameter as well as length.

Pathbuilder Transseptal Needle has two specifications of 18Ga(1.21mm) and 19Ga(1.10mm) according to its diameter. The Pathbuilder Transseptal Needle of 18Ga(1.21mm) has two curve shapes of M and L, the effective length of the M curved type includes three specifications: 710mm, 890mm and 980mm, the effective length of the L curved type includes three specifications: 710mm, 890mm and 980mm; the Pathbuilder Transseptal Needle of 19Ga(1.10mm) has two curved types of S and M, the effective length of both is 560mm. Physicians can make choices based on actual needs.

Table 1: Specification of Pathbuilder Transseptal Needle

Model	Curve	Length (mm)	Diameter (Ga/mm)
TN S 56 19	S	560	19/1.10
TN M 56 19	M	560	19/1.10

TN M 71 18	M	710	18/1.21
TN M 89 18	M	890	18/1.21
TN M 98 18	M	980	18/1.21
TN L 71 18	L	710	18/1.21
TN L 98 18	L	980	18/1.21
TN L 89 18	L	890	18/1.21

【INDICATIONS FOR USE】

PathBuilder Transseptal Needle is indicated for building a pathway from femoral vein into left atrium through atrial septum puncture.

【INDICATIONS】

Atrial fibrillation

【CONTRAINDICATIONS】

- 1) Intra-atrial septal patch;
- 2) Atrial myxoma;
- 3) Myocardial infarction within the last two weeks;
- 4) Unstable angina;
- 5) Cerebral vascular accident (CVA) within the last two weeks;
- 6) Patients who do not tolerate anticoagulation therapy;
- 7) Patients with an active infection;
- 8) Left atrial thrombus which is presence after more than three weeks of anticoagulation

【TARGET GROUP】

18-75 years old adults, men or non-pregnant women

【Accessories Required for Using Product】

Steerable introducer

PathBuilder Steerable Introducer	PathBuilder Transseptal Needle		
	Model	Diameter	Length
SIS7185	TNM9818 TNL9818	18Ga(1.21mm)	980mm
SIM7185			
SIL7185			
SIXL7185			
SIS6185	TNM8918 TNL8918	18Ga(1.21mm)	890mm
SIM6185			

Transseptal guiding introducer

PathBuilder Transseptal Guiding Introducer	PathBuilder Transseptal Needle		
	Model	Diameter	Length
TS L0 63 8	TNM7118 TNL7118	18Ga (1.21mm)	710mm
TS L1 63 8		18Ga (1.21mm)	710mm
TS L2 63 8		18Ga (1.21mm)	710mm
TS L3 63 8		18Ga (1.21mm)	710mm
TS L4 63 8		18Ga (1.21mm)	710mm

TS L0 63 85		18Ga (1.21mm)	710mm
TS L1 63 85		18Ga (1.21mm)	710mm
TS L2 63 85		18Ga (1.21mm)	710mm
TS L3 63 85		18Ga (1.21mm)	710mm
TS L4 63 85		18Ga (1.21mm)	710mm
TS L0 81 8	TNM8918	18Ga (1.21mm)	890mm
TS L1 81 8	TNL8918	18Ga (1.21mm)	890mm
TS L2 81 8		18Ga (1.21mm)	890mm
TS L3 81 8		18Ga (1.21mm)	890mm
TS L4 81 8		18Ga (1.21mm)	890mm
TS L0 81 85		18Ga (1.21mm)	890mm
TS L1 81 85		18Ga (1.21mm)	890mm
TS L2 81 85		18Ga (1.21mm)	890mm
TS L3 81 85		18Ga (1.21mm)	890mm
TS L4 81 85		18Ga (1.21mm)	890mm
TS R0 63 8	TNM7118	18Ga (1.21mm)	710mm
TS R1 63 8	TNL7118	18Ga (1.21mm)	710mm
TS R2 63 8		18Ga (1.21mm)	710mm
TS R3 63 8		18Ga (1.21mm)	710mm
TS R4 63 8		18Ga (1.21mm)	710mm
TS R0 63 85		18Ga (1.21mm)	710mm
TS R1 63 85		18Ga (1.21mm)	710mm
TS R2 63 85		18Ga (1.21mm)	710mm
TS R3 63 85		18Ga (1.21mm)	710mm
TS R4 63 85		18Ga (1.21mm)	710mm
TS R0 81 8	TNM8918	18Ga (1.21mm)	890mm
TS R1 81 8	TNL8918	18Ga (1.21mm)	890mm
TS R2 81 8		18Ga (1.21mm)	890mm
TS R3 81 8		18Ga (1.21mm)	890mm
TS R4 81 8		18Ga (1.21mm)	890mm
TS R0 81 85		18Ga (1.21mm)	890mm
TS R1 81 85		18Ga (1.21mm)	890mm
TS R2 81 85		18Ga (1.21mm)	890mm
TS R3 81 85		18Ga (1.21mm)	890mm
TS R4 81 85		18Ga (1.21mm)	890mm

【Risks and Side Effects】

Complications that may occur during use of this device include scratches of endocardium and vessel wall cause reactions like inflammation or abscess; toxic reaction; complications caused by vascular puncture (such as pseudoaneurysm, arteriovenous fistula, local hematoma, or ecchymosis and infection); radiation injury; thrombosis / air embolism; injury to coronary artery; serious coagulation reactions; pericardial tamponade; cardiac / artery perforation; pericardial effusion; cerebral embolism / cerebral apoplexy / stroke; thoracotomy to remove foreign body.

【DIRECTIONS FOR USE】

1. Prepare and assemble equipment

Preparing the equipment requires the following items:

- One transseptal sheath and dilator,
- One matching transseptal needle and a stainless steel stylet,
- One 0.032" (0.81mm), 180cm guidewire with a 3mm "J" tip,
- Syringes for aspiration and flushing,
- Sterile heparinized saline

- Flush the transseptal sheath and dilator with sterile heparinized saline and insert the dilator fully into the sheath.
- Remove the stylet from the needle and flush the needle with sterile heparinized saline. Re-insert the stylet into the needle and lock it onto the hub.
- Insert the needle assembly into the sheath/dilator. Due to the stop feature of the dilator, there will be a gap between the dilator hub and the needle pointer flange. (See Fig.1. Magnification is 2.9).


Fig. 1

- Withdraw the needle assembly until the tip of the stylet is just within the tip of the dilator.
- 2. Advance sheath/dilator assembly into superior vena cava.**
- Introduce the prepared guidewire from femoral vein (right femoral preferred) to superior vena cava. Note: 0.032 inch (0.81mm) is the maximum guidewire diameter that can be used with the transseptal dilator.
 - Insert the transseptal sheath/dilator assembly over the guidewire and advance the assembly into the superior vena cava.
- 3. Position the needle/stylet assembly inside the sheath/dilator assembly.**
- Remove the guidewire from the dilator
 - Aspirate and then flush the dilator with clean heparinized saline, ensuring that no air enters the bloodstream.

- Separate the sheath and dilator by withdrawing the dilator by a distance sufficient to accommodate the needle curve. This will facilitate passage of the needle curve through the dilator and sheath hubs.
- Insert the needle/stylet assembly into the dilator, letting it rotate freely as it advances.
- After the needle curve is advanced beyond the hemostasis valve hub of the sheath, reconnect the sheath and dilator by sliding the sheath back over the dilator slowly.(do not advance the dilator)
- Advance the needle/stylet assembly and confirm that the needle tip is inside the dilator by fluoroscopy.
- After removing the stylet, attach a syringe to the needle hub and aspirate until blood return is observed. Then discard the syringe.
- Fill the syringe and flush the needle with clean heparinized saline, ensuring that the syringe is attached to the needle hub and no air enters the bloodstream.

4. Engage the fossa ovalis

- Visualize and identify anatomic landmarks. Set the fluoroscopy unit to an appropriate angle, parallel to the plane of the mitral valve and orthogonal to the plane of the septum. This will typically be LAO.
- Adjust the needle pointer so that the needle is perpendicular to the fossa ovalis (typically between 3:00 and 5:00 o'clock, as viewed from the foot end of the patient). (See Fig.2. Magnification is 2.9).


Fig. 2

- Confirm that the needle tip is inside the dilator by fluoroscopy. Then drag the entire assembly slowly. Prevent any movement of the assembly parts relative to each other. Maintain the previous orientation of the needle pointer.
- Observe the tip of the dilator during the drag for medial (or rightward) movement, indicating the tip has engaged the fossa ovalis.

5. Puncture the fossa ovalis

- Once the correct location is confirmed, advance the needle across the interatrial septum.
- Under pressure monitoring, entry into the left atrium is confirmed when the pressure tracing shows a left atrial pressure waveform.

If there is any resistance to needle advancement, re-evaluate the anatomic landmarks.

...a distance sufficient
...needle curve through
...advances.
...th, reconnect
...nce the

Caution: if pericardial or aortic entry occurs, do not advance the dilator over the needle. If the needle has penetrated the pericardium or aorta, it must be withdrawn. Monitor vital signs closely.

6. Advance the sheath/dilator assembly

- While maintaining a fixed needle position, advance the sheath/dilator assembly over the needle.

7. Advance the sheath over the fixed dilator and needle into the left atrium.

- Withdraw the needle into the dilator until it is just inside the tip. Maintain the position of the dilator across the septum.
- With the dilator in a fixed location, advance the sheath over the dilator.

8. Withdraw the needle and the dilator

- Disconnect any attachments to the needle hub
- Immediately aspirate blood through the side arm of the sheath after withdrawing the needle and dilator slowly at the same time. The blood should be arterial blood.
- Under pressure monitoring, the pressure tracing still shows a left atrial pressure waveform. The sheath is now in place in the left atrium.

9. Advance the catheter

- Insert electrode catheter into the sheath slowly and observe extension length of catheter by fluoroscopy.

10. Remove the sheath

- After removing the catheter, re-insert the guidewire over the sheath. Then re-insert the dilator over the guidewire and keep the sheath straight. At last remove the sheath, the dilator and the guidewire as a whole.

【STERILIZATION】

PathBuilder Transseptal Needle has been sterilized by epoxy ethane (ethylene oxide) before delivery.

【STORAGE & TRANSPORTATION REQUIREMENTS】

The product shall be stored in a shady and cool, dry, clean, and well-ventilated warehouse which is in natural air circulation environment.

In transit, the product shall be protected from heavy load, direct sunlight and rain or as specified in the ordering contract.

Situation	Ambient temperature	Relative air humidity	Atmospheric pressure
Transportation	-40°C to 55°C	10% to 90%	50kPa to 106kPa
Storage	18°C to 28°C	10% to 90%	50kPa to 106kPa
Working	5°C to 40°C	10% to 80%	70kPa to 106kPa

Operation conditions: Contact of the device with the internal environment of the body at a temperature of 32°C - 42°C and a humidity of 100% (contact with blood).

【MANUFACTURING DATA】

See label statement for service life.

【PERIOD OF VALIDITY】

The sterilization period of validity is three years under specified storage conditions.

【DISPOSAL】

Dispose the devices according to local environmental requirements.


Environmental Monitoring

Class of hazard for medical wastes by reference to specification of morphological composition (in accordance with Russian Sanitary Regulations and Standards 2.1.3684-21) Class B

【SYMBOLS DESCRIPTION】

	SYMBOL	DESCRIPTION
1.		DO NOT REUSE
2.		CONSULT INSTRUCTIONS FOR USE
3.		PROTECT FROM HEAT SOURCE AND RADIATION SOURCE
4.		KEEP DRY
5.		USE BY
6.		BATCH CODE
7.		CATALOGUE NUMBER
8.		STERILIZED USING ETHYLENE OXIDE
9.		DO NOT USE IF PACKAGE IS DAMAGED
10.		QUANTITY OF PRODUCT CONTAINED 1
11.		DATE OF MANUFACTURE
12.		MANUFACTURER

13.

AUTHORISED REPRESENTATIVE IN THE EUROPEAN
COMMUNITY

14.

TEMPERATURE LIMITATION

【AFTER-SALES SERVICE】

With “providing the medical sector with high quality and efficacy medical products” as its top operational objective, Shanghai MicroPort EP MedTech Co., Ltd. (hereinafter referred to as MicroPort EP) guarantees that its products are free of defects in materials or manufacturing when the clients receive them. For other questions relating to the products, please directly consult the company.

【SOLEMN STATEMENT】

The potential adverse event listed in the user manual, collected by Shanghai MicroPort EP MedTech Co., Ltd. (MicroPort EP), is based on the information, data from clinical trials and markets of China and other countries or regions instead of real condition occurred. The potential adverse event is taken as the product features sufficiently warned and explained instead of the defects or potential risks. Subject to the current scientific level and knowledge of MicroPort EP, the product may have some undetected or unknown potential risks and that cannot be considered defects or faults of the product.

MicroPort EP will not responsible for any medical cost or direct and indirect loss caused by the improper usage (including reuse), improper model choice, maloperation or other accident caused by human, illegal purchase, risks cannot be detected by current technology or potential adverse event or risks listed in this user manual, regardless of any guarantee, contract, tort or other regulations the claim is based on.

【GENERAL STANDARDS】

Shanghai MicroPort EP MedTech Co., Ltd. declares that PathBuilder transseptal access kit with Transseptal Needle meets the requirements ISO 11070, EN ISO 13485, EN ISO 14971. A complete list of requirements is available upon request.

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Shanghai MicroPort EP MedTech Co., Ltd.

Address: Building 23&28, Lane 588, Tianxiong Rd., 201318 Shanghai, PEOPLE’S REPUBLIC OF CHINA

Postcode: 201318

Tel: +86 (21) 60969600

Fax: +86 (21) 20903925

E-mail: customerservice@everpace.com

Website: www.everpace.com

【Authorized representative the European community】

MicroPort Medical B.V.

Address: Paasheuvelweg 25, 1105BP Amsterdam, The Netherlands

Contact person: HH Chang

E-mail: cs@microport-int.com

Tel: +31 (0)20 545 0100-8 Fax: +31(0)205450109

【Authorized representative IN Russian Federation】

Synergy Consulting LLC

Address: Russian Federation, 109548, Moscow, Shosseynaya, 1, bld. 2, suite I, floor 4, office 26

E-mail: Info@synergyc.ru

**The copyright of this document is owned by Shanghai MicroPort EP MedTech Co., Ltd.
Copying is not allowed without its consent.**

A-C6H01-018, Rev.: A, Revision date: 2023-01

Instructions for use

PathBuilder transseptal access kit with Steerable Introducer

Shanghai MicroPort EP MedTech Co., Ltd.

【WARNING】

- Do not reuse this device. The biological materials and foreign bodies attached cannot be removed thoroughly after use, so any reuse of this product may lead to adverse effect of patients.
- The "Use By" date on the device package and the sterile packaging should be inspected prior to use. Do not use the device if past the "Use By" date or the packaging appears damaged.
- Do not alter this device in any way.
- Do not use excessive force to advance or withdraw the device when resistance is encountered.
- Procedures present the potential for significant x-ray exposure, which can result in acute radiation injury as well as increased risk for somatic and genetic effects, to both patients and physician due to the x-ray beam intensity and duration of the fluoroscopic imaging. Procedures should only be performed after adequate attention has been given to the potential radiation exposure associated with the procedure, and steps taken to minimize this exposure.
- When used, device shall be disposed of as medical waste according local laws and regulations.
- Only those physicians who are trained in transseptal procedures and catheter delivery systems should use this device.
- If the patient with a history of heparin induced thrombocytopenia (HIT), bivalirudin is a preferred alternative to heparin for anticoagulation during procedures.

【Packing List】

PathBuilder transseptal access kit, in variants

1. PathBuilder transseptal access kit with Steerable Introducer length 710mm, composed of:

- Steerable Introducer PathBuilder, in variants: effective length 710 mm, curve diameter 20mm; 30mm; 50mm; 65mm – 1pc.;
- Dilator length 950mm – 1pc..
- Guidewire length 180cm, diameter 0.81mm – 1pc.;
- Instructions for use

2. PathBuilder transseptal access kit with Steerable Introducer length 610mm, composed of:

- Steerable Introducer PathBuilder, in variants: effective length 610 mm, curve diameter 20mm; 30mm – 1pc.;

- Dilator length 850mm – 1pc;
- Guidewire length 180cm, diameter 0.81mm – 1pc;
- Instructions for use

【Overview】

PathBuilder™ Steerable Introducer manufactured by Shanghai MicroPort EP MedTech Co., Ltd. is composed of a steerable sheath, a dilator and a guidewire. PathBuilder™ Steerable Introducer can be inserted into various cardiovascular catheters, including guiding cardiovascular catheters into the left side of the heart.

The steerable introducer is fitted with a hemostasis valve to minimize blood loss during catheter introduction and/or exchange. A sideport with three-way stopcock is provided for air or blood aspiration, fluid infusion, blood sampling and pressure monitoring. A handle equipped with a rotating collar to deflect the tip clockwise $\geq 180^\circ$ and counter-clockwise $\geq 90^\circ$. The steerable introducer features distal vent holes to facilitate aspiration and minimize cavitation and a radiopaque tip marker to improve fluoroscopic visualization.

【PRODUCT SPECIFICATION】

PathBuilder™ Steerable Introducer can be classified into different specifications based on the curve as well as effective length.

Table 1: Specification table of Pathbuilder™ Steerable Introducer

Model		SIS7185	SIM7185	SIL7185	SIXL7185	SIS6185	SIM6185
Curve		S	M	L	XL	S	M
Diameter of curve (mm)		20	30	50	65	20	30
Steerable sheath	Total length (cm/mm)	92/920	92/920	92/920	92/920	82/820	82/820
	Effective length (cm/mm)	71/710	71/710	71/710	71/710	61/610	61/610
	Outer diameter (F/mm)	11.5/3.83	11.5/3.83	11.5/3.83	11.5/3.83	11.5/3.83	11.5/3.83
	Inner diameter (F/mm)	8.5/2.83	8.5/2.83	8.5/2.83	8.5/2.83	8.5/2.83	8.5/2.83
Dilator	Length (cm/mm)	95/950	95/950	95/950	95/950	85/850	85/850
	Distal end	Curved	Curved	Curved	Curved	Curved	Curved
Guidewire	Length (cm/mm)	180/1800	180/1800	180/1800	180/1800	180/1800	180/1800
	Outer diameter (inch/mm)	0.032/0.81	0.032/0.81	0.032/0.81	0.032/0.81	0.032/0.81	0.032/0.81

【INDICATIONS FOR USE】

PathBuilder™ Steerable Introducer is indicated for introducing various cardiovascular catheters from femoral vein into left atrium through the interatrial septum, and when used with cardiovascular catheters, for assisting catheters in reaching the target site.

【INDICATIONS】

Atrial fibrillation

【CONTRAINDICATIONS】

- 1) Intra-atrial septal patch;
- 2) Atrial myxoma;
- 3) Myocardial infarction within the last two weeks;
- 4) Unstable angina;
- 5) Cerebral vascular accident (CVA) within the last two weeks;
- 6) Patients who do not tolerate anticoagulation therapy;
- 7) Patients with an active infection;
- 8) Left atrial thrombus which is presence after more than three weeks of anticoagulation

【TARGET GROUP】

For 18-75 years old adults, men or non-pregnant women

【Risks and Side Effects】

Complications that may occur during use of this device include, but are not limited to scratches of endocardium and vessel wall cause reactions like inflammation or abscess; toxic reaction; complications caused by vascular puncture (such as pseudoaneurysm, arteriovenous fistula, local hematoma, or ecchymosis and infection); radiation injury; thrombosis / air embolism; injury to coronary artery; serious coagulation reactions; pericardial tamponade; cardiac / artery perforation; pericardial effusion; cerebral embolism / cerebral apoplexy / stroke; thoracotomy to remove foreign body.

【Accessories Required for Using Product】

- Transseptal needle (to be purchased as PathBuilder transseptal access kit with Transseptal Needle)

PathBuilder Steerable Introducer	PathBuilder Transseptal Needle		
	Model	Diameter	Length
SIS7185	TNM9818 TNL9818	18Ga(1.21mm)	980mm
SIM7185			
SIL7185			
SIXL7185			
SIS6185	TNM8918 TNL8918	18Ga(1.21mm)	890mm
SIM6185			

- Cardiovascular catheter (to be purchased separately)
Ablation or diagnostic catheters (OD: 8F/2.67mm and 8.5F/2.83mm)
- Syringe (to be purchased separately)
Syringes with 6% (Luer) taper lock fitting

【DIRECTIONS FOR USE】

1. Prepare and assemble equipment

Preparing the equipment requires the following items:

- One steerable sheath and dilator,

- One matching transseptal needle and a stainless steel stylet,
- One 0.032" (0.81mm), 180cm guidewire with a 3mm "J" tip,
- Syringes for aspiration and flushing,
- Sterile heparinized saline
- Flush the steerable sheath and dilator with sterile heparinized saline and insert the dilator fully into the sheath.
- Remove the stylet from the needle and flush the needle with sterile heparinized saline. Re-insert the stylet into the needle and lock it onto the hub.
- Insert the needle assembly into the sheath/dilator. Due to the stop feature of the dilator, there will be a gap between the dilator hub and the needle pointer flange. (See Fig.1. Magnification is 2.9).


Fig. 1

- Withdraw the needle assembly until the tip of the stylet is just within the tip of the dilator. Measure the distance from the pointer flange and the dilator hub with length measuring tool. Record this measurement for use during the procedure. (See Fig.2. Magnification is 2.9).
- Caution:** it is critical during the procedure to maintain the distance between the pointer flange and the dilator hub. This ensures that the needle assembly does not extend beyond the dilator tip until it is deployed for transseptal crossing.


Fig. 2

2. Advance sheath/dilator assembly into superior vena cava.

- Introduce the prepared guidewire from femoral vein (right femoral preferred) to superior vena cava. Note: 0.032 inch (0.81mm) is the maximum guidewire diameter that can be used with the transseptal dilator.
- Insert the transseptal sheath/dilator assembly over the guidewire and advance the assembly into the superior vena cava.

3. Position the needle/stylet assembly inside the sheath/dilator assembly.

- Remove the guidewire from the dilator.
- Aspirate and then flush the dilator with clean heparinized saline, ensuring that no air enters the bloodstream.
- Separate the sheath and dilator by withdrawing the dilator by a distance sufficient to accommodate the needle curve. This will facilitate passage of the needle curve through the dilator and sheath hubs.
- Insert the needle/stylet assembly into the dilator, letting it rotate freely as it advances.
- After the needle curve is advanced beyond the hemostasis valve hub of the sheath, reconnect the sheath and dilator by sliding the sheath back over the dilator slowly. (do not advance the dilator)
- Advance the needle/stylet assembly until the distance between pointer flange and hub of the dilator is slightly larger than d measured in part 1.
- After removing the stylet, attach a syringe to the needle hub and aspirate until blood return is observed. Then discard the syringe.
- Fill the syringe and flush the needle with clean heparinized saline, ensuring that the syringe is attached to the needle hub and no air enters the bloodstream.

4. Engage the fossa ovalis

- Visualize and identify anatomic landmarks. Set the fluoroscopy unit to an appropriate angle, parallel to the plane of the mitral valve and orthogonal to the plane of the septum. This will typically be LAO.

- Adjust the needle pointer so that the needle is perpendicular to the fossa ovalis (typical between 3:00 and 5:00 o'clock, as viewed from the foot end of the patient). See Fig.3.


Fig. 3

- Confirm that the needle tip is inside the dilator by fluoroscopy and by your previous measurements. Then drag the entire assembly slowly. Prevent any movement of the assembly parts relative to each other. Maintain the previous orientation of the needle pointer.
- Observe the tip of the dilator during the drag for medial (or rightward) movement, indicating the tip has engaged the fossa ovalis.

5. Puncture the fossa ovalis

- Once the correct location is confirmed, advance the needle across the interatrial septum.
- Connect the pressure monitor line with the needle hub. Under pressure monitoring, entry into the left atrium is confirmed when the pressure tracing shows a left atrial pressure waveform. Left atrial access can be confirmed under fluoroscopy.

If there is any resistance to needle advancement, re-evaluate the anatomic landmarks.

Caution: if pericardial or aortic entry occurs, do not advance the dilator over the needle. If the needle has penetrated the pericardium or aorta, it must be withdrawn. Monitor vital signs closely.

6. Advance the sheath/dilator assembly

- While maintaining a fixed needle position, advance the sheath/dilator assembly over the needle.

7. Advance the sheath over the fixed dilator and needle into the left atrium.

- Withdraw the needle into the dilator until it is just inside the tip. Maintain the position of the dilator across the septum.
- With the dilator in a fixed location, advance the sheath over the dilator.

8. Withdraw the needle and the dilator

- Disconnect any attachments to the needle hub.

- Immediately aspirate blood through the side arm of the sheath after withdrawing the needle and dilator slowly at the same time. Then flush it with heparinized saline, taking care to avoid air bubbles.
- The sheath is now in place in the left atrium.

9. Advance the catheter

- Insert electrode catheter into the sheath slowly and observe extension length of catheter by fluoroscopy.
- Rotate the knob of sheath handle to assist catheter in reaching the target location.

10. Remove the sheath

- After removing the catheter, re-insert the guidewire over the sheath. Then re-insert the dilator over the guidewire and keep the sheath straight. At last remove the sheath, the dilator and the guidewire as a whole.

【STERILIZATION】

Pathbuilder™ Steerable Introducer has been sterilized by epoxy ethane (ethylene oxide) before delivery.

【STORAGE & TRANSPORTATION REQUIREMENTS】

The product shall be stored in a shady and cool, dry, clean, and well-ventilated warehouse which is in natural air circulation environment.

In transit, the product shall be protected from heavy load, direct sunlight and rain or as specified in the ordering contract.

Situation	Ambient temperature	Relative air humidity	Atmospheric pressure
Transportation	-40°C to 55°C	10% to 90%	50kPa to 106kPa
Storage	18°C to 28°C	10% to 90%	50kPa to 106kPa
Working	5°C to 40°C	10% to 80%	70kPa to 106kPa

Operation conditions: Contact of the device with the internal environment of the body at a temperature of 32°C - 42°C and a humidity of 100% (contact with blood).

【PERIOD OF VALIDITY】

The sterilization period of validity is three years under specified storage conditions.

【MANUFACTURING DATA】

See label statement for service life.


【DISPOSAL】

Dispose the devices according to local environmental requirements.

Environmental Monitoring

Class of hazard for medical wastes by reference to specification of morphological composition (in accordance with Russian Sanitary Regulations and Standards 2.1.3684-21) Class B

【SYMBOLS DESCRIPTION】

SYMBOL	DESCRIPTION
1. 	DO NOT REUSE
2. 	CONSULT INSTRUCTIONS FOR USE
3. 	PROTECT FROM HEAT SOURCE AND RADIATION SOURCE
4. 	KEEP DRY
5. 	USE BY
6. 	BATCH CODE
7. 	CATALOGUE NUMBER
8. 	STERILIZED USING ETHYLENE OXIDE
9. 	DO NOT USE IF PACKAGE IS DAMAGED
10. 	QUANTITY OF PRODUCT CONTAINED 1
11. 	DATE OF MANUFACTURE
12. 	MANUFACTURER
13. 	AUTHORISED REPRESENTATIVE IN THE EUROPEAN COMMUNITY
14. 	TEMPERATURE LIMITATION

【AFTER-SALES SERVICE】

With “providing the medical sector with high quality and efficacy medical products” as its top operational objective, Shanghai MicroPort EP MedTech Co., Ltd. (hereinafter referred to as MicroPort EP) guarantees that its products are free of defects in materials or manufacturing when the clients receive them. For other questions relating to the products, please directly consult the company.

【SOLEMN STATEMENT】

Shanghai MicroPort EP MedTech Co., Ltd. (hereinafter referred to as MicroPort EP Co.) expressly states herein that its Pathbuilder™ Steerable Introducer is disposable and cannot be reused. MicroPort EP Co. will not recommend, indicate and imply in any manner the reusability of the system, and will not assume the responsibility for any accident or product damage resulting from reuse. Pathbuilder™ Steerable Introducer can be connected and used only with the compatible devices specified herein, and MicroPort EP Co. will not assume the responsibility for damage to product device, procedure failure and the like resulting from operating mistakes or any other human factors.

【GENERAL STANDARDS】

Shanghai MicroPort EP MedTech Co., Ltd. declares that PathBuilder transseptal access kit with Steerable Introducer meets the requirements ISO 11070, EN ISO 13485, EN ISO 14971. A complete list of requirements is available upon request.

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Shanghai MicroPort EP MedTech Co., Ltd.

Address: Building 23&28, Lane 588, Tianxiong Rd., 201318 Shanghai, PEOPLE'S REPUBLIC OF CHINA

Postcode: 201318

Tel: +86 (21) 60969600

Fax: +86 (21) 20903925

E-mail: customerservice@everpace.com

Website: www.everpace.com

【Authorized representative IN the European community】

MicroPort Medical B.V.

Address: Paasheuvelweg 25, 1105BP Amsterdam, The Netherlands

Contact person: HH Chang

E-mail: cs@microport-int.com

Tel: +31 (0)20 545 0100-8

Fax: +31(0)205450109

【Authorized representative IN Russian Federation】

Synergy Consulting LLC

Address: Russian Federation, 109548, Moscow, Shosseynaya, 1, bld. 2, suite I, floor 4, office 26

E-mail: Info@synergyc.ru

**The copyright of this document is owned by Shanghai MicroPort EP MedTech Co., Ltd.
Copying is not allowed without its consent.**

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