INTEGRA®

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Integra®

NeuroBalloon™ Catheter

Instructions for Use
Mode d'Emploi
Gebrauchsanleitung
Istruzioni per l'Uso
Instrucciones de Uso
Gebruiksaanwijzing

I, the undersigned

Andrea Tommasoli

certify that this copy is identical to the original.

At St Priest on October 05, 2018

Signature :

[Signature]

INTEGRA LIFESCIENCES (France) SAS
Regulatory Affairs Department
Immeuble Séquoia 2
97 Allée Alexandre Borodine
69791 SAINT PRIEST CEDEX

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
Monsieur Andrea TOMMASOLI

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
MANUFACTURER:
Integra NeuroSciences Implants (France) S.A.S.
2905 Route des Dolines
06921 SOPHIA ANTIPOLIS CEDEX, FRANCE
MADE IN FRANCE

CE
0086

BL 915009 Rev01

...
...
...
...

Notes

DESCRIPTION

The sterile, disposable NeuroBalloon™ Catheter consists of a 60 cm minimum usable length, single lumen catheter with a distal inflatable balloon attached to the catheter's tip (see figure page 2). It is designed to be inserted through a lumen presenting a French 4 minimum diameter (1.35 mm minimum). The balloon inflates with a 1 ml (1 cc) syringe which is included in the package. Upon inflation following the recommended procedure described in these instructions for use, the balloon inflates in a dumb-bell shape.

The diameter of the inflated balloon waist – secondary to a pre-inflation at 1 ml (1 cc) of air – is between 3.5 mm minimum with 0.6 ml (0.6 cc) of air and 6 mm maximum with 1 ml (1 cc) of air (when measured immediately after inflation at the atmosphere pressure). **NEVER INFLATE WITH MORE THAN 1 ML (1CC) OF AIR.**

The special dumb-bell balloon shape facilitates positioning, and maintains the balloon within the cerebral membrane. As the air is slowly injected, the sections of the balloon inflate either at the same time or in sequence (the proximal section inflates shortly before the distal section), positioning the balloon waist in the center of the prepunctured membrane opening. The membrane opening is progressively dilated to the outer diameter of the balloon waist. During the dilatation procedure, the balloon waist can be visualized through the clear membrane of the inflated proximal balloon section.

Materials in contact with body tissues or body fluids are polyether block amide, silicone elastomer and cyanoacrylate glue.

INDICATION

The NeuroBalloon™ Catheter is intended for dilatation of cerebral membrane fenestrations under direct or endoscopic visualization during intracranial procedures.

WARNING

Not for intravascular use. Do not use in rigid neuro-tissues, such as dilatation of aqueducts stenosis or thick membranes of arachnoid cysts, since the balloon is not designed to withstand high pressure. If the dilatation of the target site is not successful with 1 ml of air, the tissue may be too rigid for this instrument. Over-inflation should not be attempted as it may damage the balloon; another technique should be used.

PRECAUTIONS

- Do not expose to organic solvents.
- Have a duplicate of the sterile product available during the procedure.
- Do not attempt to pass the balloon catheter through a lumen smaller than French 4 (1.35 mm minimum).
- If a strong resistance is felt during introduction, discontinue the procedure and determine the cause of resistance before proceeding. If the cause of the resistance cannot be determined, withdraw the deflated balloon catheter under negative pressure.
- Inflate balloon with air only. Do not inflate with saline or contrast media.
- Do not inflate the balloon in excess of 1 ml (1 cc) of air. Over-inflation increases the possibility of balloon damage.
- Slowly inflate the balloon.
- Do not advance or retract the catheter unless the balloon is fully deflated under negative pressure.
- These balloon catheters are extremely delicate instruments. Extra care must be taken in their handling and use. Even though each product is inspected and tested to assure integrity during manufacture, a pin hole or leak may still occasionally develop during use due to the special construction of the device. Prior to use, inspect the balloon catheter by gently inflating with 1 ml (1 cc) of air and examine to verify balloon integrity and shape (see paragraph below: Instructions for Use – Preparation). If the balloon still does not maintain inflation during a procedure, remove it from the surgical site as indicated in step 15 of the instructions for use below, and replace it with new balloon catheter.

This product is for single use only.

Caution: This product is fragile. Resterilization and/or re-use may affect its performances and characteristics.

This product is sterilized with ethylene oxide. Do not use if the package is open or damaged. Use the device prior to the "Use Before" date on the package label.

Caution: Do not re-sterilize. Integra will not be liable for any or all damages including, but not limited to, direct, indirect, incidental, consequential or punitive damages resulting from or related to re-sterilization.

COMPLICATIONS

Complications may occur at any time during or after the procedure. Possible complications include, but are not limited to: hemorrhage, infection, elevated intracranial pressure, bradycardia, air embolism, neurological disorders, CSF leak, subdural collection, death. Intracranial procedures should not be attempted by physicians unfamiliar with the possible complications.

Risk of incident such as balloon rupture may potentially occur during procedure. This risk is minimized when following the instructions for use below.

INSTRUCTIONS FOR USE

Preparation

1. Carefully remove the balloon catheter from its protective, transparent plastic sheath before use. Caution: When pulling the catheter out of the sheath, be careful not to kink the catheter.
2. The use of the syringe provided within the package is strongly recommended. Fill the syringe with 1 ml (1cc) of air and connect it to the catheter hub.
Warning: Do not use more than 1 ml (1 cc) of air; over-inflation may damage the balloon.
3. Inflate the balloon by slowly injecting with 1 ml (1 cc) of air. The balloon should inflate into a dumb-bell shape (see figure page 2), both sections inflating either at the same time or the proximal section inflating shortly before the distal section. Carefully inspect the catheter with the balloon inflated. A balloon that leaks or inflates in a grossly asymmetric (eccentric) manner should not be used.
4. Completely deflate the balloon by pulling back on the syringe's plunger in order to collapse the balloon over the catheter tip. Disconnect the syringe from the catheter hub.
5. Gently moisten the exterior of the balloon catheter with sterile, non pyrogenic saline solution (or similar isotonic solution) to facilitate balloon introduction into the cannula or other introducing tool.

Insertion, Inflation and Withdrawal

6. Using a technique of choice, insert an introducing tool (e.g.: a neuroendoscope or a multilumen cannula) close to the membrane to be punctured.
7. Create the puncture opening in the cerebral membrane.
8. Connect the syringe unfilled to the catheter hub. Pull back on the syringe's plunger in order to collapse the balloon over the catheter tip and to remove any residual inflation of the balloon. Leave the syringe connected.
9. Carefully introduce the deflated balloon catheter through the 4 French (or bigger) lumen of the introducing tool, just proximal to the tool's distal tip. If any resistance is felt during advancement, gently remove the catheter, check for possible residual inflation of the balloon or other cause of obstruction. Never force the catheter down the lumen as this may damage the balloon.
10. Once the balloon tip clearly appears under direct visualization or on the endoscopy visualization system screen, disconnect the syringe and fill it with 1 ml (1 cc) of air. Then, reconnect the syringe to the catheter hub.
11. Carefully position the balloon within the opening in the cerebral membrane. A black marker within the balloon indicates the balloon waist. This black marker should be positioned within the membrane opening.
12. Inflate the balloon by slowly injecting 0.6 ml (0.6 cc) of air (advance the plunger to the graduation 0.4 ml (0.4 cc) of the syringe).
The balloon will inflate into a dumb-bell shape, the proximal section inflating shortly before the distal section. Upon inflation, the proximal and distal balloon sections will capture the membrane, progressively dilating the hole to the outer diameter of the balloon waist. The dilatation can be visualized through the transparent wall of the proximal balloon section.
13. Maintain inflation for 2 minutes to facilitate hemostasis of eventual microhemorrhages.
14. It should be noted that air slowly diffuses through the balloon wall. If required, inject the residual 0.4 ml (0.4 cc) left in the syringe to compensate for air loss. If necessary, maintain the inflation for 2 more minutes.
Warning: Do not try to inject more than 1 ml (1 cc) of air as over-inflation may damage the balloon.
15. Deflate the balloon by completely pulling back on the syringe's plunger and, without disconnecting the syringe, gently remove the catheter from the introducing tool under negative pressure.

HOW SUPPLIED

Integra NeuroBalloon™ Catheter is supplied sterile and non-pyrogenic in a double-wrap packaging.

Disposal

After patient use, the system must be handled as biohazardous material and disposed of in accordance with applicable Federal, State, local, or international environmental requirements following facility protocols.

PRODUCT INFORMATION DISCLOSURE

INTEGRA HAS EXERCISED REASONABLE CARE IN THE CHOICE OF MATERIALS AND MANUFACTURE OF THIS PRODUCT. INTEGRA EXCLUDES ALL WARRANTIES, WHETHER EXPRESSED OR IMPLIED BY OPERATION OF LAW OR OTHERWISE, INCLUDING BUT NOT LIMITED TO, ANY IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. INTEGRA SHALL NOT BE LIABLE FOR ANY INCIDENTAL OR CONSEQUENTIAL LOSS, DAMAGE, OR EXPENSE, DIRECTLY OR INDIRECTLY ARISING FROM USE OF THIS PRODUCT. INTEGRA NEITHER ASSUMES NOR AUTHORIZES ANY PERSON TO ASSUME FOR IT ANY OTHER OR ADDITIONAL LIABILITY OR RESPONSIBILITY IN CONNECTION WITH THESE PRODUCTS. INTEGRA INTENDS THAT THIS DEVICE SHOULD BE USED ONLY BY PHYSICIANS WITH EDUCATIONAL AND TRAINING BACKGROUND ENABLING THE PROPER USE OF THE DEVICE.

Returned Goods Policy

Authorization from customer service must be obtained prior to returning product.

Sterile products must be returned in unopened packages with manufacturer's seals intact to be accepted for replacement or credit unless returned due to a complaint or product defect.

Determination of a product defect will be made by Integra.

Credit will be issued for goods returned prior to 90 days from ship date and may be subject to a restocking charge of up to 20%. This assumes that the product returned is not damaged and can be verified to have not been used or opened.

Product Order Information

All products can be ordered through your Integra NeuroSpecialist or customer service representative or by contacting:

Integra LifeSciences Corporation

311 Enterprise Drive
Plainsboro, NJ 08536, U.S.A.
Telephone: 1-800-654-2873
Outside the US: +1-609-275-0500
Fax: +1-609-275-5363

Integra NeuroSciences Implants (France) S.A.S.

2905 Route des Dailines
06921 Sophia Antipolis CEDEX, France
Telephone: +33 (0) 4 93 95 56 00
Fax: +33 (0) 4 93 65 40 20

Integra LifeSciences Services (France)

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Date of revision: March 2016

SYMBOLS/SYMBOLS/SYMBOLS/SYMBOLS/SYMBOLS/SYMBOLS

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U.S. Patent 5,718,712

Integra and the Integra logo are registered trademarks of Integra LifeSciences Corporation.
NeuroBalloon is a trademark of Integra LifeSciences Corporation or its subsidiaries
in the United States and/or other countries.

© 2016 Integra LifeSciences Corporation. All Rights Reserved.
<http://www.integralife.com>

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