APOSTILLE

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1. Country: United States of America
Pays | País:

This public document Le présent acte public | El presente documento público
JENNIFER KESTER

2. has been signed by
a été signé par | ha sido firmado por
Notary Public, State of Colorado

3. acting in the capacity of
agissant en qualité de | quien actúa en calidad de
JENNIFER KESTER, Notary Public State of Colorado

4. bears the seal / stamp of
est revêtu du sceau / timbre de | y está revestido del sello / timbre de

Certified Attesté | Certificado

5. at Denver, Colorado
à | en

6. the 3rd day of February, 2021
le | el día

7. by Jena Griswold, Secretary of State of the State of Colorado
par | por

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9. Seal / stamp
Sceau / timbre | Sello / timbre

10. Signature:
Signature | Firma

Jena Griswold

This Apostille only certifies the authenticity of the signature and the capacity of the person who has signed the public document, and, where appropriate, the identity of the seal or stamp which the public document bears. This Apostille does not certify the content of the document for which it was issued. This Apostille is not valid for use anywhere within the United States of America, its territories or possessions. To verify the issuance of this Apostille, see:

www.sos.state.co.us/auth

This certificate does not constitute an Apostille under the Hague Convention of 5 October 1961, when it is presented in a country which is not a party to the Convention. In such cases, the certificate should be presented to the consular section of the mission representing the country.

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
To whom it may concern

Subject: IFU

The company Terumo BCT Inc., having its registered office at 10811 W. Collins Ave. Lakewood, Colorado, 80215, USA, hereby declare that the attached IFU with addition is true.

Best regards,


Janet G Johnson

Terumo BCT
10811 W. Collins Ave
Lakewood, CO 80215
USA

Vice President, Global Regulatory Affairs

Subscribed and sworn to before me this 13th of January 2021 in the County of Jefferson, State of Colorado, USA.

Seal:

Notary public

Terumo BCT, Inc.
10811 West Collins Ave.
Lakewood, Colorado 80215-4440
USA
SA Phone: 1.877.339.4228
Phone: +1.303.231.4357
Fax: +1.303.542.5215

Terumo BCT Europe N.V.
Europe, Middle East and Africa
Ikarostraat 41
1930 Zaventem
Belgium
Phone: +32.2.715.0590
Fax: +32.2.721.0770

Terumo BCT Asia Pte. Ltd.
89 Science Park Drive
#04-25 (Lobby B)
The Rutherford
Singapore 118261
Phone: +65.6715.3778
Fax: +65.6774.1419

Terumo BCT Latin America S.A
La Pampa 1517-12th Floor
C1428DZE
Buenos Aires
Argentina
Phone: +54.11.5530.5200
Fax: +54.11.5530.5201

Terumo BCT Japan, Inc.
Tokyo Opera City Tower 49F,
3-20-2, Nishi-Shinjuku,
Shinjuku-ku, Tokyo 163-1450,
Japan
Phone: +81.3.6743.7890
Fax: +81.3.6743.9800

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Description

The Mirasol® Pathogen Reduction Technology (PRT) System consists of the Mirasol Illuminator, the Mirasol Manager Software (optional), and a Mirasol Disposable Kit. The sterile disposable kit contains an illumination bag, a storage bag, and riboflavin solution.

The Mirasol PRT system is for use with Whole Blood. The treatment procedure involves mixing the Whole Blood product with riboflavin solution in the illumination bag and illuminating the mixture with a calculated light energy dose in the Illuminator. The Mirasol-treated Whole Blood product is then transferred from the illumination bag into the storage bag.

Intended Use

The Mirasol Pathogen Reduction Technology (PRT) System for Whole Blood is an extracorporeal blood treatment process that uses ultraviolet light energy and riboflavin (vitamin B₂) to reduce the pathogen load and to inactivate white blood cells in Whole Blood for transfusion.

Indications for Use

- Mirasol PRT system treatment can be used instead of gamma irradiation to prevent transfusion-associated graft-versus-host disease (TA-GVHD).
- Mirasol PRT system treatment has demonstrated the ability to reduce the occurrence of transfusion-related transmission of malaria.
- For additional information regarding the Mirasol-treated Whole Blood product, including in vitro and in vivo data, refer to the Physician's Information Booklet (Terumo BCT part number 777714-013).

Note: Laboratory studies of photochemical treatment processing with the Mirasol PRT system for Whole Blood have shown a reduction in the levels of certain viruses, bacteria, and parasites, as well as effective inactivation of residual white blood cells. However, no pathogen reduction process has been shown to eliminate all pathogens or excludes the possibility of adverse events due to contaminating white blood cells.

Contraindications

There are no known contraindications for use of this product.

Warnings

1. **Do Not Reuse/Not for Reuse:** Terumo BCT, Inc. products bearing the "Do Not Reuse" symbol are intended for single use only and are not intended to be reused or re-sterilized in any manner. Terumo BCT cannot ensure the functionality or sterility of the product if it is reused or re-sterilized.

Reuse of a single-use product could result in:

- Product performance issues due to a loss of product integrity, including but not limited to the following:
 - Fluid leaks
 - Parts that are warped or deformed
 - Plastics that are brittle and discolored
 - Filters that have reduced filtration capabilities
- Exposure to excessive ethylene oxide (EO) residuals

- Viral infections such as hepatitis or human immunodeficiency virus (HIV)
- Bacterial infections
- Cross-contamination

Any of these risks could result in serious injury or death. These risks are shared by product users, donors, patients, and recipients of end products of the device.

2. This product contains phthalates (DEHP). Whole blood donors are not exposed to DEHP; the potential health risk to apheresis donors is low, because the time-averaged DEHP dose exposure is very low. Patient groups that include pregnant or nursing women and children are considered to be the most at risk to potential harmful effects of exposure to DEHP. However, regulatory bodies have noted that the benefit of doing a needed procedure is far greater than the risk associated with exposure to DEHP. It is the responsibility of the treating physician to balance this risk for the patient.
3. Only Whole Blood products collected in either CPD anticoagulant or CPDA-1 anticoagulant have been evaluated for use with the Mirasol PRT system. Do not use Whole Blood products collected in any other anticoagulant.

Cautions

For additional information regarding the Mirasol-treated Whole Blood product, including in vitro and in vivo data, refer to the Physician's Information Booklet (Terumo BCT part number 777714-013).

General

1. Each operator should be thoroughly familiar with these Instructions for Use and the Mirasol Illuminator operator's manual.
2. Use aseptic technique throughout all procedures and follow standard blood establishment practices.

Disposable Kit

1. Use the Mirasol Whole Blood Disposable Kit described here only for treatment of Whole Blood. Contact your Terumo BCT representative for information about other Mirasol products.
2. Use only Disposable Kits manufactured for the Mirasol PRT system to ensure compatibility and adherence with operational specifications.
3. Do not use the disposable kit if the kit label, the bag label, or the riboflavin solution bag label has passed the printed expiration date (YYYY-MM-DD).
4. Carefully inspect for leaks when using the disposable kit. If a leak is observed, manage the product according to your facility's standard operating procedure (SOP).
5. The Illumination/Storage Set is sterilized with ethylene oxide (EO), which in rare cases has been reported to cause recipient reactions.
6. The Mirasol Illumination Bag has not been qualified for centrifugation.

Riboflavin Solution

1. Riboflavin solution is sensitive to light. Add the riboflavin solution to the Whole Blood product within 1 hour of removing the riboflavin solution bag from the foil pouch.
2. Riboflavin solution is not intended for direct intravenous infusion.
3. Do not add more than one bag of riboflavin solution to any one product.

Blood Components

1. All Whole Blood products must be illuminated according to the treatment time constraints listed in Table 1 in order to ensure adequate product quality after treatment and storage.
2. After adding riboflavin solution and until the time of transfusion, protect the Whole Blood product from exposure to direct sunlight or strong ambient light sources. See Table 3 for storage specifications.
3. Mirasol PRT system treatment and storage of Whole Blood result in levels of extracellular potassium that are higher than those in non-treated Whole Blood units. Considerations should be made for transfusion recipients who may be susceptible to hyperkalemia.

Disposable Kit Overview

The Mirasol Whole Blood Disposable Kit consists of:

- Clear outer packaging with a white peel strip for opening
- One opaque foil pouch containing one bag of riboflavin solution (500 μ M riboflavin in normal saline) (see Figure 1)
- One illumination/storage set with an illumination bag, a storage bag, and attached tubing and frangible connectors (see Figure 1)


Figure 1: Mirasol Whole Blood Disposable Kit

- 1 Inlet line to the illumination bag with a luer and a blue end cap
- 2 Sealed inlet line
- 3 Drain line from the illumination bag to the storage bag
- 4 Vent line with a blue end cap on the storage bag
- 5 Storage bag
- 6 Pre-attached label
- 7 Set type label

- 8 Spike port
- 9 Frangible connector on the spike port
- 10 Frangible connector on the drain line from the illumination bag to the storage bag
- 11 Illumination bag
- 12 Opaque foil pouch
- 13 Riboflavin solution bag
- 14 Drain line for the riboflavin solution

Note: Callouts ① through ④ are reference points for the steps in "Treating the Whole Blood Product with the Mirasol PRT System" beginning on page 7.

Ancillary Equipment

The following ancillary equipment is needed to complete a Mirasol treatment procedure:

- Tubing sealer device
- Sterile connection device
- Hematocrit measurement device

Use this equipment according to your facility's SOP.

Incoming Product Specifications

Table 1: Incoming product specifications

Characteristic	Specification
Anticoagulant	CPD or CPDA-1
Starting Whole Blood volume	483 mL to 541 mL, including anticoagulant and donor blood This acceptable range is calculated as follows: (450 mL $\pm$ 5%/-10% donor blood) + (63 mL $\pm$ 5 mL AC)
Treatment time constraints	Whole Blood must be treated within 24 hours of collection.
Illumination hematocrit (after riboflavin solution has been added)	30% to 45%
Storage temperature before Mirasol PRT system treatment	Follow your facility's SOP for Whole Blood products.

Whole Blood donation volume for Mirasol PRT system treatment is specified as 450 mL $\pm$ 5%/-10%. The addition of anticoagulant increases the volume by 63 mL $\pm$ 5 mL. Prior to Mirasol PRT system treatment, 35 mL $\pm$ 5 mL riboflavin solution is added to the Whole Blood unit. The following table shows the low limit and high limit of possible volumes and the weight of each volume.

Table 2: Whole Blood volume and weight specifications

	Low Limit	Target	High Limit
Volume/weight* of Whole Blood collection (450 mL +5%/-10%)	405 mL/426 g	450 mL/474 g	473 mL/502 g
Volume/weight* of Whole Blood unit including anticoagulant (63 mL $\pm$ 5 mL)	463 mL/485 g	513 mL/539 g	541 mL/572 g
Volume/weight* of Whole Blood unit including anticoagulant and riboflavin solution (35 mL $\pm$ 5 mL)	493 mL/516 g	548 mL/575 g	581 mL/613 g
*Weight does not include disposable set tare weight.			

Instructions for Use

Treating the Whole Blood Product with the Mirasol PRT System

Caution: Use the Mirasol Whole Blood Disposable Kit described here only for treatment of Whole Blood. Different Mirasol disposable kits are available for treatment of platelets and plasma.

1. Review the treatment time constraints listed in Table 1 and the storage conditions listed in Table 3.
2. Just prior to use, open the white peel strip on the outer packaging of the disposable kit. Remove and inspect the illumination/storage set.

Caution: Do not use the disposable kit under any of the following circumstances:

- a. The blue end cap on either the inlet line (⊙) or the storage bag line (⊗) is not attached.
- b. The integrity of the illumination/storage set is compromised for any reason.
- c. The frangible connectors are broken.
- d. The tubing is severely kinked.
- e. There is debris on the outer surface of the illumination bag. Do not clean with disinfection chemicals.
- f. The barcodes are illegible.
- g. The pre-attached label is missing or not securely attached.

Note: In general, you should use the riboflavin solution bag with the illumination/storage set from the same disposable kit. However, if any of the circumstances listed in items c through g in the **Caution** above occurs after the riboflavin solution has been added to the Whole Blood product, transfer the Whole Blood product to a new illumination/storage set from the same lot number. Discard the original illumination/storage set and the unused riboflavin solution.

3. Permanently seal and disconnect the vent line on the storage bag (Line ④) using a tubing sealer device. Seal the line as near to the storage bag as possible, as shown in Figure 2.


Figure 2: Removing the vent line

4. Before connecting the whole blood product bag to the illumination bag, verify that the Whole Blood product contains only CPD anticoagulant or CPDA-1 anticoagulant. If the Whole Blood product does not contain either CPD anticoagulant or CPDA-1 anticoagulant, do not continue with the Mirasol process. Only Whole Blood products collected in either CPD anticoagulant or CPDA-1 anticoagulant have been evaluated for use with the Mirasol PRT system.

5. Using a sterile connection device, connect the whole blood product bag to the illumination bag inlet line (Line ⑥), and remove the part of the line that has a luer and a blue end cap, as shown in Figure 3. Do not open the sterile connection at this time. (You will transfer the Whole Blood and the riboflavin solution to the illumination bag at the same time in step 8.)


Figure 3: Connecting the whole blood product bag to the illumination bag inlet line

6. Transfer an identification label from the donor whole blood product bag to the storage bag according to your facility's standard operating procedure (SOP), using the space on the pre-attached label. Additional labels may be applied over the pre-attached label on the storage bag.

Note: Labels may not be applied to the illumination bag before illumination. The amount of UV energy delivered to the Whole Blood product is calculated using the surface area of the non-labeled illumination bag, among other factors.

Note: The following set type label should be present on the illumination bag in the position shown in Figure 3, Item 7:

7. Just before you are ready to transfer the riboflavin solution to the illumination bag, open the foil pouch, remove the riboflavin solution bag, and inspect the bag and solution.

Caution: Do not use the riboflavin solution if:

- The integrity of the riboflavin solution bag is compromised for any reason.
- Any crystals or particulates are visible in the riboflavin solution.

Note: In general, you should use the riboflavin solution bag with the illumination/storage set from the same disposable kit. However, if either of the two circumstances listed in the **Caution** above occurs, use a new riboflavin solution bag of the same lot number from another disposable kit and discard the unused illumination/storage set.

8. Transfer the Whole Blood product and the riboflavin solution into the illumination bag at the same time by performing these steps:
 - a. Using a sterile connection device, connect the tubing on the riboflavin solution bag (Line E) to the sealed inlet line on the illumination bag (Line D). See Figure 4.
 - b. Open both sterile connections and drain the Whole Blood product and the riboflavin solution into the illumination bag.


Figure 4: Transferring the Whole Blood product and the riboflavin solution into the illumination bag

- c. Permanently seal and disconnect the inlet line between the empty whole blood product bag and the Y connector, leaving the riboflavin solution bag attached. See Figure 5.


Figure 5: Disconnecting the whole blood product bag

9. Express residual air from the illumination bag into the empty riboflavin solution bag by performing these steps:
- Hold the illumination bag so that the ports are facing up, as shown in Figure 6.
 - Gently apply pressure to the illumination bag to move the residual air and Whole Blood up through the inlet line, allowing the Whole Blood to fill the riboflavin solution bag approximately half full.
 - Close the pinch clamp, and then gently mix the Whole Blood in the riboflavin solution bag to ensure that the solution is homogeneous. Gently tap the ports on the riboflavin solution bag to remove any solution that may have been trapped during mixing. Do not shake or vigorously mix the Whole Blood in the riboflavin solution bag, as this can cause the formation of air bubbles and possible hemolysis.
 - Open the pinch clamp to drain the Whole Blood from the riboflavin solution bag and the inlet line into the illumination bag, closing the pinch clamp on the inlet line before the air interface reaches the illumination bag port. Gently mix the Whole Blood and the riboflavin solution to obtain a homogeneous solution.

Note: It is very important that the product be gently but thoroughly mixed to achieve an accurate hematocrit measurement.

- e. Open the pinch clamp on the inlet line and hold the illumination bag so that the ports are facing up, as shown in Figure 6. Gently apply pressure to the illumination bag to move the residual air and Whole Blood up through the inlet line until it reaches halfway between the Y connector and the riboflavin solution bag. Close the pinch clamp on the inlet line.


Figure 6: Expressing residual air from the illumination/storage set

Residual air should now be removed from the illumination bag. The following figure shows acceptable and unacceptable amounts of residual air.


Figure 7: Residual air in the illumination bag

- 1** Acceptable amount of residual air
- 2** Too much residual air

10. Permanently seal the inlet line just below the riboflavin solution bag using a tubing sealer device (see Figure 8).
11. Disconnect and discard the riboflavin solution bag.
12. Permanently seal the inlet line as near as possible to the illumination bag, using a tubing sealer device. The sealed inlet line must be less than 2.5 cm (1 in) long for proper fit into the Illuminator. Disconnect the tubing segment between the illumination bag and the riboflavin solution bag. (This segment is approximately 18 cm to 30 cm [8 in to 12 in] long. See Figure 8.) This segment will be used to collect a hematocrit sample.


Figure 8: Sealing and disconnecting the riboflavin solution bag and the hematocrit sample segment

13. Extract a hematocrit sample from the segment of inlet line tubing that was disconnected.
 - Extract the sample using, for example, a plain tube or a microcentrifuge tube; do not use a tube that contains EDTA or AC, which will dilute the sample and change its hematocrit.
 - Follow your facility's SOP first to mix the sample so that it is homogeneous, and then to measure the hematocrit.

Note: Always use the hematocrit of the pre-illumination Whole Blood (including riboflavin solution). If the hematocrit of the blood donor or of the collected Whole Blood is used, the blood product will receive more light energy than necessary, and the quality of the Whole Blood may be decreased.

14. Move the pinch clamp on the line between the storage bag and the illumination bag (Line (a)) so it is as near as possible to the storage bag (see Figure 8).
15. Perform the illumination procedure according to the instructions in the Mirasol PRT System for Whole Blood Illuminator operator's manual.
16. After the illumination procedure, remove the illumination/storage set from the Illuminator, hang the illumination bag, and break the frangible connector on the outlet port of the line that connects the illumination bag to the storage bag (Line (a)). See Figure 9.

Note: The broken part of the frangible connector will remain in the illumination bag before and after the Mirasol-treated Whole Blood is transferred to the storage bag.

17. Drain the treated Whole Blood product into the storage bag.



Figure 9: Draining the Whole Blood product

18. Express residual air from the storage bag by performing these steps:


- a. Hold the storage bag so that the ports are facing up, as shown in Figure 10.
- b. Gently apply pressure to the storage bag to move the residual air and Whole Blood up through the storage bag inlet line (Line B) until it just reaches the illumination bag.
- c. Close the pinch clamp on the inlet line to hold the blood in the tubing.

- 1 Illumination bag
- 2 Drain line from the illumination bag to the storage bag
- 3 Pinch clamp
- 4 Storage bag spike port
- 5 Storage bag

Figure 10: Expressing residual air from the storage bag

19. Permanently seal the storage bag inlet line (Line Ⓐ) near the illumination bag, and then disconnect and discard the illumination bag (see Figure 11).

- 1 Disconnected illumination bag
- 2 Drain line from the illumination bag to the storage bag
- 3 Sealed storage bag inlet line
- 4 Pinch clamp
- 5 Storage bag

Figure 11: Sealing and disconnecting the illumination bag

20. Remove the pinch clamp from the storage bag inlet line (Line Ⓐ), and make tubing segments according to your facility's SOP.
21. Label the product and release it to inventory according to your facility's SOP.
- Until the time of transfusion, protect the treated Whole Blood product from exposure to direct sunlight or strong ambient light sources.
 - Follow the Whole Blood product storage specifications listed in Table 3.

Storage Conditions

Table 3: Mirasol-treated Whole Blood storage conditions and quality*

Characteristic	Specification
Time Between Treatment and Storage	Place the treated Whole Blood product into storage within 45 minutes. Until the time of transfusion, protect the treated Whole Blood product from exposure to direct sunlight or strong ambient light sources.
Maximum Storage Period	- Whole Blood units that are held at controlled room temperature or below prior to treatment on the Mirasol PRT system within 24 hours of collection have been verified for storage for up to 14 days at 1 °C to 6 °C. - Whole Blood units that are held at controlled room temperature or below prior to treatment on the Mirasol PRT system within 8 hours of collection have been verified for storage for up to 21 days at 1 °C to 6 °C.
Storage Temperature	1 °C to 6 °C
Transport Temperature	1 °C to 10 °C
Leukoreduction	Optional. Follow your facility's SOP.
*The Mirasol PRT system for Whole Blood has been verified for these product storage conditions. Certain conditions outside of these ranges may be acceptable but should be evaluated by the user as required.	

Mirasol Whole Blood Disposable Kit Specifications

Table 4: Mirasol Whole Blood disposable kit

Characteristic	Specification
Mirasol Whole Blood Disposable Kit	Each disposable kit consists of one illumination/storage set and one 35 mL bag of 500 µM riboflavin in 0.9% sodium chloride solution.
Illumination/storage set	Citrate vinyl bag connected to a PVC universal blood storage bag.
Illumination/storage set sterility	The illumination/storage set has a sterile and non-pyrogenic fluid pathway, sterilized by ethylene oxide (EO).
Illumination/storage set label	The pre-attached label covers a 110 cm ² (17 in ²) surface area of the storage bag. Additional labels may be applied over the pre-attached label. Note: Do not apply any labels to the illumination bag prior to illumination.
Tubing dimensions for sterile connections and seals	Outer diameter = 4.06 mm ± 0.08 mm Wall thickness = 0.56 mm ± 0.05 mm

Table 4: Mirasol Whole Blood disposable kit (continued)



Characteristic	Specification
Riboflavin solution content	500 µM riboflavin in 0.9% sodium chloride solution pH 3.5 to 5.1 Riboflavin Material Safety Data Sheets are available upon request.
Riboflavin solution sterility	The riboflavin solution is sterile and non-pyrogenic, sterilized by steam.
Riboflavin solution volume	35 mL ± 5 mL
Riboflavin solution raw materials	Riboflavin: USP, EP Sodium chloride: USP, EP Hydrochloric acid: USP, EP Water for injection: USP, EP
Disposable kit storage conditions	Long-term storage range: 2 °C to 30 °C (36 °F to 86 °F) Storage humidity: Up to 75% relative humidity (RH), non-condensing Permitted excursions: Up to 60 °C (140 °F) for up to 2 weeks
Case quantity	25 kits per case
Case size	24.9 cm × 35.8 cm × 35.8 cm
Case weight	6.0 kg

Symbols and Certification

Table 5: Symbols and certification

Symbol	Description
	Indicates the product was manufactured in accordance with Annex II of the European Council Directive 93/42/EEC, as amended.
	Indicates the Authorized representative in the European Community.
	Indicates the medical device manufacturer.
	Indicates the package contents.
	Indicates that the user should consult the instructions for use.
	Indicates the product quantity when the quantity is placed in the square.

Table 5: Symbols and certification (continued)

Symbol	Description
	Indicates the manufacturer's catalog number so that the medical device can be identified.
	Indicates the manufacturer's batch code so that the batch or lot can be identified.
	Indicates that the product should not be used if the package is damaged.
	Indicates the date after which the medical device is not to be used.
	Indicates a non-pyrogenic fluid pathway.
	Indicates the presence of a sterile fluid path. The method of sterilization is ethylene oxide (EO).
	Indicates the presence of a sterile fluid path. The method of sterilization is steam.
	Indicates the range of temperatures to which the device can be safely exposed.
	Indicates that the product should not be vented. Do not remove the end caps from the set.
	Indicates that the product contains phthalates, specifically Di(2-ethylhexyl)phthalate (DEHP).
	Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure.
	Indicates a medical device that is not to be resterilized.
	Indicates that the product packaging complies with European Directive 94/62/EC for packaging and packaging waste.
	Indicates the date when the medical device was manufactured (or sterilization date, if the product is sterile).

Return of Used Product

If for any reason the product must be returned to Terumo BCT, Inc., a returned goods authorization (RGA number) is required prior to shipping.

It is the responsibility of the healthcare institution to adequately prepare and identify the product for return shipment.

Shipping instructions, including an RGA number and instructions for cleaning, packaging, and labeling, may be obtained from the Terumo BCT Quality Assurance Department. Please contact your local representative for information regarding returned goods and product complaints.

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This page intentionally left blank.

Terumo BCT, Inc
10811 W. Collins Avenue
Lakewood, Colorado 80215
USA
USA Phone: +1.877.339.4228
Phone: +1.303.231.4357
USA Fax: +1.866.715.6768
Fax: +1.303.542.5215

Terumo BCT Europe N.V
Ikarostraat 41
1930 Zaventem
Belgium
Phone: +32.2.715.05.90
Fax: +32.2.721.07.70

TERUMOBCT.COM

©2015, 2019 Terumo BCT, Inc

Terumo BCT Biotechnologies, LLC
10811 W Collins Avenue
Lakewood, Colorado 80215
USA
USA Phone: +1 877.339.4228
Phone: +1 303.231.4357
USA Fax: +1.866.715.6768
Fax: +1.303.542.5215

2019-01

Part No 777714-015

TERUMO BCT

Addition to IFU

Mirasol® Whole Blood Disposable Set

2021

Addition to IFU

Table of content

1	NAME OF MEDICAL DEVICE	3
1.1	CONTENTS.....	3
2	INTENDED USE	3
3	POSSIBLE SIDE EFFECTS	3
4	OPERATION CONDITIONS	4
5	INFORMATION ON MANUFACTURERS OF THE MEDICAL DEVICE.....	4
6	MEDICAL DEVICE CLASSIFICATION	5
6.1	CLASS OF ELECTRICAL SAFETY, LASER SAFETY AND OTHER TYPES OF DANGER.....	5
6.2	TYPE OF CONTACT WITH A HUMAN BODY	5
7	PRINCIPLE OF MEDICAL DEVICE OPERATION	5
7.1	RIBOFLAVIN.....	5
7.2	OPERATION AS PART OF THE SYSTEM	6
7.3	SCOPE OF OPERATION	8
8	OPPORTUNITIES AND WAYS OF INTEGRATION WITH OTHER MEDICAL DEVICES.....	8
9	PACKAGE AND LABELING.....	8
9.1	INFORMATION OF MARKING AND PACKAGING	8
9.2	PACKAGING.....	10
10	MAIN TECHNICAL CHARACTERISTICS / PARAMETERS / PROPERTIES OF A PRODUCT AND ACCESSORIES, WHICH ARE GUARANTEED BY THE MANUFACTURER (INCLUDING MATERIALS OF PRODUCTS).....	10
10.1	MATERIALS OF WHICH THE MEDICAL DEVICE OR ITS MAIN COMPONENTS ARE MADE	10
10.2	LIST OF MATERIALS OF ANIMAL AND (OR) HUMAN ORIGIN	11
10.3	LIST OF DRUG FOR MEDICAL USE CONTAINED IN THE MEDICAL DEVICE.....	11
10.4	RATE TECHNICAL CHARACTERISTICS	11
11	TRANSPORT AND STORAGE CONDITIONS	13
12	SHELF LIFE	13
13	REQUIREMENTS TO ENVIRONMENT PROTECTION.....	14
13.1	PROCEDURE AND CONDITIONS FOR THE DISPOSAL OF A MEDICAL DEVICE	14
14	METHODS AND MEANS OF DISINFECTION AND PRESTERILIZING CLEANING.....	14
15	STERILITY AND STERILIZATION METHOD	14
16	THE LIST OF NATIONAL AND INTERNATIONAL REGULATORY DOCUMENTS, STANDARDS TO WHICH THE MEDICAL DEVICE COMPLIES	14
17	REQUIREMENTS FOR MAINTENANCE AND REPAIR OF THE MEDICAL DEVICE	15
18	WARRANTY LIABILITY	15
19	CLAIM.....	16

1 NAME OF MEDICAL DEVICE

Mirasol® Whole Blood Disposable Set

Further in the text defines as «Mirasol WB», «Mirasol kit», «Mirasol set», «kit», «set», «medical device».

1.1 CONTENTS

Mirasol® Whole Blood Disposable Set:

- I. 500µM Riboflavin in 0.9% Sodium Chloride solution bag, 35mL
- II. Whole blood Illumination/storage set:
 - 1. Illumination bag
 - 2. Storage bag

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- 1. Local reactions related to needle insertion:
 - a. Vessel injuries: haematomas, arterial puncture, thrombophlebitis.
 - b. Nerve injuries: injury of nerve, injury of nerve by haematoma.
 - c. Other complications: tendon injury, allergic reaction (local), infection (local)
- 2. General reactions:
 - a. Vasovagal reaction: immediate type, delayed type.
- 3. Rare, significant complications :
 - a. Related to vessel injury: brachial artery pseudoaneurysm, arteriovenous fistula, compartment syndrome, axillary vein thrombosis.
 - b. Accidents: accidents or injuries related to vasovagal syncope, other kinds of accidents.
 - c. Cardiovascular reactions: angina pectoris, myocardial infarction, cerebral ischaemia.

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The Mirasol® PRT System (Russian registration certificate № ФС3 2009/05465), with which the set is applied, uses a unique process to photochemically inhibit the replication of DNA and RNA to reduce the pathogen load and inactivate WBC. The ability of riboflavin to act as a photosensitizer results in selective damage to nucleic acids upon exposure to UV light, without binding to cells and proteins. The therapeutic benefit of Mirasol treatment derives from the blood product itself, not from the Mirasol treatment. There are no intended pharmacologic functions directly associated with the device or its components. The Mirasol PRT System increases the safety of using blood products.

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The Mirasol PRT System uses riboflavin, as a photosensitizer, at a final concentration of approximately 50 µmol/L. Riboflavin associates with genomic nucleic acids and mediates an oxygen-independent electron transfer process upon exposure to UV light. There are no intended pharmacologic functions associated with the device or its components. Riboflavin does not act as a medicinal product in this application.

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The Mirasol system uses riboflavin at a final concentration of approximately 50 $\mu\text{mol/l}$, in combination with UV light from a fluorescent lamp with the energy centered at 313 nm. The majority (99%) of the total energy reaching the product through the plastic of the Illumination bag is in the UVB (280–315 nm) and UVA (315–400 nm) spectral regions. The peak wavelength chosen for the Mirasol process preferentially targets riboflavin associated with nucleic acids, directing the damage specifically to these nucleic acids. This process involves an oxygen-independent electron transfer leading to the modification of nucleic acids, primarily on guanine residues, and the conversion of riboflavin to its primary photoproduct lumichrome. Compared to the use of UV light alone, which causes *reversible* nucleic acid damage, damage induced by riboflavin is irreversible since replication and repair processes are impaired due to the guanine base modification. The numbers of lesions occur at a frequency of approximately one in every 350 base pairs.

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- Equipment: Mirasol Illuminator and tool for carrying out a light calibration check (Radiometer). The Illuminator delivers light and agitation for the Mirasol process.
- Disposables: A disposable kit is available for treatment of Whole Blood

- Mirasol Whole Blood Disposable Kit
- **Software:** An accessory data management software (Mirasol Manager is available that interfaces with one or multiple Illuminators to capture and manage data and produce reports. This software also includes a data export tool.

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Optional Mirasol Manager System Components:

- Barcode reader (distributed accessories)
- Electronic scale(distributed accessories)
- Report Printer (Not distributed by Terumo BCT)
- Process Label Printer (Not distributed by Terumo BCT)

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The Mirasol® Whole Blood is directly compatible with the Mirasol® Pathogen Reduction Technology (PRT) System (RC № ФСЗ 2009/05465). There are no other known compatible medical devices.

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



Devices supplied in packing marked with the following information:

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Symbols:

Symbol	Description
	Indicates that the product was manufactured in accordance with Annex II of the European Council Directive, 93/42/EEC, as amended.
	Indicates the European Authorized Representative of the product when accompanied by the name of the European Authorized Representative.
	Indicates the manufacturer of the product when accompanied by the name of the manufacturer
	Date of manufacture
	Indicates the package contents
	Indicates that the user should consult the instructions for use.
	Indicates the product quantity when the quantity is placed in the square.
	Indicates the product catalog number when accompanied by a number.
	Indicates the lot of the product when accompanied by a number.
	Indicates that the user should not use the product if the package is damaged.
	Indicates the expiration date of the product when accompanied by a specific date.
	Indicates a non-pyrogenic fluid pathway.
	Indicates a sterile fluid pathway using ethylene oxide (EO).
	Indicates a sterile fluid pathway using steam or dry heat.
	Indicates the temperature limitations for product storage when accompanied by temperature(s).
	Do not vent the product. Do not remove the end caps from the set.
	Indicates that the product contains phthalates, specifically Di(2-ethylhexyl)phthalate (DEHP).
	Indicates that the product should not be reused.

	Indicates a medical device should not be re-sterilized.
	Indicates that the product packaging complies with European Directive 94/62/EC for packaging and packaging waste.
	Time & temperature sensitive
	Keep away from sunlight
	Handle with care
	This side up
	Fragile
	Keep dry

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	<u>Whole blood Illumination/storage set:</u> <u>Storage bag:</u> PVC Container holder: PP <u>Illumination bag:</u> PVC Polycarbonate <u>Tubes:</u> PVC <u>Y-type connector:</u> PETG <u>White clamp:</u> PP/Polyalkylene oxide block copolymer colored white <u>Blue cap:</u> PVC colored blue
--	---

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General Specifications	
Bag Capacity	• Illumination bag: maximum volume 613 mL (including riboflavin solution) • Storage bag: maximum volume 613 mL (including riboflavin solution) • Riboflavin solution bag: maximum volume 50 mL (Riboflavin solution volume 35 mL $\pm$ 5 mL)
Illumination bag dimensions	• Length: 239 mm • Width: 177 mm
Storage bag dimensions	▪ Length: 238 mm ▪ Width: 130 mm
Riboflavin Solution bag	▪ Length: 98 mm ▪ Width: 75 mm

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Transportation conditions

Parameter	Value
Humidity	35% to 75% relative humidity (RH), non-condensing Permitted humidity excursions: – from 75% ($\pm 5\%$) up to 72 hours – 15% to 35% for up to 6 months
Transport temperature	2 – 30 °C (protect from freezing), with permitted excursions up to 60 °C for up to 2 weeks Permitted temperature excursions: 1 °C to 10 °C

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13.1 Procedure and conditions for the disposal of a medical device

The Mirasol® WB Set should be disposed of in accordance with local regulations:

Class A of medical waste in accordance with SanPin 2.1.3684-21 for unused medical device

Class B of medical waste in accordance with SanPin 2.1.3684-21 for used medical device.

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List of applicable standards used:

- Council Directive 93/42/EEC
- Council Directive 94/62/EC concerning Packaging and Packaging Waste Directive
- European Commission Regulation 722/2012
- Guidelines on Medical Devices MEDDEV.2.7.1 Rev.4, June 2016 Clinical Evaluation: Guide for Manufacturers and Notified Bodies Under Directives 93/42/EEC and 90/385/EEC
- ICH Q1A (R2):2003 Stability Testing of New Drug Substances and Products, Section 2.2.7.3 Guidelines for Semi-Permeable Containers, Climatic Zones I-II.
- ICH-Q1B:1996 Photostability Testing of New Drug Substances and Products.
- EN ISO 3826-1:2017 Plastics Collapsible Containers for Human Blood and Blood Components - Part 1: Conventional containers.
- EN ISO 3826-3:2013 Plastics Collapsible Containers for Human Blood and Blood Components - Part 3: Blood bag systems with integrated features.
- EN ISO 13485:2012/AC:2012 Medical Devices -- Quality Management Systems -- Requirements for Regulatory Purposes
- EN ISO 14971:2012 Medical devices – Application of risk management to medical devices.
- ISO 80369-7: 2016, Small-bore connectors for liquids and gases in Healthcare applications - Part 7: Connectors for intravascular or hypodermic applications

GHTF/SG5/N2R8:2007

- EN ISO 11135-1:201 Medical devices - Validation and routine control of ethylene oxide sterilization. Half Cycle Method (Method C).
- EN ISO 17665-1:2006 Sterilization of health care products – Moist Heat – Part 1: Requirements of the development, validation and routine control of a sterilization process for medical devices.
- EN 556-1:2001, COR 2006, Sterilization of medical devices – requirements for medical devices to be designated “sterile” – Part 1: Requirements for Terminally sterilized medical devices.
- ISO 14644-1:2015, Cleanrooms and associated controlled environments – Part 1: Classification of air cleanliness by particle concentration
- ISO 14644-2:2015, Cleanrooms and associated controlled environments – Part 2: Specifications for testing and monitoring to prove continued compliance with ISO 14644-1
- EN ISO 14698-1:2014, Clean Rooms and Associated Controlled Environments - Part 1: Biocontamination Control – General Principles and Methods
- EN 1041:2008, A1:2013 Information supplied by the manufacturer with medical devices.
- EN ISO 15223-1:2016 Medical devices – Symbols to be used with medical device labels, labeling and information to be supplied. Part 1: General requirements.
- EN 15986: 2011 Symbol for use in the labeling of medical devices – requirements for labeling of medical devices containing phthalates.
- ISTA 2A:2011 Partial Simulation Performance Test; Packaged Products 150 lbs. or less.
- ISTA 3E:2017 Partial Simulation Performance Test; Unified loads of Same Product
- EN ISO 11607-1:2017 Packaging for terminally sterilized medical devices.
- European Pharmacopeia, 98th Edition, 2017
- USP 37/NF 32 (2014) United States Pharmacopeia
- ANSI/AAMI ST72: 2011 Bacterial endotoxin - Test methods, routine monitoring and alternatives to batch testing.
- EN ISO 10993-1:2018 Biological Evaluation of Medical Devices -- Part 1: Evaluation and testing.
- EN ISO 10993-7:2008, COR 2009 Biological evaluation of medical devices -- Part 7: Ethylene oxide sterilization residuals.
- EN ISO 15747:2018 Plastic containers for intravenous injections.

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LIMITED WARRANTY. Warranty information applicable to the equipment or disposables specified on this invoice is set forth in the price list and any contract for such equipment or disposables and shall be provided by Terumo BCT upon request. For a period of ninety (90) days from the date of the invoice, Terumo BCT will repair or replace Terumo BCT spare parts or optional feature components that Terumo BCT finds defective in material or workmanship. Terumo BCT is not responsible for, and shall have no obligation with respect to any failure caused by normal wear and tear, misuse, unauthorized alterations, accident, neglect, use of non-standard accessory attachments, and/or improper customer maintenance. Terumo BCT does not verify the safety or efficacy of, and makes no warranties, express or implied, with respect to any non-Terumo BCT components or non-Terumo BCT spare parts included at Purchaser's request or any components or products used other than as expressly intended by their manufacturer. Failure to perform preventive maintenance procedures will invalidate this warranty. Terumo BCT is not responsible for any events related to non-Terumo BCT parts installed in the equipment. **THIS LIMITED WARRANTY CONTAINS THE PURCHASER'S EXCLUSIVE WARRANTIES AND REMEDIES.**

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Representative office in Russia:

Terumo Rus LLC

123112, Moscow, Testovskaya st, 10, floor/flat/room 13/1/5

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This Certificate is not valid for use anywhere within the United States of America, its territories or possessions.

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10. Signature:
Signature | Firma

Jena Griswold

This Apostille only certifies the authenticity of the signature and the capacity of the person who has signed the public document, and, where appropriate, the identity of the seal or stamp which the public document bears. This Apostille does not certify the content of the document for which it was issued. This Apostille is not valid for use anywhere within the United States of America, its territories or possessions. To verify the issuance of this Apostille, see: www.sos.state.co.us/auth

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
Best regards,

Corrine Lucas
Director Regulatory Affairs,
Compliance and Systems

Subscribed and sworn to before me this 21st (day) in the County of Jefferson, State of Colorado, USA.
01 of February 2022

Seal:

Notary public

Terumo BCT, Inc.
10811 West Collins Ave.
Lakewood, Colorado 80215-4440
USA
USA Phone: 1.877.339.4228
Phone: +1.303.231.4357
Fax: +1.303.542.5215

TERUMOBCT.COM

Terumo BCT Europe N.V.
Europe, Middle East and Africa
Barossaan 41
1930 Zaventem
Belgium
Phone: +32.2.715.0590
Fax: +32.2.721.0770

Terumo BCT Asia Pte. Ltd.
89 Science Park Drive
#04-25 (Lobby B)
The Rutherford
Singapore 118251
Phone: +65.6715.3778
Fax: +65.6774.1419

Terumo BCT Latin America S.A.
La Pampa 1517-12th Floor
C1428QZE
Buenos Aires
Argentina
Phone: +54.11.5530.5200
Fax: +54.11.5530.5201

Terumo BCT Japan, Inc.
Tokyo Opera City Tower 49F,
3-20-2, Nishi-Shinjuku,
Shinjuku-ku, Tokyo 163-1480,
Japan
Phone: +81.3.6743.7899
Fax: +81.3.6743.9850

001-XXXX-XXXX

Addition to IFU

Table of content

1	NAME OF MEDICAL DEVICE.....	3
1.1	CONTENTS.....	3
2	INTENDED USE	3
3	POSSIBLE SIDE EFFECTS	3
4	OPERATION CONDITIONS	4
5	INFORMATION ON MANUFACTURERS OF THE MEDICAL DEVICE.....	4
6	MEDICAL DEVICE CLASSIFICATION	5
6.1	CLASS OF ELECTRICAL SAFETY, LASER SAFETY AND OTHER TYPES OF DANGER	5
6.2	TYPE OF CONTACT WITH A HUMAN BODY	5
7	PRINCIPLE OF MEDICAL DEVICE OPERATION.....	5
7.1	RIBOFLAVIN	5
7.2	OPERATION AS PART OF THE SYSTEM.....	6
7.3	SCOPE OF OPERATION	8
8	OPPORTUNITIES AND WAYS OF INTEGRATION WITH OTHER MEDICAL DEVICES.....	8
9	PACKAGE AND LABELING	8
9.1	INFORMATION OF MARKING AND PACKAGING.....	8
9.2	PACKAGING	10
10	MAIN TECHNICAL CHARACTERISTICS / PARAMETERS / PROPERTIES OF A PRODUCT AND ACCESSORIES, WHICH ARE GUARANTEED BY THE MANUFACTURER (INCLUDING MATERIALS OF PRODUCTS).	10
10.1	MATERIALS OF WHICH THE MEDICAL DEVICE OR ITS MAIN COMPONENTS ARE MADE.....	10
10.2	LIST OF MATERIALS OF ANIMAL AND (OR) HUMAN ORIGIN	11
10.3	LIST OF DRUG FOR MEDICAL USE CONTAINED IN THE MEDICAL DEVICE.	11
10.4	RATE TECHNICAL CHARACTERISTICS	11
11	TRANSPORT AND STORAGE CONDITIONS	13
12	SHELF LIFE	13
13	REQUIREMENTS TO ENVIRONMENT PROTECTION	14
13.1	PROCEDURE AND CONDITIONS FOR THE DISPOSAL OF A MEDICAL DEVICE.....	14
14	METHODS AND MEANS OF DISINFECTION AND PRESTERILIZING CLEANING	14
15	STERILITY AND STERILIZATION METHOD	14
16	THE LIST OF NATIONAL AND INTERNATIONAL REGULATORY DOCUMENTS, STANDARDS TO WHICH THE MEDICAL DEVICE COMPLIES	14
17	REQUIREMENTS FOR MAINTENANCE AND REPAIR OF THE MEDICAL DEVICE	15
18	WARRANTY LIABILITY	15
19	CLAIM	16

1 NAME OF MEDICAL DEVICE

Mirasol® Whole Blood Disposable Set

Further in the text defines as «Mirasol WB», «Mirasol kit», «Mirasol set», «kit», «set», «medical device».

1.1 CONTENTS

Mirasol® Whole Blood Disposable Set:

- I. 500µM Riboflavin in 0.9% Sodium Chloride solution bag, 35mL
- II. Whole blood Illumination/storage set:
 1. Illumination bag
 2. Storage bag

Further in the text the bag defines as «bag», «pouch».

2 INTENDED USE

Intended for extracorporeal processing of blood using the energy of ultraviolet light and riboflavin (vitamin B2), carried out in order to reduce the pathogenic load and inactivate leukocytes in whole blood intended for transfusion.

3 POSSIBLE SIDE EFFECTS

Known adverse reactions related to blood collection

1. Local reactions related to needle insertion:
 - a. Vessel injuries: haematomas, arterial puncture, thrombophlebitis.
 - b. Nerve injuries: injury of nerve, injury of nerve by haematoma.
 - c. Other complications: tendon injury, allergic reaction (local), infection (local)
2. General reactions:
 - a. Vasovagal reaction: immediate type, delayed type.
3. Rare, significant complications:
 - a. Related to vessel injury: brachial artery pseudoaneurysm, arteriovenous fistula, compartment syndrome, axillary vein thrombosis.
 - b. Accidents: accidents or injuries related to vasovagal syncope, other kinds of accidents.
 - c. Cardiovascular reactions: angina pectoris, myocardial infarction, cerebral ischaemia.

Source: Guide to the preparation, use and quality assurance of blood components (Recommendation No. R (95) 15, 17th Edition) published by Directorate for the Quality of Medicines & HealthCare of the Council of Europe (EDQM).

4 OPERATION CONDITIONS

The collection of whole blood using the Mirasol Set typically occurs at room temperature. The processing of whole blood occurs in a laboratory setting in further environmental conditions:

Specification Type	Value
Relative Humidity	10 to 80%, non-condensing
Operating temperature	18 °C to 27 °C
Altitude	Usually < 2000m

The patient's body temperature is 32° C – 42° C, hence the operating conditions are:

Parameter	Value
Humidity	Up to 100%
Operating temperature	From 32 °C to 42 °C
Pressure	700 – 1060 gPa

5 INFORMATION ON MANUFACTURERS OF THE MEDICAL DEVICE

Legal Manufacturer

Terumo BCT, Inc.
10811 W. Collins Avenue Lakewood, Colorado 80215 USA Tel: +1 303 231 4357
Fax: +1 303 542 5215

Developer

Terumo BCT, Inc.
10811 W. Collins Avenue Lakewood, Colorado 80215 USA Tel: +1 303 231 4357
Fax: +1 303 542 5215

Manufacturing Facility

Terumo BCT, Inc.
10811 W. Collins Avenue Lakewood, Colorado 80215 USA Tel: +1 303 231 4357
Fax: +1 303 542 5215

6 MEDICAL DEVICE CLASSIFICATION

Mirasol® Whole Blood Disposable Set is a CE-marked class IIb medical device complying with the Medical Device Directive (93/42/EEC) as amended

Potential risk class - 3

Sterility - sterile device

Invasiveness - non-invasive device

Frequency rate of application: single-use device

6.1 Class of electrical safety, laser safety and other types of danger

Not applicable

6.2 Type of contact with a human body

500µM riboflavin solution in 0.9% sodium chloride solution, 35 ml; Whole Blood Storage

Container - Long-term contact (24 hours to 30 days) with blood.

Irradiation container – Short-term contact (less than 24 hours) with blood.

Packaging - Short-term indirect contact (less than 24 hours) with the body

7 PRINCIPLE OF MEDICAL DEVICE OPERATION

The Mirasol® PRT System (Russian registration certificate № ФЦ3 2009/05465), with which the set is applied, uses a unique process to photochemically inhibit the replication of DNA and RNA to reduce the pathogen load and inactivate WBC. The ability of riboflavin to act as a photosensitizer results in selective damage to nucleic acids upon exposure to UV light, without binding to cells and proteins. The therapeutic benefit of Mirasol treatment derives from the blood product itself, not from the Mirasol treatment. There are no intended pharmacologic functions directly associated with the device or its components. The Mirasol PRT System increases the safety of using blood products.

The simple, effective process involves transferring the platelet concentrate to the Illumination/Storage bag and adding 35 ± 5 mL of riboflavin solution (containing 6.5 mg of riboflavin). The product is then placed into the Mirasol illuminator where UV light is delivered, depending on the blood product type and weight, to reach the target energy dose.

7.1 Riboflavin

Riboflavin, or vitamin B2, occurs naturally in foodstuffs such as cheese, almonds, and red meat and has a long history of safe use as a nutritional supplement and a permitted food coloring. Many studies and data from literature demonstrate the safety of riboflavin.

The Mirasol PRT System uses riboflavin, as a photosensitizer, at a final concentration of approximately 50 µmol/L. Riboflavin associates with genomic nucleic acids and mediates an oxygen-independent electron transfer process upon exposure to UV light. There are no intended pharmacologic functions associated with the device or its components. Riboflavin does not act as a medicinal product in this application

Compared with the use of UV light alone, which mainly causes reversible nucleic acid dimerization, damage induced by the combination of UV light and riboflavin is irreversible since the replication and repair processes are impaired due to the guanine base modifications. Thereby DNA and RNA replication are inhibited.

7.2 Operation as part of the system

The Mirasol PRT System was developed to improve safety of transfused blood and blood components by reducing the possibility of a transfusion transmitted infection (TTI) and of immune reactions resulting from transfusion. Nucleic acid bound riboflavin molecules are preferentially activated by UV light. The riboflavin is reduced upon activation, thereby modifying the nucleic acids (primarily guanine) and impairing replication. The therapeutic benefit of Mirasol treatment derives from the blood product itself, not from the Mirasol treatment. There are no intended pharmacologic functions directly associated with the device or its components.

The Mirasol PRT System for Whole Blood uses riboflavin (vitamin B2) and UV light to induce damage to nucleic acids that are present in viruses, bacteria, parasites and white blood cells. The ability of riboflavin to act as a photosensitizer results in selective damage to nucleic acids upon exposure to light, without binding to cells and proteins.

The Mirasol system uses riboflavin at a final concentration of approximately 50 µmol/L, in combination with UV light from a fluorescent lamp with the energy centered at 313 nm. The majority (99%) of the total energy reaching the product through the plastic of the Illumination bag is in the UVB (280–315 nm) and UVA (315–400 nm) spectral regions. The peak wavelength chosen for the Mirasol process preferentially targets riboflavin associated with nucleic acids, directing the damage specifically to these nucleic acids. This process involves an oxygen-independent electron transfer leading to the modification of nucleic acids, primarily on guanine residues, and the conversion of riboflavin to its primary photoproduct lumichrome. Compared to the use of UV light alone, which causes reversible nucleic acid damage, damage induced by riboflavin is irreversible since replication and repair processes are impaired due to the guanine base modification. The numbers of lesions occur at a frequency of approximately one in every 350 base pairs.

Nucleic acid damage may be less frequent in mitochondrial DNA. The light source used by the process does not emit light in spectral regions where cytochrome C, flavin mononucleotides (FMN) and flavin adenine dinucleotides (FAD) absorb, explaining why platelet mitochondria retain their normal structure and function following Mirasol treatment and subsequent storage.

The treatment procedure involves mixing a blood product with Riboflavin Solution in the Illumination/Storage Bag and illuminating the mixture with a calculated light energy dose in the Illuminator. The Mirasol System is comprised of three main elements

- **Equipment:** Mirasol Illuminator and tool for carrying out a light calibration check (Radiometer). The Illuminator delivers light and agitation for the Mirasol process.
- **Disposables:** A disposable kit is available for treatment of Whole Blood

- Mirasol Whole Blood Disposable Kit
- **Software** An accessory data management software (Mirasol Manager is available that interfaces with one or multiple Illuminators to capture and manage data and produce reports. This software also includes a data export tool.

Depending on the Mirasol treatment, a specific Mirasol Disposable Tubing Kit is used.

Additional details for the various system components follow below. Optional accessories that are also distributed by Terumo BCT and that are intended to facilitate the Mirasol system's operation in routine include a barcode reader and an electronic scale. An example of a full system set-up (including standard and optional components) is shown in Figure 1.1.


Figure 1.1.

Mirasol System Standard Components:

- Mirasol Illuminator
- Mirasol Manager-Integrated Data Management and Reporting Software
- Mirasol Whole Blood Disposable Set

Optional Mirasol Manager System Components:

- Barcode reader (distributed accessories)
- Electronic scale (distributed accessories)
- Report Printer (Not distributed by Terumo BCT)
- Process Label Printer (Not distributed by Terumo BCT)

7.3 Scope of operation

Processing of whole blood, component storage and transfusion after processing.


Figure 1.2.

8 OPPORTUNITIES AND WAYS OF INTEGRATION WITH OTHER MEDICAL DEVICES

The Mirasol® Whole Blood is directly compatible with the Mirasol® Pathogen Reduction Technology (PRT) System (RC № ΦC3 2009/05465). There are no other known compatible medical devices.

The Mirasol® Whole Blood Set is validated for use with Mirasol PRT System. It allows the processing (illumination) of whole blood, after addition of riboflavin solution, and storage of the whole blood in the attached blood bags after processing. The Mirasol Illuminator is intended to provide light energy as part of the Mirasol System.

9 PACKAGE AND LABELING

9.1 Information of marking and packaging

Labels supplied in packing marked with the following information:

- device name and reference number;

- LOT of the product;
- logo and address of the manufacturer;
- Information about sterile;
- Number of products;
- expiry date;
- and other information and symbols.

Symbols:

Symbol	Description
	Indicates that the product was manufactured in accordance with Annex II of the European Council Directive, 93/42/EEC, as amended.
	Indicates the European Authorized Representative of the product when accompanied by the name of the European Authorized Representative
	Indicates the manufacturer of the product when accompanied by the name of the manufacturer
	Date of manufacture
	Indicates the package contents
	Indicates that the user should consult the instructions for use.
	Indicates the product quantity when the quantity is placed in the square.
	Indicates the product catalog number when accompanied by a number.
	Indicates the lot of the product when accompanied by a number.
	Indicates that the user should not use the product if the package is damaged.
	Indicates the expiration date of the product when accompanied by a specific date.
	Indicates a non-pyrogenic fluid pathway.
	Indicates a sterile fluid pathway using ethylene oxide (EO).
	Indicates a sterile fluid pathway using steam or dry heat
	Indicates the temperature limitations for product storage when accompanied by temperature(s).
	Do not vent the product. Do not remove the end caps from the set.
	Indicates that the product contains phthalates, specifically Di(2-ethylhexyl)phthalate (DEHP).
	Indicates that the product should not be reused

	Indicates a medical device should not be re-sterilized
	Indicates that the product packaging complies with European Directive 94/62/EC for packaging and packaging waste.
	Time & temperature sensitive
	Keep away from sunlight
	Handle with care
	This side up
	Fragile
	Keep dry

9.2 Packaging

Each group carton contains 25 Mirasol WB sets per case. The shipping carton also contains one instruction for use.

Group packaging is put into a shipping box, while all sides are protected by foam plastic and refrigerant.

10 MAIN TECHNICAL CHARACTERISTICS / PARAMETERS / PROPERTIES OF A PRODUCT AND ACCESSORIES, WHICH ARE GUARANTEED BY THE MANUFACTURER (INCLUDING MATERIALS OF PRODUCTS).

10.1 Materials of which the medical device or its main components are made

Name	Materials
Mirasol Whole Blood Disposable Set:	500µM Riboflavin in 0.9% Sodium Chloride solution bag, 35mL.
I. 500µM Riboflavin in 0.9% Sodium Chloride solution bag, 35mL	Bag: PVC
II. Whole blood	Riboflavin solution:
Illumination/storage set:	<1% - Riboflavin 6,58 mg (CAS№ 83-88-5)
1. Illumination bag	<1% - Sodium Chloride 315 mg (CAS№ 7647-14-5)
2. Storage bag	<1% - Hydrochloric acid (CAS№ 7647-01-0)
	>97% - Injection water 35 ml (CAS№ 7732-18-5)

	<u>Whole blood Illumination/storage set</u> <u>Storage bag:</u> PVC Container holder: PP <u>Illumination bag:</u> PVC Polycarbonate <u>Tubes:</u> PVC <u>Y-type connector:</u> PETG <u>White clamp:</u> PP/Polyalkylene oxide block copolymer colored white <u>Blue cap:</u> PVC colored blue
--	--

10.2 List of materials of animal and (or) human origin

Not applicable

10.3 List of drug for medical use contained in the medical device.

The device contains a drug «Riboflavin» that is registered in the Russian Federation (Registration number: JICP-002944/07 from 01.10.2007)

10.4 Rate technical characteristics

The permissible deviation of the parameters is $\pm 10\%$, unless otherwise specified.

General Specifications	
Bag Capacity	- Illumination bag: maximum volume 613 mL (including riboflavin solution) - Storage bag: maximum volume 613 mL (including riboflavin solution) - Riboflavin solution bag: maximum volume 50 mL (Riboflavin solution volume 35 mL $\pm$ 5 mL.)
Illumination bag dimensions	- Length: 239 mm - Width: 177 mm
Storage bag dimensions	- Length: 238 mm - Width: 130 mm
Riboflavin Solution bag	- Length: 98 mm - Width: 75 mm

Tensile properties (extension strength) of connections between bags and tubings	20 N applied to the tubing for at least 15 seconds The tensile force shall be applied at a right angle (90°) to the connection as well as along the longitudinal axis of the tubing at a temperature of 23 ± 5 °C
Tubing	- Outer diameter of the tubes $4,06 \pm 0,8$ mm - Wall thickness $0,56 \text{ mm} \pm 0,05$ mm
Inlet line to the illumination bag with a luer and a blue end cap	- Luer type Male - Blue cap – Luer type female
Y connector dimensions (mm)	- Outer diameter of the tubes: 6 mm - Thickness of the tube: 0,75 mm - Length: 21,50 mm
Spike port dimensions (mm)	- Outer diameter of the tubes: 6,7 mm - Thickness of the tube: 6,90 mm - Length: 35,25 mm
Storage bag spike port dimensions (mm)	- Outer diameter at base: 8 mm - Outer diameter at tube connection: 7,5 mm - Length: 18 mm
Riboflavin solution bag spike port dimensions (mm)	- Length: 30,63 mm - Width: 11,5 mm - Thickness: 2,76 mm
Pinch clamp dimensions in closed position, (mm)	- Length: 24,5 mm - Width: 12,75 mm - Height: 10,45 mm
Luer cap dimensions (mm)	- Diameter: 10,75 mm - Height: 9,75 mm
Weight of the Luer cap, g	0,41 g
The diameter of the fitting on the bags	- Min 6,64 mm - Max 6,75 mm
Vent line with a blue end cap on the storage bag dimensions, (mm)	- Length: 135,63 - Outer diameter: 6,50 - Inner diameter: 4,20
Frangible connector dimensions, (mm)	- Outer diameter of the tubes: 9,25 mm - Thickness of the tube: 11,27 mm - Length: 41,5 mm
Bag holder	- Length: 260 mm - Width: 180 mm
Mass	- Illumination bag: 41,97 g - Storage bag: 29,18 g - Riboflavin Solution bag: 48,96 g - Mirasol Whole Blood Disposable Set net weight: 129 g

Tubing specifications of Mirasol WB

Item	Wall thickness (mm)	Outer diameter (mm)	Inner diameter (mm)
Inlet lines (1, 2, 3, 4)	0,56 mm	4,06 mm	2,94 mm
Riboflavin solution bag tubing (14)	0,56 mm	4,06 mm	2,94 mm

11 TRANSPORT AND STORAGE CONDITIONS

Storage conditions of Mirasol-treated whole blood

Use blood bag sets within 28 days after you open the outer aluminum foil pouch. To store unused blood bag sets, return them to the outer aluminum foil pouch and re-close the pouch with tape or a clip. Once you open the transparent packing wrap, you must use the blood bag set within 7 days, not exceeding 28 days from when you opened the outer aluminum foil pouch.

- Storage humidity: 35% to 75% relative humidity (RH), non-condensing
- Protect from freezing. Avoid excessive heat.

Storage conditions of the disposable kits

Shipping and storage conditions of the Mirasol Whole Blood disposable kits are:

- Long-term storage range: 2 – 30 °C (protect from freezing), with permitted excursions up to 60 °C for up to 2 weeks;
- Relative humidity conditions of up to 75% (non-condensing).

Transportation conditions

Parameter	Value
Humidity	35% to 75% relative humidity (RH), non-condensing Permitted humidity excursions: – from 75% (± 5%) up to 72 hours – 15% to 35% for up to 6 months
Transport temperature	2 – 30 °C (protect from freezing), with permitted excursions up to 60 °C for up to 2 weeks Permitted temperature excursions: 1 °C to 10 °C

Avoid excessive heat (> 40 °C). Protect from freezing (< 0 °C).

12 SHELF LIFE

The product expires 3 years from date of manufacture. Refer to product labelling for expiry date.

13 REQUIREMENTS TO ENVIRONMENT PROTECTION

Medical device does not influence on the environment.

13.1 Procedure and conditions for the disposal of medical device

The Mirasol® WB Set should be disposed of in accordance with local regulations:
Class A of medical waste in accordance with SanPin 2.1.3684-21 for unused medical device
Class B of medical waste in accordance with SanPin 2.1.3684-21 for used medical device.

14 METHODS AND MEANS OF DISINFECTION AND PRESTERILIZING CLEANING

Not applicable. The Mirasol® WB is provided sterile and does not need additional sterilization.

15 STERILITY AND STERILIZATION METHOD

The riboflavin solution is steam sterilized. The irradiation / storage kit is sterilized with ethylene oxide, including a sterile, pyrogen-free fluid channel.

16 THE LIST OF NATIONAL AND INTERNATIONAL REGULATORY DOCUMENTS, STANDARDS TO WHICH THE MEDICAL DEVICE COMPLIES

List of applicable standards used:

- Council Directive 93/42/EEC
- Council Directive 94/62/EC concerning Packaging and Packaging Waste Directive
- European Commission Regulation 722/2012
- Guidelines on Medical Devices MEDDEV.2.7.1 Rev.4, June 2016 Clinical Evaluation: Guide for Manufacturers and Notified Bodies Under Directives 93/42/EEC and 90/385/EEC
- ICH Q1A (R2):2003 Stability Testing of New Drug Substances and Products, Section 2.2.7.3 Guidelines for Semi-Permeable Containers, Climatic Zones I-II.
- ICH-Q1B:1996 Photostability Testing of New Drug Substances and Products.
- EN ISO 3826-1:2017 Plastics Collapsible Containers for Human Blood and Blood Components - Part 1: Conventional containers.
- EN ISO 3826-3:2013 Plastics Collapsible Containers for Human Blood and Blood Components - Part 3: Blood bag systems with integrated features.
- EN ISO 13485:2012/AC:2012 Medical Devices -- Quality Management Systems -- Requirements for Regulatory Purposes
- EN ISO 14971:2012 Medical devices -- Application of risk management to medical devices.
- ISO 80369-7: 2016, Small-bore connectors for liquids and gases in Healthcare applications - Part 7: Connectors for intravascular or hypodermic applications.

GHTE/SG5/N2R8:2007

- EN ISO 11135-1:2011 Medical devices – Validation and routine control of ethylene oxide sterilization, Half Cycle Method (Method C).
- EN ISO 17665-1:2006 Sterilization of health care products – Moist Heat – Part 1: Requirements of the development, validation and routine control of a sterilization process for medical devices.
- EN 556-1:2001, COR 2006 Sterilization of medical devices – requirements for medical devices to be designated “sterile” – Part 1: Requirements for Terminally sterilized medical devices.
- ISO 14644-1:2015, Cleanrooms and associated controlled environments – Part 1: Classification of air cleanliness by particle concentration
- ISO 14644-2:2015, Cleanrooms and associated controlled environments – Part 2: Specifications for testing and monitoring to prove continued compliance with ISO 14644-1
- EN ISO 14698-1:2014, Clean Rooms and Associated Controlled Environments - Part 1: Biocontamination Control – General Principles and Methods
- EN 1041:2008, A1:2013 Information supplied by the manufacturer with medical devices.
- EN ISO 15223-1:2016 Medical devices – Symbols to be used with medical device labels, labeling and information to be supplied. Part 1: General requirements.
- EN 15986: 2011 Symbol for use in the labeling of medical devices – requirements for labeling of medical devices containing phthalates.
- ISTA 2A:2011 Partial Simulation Performance Test; Packaged Products 150 lbs. or less.
- ISTA 3E:2017 Partial Simulation Performance Test; Unified loads of Same Product
- EN ISO 11607-1:2017 Packaging for terminally sterilized medical devices.
- European Pharmacopoeia, 98th Edition, 2017
- USP 37/NF 32 (2014) United States Pharmacopoeia
- ANSI/AAMI ST72: 2011 Bacterial endotoxin - Test methods, routine monitoring and alternatives to batch testing.
- EN ISO 10993-1:2018 Biological Evaluation of Medical Devices -- Part 1: Evaluation and testing.
- EN ISO 10993-7:2008, COR 2009 Biological evaluation of medical devices -- Part 7: Ethylene oxide sterilization residuals.
- EN ISO 15747:2018 Plastic containers for intravenous injections.

17 REQUIREMENTS FOR MAINTENANCE AND REPAIR OF THE MEDICAL DEVICE

The Mirasol® WB Set is a single-use medical device and it does not need technical maintenance or repair.

18 WARRANTY LIABILITY

Warranty shelf life: 3 years

When purchasing a Terumo BCT product, the general terms and conditions are applicable. This includes the following warranty statement:

LIMITED WARRANTY Warranty information applicable to the equipment or disposables specified on this invoice is set forth in the price list and any contract for such equipment or disposables and shall be provided by Terumo BCT upon request. For a period of ninety (90) days from the date of the invoice, Terumo BCT will repair or replace Terumo BCT spare parts or optional feature components that Terumo BCT finds defective in material or workmanship. Terumo BCT is not responsible for, and shall have no obligation with respect to any failure caused by normal wear and tear, misuse, unauthorized alterations, accident, neglect, use of non-standard accessory attachments, and/or improper customer maintenance. Terumo BCT does not verify the safety or efficacy of, and makes no warranties, express or implied, with respect to any non-Terumo BCT components or non-Terumo BCT spare parts included at Purchaser's request or any components or products used other than as expressly intended by their manufacturer. Failure to perform preventive maintenance procedures will invalidate this warranty. Terumo BCT is not responsible for any events related to non-Terumo BCT parts installed in the equipment. THIS LIMITED WARRANTY CONTAINS THE PURCHASER'S EXCLUSIVE WARRANTIES AND REMEDIES.

19 CLAIM

When purchasing a Terumo BCT product, the general terms and conditions are applicable. This includes the following liability statement:

Limited liability: To the extent permitted under applicable law, Purchaser's exclusive remedy with respect to any claim, whether in contract, in tort, under statute, under warranty or otherwise shall be limited solely to the refund of the purchase price or replacement of all Products shown to be other than as warranted. Terumo BCT's liability for losses which do not affect the Product (e.g. indirect, incidental, consequential losses, damages caused by delay) shall be excluded. To the extent the foregoing limitation is not permissible under applicable law, the applicable legal limitation of liability under statute or otherwise if there shall be one, shall apply.

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