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CERTIFICATE

中国国际贸易促进委员会暨中国国际商会
China Council for the Promotion of International Trade is China Chamber of International Commerce

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China Council for the Promotion of International Trade
China Chamber of International Commerce

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THIS IS TO CERTIFY THAT: the seal of HUZHOUE MEIQI MEDICAL INSTRUMENTS CO., LTD on the annexed DOCUMENT is genuine.

China Council for the Promotion
of International Trade

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书查询网址 Website for verifying the certificate: <http://www.rzccpit.com/validate.html>

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"I'm certify"
Zhang Feng
General Manager
Huzhou Meiqi Medical Instruments Co., Ltd.
Zhang Feng

" 08 " December 2023
2023.12.08

Stamp / Signature _

湖州美奇医疗器械有限公司
HUZHOU MEIQI MEDICAL INSTRUMENTS CO., LTD.

MANUAL

Continuous Glucose Monitoring System Transmitter

Dr. Briander RGMS-II

2023

1. Description of the medical device

1.1. Name of medical product

Continuous Glucose Monitoring System Transmitter Dr Brinner RGMS-II. (hereinafter referred to as transmitter, product).

1.2. Composition of the medical device/equipment

The "Transmitter for continuous glucose monitoring system Dr. Brinner RGMS-II" includes:

Option 1:

- 1) Transmitter - 1 pc.
- 2) Applicator - 1 pc.
- 3) Separator - 1 pc.
- 4) Instruction manual - 1 pc.

Option 2:

- 1) Transmitter - 1 pc.
- 2) Instruction manual - 1 pc.

1.3. Purpose of the medical device

The transmitter, as part of the Continuous Glucose Monitoring System (CGMS), is designed to collect, process according to a given algorithm and wirelessly transmit signals from an implanted sensor about the glucose content in the interstitial fluid in adults diagnosed with diabetes mellitus and children diagnosed with diabetes mellitus over 2 years of age. The established glycemic range is 1–40 mmol /L (18–720 mg/ dL). The received signal is wirelessly transmitted to the application using a transmitter.

1.4. Application area

Pediatrics. Pediatric endocrinology. Endocrinology. Therapy. Use as prescribed by a pediatrician, endocrinologist, or general practitioner in medical institutions and at home in accordance with the Instructions for Use.

1.5. Indications for use

Diabetes mellitus, prediabetes.

1.6. Contraindications for use

Use in patients is possible only on the recommendation of the attending physician. When wearing the sensor, you must strictly follow the instructions.

Do not use in the presence of the following factors:

- Allergic reaction to product materials.

- Trophic disorders and skin diseases in the areas where the sensor is installed and fixed.
- Self-administration without supervision by patients who suffer from alcoholism, drug addiction, severe mental disorders (e.g., depression, schizophrenia), are unconscious, are unable to understand or master the use of the device, or have severe hearing or vision impairments that may lead to inappropriate use of the cerebrospinal fluid. .

1.7. Possible side effects

Possible side effects are associated with the installation of a sensor, the signal from which is processed by the transmitter. Most likely possible side effects:

1) The appearance of bruising or capillary bleeding. Solution:

- Press the sensor and transmitter firmly for a few minutes until the bleeding stops. You can continue to use the system and monitor your blood sugar levels. If there is a large difference between the system readings and the finger prick meter blood glucose results , you should consider replacing the sensor and installing it in a different location.
- If after applying pressure on the sensor for 5 minutes the bleeding does not stop, it is necessary to remove the sensor and apply a pressure bandage to the wound. If bleeding does not stop within 10 minutes, seek medical attention.

If the patient does not have a known bleeding disorder, place the new sensor on a different area of the skin. If the patient has known bleeding disorders, the decision about the possibility of using LUMG is made after consultation with a doctor, taking into account the individual characteristics of the patient.

2) Appearance of skin irritation or redness

Some users may be hypersensitive to components in the adhesive that adheres the sensor and transmitter to the skin.

If you notice significant skin irritation around or underneath the sensor, remove the sensor and stop using the system. The decision on the possibility of using CGMS is made after consultation with a doctor, taking into account the individual characteristics of the patient.

1.8. System Features

1) Warnings

Interstitial fluid glucose measurements obtained by CGMS are intended to complement, but not replace, capillary blood glucose measurements using a glucometer . They give an idea of the trends in glucose levels and indicate when it is necessary to check it in capillary blood with a glucometer . It is not allowed to adjust the treatment process solely on the basis of information received by the CGMS.

When checking your glucose levels and entering a calibration or reference value, it is important to note:

- Calibration of the CGMS is not mandatory, but is recommended to achieve the best measurement accuracy.
 - The same properly registered meter and test strips with the same lot number should be used throughout the sensor's life cycle.
 - Monitoring and calibration (fixation of reference values of capillary blood glucose) is recommended to be carried out in accordance with the calibration instructions (see clause 5.9.). Calibration must be done strictly in the morning on an empty stomach, when blood glucose levels are most stable.
 - After measuring the capillary blood glucose level, this value should be entered into the MQ CGM application of the user's smartphone immediately, but no later than 2 minutes after drawing blood from a finger prick. If more than 3 minutes have passed since the measurement, it is necessary to analyze the blood again with a glucometer .
 - If you accidentally made a mistake when entering the glucose value in the application, you need to re-test your blood with a glucometer and enter the calibration again, in accordance with the instructions.
 - The monitoring range corresponds to the range of capillary blood glucose levels: 1–40 mmol /L (18–720 mg/ dL).
- 2) The patient, while using the continuous glucose monitoring system, can record information about diet, sports exercises, medications, and insulin in the smartphone application.
 - 3) The sensor and transmitter can be wetted, washed, and you can swim with it, but you cannot keep it under water for more than 30 minutes. After washing, be sure to blot the patch with a towel (to remove excess moisture). Also, if possible, you should dry the LUMI and patch with a hairdryer on low or no heat to remove moisture between the sensor and the skin to prevent moisture from contacting the sensor electrodes.

4) Measurement failure in contact with water: important requirements

If you notice a high, abnormal glucose level, it may be caused by water or sweat getting on the electrodes, which may cause a temporary disruption in the transmission of the electrical signal expressing the glucose level to the sensor. In this case, it is necessary to dry it with a hairdryer at low temperature or without heating the air so that the skin under the sensor dries. After drying, the functioning of the sensor is completely restored. In a situation where water or sweat on the sensor electrodes causes a temporary measurement failure, the monitoring data cannot be considered valid until the system has dried and returned to normal operation.

5) Potential consumers


Adult patients diagnosed with diabetes mellitus and children diagnosed with diabetes mellitus, over 2 years of age (aged 2 to 17 years inclusive, under the supervision of a caregiver at least 18 years old), prediabetes .

A caregiver, legal representative or attending physician can independently monitor the patient's/user's glycemic level through remote access and, if necessary (hypoglycemia/hyperglycemia), provide assistance, track and interpret the results of


measurements obtained by the system, or provide assistance and support to the child in obtaining and interpreting the results measurements.

6) Operating principle

The patient determines the location of the CGMS: on the front side of the shoulder (the area between the shoulder and elbow joints, in a place where there are no muscles and there is a sufficient layer of fat), or on the stomach, at a distance of 2–7 cm from the navel to the left or right, using The applicator places the system (consisting of a sensor and transmitter) on the skin by inserting the sensor electrodes under the skin.

The installed system determines the glucose concentration in the interstitial fluid using electrochemical methods. The glucose oxidase involved in the process is applied to the sensor electrode, also called a biosensor. When glucose from the interstitial fluid enters the enzyme layer through the polymer membrane, it is catalyzed by the enzyme to produce an equivalent amount of hydrogen peroxide. Next, hydrogen peroxide interacts with the electrode. During electrochemical oxidation, the intensity of the generated electrical signal is proportional to the concentration of glucose in the physiological environment in which the sensor is located.

The function of the continuous glucose monitoring system transmitter is to apply constant voltage to the sensor electrode and generate microcurrent, and send a signal to the continuous blood glucose monitoring application. The function of the continuous blood glucose monitoring application is to receive and process the signal sent by the transmitter, convert it into data about the current interstitial fluid glucose level and display it on the smartphone interface.

7) Advantages

- Detection of hidden increases and decreases in blood sugar, especially postprandial (reactive) hyperglycemia, nocturnal asymptomatic hypoglycemia, the dawn phenomenon, the Somogyi effect, which are difficult to detect using traditional monitoring methods.
- Detect changes in blood sugar levels related to factors such as types of foods, types of exercise, treatment plans, mental factors, etc.
- Additional data for creating individual treatment plans.
- Improving patient compliance with treatment regimen.
- Providing visual aids for learning about living with diabetes.

8) Classification of medical device and safety requirements

Table 1 – Information on the classification of medical device and safety requirements

Parameter	Meaning
Potential risk class of use	Class IIb
Protection against electric shock according to EN 60601-1	Class II
Classification of working parts according to EN 60601-1	Has a BF type working part
Electromagnetic compatibility classification according to CISPR 11	Group 1 Class B
Protection against hazardous ingress of water or solid particles into the transmitter attached to the sensor.	IPX7
Operating mode	Designed for continuous operation. Each transmitter's battery life is limited to 60 days or 4 sensor cycles.

2. Description of the main functional elements and components of the medical device

1.1. Transmitter

The transmitter is a device with a battery and a near field communication (Bluetooth) radio transmitter that connects to the sensor of the continuous glucose monitoring system Dr.Brinner RGMS-II, duly registered. Transmitter functions:

- calculation and transmission of a signal about the current level of glucose in the blood to the receiving device;

- applying a constant voltage to create a microvoltage on the sensor electrode, which initiates a reaction with glucose in the interstitial fluid and generates a readable signal proportional to the glucose content. The power supply in the transmitter is a non-removable 3 V battery.

Battery replacement: not provided.

The service life of the transmitter is 60 days, or 4 cycles of sensor operation.

According to the manufacturer's catalog, the transmitter has a REF number (article) MQ-TS II or MQ-TS II (1), depending on the configuration.

Figure 1. Transmitter appearance, top view


Figure 2. Technical drawing of the transmitter structure

Table 2 – Transmitter structure

No.	Name	Function
1	Printed circuit board	Data processing
2	Connector	Connects sensor and transmitter
3	Frame	Protects your device

1.2. Applicator

An applicator is a device for implanting sensor electrodes connected to a transmitter under the skin. The applicator must be secured on top of the transmitter attached to the sensor in its sterile packaging. After this, the applicator is placed on the surface of the skin at the selected implantation site, the user presses a button - a spring is released and the applicator implants a sensor with a transmitter.

Figure 3. Applicator appearance


Figure 4. Applicator structure

Table 3 - Continuous Glucose Monitoring System Applicator

No.	Name	Function
1	Limiter	Control of insertion depth
2	Fixing block	Spring fixation
3	Spring	Sensor implantation
4	Pusher	Transmitter pusher
5	Side button	Spring release
6	Frame	Protects your device
7	Stub	Fixing the mechanism in the body

1.3. Delimiter

The separator is designed to detach the transmitter from the used sensor. The separator contains a lever that separates the transmitter and sensor housings.


Figure 5. External view of the separator, side view

1.4. MQ CGM Continuous Glucose Monitoring System Application,

The MQ CGM glucose monitoring system application is downloaded from the website brinner-tech.ru and installed on a smartphone equipped with a Bluetooth near-field communication channel. The application reads the data transmitted by the transmitter about the level of glucose in the interstitial fluid, visualizes and accumulates them. The application allows you to monitor the effect of medications, diet, exercise, and insulin administration on the glucose content of the patient's interstitial fluid. Providing this information to the physician helps to timely adjust the clinical treatment plan.

Software function: download data, store continuous glucose monitoring results, enter calibration reference values, output result report, etc.

According to the manufacturer's catalog, the application has a REF number (article) MQ-AP II.

1.5. Specifications

Table 4 - Technical specifications

Parameter	Meaning
MQ-TS II Continuous Glucose Monitoring System Transmitter	
Width * Length * Height, mm:	25.6*29.6*8.2 +/- 15%
Weight, g :	5 +/- 15%
Battery capacity, mAh :	70
Battery life:	up to 60 days or 4 sensor operation cycles
Communication method:	Bluetooth
Data transmission and reception distance, m:	1
Operating frequency of Bluetooth transmitter, MHz:	2402-2480
Maximum Bluetooth transmission power , dBm :	7
Continuous Glucose Monitoring System Applicator	
Width * Length * Height, mm:	54.5*60.5*93.5 +/- 15%
Weight, g :	99.5 +/- 15%
Injection force:	>20N, 50N~80N
Depth of insertion of sensor electrodes, mm:	4.5
Delimiter	
Width * Length * Height, mm:	40*30*36 +/- 15%
MQ CGM Continuous Glucose Monitoring System Application	
Monitoring	Glucose monitoring
Calibration	Entering a blood glucose value to calibrate the system

System parameters of sensor and transmitter	
Waterproof (when connecting the transmitter to the sensor):	IPX7
Tracking results update cycle, min:	3
Number of measurements during the update cycle:	18
Electrode voltage:	dynamic
Upper limit of glucose value, mmol /l:	40.0
Lower limit of glucose value, mmol /l:	1.0
Linear correlation coefficient of the current measured by one blood glucose sensor relative to the current of the glucose solution	$r > 0.98$
Linear correlation coefficient of data display	$r > 0.99$

Clinical trials

Continuous Glucose Monitoring System Features Dr.Brinner RGMS-II

Table 5 - Measurement accuracy

BG < 80 mg/ dl (4.4 mmol /l)		BG > 80 mg/ dl (4.4 mmol /l)		General	
Deviation 20 mg/ dl	91.4%	Deviation 20%	88.7%	Deviation 20/20%	87.6%

Finding glucose measurements in clinically significant zones: **percentage (%) of values in zones A and B on the Parkes matched error grid : 99.9%.**

MARD - 9.26% (lower limit of 95% CI - 8.25%, upper limit of 95% CI - 10.27%).

Patients (n=60) during a clinical study in the Russian Federation were surveyed twice in the first and second half of the system study, all completed questionnaires were analyzed.

The survey results are presented below:

Table 6 – Results of surveys of patient groups during clinical trials

Patients	Questionnaire item	Share (%)	95% confidence interval for proportion (%)
Adults	Did not experience any physical pain when installing the sensor	66.67	47.19; 82.71
	Very little physical pain during installation	30.00	14.73; 49.4
Children	Did not experience any physical pain when installing the sensor	36.84	24.45; 50.66

	Very little physical pain during installation	29.82	18.43; 43.40
Adults	Doesn't interfere at all, there was no physical pain after installation	80.00	61.43; 92.29
Children	Doesn't interfere at all, there was no physical pain after installation	71.93	58.46; 83.03
Adults	Does not limit taking a hygienic shower	73.33	54.11; 87.72
	Slightly limits taking a hygienic shower	26.67	12.28; 45.89
Children	Does not limit taking a hygienic shower	77.19	64.16; 87.26
	Slightly limits taking a hygienic shower	22.81	12.74; 35.84
Adults	Does not limit moderate physical activity	91.95	84.12; 96.70
Children	Does not limit moderate physical activity	89.47	78.48; 96.04

When analyzing patient questionnaires to the question of assessing the level of physical pain when installing the sensor of a registered continuous glucose monitoring system, the following results were obtained. "They did not experience pain" or "The pain was characterized as very mild" – in 77.02%, which can be regarded as good tolerance of the installation of the system by patients.

When analyzing patient questionnaires to the question "Assess the ease of use of the studied continuous glucose monitoring system," the following results were obtained. 75.86% of patients considered the system convenient and believed that it did not limit the ability to take a shower, while there were not a single patient who thought that the system significantly limited taking a hygienic shower; and only 24.14% of patients reported that the system slightly limited this ability, which can be regarded as a good indicator of the usability of the studied system.

When analyzing patient questionnaires to the question "Assess the usability of the application of the continuous glucose monitoring system being studied," the following results were obtained. The smartphone application was considered convenient by 81.61% of patients who answered that the statement was true or mostly true, which can be regarded as a good indicator of the ease of use of the application of the system under study.

When analyzing patient questionnaires to the question "Assess possible restrictions when taking a hygienic shower using the studied continuous glucose monitoring system," the following results were obtained. The smartphone application was considered convenient by 81.61% of patients who answered that the statement was true or mostly true, which can be regarded as a good indicator of the ease of use of the application of the system under study.

When analyzing patient questionnaires, 91.95% of patients answered that the system under study does not limit moderate physical activity.

1.6. List and description of medical device materials, including those that come into direct or indirect contact with the human body

Table 5 - Continuous Glucose Monitoring System Transmitter Components

No .	The name of detail	Artie cool for the spare part	Specification or characteristic	Manufacturer	Type of contact with a person
				–	

1	Frame	960702	Acrylonitrile butadiene styrene (ABS), BASF GP-22	Guangdong Zhuyuan Molding Co.,Ltd	Brief contact with intact skin
2	Printed circuit board consisting of:	960697	BL530		No contact
	CPU			Texas Instruments	
	Chips			Texas Instruments	
	Operational amplifier			Texas Instruments	
	Transmitter module			Oviphone Technology Limited	
3	Connector	960696	Brass		No contact
4	Battery			Panasonic	No contact

The applicator body and separator are made of acrylonitrile butadiene styrene (ABS), manufactured by LG CHEM LTD, Republic of Korea, and come into short-term contact with intact skin.

1.7. Information on medicinal products and materials of animal/human origin

The product does not contain drugs or pharmaceutical substances.

No materials of animal/human origin are used in the manufacture of the product.

3. Products intended for use in combination with MI

The transmitter, when used as intended, is part of a continuous glucose monitoring system consisting of a transmitter and a disposable sensor.

2.1. Continuous Glucose Monitoring System Sensor Dr.Brinner RGMS-II

Continuous Glucose Monitoring System Sensor Dr.Brinner RGMS-II consists of electrodes implanted under the skin, a sensor base, and a patch.

Figure 6. General view of the sensor secured in an opened individual package

The subcutaneous part of the sensor consists of three electrodes, which are shaped like microneedles. The sensor electrodes are coated with a polymer material.


Figure 7. Appearance of the back of the sensor, electrodes and medical patch applied to the base of the sensor

Glucose in the interstitial fluid penetrates the polymer membrane into the reagent layer and produces hydrogen peroxide. Next, electrochemical oxidation of hydrogen peroxide occurs on the surface of the electrode, during which the sensor current strength and the parameters of the physiological environment have a certain proportional dependence on the glucose concentration.

The electrodes should be positioned vertically and perpendicular to the surface of the skin and the plane of the base with the medical patch. The electrodes penetrate the skin and fix the sensor. The outer surface of contact with human tissue is a biocompatible polymer film.

According to the manufacturer's catalog, the sensor has a REF number (article) MQ-GS II.

2.2. Smartphone

The application is installed on the smartphone.

Table 6 - Smartphone requirements:

Software name:	MQ CGM (Android) / MQ CGM (iOS)
Version:	1.0.2.2 / 1.0.3.2
Operating Environment:	Android 5.0 and later; iOS 10.0 and later
Permission:	1334*750 and higher
Processor core frequency:	1.0 GHz
RAM:	4 GB
Working volume on the drive:	At least 16 GB, the program takes up 500 MB
Interface:	Touch screen
Communication method:	Bluetooth
Data transfer method:	HTTP

4. Precautionary measures

- Patients with possible allergic reactions should monitor their condition closely and remove the sensor from the transmitter if a severe allergic reaction occurs.

- Patients with an installed sensor and transmitter should not be in strong magnetic fields, and X-ray, MRI and CT scans should not be performed.
- Use with caution in patients with trophic disorders and skin diseases.

Patients installing the sensor on their arm should avoid hitting the door frame, moving when wiping with a towel, should not lie on the sensor, lean on it, sleep on it, sit on it, pull it with an elastic band or belt, hit it, cling to clothing and allow other physical impact or pressure: this can lead to abnormal reactions on the skin and incorrect operation of the sensor.

Try to avoid strenuous exercise or frequent squeezing or twisting of the skin near the sensor while you sleep, as this may cause the sensor to come off prematurely or become damaged.

5. Operating procedure

The transmitter is used as part of a continuous glucose monitoring system (CGMS), which is formed from a sensor and products specified in paragraph 3. Before use, make sure that you have all the products required for use.

4.1. terms of Use

Ambient temperature from plus 10 °C to plus 42 °C.

Relative humidity from 5% to 80%.

Atmospheric pressure from 70 kPa to 106 kPa.

4.2. Preliminary steps for installing the system


- 1) Make sure you have alcohol wipes, sensor, transmitter, applicator, blood glucose meter .
- 2) Before you start working with CHMG, make sure that mobile communications, Bluetooth and geolocation are turned on on your phone.
- 3) Download the MQ CGM application from brinner-tech.ru*, following the instructions in the application.

*Note: If the site is being updated, please download the MQ CGM app from the app store: App Store for iPhone and Play mobile smartphones Market for Android respectively.

Allow all requested access on the phone. Including:

- 1) device memory;
- 2) bluetooth connection ;
- 3) geolocation ;
- 4) notifications;
- 5) battery;
- 6) camera.

After installing the application, additionally go to “Settings” - “Applications” on your phone and make sure that all the settings are given to the application, including the setting not to close the application when the battery is low and all possible permissions.

4.3. Please Note: Application Permissions

The MQ CGM app is a medical app that is critical to your health. Make sure that the application is set to the maximum phone settings, otherwise the application may work with errors and interruptions. The mobile phone may close the application spontaneously. It is recommended to open the application on your phone from time to time and ensure that it is active.

When installing an iPhone app, go to your phone's Settings, Applications section, and grant the app all the requested permissions.

When installing an Android app, go to the Permissions section of the downloaded app and grant it all requested permissions, including:

- Battery consumption without restrictions;
- Permissions without restrictions;
- Bluetooth without restrictions;
- Internet without restrictions;
- Notifications without restrictions;
- Camera without restrictions;
- Geolocation without restrictions;
- Unlimited storage and memory;
- Autostart function without restrictions.

4.4. Continuous Glucose Monitoring System Sensor Operating Instructions

- 1) It is recommended to install the sensor in the morning on an empty stomach. Implantation

and warming up of the sensor must be carried out in the morning on an empty stomach, before meals.

- 2) Determine where to install the system


Recommended areas for sensor placement include the front of the shoulder and abdomen.

- If you have a thin build and lack of subcutaneous fat on your arm, it is recommended to place


it on your stomach. In patients with a small amount of subcutaneous fat in the upper arms and abdomen, the sensor electrodes may become embedded in muscle tissue, causing capillary bleeding, pain, and sensor malfunction. Carefully select the appropriate sensor installation area in this situation. For children, it is recommended to place it on the back of the shoulder, on the skin fold of the triceps.

- If you have a normal build, place it on your arm or stomach.
- If you are overweight, it is recommended to place it on your arm.


- 3) Shoulder mounting

The ideal place to place the sensor for adults and older children is the front of the shoulder, along the chest line, where the deltoid muscle ends. In this place, there is the lowest probability that external influence on the sensor will occur during sleep, movement and in other situations.


It is recommended to install the system on a non-dominant hand, that is, for right-handers - not the left, for left-handers - on the right.

- The most preferred placement location, #1 in the picture below, is where the deltoid ends and the biceps begins, on the front of the shoulder, in line with the chest.
- The second most preferred placement site, #2 in the picture below, is where the deltoid ends and the brachialis begins, on the lateral side of the shoulder,
- The third preferred placement location, #3 in the picture below, is the triceps skin fold location on the back of the upper arm. For young children with insufficient subcutaneous fat, installation in this area is recommended.


If you are installing the system yourself on the front side of the shoulder, stand up straight, clasp your chosen arm around your stomach, press it tightly to your body so that the shoulder muscles are tightly pressed to the body and at the same time relaxed.

In this position, install the CHMG with your other hand.


If you are receiving assistance during installation, place your chosen hand on your waist with your palms facing out.

4) Installation on the stomach


The recommended area for self-installation is on the abdomen - 2–7 cm from the navel to the left or right and more than 7 cm from the site of implantation of the insulin pump cannula. Avoid the belt area as pressure on the sensor may cause malfunction, deformation of the electrodes, or peeling off.

For best results, avoid excessive movements, which can weaken the sensor and medical tape. Avoid accidentally detaching the sensor.

5) Prepare the installation site

Use an alcohol wipe to clean the selected area, then allow the alcohol to dry completely before installing. Use only alcohol for cleaning; do not use iodine or cosmetics: they will affect the stability of the patch. Before installing, wait until the alcohol has completely dried, otherwise the patch will not adhere as well.


6) Add a sensor in the application

Scan the QR code on the sensor packaging or enter the sensor number manually. After successfully entering the code, click the


“Next” button.

7) Tear off the film of the sterile sensor packaging


8) Snap the transmitter to the sensor


9) Make sure the connection is correct

Pay attention to the matching connections of the transmitter and the gold-colored sensor, as well as the installation order.

Place the transmitter on top of the sensor, first with the side opposite the gold-colored metal connection. It is recommended that the package containing the transmitter and sensor be placed on a hard surface.


After this, press the transmitter firmly . The connection will be established correctly after clicking.

Please note that the transmitter must be pressed tightly and evenly to the sensor; if the transmitter is pressed loosely or diagonally, the connection was incorrect, disconnect the transmitter and connect again.

10) Connect the transmitter to the receiving device - a smartphone via the application

Remove the transmitter and check the code on the transmitter housing. It is located on the inside of the transmitter itself.


For ease of entry, the last three digits of the code are duplicated on the side of the transmitter body. Enter these last three digits of the code into the MQ CGM application following the application prompts.


11) Preparing the applicator

The following shows the operating state of the applicator in which it can be used to install the SIMG.

The inside of the applicator should not be pressed into the applicator. In the illustration below, the applicator is cocked and not ready for use. You need to press a key on the applicator to bring it into working condition, or replace the applicator.

Make sure that the oval of the applicator is oriented in the same way as the oval of the transmitter body, and the inscriptions on the transmitter and applicator are oriented the same way.

Install the applicator onto the transmitter .



Press the applicator onto the transmitter until it clicks.

Lift the applicator together with the sensor and transmitter strictly vertically.

The system will be inside the applicator. The applicator is ready for installation.


Hold the applicator with your thumb on the applicator button.


12) Position when installing the system with an applicator

When installing the system on the shoulder, you need to position your arm in a certain way so that it is relaxed and the installation site is accessible.



When installing yourself, stand up straight, clasp your chosen arm around your stomach, press it tightly to your body so that the shoulder muscles are tightly pressed to the body and at the same time relaxed.

Press the applicator firmly onto the installation site, previously treated with alcohol, holding it strictly perpendicular to the surface of the skin and pressing it against the body. The applicator should be pressed as tightly as possible into the body, the installation site should be even and taut.


Press the applicator firmly onto the installation site, previously treated with alcohol, holding it strictly perpendicular to the surface of the skin and pressing it against the body. The applicator should be pressed as tightly as possible into the body, the installation site should be even and taut.

When installing with assistance, lower your chosen arm along your body, palms inward. In this case, the hand should be relaxed.

Grab your hand from behind so that the skin on your hand becomes taut and the installation site is level. Press the applicator firmly onto the installation site, previously treated with alcohol, holding it strictly perpendicular to the surface of the skin and pressing it against the body. After pressing the applicator, with your hand, with which you clasped the patient's arm from behind, press down on the arm so that the place where the applicator is installed becomes flat.

Please note that when installing the system for children whose hand diameter is quite small, it is especially important, after pressing the applicator to the hand, with the hand you clasped the hand from behind, press on the hand so that the installation site of the applicator becomes flat, otherwise after installation the sensor patch will not be tight fit to the hand.


13) Installation position: important requirements

- If the applicator is not applied perpendicularly to the skin, but at an angle, the sensor may not install correctly and show less accurate measurement results.
- The applicator button is quite tight. It is recommended to hold your thumb or index finger on the button when installing it so that it is comfortable to press. Do not move the applicator while you press the button - make sure that your hand does not shake or move when installing.
- Movement, displacement and vibration of the applicator during installation is the main reason that can lead to damage or incorrect installation of the sensor!


14) Installing the system with an applicator

Press the applicator button, a pop will be heard, this means that the sensor and transmitter are installed.

Hold the applicator in place for 5-10 seconds to allow the patch to adhere.

After waiting, slowly, strictly perpendicularly, move the applicator away from the skin with the sensor and transmitter installed.


Gently smooth the patch by pressing it firmly against the skin with your fingers and press firmly on the transmitter to ensure it is securely attached.

15) Using a retention patch

If necessary, use additional adhesive to secure the sensor and transmitter in place by applying it to the skin over the transmitter.

16) **Retention patch: Important requirements**

It is recommended to use an additional patch only when the patch built into the sensor has become loose during use. It is not recommended to use an additional patch from the very beginning of use to prevent irritation and redness of the skin under the patch.

4.5. **Directions for Use: Important Requirements**

- Keep your mobile phone turned on no more than 1 meter from the transmitter to ensure uninterrupted data transmission.
- It is prohibited to remove the sensor from the skin and reinstall it.
- It is prohibited to disconnect the sensor and transmitter during operation.
- Bluetooth connection, use another mobile phone.
- The transmitter has a shelf life of 4 years. Use the product before the expiration date.
- The transmitter consumes power while in use. When the transmitter battery is low, the app will show a notification. The transmitter is designed for 4 full sensor operation cycles of 14 days. When installing the next sensor, use a new transmitter.
- Please note that the sensor must be installed in a place where there is a sufficient layer of subcutaneous fat directly under the skin.

4.6. **Connecting the transmitter with a smartphone**

Operating frequency of Bluetooth transmitter: 2402–2480 MHz;

Maximum Bluetooth transmission power : 7 dBm .

4.7. **Testing with a glucometer**

Always **use the same meter** for calibration . The use of another glucometer is not permitted. It is prohibited to use or alternate readings from two different glucometers .

4.8. **Removing the sensor, preparing to use the transmitter with the new sensor**

- 1) After 14 days of use, remove the retaining patch (if it was attached to the system for additional fixation) around the entire perimeter around the SIMG, then remove the medical

sensor patch around the entire perimeter of the sensor body, without moving the sensor itself along the surface of the skin.

- 2) Remove the sensor with the transmitter from the skin strictly perpendicular to the surface of the skin.
- 3) Place the sensor with the transmitter electrodes on a hard surface and press firmly, crushing the sensor electrodes.
- 4) To disinfect after removing the sensor, wipe the surface of the skin where the sensor was installed with an alcohol wipe.
- 5) Apply adhesive tape as needed to protect and heal the sensor site.
- 6) Separate the transmitter from the sensor using a separator: place the separator on the top of the transmitter so that the top characters line up. Using the adjustable part of the separator, hook the transmitter by the recess and press it inward to disconnect from the sensor.
- 7) Dispose of the sensor according to local guidelines.
- 8) Use the transmitter for the next sensor until the battery is completely depleted, then dispose of it according to local guidelines.
- 9) Rotate sensor locations regularly to avoid discomfort or skin irritation. The areas of skin successively used to install the sensor must be separated by at least 5 cm.

10) Please note: Implantation site after sensor removal

After removing the sensor, 3 noticeable points will remain on the skin through which penetration and contact of the sensor electrodes with the interstitial fluid occurred. Electrode penetration points are normal and a standard feature of all existing continuous glucose monitoring systems. After removing the sensor, it is recommended to treat the removal site with an alcohol wipe and apply a bactericidal patch. After a few days, the penetration points of the electrodes will tighten and heal.

4.9. System Callibration

Once you have installed CHMG Dr.Brinner , in the already installed MQ CGM application on your smartphone, the system will warm up for an hour (60 minutes), after which the system will start working.

1) Please note: abstain from food while the system is warming up

Warming up the system before starting operation takes 1 hour. During this time, you should not eat, take medications, or exercise to ensure your glucose levels remain as stable as possible while the system warms up.

After the sensor warms up, you will see the first readings, including on the graph of your application.

2) Please Note: Achieving Best Results Using Calibration

For best results and maximum measurement accuracy by the system, the patient can calibrate every day, but not more than once a day, strictly in the morning on an empty stomach. To determine the need for calibration, you need to check your blood with a glucometer in the morning, before meals and treatment, and compare the readings with the glucometer. If the discrepancy between the glucometer and SNMG readings is more than 1 mmol /l in the range up to 10 mmol /l, or more than 2 mmol /l in the range above 10 mmol /l, it is recommended to perform calibration. In this case, calibration should be carried out strictly in the morning on an empty stomach.

3) Please note: mandatory calibration requirements

- Calibration cannot be performed more than once a day.
- Calibration cannot be performed on an empty stomach. If you forget to calibrate in the morning on an empty stomach, there is no need to do it later. Reschedule calibration to the next day.
- meter to measure your glucose and record only its readings.
- Do not change and use the same meter for the entire 14 days of sensor use.
- Enter the exact value you obtained when measuring your blood. Do not round off your meter readings. Record exactly the meter reading that you see on the meter screen.
- After measuring the blood glucose level with a finger stick glucometer, no later than 3 minutes after drawing blood, enter the blood glucose level into the MQ CGM application of your smartphone.

- If more than 3 minutes have passed since the blood was drawn, you must analyze the blood again with a glucometer and enter the updated value during calibration.
- If you make a mistake when entering your glucose value, you will need to retest your blood with the meter and enter the correct calibration value again. Please note that calibration must only be performed on an empty stomach.

4) Please note: special calibration cases, slight discrepancies with the meter

If you notice the glucose level value of the Dr.Brinner system is below 3.9 mmol /L or above 15 mmol /L, it is recommended to check the readings on a glucometer .

- If the discrepancy between the glucometer and SNMG readings is more than 1 mmol /L in the range up to 10 mmol /L, or more than 2 mmol /L in the range above 10 mmol /L, it is recommended to perform double calibration. Calibrate as soon as convenient, but only on an empty stomach, at least 3 hours after eating.
- Calibrate as usual the next morning on an empty stomach.

5) Please note: special calibration cases, monitoring values too high or too low

If the system consistently shows the same results, too high or too low, this can happen in three cases: either due to the system becoming detached, or due to water getting on the sensor electrodes, or due to damage to the sensor electrodes.

- Check to see if the patch has loosened or the system has come loose. Add or replace an additional patch, pressing the system tightly to the body.
- Check to see if any water has gotten under the patch. Remove the extra patch. Dry the sensor using a cold hair dryer. Install an additional patch.

After this, a double calibration is required:

- Calibrate as soon as convenient, but only on an empty stomach, at least 3 hours after eating.
- Calibrate as usual the next morning on an empty stomach.

6) Please note: special calibration cases, moisture on the electrodes

If you notice a high, abnormal glucose level, it may be caused by water or sweat getting on the electrodes, which may cause a temporary disruption in the transmission of the electrical signal expressing the glucose level to the sensor. In this case, it is necessary to dry it with a hairdryer at low temperature or without heating the air so that the skin under the sensor dries.

After this, a double calibration is required:

- Calibrate as soon as convenient, but only on an empty stomach, at least 3 hours after eating.
- Calibrate as usual the next morning on an empty stomach.

After drying and calibration, the functioning of the sensor is completely restored.

6. Maintenance

The transmitter is a precision device. If it fails, it can be returned to the manufacturer's representative at an authorized service center to account for defects (if any) or disposed of if it can be replaced with a new transmitter.

7. Cleaning the transmitter connected to the sensor

- Clean the outside of the transmitter with a mild detergent and a soft, damp cloth. To wipe dry, use another cloth. When cleaning the transmitter, hold it steady with your other hand.
- Wipe the transmitter with a disinfectant alcohol wipe.
- Do not use solvents, nail polish remover or thinner to clean the exterior surface.
- Keep the transmitter in a dry place when not in use.
- Do not use any lubricants.

- Do not immerse the transmitter in water or other liquids.

8. Errors, adverse events and solutions

The continuous glucose monitoring system is a precision electronic device. In rare cases, the application may show various errors.

Four simple steps for the user for any situation if the system does not work correctly, does not show notifications, does not connect to the transmitter:

- Exit the application completely, turn the application off and on again, wait a few minutes;
- Turn Bluetooth off and on again on your phone, wait a few minutes;
- Restart your phone, turn on the application, wait a few minutes;
- Go to the phone's system settings in the "Settings" - "Applications" section, select MQ CGM from the list of applications. Make sure the app has maximum permissions on your phone, otherwise the app may experience errors and interruptions, including:
 - Battery consumption and optimization without restrictions;
 - Permissions without restrictions;
 - Bluetooth without restrictions;
 - Internet use without restrictions;
 - Sound, light and vibration notifications without restrictions;
 - Camera without restrictions;
 - Geolocation without restrictions;
 - Storage, files and memory without restrictions;
 - Devices nearby without restrictions;
 - Unlimited use time;
 - Autostart function without restrictions;
 - Autostart function without restrictions;
 - The function of placing the application in the background is disabled;
 - The function to close an inactive application is disabled.

Please note: Application settings on Android mobile phones from different manufacturers may vary. The MQ CGM application is a medical application that is critical for your health, the same approach to settings with any mobile phones is the need to provide maximum permissions and remove all restrictions when using the application.

Please note: the operating system of Apple mobile phones tends to put applications into the background to save battery power and mobile phone resources. When using the system with an

Apple mobile phone , you must always keep the application open, do not close or minimize the application, and, if possible, close other applications that are not in use.

Below are the main errors and how to solve them.

1) Sensor error

The transmitter is out of range of the signal, or the application is turned off, or the application cannot connect to the transmitter in sleep mode due to phone settings, or the Internet, Bluetooth , geolocation , access to the battery, and phone memory are disabled. Solution:

- Make sure your phone is no more than 1 meter away from the transmitter.
- Reboot your phone. Close the application and then open it again.
- Make sure Internet, Bluetooth , geolocation are turned on .
- Make sure that the application has permission to access the Internet, Bluetooth , geolocation , memory, battery.
- Make sure the sensor is securely attached to your body.
- The sensor may be broken; If this reminder does not disappear within 6 hours, please replace the sensor.

2) Low charge

The battery inside the transmitter is not enough to charge, replace the transmitter with a new one, and the sensor also needs to be replaced.

3) Voltage error

Check that the sensor is securely installed and the installation location is correct; if necessary, replace the sensor with a new one. Check to see if the sensor is wet or if water or sweat has come into contact with the sensor electrodes, and dry the sensor with a cold hairdryer.

4) Test value too high

Your glucose level is above the limit. Once your blood glucose level drops to the set value, the reminder will disappear.

In addition, the error may be caused by incorrect installation of the sensor, water on the electrodes, or damage to the sensor due to impact, pressure, or displacement. If water or sweat gets on the electrodes, the problem will be eliminated after drying the sensor with a cold hairdryer.

5) Test value too low

Your glucose level is above the limit. Once your blood glucose level drops to the set value, the reminder will disappear.

In addition, the error may be caused by incorrect installation of the sensor, water on the electrodes, or damage to the sensor due to impact, pressure, or displacement. If water or sweat gets on the electrodes, the problem will be eliminated after drying the sensor with a cold hairdryer.

6) Scan error, try again

The phone was unable to scan the QR code from the sensor packaging using the app. The Internet, Bluetooth, access to the battery or memory, geolocation or camera may be disabled.

- Reboot your phone.
 - Close and then reopen the application.
 - Make sure the Internet, Bluetooth, geolocation, camera are turned on.
 - Make sure that the application has permission to access the Internet, Bluetooth, geolocation, camera, battery, memory.
- 7) Error communicating with the server, synchronize again / Network error, make sure your phone is connected to the Internet.**

The application cannot connect to the network. Make sure that the smartphone is connected to a mobile network or Wi-Fi network, the application has permission to access the Internet, battery.

8) Transmitter error

- Reboot your phone.
- Close and then reopen the application.
- Make sure Internet, Bluetooth, geolocation are turned on. Make sure that the application has permission to access the Internet, Bluetooth, geolocation, battery, memory.

9) Monitoring termination error

The transmitter is out of signal range. Make sure your phone is no more than 1 meter away from the transmitter. Reboot your phone.

10) The connection is not established

- The transmitter is out of signal range, or the application is turned off, or the application cannot connect to the transmitter in sleep mode due to phone settings, or the Internet, Bluetooth, geolocation are disabled.
- Make sure your phone is no more than 1 meter away from the transmitter.
- Reboot your phone. Turn the application off and on again.
- Make sure Internet, Bluetooth, geolocation are turned on.
- Make sure the application has permission to access the Internet, Bluetooth, geolocation.
- Make sure the sensor is securely attached to your body.

11) Itchy skin due to sensor adhesive tape or electrodes

If the skin itching is caused by the electrode patch or microneedles, don't worry, remove and discard the sensor, the skin will repair itself. Try applying the new sensor to a different area of the skin; if the itching continues, seek medical advice.

12) Constant pain at the implantation site, bruising or capillary bleeding.

Some itching may occur, especially during the first few days of wearing the sensor. If the sensor site is very painful, or has a bruise or bruise, the electrodes may have become embedded in a muscle, or nerve, or vessel, or scar, or other skin lesion during implantation. You must remove the used sensor, discard it, and install a new sensor in another location. Please note that the sensor must be installed in a place where there is no muscle directly under the skin and there is a sufficient layer of subcutaneous fat.

9. Labeling, delivery set, packaging

8.1. Transmitter markings

On the outside of the transmitter is indicated:

- Trademark symbol (at the discretion of the manufacturer)

On the inside of the transmitter it is indicated:

- product name
- Date of manufacture
- Model number REF MQ-TS II or REF MQ - TS II (1)
- Serial number
- Symbol "Attention! Refer to user manual"
- "Read the operating instructions" symbol

Table 7 - Symbols used for marking

Symbol	Transcript	Symbol	Transcript
	Do not use if packaging is damaged		Humidity range
	Do not reuse (Do not reuse)		Atmospheric pressure limitation
	Working part type BF		Manufacturer
	Storage temperature		Batch code

	Fragile, handle with care		Date of manufacture
	Do not throw away with household waste		Use before
	Carefully!		Serial number
IPX7	Waterproof level		Non-pyrogenic
	Read the instruction manual		Avoid exposure to sunlight
	Non-ionizing radiation		Keep away from moisture
	Class 2 equipment		

10. Transportation and storage requirements

Products are transported by all types of transport in covered vehicles in accordance with the rules for the transportation of goods in force for this type of transport in accordance with the transportation conditions.

9.1. Transportation conditions

- Transportation conditions: +2°C - +40°C.
- Relative humidity: 5% to 80%.
- Atmospheric pressure: from 70 kPa to 106 kPa.

9.2. Storage conditions

- Storage conditions: +2°C - +40°C. Store in a dark, dry place. Avoid exposure to direct sunlight and rain.
- Shelf life: 4 years.
- Operating conditions: +5°C - +42°C.
- Service life: 60 days or 4 sensor operation cycles.

11. Environmental protection requirements, including requirements for safe destruction and disposal

The transmitter is an electrical medical device; its disposal together with municipal solid waste is prohibited! Use a separate waste disposal system and dispose of the transmitter with electrical products.

The packaging is epidemiologically safe and can be disposed of together with municipal solid waste.

12. Manufacturer's warranty

11.1. Warranty

- Warranty shelf life of the SNMG transmitter: 4 years from the date of production.
- Warranty shelf life of the applicator: 24 months.
- Warranty shelf life of the separator: 24 months.

11.2. Operating data

- Service life of the SNMG transmitter: up to 60 days or 4 cycles of sensor operation.
- Applicator service life: 24 months (until expiration date).
- Separator service life: 24 months (until expiration date).

13. List of international standards applied to ensure the safety, effectiveness and quality of a medical device

- National standard ISO 10993-5-2011 "Evaluation of the biological effect of medical devices. Part 5. Cytotoxicity studies : methods in vitro " by MTT method.
- Interstate standard GOST ISO 10993-10-2011 "Assessment of the biological effect of medical devices. Part 10. Studies of irritant and sensitizing effects."
- National standard GOST R ISO 10993-3-2010 "Assessment of the biological effect of medical devices. Part 3. Study of genotoxicity , carcinogenicity , and toxic effects on reproductive function."
- National standard GOST R ISO 10993-6-2009 "Assessment of the biological effect of medical devices. Part 6. Studies of local effects after implantation."
- IEC 60601-1:2012 "Medical electrical apparatus. Part 1. General requirements for basic safety and essential characteristics."
- ETSI EN 300 328-2016: Broadband transmission systems; Data transmission equipment operating in 2.4 GHz ISM groups.
- IEC 623204-2022 "Medical devices. Software. Life cycle processes".

12.1. List of certificates, declarations of conformity with national and international standards

The production quality management system complies with the requirements of ISO 13485:2016, which is confirmed by the following documents:

- Certificate No. SX2112917-1.

14. Main results of risk analysis and management

A list of risks identified during the risk analysis process and a description of how to manage these risks in order to reduce them to an acceptable level are described in the document "Risk Management File".

15. MQ CGM Software Details

14.1. Mobile application launch process

- 1) Launch the MQ CGM mobile app from your phone.
- 2) Grant permissions for quality use of the application:

For iPhone : Go to your phone 's Settings, Applications section, and grant the app all the requested permissions.

For Android : Go to the Rights or Permissions section of the downloaded app and grant the app maximum access for all requested permissions, including:

- Battery consumption without restrictions;
- Permissions without restrictions;
- Bluetooth without restrictions;
- Internet without restrictions;
- Notifications without restrictions;
- Camera without restrictions;
- Geolocation without restrictions;
- Unlimited storage and memory;
- Autostart function.

The MQ CGM app is a medical app that is critical to your health. Make sure that the application is set to the maximum phone settings, otherwise the application may work with errors and interruptions. The mobile phone may close the application spontaneously. It is recommended to open the application on your phone from time to time and ensure that it is active.

14.2. Login and registration

To register a new user, click "Subscribe now", in the window that opens, enter the following data:

- username in Russian (from 6-20 characters);
- password (in Latin letters, 6 characters or more);
- repeat password entry;
- Secret Question;
- answer to the secret question.

Next, click "Subscribe".

For already authorized users

- Login:

- Enter your username and password.
- Log in and enter the tracking interface.
- On the main page of the application, click “Start Tracking”.
- The system will ask you to enter the sensor serial number.
- Scan the QR code on the sensor packaging or enter it manually.
- After successfully entering the code, click the “next” button.
- In the window that opens, you will see instructions for installing the system; read them carefully, clicking the “Next” button each time.
- The system will search for a transmitter.
- To confirm the transmitter number, click on the displayed code on the phone screen.
- Carefully read the instructions for preparing the system.
- Once you have successfully installed the system, the app will countdown the heating time. Click Finish.
- The device will initialize within an hour.
- After an hour, you will receive a push notification from the application that initialization has been successfully completed:
- On the application screen that appears, calibrate using the current glucose value measured by the meter , click “Confirm”.
- The system will update the monitoring value within 3 minutes.

14.3. Settings

Go to the main screen of the application, switch to “Settings”.

- Set your blood glucose lower and upper limits (target range). Get a reminder from the app when your blood glucose levels are outside your target range.
- Enter information about insulin injections, exercise, diet, and medications.

14.4. FAQ (Possible problems during use and corresponding solutions)

In most cases of incorrect operation of the application, you need to make sure that on your mobile phone the application is given access to the Internet, Bluetooth , geolocation , battery, and phone memory. The notification and sound reminder function is enabled in the application settings.

Go to the phone settings, application section and make sure that the MQ CGM application is given access to the Internet, Bluetooth , geolocation , battery, phone memory with the highest permissions.

Minimize the number of applications running in the background and leave the MQ CGM application open.

- 1) The MQ CGM application was unable to detect the continuous glucose monitoring system after installing the sensor in the transmitter.
 - Step 1: Make sure that on your mobile phone the connected application is given access to the Internet, Bluetooth , geolocation , battery, phone memory.
 - Step 2: Search for the device via Bluetooth or try searching for the device using another mobile device. If the transmitter cannot be detected, please separate the transmitter from the sensor using a separating tool and install again.
 - Step 3: If the transmitter cannot be detected after step 2, please replace it with a new one.
- 2) What to do if the application shows that the device is connecting or disconnecting.

- Step 1: Make sure that the application is connected to the mobile phone and has access to the Internet, Bluetooth , geolocation , battery, and phone memory.
- Step 2: Log out of the background app, then log back in and wait 3 minutes.
- Step 3: If the device still doesn't work, turn off your phone and restart it. After rebooting, connect Bluetooth and geolocation , and then re-enter the application.

- 3) What should I do if the user mistakenly entered the wrong calibration?

If three minutes or less have passed since the meter was calibrated , enter the correct calibration again; If more than three minutes have passed, you must retest your blood with the meter and enter the correct calibration value again. Please note that calibration must only be performed on an empty stomach.

- 4) The glucose data in the app is not updated on time.

- Step 1: Make sure that the application is connected to the mobile phone and has access to the Internet, Bluetooth , geolocation , battery, and phone memory.
- Step 2: Make sure the system and the phone on which the app is installed are less than 2 meters away, and check the phone to make sure the Bluetooth function is turned on.

- 5) Why doesn't the high or low glucose alarm go off when it happens.

- Make sure that the application is connected to the mobile phone and has access to the Internet, Bluetooth , geolocation , battery, and phone memory.
- Enable the notification and sound reminder function in the application settings.
- Go to the phone settings, application section and make sure that the MQ CGM application is given access to the Internet, Bluetooth , geolocation , battery, phone memory with the highest permissions.
- Minimize the number of applications running in the background and leave the MQ CGM application open.

- 6) What should I do if I accidentally click end wearing session in the app?

Click "my" in the bottom right corner of the application, find data history, find the current session and click continue monitoring.

- 7) The app shows low battery during a session.

You can wear the sensor for the rest of the session and use a new transmitter for the next sensor. ,
Additional information about solving common problems and errors can be found on the website
www.brinner-tech.ru

16. Information about the manufacturer of the medical device and place of production

Name: Huzhou Meiqi Medical Instruments Co., Ltd.

Address: Building 5, NO. 58 Yaochi Road, Huzhou, 313000 Zhejiang, P.R. China

17. Authorized representative of the manufacturer

Name of company

Full: Limited Liability Company "Prodrom", Russia

Abbreviated: Prodrom LLC, Russia

Legal address: 364028, Grozny, lane. Reznoy, building 4, office 11

Postal/actual address: 123060, Moscow, st. People's Militia 43, bldg./p. 2, office 8

Phone: +79996138730

Email address: info@brinner.ru, website: www.brinner-tech.ru

Appendix A. Declaration of Electromagnetic Compatibility

Please install and use the device in accordance with EMC instructions.

The continuous glucose monitoring system is a portable radio frequency medical electrical device, please keep it away from other medical electrical equipment.

The Continuous Glucose Monitoring System cannot connect or interface with any device other than a smartphone with the proprietary app installed.

1) Basic work scenario

- Working conditions. Normal transmitter operating range: 260-290 nA.

2) Test method

The glucose emitter is connected to the analog sensor, connected to the glucose emitter. During the countdown, you can check the data update via the smartphone interface, the correct value range is 260nA-290nA.

3) EMC

Guidance and manufacturer's declaration - electromagnetic emissions		
LMWIs are intended for use in the electromagnetic environment specified below. The purchaser or user should ensure that they are used in the specified electromagnetic environment.		
Electromagnetic Emission Test	Correspondence	Electromagnetic environment - instructions
Radio interference according to CISPR 11	I Group	SNMGs use RF energy only to perform internal functions. RF emissions are low and are not likely to cause interference with nearby electronic equipment. LUMGs must emit electromagnetic energy to perform their primary function. Possible impact on nearby electronic equipment.
Radio interference according to CISPR 11	B Category	
Harmonic current components according to IEC 61000-3-2	N.A.	
Voltage fluctuations and flicker according to IEC 61000-3-3	N.A.	

4) EMI

Guidance and Manufacturer's Declaration - Immunity			
system is intended for use in the electromagnetic environment specified below. The purchaser or user should ensure that they are used in the specified electromagnetic environment.			
Immunity test	Test level according to IEC 60601	Compliance level	Electromagnetic environment - instructions
Electrostatic discharge (ESD) according to IEC 61000-4-2	±6 kV - contact discharge		The floor in the room is made of wood, concrete or ceramic tiles. For floors covered with synthetic material, the relative air humidity is at least 30%.
	±8 kV - air discharge		
Nanosecond impulse noise according to IEC 61000-4-4	±2 kV - for power supply lines ±1 kV - for input/output lines		Electrical power quality on the grid according to typical conditions in a commercial or hospital environment

High energy microsecond impulse noise according to IEC 61000-4-5	± 1 kV when applying interference in a wire-to-wire diagram ± 2 kV when applying interference in a wire-to-ground diagram		The electrical power quality of the electrical network should be consistent with that of a typical commercial or hospital environment.
Voltage dips, short-term interruptions and voltage changes in the input power lines according to IEC 61000-4-11	$< 5\% U_{H1}$ (voltage dip $> 95\% U_{H1}$) for 0.5 periods $40\% U_{H1}$ (voltage dip $60\% U_{H1}$) for 5 periods $70\% U_{H1}$ (voltage dip $30\% U_{H1}$) for 25 periods $< 5\% U_{H1}$ (voltage dip $> 95\% U_{H1}$) for 5 s		The quality of electrical power on the grid is consistent with typical conditions in a commercial or hospital environment. If the user of the [ME EQUIPMENT or ME SYSTEM] needs to ensure continuous operation under conditions of possible interruptions in the mains voltage, it is recommended that the [ME EQUIPMENT or ME SYSTEM] be powered from an uninterruptible power supply or battery
Power frequency magnetic field (50/60 Hz) according to IEC 61000-4-8	3 A/m		Power frequency magnetic field levels should be consistent with typical commercial or hospital environments
Note - U_{H1} is the voltage level of the electrical network before the test impact is applied.			

Guidance and Manufacturer's Declaration - Immunity

system is intended for use in the electromagnetic environment specified below. The purchaser or user should ensure that they are used in the specified electromagnetic environment.

Immunity test	Test level according to IEC 60601	Compliance level	Electromagnetic environment - instructions
---------------	-----------------------------------	------------------	--

Conducted interference induced by radio frequency electromagnetic fields according to IEC 61000-4-6	3 Vrms in the 150 kHz to 80 MHz band outside the frequency ranges designated for industrial, scientific and medical high frequency (ISM RF) devices	3 V	The distance between the mobile radiotelephone communication systems used and any element of the Continuous Monitoring System, including cables, must be no less than the recommended spatial separation, which is calculated in accordance with the expression below in relation to the frequency of the transmitter.
RF electromagnetic field according to IEC 61000-4-3	3 V/m in the band from 80 MHz to 2.5 GHz	3 V/m	Recommended spacing: $d = 1.2 \sqrt{P}$ $d = 1.2 \sqrt{P} \quad 80 \text{ to } 800 \text{ MHz}$ $d = 2.3 \sqrt{P} \quad \text{from } 800 \text{ MHz to } 2.5 \text{ GHz}$ Note: P rated maximum output power of the transmitter, W, set by the manufacturer.; recommended spatial spacing, m. The field strength of radio waves from fixed radio transmitters, based on observations of the electromagnetic environment, should be lower than the compliance level in each frequency band. Interference may occur in the vicinity of equipment marked with
Note 1: At frequencies of 80 and 800 MHz, higher field strengths apply. Note 2: Expressions do not apply in all cases. The propagation of electromagnetic waves is affected by absorption or reflection from structures, objects and people.			
The field strength of radio waves from fixed radio transmitters, such as base stations of radiotelephone networks (cellular/wireless), and land mobile radio stations, amateur radio stations, AM and FM radio broadcast transmitters, television transmitters, cannot be determined by calculation with sufficient accuracy. To do this, practical measurements of the field strength must be carried out. If the measured values at the location of the continuous monitoring system exceed the applicable compliance levels, then the operation of the LMWH should be monitored to verify that it is functioning properly. If monitoring reveals a deviation from normal functioning, additional measures may need to be taken, such as reorienting or relocating the LPG. Outside the band from 150 kHz to 80 MHz, the field strength should be considered less than 3 V/m.			

5) Recommended separation distances between portable and mobile RF communications equipment and LUMS

Recommended separation distances between portable and mobile RF communications equipment and LUMS

The Continuous Glucose Monitoring System is intended for use in radiofrequency environments exposed to a controlled electromagnetic field. Based on the maximum output power of the communication equipment, the buyer or user can use the following method to avoid electromagnetic interference: Maintain a minimum distance between the mobile RF communication equipment (transmitter) and the continuous glucose monitoring system.

Rated maximum output power of the transmitter, W	Spatial separation d depending on transmitter frequency (m)		
	150KHz~80MHz $d=1.2VP$	80 MHz ~800 MHz $d=1.2VP$	800 MHz ~2.5 GHz $d=2.3 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.27
100	12	12	23

When determining the recommended spacing values d for transmitters with a rated maximum output power not listed in the table, substitute the rated maximum output power p in watts specified in the transmitter manufacturer's documentation into the given expressions.

Note 1: At frequencies of 80 and 800 MHz, higher field strengths apply.

Note 2: The expressions given do not apply in all cases. The propagation of electromagnetic waves is affected by absorption or reflection from structures, objects and people.

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