BE IT KNOWN that I, Laura Anne Pawley, Mezzanine Floor, Epsom Square, Epsom, Surrey KT19 8AG;


CERTIFY ONLY that Brian Michael Howes, a United Kingdom Citizen, who is well known to me, who, is duly authorised by Bryan Burton of Reckitt Benckiser Healthcare (UK) Limited ("the Company") to represent them in this matter, has today caused the annexed:

- Technical Documentation of the Medical Device – Contex® Classic natural rubber latex condoms
- Maintenance Documentation of the Medical Device – Contex® Classic natural rubber latex condoms
- Instruction for Consumers – Contex® Classic natural rubber latex condoms

For use in Russia to be produced to me and that he has represented to me on behalf of the Company that the said documents are true electronic copies of the documents.

SIGNED and sealed at Mezzanine Floor, Epsom Square, Epsom, Surrey KT19 8AG on the 15th of November 2023.

Protocol No: 2023-6781

Laura Anne Pawley
Notary Public

Dear Sir / Madam

I Bryan Burton of Reckitt Benckiser Healthcare (UK) Limited CERTIFY that the attached bundle of documents containing;

- TECHNICAL DOCUMENTATION OF THE MEDICAL DEVICE - Context® Classic natural rubber latex condoms
- MAINTENANCE DOCUMENTATION OF THE MEDICAL DEVICE - Context® Classic natural rubber latex condoms
- INSTRUCTION FOR CONSUMERS - Context® Classic natural rubber latex condoms

is a true electronic copy document for use in The Russian Federation.

Name: Bryan Burton

Signature: *Bryan Burton* Electronically signed by: Bryan Burton
Reason: I approve this document
Date: Nov 14, 2023 09:18 GMT

Date: 14-Nov-2023

Reckitt Benckiser Healthcare (UK)
Limited

140, Old Street, London, EC1A 3DF
United Kingdom
T: +44 (0)20 7654 3800
F: +44 (0)20 7654 3801
E: info@russia.reckitt.com

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Approved by: Andrew Brownhill

Global Regulatory Affairs Senior Manager

Date: 13-Nov-2023

Signature: *Andrew Brownhill*
Electronically signed by: Andrew Brownhill
Reason: I approve this document
Date: Nov 13, 2023 10:00 GMT

TECHNICAL DOCUMENTATION OF THE MEDICAL DEVICE

Contex® Classic natural rubber latex condoms

Reckitt Benckiser Healthcare (UK) Ltd., UK

Version 1.3

2023

NAME OF THE MEDICAL DEVICE

Contex® Classic natural rubber latex condoms in variations:

1. Contex® Classic natural rubber latex condoms 3 pcs. per pack.
2. Contex® Classic natural rubber latex condoms 6 pcs. per pack.
3. Contex® Classic natural rubber latex condoms 12 pcs. per pack.
4. Contex® Classic natural rubber latex condoms 18 pcs. per pack.
5. Contex® Classic natural rubber latex condoms 30 pcs. per pack.

PURPOSE OF THE MEDICAL DEVICE

Condoms are used as a barrier contraception to prevent unwanted pregnancy and to protect against HIV and other sexually transmitted diseases.

POTENTIAL CONSUMERS OF THE MEDICAL DEVICE

Medical Device is intended to be used by laypersons.

MEDICAL DEVICE POTENTIAL RISK CLASS: 2b.

STERILITY

The device is non-sterile.

MEDICAL DEVICE DESCRIPTION

Contex® Classic natural rubber latex condoms are colorless, smooth, thin, straight walled, teat ended, in silicone lubricant. Made from natural latex.

TECHNICAL CHARACTERISTICS

Technical parameter	Specification	Test Method used by the manufacturer	Frequency
Condom length, mm	180 ±10	ISO 4074:2015 clause 9.3.1 Annex D	Every batch
Condom width, mm (measured at 30 mm from the bead)	52 ±2	ISO 4074:2015 clause 9.3.2 Annex E	Every batch
Teat length, mm	Min 11	Appendix 1	Once a month
Thickness, mm	0.070+0.004/-0.008	ISO 4074:2015 clause 9.3.3 Annex F	Every batch
Bead thickness, mm	1-2	Appendix 2	Every batch
Condom weight (after removed from the primary package, W2), g	1.60 ± 0.5	Appendix 3	Once a month
Elongation, %	Median ≥ 700	ISO 4074:2015 Annex J	Once a month
Elongation, % (aged)	Median ≥ 600	ISO 4074:2015 Annex J	Once a month
Tensile strength, mPa	Median ≥ 20	ISO 4074:2015 Annex J	Once a month
Tensile strength, mPa (aged)	Median ≥ 17	ISO 4074:2015 Annex J	Once a month
Freedom from holes	Impermeable	ISO 4074:2015 clause 12 Annex M	Every batch

RAW MATERIALS (% material is based on theoretical value)

Material	Nominal content of the material, % $\pm$ 10%	Trade name
Basic material: - Compounded Natural Rubber latex	76.92	NR Latex Concentrate CL60 HA
Dusting material: - Calcium Carbonate - Magnesium Carbonate	2.41 0.15	Calcium carbonate light Bluestar Extra Light Magnesium Carbonate
Lubricant: - Silicone - Odour masker	20.46* – 20.50** 0.02* - 0.05**	BRB Silicone Oil DM 350 EP*/ Angeläy ® 350CS** Zeus Odour masking agent FLA70281*/ Transatak TAK-181256**

Raw Material	Additional Information
Compounded Natural Rubber latex	NR Latex Concentrate CL60 HA is a trade name of natural latex raw material used in compounded NRL.
Silicone	*BRB Silicone Oil DM 350 EP and **Angeläy ® 350CS are trade names of the silicone.
Odour masker	Zeus Odour masking agent FLA70281*/ Transatak TAK-181256** are trade names of the odour masker.

PACKAGING INFORMATION

Primary packaging

Package dimensions:

Individual foil: 60 $\pm$ 2 mm x 60 $\pm$ 2 mm.

Each condom is packaged in primary packaging of multi-layer material (foil).

The table below shows the materials used in the manufacture of primary packaging

Contex® Classic	Frontal foil layer: PET/AL/LDPE Back foil layer: PET/AL/LDPE
-----------------	---

From the outer to the inner foil layer:

PET = Polyethylene terephthalate (Polyethyl benzene-1,4-dicarboxylate)

AL = Aluminium

LDPE = Low Density Polyethylene

The condom in the primary packaging is packed into the secondary one – carton pack – of the quantity indicated in the secondary package.

Secondary packaging

Secondary packaging is a carton pack wrapped with the transparent tape where applicable. The quantity of condoms is indicated in the pack.

Below is information about the overall dimensions of secondary packages for all medical device variations.

1. Contex® Classic natural rubber latex condoms 3 pcs. per pack.

Secondary packaging specification

Dimensions		
	Value	Tolerance
Length, mm	64,5	+/- 1
Width, mm	19,5	+/- 1
Height, mm	74	+/- 1

2. Contex® Classic natural rubber latex condoms 6 pcs. per pack.

Secondary packaging specification

Dimensions		
	Value	Tolerance
Length, mm	63	+/- 1
Width, mm	21	+/- 1
Height, mm	124	+/- 1

3. Contex® Classic natural rubber latex condoms 12 pcs. per pack.

Secondary packaging specification

Dimensions		
	Value	Tolerance
Length, mm	63	+/- 1
Width, mm	37	+/- 1
Height, mm	124	+/- 1

4. Contex® Classic natural rubber latex condoms 18 pcs. per pack.

Secondary packaging specification

Dimensions		
	Value	Tolerance
Length, mm	82	+/- 1
Width, mm	51	+/- 1
Height, mm	124	+/- 1

5. Contex® Classic natural rubber latex condoms 30 pcs. per pack.

Secondary packaging specification

Dimensions		
	Value	Tolerance
Length, mm	95	+/- 1
Width, mm	64	+/- 1
Height, mm	124	+/- 1

PACKAGE CONTENTS OF THE MEDICAL DEVICE

The medical device is packed in a secondary / consumer package of 3 pcs. a pack, 6 pcs. a pack, 12 pcs. a pack, 18 pcs. a pack or 30 pcs. a pack. It is accompanied by the instruction for consumers that can be enclosed in the package or printed in the inside and / or outside of the package.

LABELLING INFORMATION

Each individual package is labelled with:

- trademark;
- lot number (LOT);
- expiry date (EXP)

Symbols used for labelling and what they mean:

	Expiry date
	Lot number

Additional information can be printed on the individual package including reference to international standards.

The following information is applied to the outer surface of the secondary/consumer package:

- condom description;
- quantity of condoms in consumer package;
- nominal condom length;
- nominal condom width;
- indication that condom contains lubricant;
- indication that the lubricant does not contain spermicide;
- - indication that the condom contains natural latex;
- storage guidelines;
- name and address of the manufacturer;
- country of origin (if applicable);
- - lot number;
- expiry date (EXP)

Symbols used for labelling and what they mean:

	Expiry date
	Lot number
	Contains natural latex
	For single use only
	Caution!

Additional information can be printed on the outer surface of the secondary/ consumer package including reference to international standards.

Outer or inner surface of consumer packages and/or package insert contains prescribing information, which includes:

- reminder to handle the condoms carefully, including when removing them out of the package;
- indications, contraindications, precautions for use, possible side effects;
- method of application;
- necessity to remove condom right after sexual intercourse;
- necessity to exclude petroleum- and oil-based lubricants;
- necessity to use safe additional lubricant if needed;
- necessity to consult a specialist on compatibility with topical medications;
- disposal instructions for used condoms;
- the condom is a disposable product.

INDICATIONS

- prevention of unwanted pregnancy;
- protection from HIV and other sexually transmitted diseases.

CONTRAINDICATIONS

Hypersensitivity or allergy to any of the components in condom composition.

PRECAUTIONS FOR USE

Put a condom on before the penis touches the other person. This helps to prevent pregnancy and the possibility of catching sexually transmitted infections. Check the expiry date on the condom wrapper before you use it. Do not use the condom after the expiration date. If the individual container is obviously damaged or deterioration in the properties of the condom is detected discard that condom and use a new one from an undamaged package.

Use the lubricant recommended for use with condoms. Petroleum- and oil-based lubricants (petroleum jelly, baby oil and some pessaries) and some topical medicines that are applied to the penis or vagina can damage condoms.

Stop use if you feel discomfort or irritation while using the condom.

Consult a specialist:

- about the compatibility of condoms with other contraceptives;
- if you are using topical medicines applied to penis or vagina;
- irritation/discomfort persists after use;
- you are unsure if you are allergic to latex;
- a condom leaks or bursts during use (do this as soon as possible and within 72 hours).

For anal sex, if it is necessary use additional lubrication on the outside of the condom.

Condoms may present a choking hazard. Keep out of reach of children.

If you use a condom for oral sex first, you should use a new condom for any other type of sex that follows. Use a condom only once. Reuse may increase risk of failure or infection.

None of the contraception methods can guarantee 100% protection against unwanted pregnancy, HIV and other sexually transmitted diseases.

POSSIBLE SIDE EFFECTS

Condoms are made of natural rubber latex, which, like any component of the condom, might cause allergic reactions including anaphylactic shock if user is allergic to latex.

APPLICATION

Tear the package at the toothed edge and handle the condom carefully as it can be damaged by fingernails and sharp objects.

Make sure that the condom is not inside out. Having squeezed the teat place the condom on the top of the erect penis. Roll down to its base. If the condom doesn't unroll, it may be on backwards, discard and use a new one. Stop and check if you feel the condom slipping off or tightening excessively onto the penis because this might lead to breakage.

Soon after ejaculation pull the penis out holding the condom firmly at the base of the penis and carefully take the condom off.

Dispose in a bin. Do not flush the condom down the toilet.

TRANSPORTATION

There are no special conditions required for transportation.

EXPLOITATION CONDITONS

There are no special conditions required when the medical device is used for its purpose and as indicated in the labelling.

STORAGE AND SHELF LIFE

Store the product in a dry place at 0-25 °C. Avoid mechanical impact and direct sunlight.

The condom must not come into contact with petroleum products, phenols and their derivatives, oils of petroleum origin, gasoline, kerosene and other organic substances.

Shelf life - 5 years.

RISK MANAGEMENT SUMMARY

In accordance with the BS EN ISO 14971, Medical devices - Application of risk management to medical devices, the international standard for medical device risk management an assessment of risk relating to the use of medical device has been conducted.

INFORMATION ABOUT BIOCOMPATIBILITY

Test	Standard	Biological effect
In vitro cytotoxicity testing	ISO 10993-5*	Toxic at 1:4 dilution
Maximization Sensitization	ISO 10993-10	The test articles showed no evidence of causing delayed dermal contact sensitization
Irritation	ISO 10993-10	Test articles was considered a non-irritant to vaginal tissue

Acute systemic toxicity	ISO 10993-11	There was no mortality or evidence of systemic toxicity from test article
-------------------------	--------------	---

*ISO 10993-5 suggests that a response grade of moderate or severe should be considered a cytotoxic effect. The typical condom result has been found not to be associated with any problem in the regular use of condoms by male or the female partner. If the results of other testing (typically sensitization and irritation tests) and associated evaluations are seen as acceptable then the cytotoxicity result should be considered in context and the condom considered to meet the requirements of ISO 10993. The results indicate that Context® Classic natural rubber latex can be considered to be potentially biologically safe.

MANUFACTURING PROCESS

See Appendix 4.

MEDICINAL OR PHARMACEUTICAL SUBSTANCE CONTAINED IN THE MEDICAL DEVICE

Not applicable.

LIST OF MATERIALS OF ANIMAL AND (OR) HUMAN ORIGIN

Compounded natural rubber latex contains casein (0.02% ± 10% of the contents of it, functions as a stabilizer) that is considered as a derivative from tissue/materials of animal origin as it is extracted from milk.

Casein (milk protein) is a natural by-product of milk. It has no known toxic effects although it may cause a reaction in individuals with an allergy to milk. The exposure from the condoms is considered not to pose any additional risk to users of the condoms.

The device does not contain any human blood derivatives.

STERILIZATION

The medical device is presented as a non-sterile device and is not intended to be sterilized thereafter.

SOFTWARE DESIGN, DEVELOPMENT AND VALIDATION

Not applicable.

MAINTENANCE AND REPAIR

Not applicable, as the product should be disposed of after a single use by a user.

GUIDELIENSS FOR DESTRUCTION AND DISPOSAL

The devices should be used in accordance with the instructions for use. In household use, the product is disposed of as household waste. Do not put it down the drain. In all other cases, the condoms are disposed of in accordance with local requirements.

ENVIRONMENT SAFETY REQUIREMENTS

With proper usage, storage and disposal Context® Classic natural rubber latex condoms have no negative impact on the environment.

Test Method used by the manufacturer

Appendix 1: Test method for condom teat length determination (D8291875)

- Sample preparation:
 - Dipped condom

Dry dipped condom before rolling will be used for this test. Stretch the condom slightly twice to smooth out the wrinkles.

- Condom in container
Move the condom inside the individual container such that it is away from the area where the container is to be torn. Tear the container and remove the condom. Do not use scissors or other sharp instruments to open the package.
Unroll the condom, stretch it slightly twice to smooth out the wrinkles. Lubricants may be removed by washing with a suitable solvent. If a solvent is used to remove the lubricant an adequate period of drying should be allowed before testing.
- Put the condom over the mandrel, having a rounded end with a diameter of (25 ± 2) mm with a scale divided into millimetres and dimensions with the zero beginning at the rounded end, and let it hang freely, stretched only by its own mass.
- Mark a point at the junction of condom teat and the shoulder of the condom. See Figure 1 and Picture 1.


Figure 1. Length of condom teat

Picture 1. A point marked at the junction of condom teat and the shoulder.

- Remove the marked condom out from the mandrel. Lay it flat and put the ruler on top and measure the distance between the marked point and teat end by using a calibrated ruler and record the result.

Appendix 2: Test method for measurement of condom bead thickness (LAB/LT-DMS/16, D8291143)


- Take a condom out of the pack or use a dried condom and spread it out gently.
- Measure the thickness of the bead by thickness gauge with the nearest reading of 0.01 mm.
- Repeat the measurement another 2 times.
- Average condom bead thickness for individual condom will be calculated.

Appendix 3: Test method for condom weight (after removed from the primary package) (D8291956)

- Weigh individual foiled condom and record the weight (W1).
- Weigh a dry and clean petri dish and tare zero.
- Move the condom inside the individual package such that it is away from the area where the package is to be torn. Tear the package. Do not use scissors or other sharp instruments to open the package.

- Remove the condom carefully out of the opened foil package by using a dry and clean forceps and place it on the petri dish. Record the weight of condom (W2).
- Ensure that petri dish and forceps are dry and clean before testing next sample.
- Weight of foil package and remaining lubricant and powder after condom removal (W3) = W1 - W2.

Appendix 4

OVERVIEW PROCESS CONTROL FLOWCHART

TF Contex Classic_LETO_new pack size variation_RTQ2_KH reviewed+ amended from LRA

Final Audit Report

2023-11-13

Created:	2023-11-13
By:	Beth Hefferan (Beth.Hefferan@reckitt.com)
Status:	Signed
Transaction ID:	CBJCHBCAA3AA5hMvhFiBxjo4pmP31NEhYoBPZiyFUCN0

"TF Contex Classic_LETO_new pack size variation_RTQ2_KH reviewed+ amended from LRA" History

-  Document created by Beth Hefferan (Beth.Hefferan@reckitt.com)
2023-11-13 - 09:51:16 GMT- IP address: 147.161.225.84
-  Document emailed to Andrew Brownhill (Andrew.Brownhill@reckitt.com) for signature
2023-11-13 - 09:52:02 GMT
-  Email viewed by Andrew Brownhill (Andrew.Brownhill@reckitt.com)
2023-11-13 - 10:02:33 GMT- IP address: 147.161.237.27
-  Document e-signed by Andrew Brownhill (Andrew.Brownhill@reckitt.com)
Signing reason: I approve this document.
Signature Date: 2023-11-13 - 10:06:45 GMT - Time Source: server- IP address: 147.161.237.27
-  Agreement completed.
2023-11-13 - 10:06:45 GMT
-  Andrew Brownhill (Andrew.Brownhill@reckitt.com) authenticated with Adobe Acrobat Sign.
2023-11-13 - 10:06:45 GMT

Powered by
Adobe
Acrobat Sign

Approved by: Andrew Brownhill

Global Regulatory Affairs Senior Manager

Date: 13-Nov-2023

Signature: *Andrew Brownhill*

Electronically signed by: Andrew
Brownhill
Reason: I approve this document
Date: Nov 13, 2023 10:02 GMT

MAINTENANCE DOCUMENTATION OF THE MEDICAL DEVICE

Contex® Classic natural rubber latex condoms

Reckitt Benckiser Healthcare (UK) Ltd., UK

Version 1.3

NAME OF THE MEDICAL DEVICE

Contex® Classic natural rubber latex condoms in variations:

1. Contex® Classic natural rubber latex condoms 3 pcs. per pack.
2. Contex® Classic natural rubber latex condoms 6 pcs. per pack.
3. Contex® Classic natural rubber latex condoms 12 pcs. per pack.
4. Contex® Classic natural rubber latex condoms 18 pcs. per pack.
5. Contex® Classic natural rubber latex condoms 30 pcs. per pack.

MANUFACTURER

Reckitt Benckiser Healthcare (UK) Ltd., 103-105 Bath Road, Slough, Berkshire, SL1 3 UH, UK.

MANUFACTURING SITE

SSL Manufacturing (Thailand) Ltd., Wellgrow Industrial Estate 96 & 100, Moo 5, Bangna-Trad Road KM.36, Bangsamak, Bangpakong, Chachoengsao, 24180, Thailand;

AUTHORIZED REPRESENTATIVE OF THE MANUFACTURER

Reckitt Benckiser Healthcare LLC, Floor 3, 4 Shlyuzavaya emb., Moscow, 115114, Russia.
Phone: +7 (495) 961 25 65; fax: +7 (495) 961 25 66

PURPOSE OF THE MEDICAL DEVICE

Condoms are used as a barrier contraception to prevent unwanted pregnancy and to protect against HIV and other sexually transmitted diseases.

POTENTIAL CONSUMERS OF THE MEDICAL DEVICE

Medical Device is intended to be used by laypersons.

MEDICAL DEVICE DESCRIPTION

Contex® Classic natural rubber latex condoms are colorless, smooth, thin, straight walled, teat ended, in silicone lubricant. Made from natural latex.

TECHNICAL CHARACTERISTICS

Technical parameter	Specification
Condom length, mm	180 ± 10
Condom width, mm (measured at 30 mm from the bead)	52 ± 2
Teat length, mm	Min 11
Thickness, mm	0.070+0.004/-0.008
Bead thickness, mm	1-2
Condom weight (after removed from the primary package, W2), g	1.60 ± 0.5
Elongation, %	Median ≥ 700
Elongation, % (aged)	Median ≥ 600
Tensile strength, mPa	Median ≥ 20
Tensile strength, mPa (aged)	Median ≥ 17

Freedom from holes	Impermeable
--------------------	-------------

INDICATIONS

- prevention of unwanted pregnancy;
- protection from HIV and other sexually transmitted diseases.

CONTRAINDICATIONS

Hypersensitivity or allergy to any of the components in condom composition.

PRECAUTIONS FOR USE

Put a condom on before the penis touches the other person. This helps to prevent pregnancy and the possibility of catching sexually transmitted infections. Check the expiry date on the condom wrapper before you use it. Do not use the condom after the expiration date. If the individual container is obviously damaged or deterioration in the properties of the condom is detected discard that condom and use a new one from an undamaged package. Use the lubricant recommended for use with condoms. Petroleum- and oil-based lubricants (petroleum jelly, baby oil and some pessaries) and some topical medicines that are applied to the penis or vagina can damage condoms.

Stop use if you feel discomfort or irritation while using the condom.

Consult a specialist if:

- about the compatibility of condoms with other contraceptives;
- if you are using topical medicines applied to penis or vagina;
- irritation/discomfort persists after use;
- you are unsure if you are allergic to latex;
- a condom leaks or bursts during use (do this as soon as possible and within 72 hours).

For anal sex, if it is necessary use additional lubrication on the outside of the condom.

Condoms may present a choking hazard. Keep out of reach of children.

If you use a condom for oral sex first, you should use a new condom for any other type of sex that follows. Use a condom only once. Reuse may increase risk of failure or infection.

None of the contraception methods can guarantee 100% protection against unwanted pregnancy, HIV and other sexually transmitted diseases.

POSSIBLE SIDE EFFECTS

Condoms are made of natural rubber latex, which, like any component of the condom, might cause allergic reactions including anaphylactic shock if user is allergic to latex.

APPLICATION

Tear the package at the toothed edge and handle the condom carefully as it can be damaged by fingernails and sharp objects.

Make sure that the condom is not inside out. Having squeezed the teat place the condom on the top of the erect penis. Roll down to its base. If the condom doesn't unroll, it may be on backwards, discard and use a new one. Stop and check if you feel the condom slipping off or tightening excessively onto the penis because this might lead to breakage.

Soon after ejaculation pull the penis out holding the condom firmly at the base of the penis and carefully take the condom off.

Dispose in a bin. Do not flush the condom down the toilet.

STORAGE CONDITIONS AND SHELF LIFE

Store the product in a dry place at 0-25°C. Avoid mechanical impact and direct sunlight.

The condom must not come into contact with petroleum products, phenols and their derivatives, oils of petroleum origin, gasoline, kerosene and other organic substances.

Shelf life - 5 years.

MEDICINAL OR PHARMACEUTICAL SUBSTANCE CONTAINED IN THE MEDICAL DEVICE

Not applicable.

LIST OF MATERIALS OF ANIMAL AND (OR) HUMAN ORIGIN

Compounded natural rubber latex contains casein (0.02% ± 10% of the contents of it), functions as a stabilizer that is considered as a derivative from tissue/materials of animal origin as it is extracted from milk.

Casein (milk protein) is a natural by-product of milk. It has no known toxic effects although it may cause a reaction in individuals with an allergy to milk. The exposure from the condoms is considered not to pose any additional risk to users of the condoms.

The device does not contain any human blood derivatives.

STERILIZATION

The medical device is presented as a non-sterile device and is not intended to be sterilized thereafter.

USAGE WITH OTHER MEDICAL DEVICES

Condoms can be used on their own to prevent from STIs and pregnancy and aren't intended to be used in conjunction with other medical devices.

However, consumers can use condoms with condom compatible lubricants if needed.

INSTALLATION, MAINTENANCE, CALIBRATION

Not applicable.

CLEANING, DISINFECTION, STERILISATION OF THE MEDICAL DEVICE

Not applicable.

INTERNATIONAL STANDARDS

ISO 13485

ISO 14971

ISO 4074

PRECAUTIONS TO BE TAKEN BY USER RELATED TO THE DISPOSAL OF THE DEVICE

Context[®] Classic natural rubber latex condoms should be disposed in a bin. The condoms should not be flushed the condom down the toilet.

DATE OF THE ISSUE / THE LATEST REVISION OF THE MAINTENANCE DOCUMENTATION

11/2022.

GUIDELIENSS FOR DESTRUCTION AND DISPOSAL

The devices should be used in accordance with the instructions for use.

In household use, the product is disposed of as household waste. Do not put it down the drain.

In all other cases, the condoms are disposed of in accordance with local requirements.

ENVIRONMENT SAFETY REQUIREMENTS

With proper usage, storage and disposal Context[®] Classic natural rubber latex condoms have no negative impact on the environment.

MD_Context Classic_LETO_new pack size variation_RTQ2_KH reviewed+ LRA amended

Final Audit Report

2023-11-13

Created:	2023-11-13
By:	Beth Hefferan (Beth.Hefferan@reckitt.com)
Status:	Signed
Transaction ID:	CBJCHBCAABAAqJn9UQI1pUq6TLC-woVs7c1yKfUVnQ8o

"MD_Context Classic_LETO_new pack size variation_RTQ2_KH reviewed+ LRA amended" History

-  Document created by Beth Hefferan (Beth.Hefferan@reckitt.com)
2023-11-13 - 09:47:59 GMT - IP address: 147.161.225.84
-  Document emailed to Andrew Brownhill (Andrew.Brownhill@reckitt.com) for signature
2023-11-13 - 09:48:42 GMT
-  Email viewed by Andrew Brownhill (Andrew.Brownhill@reckitt.com)
2023-11-13 - 09:58:17 GMT - IP address: 147.161.237.27
-  Andrew Brownhill (Andrew.Brownhill@reckitt.com) authenticated with Adobe Acrobat Sign.
2023-11-13 - 10:02:17 GMT
-  Document e-signed by Andrew Brownhill (Andrew.Brownhill@reckitt.com)
Signing reason: I approve this document.
Signature Date: 2023-11-13 - 10:02:18 GMT - Time Source: server- IP address: 147.161.237.27
-  Agreement completed.
2023-11-13 - 10:02:18 GMT

Powered by
Adobe
Acrobat Sign

Approved by: Andrew Brownhill

Global Regulatory Affairs Senior Manager

Date: 13-Nov-2023

Signature: *Andrew Brownhill*

Electronically signed by: Andrew
Brownhill
Reason: I approve this document.
Date: Nov 13, 2023 09:57 GMT

INSTRUCTION FOR CONSUMERS
Contex® Classic natural rubber latex condoms

Reckitt Benckiser Healthcare (UK) Ltd., UK

Version 1.3

2023

Read the instruction carefully before use, especially if condoms are used for oral or anal sex.

Put a condom on before the penis touches the other person. This helps to prevent pregnancy and the possibility of catching sexually transmitted infections.

ATTENTION: Check the expiry date on the condom wrapper before you use it. If the individual container is obviously damaged or deterioration in the properties of the condom is detected discard that condom and use a new one from an undamaged package.

PURPOSE

Condoms are used as a barrier contraception to prevent unwanted pregnancy and to protect against HIV and other sexually transmitted diseases.

APPLICATION

1. Tear the package at the toothed edge and handle the condom carefully as it can be damaged by fingernails and sharp objects.

2. Make sure that the condom is not inside out. Having squeezed the teat place the condom on the top of the erect penis.

3. Roll down to its base. If the condom doesn't unroll, it may be on backwards, discard and use a new one. Stop and check if you feel the condom slipping off or tightening excessively onto the penis because this might lead to breakage.

Soon after ejaculation pull the penis out holding the condom firmly at the base of the penis and carefully take the condom off.

4. Dispose in a bin. Do not flush the condom down the toilet.

For anal sex or if it is necessary use additional lubrication on the outside of the condom.

Condoms may present a choking hazard. Keep out of reach of children.

If you use a condom for oral sex first, you should use a new condom for any other type of sex that follows.

Use a condom only once. Reuse may increase risk of failure or infection.

CONTRAINDICATIONS

Hypersensitivity or allergy to any of the components in condom composition.

POSSIBLE SIDE EFFECTS

Condoms are made of natural rubber latex, which might cause allergic reactions including anaphylactic shock if user is allergic to latex.

PRECAUTIONS FOR USE

Use the lubricant recommended for use with condoms. Petroleum- and oil-based lubricants (petroleum jelly, baby oil), some pessaries and some topical medicines that are applied to the penis or vagina can damage condoms.

Stop use if you feel discomfort or irritation while using the condom.

Consult a specialist:

- about the compatibility of condoms with other contraceptives;
- if you are using topical medicines applied to penis or vagina;
- if irritation/discomfort persists after use;
- if you are unsure if you are allergic to latex;
- if a condom leaks or bursts during use (do this as soon as possible and within 72 hours).

None of the contraception methods can guarantee 100% protection against unwanted pregnancy, HIV and other sexually transmitted diseases.

Date of the latest revision of the instruction: 11/2022.

IFU_ Contex Classic_LETO_new pack size variation_RTQ2_LRA amended


Final Audit Report

2023-11-13


Created:	2023-11-13
By:	Beth Hefferan (Beth.Hefferan@reckitt.com)
Status:	Signed
Transaction ID:	CBJCHBCAABAAICzmzw91XRYmnWQjTcHacu7O3H-7OsgS

"IFU_ Contex Classic_LETO_new pack size variation_RTQ2_LRA amended" History

-  Document created by Beth Hefferan (Beth.Hefferan@reckitt.com)
2023-11-13 - 09:46:59 GMT- IP address: 147.161.225.84
-  Document emailed to Andrew Brownhill (Andrew.Brownhill@reckitt.com) for signature
2023-11-13 - 09:47:41 GMT
-  Email viewed by Andrew Brownhill (Andrew.Brownhill@reckitt.com)
2023-11-13 - 09:56:28 GMT- IP address: 147.161.237.27
-  Andrew Brownhill (Andrew.Brownhill@reckitt.com) authenticated with Adobe Acrobat Sign.
2023-11-13 - 09:57:47 GMT
-  Document e-signed by Andrew Brownhill (Andrew.Brownhill@reckitt.com)
Signing reason: I approve this document.
Signature Date: 2023-11-13 - 09:57:48 GMT - Time Source: server- IP address: 147.161.237.27
-  Agreement completed.
2023-11-13 - 09:57:48 GMT

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