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APOSTILLE

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1. Ülke/Country/Pays/Staat TÜRKİYE - LA TURQUIE

İşbu resmî belge/This public document/Le présent acte public/Dieses zeugnis wurde

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
9. Mühür - Damga/Seal-stamp/Sceau-
timbre/Siegel-Stempel

10. İmza/Signature/Signature/Unterschrift:

2


Handwritten signature of İsmail KARA.

KAPSAM SAĞLIK ÜRÜNLERİ İNŞAAT TAAHHÜT TURİZM VE TİCARET LİMİTED ŞİRKETİ is a registered member of our chamber since 25-12-2013, under reg. No: 24314

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INSTRUCTIONS FOR USE

for medical device

EVA infusion system (bag) for administration of intravenous solutions
and parenteral nutrition

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1. Medical device name

EVA infusion system (bag) for the administration of intravenous solutions and parenteral nutrition, versions.

EVA infusion system (bag) for the administration of intravenous solutions and parenteral nutrition, versions.

1. EVA Infusion system (transparent bag) for the administration of intravenous solutions and parenteral nutrition, contents:

1.1. Transparent EVA bag for the administration of intravenous solutions and parenteral nutrition, volume: 150 ml. - 1 pc.

- Luer-lock stopper with filter - 1 pc.

- Instructions for use - 1 pc.

1.2. Transparent EVA bag for the administration of intravenous solutions and parenteral nutrition, volume: 250 ml. - 1 pc.

- Luer-lock stopper with filter - 1 pc.

- Instructions for use - 1 pc.

1.3. Transparent EVA bag for the administration of intravenous solutions and parenteral nutrition, volume: 500 ml. - 1 pc.

- Luer-lock stopper with filter - 1 pc.

- Instructions for use - 1 pc.

1.4. Transparent EVA bag for the administration of intravenous solutions and parenteral nutrition, volume: 1000 ml. - 1 pc.

- Luer-lock stopper with filter - 1 pc.

- Instructions for use - 1 pc.

1.5. Transparent EVA bag for the administration of intravenous solutions and parenteral nutrition, volume: 2000 ml. - 1 pc.

- Luer-lock stopper with filter - 1 pc.

- Instructions for use - 1 pc.

1.6. Transparent EVA bag for the administration of intravenous solutions and parenteral nutrition, volume: 3000 ml. - 1 pc.

- Luer-lock stopper with filter - 1 pc.

- Instructions for use - 1 pc.

1.7. Transparent EVA bag for the administration of intravenous solutions and parenteral nutrition, volume: 4000 ml. - 1 pc.

- Luer-lock stopper with filter - 1 pc.

- Instructions for use - 1 pc.

2. EVA Infusion system (light-protected bag) for the administration of intravenous solutions and parenteral nutrition, contents:

2.1. Light-protected EVA bag for the administration of intravenous solutions and parenteral nutrition, volume: 150 ml. - 1 pc.

- Luer-lock stopper with filter - 1 pc.

- Instructions for use - 1 pc.

2.2. Light-protected EVA bag for the administration of intravenous solutions and parenteral nutrition, volume: 250 ml. - 1 pc.

- Luer-lock stopper with filter - 1 pc.

- Instructions for use - 1 pc.

2.3. Light-protected EVA bag for the administration of intravenous solutions and parenteral nutrition, volume: 500 ml. - 1 pc.

- Luer-lock stopper with filter - 1 pc.

- Instructions for use - 1 pc.

2.4. Light-protected EVA bag for the administration of intravenous solutions and parenteral nutrition, volume: 1000 ml. - 1 pc.

- Luer-lock stopper with filter - 1 pc.

- Instructions for use - 1 pc.

2.5. Light-protected EVA bag for the administration of intravenous solutions and parenteral nutrition, volume: 2000 ml. - 1 pc.

- Luer-lock stopper with filter - 1 pc.

- Instructions for use - 1 pc.

2.6. Light-protected EVA bag for the administration of intravenous solutions and parenteral nutrition, volume: 3000 ml. - 1 pc.

- Luer-lock stopper with filter - 1 pc.

- Instructions for use - 1 pc.

2.7. Light-protected EVA bag for the administration of intravenous solutions and parenteral nutrition, volume: 4000 ml. - 1 pc.

- Luer-lock stopper with filter - 1 pc.

- Instructions for use - 1 pc.

(hereinafter referred to as EVA infusion system (bag) for the administration of intravenous solutions and parenteral nutrition, EVA bag for the administration of intravenous solutions and parenteral nutrition, EVA bag, medical device, product).

2. Information about the developer, manufacturer, and production sites

Medical product manufacturer:

Kapsam Saglik Urunleri Insaat Taahhut Turizm Ve Ticaret Limited Sirketi (KAPSAM HEALTH PRODUCTS).

Address - CUMHURIYET MH. ISMET INONU CD. NO: 57 CAYIROVA / Gebze - Kocaeli - Turkey.

Medical Device Developer:

Kapsam Saglik Urunleri Insaat Taahhut Turizm Ve Ticaret Limited Sirketi (KAPSAM HEALTH PRODUCTS).

Address - CUMHURIYET MH. ISMET INONU CD. NO: 57 CAYIROVA / Gebze - Kocaeli - Turkey.

Device production site:

Kapsam Saglik Urunleri Insaat Taahhut Turizm Ve Ticaret Limited Sirketi (KAPSAM HEALTH PRODUCTS).

Address - CUMHURIYET MH. ISMET INONU CD. NO: 57 CAYIROVA / Gebze - Kocaeli - Turkey.

3. Information about the authorized representative in the Russian Federation

MEDSERVICE Limited Liability Company (MEDSERVICE LLC)

127299, Moscow, internal territory of the Kopteevo municipality, Bolshaya Akademicheskaya str., 4, office 5A

4. Intended use of the medical device

EVA infusion system (bag) is designed for mixing and administration of intravenous solutions and parenteral nutrition.

5. Application area

The medical device is used in outpatient or inpatient settings.

6. Indications for use of the medical device

The medical device is used in a situation of prolonged prohibition of food intake and with increased metabolism, when the patient's nutrient solution must be fed parenterally and infusion of metered-mixed drug solutions is required.


7. Contraindications for the medical device

This medical device has no documented contra-indications when used in accordance with the operational documentation.

8. Possible side effects

Individual intolerance to the medical device components.

9. Method of using the medical device

The EVA bag has 3 ports:

A - Luer-Lock inlet port;

B - Outlet port (screw-on);

C - Injection port with cover.

A - Luer-Lock inlet port

The port is used with an automated compounding device to transfer sterile solutions.

Step 1

Insert the luer-lock connector from the transfusion system into this port. When the bag is filled, close the safe lock cap.

Step 2


Disconnect the infusion system's luer-lock from the port. Close the port by attaching an additional stopper, that is sterile and packaged with the light-protected/transparent EVA bag.

B - Output port (threaded)

Step 1:
Detach the cap by
unscrewing it.

Step 2:
Insert the IV tubing
spike into the port of
the bag.



Step 3:

Administer the solution into the patient's body appropriately.

C - Injection port with a cap.

This port is used for manual administration of sterile solutions.

Step 1:

Detach / break the safe-lock cap.

Step 2:

Inject the solution into the EVA bag manually, with the use of an appropriate syringe.

Step 3:

When the injection process is finished, close the safe-lock cap.

10. Information about potential consumers

The medical device must be used by qualified, specially trained personnel.


11. Product description

EVA infusion system (bag) for the administration of intravenous solutions and parenteral nutrition consists of:

- a transparent/light-protected EVA bag consisting of: 25-28 mm/40-43 mm/110 mm tube, exit port (screw-type), injection port with lid, inlet port with an internal Luer-lock fitting, clamp, stopper for Luer-lock with filter;

The bag is filled with a nutrient solution or medication through the inlet port. If necessary, liquid medication can be added to the nutrient solution bag through a closed injection port. The bag is then suspended from the infusion rack, after which it is connected to the single-use IV system using the outlet port (screw-type), and the IV system flow regulator is opened to release air. The system is connected to a peripheral venous or central venous catheter, and parenteral administration of a nutrient solution to the patient is performed.


Available sizes of light-protected/transparent EVA bag for the administration of intravenous solutions and parenteral nutrition: 150 ml, 250 ml, 500 ml, 1,000 ml, 2,000 ml, 3000 ml, 4,000 ml

No.	Component name
1	Transparent bag for total parenteral nutrition
2	25-28 mm tube

No.	Component name
3	40-43 mm tube
4	110 mm tube
5	Outlet port (screw-type)
6	Injection port with cover
7	Inlet port with a an internal Luer-lock fitting
8	Clamp
9	Plug for Luer-lock stopper with filter

Available sizes of transparent EVA bags for administration of intravenous solutions and parenteral nutrition; 150 ml, 250 ml, 500 ml, 1,000 ml, 2,000 ml, 3000 ml, 4,000 ml

No.	Component name
1	Transparent EVA bag for intravenous solutions and parenteral nutrition
2	25-28 mm tube
3	40-43 mm tube

No.	Component name
4	110 mm tube
5	Outlet port (screw-type)
6	Injection port with cover
7	Inlet port with a an internal Luer-lock fitting
8	Clamp
9	Plug for Luer-lock stopper with filter

13. Technical specifications of the medical device

Transparent polymeric bag for parenteral nutrition solution, 150 ml, 250 ml, 500 ml, 1,000 ml, 2,000 ml, 3,000 ml, 4,000 ml

Volume	Maximum volume	Weight	Overall size
150 ml	380 ml	38.2 g	120 x 195 mm
250 ml	380 ml	38.8 g	120 x 195 mm
500 ml	680 ml	43.7 g	120 x 260 mm
1,000 ml	1,300 ml	52.4 g	160 x 280 mm
2,000 ml	2,450 ml	81.0 g	204 x 360 mm
3,000 ml	3,800 ml	104.9 g	245 x 380 mm
4,000 ml	4,900 ml	112.4 g	260 x 420 mm

Note: Permissible deviation from nominal values $\pm 5\%$

Polymeric light-protected bag for parenteral nutrition solution, 150 ml, 250 ml, 500 ml, 1,000 ml, 2,000 ml, 3,000 ml, 4,000 ml

Volume	Maximum volume	Weight	Overall size
150 ml	380 ml	35.4 g	120 x 195 mm
250 ml	380 ml	38.5 g	120 x 195 mm
500 ml	680 ml	43.7 g	120 x 260 mm

1,000 ml	1,300 ml	52.8 g	160 x 280 mm
2,000 ml	2,450 ml	80.6 g	204 x 360 mm
3,000 ml	3,800 ml	108.2 g	245 x 380 mm
4,000 ml	4,900 ml	110.0 g	260 x 420 mm

Note: Permissible deviation from nominal values $\pm 5\%$

The EVA bag for administration of intravenous solutions and parenteral nutrition has a hole for hanging from the rack; when the bag is filled to its rated capacity and suspended, and when a static force of 15 N is applied for 24 hours, no sagging or tearing of the bag occurs.

Graduation of measuring instruments:

The graduation of the bag and the units of measurement on the bag are clear and accurate. When the bag is filled to the nominal capacity mark in the suspended state, the permissible error of graduation is not more than $\pm 10\%$ of the specified value.

Volume between release lines with numbers:

- 50 ml for 150 ml and 250 ml bags
- 100 ml for 500 ml bags
- 200 ml for 1,000 ml bags
- 250 ml for 2,000 ml and 3,000 ml bags
- 500 ml for 4,000 ml bags.

Transparent EVA bags for administration of intravenous solutions and parenteral nutrition

150 ml, 250 ml, 500 ml, 1,000 ml, 2,000 ml, 3,000 ml, 4,000 ml

The bag is soft, smooth and without obvious foreign inclusions. With normal or corrected vision, the passage of bubbles through it as well as the interface between the liquid and the air can be seen. The heat-sealing line is transparent and uniform.

25-28 mm tube

The tube is transparent, with normal or corrected vision the passage of bubbles through it, as well as the interface between the liquid and the air, are visible.

Tube length from 25 to 28 mm; tube wall thickness from 0.90 to 1.10 mm; outer diameter from 8.14 to 8.30 mm.

40-43 mm tube

The tube is transparent, with normal or corrected vision the passage of bubbles through it, as well as the interface between the liquid and the air, are visible.

Tube length from 40 to 43 mm; tube wall thickness from 0.90 to 1.10 mm; outer diameter from 7.70 to 8.14 mm.

110 mm tube

The tube is transparent, with normal or corrected vision the passage of bubbles through it, as well as the interface between the liquid and the air, are visible.

Tube length from 110 to 120 mm; tube wall thickness from 0.97 to 1.12 mm; outer diameter from 6.55 to 6.75 mm.

Exit port (twist off)

After piercing the hole with an IV system needle, the bag is not damaged. There is no leakage of fluid from the hole during operation.

Injection port with cover

After puncturing the hole with the injection needle, the bag is not damaged. During operation there is no fluid leakage from the hole.

Inlet port with a an internal "Luer-screw connection" fitting

The open tube port has a locking connector.

Fitting dimensions with internal "Luer-screw connection":

Outer diameter of the tube inlet	5.2 ± 0.1 mm
Inner diameter of the tube inlet	4.025 ± 0.125 mm
Inlet tube cone length	13.00 (-0/+0.50) mm
Length from the tube inlet to the Luer-lock	13.00 (-0/+0.50) mm
External diameter of the Luer-lock	11.20(±0.1) mm
Inner diameter of the Luer-lock	7.0(±0.1) mm
Maximum inside diameter of the male thread	7.10 (-0/+0.10) mm
Outer diameter without reducer	10.00(+0/-0.10) mm
Cone	6°

Clamp

The surface is clean, free of dirt and processing defects.

When closed, the clamp blocks the flow, and when not closed, in the open position, it allows the flow to move freely.

When the air-filled bag clamp is closed, there is no air leakage at 3 bar.

Luer-lock stopper with filter

The cone-shaped connector with an attached filter and threaded Luer-Lock has bumps on the side for better grip.

Plug dimensions:

Outer diameter of the outer cone:	11.05 ± 0.65 mm
Inner diameter of the outer cone:	6.15 ± 0.75 mm
Outer length of the cap	21,00 (-0/+0.10) mm
Inner cone length	19,90 (-0/+0.10) mm
Cone	6%

The pore size of the filter is 3 µm. The effective filtration area is 35.25 mm². The bacterial filtration efficiency is 99.95%.

Light-protected EVA bags for administration of intravenous solutions and parenteral nutrition
150 ml, 250 ml, 500 ml, 1,000 ml, 2,000 ml, 3,000 ml, 4,000 ml

The bag is soft, smooth and without obvious foreign inclusions. With normal or corrected vision, the passage of bubbles through it as well as the interface between the liquid and the air can be seen. The heat-sealing line is transparent and uniform.

The bag provides 89-99% protection from light in the range of 290 to 450 nm, blocking the harmful effects of daylight.

25-28 mm tube

The tube is yellow, with normal or corrected vision the passage of bubbles through it, as well as the interface between the liquid and the air can be seen.

Tube length from 25 to 28 mm; tube wall thickness from 0.90 to 1.10 mm; outer diameter from 8.14 to 8.30 mm.

40-43 mm tube

The tube is yellow, with normal or corrected vision the passage of bubbles through it, as well as the interface between the liquid and the air can be seen.

Tube length from 40 to 43 mm; tube wall thickness from 0.90 to 1.10 mm; outer diameter from 7.70 to 8.14 mm.

110 mm tube

The tube is yellow, with normal or corrected vision the passage of bubbles through it, as well as the interface between the liquid and the air can be seen.

Tube length from 105 to 120 mm; tube wall thickness from 0.97 to 1.12 mm; outer diameter from 6.55 to 6.75 mm.

Outlet port (screw-type)

After piercing the hole with an IV system needle, the bag is not damaged. During operation there is no fluid leakage from the hole.

Injection port with cover

After puncturing the hole with the injection needle, the bag is not damaged. During operation there is no fluid leakage from the hole.

Inlet port with a an internal "Luer-screw connection" fitting

The open tube port has a locking connector.

Fitting dimensions with internal "Luer-screw connection":

Outer diameter of the tube inlet	5.2 ± 0.1 mm
Inner diameter of the tube inlet	4.025 ± 0.125 mm
Inlet tube cone length	13.00 (-0/+0.50) mm
Length from tube inlet to Luer-lock	13.00 (-0/+0.50) mm
Outer diameter of the Luer-lock	11.20(±0.1) mm
Inner diameter of the Luer-lock	7.0(±0.1) mm
Maximum inside diameter of the male thread	7.10 (-0/+0.10) mm
Outer diameter without reducer	10.00(+0/-0.10) mm
Cone	6%

Clamp

The surface is clean, free of dirt and processing defects.

When closed, the clamp blocks the flow, and when not closed, in the open position, it allows the flow to move freely.

When the air-filled bag clamp is closed, there is no air leakage at 3 bar.

Luer-lock stopper with filter

The cone-shaped connector with an attached filter and threaded Luer-Lock has bumps on the side for better grip.

Plug dimensions:

Outer diameter of the outer cone:	11.05 ± 0.65 mm
Inner diameter of the outer cone:	6.15 ± 0.75 mm
Outer length of the cap	21,00 (-0/+0.10) mm
Inner cone length	19,90 (-0/+0.10) mm
Cone	6%

The pore size of the filter is 3 μ m. The effective filtration area is 35.25 mm². The bacterial filtration efficiency is 99.95%.

General requirements

Bioburden

Total number of bacteria not more than 10 CFU before packaging in bulk product; coliform bacteria should be 0 CFU.

Leakage

If an air pressure of 2 bar is maintained inside the bag and it is held under water for 2 minutes, no air leakage occurs.

Joint strength

Bag component joints can withstand a static load of 15 N for 2 minutes without dropping or breaking.

Resistance force at fluid outlet

When the bag is filled to its nominal capacity, the liquid is drained without applying pressure or air exchange.

14. Medical device materials

Name of component	Material name	Material class	Form of contact with the patient	Manufacturer of material
Light-protected bag for total parenteral nutrition	Ethylene vinyl acetate	EVA4500	Short-term Indirect contact with the blood circulation system through the solution.	RAUMEDIC AG (Germany)
	Yellow dye	93055 yellow transparent		RAUMEDIC AG (Germany)
Transparent bag for total parenteral nutrition	Ethylene vinyl acetate	EVA4500	Short-term Indirect contact with the blood circulation system through the solution.	RAUMEDIC AG (Germany)

Name of component	Material name	Material class	Form of contact with the patient	Manufacturer of material
25-28 mm tube (Light-protected bag for total parenteral nutrition)	Ethylene vinyl acetate + PVC or Polyurethane	EVA4500 + PVC NO DOP / NEUTRO	Short-term Indirect contact with the blood circulation system through the solution.	RAUMEDIC AG (Germany) + MEDIFLY S.R.L. (Italy)
	Yellow dye	93055 yellow transparent		RAUMEDIC AG (Germany)
25-28 mm tube (Transparent bag for total parenteral nutrition)	Ethylene vinyl acetate + PVC or Polyurethane	EVA MEDICAL + PVC NO DOP / NEUTRO	Short-term Indirect contact with the blood circulation system through the solution.	RAUMEDIC AG (Germany)+ MEDIFLY S.R.L.
40-43 mm tube (Light-protected bag for total parenteral nutrition)	Ethylene vinyl acetate + PVC or Polyurethane	EVA MEDICAL + PVC NO DOP / AMBRATO	Short-term Indirect contact with the blood circulation system through the solution.	RAUMEDIC AG (Germany) + MEDIFLY S.R.L. (Italy)
	Yellow dye	93055 yellow transparent		RAUMEDIC AG (Germany)
40-43 mm tube (Transparent bag for total parenteral nutrition)	Ethylene vinyl acetate + PVC or Polyurethane	EVA MEDICAL + PVC NO DOP / NEUTRO	Short-term Indirect contact with the blood circulation system through the solution.	RAUMEDIC AG (Germany) + MEDIFLY S.R.L. (Italy)

Name of component	Material name	Material class	Form of contact with the patient	Manufacturer of material
110 mm tube (Light-protected bag for total parenteral nutrition)	PVC	AM14TW502	Short-term Indirect contact with the blood circulation system through the solution.	TPV Compound Srl (Italy)
110 mm tube (Transparent bag for total parenteral nutrition)	PVC	AM69M	Short-term Indirect contact with the blood circulation system through the solution.	TPV Compound Srl (Italy)
Outlet port (screw-type)	PVC	RB 3/D NDG	Short-term Indirect contact with the blood circulation system through the solution.	MODENPLAST MEDICAL SRL (Italy)

Name of component	Material name	Material class	Form of contact with the patient	Manufacturer of material
Injection port with cover	Polypropylene + Polyisoprene + Acrylonitrile-butadiene styrene	MOPLEN HP 456H + 70D16D3 + TERLUX 2802 HD	Short-term Indirect contact with the blood circulation system through the solution.	Basell Sales & Marketing Company BV (Netherlands) + Datwyler Pharma Packaging International NV (Belgium) + INEOS Styrolution Group GmbH (Germany)
Inlet port with an internal Luer-lock fitting	Acrylonitrile-butadiene styrene	LUSTRAN 348	Short-term Indirect contact with the blood circulation system through the solution.	INEOS Styrolution Group GmbH (Germany)
Clamp	Polypropylene	K-RESIN KP 03	No contact	INEOS Styrolution Group GmbH (Germany)

Name of component	Material name	Material class	Form of contact with the patient	Manufacturer of material
Luer-lock stopper with filter	Acrylonitrile-butadiene styrene + Polypropylene + Acrylic copolymer	LUSTRAN 348 + PPC 10712 + VERSAPOR 3000R	No contact	INEOS Styrolution Group GmbH (Germany) + TOTALENERGIE S BRANCHE RAFFINAGE CHIMIE + Pall Corporation (France)

15. Description of accessories, medical or non-medical items intended to be used in conjunction with the claimed medical devices

- For the input port with a Luer-Lock connector - any duly registered automated compounding device for transfusion of sterile solutions, with a large port, that is compatible with the female threads, with the male end for the fill/inlet port (i.e., KS-124)
- For the exit port (threaded) - any duly registered sterile transfusion system with a 5-5.5 mm needle tip (must be spike-structured) for the exit port (threaded) (e.g., KS-IV501)
- For the injection port, with a cap - Any duly registered syringe with an injection port needle (any needle under 16G)

16. List of international and national regulatory documents and standards to which the medical device complies

DOCUMENT NO.	DOCUMENT NAME
Directive 93/42 CE	Council Directive 93/42/EEC of June 14, 1993 on medical devices
EN ISO 13485:2016	Medical devices. Quality Management Systems. System requirements for regulatory purposes

DOCUMENT NO.	DOCUMENT NAME
EN ISO 14971:2019	Medical devices. Applying risk management to medical devices
EN 1041: 2008 +A1:2013	Medical devices. Information provided by the manufacturer
EN ISO 10993-1:2020	Medical devices. Biological Action Evaluation. Part 1: Evaluation and testing within a risk management process
EN ISO 10993-4:2017	Medical devices. Biological evaluation of medical devices. Part 4. Testing of devices that interact with blood
EN ISO 10993-5:2009	Medical devices. Biological evaluation of medical devices. Part 5. Tests for in vitro cytotoxicity
EN ISO 10993-7:2008/AC:2009	Medical devices. Biological evaluation of medical devices. Part 7. Residual ethylene oxide content after sterilization
EN ISO 10993-10:2013	Biological evaluation of medical devices. Part 10. Tests for irritation and skin sensitization
EN ISO 10993-11:2018	Medical devices. Biological evaluation of medical devices. Part 11. Studies of general toxic effects
EN ISO 15223-1:2021	Medical devices. Symbols used in labeling medical devices, on labels and in accompanying documentation. Part 1. Basic requirements
EN ISO 11607-1:2020	Packing of medical devices subject to final sterilization. Part 1. Requirements for materials, barrier systems for sterilization and packaging systems
EN ISO 11607-2:2020	Packaging of medical devices that are sterilized at the final stage of production. Part 2. Evaluation requirements for forming, sealing and assembly processes
EN ISO 11135:2014/A1:2019	Sterilization of health care products. Ethylene oxide. Requirements for the development, validation and routine control of a sterilization process for medical devices.
EN ISO 11737-1:2018	Sterilization of medical products. Microbiological methods. Part 1. Determination of microbial populations on products
EN ISO 14644-1:2015	Clean rooms and related controlled environments. Part 1. Classification of air purity based on particle concentration
EN ISO 14644-2:2015	Clean rooms and related controlled environments. Part 2. Routine control to confirm continued compliance with air cleanliness based on particle concentration

DOCUMENT NO.	DOCUMENT NAME
EN ISO 14937:2009	Sterilization of medical products. General requirements for characterization of a sterilizing agent and the development, validation and routine control of a sterilization process for medical devices
MEDDEV.2.7.1 Rev.4	Clinical Assessment: Guidelines for manufacturers and Notified Bodies
MEDDEV.2.12/1 Rev.8	Manufacturer's incident report form. System of supervision of medical devices
ГОСТ ISO 11607-1-2018	Packing of medical devices subject to final sterilization, Part 1. Requirements for materials, barrier systems for sterilization and packaging systems
ГОСТ ISO 14971-2021	Medical devices. Applying risk management to medical devices
ГОСТ Р ISO 15223-1-2020	National standard of the Russian Federation. Medical devices. Symbols used in labeling, packaging and documentation, associated with medical devices
ГОСТ ISO 10993-1-2021	Medical devices. Biological evaluation of medical devices
ГОСТ Р 52770-2016	Medical devices. Security requirements. Methods of sanitary-chemical and toxicological tests
ГОСТ 31214-2016	Medical devices. Requirements for the materials and documentation provided for toxicological, sanitary-chemical, sterility and pyrogenicity tests
ГОСТ 31870-2012	Drinking water. Determination of element content by atomic spectrometry methods
ГОСТ 31209-2003	Containers for blood and blood components. Chemical and biological security requirements and testing methods.
ГОСТ Р 55227-2012	Water. Methods of detecting formaldehyde levels
MYK 4.1.3166-14	Gas chromatographic determination of hexane, heptane, acetaldehyde, acetone, methyl acetate, ethyl acetate, methanol, isopropanol, acrylonitrile, n-propanol, n-propyl acetate, butyl acetate, isobutanol, n-butanol, benzene, toluene, ethylbenzene, m-, o- and p- xylenes, isopropylbenzene, styrene, alpha-methylstyrene in water and aqueous extracts from materials of various compositions
MYK 4.1.3169-14	Gas chromatographic determination of dimethyl phthalate, dimethyl terephthalate, diethyl phthalate, dibutyl phthalate,

DOCUMENT NO.	DOCUMENT NAME
	butyl benzyl phthalate, bis(2-ethylhexyl) phthalate and dioctyl phthalate in water and aqueous extracts from materials of various compositions
Государственная Фармакопея Российской Федерации. ОФС.1.2.4.0005.XV	Pyrogenicity
OCT 42-513-99	Complex single use exfusion, infusion and transfusion devices. Chemical and biological security requirements. Clinical testing methods in the environment of the manufacturer's facilities
ГОСТ Р ISO 10993-2-2009	Medical devices. Biological evaluation of medical devices. Part 2. Animal welfare requirements
ГОСТ ISO 10993-4-2020	Medical devices. Biological evaluation of medical devices. Part 4. Testing of devices that interact with blood
ГОСТ ISO 10993-5-2011	Medical devices. Biological evaluation of medical devices. Part 5. Tests for in vitro cytotoxicity
ГОСТ ISO 10993-7-2016	Medical devices. Biological evaluation of medical devices. Part 7. Residual ethylene oxide content after sterilization
ГОСТ ISO 10993-10-2011	Medical devices. Biological evaluation of medical devices. Part 10. Tests for irritation and skin sensitization
ГОСТ ISO 10993-11-2021	Medical devices. Biological evaluation of medical devices. Part 11. Studies of general toxic effects
ГОСТ ISO 10993-12-2015	Medical devices. Biological evaluation of medical devices. Part 12. Sample preparation and reference materials

17. Information about marking

Labels are in Russian, English and Turkish.

Primary packaging labeling

- Medical device name;
- Name and address of the manufacturing company;
- "1984 - Reference to a notified body with CE marking" symbol
- "Catalogue number" symbol
- "Expiration date" symbol
- "Caution" symbol

- "Refer to instructions for use" symbol
- "Sterilized with ethylene oxide" symbol
- "Single-use product" symbol
- "Do Not Re-sterilize" symbol
- "Batch Code" symbol
- "Keep away from sunlight" symbol
- "Do not use if package is damaged" symbol
- "Does not contain DEHP" symbol
- "Product made without the use of natural rubber latex" symbol
- "Store at temperatures between 0 and 30 degrees" symbol
- "Filling/storing volume per meter of tube length" symbol
- "Humidity limit" symbol
- Nominal volume;
- Batch/Lot Code;
- Shelf life.

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Purpose of usage: This product is used for mixing and administration of Total Parenteral Nutrition (TPN) solutions.

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
YYYY/AA

KPSM.LC.FR.17-1015015

Layout of additional marking of primary packaging in Russian for circulation in the Russian Federation

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
Secondary packaging labeling (cardboard box):

Medical device name

- Name and address of the manufacturing company
- "1984 - Reference to a notified body with CE marking" symbol
- "Catalogue number" symbol
- Catalog number;

- "Reuse Prohibited" symbol
- "Do Not Re-sterilize" symbol
- "Ethylene oxide sterilization" symbol
- "Attention! Refer to the instructions for use" symbol
- "Do not use if package is damaged" symbol
- "Batch Code" symbol
- Symbol "Use up to ..."
- Use up to
- "Keep away from sunlight" symbol
- "Temperature Range Limit" symbol
- "Manufacturer" symbol
- "Does not contain DEHP" symbol
- "Product made without the use of natural rubber latex" symbol
- "Quantity in Box" symbol
- Quantity in box
- "Batch/lot code/" symbol
- Batch/lot code
- "Humidity limit" symbol
- "Keep away from sunlight" symbol
- "See instructions for use" symbol

Secondary packaging labeling (cardboard box)

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
Symbols used on labels:

Symbol	Description
	Manufacturer
	CE marking
	Catalog number
	Use up to...
	Attention! Refer to the instructions for use
	See instructions for use
	Ethylene oxide sterilization
	Reuse Prohibited
	Batch code
	Keep away from sunlight
	Protect from moisture
	Do not use if package is damaged
	Does not contain diethylhexyl phthalate
	Latex-free natural rubber
	Do not re-sterilize

Symbol	Description
	Manufacturer
	Temperature Limit
	Humidity range
	Apyrogenic

18. Warnings and Precautions

Read the manufacturer's instructions for use before application.

Prior to opening the package, see expiration date.

Do not use the product if the packaging is damaged or opened.

Always use aseptic methods when working with the products.

The bag is disposable, do not reuse or re-sterilize.

Products contaminated with cytotoxic drugs must be disposed of in a container for cytotoxic waste.

Store EVA bags away from light.

Use an infusion line to administer the prepared medicine to the patient.

This product should be used only by qualified personnel or those with special training.

If the physician has prescribed light- or UV-sensitive solutions, use a light-protected EVA bag (yellow).

The storage time of solutions prepared in EVA bags should be determined according to the information sheet prepared by the drug manufacturer.

The product should have protective covers that should not be reused, as there is a risk of contamination if they are removed.

The dimensions shown on the bag are indicative and cannot be used to determine the exact volume of solutions.

Medical device inlets to be used with the EVA bag must have a luer lock (6%).

If an EVA bag is found to be leaking, contaminated, perforated, torn or if the contents is missing, discard the EVA bag in a cytotoxic waste container and use a new sterile EVA bag.

Ports of any medical devices used with the EVA bags must have the structure of a spike.

Never put the needle on the needleless connector.

The device should not be used in blood or blood product transfusion.

Due to the risk of intravenous occlusion, the device should be replaced every 24 hours.

Keep the source vessel at a higher level than the catheter and connection points.

19. Packaging information

Primary packaging:

EVA bag for the administration of intravenous solutions and parenteral nutrition and Luer-Lock stopper with filter in sterile packaging, are placed and sealed inside of sterile packaging made from medical craft paper and laminated film, on the packaging line under aseptic conditions.

Initial packaging (laminated film and medical craft paper) is manufactured in accordance with the EN 868-5 and ISO 11607 standards, medical sterilization pouch complies with the requirements of the council directive 93/42/EEC concerning medical devices, with the 2007/47/EC amendment.

Packaging materials and construction provide:

- High barrier packaging for sterile medical devices;
- Maintaining of the sterility during the entire shelf life of the device under the usual storage and transportation conditions;
- Impossibility of second-hand packaging: sterility indicator for the medical device.

In the middle of both sides of the craft paper packaging, there are grooves that facilitate pre-sterilization packaging and opening.

Sterilization method declared by the manufacturer of the primary packaging, must be the ethylene oxide method.

None of the primary packaging materials contain chemical, toxic components or latex.

Size of the primary packaging:

- 200 x 300 mm (contains transparent/light-protected EVA bag for the administration of intravenous solutions and parenteral nutrition, volume - 150, 250, 500 and 1000 ml, and a Luer-Lock stopper with filter);

- 300 x 320 mm (contains transparent/light-protected EVA bag for the administration of intravenous solutions and parenteral nutrition, volume - 2000 ml, and a Luer-Lock stopper with filter);

- 400 x 400 mm (contains transparent/light-protected EVA bag for the administration of intravenous solutions and parenteral nutrition, volume - 3000 ml and 4000 ml and a Luer-Lock stopper with filter);

Transport packaging (cardboard box):

Bags in the primary packaging are stacked in a shipping carton (double wall corrugated cardboard) in an amount of 50 or 100 pieces. In order to prevent violation of integrity of the transport packaging, it is covered with tape.

Information about the materials used in the production of packaging

Component name	Name of material	Material class	Manufacturer of material
Primary packaging and packing of the plug	Medical bag for EO sterilization	BOPET-PP - Medical Kraft Paper	SUDEMED and PMS
Secondary packaging	Cardboard	Double wall corrugated cardboard	Marmara koli ve ambalaj san.

Characteristics of the primary packaging material, that undergoes final sterilization.

Parameter	Value
Material density	60 - 70 g/m ²
Material thickness	75-85 µm.
Breathability	11,4 µm/(Pas)
Peeling force	≥ 1.2 N (EN868-5)
MD tensile strength limit	6.49 kN/m
CD tensile strength limit	3.4 kN/m
Elongation of the packaging material at rupture	-
Tensile strength	350 kPa

MD tensile strength	600 mH
CD tensile strength	650 mH
Seam width	≥ 3 mm
pH	5,0 - 8,0 (ISO 6588-2)
Chloride levels	$< 0,05\%$ (ISO 9197)
Sulfate level	$< 0,25\%$ (ISO 9198)
Material thickness (film)	65 ± 5 μ m
Heat seal (film)	110 ± 5 °C
MD tensile strength (film)	490 ± 49 mH
CD tensile strength (film)	580 ± 49 mH

Dimensions and weights of products in transport packaging (double wall corrugated cardboard) are shown in the tables below:

EVA infusion system (transparent bag) for administration of intravenous solutions and parenteral nutrition

Volume	Quantity	*Box size	*Box weight
150 ml	100 pc.	60 x 40 x 28 cm	4.55 kg
250 ml	100 pc.	60 x 40 x 28 cm	5.05 kg
500 ml	100 pc.	60 x 40 x 28 cm	5.38 kg
1,000 ml	100 pc.	60 x 40 x 28 cm	5.95 kg
2,000 ml	50 pc.	60 x 40 x 28 cm	4.85 kg
3,000 ml	50 pc.	60 x 40 x 28 cm	5.90 kg
4,000 ml	50 pc.	60 x 40 x 28 cm	5.98 kg

*Note: The permissible deviation of the values is $\pm 5\%$.

EVA infusion system (light-protected bag) for administration of intravenous solutions and parenteral nutrition

Volume	Quantity	*Box size	*Box weight
150 ml	100 pc.	60 x 40 x 28 cm	4.55 kg
250 ml	100 pc.	60 x 40 x 28 cm	5.05 kg
500 ml	100 pc.	60 x 40 x 28 cm	5.38 kg
1,000 ml	100 pc.	60 x 40 x 28 cm	5.95 kg
2,000 ml	50 pc.	60 x 40 x 28 cm	4.85 kg
3,000 ml	50 pc.	60 x 40 x 28 cm	5.90 kg
4,000 ml	50 pc.	60 x 40 x 28 cm	5.98 kg

*Note: Permissible deviation of values $\pm 5\%$

20. Information about sterilization

EtO sterilization (ethylene oxide). In accordance with the 93/42/EEC Medical Device Directive and the EN ISO 11135-1 sterilization standard for medical devices, EN ISO 14937 standard for sterilization of medical devices, EN ISO 15223-1 labelling standard; the label is marked as defined above.

See sterilization inspection report.

Information about pyrogenicity: the product is apyrogenic.

21. The presence in the product of drugs, materials of animal and/or human origin

Medicinal substances, there are no materials of animal and/or human origin in the composition of the product.

22. Form and duration of contact with the patient.

Indirect contact with the blood circulation system. Contact with the patient's body through solution or parenteral nutrition in the device.

Short-term contact (A category)

23. Environmental requirements

During storage, the medical device does not have a negative impact on the patient or the environment.

24. Operation conditions:

The medical device is intended for use in medical facilities and at home

Operation conditions:

Temperature range: 0 to 27°C;

Relative humidity range: 20 to 60%;

Atmospheric pressure range: 1 Atm Pressure.

25. Storage conditions

This medical device should not be exposed to direct sunlight and should be protected from high temperature and humidity

Storage conditions:

Temperature range: 0 to + 30°C;

Relative humidity range: 20 to 60%;

Atmospheric pressure range: 500 – 1,060 hPa.

26. Shipping requirements

This medical device must be transported in the manufacturer's packaging by all types of covered vehicles in accordance with the shipping conditions applicable to each mode of transportation.

Shipping conditions:

Temperature range: 0 to 30°C;

Relative Humidity Range: 0 - 70%;

Atmospheric pressure range: 500 – 1,060 hPa.

27. Shelf life

The shelf life is 5 years.

Shelf life confirmed by the manufacturer (see validation report on stability)

The sterility of the contents must be maintained throughout the shelf life of the product, subject to storage conditions and the integrity of the packaging.

The shelf life is shown on the package label.

28. Manufacturer's warranty

The manufacturer guarantees the quality of the medical device if the requirements for transportation, storage and use of the product set forth in the technical documentation are met.

Do not use an expired medical product.

Product warranty period: 5 years from the date of sterilization.

29. Maintenance

This medical device is a single-use, maintenance-free device.

30. Requirements for safe destruction and disposal

Used products, as well as failed unused products are disposed of accordance with the hazard class specified in SanPiN 2.1.3684-21.

Disposal must be carried out in accordance with the requirements of SanPiN 2.1.3684-21:

- for used items as class B "epidemic hazardous" waste;
- for packaging materials and items with expired shelf life, which were not used for their intended purpose, as waste of class A "epidemiologically safe waste, the composition of which is close to MSW".
- the items are disposed of in accordance with the medical facility regulations.

31. Contact information (complaints)

For all questions related to the circulation of the medical device in the Russian Federation, you should contact the Authorized Representative of the manufacturer:

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