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arthroscope in direct contact with patient tissue or combustible materials, as burns may result. Lower the light source output when working in close proximity to the object.

- This device has not been evaluated for safety and compatibility in the MR environment.

- The arthroscope light post and adaptor may exceed temperatures of 106° F (41° C). Arthroscopes should not be placed on the patient or on combustible materials, as burns may result.

- To prevent potential safety hazard to the patient caused by accidental loss of function of the device (i.e., front end damage by surgical instruments) it is recommended to have an additional sterile "stand-by" device during surgical procedures.

- When using RF surgical equipment, keep the working part of the active electrode in the field of view to avoid accidental RF burns.

Avoid contact with metal parts of the arthroscope and other conductive accessories by ensuring, before activation of the RF output, that the active electrode is at a sufficient distance from the tip of the arthroscope. Ensure that only medical electrical equipment that complies with IEC 60601-1 (and its relative particular standards) is connected to, or used in conjunction with, the arthroscope.

- Video arthroscope and camera head MUST be detached prior to cleaning and sterilization.

- Failure to detach them from each other prior to cleaning and sterilization will compromise sterility and constitute a patient hazard.

- Do not modify this equipment. Attempts to modify this equipment will void the warranty.

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Precautions

U.S. Federal law restricts this device to sale by or on the order of a physician.

- Prior to each use, inspect the device to ensure it is functioning properly and not damaged. Do not use a damaged device.

- To avoid fogging during surgery, the rear portion of the arthroscope must be entirely free of moisture before attachment to either the camera head or camera coupler.

- The aluminum oxide powder (REF 7205548) is shipped non-sterile and should only be used during the cleaning process of autoclavable arthroscopes prior to sterilization.

- Do not use any ultrasonic cleaning methods. The energy transmitted through fluid cavitation will damage seals and optical surfaces and will void the warranty.

- Only arthroscopes marked "Autoclavable" can be steam autoclaved.

With the exception of these arthroscopes which have been designed specifically for autoclavability, arthroscopes cannot be steam autoclaved. Steam

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autoclaving arthroscopes not specifically designed for this purpose would permanently damage the optical components and void the warranty. Do not expose these arthroscopes to temperatures exceeding 140° F (60° C).

- Before autoclaving, ensure that exposed optical surfaces of the arthroscope (distal tip and proximal window) are thoroughly clean and free of residue.

- The extreme heat from steam autoclaving will cause foreign materials left on these optical surfaces to harden and discolor. Any foreign matter remaining on the distal tip will reduce image quality and may lead to a completely obscured picture.

- Foreign matter remaining on the fiber surface of the light post after cleaning may tend to burn and discolor the surface when exposed to a high intensity light source.

- Allow the arthroscope to cool to room temperature after it has been autoclaved.

- Caution must be taken to dry-out mating portions of both the arthroscope and camera head; otherwise trapped moisture may condense on glass surfaces, blurring the image.

- Arthroscopes are delicate instruments and must be handled carefully.

Isolate the arthroscope from other instruments immediately after the case to prevent the arthroscope from being damaged from contact with other instruments.

- In some cases, a specific filter must be put between the arthroscope and camera head to prevent camera head damage by high-power laser radiation. Contact your laser supplier for details.

- To prevent scope damage by high-power laser radiation, always ensure that the laser delivery fiber is visualized through the arthroscope and is not directed at the arthroscope before energizing the laser.

- Any mechanical manipulation of the direct view arthroscope eyepiece may result in seal breakage; therefore, do not attempt to remove the eyepiece.

- When arthroscopes are used with laser equipment, appropriate filtering spectacles must be worn by the operating team.


- To comply with IEC 60601-2-18 standard requirements for BF-type and CF-type equipment, the light cable and camera head must provide electrical isolation from the light source and camera control unit, respectively.

Safety Precautions

- Personal Protective Equipment (PPE) should be worn when handling or working with contaminated or potentially contaminated devices. PPE includes: gown, mask, goggles or face shield, and shoe covers.

- Universal precautions are standards of infection control practices designed to reduce the risk of transmission of bloodborne infections.

- Precautions should be observed by all hospital personnel that work with contaminated or potentially contaminated devices. Exercise caution

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For VideoArthroscopes

• Check that the O-ring located at the base of the mounting threads is in place and undamaged. If damaged or missing, replace it with a new O-ring which can be ordered through Smith & Nephew Customer Service.

• Check that the proximal window is clean.

• Screw the VideoArthroscope onto the desired camera head in a clockwise direction, hand tightening only, until the VideoArthroscope base is firmly seated against the camera flange.

• Attachment of arthroscope accessories, such as cannulas, is generally accomplished via the "J" lock. To lock, simply push down and twist clockwise.

CAUTION: To avoid fogging during surgery, the rear portion of the arthroscope must be entirely free of moisture before attachment to either the camera head or camera coupler

Maintenance

We recommend that you inspect the endoscope carefully before and after the case for possible signs of damage. This will allow immediate detection of minor damage which, if repaired immediately, will increase the life of your investment.

First, check the image quality of the arthroscope by looking through it, if it is a direct-view arthroscope, or by viewing the monitor if it is a VideoArthroscope. Impaired image quality can result from any of the following conditions:

• Check the distal and proximal lenses of the arthroscope for cracked or scratched lenses.

• Inspect the shaft of the arthroscope to detect any noticeable dents. Dents along the shaft might cause lenses inside the tube to crack.

• A foggy or cloudy image can be the result of moisture entering the optical system or cleanliness of exterior surfaces. Check the surface cleanliness of the distal and proximal lenses. When viewing reflected light, the surfaces should appear smooth and shiny. For autoclavable arthroscopes, specific instructions to remove deposits are provided below. For other arthroscopes, please follow the "Recommended Cleaning Steps."

As a second step, check the illumination system of the arthroscope. Reduced brightness can result from fiber damage:

• Check for fiber optic damage in the arthroscope by holding the distal end of the arthroscope toward a light and observing the light post on the hub. The center of the light post should appear clear or white.

Noticeable black spots indicate serious damage to the fiber optic bundle in the arthroscope. This will affect light transmission and the brightness of the image viewed on the monitor.

• Check the light cable for damaged fiber optics by holding one end of the cable toward a light and observing the other end. Broken fibers will appear as black spots in the light field. A damaged light

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cable will affect its ability to transmit light and the brightness of the image viewed on the monitor.

Note: VideoArthroscopes are designed to be used directly coupled to an appropriate camera head, and therefore cannot be used as direct-view arthroscopes. When viewing through a VideoArthroscope, you may see small particles that have been magnified through the optical system but that will in no way affect the image displayed on the monitor when properly viewed in conjunction with the camera head.

Service Life Checklist

Perform a performance check for the device before each use case begins.

Check the following items to ensure that the device functions properly and the device is ready to use.

- Confirm that the labeling information on the device is still readable.
- Check for any damage to the lens at the proximal and distal end of the device.
- Check for any damage to the device (cracked, damaged, bent misplaced and/or missing).
- Confirm an image can be visually seen through the endoscope. For direct view scopes this can be done using the device only, for video scopes a camera head/CCU (camera control unit) will be needed to aid the visualization by displaying on a monitor.
- Confirm the focusing mechanism can be properly turned and focus can be adjusted (video scopes only).
- Confirm the direction of view can be adjusted (video scopes only).
- Confirm no moisture is observed on the internal side of the lens.

Cleaning Instructions

 **WARNING:** This product is shipped non-sterile. It must be disassembled, cleaned, and sterilized before the first use. It must be disassembled, cleaned, and sterilized before every subsequent use.

The Smith & Nephew line of rigid Arthroscopes should be soaked and cleaned for five minutes in an enzymatic, neutral pH cleaner recommended for reprocessing of medical endoscopes, prepared according to the manufacturer's instructions. Any light post adaptors must be removed prior to cleaning. Thoroughly clean all surfaces of the arthroscope and adaptor, allowing for any recommended soak time.

Thoroughly rinse the arthroscope and adaptor for two minutes with water to completely remove the

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cleaning solution. Dry and reassemble the adaptor to the arthroscope.

The arthroscope has three optical surfaces that must be thoroughly cleaned and checked routinely to ensure both maximum transmission of illumination and a high-quality image. These are the distal tip, the proximal window or eyepiece, and the fiber optic light post.

Recommended Manual Cleaning Steps

1. Rinse in cold utility water $<113^{\circ}\text{F}$ (45°C) for a minimum of 2 minutes to remove visible debris and to prevent coagulation of blood.

2. Immerse and soak for a minimum of five (5) minutes in a neutral pH enzymatic detergent*.

3. With the device completely immersed in the cleaning solution*, remove additional soil from the surfaces and from the challenging design features (i.e. crevices) using common hospital cleaning tools. Move all moveable parts and remove visible debris using a brush.

a. Scrub crevices with a brush.

b. Scrub the surfaces of the device with a toothbrush styled brush.

4. Rinse the device thoroughly for a minimum of 2 minutes using a sufficient volume of critical water that will completely immerse the device. Move any movable features of the device during rinsing.

5. Inspect the device using lighted magnification. If soil is observed repeat cleaning steps.

CAUTION: Do not use any ultrasonic cleaning methods. The energy transmitted through fluid cavitation will damage seals and optical surfaces and will void the warranty.

Automated Cleaning Steps

1. Minimum cycle parameters: a. Five (5) minute cold prewash $<113^{\circ}\text{F}$ (45°C).

- b. Five (5) minute neutral enzymatic detergent wash*.

- c. Five (5) minute neutral detergent wash*.

- d. Two (2) minute rinse at 113°F (45°C)

Note: For the enzymatic and neutral detergent wash cycles, follow the detergent manufacturer's recommended instructions for the washer that is used for cleaning.

Note: For the enzymatic and neutral detergent wash cycles, follow the detergent manufacturer's recommended instructions for the washer that is used for cleaning.

* Enzymatic detergents with a neutral pH range between 6.0 and 8.0 are recommended. Detergents with a pH outside this range can have an adverse effect or be damaging to some medical devices and containment devices. Enzymatic detergents aid in the removal of organic soil such as blood. Detergents should be used at the concentration and temperature recommended by

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the detergent manufacturer. 1 to 4 mL/L of Prolystica 2X Concentrate Enzymatic Presoak and cleaner (neutral pH enzymatic detergents manufactured by STERIS Corporation) or equivalent was used in the validation of the cleaning processes for Smith & Nephew Endoscopy devices. Composition of «Prolystica 2X Concentrate Enzymatic Presoak and cleaner»

Name	Product Identifier	%
Alcohols, C9-11, ethoxylated surfactant	(CAS No) 68439-46-3	1 - 5
Glycerin	(CAS No) 56-81-5	1 - 5
Citric acid	(CAS No) 77-92-9 (REACH No) 01-2119457026-42-0067	1 - 5
N,N-Dimethyloctadecylamine oxide	(CAS No) 2571-88-2	1 - 2
Subtilisins (proteolytic enzymes)	(CAS No) 9014-01-1	0.1 - 1

Thermal Disinfection

1. The parameters for thermal disinfection would be 1 min at 195.8° F (91° C).
2. Inspect the device and check the service life checklist prior to use.

Removing Deposits – Autoclavable Scopes Only

The aluminum oxide powder (REF 7205548) is intended to be used, when required, to remove deposits on optical surfaces thus maintaining optical integrity and allowing the scopes to perform as per their intended uses.

Note: Cleaning with aluminum oxide powder should only be performed prior to sterilization when the image as viewed through the arthroscope is cloudy and not as part of your routine cleaning procedures.

To remove deposits, moisten a cotton-tipped swab with water. Dip the moistened swab into the aluminum oxide powder. Position the swab tip against the optical surface of the scope to be cleaned. While applying gentle pressure to the swab, scrub the window to remove any deposits. For smaller windows, rotate the swab instead of scrubbing. Rinse the optical surface with tap water and brush with a soft bristle brush while under running water to remove the aluminum oxide powder.

Clean optical surfaces as instructed in the "Recommended Cleaning Steps".

CAUTION: The extreme heat from steam autoclaving will cause foreign materials left on these optical surfaces to harden and discolor. Any foreign matter remaining on the distal tip will reduce image quality and may lead to a completely obscured picture.

CAUTION: Foreign matter remaining on the fiber surface of the light post after cleaning may tend to burn and discolor the surface when exposed to a high intensity light source.

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Sterilization

WARNING: This product is shipped non-sterile. It must be disassembled, cleaned, and sterilized before the first use.

It must be disassembled, cleaned, and sterilized before every subsequent use.

READ ALL WARNINGS AND CAUTION BEFORE STERILIZING.

Note:

- Disinfection is not a substitute for sterilization.
- Always inspect for damage and debris prior to sterilization.

- Select the appropriate sterilization method for the equipment. For sterilization of light guides, couplers and camera heads refer to the "Sterilization" section of the respective product's Instructions for Use.

- Smith & Nephew arthroscopes should be sterilized in containers which are designed to secure the instruments in place while allowing proper circulation of the sterilant medium.

- These parameters have been validated to ensure sterility. Sterilizer functioning should be monitored at regular intervals with biological indicators to ensure products have been subjected to sterilization conditions.

- VideoArthroscopes are designed to be used with solid-state video cameras. These scopes thread directly onto the camera head with no additional coupler components.

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CAUTION: Allow the arthroscope to cool to room temperature after it has been autoclaved.

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Ordering Information

Adverse Reaction Reporting

If, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and to your national authority and complaints@smith-nephew.com

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Warranty

Smith & Nephew products are guaranteed to be free from defects in material and workmanship for the warranty period for a particular product, beginning from date of invoice. Refer to the current Smith & Nephew Product Catalog or contact Smith & Nephew Customer Service for specific warranty information.

This limited warranty is restricted to repair or replacement by Smith & Nephew, at its option, of any product found to be defective during the warranty period. Damage inflicted to a product by the user that causes it to be unsuitable for refurbishment may result in additional charges, regardless of warranty status. All warranties apply to the original buyer only. In no event shall Smith & Nephew be liable for any anticipated profits, consequential damages or loss of time incurred by the buyer with the purchase or use of any product. **NO OTHER WARRANTY, EXPRESSED OR IMPLIED, IS GIVEN.**

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On the territory of the Russian Federation, it is recommended to use the following sterilization parameters:

Method	Cycle
Pre-vacuum steam	Temperature: 132°C+2°C /270°F+35,6 °F* Exposure time**: 4 minutes Drying time**: 30 minutes Steam pressure in the sterilization chamber, MPa (kgf / cm ²): 0,20 (2,0) ± 0,02 (0,2)*

*in accordance with the information of the MU-287-113 guidelines for DISINFECTION, PRE-STERILIZATION CLEANING AND STERILIZATION OF MEDICAL DEVICES (APPROVED by the DEPARTMENT of GOSSANEPIDNADZOR-Ministry of HEALTH of the Russian Federation on DECEMBER 30, 1998)

** The sterilization exposure time corresponding to the sterilization temperature value, taking into account the maximum temperature deviations in the loaded sterilization chamber, must be indicated in the operating manual for a specific sterilizer model.

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

Interpretation of symbols used on the package.

	Non-steril
	Serial Number
	Catalog number
	Batch code
	Temperature limitation
	Humidity limitation
	Atmospheric pressure

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The existing Smith & Nephew arthroscopes are the comparable with video systems LENS (registration certificate № RZN 2019/9141) и 4K LENS (registration certificate № RZN 2021/15938). To check the compatibility of arthroscopes with other Smith & Nephew systems and devices please contact authorized Smith & Nephew representative or distributor.

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The subject Smith & Nephew products do not incorporate materials of animal and/or human origin, nor are they present in the manufacturing process or present in the final device.

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The subject Smith & Nephew products do not incorporate a medicinal substance, nor is it present in the manufacturing process or present in the final device, therefore this section does not apply.

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The packaging provides protection against mechanical and climatic factors during transportation and storage under the following conditions:

Transport and/or Storage

Ambient Temperature	-29...60°C (-20...140°F)
Relative Humidity	15 ...85%
Atmospheric Pressure	70...106kPa
Operating Altitude	3000 meters maximum

Used devices must be transported to the processing location (i.e. Sterile Processing) in closed or covered containers to prevent contamination risks.

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The shelf life of the product is not limited.

The following products can withstand up to 300 sterilization cycles:

1. Arthroscope, viewing angle 0°, diameter 4.0 mm, length 160 mm, camera mounting type – eyepiece cup;
2. Arthroscope, viewing angle 30°, diameter 4.0 mm, length 160 mm, camera mounting type – eyepiece cup
3. Arthroscope, viewing angle 70°, diameter 4.0 mm, length 160 mm, camera mounting type – eyepiece cup, (RND TIP)
4. Arthroscope, viewing angle 45°, diameter 4.0 mm, length 160 mm, camera mounting type – eyepiece cup.

The service life of the product is at least 1 year.

The following products can withstand up to 100 sterilization cycles:

1. Arthroscope, viewing angle 30°, diameter 4.0 mm, length 160 mm, camera mounting type C-Mount;
2. Arthroscope, viewing angle 70°, diameter 4.0 mm, length 160 mm, camera mounting type C-Mount, (RND TIP).

The service life of the product is at least 4 months.

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